

Since 1980



KOSDAQ
LISTED COMPANY



Smart-Lab™
On Mobile



- ✓ CE Certified
- ✓ Excellent Service Quality - Certified
- ✓ GMP Certified for Medical Instruments
- ✓ ISO 9001 Certified
- ✓ ISO 13485 Certified for Medical Instruments
- ✓ KOLAS Accredited KC23-433
- ✓ UL Certified
- ✓ Patented(No. 10-0554098) Jog Shuttle Control System
- ✓ PL Insurance, Validation(IQ & OQ) Service

DAIHAN-Sci.com Shopping-mall
ALLforLAB.com Open-market

DAIHANmedical.com
DAIHANused.com
LabRental.co.kr
Lab-Service.co.kr
Pipettermall.com

Calibration.kr (ml & g)

검색 대 한 과 학

2026 WORLDWIDE DAIHAN™ Instruments for Laboratory and Medical use



Since 1980



DAIHAN
SCIENTIFIC
GROUP

Your Neighborhood Friend
SOUTH KOREA

Sales & Service
- by Distributors



Ordering & Inquiring

Easy, Fast, Comfortable & Economical!



via Phone

DAIHAN's Official Distributor Nearest You
or DAIHAN Call Center
080-008-3000 or 070-5020-1516
Int'l ; +82-70-5020-9000



via Fax

DAIHAN's Official Distributor Nearest You
or DAIHAN Fax ; 033-737-7575
Int'l ; +82-2-963-5231



via Website

DAIHAN-Sci.com / ALLforLAB.com
DAIHANmedical.com
DAIHANused.com
LabRental.co.kr
Lab-Service.co.kr
Pipettermall.com
Calibration.kr (ml & g)



via Mail

dh@daihansci.co.kr
Int'l ; info@daihan-sci.com

Comprehensive Solution

in Laboratory, Medical & Medical-Lab, and Industry MRO



DAIHAN General Catalog



Time-honored Business, **DAIHAN Scientific Group** manufactures a wide range of "DAIHAN™" Brand Premium Laboratory / Medical Instruments, and has the largest comprehensive supply system with over 30,000 items from 300 suppliers specialized in Laboratory, Medical & Medical Lab and Industrial MRO products. "DAIHAN Scientific Catalog" contains about 18,000 Items that you might surely want for your satisfaction in your field. Based on our own SCM Systems and WMS Wonderful Warehouse, we have realized a "One-day Delivery Service" and have the Highest Customer Satisfaction in the market. It will give people a much Better Quality of Life. This is our mission. We have excellent 120 global partners all over the world, and an On-Line shopping mall "DAIHAN-Sci.com" & "ALLforLAB.com" where you can easily get all products information anytime, anywhere. Wherever you are, the DAIHAN worldwide organization will always be with you for human health, safety and best service as "Your Neighborhood Friend".

Thank you very much! 23 Mar. 2026

Shu, Hyun Jeong
President, CEO

Shu, Chong Ku
Chairman



A

Anemometers / Airflow Meters	206 ~ 207
Autoclaves / Steam Sterilizers	299 ~ 304
Bench-top Front Door.....	299
Front Door.....	304
Top-loading Vertical New	300 ~ 303
Accessories.....	303, 307
Autoclaving Deodorant.....	303
Sterilization Drapes New	307
Sterilization Indicators New	307

B

Bag / Sample Mixer	66
Balances	13 ~ 46
[d] 0.01mg Analytical Semi-micro.....	14 ~ 15
[d] 0.1~1mg Analytical New	15 ~ 19
Lab Precision New	20 ~ 25
Lab-General.....	26 ~ 29
Medical / Healthcare.....	35 ~ 37
Baby Scales.....	35
Mobile Chair-type Scale New	37
Step-On Body Scales.....	36
Moisture Analyzers New	38 ~ 39
Pocket Balances.....	40
Scales.....	30 ~ 35
Multi-Function.....	30 ~ 31
High Capacity.....	32 ~ 33
Pet Scales.....	34
Kitchen Scale New	35
Weight.....	41 ~ 42
"E1"-Class.....	41
"F1"- & "M1"-Class.....	42
Weight Calibration & Recalibration.....	43
Accessories.....	43 ~ 46
: Weighing Bottles, Brushes, Dishes, Forceps, Paper, Scoops, Spatulas, Spoons	
Bath	47 ~ 55
Basic 100°C.....	48
Basic 100°C Circulation(Internal).....	49
Basic 250°C Oil Bath.....	50
Cold Trap Baths.....	51
Extraction Water Baths.....	52, 124
Magnetic Stirring Baths.....	53
Shaking Water Baths.....	54
Viscosity Baths.....	55, 342
Bath Circulators	56 ~ 65
External, Chiller -30°C/-20°C~R/T New	56 ~ 57, 117 ~ 118
External, Heated/Refrigerated -30°C/-20°C~200°C New	58 ~ 59, 270 ~ 271
External / Internal, Heated +250°C.....	60
External / Internal, Heated / Refrigerated -35°C/-25°C~150°C.....	61 ~ 62, 119
Accessories.....	63 ~ 65
Circulation Hoses New	63 ~ 64
Cooling / Heating Transfer Fluids.....	65, 155, 269
Biological Safety Cabinet(BSC)	68 ~ 69
BioReactor / Fermentor <i>Coming Soon!</i>	252 ~ 255
Blenders	66 ~ 67
Bag / Sample Mixer.....	66
Blenders with Grinders New	67

C

Cabinets & Enclosures	68 ~ 78
Biological Safety Cabinet(BSC).....	68 ~ 69
Clean Benches / PCR Workstation.....	70 ~ 71
PC-type, 3-Side Viewing, Fabricated New	70
Steel-type.....	71
Ductless Filtering Reagent Cabinets.....	72 ~ 73
Powder-coated Steel Chamber New	72
PP Chamber.....	73
Safety Cabinets.....	74 ~ 78
Acids / Corrosives New	74 ~ 75
Flammables New	76 ~ 78
Vapor Purifier Box, Ductless Filtering.....	78
Carts	78 ~ 83
Plastic Carts.....	79 ~ 80
ABS Assembly Carts.....	79
ABS Drawer Dressing Carts.....	79 ~ 80
Stainless-steel Carts.....	80 ~ 83
Bottle Rack Carts.....	82 ~ 83
Dressing Carts New	80 ~ 81
Drying Cart with 3 Shelves and 24 Pegs.....	83
Laundry Bag Carts.....	82

<Continued from **Carts**>

Storage Carts with Curtain.....	83
Utility Carts.....	81 ~ 82
Steel Carts for Gas Cylinder.....	83
Centrifuges	84 ~ 97
0.2~2.0ml(50ml) Microcentrifuges.....	84 ~ 87
6-Place 5,500 rpm.....	84
12-Place 13,500rpm.....	85 ~ 86
24-Place 17,000rpm.....	87
24-Place Refrigerated, 17,000rpm.....	87
Clinical Centrifuges.....	88 ~ 90
Microhematocrit(Capillaries) 12,000rpm.....	88
1.5~50ml 4,000rpm.....	90
7 · 10 · 15ml 4,000rpm.....	89
General Purpose Centrifuges.....	91 ~ 92
0.2~100ml 16,500rpm Hi-Speed.....	91
1.5~500ml 4,000rpm Hi-Volume.....	92
Refrigerated Centrifuges.....	93 ~ 97
0.2~100ml 18,000rpm.....	93
0.2~250ml 12,000rpm New	94 ~ 95
1.5~250ml 4,300rpm.....	96
1.5~1,000ml 5,000rpm.....	97
Chillers / Cooling Circulators -30°C/-20°C~R/T New	56 ~ 57, 117 ~ 118
Clean Benches / PCR Workstation	70 ~ 71
PC-type, 3-Side Viewing, Fabricated New	70
Steel-type.....	71
Cryogenic LN₂ Storage Tanks	98 ~ 99
Canister-types 10~50Lit.....	98
Rack-types 47~95Lit.....	99

D

Desiccator Cabinets, PMMA	100 ~ 102
Auto-Dry.....	101
Gas Exchangeable.....	100
General Purpose.....	100
UV Protection Auto-Dry.....	102
Vacuum.....	102
Accessory; Desiccant Silica-gel.....	102
Dispensers, Wisdoser™	103 ~ 104
Application / Chemical Compatibility List.....	103
HF(Hydrofluoric) Acid Dispenser Kits.....	104
Standard Dispenser Kits.....	104
Dry Heat Sterilizers	305
Dryers	105 ~ 106
Drying Chambers.....	105
Labware Dryer.....	106

E

Electrophoresis	107 ~ 109
Gel Documentation Systems.....	107 ~ 108
UV Transilluminators.....	109
Evaporators	110 ~ 122
Rotary Vacuum Evaporators.....	110 ~ 116
0.1~5Lit, Vertical-type.....	110 ~ 113
0.5~2Lit, Vertical & Diagonal-type.....	114
5~50 Lit, Large Capacity, Vertical -type.....	115 ~ 116
Accessories.....	117 ~ 122
Bath Circulators, Heated / Refrigerated.....	61 ~ 62, 119
Chillers / Cooling Circulators New	56 ~ 57, 117 ~ 118
Evaporating Flasks & Adapters New	121 ~ 122
Vacuum Pumps for Evaporators.....	120
Extraction	123 ~ 124
Extraction Water Baths.....	52, 124
Heating Mantles with 3-& 6-Place.....	123, 168

F

Fat Analyzer New	126
Fiber Analyzer New	127
Fermentor / BioReactor <i>Coming Soon!</i>	252 ~ 255
Food Analysis New	125 ~ 130
Fat Analyzer.....	126
Fiber Analyzer.....	127
Kjeldahl Digestion System.....	128
Kjeldahl Nitrogen Analyzers.....	129 ~ 130
Freeze Dry System	131 ~ 133
8 & 16Lit Freeze Dryers.....	131 ~ 132
Accessories.....	132 ~ 133
Dry Chambers.....	132 ~ 133
Manifold, 12-port T-type.....	132
Freezers	273 ~ 281
Bio Freezers -35°C.....	281
Ultra Low Temp Freezers.....	273 ~ 280

<Continued from **Freezers**>

Energy Saving SMART ULT Freezer ReNew	273
Energy Saving ULT Freezing System.....	274
-95°C Upright-type ReNew	275
-86°C Personal Upright & Underbench.....	277
-86°C Upright-type ReNew	276
-45°C Upright-type.....	278
Accessories; Racks, Boxes, etc.....	279 ~ 280
Fume Hoods	134 ~ 138
Ducted Fume Hoods.....	137 ~ 138
Acids / Corrosives Resistance, PP.....	137
Multiuse, Aluminum Frame.....	138
Ductless Filtering Fume Hoods.....	135 ~ 136
Fume / Gas Purifiers	139 ~ 141
Dust Collector of Fine Dust · Fume · Bacteria.....	141
Fume Multi Gas Purifiers.....	139
Toxic Gas Multi-Purifiers.....	140
Vapor(Fume & Gas) Purifier Box.....	141
Furnaces	142 ~ 147
Dental Furnaces.....	142
1100°C, 0.4Lit.....	142
1650°C, 0.7 & 1.8Lit.....	142
Muffle Furnaces, 1000/1200°C.....	143 ~ 144
Muffle Furnaces, High Temp 1450/1650°C.....	145 ~ 146
Tube Furnaces, 1000°C.....	147
Furniture	148 ~ 151
Assembly Center Tables.....	149
Assembly Side & Corner Tables.....	150
Assembly Sink Tables.....	148
Chairs.....	151

G

Gas Detectors / Monitors	208 ~ 209
Gel Documentation Systems	107 ~ 108
Glove Boxes	152 ~ 154
10mm PMMA Complete Sets.....	152
15mm PMMA Heavy Duty Complete Sets.....	153
40mm PMMA Vacuum Complete Set.....	154
Accessory; Gloves For Glove Box.....	153

H

Heating Circulators => See "Bath Circulators"	58 ~ 65
Heating Blocks / Dry Bath Incubators	156 ~ 159
Heating, Single & Dual Block.....	156 ~ 157
Heating & Cooling.....	159
Heating · Cooling · Mixing.....	158
Heating & Mixing.....	158
Heating Mantles	160 ~ 170
Controllers.....	169
For Beakers · Flasks · Reaction Vessels.....	160 ~ 167
For Soxhlet / Extraction, 3 & 6 Places.....	123, 168
Heating Tapes & Custom-made.....	170
Homogenizers	171 ~ 174
Direct Control.....	171
Remote Digital Control.....	172
Tools & Stands.....	173
Ultrasonic Homogenizer New	174, 338
Homogenizing Stirrers	175
Hotplates	176 ~ 179
380°C, Analog & Digital.....	178
500°C, Solid Ceramic Glass.....	177
550°C, Hi-Light.....	176
Large Plate, 350°C.....	179
Hotplate Stirrers	180 ~ 182, 313 ~ 315
380°C, Analog & Digital.....	181, 314
500°C, Solid Ceramic Glass.....	180, 313
Multi Point 3- & 6-Place, 350°C.....	182, 315

I

Ice Makers	183
Flake-type.....	183
Snow-type.....	183
Incubators	184 ~ 203
CO ₂ Incubators.....	185 ~ 186
General Low Temp Mini.....	187
General Purpose.....	188 ~ 193
Forced-air.....	188 ~ 189, 192
Gravity-air.....	190 ~ 192
Multi-Room.....	193
Heating Blocks / Dry Bath Incubators.....	154 ~ 158
Humidity / Temp Test(Environmental) Chambers.....	194 ~ 196
Low Temp(B.O.D) Incubators.....	197 ~ 198
Plant Growth Chambers.....	199

<Continued from Incubators>

Shaking Incubators	200 ~ 203
Bench Top-type	200
Front Door-type	201
Mini, CO ₂ & Low Temp	185, 187
Multi-stack Refrigerated	203
Top Door-type	202

J

Jar tester / 6-Paddle Stirrer	204, 324
-------------------------------	----------

K

Kjeldahl Digestion System New	128
Kjeldahl Nitrogen Analyzers New	129 ~ 130

L

LN ₂ Storage Tanks	98 ~ 99
Canister-types 10~50Lit	98
Rack-types 47~95Lit	99

M

Magnetic Stirrers	309 ~ 312
120 × 120 Mini & 180 × 180 Midi	309
Multi-point, 4 & 8-Point	312
Premium Solid Ceramic Glass Plate	311
Standard Analog & Digital	310
Magnetic Stirrer Bars ⇒ See "Stirrer Bars"	330 ~ 332
Magnifiers	205
Fluorescent Light Magnifiers	205
LED Light Magnifiers	205
Manometers	210
Measuring Apparatus	206 ~ 229
Anemometers / Airflow Meters	206 ~ 207
Gas Detectors / Monitors	208 ~ 209
Manometers	210
Salinity Meter	210
Sound Level(Noise) Meters	211
Tachometers	212
Thermo-Hygrometers	213 ~ 217
Thermometers	218 ~ 229
Mills	230 ~ 232
Ball Mills	230 ~ 231
Balls, Alumina · Porcelain · STS · Zirconia	231
Mill Pots, HDPE & Porcelain	231
Cutting Mill & Pin Mill	232
Mixers	287 ~ 288
Tube Mixers	287
Vortex Mixers	288
Moisture Analyzers New	38 ~ 39

O

Ovens / Lab Drying	233 ~ 245
General Purpose	234 ~ 241
Forced-air R/T~250°C	234 ~ 237
Gravity-air R/T~230°C	238 ~ 241
High Clean Air 60°C~250°C	242
Industrial / Large-capacity R/T~250°C	243
Vacuum-Gas Exchangeable R/T~200°C	244 ~ 245
Pumps for Vacuum Ovens	250

P

PCR Workstation / Clean Benches	70 ~ 71
PC-type, 3-Side Viewing, Fabricated New	70
Steel-type	71
Pipettors, WisdPet™	246
Macropipettors	246
Micropipettors	246
Plasma Sterilizers	306
Pumps	247 ~ 250
Vacuum-Aspirators	247
Vacuum-Diaphragm	248 ~ 249
Vacuum-Oil Sealed	250

R

Reactor	251 ~ 268
BioReactor / Fermentor Coming Soon!	252 ~ 255
High Pressure Stainless-steel Reactors	256 ~ 258
Mini Autoclave Reactors of STS	256 ~ 257
Hydrothermal Synthesis Reactors of STS & PTFE	258
Vacuum Glass Reactors ReNew	259 ~ 266
Classic Round Flask Reactors, 0.25~50Lit	259

<Continued from Reactor>

Cylindrical Single Wall Reactors, 0.1~50Lit	260 ~ 261
Double Wall Reactors-Self Standing, 0.25~20Lit	262
Double Wall Reactors-Standard 1~100Lit	263 ~ 266
Autoflitting & Tilting, 20~100Lit	266
Tilting, 20~100Lit	265
Vacuum Stainless-steel Double Wall Reactors ReNew	267 ~ 268
Double Wall Reactors-Self Standing, 1~10Lit	267
Double Wall Reactors-Standard 3~100Lit	268

Reactor Circulators

External, Heated/Refrigerated	58 ~ 62, 269 ~ 271
-30°C/-20°C~200°C New	58 ~ 59, 270 ~ 271
External / Internal, Heated +250°C	60
External / Internal, Heated / Refrigerated	61 ~ 62, 119
-35°C/-25°C~150°C	61 ~ 62, 119
Accessories	63 ~ 65
Circulation Hoses New	63 ~ 64
Cooling / Heating Transfer Fluids	65, 155, 269

Refrigerators & Freezers

Freezers, Bio -35°C	281
Freezers, Ultra Low Temp	273 ~ 280
Energy Saving SMART ULT Freezer ReNew	273
Energy Saving ULT Freezing System	274
-95°C Upright-type ReNew	275
-86°C Personal Upright & Underbench	277
-86°C Upright-type ReNew	276
-45°C Upright-type	278
Accessories, Racks/Boxes/etc.	279 ~ 280
Refrigerators	282 ~ 285
0°C~10°C Laboratory Refrigerators	282
2~10°C Blood Bank Refrigerators	283
2~10°C Pharm. Refrigerators	284 ~ 285
Rockers	293
Compact Tube Rockers	293
Platform Rockers	293
Rotator	294

S

Safety Cabinets	74 ~ 78
Acids / Corrosives New	74 ~ 75
Flammables New	76 ~ 78
Vapor Purifier Box, Ductless Filtering	78
Salinity Meter	210
Sample / Bag Mixer	66
Shakers	286 ~ 295
Mixers	287 ~ 288
Tube Mixers	287
Vortex Mixers	288
Orbital Motion Shakers	289 ~ 292
Microplate Mixer	289
Orbital Shakers	290 ~ 292
Reciprocating Motion Shakers	291 ~ 292
Rockers	293
Compact Tube Rockers	293
Platform Rockers	293
Rotator	294
Sieve Shakers	295
Sieve	296 ~ 297
φ 200 × h50mm All Stainless-steel, Mesh Holes	296
φ 203 × h41mm Brass-Frame, STS Mesh Holes	297
Sieve Shakers	295
Electromagnetic 3D Movement	295
Vibrating	295
Sound Level(Noise) Meters	211
Sterilizers	298 ~ 307
Autoclaves / Steam Sterilizers	299 ~ 304
Bench-top Front Door	299
Front Door	304
Top-loading Vertical New	300 ~ 303
Dry Heat Sterilizers	305
Plasma Sterilizers	306
Accessories	303, 307
Autoclaving Deodorant	303
Sterilization Drapes New	307
Sterilization Indicators New	307
Stirrer	308 ~ 330
Magnetic Stirrers	309 ~ 312
120 × 120 Mini & 180 × 180 Midi	309
Multi-point, 4 & 8-Point	312
Premium Solid Ceramic Glass Plate	311
Standard Analog & Digital	310
Magnetic with Heating / Hotplate Stirrers	180 ~ 182, 313 ~ 315

<Continued from Stirrer>

380°C, Analog & Digital	181, 314
500°C, Solid Ceramic Glass	180, 313
Multi Point 3- & 6-Place, 350°C	182, 315
Overhead Stirrer	316 ~ 324
6-Paddle Stirrer / Jar tester	204, 324
High-Speed	317 ~ 319
High-Torque	320 ~ 323
Selection Guide	316
Overhead Stirrer Accessories	325 ~ 330
Flexible Couplings	326
Stands · Rods · Clamp Holder	330
Stir-Guides / Seals	325
Stirring Shafts / Impellers	326 ~ 329
Stirrer Bars	331 ~ 333
Cross-type	332
Elliptical-type	332
Micro	331
Octagonal + Oval Mixed Sets	332
Octagonal-type & Sets	331
Oval-type	332
Pivot Ring-type	331
Retrievers	333
Wing-type	333

T

Tachometers	212
Thermo-Hygrometers	213 ~ 217
Desktop / Wallmount-type New	214 ~ 216
Handy-type	213 ~ 214
Temp & Humidity Data Loggers	216
WBGT Monitors	217
Thermometers	218 ~ 229
Infra-Red	218 ~ 219
K-Type	220 ~ 222
NTC-Type	223 ~ 227
Pt100	228
Temp Data Loggers New	228 ~ 229
Timers / Stopwatches	334
3-/4- Multi Timers	334
Basic Digital Timers	334
Timer / Stopwatch	334

U

Ultrasonic Cleaners	335 ~ 337
Analog Cleaners	335
Digital Cleaners	336
Large Capacity, Remote Control	337
Ultrasonic Cleaning Detergents	337
Ultrasonic Homogenizer New	174, 338
UV Intensity Meter	109
UV Transilluminators	109

V

Vacuum Manifolds	339 ~ 340
Chromatography Solid-phase Extraction	340
Glass Manifolds, Glass Stopcock Type	339
Glass Manifolds, PTFE Needle Valve Type	340
STS Manifolds, Filter Holder Funnel Type	339
Viscosity	341 ~ 345
Capillary Viscometer Bath	55, 343
Capillary Viscometers New	341 ~ 342
Cannon-Fenske	341 ~ 342
Ubbelohde	341 ~ 342
Rotary Digital Viscometers	344 ~ 345
Multi-speed	345
Standard	344

W

Washers, Glass/Labware Washers	346 ~ 347
Water Purification	348 ~ 352
Classic Water Stills	348 ~ 349
Water Purification Systems	350 ~ 352
Reverse Osmosis (RO)	350
Ultra Pure (UP)	351
Ultra Pure(UP) & Pure(RO)	352
WBGT Monitors	217
Weight Calibration & Recalibration	43

Since 1980

**DAIHAN
SCIENTIFIC
GROUP**
Your Neighborhood Friend

DAIHAN-Sci.com Shopping-mall
ALLforLAB.com Open-market

2026

DAIHAN™ Lab-Safety

- ▶ Biological Safety Cabinet(BSC) P. 68
- ▶ Clean Bench / PCR Workstation **NEW** P. 70
- ▶ Ductless Filtering Fume Hood P. 135
- ▶ Ductless Filtering Reagent Storage Cabinet, Steel **NEW** ... P. 72
- ▶ Filtering Corrosives/Reagent Storage Cabinet, PP P. 73
- ▶ Safety Cabinet - Acids / Corrosives **NEW** P. 74
- Flammables **NEW** P. 76



Call Center ; 080-008-3000
International Tel ; +82-70-5020-9000

Fax ; 033-737-7575
Fax ; +82-2-963-5231

E-mail ; dh@daihansci.co.kr
E-mail ; info@daihan-sci.com

(1) 부피 질량 검교정 서비스

(2) Digi-Burets, Motorized

(3) Dispensers

(4) Pipets & Aids

(5) Pipettors

(6) Pipette Tips

“대한과학” 부피(ml)와 질량(g)

검교정서비스

with Best Price

검색

검교정.kr

Volume and Weight (Mass) Calibration & Recalibration Services



공인; KC23-433
KS Q ISO/IEC 17025

교정서비스 신청하기 >

Direct Call; 080-858-3002

E-mail; cal@dhsci.co.kr

부피(ml) 검교정; 모든 부피계

Burets, Cylinders, Dispensers, Flasks-Volumetric,
Pipets, Pipettors (Micro · Macro), Titrators, etc.



질량(g) 검교정; 모든 질량계

Balances, Scales,
Weight, etc.



“대한과학” 피펫터 & 디스펜서 판매

Pipetter Mall

검색

Pipettermall.com

Direct Call; 080-858-3002

E-mail; pipet@dhsci.co.kr

▶ Almost all Pipettors

; Brand™, DLAB™, Eppendorf™, Gilson™,
Microlit™, Sartorius™, ThermoFisher,
VITLAB™, WisdPet™, and WITOPET™

▶ Dispensers

; Brand™, Microlit™, Wisdose™,
and witeg LABMAX™



대한과학 교정센터 / 피펫터몰

강원특별자치도 원주시 지정면 신평석화로 326

Direct Call; 080-858-3002, E-mail; cal@dhsci.co.kr / pipet@dhsci.co.kr

대한과학그룹

LAB SERVICE

Total Laboratory Solution

• Lab-service.co.kr • Labas.kr

검색 랩서비스

검색 랩 A/S



제조사에 구매없는
모든 실험장비의 유지보수

One Stop Service



전문가에 의한
안전한 실험장비

이전 설치 서비스



완벽한 절차의
적격성 평가

IQ & OQ Service



연구실험실 안전진단 포함
장비의 정기점검 및 유지보수

Total Solution

Full Service Network! Covering the Whole Nation

서비스 본부

Direct Call ; 080-858-3001(수신자 부담)

E-mail ; service@daihansci.co.kr

전국 Service Network

광주 ; 062-525-1155 (광주과학)

대구 ; 053-381-8445 (대한과학대구)

부산 ; 051-247-2875 (대한과학부산)

전주 ; 063-277-8901 (대한과학전북)

진주 ; 055-745-8622 (선일과학)

제주 ; 064-721-7749 (금성하이텍)

청주 ; 043-234-6200 (충북과학)



대한과학그룹 Total Laboratory Solution

고객만족을 최우선 정책으로 생각하는 대한과학그룹이 부담없는 연구실 구축을 목표로
실험장비 종합렌탈 서비스와
중고실험기기 판매서비스를 통해 다시한번 나아갑니다.

LAB RENTAL

검색 대한렌탈

검색 랩렌탈

부담없는 사용료 & 최상의 품질
할부구매식 인수렌탈
지속적인 유지보수 관리 서비스
필요시사용 일시렌탈

Direct Call & E-mail.

☎ 080-858-3001 (수신자부담)
✉ rental@daihansci.co.kr

🌐 DAIHANRental.co.kr / LabRental.co.kr



DAIHAN Used

검색 중고실험기기

검색 대한중고

구입후 1년간 무상보증
효율적인 가격 & 세심한 성능진단, 최상의 품질
구매, 설치, 사용법 안내까지 한번에
< One Stop 서비스 >
(주)대한과학 기술영업팀의 빠르고 정확한 A/S

Direct Call & E-mail.

☎ 080-008-9009 (수신자부담)
✉ used@daihansci.co.kr

🌐 DAIHANUsed.com / UsedLab.kr



2026 DAIHAN Scientific Laboratory & Medical Instruments “Overview”

All DAIHAN™ Instruments have a **3-Year Global Warranty**

• Balance



Analytical

P.14



Lab Precision

P.21



Moisture Analyzer

P.38

• Bath



Basic Bath, 100℃~250℃

P.48



Circulation & Shaking Water Bath

P.49

• Bath



Cold Trap Bath

P.51



Extraction Water Bath

P.52



Viscosity Bath

P.55

• Bath Circulator



External, Chiller -30℃ & -20℃

P.56



External, Heated / Refrigerated

P.58

• Bath Circulator



External / Internal, Heated / Refrigerated

P.61

• Blender



Bag / Sample Mixer

P.66

• Cabinet & Enclosure



Biological Safety Cabinet (BSC) Class II Type A2

P.68



Clean Bench / PCR Workstation

P.70



Filtering Reagent Storage Cabinet

P.72

• Cabinet & Enclosure



Safety Cabinet

P.74

• Cart



ABS and Stainless-steel Cart

P.79

• Centrifuge



0.2~50ml Micro

P.87



Clinical Centrifuge

P.88



Refrigerated, Hi-Volume

P.94

• Cryogenic



Liquid Nitrogen (LN₂) Storage System

P.98

• Desiccator Cabinet



Desiccator Cabinet

P.100


• Dispenser



Wisdose™ Bottle-top Standard & HF

P.104

• Dryer



Labware Dryer, 39 Places

P.106

• Drying Chamber



Drying Chamber

P.105

<Continued on next page>

2026 DAIHAN Scientific Laboratory & Medical Instruments “Overview”

All DAIHAN™ Instruments have a **3-Year Global Warranty**

• Electrophoresis



Built-in Transilluminator
Gel Documentation System
P.108



UV Transilluminator
P.109

• Evaporator



Rotary Vacuum Evaporator
P.110

• Food Analyzer



New
Fat Analyzer
P.126



New
Fiber Analyzer
P.127

• Food Analyzer



New
Kjeldahl Digestion
P.128



New
Kjeldahl Nitrogen Analyzer
P.129

• Freeze Dry System



Freeze Dryer, 3 & 6Lit / day
P.131

• Fume Hood



Vapors Safeguard
Ductless Filtered Fume Hood
P.135



Multiuse
Acids/Corrosives
Ducted Fume Hood
P.137

• Fume Hood



P.139
Fume Purifier

• Furnace



P.142
Dental Furnace,
1100 °C / 1650 °C




Up to 1650°C
Muffle Furnace,
1000 °C~1650 °C
P.143

• Furniture



P.148
Assembly-type
Laboratory Furniture

• Glove Box



P.152
PMMA Vacuum
Glove Box

• Heating Block / Dry Bath Incubator



P.156
Heating Block




P.158
Heating · Cooling · Mixing

• Heating Mantle



P.160
Heating Mantle



P.168
Multi Heating Mantle

• Homogenizer



P.171
Analog & Digital

• Homogenizer



New
P.174
Ultrasonic Homogenizer

• Hotplate



P.176
550 °C High Temp.
Hi-Light Hotplate



P.177
500 °C & 380 °C Hotplate



Remote Control
Built-in Control
P.179
350 °C Large Hotplate

• Hotplate Stirrer



P.180
500 °C & 380 °C
Hotplate Stirrer

<Continued on next page>

2026 DAIHAN Scientific Laboratory & Medical Instruments “Overview”

All DAIHAN™ Instruments have a **3-Year Global Warranty**

• Hotplate Stirrer



P.182

Multi Hotplate Stirrer
3- / 6-Place

• Ice Maker



P.183

Ice Maker
Snow & Flake Type

• Incubator



P.186

CO₂ Incubator



P.187

General Mini Low
Temp / Shaking, 20L



P.188

General Purpose
Forced- / Gravity-Air

• Incubator



P.195

Humidity / Temp.
Chamber



P.197

Low Temp (B.O.D)
Incubator



P.199

Plant Growth Chamber



P.202

Shaking Incubator



P.204

6-Paddle Stirrer

• Magnifier



P.205

Light Magnifier
Fluorescent / LED



P.206

Anemometer and
Gas Detector



P.210

Manometer and
Sound Level Meter



P.213

Thermo-Hygrometer,
WBGT Monitor



P.218

Thermometer, Infra-Red ·
K-type · NTC-type · Pt100

• Mill



P.230

Ball Mill



P.232

Cutting Mill & Pin Mill

• Oven



P.234

General Oven
Forced- / Gravity-Air



P.242

High Clean Air Oven



P.243

Industrial Forced-Air
Oven

• Oven



P.244

Vacuum Oven

• Pipetter



P.246

WisdPet™ Adjustable
Micro & Macro pipetter

• Pump



P.248

Vacuum Pump

• Reactor



P.252

BioReactor / Fermentor



P.256

Pressure ; Stainless-steel

<Continued on next page>

2026 DAIHAN Scientific Laboratory & Medical Instruments “Overview”

All DAIHAN™ Instruments have a **3-Year Global Warranty**

● Reactor



P.265

Vacuum ; Glass & Stainless steel

● Refrigerator & ULT Freezer



P.272

Energy Saving & Smart Control
INVERTER COMPRESSOR SYSTEM™

● Shaker



P.287


Mixer, Tube Mixer



P.288

Mixer, Vortex Mixer

● Shaker



P.289

Orbital Microplate Mixer



P.291

Orbital / Reciprocating Shaker



P.293

Rocker



P.294

Rotator



P.295

Sieve Shaker

● Sterilizer



P.299

Autoclave, 25 & 40Lit
Bench-top Front Door



P.300

Autoclave, 47~100Lit
Toploading / Vertical



P.304

Autoclave, 90Lit
Front Door



P.305

Dry Heat Sterilizer
17 · 32 · 50Lit



P.306

Plasma Sterilizer
18Lit

● Stirrer, Magnetic



P.309

Magnetic, Digital
Mini & Midi



P.311

Magnetic,
High-Performance



P.312

Magnetic, Multi-Point

● Stirrer, Overhead



P.317

Overhead,
Hi-Speed & Hi-Torque



P.321

Overhead, Hi-Torque Ultra
~350,000mPas

● Ultrasonic Cleaner



P.335

Ultrasonic Cleaner



P.344

Rotary Viscometer

● Washer



P.346

Glass / Labware Washer

● Water Purification



P.349

Classic Water Still



P.352

Water Purification System
UP & RO

DAIHAN Scientific

BALANCES

Balances	P. 14~29, 40
Analytical, Semi-Micro	P. 14~15
Analytical	P. 15~19
Lab Precision	P. 20~25
Lab-General	P. 26~29
Pocket	P. 40
Medical / Healthcare Scales	P. 35~37
Moisture Analyzers	P. 38~39
Scales	P. 30~35
Weight	P. 41~43
Accessories - Weighing Bottles, Dishes, Spatulas &c.	P. 43~46



Analytical Balance
"WBA-220A"
 [d] 0.1mg, max 220g



Lab Precision
"BAL-310DS"
 [d] 1mg, max 310g



Moisture Analyzer
"BAL-110M"
 [d] 1mg, max 110g

질량(g) 검교정 ; 모든 질량계
Weight Calibration & Recalibration

“대한과학은 KOLAS 인증을 획득한
 국제공인 부피 및 질량 교정 기관입니다.”

KS Q ISO/IEC 17025
 공인 ; KC23-433

대한과학 교정센터
 Direct Call : 080-858-3002
 E-mail : call@dhsc.co.kr

Kern® Premium [d] 0.01mg, max 220g Hi-professional Analytical Semi-Micro Balance "ABP5M", with the Latest Single-Cell Generation Extremely Rapid, Internal Calibration, Bright OLED Display, with Multi-function Weighing Plate, GLP/ISO Record Keeping

고급형 고정밀 세미 마이크로-분석/화학천평, 내부자동보정

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

Single-cell advanced technology

- Fully Automatic Manufactured Weighing Cell from One Piece of Material
- Stable Temperature Behavior
- Short Stabilization Time ; steady weight values within Approx. 2 s! (models with [d] = 0,1 mg), approx. 8 s! (models with [d] = 0,01 mg)
- Shock Proof Construction
- High Corner Load Performance

Specifications

- Luminescent OLED Display ;** Height 12 mm, Bright with High Contrast, for Easy Reading of the Weight, even in Poor Lighting Conditions
- Weighing Plate Dimensions ;** $\phi 91$ mm
- Overall Dimensions ;** 213×407×344 mm (with Draught Shield)
- Net Weight ;** approx. 8 kg
- Ambient Temperature ;** 10°C~30°C

"Auto Int-CAL" "ABP.5M"

3-Year Warranty



Shown with Min. 0.00001g



A incl. Multi-function Weighing Plate

Features



Internal Calibration



Data Interface RS 232



USB Data Interface PC or Printer



GLP/ISO record keeping



Piece Counting



Net-total Weighing



% Determination



Weighing Unit



Animal Weighing Program ; Strong Vibrations are Filtered Out.



Suspended Weighing



Universal Adapter



Single Cell Tech

* **Navigation Pad** for Super Quick Navigating through the Menus

* The Minimum Weight of Sample can be manually Stored in the Device or Automatically Calculated.

* **Simple Recipe Weighing and Documenting** with a Combined Tare/Print Function. In addition, the Ingredients for the Recipe are Numbered Automatically and Printed out with Their Corresponding Number and Nominal Weight

* **Individual User Settings for up to 10 users can be saved.** * Protective working cover included with delivery.

* **U.S. FDA 21 Part 11 ;** assists you in data integrity in accordance with U.S. Pat. FDA 21 Part 11

* Automatic data output to the PC/printer each time the balance is steady

* **A Multi-function weighing plate included ;** minimises the effect of currents of air in the weighing space and therefore significantly improves the stabilization time and repeatability. In addition samples, sample paper, PCR containers, micro centrifuge tubes and many other items which protrude can be easily fixed in place and weighed easily

Specifications & Ordering information

<L1>

Cat. No.	Weighing -range Max g	Read -out d mg	Verif. value e mg	Reproducibility mg	Linearity mg	Price(W)
KE.ABP100.5M	135	0.01	1	0.05	± 0.2	
KE.ABP200.5M	220	0.01	1	0.05	± 0.2	
KE.ABP200.5DM	102/220	0.01/0.1	1	0.05/0.1	± 0.1/0.2	

- KE.ABPA01 (1) **Optional Draft Shield Rear Panel with Integrated Ionizer** to neutralise electrostatic charge
- KE.YKS.01 (2) **Optional Statistics Printer**, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50 mm, Adapter 230V/50Hz
- KE.YKB.01 (3) **Optional Thermal Printer**, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165×140×50 mm, Adapter 230V/50Hz
- KE.MLB.A05 **RS-232 Interface Cable**, PC or Thermal Printer, for "YKS.01", "YKB.01", 1.5m

Ea. Price, ₩

(1) **Optional, KE.ABPA01**
Draft Shield Rear Panel
with Integrated Ionizer
정전기 방지(중화)용



(2) Statistics Printer
KERN YKS.01



(3) Thermal Printer
KERN YKB.01

$\phi 75 \cdot 90 \cdot 102$ mm Reusable Aluminum Weighing Dish
Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

■ Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005 (1) $\phi 75$ Alumi Weigh-Dish, $\phi 75 \times h8$ mm, 4.5g, 30ml
- DH.Dis8006 (2) $\phi 90$ Alumi Weigh-Dish, $\phi 90 \times h4$ mm, 3.05g, 25ml **New**
- KE.MLB.A01 , $\phi 90 \times h7$ mm, 3.7g, 40ml
- DH.Dis8007 (3) $\phi 102$ Alumi Weigh-Dish, $\phi 102 \times h8$ mm, 3.8g, 57ml

<L3>

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	



(2) $\phi 90$ mm Reusable Al. Dish

Articles

Cat. No

Description

Kern® [d] 0.01/0.1mg, max 220g Semi-Micro & Precise Analytical Balance "ABT", Single-Cell Weighing System
With Internal Calibration, Multi-Function ; ex. Density-Measure, 고정밀-분석/화학천평, 내부자동보정

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

Features

Single-cell Advanced Technology

- Fully Automatic Manufacture Weighing Cell from One Piece of Material
- Stable Temperature Behavior
- Short Stabilization Time
- Shock Proof Construction
- High Corner Load Performance

Specifications

- "Int-CAL" ; Automatic Internal Adjustment
- Fast Display** ; Short Stabilization Time
- Large LCD Display** ; Digit Height 14mm
- Identification Number** ; 4 Digits, Printed on Calibration Protocol freely Programmable
- Weighing Range** ; 0.01mg~101g (ABT100.5M Model)
0.01/0.1mg ~ 82/220g (ABT220.5DM Model)
0.1mg ~ 120g (ABT120.4M Model)
0.1mg ~ 220g (ABT220.4M Model)
- Overall Dimensions** ; 217 × 356 × h338mm
- Weighing Plate** ; ϕ 80 mm
- Weighing Space** ; 168 × 172 × h223 mm
- Net Weight** ; approx. 7kg
- Permissible Ambient Temperature** ; 10°C / 30°C

"Auto Int-CAL" "ABT"

3-Year Warranty



Large LCD Display
Digit Height 14 mm and
Capacity Display



Internal Calibration



Data Interface RS 232



GLP/ISO record keeping



Piece Counting



Net-total Weighing



% Determination



Weighing Unit



Suspended Weighing



AC 230V Adapter



Single Cell Tech

ALL for the Best BALANCES
-- via "KERN" since 1844

* **Simple Adding Function** with Auto-tare and Automatic Subsequent Data Output after each Weighing Procedure

* **Dosing assistance** High-Stability Mode and other Filter Settings can be Selected

* **Simple Recipe-weighing/Documenting** with a Combined Tare/Print Function. In addition, the Mixing Ingredients for the Recipe are Numbered Automatically and Printed out with its Corresponding Number and Weight Value

* **Identification number** of Balance is Printed on Calibration Protocol

Specifications & Ordering information


<L1>

Cat. No.	Weighing -Range Max g	Read -out d mg	Verif. value e mg	Reproducibility mg	Linearity mg	Price(₩)
1) 0.01mg, max 200g Semi-micro						
KE.ABT100.5M New	101	0.01	1	0.05	± 0.15	
KE.ABT220.5DM	82 220	0.01 0.1	1	0.05 0.1	± 0.1 ± 0.2	
2) 0.1mg, max 200g Hi-Quality						
KE.ABT120.4M	120	0.1	1	0.1	± 0.2	
KE.ABT220.4M	220	0.1	1	0.1	± 0.2	

Optional Accessories

- KE.YKS.01 (1) **Statistics Printer**, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165 × 140 × 50 mm, Adapter 230V/50Hz
- KE.YKB.01 (2) **Thermal Printer**, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165 × 140 × 50 mm, Adapter 230V/50Hz
- KE.YDB.03 (3) **Universal Density Determination Set**, Including Beaker, Thermometer, Sample Holder, Adapter, Equalizing Weights, Volume Plummet (STS, 20g), Universal Submersible Basket for Solid Bodies which Float and Sink
- KE.770.926 **RS-232 Interface Cable**, PC or Thermal Printer "YKS.01", "YKB.01", 1.5m for ABS-/ABJ-/ABT (Year of Production 2000~2012)- Model

Ea. Price, ₩



(1)
Statistics Printer
KERN YKS.01

(2)
Thermal Printer
KERN YKB.01



(3) Set for Density Determination of Liquids and Solids, Density is Indicated directly on the Display KERN YDB.03

ϕ 75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

■ Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005 (1) ϕ 75 Alumi Weigh-Dish, ϕ 75 × h8mm, 4.5g, 30ml
- DH.Dis8006 (2) ϕ 90 Alumi Weigh-Dish, ϕ 90 × h4mm, 3.05g, 25ml **New**
- KE.MLB.A01 , ϕ 90 × h7mm, 3.7g, 40ml
- DH.Dis8007 (3) ϕ 102 Alumi Weigh-Dish, ϕ 102 × h8mm, 3.8g, 57ml

<L3>

Pk. Price, ₩

100ea
50ea
80ea
50ea



(2) ϕ 90mm Reusable Al. Dish

DAIHAN™ [d] 0.1mg, max 220g Calibration Certified Standard Analytical Balance, ϕ 80 · 90mm Weighing Plate

A. Ext-CAL "WBA-220", B. Auto Int-CAL "WBA-220A", with Glass Draft Shield, Back Light LCD, Counting Function, Various Weight Mode

"Ext-CAL 외부 보정형" & "Int-CAL 자동 보정형" 표준 분석/화학 천평, 계수계 겸용, 작업보호커버 포함, 보정 성적서 포함

"Ext-CAL" "WBA-220"
[d] 0.1mg, max 220g



Calibration
Certificate 포함



3-Year Warranty

A. Ext-CAL "WBA-220"
with Calibration Certificate

Includes 200g
Calibration Weight



- Supplied by DAIHAN Korea
- CE Certified
- With Calibration Certificate(보정 성적서 포함)
- PL(Product Liability) Insurance
- 3-Years Warranty
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Supplied with ; - ϕ 80/90mm Weighing Plate and Glass Draft Shield
- AC 230V Adapter for DC Power 15V / 500mA(WBA-220A)
- 200g Calibration Weight (Ext-CAL Model, WBA-220)

* Various Weight Mode (g/ct/%/GN/dwt/etc)

* Numeric Linearity Calibration

* Appearance Parameters

* Electro Magnetic Base

* Weighing Response Time Adjustable

* Built in RS232C for Printer

Specifications & Ordering information

<L1>

Spec.	Model	A. "WBA-220" Ext-CAL	B. "WBA-220A" Auto Int-CAL
Read-Out [d]. mg		0.1mg	
Weighing Range. g		max 220g	
Weighing Plate. mm		ϕ 80mm	ϕ 90mm
Reproducibility. mg		0.2mg	
Linearity. mg		0.3mg	0.2mg
Ambient Temperature		5.0°C~35°C	
Stabilization Time		≤ 3 Seconds (Typical at Fast)	
Weight Mode		g / ct	g / ct / % / GN / dwt / etc.
Counting Mode		Counting & Weighing Mode Available	
Calibration		External Calibration, Included 200g Calibration Weight	Internal Calibration, AUTO
Display		Digital LCD with Back-Light Function	
Dimension & Net Weight		190 × 240 × h265mm / 6.3kg	223 × 345 × h331mm / 5.7kg
Packing Size & Gross Weight		290 × 420 × h490mm / 8.3kg	390 × 510 × h520mm / 10.1kg
Electrical Reqts.		AC 220V, 50 / 60Hz	DC 15V, with AC 230V Adapter, 50/60Hz
Cat. No.		DH.WBA0220	DH.WBA0220A
Price, Each. ₩			

"Auto Int-CAL" "WBA-220A"
[d] 0.1mg, max 220g



Calibration
Certificate 포함



3-Year Warranty

B. Int-CAL "WBA-220A"
with Calibration Certificate

Optional Thermal Printers ; Specify the Balance Model for Printer Installation

* 프린터 구매 시 바란스 모델을 알려주어야 합니다.

Ea. Price, ₩

SL.WBA.100P (1) Thermal Printer, Paper Roll, "WBA-100"
Included Cable Module, 1 × Paper Roll

DH.WBA0084 Thermal Paper Package, 5 × Paper Roll, for "WBA-100"

SL.WBA.420P (2) Thermal Printer, Label Paper Roll, "WBA-420"
Included Cable Module, 1 × Label Paper Roll

DH.WBA0771 Thermal Label Paper Package, 5 × Label Paper Roll, for "WBA-420"



(1) "WBA-100" RS-232C
Thermal Printer



(2) "WBA-420" RS-232C
Thermal Label Printer

Kern® [d] 0.1mg, max 210g Popular Analytical Balance Ext-CAL "ADB" & Auto Int-CAL "ADJ"

With Glass Draft Shield, Backlit LCD, 인기형 표준분석/화학천평, 외부 보정형 & 자동 보정형

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

Features

Specification

- Large backlit LCD Display ; Digit Height 16 mm
- Weighing Plate Dimensions ; ϕ 90 mm
- Overall Dimensions ; 230×310×330 mm (with Draught Shield)
- Net Weight ; approx. 4.4 kg
- Ambient Temperature ; 10°C~30°C

"Ext-CAL" "ADB"
"Auto Int-CAL" "ADJ"

NEW



3-Year Warranty



External Calibration



Internal Calibration



Data Interface RS 232



GLP/ISO record keeping



Piece Counting



% Determination



Weighing Unit



Universal Adapter



Electromagnetic Force Compensation

* ADB ; Adjusting Program CAL for Quick Setting of the Balance Accuracy Using an External Test Weight.

* ADJ ; Automatic Internal Adjustment in the Case of a Change in Temperature $\geq 2^\circ\text{C}$ or Timecontrolled Every 3 Hours

* Level Indicator and Foot Screws to Level the Balance Precisely, Fitted as Standard, to give the Most Accurate Weighing Results.

* Large Glass Draught Shield with 3 Sliding Doors for Easy Access to the Items being Weighed.

* Compact Size, Practical for Small Spaces. * Simple and Convenient 6-key Operation.

Specifications & Ordering information

<L1>

Cat. No.	Weighing -Range Max g	Read -out d mg	Reproducibility mg	Linearity mg	Weighing Plate	Price(W)
1) Ext-CAL						
KE.ADB100.4	120	0.1	0.2	± 0.4	ϕ 90 mm	
KE.ADB200.4	210	0.1	0.2	± 0.4	ϕ 90 mm	
2) Auto Int-CAL						
KE.ADJ200.4	210	0.1	0.2	± 0.4	ϕ 90 mm	

Optional Accessories

- KE.YKS.01 (1) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50mm, Adapter 230V/50Hz
- KE.YKB.01 (2) Thermal Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165×140×50mm, Adapter 230V/50Hz
- KE.MLB.A05 RS-232 Interface Cable, PC or Thermal Printer, for "YKS.01" / "YKB.01", 1.5m

Ea. Price, ₩

(1)
Statistics Printer
KERN YKS.01(2)
Thermal Printer
KERN YKB.01ALL for the Best BALANCES
--- via "KERN" since 1844 ϕ 75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

<L3>

■ Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005 (1) ϕ 75 Alumi Weigh-Dish, ϕ 75×h8mm, 4.5g, 30ml
- DH.Dis8006 (2) ϕ 90 Alumi Weigh-Dish, ϕ 90×h4mm, 3.05g, 25ml **New**
- KE.MLB.A01 , ϕ 90×h7mm, 3.7g, 40ml
- DH.Dis8007 (3) ϕ 102 Alumi Weigh-Dish, ϕ 102×h8mm, 3.8g, 57ml

Pk. Price, ₩

100ea

50ea

80ea

50ea

(2) ϕ 90mm Reusable Al. Dish

Kern® [d] 0.1mg, max 320g High-Standard Analytical Balance Ext-CAL "ACS" & Auto Int-CAL "ACJ"

Multi-function ; ex. Density-measure, and Glass Draft Shield, 표준분석/화학천평, 외부 보정형 & 자동 보정형

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

Single-cell Advanced Technology

- Fully Automatic Manufacture Weighing Cell from One Piece of Material
- Stable Temperature Behavior
- Short Stabilization Time
- Shock Proof Construction
- High Corner Load Performance

Specifications

- ACS** ; Adjusting Program CAL for Quick Setting of the Balance Accuracy Using an External Test Weight.
- ACJ** ; Automatic Internal Adjustment in the Case of a Temperature Change > 2°C and Time-controlled Every 4 Hours
- Simple Recipe Weighing and Documenting with a Combined Tare / Print Function
- Automatic Data Output** to the PC / Printer
- Weighing Range** ; 0.1mg ~ 120 / 220 / 320g
- Overall Dimensions** ; 210 × 340 × h325mm
- Weighing Plate** ; $\phi 91$ mm Stainless-steel
- Weighing Space** ; 174 × 162 × h227mm
- Net Weight** ; Approx. 6 kg
- Permissible Ambient Temperature** ; 10°C / 30°C

"Ext-CAL" "ACS"

"Auto Int-CAL" "ACJ"

NEW

3-Year Warranty

CE



Capacity Display
An Ascending Bar Graph Displays the Remaining Available Weighing Range

Features

- External Calibration**
- Internal Calibration**
- GLP/ISO record keeping**
- Piece Counting**
- Net-total Weighing**
- % Determination**
- Weighing Unit**
- Weighing with Tolerance Range**
- AC 230V Adapter**
- Single Cell Tech**

Carat switch-over from g to ct is possible.

Model	d	Max
KERN	ct	ct
ACJ 100.4M	0.001	600
ACJ 200.4M	0.001	1100
ACJ 300.4M	0.001	1600

The carat display can also be verified.

* **Metal Housing** robust and Sturdy

* **Large LCD Display** Digit height 14mm

* **Identification number** 4 numeric Positions, 2 key Programming method

Specifications & Ordering information

<L1>


Cat. No.	Weighing Range Max g	Read-out d mg	Verif. value e mg	Reproducibility mg	Linearity mg	Price(W)
1) Ext-CAL						
KE.ACS100.4M	120	0.1	-	0.2	± 0.3	
KE.ACS200.4M	220	0.1	-	0.2	± 0.3	
KE.ACS300.4M	320	0.1	-	0.2	± 0.3	
2) Auto Int-CAL						
KE.ACJ100.4M	120	0.1	1	0.2	± 0.3	
KE.ACJ200.4M	220	0.1	1	0.2	± 0.3	
KE.ACJ300.4M	320	0.1	1	0.2	± 0.3	

Optional Accessories

- KE.YKS.01** (1) **Statistics Printer**, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165 × 140 × 50mm, Adapter 230V/50Hz
- KE.YKB.01** (2) **Thermal Printer**, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165 × 140 × 50mm, Adapter 230V/50Hz
- KE.YDB.03** (3) **Universal Density Determination Set**
Including Beaker, Thermometer, Sample Holder, Adapter, Equalizing Weights, Volume Plummets (stainless-steel, 20g), Universal Submersible Basket for Solid Bodies which Float and Sink
- KE.ACS.A01** **RS-232 Interface Cable**, PC or Statistics Printer "YKS.01", "YKB.01", 1.5m for ACS / ACJ Model

Ea. Price, ₩

ALL for the Best BALANCES
--- via "KERN" since 1844



$\phi 75 \cdot 90 \cdot 102$ mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평판 디쉬

■ Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005** (1) $\phi 75$ Alumi Weigh-Dish, $\phi 75 \times h8$ mm, 4.5g, 30ml
- DH.Dis8006** (2) $\phi 90$ Alumi Weigh-Dish, $\phi 90 \times h4$ mm, 3.05g, 25ml **NEW**
- KE.MLB.A01** , $\phi 90 \times h7$ mm, 3.7g, 40ml
- DH.Dis8007** (3) $\phi 102$ Alumi Weigh-Dish, $\phi 102 \times h8$ mm, 3.8g, 57ml

<L3>

Pk. Price, ₩

100ea
50ea
80ea
50ea



(2) $\phi 90$ mm Reusable Al. Dish

Kern® Premium [d] 0.1mg, max 320g Hi-professional Analytical Balance "ABP.4M", with the Latest Single-Cell Generation
Extremely Rapid, Internal Calibration, Bright OLED Display, with GLP/ISO Record Keeping
고급형-분석/화학천평, 내부자동보정

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

Features



Internal Calibration



Data Interface RS 232



USB Data Interface
PC or Printer



GLP/ISO record keeping



Piece Counting



Net-total Weighing



% Determination



Weighing Unit



Animal Weighing Program ;
Strong Vibrations are Filtered Out.



Suspended Weighing



Universal Adapter



Single Cell Tech

Single-cell advanced technology

- Fully Automatic Manufactured Weighing Cell from One Piece of Material
- Stable Temperature Behavior
- Short Stabilization Time ; steady weight values within Approx. 2 s! (models with [d] = 0,1 mg).
- Shock Proof Construction
- High Corner Load Performance

Specifications

- Luminescent OLED Display ;** Height 12 mm, Bright with High Contrast, for Easy Reading of the Weight, even in Poor Lighting Conditions
- Weighing Plate Dimensions ;** ϕ 91 mm
- Overall Dimensions ;** 213×407×344 mm (with Draught Shield)
- Net Weight ;** approx. 8 kg
- Ambient Temperature ;** 10°C~30°C

"Auto Int-CAL" "ABP.4M"

3-Year Warranty



* **Navigation Pad** for Super Quick Navigating through the Menus

* The Minimum Weight of Sample can be manually Stored in the Device or Automatically Calculated.

* **Simple Recipe Weighing and Documenting** with a Combined Tare/Print Function. In addition, the Ingredients for the Recipe are Numbered Automatically and Printed out with Their Corresponding Number and Nominal Weight

* **Individual User Settings for up to 10 users can be saved.**

* **U.S. FDA 21 Part 11 ;** assists you in data integrity in accordance with U.S. Pat. FDA 21 Part 11

* Automatic data output to the PC/printer each time the balance is steady

* Protective working cover included with delivery.

Specifications & Ordering information

<L1>

Cat. No.	Weighing -range Max g	Read -out d mg	Verif. value e mg	Reproducibility mg	Linearity mg	Price(₩)
KE.ABP100.4M	120	0.1	1	0.1	±0.2	
KE.ABP200.4M	220	0.1	1	0.1	±0.2	
KE.ABP300.4M	320	0.1	1	0.2	±0.3	

Optional Accessories

Ea.	Price, ₩
KE.ABPA.01 (1) Draft Shield Rear Panel with Integrated Ionizer to neutralise electrostatic charge	
KE.YKS.01 (2) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50 mm, Adapter 230V/50Hz	
KE.YKB.01 (3) Thermal Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165×140×50 mm, Adapter 230V/50Hz	
KE.MLB.A05 RS-232 Interface Cable, PC or Thermal Printer, for "YKS.01", "YKB.01", 1.5m	



(1) Optional, KE.ABPA.01
Draft Shield Rear Panel
with Integrated Ionizer
정전기 방지(중화)용



(2)
Statistics Printer
KERN YKS.01



(3)
Thermal Printer
KERN YKB.01

ϕ 75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

■ Ideal for General Weighing and Moisture Drying Application.

Pk.	Price, ₩
DH.Dis8005 (1) ϕ 75 Alumi Weigh-Dish, ϕ 75×h8mm, 4.5g, 30ml	100ea
DH.Dis8006 (2) ϕ 90 Alumi Weigh-Dish, ϕ 90×h4mm, 3.05g, 25ml New	50ea
KE.MLB.A01 , ϕ 90×h7mm, 3.7g, 40ml	80ea
DH.Dis8007 (3) ϕ 102 Alumi Weigh-Dish, ϕ 102×h8mm, 3.8g, 57ml	50ea



(2) ϕ 90mm Reusable Al. Dish

DAIHAN™ [d] 1mg, max 320 & 620g Calibration Certificated Hi-Precision Lab Balance, $\phi 90 \cdot 110 \cdot 128$ mm Weighing Plate
Ext-CAL "WBA-320 & 620", Auto Int-CAL "WBA-620A", with Glass Draft Shield, Super Size Back Light LCD, Counting Function, Various Weight Mode
"Ext-CAL 외부 보정형" & "Int-CAL 자동 보정형" 다기능 랩 바런스, g · 계수 · % 측정, 작업보호커버 포함, 보정 성적서 포함

"Ext-CAL" "WBA-320"
[d] 1mg, max 320g / 620g



Calibration
Certificate 포함



3-Year Warranty

Includes 200g or 500g
Calibration Weight

A. Ext-CAL "WBA-320" & B. "WBA-620"
with Calibration Certificate

"Auto Int-CAL" "WBA-620A"
[d] 1mg, max 620g



Calibration
Certificate 포함



3-Year Warranty

C. Int-CAL "WBA-620A"
with Calibration Certificate

- Supplied by DAIHAN Korea ■ CE Certified ■ PL Insurance ■ 3-Years Warranty
- With Calibration Certificate(보정 성적서 포함), Numeric Linearity Calibration
- With Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-Information, and Traceable Data Base System
- Supplied with ; - Glass Draft Shield
- 200g or 500g Calibration Weight (WBA-320 ; 200g, WBA-620 ; 500g)
- AC Adapter for DC Power 12V/500mA(WBA-320), 15V/300mA(WBA-620/620A)
- * Various Weight Mode (g/ct/GN/dwt/etc.) and Percent(%) Weighing Mode
- * Numeric Linearity Calibration * Electro Magnetic Base * Display Present Temperature
- * Four Level Anti-shock Adjustable Function * Time & Date Setting for Recording Data
- * Stabilization Time(Typical at Fast) ≤ 3 Seconds * Built in RS232C port for Printer
- * Default and Auto Lock Present Working Mode * Weighing Response Time Adjustable

Specifications & Ordering information

<L1>

Spec. \ Model	A. “WBA-320” Ext-CAL, 320g	B. “WBA-620” Ext-CAL, 620g	C. “WBA-620A” Auto Int-CAL, 620g
Read-Out [d]. g	0.001g		
Weighing Range. g	max 320g	max 620g	
Weighing Plate. mm	φ 90mm	φ 128mm	φ 110mm
Reproducibility. g	0.001g	0.002g	
Linearity. g	0.002g	0.003g	
Ambient Temp	5.0℃~35℃		
Weight Mode	% Mode and g-/ct-/GN-/dwt-/etc.		
Counting Mode	Counting & Weighing Mode Available		
Calibration	External Calibration		Internal Calibration
	Included 200g Calibration Weight	Included 500g Calibration Weight	AUTO
Display	Digital Super Size LCD with Back-Light Function		
Dimension & Net Weight	223×345×h331mm 3.6kg	223×345×h331mm 5.7kg	
Packing Size & Gross Weight	290×385×h365mm 5.1kg	390×510×h520mm 10.1kg	
Electrical Reqts.	DC 12V, with AC 230V Adapter, 50/60Hz	DC 15V, with AC 230V Adapter, 50/60Hz	
Cat. No.	DH.WBA0320	DH.WBA0620	DH.WBA0620A
Price, Each. ₩			

Optional Thermal Printers ; Specify the Balance Model for Printer Installation

* 프린터 구매 시 바런스 모델을 알려주어야 합니다.

Ea. Price, ₩

- SL.WBA.100P (1) Thermal Printer, Paper Roll, "WBA-100"
Included Cable Module, 1×Paper Roll
- DH.WBA0084 Thermal Paper Package, 5×Paper Roll, for "WBA-100"
- SL.WBA.420P (2) Thermal Printer, Label Paper Roll, "WBA-420"
Included Cable Module, 1×Label Paper Roll
- DH.WBA0771 Thermal Label Paper Package, 5×Label Paper Roll, for "WBA-420"



(1) "WBA-100" RS-232C
Thermal Printer



(2) "WBA-420" RS-232C
Thermal Label Printer



(2) $\phi 90$ mm Reusable Al. Dish

$\phi 75 \cdot 90 \cdot 102$ mm Reusable Aluminum Weighing Dish Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

<L3>

■ Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005 (1) $\phi 75$ Alumi Weigh-Dish, $\phi 75 \times h8$ mm, 4.5g, 30ml
- DH.Dis8006 (2) $\phi 90$ Alumi Weigh-Dish, $\phi 90 \times h4$ mm, 3.05g, 25ml **New**
- KE.MLB.A01 , $\phi 90 \times h7$ mm, 3.7g, 40ml
- DH.Dis8007 (3) $\phi 102$ Alumi Weigh-Dish, $\phi 102 \times h8$ mm, 3.8g, 57ml

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	

Articles

Cat. No

Description

DAIHAN™ [d] 1 & 10mg, max 310 & 3,100g High-Precision Lab Balance "BAL-D", RS-232C for Printer

With Counting Function, Various Weight Mode, DC Adaptor, Standard- & Shield-type, 정밀 랩 바란스, 계수계 겸용

NEW



3-Year Warranty

(1) Standard-type "BAL-3100D"
[d] 10mg, max 3100g

3-Year Warranty

(2) PC Shield-type "BAL-310DS"
[d] 1mg, max 310g(3) "WBA-100" RS-232C
Thermal Printer(4) "WBA-420" RS-232C
Thermal Label Printer

- Supplied by DAIHAN Korea ■ 3-Year Warranty
- Digital Large LCD Display and Ergonomic Design
- Easy Operation Weighing Mode ; g, PCS, %, CT

- * Material ; Stainless-steel Weighing Plate and ABS Resin Body
- * Stabilization Time ; 3sec
- * Permissible Ambient Temp 5℃~35℃
- * RS-232C for Printer and Level-Gauge integrated
- * Rechargeable Battery integrated, and with AC 230V Adapter for DC Power 12V/300mA

(1) Standard Model, Without Draft Shield

Specifications & Ordering Information

<L1>

Spec.	Model	"BAL-310D"	"BAL-3100D"
Weighing Range. g		max 310g	max 3100g
Read-Out [d]. g		0.001g <1mg>	0.01g <10mg>
Reproducibility. g		0.002g	0.02g
Linearity. g		0.002g	0.02g
Weighing Mode		g, PCS, %, CT	
Display		Digital LCD with Back-Light Function(110×h35mm)	
Weighing Plate. mm		120×120mm	160×170mm
Dimension & Weight		170×250×h65mm, 1.4kg	
Packing & Gross Weight		210×315×h115mm, 1.6kg	
Power Supply		Adaptor DC 12V/300mA, 220V/60Hz	
Cat. No.		DH.Bal8052	DH.Bal8053
Price, Each. ₩			

(2) With Clear Polycarbonate Draft Shield ; 4면 투명·탈부착형커버형 실드

- * With Removable Powder Coated steel Cover
- * Draft Shield Size ; 156×162×h120mm
- * Removable Cover with Holes(30mm) for Open

Specifications & Ordering Information

Spec.	Model	"BAL-310DS"	"BAL-3100DS"
Weighing Range. g		max 310g	max 3100g
Read-Out [d]. g		0.001g <1mg>	0.01g <10mg>
Reproducibility. g		0.002g	0.02g
Linearity. g		0.002g	0.02g
Weighing Mode		g, PCS, %, CT	
Display		Digital LCD with Back-Light Function(110×h35mm)	
Weighing Plate. mm		120×120mm	160×170mm
Dimension & Weight		170×250×h185mm, 1.7kg	
Packing & Gross Weight		210×315×h115mm, 2.5kg	
Power Supply		Adaptor DC 12V/300mA, 220V/60Hz	
Cat. No.		DH.Bal8054	DH.Bal8055
Price, Each. ₩			

Optional Thermal Printers ; Specify the Balance Model for Printer Installation

* 프린터 구매 시 바란스 모델을 알려주어야 합니다.

		Ea. Price, ₩
SL.WBA.100P	(3) Thermal Printer, Paper Roll, "WBA-100" Included Cable Module, 1× Paper Roll	
DH.WBA0084	Thermal Paper Package, 5× Paper Roll, for "WBA-100"	
SL.WBA.420P	(4) Thermal Printer, Label Paper Roll, "WBA-420" Included Cable Module, 1× Label Paper Roll	
DH.WBA0771	Thermal Label Paper Package, 5× Label Paper Roll, for "WBA-420"	

(2) ϕ 90mm Reusable Al. Dish ϕ 75 · 90 · 102mm Reusable Aluminum Weighing Dish
Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

<L3>

- Ideal for General Weighing and Moisture Drying Application.

DH.Dis8005	(1) ϕ 75 Alumi Weigh-Dish, ϕ 75×h8mm, 4.5g, 30ml	
DH.Dis8006	(2) ϕ 90 Alumi Weigh-Dish, ϕ 90×h4mm, 3.05g, 25ml	New
KE.MLB.A01	, ϕ 90×h7mm, 3.7g, 40ml	
DH.Dis8007	(3) ϕ 102 Alumi Weigh-Dish, ϕ 102×h8mm, 3.8g, 57ml	

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	

Kern® [d] 1 & 10mg, max 350 & 3,500g High-Resolution Lab Balance “PCD”, with Remote Display - Cable Length 1.2m
Ideal for Working in Fume Hoods · Glove Boxes for Toxic · Volatile · Contaminated Substances, with Counting & PRE-TARE Func.

원격 계량 정밀 바란스, 계수계 겸용, 동물계량 가능

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

2-Year Warranty



“PCD300.3” [d] 1mg, max 350g
with Draft Shield



“PCD3000.2” [d] 10mg, max 3,500g
with Stainless-steel Plate

Specification

• Backlight LCD Display, Digit Height 21mm

• Weighing Plate Dimensions

[A] ; ϕ 105mm [B] ; ϕ 160 × 160mm

• Weighing Plate Material

[A] ; Plastic, with Conductive Lacquer

[B] ; Stainless-steel

• Overall Dimensions 165 × 280 × h141mm

• Cable Length of Display Device approx. 1.2m

• Net Weight Approx. 1.2kg

• Ambient-air Temperature 5°C ~ 35°C

■ High-Resolution Precision Balance with Removable Display for Maximum Flexibility

■ Laboratory Balance with Separate Platform ; Ideal When Working in Glove Bag or Fume Cupboard. Particularly Practical for Weighing Toxic, Volatile or Contaminated Substances

* PRE-TARE Function for Manual Subtraction of a Known Container Weight, Useful for Checking Fill-Levels

* Freely Programmable Weighing Unit, e.g. Display Direct in Special Units such as Length of Thread g/m, Paper Weight g/m², or Similar

* Level Indicator and Levelling Feet to Level the Balance Precisely, Fitted as Standard, to Give the Most Accurate Weighing Results

* Draft Shield Standard for “PCD300.3”, Weighing Space 146 × 146 × h80mm

Specifications & Ordering information

<L1>

Cat. No.	Weighing -Range Max g	Read -out d g	Reproducibility g	Linearity g	Min. piece weight [Counting] g/ piece	Weighing Plate	Price(W)
KE.PCD300.3 New	350	0.001	0.002	± 0.005	0.002	[A] ϕ 105mm	
KE.PCD3000.2	3500	0.01	0.02	± 0.05	0.02	[B] 160 × 160mm	

Optional Accessories

Ea.	Price, ₩
KE.PCD.A05	Protective Working Cover, Over the Display Device
KE.PCD.A03	(1) Stand to Elevate Display Device, h250mm
KE.PCD.A04	Rechargeable Battery Pack Internal Operating Time up to 48 h, Charging Time Approx. 8h.
KE.YKI.01	RS-232/Ethernet Adapter
KE.YKS.01	(2) Statistics Printer, “GLP” Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165 × 140 × 50mm, Adapter 230V/50Hz
KE.YKB.01	(3) Thermal Printer, “GLP” Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165 × 140 × 50mm, Adapter 230V/50Hz
KE.572.926	Interface Cable, PC or Thermal Printer, “YKS.01”, “YKB.01”, L1.5m

Features



External Calibration



Data Interface RS 232



GLP/ISO record keeping



Piece Counting



Net-total Weighing



% Determination



Animal Weighing Program ;
Strong Vibrations are Filtered Out.



Weighing Unit



Battery Operation
9V Block



AC 230V Adapter



Strain Gauges

OPTION



Rechargeable battery



(1) Stand, KERN PCD.A03



(2) Statistics Printer
KERN YKS.01



(3) Thermal Printer
KERN YKB.01

Articles

Cat. No

Description

Kern® [d] 1 & 10mg, max 360 & 3,600g Standard Precision Lab Balance "PCB", Compact, with Pcs. Counting
With Pre-Tare & Recipe Function, Freely Programmable Weighing Unit, 표준형 정밀 바란스, 계수계 겸용

■ Supplied by Kern® in Germany

2-Year Warranty



"PCB300.3" [d] 1mg, max 360g



"PCB3000.2" [d] 10mg, max 3,600g

Specification

- Backlight LCD Display, Digit Height 21mm
- Weighing Plate Dimensions
[A]; $\phi 81\text{mm}$ [B]; $130 \times 130\text{mm}$
- Weighing Plate Material
[A]; Plastic, with Conductive Lacquer
[B]; Stainless-steel

- Overall Dimensions ; w163×d245×h65mm
- Permissible Ambient Temperature ; -10℃~40℃
- Batteries Operation ; 4×1.5V AA Block Standard, Operating Time upto 20h

■ Ideal for Standard Weighing Performance in Laboratory, with Compact Size

■ CE & ISO 9001 Certified

■ External Calibration

- * Pre-Tare Function for Manual Subtraction of a Known Container Weight, Useful for Checking Fill-Levels
- * Freely Programmable Weighing Unit, e.g. Display Direct in Special Units such as Length of Thread g/m, Paper Weight g/m²
- * With the Recipe Function you can Weigh the Different Ingredients of a Mixture, also Call up the Total Weight of all the Ingredients
- * Plus & Minus Weighing Procedures, e.g Unit Weight Control, Fabrication Control etc.
- * Includes Ring-Shaped Draught Shield Standard (Dimension ; h90× ϕ 40mm), only for KE.PCB300.3
- * Includes Protective Working Cover. * Auto-Off Function to Preserve the Battery

Features



External Calibration



Data Interface KUP



GLP/ISO record keeping



Piece Counting



Net-total Weighing



% Determination



Weighing Unit



Animal Weighing Program ;
Strong Vibrations are Filtered Out.



Suspended Weighing



Battery Operation
4×1.5V AA



Universal Adapter

Specifications & Ordering information


<L1>

Cat. No.	Weighing Capacity Max g	Readability d g	Reproducibility g	Linearity g	Dimensions & Weight (W×D×H, g)	Weighing Plate	Price(W)
KE.PCB300.3	360	0.001	0.002	± 0.005	163×245×65mm, 0.85kg	[A] $\phi 81\text{mm}$, Plastic with Conductive Lacquer	
KE.PCB3000.2	3600	0.01	0.02	± 0.05	163×245×65mm, 1.6kg	[B] $130 \times 130\text{mm}$, Stainless Steel	

Optional Accessories

- KE.YKS.01 (1) Statistics Printer, "GLP" capable, 0.75 lines/sec, max 40 Characters/line,
with 112mm Paper & RS232, 165×140×50cm, Adapter 230V/50Hz
- KE.YKB.01 (2) Thermal Printer, Capable, 0.75 lines/sec
max 40 Characters/line, with 112mm Paper
- KE.KUP.01 Interface Adapter with Cable, KUP to RS-232 Standard, 1.5m
- KE.KUP.03 , KUP to USB-Device, 1.5m
- KE.YKR.01 Rechargeable Battery Pack, Internal, Operating up to 24h, Charging approx 8h

Ea. Price, ₩



(1)
Statistics Printer
KERN YKS.01

(2)
Thermal Printer
KERN YKB.01

$\phi 75 \cdot 90 \cdot 102\text{mm}$ Reusable Aluminum Weighing Dish Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

■ Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005 (1) $\phi 75$ Alumi Weigh-Dish, $\phi 75 \times h8\text{mm}$, 4.5g, 30ml
- DH.Dis8006 (2) $\phi 90$ Alumi Weigh-Dish, $\phi 90 \times h4\text{mm}$, 3.05g, 25ml **New**
- KE.MLB.A01 , $\phi 90 \times h7\text{mm}$, 3.7g, 40ml
- DH.Dis8007 (3) $\phi 102$ Alumi Weigh-Dish, $\phi 102 \times h8\text{mm}$, 3.8g, 57ml

<L3>

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	



(2) $\phi 90\text{mm}$ Reusable Al. Dish

Kern® [d] 1 & 10mg, max 300 & 3,000g Hi-Function Precision Lab-Balance "EWJ", with Draft Shield On 300Model With Auto Internal Adjustment, Plate Size ϕ 80 & 120mm, 220×315×h105mm, "Int-CAL 자동 보정형" 정밀 바란스, 계수계 겸용

■ Supplied by Kern® in Germany ■ CE & ISO 9001 Certified

"Auto Int-CAL" "EWJ"

NEW

3-Year Warranty



Draft Shield Included On Model "EWJ300.3" But, Not on 3000 Model. Draft Shield는 300모델에 한함

"EWJ300.3" [d] 1mg, max 300g with Draft Shield

Specification

- Backlight LCD Display, Digit Height 16.5 mm
- Weighing Plate Dimensions, Stainless-steel, ϕ 80 & 120mm
- Ambient-air Temperature 15°C~35°C

- High-quality Precision Balance with Automatic Internal Adjustment
- Draft Shield ; Standard on "Model EWJ300.3", But not on "Model EWJ3000.2"

- * Protective Working Cover ; Included with Delivery
- * Ideal for Mobile Applications which Require Verification, such as Ambulatory Gold and Jewelry Purchasing
- * Automatic Internal Adjustment Time-Controlled Every 4 Hours. Guarantees High Degree of Accuracy and Makes the Balance Independent of its Location
- * Draft Shield with Removable Cover and Opening for Using a Pipette
- * RS-232 and USB Interfaces for Transferring Weighing Data to The PC, Printer, USB Sticks etc.

Specifications & Ordering information

<L1>

Cat. No.	Weighing Range Max g	Read-out d g	Linearity g	Overall Dimension & Weight (W×D×H, kg)	Weighing Plate	Price(W)
KE.EWJ300.3	300	0.001	± 0.005	220×315×h105mm, 2kg	ϕ 80 mm	
KE.EWJ3000.2	3000	0.01	± 0.05		ϕ 120 mm	

Optional Accessories

- KE.KFB.A01 Rechargeable Battery Pack Internal
Operating time up to 15 h without backlight, charging time approx. 4h.
- KE.DBS.A02 USB Interface Kit, USB Cable, Driver, Software Balance Connection
- KE.YKS.01 (1) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50 mm, Adapter 230V/50Hz
- KE.YKB.01 (2) Thermal Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165×140×50 mm, Adapter 230V/50Hz
- KE.CFS.A01 RS-232 Interface Cable, PC or Thermal Printer, "YKS.01" / "YKB.01", 1.5m for PFB/EWJ-Model

Ea. Price, ₩



(1) Statistics Printer KERN YKS.01
(2) Thermal Printer KERN YKB.01

ϕ 75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평판 디쉬

- Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005 (1) ϕ 75 Alumi Weigh-Dish, ϕ 75×h8mm, 4.5g, 30ml
- DH.Dis8006 (2) ϕ 90 Alumi Weigh-Dish, ϕ 90×h4mm, 3.05g, 25ml
- KE.MLB.A01 , ϕ 90×h7mm, 3.7g, 40ml
- DH.Dis8007 (3) ϕ 102 Alumi Weigh-Dish, ϕ 102×h8mm, 3.8g, 57ml

<L3>

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	



(2) ϕ 90mm Reusable Al. Dish

Kern® [d] 1 & 10mg, max 620 & 3,200g Premium Precision Lab Balance Ext-CAL "PNS" & Auto Int-CAL "PNJ"

With Robust Tuning Fork Weighing System, External & Auto Internal Calibration, 정밀 바란스, 외부 보정형 & 자동 보정형, 계수계 겸용, 작업보호커버 포함

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

3-Year Warranty CE

"Ext-CAL" "PNS"
"Auto Int-CAL" "PNJ"[d] 1mg, 0.001~620g
with Large Glass Draught Shield

[d] 10mg, 0.01~3,200g

Specification

- Weighing Plate Dimensions, ϕ 140mm
- Overall Dimensions, 202 × 293 × h266mm (with Draught Shield)
- Net Weight, Approx. 2.2 kg
- Ambient Temperature, 5 °C~40 °C

Specification

- Weighing Plate Dimensions, w190 × d190mm
- Overall Dimensions, 195 × 290 × h90mm
- Net Weight, Approx. 2.8kg
- Ambient Temperature 5 °C~40 °C

Features



External Calibration



Internal Calibration



Data Interface RS 232



GLP/ISO record keeping



Piece Counting



% Determination



Weighing Unit



Weighing with Tolerance Range



Universal Adapter

Weighing Principle
; Tuning Fork

A resonating body is electromagnetically excited, causing it to oscillate.

* High-quality Tuning Fork Measuring System for Steady Weight Values and Continuous Weighing.

* PNS ; Adjusting Program CAL for Quick Setting of the Balance Accuracy Using an External Test Weight.

* PNJ ; Automatic Internal Adjustment, Guarantees High Degree of Accuracy and Makes the Balance Independent of Its Location of Use.

* Capacity Display ; A Bar Lights up to show How Much of the Weighing Capacity is still available.

* Precise Counting ; The Automatic Reference Weight Optimisation of Reference Weight Gradually Improves the Average Piece Weight Value

* Protective working cover included with delivery

Specifications & Ordering information

<L1>

Cat. No.	Weighing-Range Max g	Read-out d g	Reproducibility g	Linearity g	Weighing Plate	Price(₩)
1) Ext-CAL						
KE.PNS600.3	620	0.001	0.001	± 0.004	ϕ 140 mm	
KE.PNS3000.2	3200	0.01	0.01	± 0.02	190 × 190 mm	
2) Auto Int-CAL NEW						
KE.PNJ600.3	620	0.001	0.001	± 0.004	ϕ 140 mm	
KE.PNJ3000.2	3200	0.01	0.01	± 0.02	190 × 190 mm	

Optional Accessories

Ea.	Price, ₩
KE.YKI.02	RS-232/Bluetooth Adaptor, to connect to Bluetooth Capable Devices, Such as Bluetooth Printers, Tablets, Laptops, Smartphones
KE.YKS.01	(1) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165 × 140 × 50 mm, Adapter 230V/50Hz
KE.YKB.01	(2) Thermal Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165 × 140 × 50 mm, Adapter 230V/50Hz
KE.MLB.A05	RS-232 Interface Cable, PC or Thermal Printer, for "YKS.01", "YKB.01", 1.5m

(1) Statistics Printer
KERN YKS.01(2) Thermal Printer
KERN YKB.01 ϕ 75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

■ Ideal for General Weighing and Moisture Drying Application.

DH.Dis8005	(1) ϕ 75 Alumi Weigh-Dish, ϕ 75 × h8mm, 4.5g, 30ml
DH.Dis8006	(2) ϕ 90 Alumi Weigh-Dish, ϕ 90 × h4mm, 3.05g, 25ml New
KE.MLB.A01	, ϕ 90 × h7mm, 3.7g, 40ml
DH.Dis8007	(3) ϕ 102 Alumi Weigh-Dish, ϕ 102 × h8mm, 3.8g, 57ml

<L3>

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	

(2) ϕ 90mm Reusable Al. Dish

Kern® [d] 0.01g, max 2,000g Convenient Lab Balance "EMB", Underfloor Weighing Model, ϕ 150mm Weighing Plate
With Tare & Auto-Off Function, Hook for Underfloor Weighing, 다용도랩 바란스, 언더웨이잉 기능형



2-Year Warranty



Simple 2-Key Operation

"EMB" Under Weighing Model



(1) Statistics Printer (2) Thermal Printer

- Supplied by Kern® in Germany
- CE & ISO 9001 Certified
- Ideal for Economical Precise Weighing Performance in Laboratory
- External Calibration

- * Hook for Underfloor Weighing ; Objects which, because of their Size or Shape, cannot be Put on the Scale, can be Weighed by Means of Underfloor Weighing
- * Simple and Convenient 2-Key Operation
- * Tare Function ; Used for Setting Weight to Zero, Facilitates Formulation Work
- * Auto-Off Function to preserve the battery
- * Batteries Operation ; 9V Block Standard
- * Large LCD Display, Digit Height 15mm
- * Permissible Ambient Temperature ; 5°C~35°C

Specifications & Ordering information

<L1>

Cat. No.	Weighing Capacity Max g	Readability d g	Reprodu- -cibility g	Linearity g	Dimensions & Weight (W×D×H, g)	Weighing Plate (ϕ)	Price, ₩
KE.EMB2000.2	2,000g	0.01g	0.01g	±0.05g	170×240×52mm, 0.75kg	ϕ 150, Plastic	

Optional Printers

Ea. Price, ₩

- KE.YKS.01 (1) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50mm, Adapter 230V/50Hz
- KE.YKB.01 (2) Thermal Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165×140×50mm, Adapter 230V/50Hz
- KE.572.926 Interface Cable, PC or Thermal Printer, "YKS.01, YKB.01", L1.5m

Kern® [d] 0.01g, max 2,000g Reference Lab Balance "EMB.V", Density Determination Model(Density Part ; Option)
With Self-Explanatory Graphic 4 steps of Control Panel, 비중 측정형 랩 바란스, 비중측정 부품은 별매



2-Year Warranty

"EMB.V" Density Model



(1) Statistics Printer (2) Thermal Printer

- Supplied by Kern® in Germany
- CE & ISO 9001 Certified
- Ideal for Density Determination & Precision Weighing Performance in Education
- For Density Determination, Please Order the Following Optional Set (KE.YDB02) at the Same Time!
- Also with Carat Weighing Unit ; [d]0.05ct, max 10,000ct
- External Calibration

- * The 4 steps are Carried Out from Left to Right ;
- ① Tare Function ② Select Density Determination Mode
- ③ Weighing of Samples / Plummets in Air
- ④ Weighing of Samples / Plummets in Liquid
- * Thanks to Self-Explanatory Graphic 4 steps Control Panel, the Density of Solids and Liquids can be Determined in Seconds
- * Large LCD Display, Digit Height 15mm
- * Permissible Ambient Temperature ; 5°C~35°C
- * Batteries Operation ; 9V Block Standard, Operating Time up to 12h
- * Auto-Off Function to Preserve the Battery
- * Hook for Underfloor Weighing ; Objects which, because of their Size or Shape, cannot be Put on the Scale, can be Weighed by Means of Underfloor Weighing.



Specifications & Ordering information

<L1>

Cat. No.	Weighing Capacity Max g	Readability d g	Reprodu- -cibility g	Linearity g	Dimensions & Weight (W×D×H, g)	Weighing Plate (ϕ)	Price, ₩
KE.EMB2000.2V	2,000g	0.01g	0.02g	±0.05g	175×245×54mm, 0.85kg	ϕ 150, Plastic	

Optional Printers

Ea. Price, ₩

- KE.YKS.01 (1) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50mm, Adapter 230V/50Hz
- KE.YKB.01 (2) Thermal Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165×140×50mm, Adapter 230V/50Hz
- KE.572.926 Interface Cable, PC or Thermal Printer, "YKS.01, YKB.01", L1.5m

Optional Set for Density Determination (Balance not Included)

- Note ; Balance and Appropriate Set for Density Determination should be Ordered at the same time!
- * Complete Set with ; Weighing Plate, Beaker(H135× ϕ 100mm), Sample Holder, Plummets

Ea. Price, ₩

KE.YDB02 Optional ; Set for Density Determination

※ Figure ; Optional Set for
Density Determination Installed



Articles

Cat. No

Description

DAIHAN™ [d] 0.01g, max 3,200 & 6,200g Calibration Certified Hi-Standard Lab Balance “WBA-3200 & 6200”, 168×168mm Plate
With Super Size Back Light LCD, Counting Function, 2kg Cali. Weight, Various Weight Mode
표준 다용도 랩 바란스, g · 계수 · % 측정, 보정용 표준분동 포함, 4각 플레이트, 보정 성적서 포함



[d] 0.01g, max 3,200g



Calibration Certificate 포함

3-Year Warranty

Includes 2kg(2 × 1kg) Calibration Weight
(Only “WBA-3200” model)

Hi-Standard Lab Balance, “WBA-3200”



(1) Glass Draft Shield (Optional)



(2) “WBA-100” RS-232C Thermal Printer



(3) “WBA-420” RS-232C Thermal Label Printer

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- With Calibration Certificate(보정 성적서 포함)
- With Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-Information, and Traceable Data Base System
- Various Weight Mode (g/ct/GN/dwt/etc.) and Percent(%) Weighing Mode
- * Numeric Linearity Calibration
- * Includes AC 230V Adapter for DC Power 12V/500mA
- * Four Level Anti-shock Adjustable Function
- * Time & Date Setting for Recording Data
- * Weighing Response Time Adjustable
- * Default and Auto Lock Present Working Mode
- * Display Present Temperature
- * Built in RS232C for Printer (SL.WBA.100P, SL.WBA.420P)
- * Electro Magnetic Base
- * Glass Draft Shield(Optional)

Specifications & Ordering information

<L1>

Spec.	Model	“WBA-3200”	“WBA-6200”
Read-Out [d]. g		0.01g	
Weighing Range. g		max 3,200g	max 6,200g
Weighing Plate. mm		168 × 168 mm	
Reproducibility. g		0.01g	0.02g
Linearity. g		0.02g	0.03g
Ambient Temperature		5.0℃~35℃	
Stabilization Time		≤ 3 Seconds (Typically Fast)	
Weight Mode		% Mode and g-/ct-/GN-/dwt-/etc.	
Counting Mode		Counting & Weighing Mode Available	
Calibration		External Calibration	
		Includes 2 × 1kg Calibration Weight	None
Display		Digital Super Size LCD with Back-Light Function	
Dimension & Net Weight		195 × 295 × h90mm / 2.6kg	
Packing Size & Gross Weight		243 × 377 × h163mm / 3.2kg	
Electrical Reqts.		DC 12V, with AC 230V Adapter, 50/60Hz	
Cat. No.		DH.WBA3200	DH.WBA6200
Price, Each. ₩			

Optional Glass Draft Shield

- * Target Model Required for Shielding Installation
- * 설치를 위해 사용 중인 발란스를 제조사로 보내야 합니다.

DH.WBA0003 (1) Glass Draft Shield, 3 Slide Doors(Left, Right, Top),
175 × 195 × h230mm

Ea. Price, ₩

Optional Thermal Printers ; Specify the Balance Model for Printer Installation

- * 프린터 구매 시 바란스 모델을 알려주어야 합니다.

SL.WBA.100P (2) Thermal Printer, Paper Roll, “WBA-100”
Included Cable Module, 1 × Paper Roll
DH.WBA0084 Thermal Paper Package, 5 × Paper Roll, for “WBA-100”
SL.WBA.420P (3) Thermal Printer, Label Paper Roll, “WBA-420”
Included Cable Module, 1 × Label Paper Roll
DH.WBA0071 Thermal Label Paper Package, 5 × Label Paper Roll, for “WBA-420”

Ea. Price, ₩



(2) φ 90mm Reusable Al. Dish

φ 75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

<L3>

- Ideal for General Weighing and Moisture Drying Application.

DH.Dis8005 (1) φ 75 Alumi Weigh-Dish, φ 75 × h8mm, 4.5g, 30ml
DH.Dis8006 (2) φ 90 Alumi Weigh-Dish, φ 90 × h4mm, 3.05g, 25ml New
KE.MLB.A01 , φ 90 × h7mm, 3.7g, 40ml
DH.Dis8007 (3) φ 102 Alumi Weigh-Dish, φ 102 × h8mm, 3.8g, 57ml

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	

DAIHAN™ [d] 0.01 & 0.1g, max 610~6,100g Multi-use Balance “BA-Series”, DC & AC Adapter

With g · Metric · Pcs · % · Carat · CHK · Bargraph · Density Test-Function, 다용도 바런스, 무게 · 계수 · 퍼센트 · 캐럿 · 비교계량 · 비중측정 & 동물계량 기능

■ Supplied by DAIHAN Korea ■ CE Certified ■ 3-Year Warranty ■ Digital Large LCD Display and Ergonomic Design

* Metric Function ; Animal Weighing, Bar-Graph Display, Comparison Weighing CHK, Counting 10~100 Pcs., Density Determination, Percent % and Sub-Scale Dual Display

* 5 × Touch Key ; 1/2D*, MODE, ZERO, TARE, ENTER, and PRINT

* 3sec. Stabilization Time, and Permissible Ambient Temp 5℃~35℃

* Rechargeable Battery integrated, and with AC 230V Adapter for DC Power 12V/300mA

* Easy Operation Weighing Mode ; Gram, %, and Carat

* Material ; Stainless-steel Weighing Plate and ABS Resin Body

* RS-232C Interface and Level-Gauge integrated

(1) [d] 0.01g, max 610g Models

Weighing Plate ; 120 × 130mm

CE



3-Year Warranty

Specifications & Ordering information

<L1>

Spec.	Model	(1) “BA610”	(2) “BA3.1K”	(3) “BA6.1K”
Weighing Range max		610g	3,100g	6,100g
Read-Out [d] & Reproducibility. g		0.01g	0.1g	
Linearity. g		0.01g	0.1g	
Weighing Mode		Gram, %, Carat		
Weighing Plate. mm		120×130mm	160×170mm Stainless-steel	
Color		Silver-color	Light Grey-color	Silver-color
Dimension		170×250×h65mm		
Weight		0.8kg	1.0kg	
Packing		210×315×h115mm		
Gross Weight		1.7kg	1.9kg	
DC Power & AC Adapter		with Rechargeable Battery & AC 220 Adapter		
Cat. No.		DH.Bal8013	DH.Bal8016	DH.Bal8017
Price, Each. ₩				

(2) & (3) [d] 0.1g, max 3100g & 6100g Models

Weighing Plate ; 160 × 170mm

CE



3-Year Warranty

Optional Thermal Printers ; Specify the Balance Model for Printer Installation

* 프린터 구매 시 바런스 모델을 알려주어야 합니다.

Ea. Price, ₩

SL.WBA.100P	(4) Thermal Printer, Paper Roll, “WBA-100” Included Cable Module, 1 × Paper Roll	
DH.WBA0084	Thermal Paper Package, 5 × Paper Roll, for “WBA-100”	
SL.WBA.420P	(5) Thermal Printer, Label Paper Roll, “WBA-420” Included Cable Module, 1 × Label Paper Roll	
DH.WBA0771	Thermal Label Paper Package, 5 × Label Paper Roll, for “WBA-420”	



(4) “WBA-100” RS-232C Thermal Printer



(5) “WBA-420” RS-232C Thermal Label Printer



(2) φ90mm Reusable Al. Dish

φ75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

<L3>

■ Ideal for General Weighing and Moisture Drying Application.

		Pk.	Price, ₩
DH.Dis8005	(1) φ75 Alumi Weigh-Dish, φ75 × h8mm, 4.5g, 30ml	100ea	
DH.Dis8006	(2) φ90 Alumi Weigh-Dish, φ90 × h4mm, 3.05g, 25ml	50ea	
KE.MLB.A01	, φ90 × h7mm, 3.7g, 40ml	80ea	
DH.Dis8007	(3) φ102 Alumi Weigh-Dish, φ102 × h8mm, 3.8g, 57ml	50ea	

Kern® [d] 0.05 & 0.1g, max 12 · 20 · 24kg Precision Heavy-Duty Balance "572", 180×310×h90mm

With Piece Count & Large LCD Display, 대용량 정밀 바란스, 계수계 겸용

■ Supplied by Kern® in Germany ■ CE & ISO 9001 Certified



■ All-round ; as Laboratory, Counting, or Checking

■ Thanks to the Many Typical Laboratory Functions, Ideal for Day-to-day Work in the Laboratory

* With Recipe Function & Percentage Determination & GLP Record Keeping

* With Piece-Counting, Vibrationfree Weighing and the Large Weighing Capacities also Make these Balances Ideal for all Industrial Applications

* Freely Programmable Weighing Unit, e.g. Display Direct in Special Units such as Length of Thread g/m, Paper Weight g/m², or Similar

* The Robust Aluminium Diecast Housing Maintains the Stability, Protects the Weighing Technology Elements

* Protective Working Cover Included


Specifications & Ordering information

Cat. No.	Weighing-Range Max g	Read-out d g	Repro-ducibility g	Linea- rity g	Overall Dimension & Weight (W×D×H, kg)	Weighing Plate	Price(W)
KE.572.45	12000	0.05	0.05	±0.15	130×310×h90mm, 2.7kg	160×200mm	
KE.572.55	20000	0.05	0.1	±0.15	130×310×h90mm, 2.7kg	160×200mm	
KE.572.57	24000	0.1	0.1	±0.3	130×310×h90mm, 2.7kg	160×200mm	

Optional Accessories

KE.YKS.01	(1) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50 mm, Adapter 230V/50Hz
KE.YKB.01	(2) Thermal Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper, 165×140×50 mm, Adapter 230V/50Hz
KE.KUP.01	Interface Adapter with Cable, KUP to RS-232 Standard, 1.5m
KE.KUP.03	, KUP to USB-Device, 1.5m
KE.YKR.01	Rechargeable Battery Pack, Internal, Operating up to 24h, Charging approx 8h

Ea.	Price, ₩



(1) Statistics Printer
KERN YKS.01

(2) Thermal Printer
KERN YKB.01

φ 75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

■ Ideal for General Weighing and Moisture Drying Application.

DH.Dis8005	(1) φ 75 Alumi Weigh-Dish, φ 75×h8mm, 4.5g, 30ml
DH.Dis8006	(2) φ 90 Alumi Weigh-Dish, φ 90×h4mm, 3.05g, 25ml New
KE.MLB.A01	, φ 90×h7mm, 3.7g, 40ml
DH.Dis8007	(3) φ 102 Alumi Weigh-Dish, φ 102×h8mm, 3.8g, 57ml

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	



(2) φ 90mm Reusable Al. Dish

Features

- External Calibration
- Data Interface KUP
- GLP/ISO record keeping
- Piece Counting
- Net-total Weighing
- % Determination
- Weighing Unit
- Animal Weighing Program ; Strong Vibrations are Filtered Out.
- Suspended Weighing
- AC 230V Adapter
- Strain Gauges
- Rechargeable Battery Pack

Specifications

- * **Large Backlit LCD Display**, Digit Height 18 mm
- * **Weighing Plate Dimensions** ; 160×200mm, Stainless-steel
- * **Permissible Ambient Temperature** ; 10°C~40°C

ALL for the Best BALANCES
--- via "KERN" since 1844


Kern® [d] 0.1 & 0.2g, max 36 & 65kg **Premium** High Precision Lab-Scale "CKE", Stainless-steel Weigh-Plate 340×240mm
With Laboratory- & Self Explanatory Counting-Function, 프리미엄 정밀 대용량 저울-스케일, 계수기능

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

"CKE" Scale

3-Year Warranty CE



The 4 steps are Carried Out from Left to Right

350×390×h120mm

Features



External Calibration



Data Interface KUP



Piece Counting



Totalising



% Determination



Weighing Unit



Weighing with Tolerance Range



Rechargeable Battery



AC 230V Adapter



Strain Gauges

Specification

- Backlit LCD Display ; Digit Height 25mm
- Batteries Operation ; 6×1.5V Size C, Operating Time up to 40h
- Permissible Ambient Temperature ; 10℃~40℃

- Ideal for High Capacity Precision Weighing & Piece Counting Performance in Laboratory
- Thanks to Self-explanatory Graphic Control Panel, Easy to Use, Ideal for Untrained Users

* The 4 steps are Carried Out from Left to Right ;

- ① Place the Empty Container onto the Weighing Plate and Tare by Pressing the TARE Key
- ② Place the Reference Quantity for the Goods to be Counted into the Container (5, 10 or 20 pieces)
- ③ Confirm the Selected Reference Quantity by Pressing the Key (5, 10 or 20)
- ④ Pour in the Goods to be Counted. The Number of Pieces will Immediately be Shown in the Display

* **Two Balance in One** ; Switching from Counting Mode to Weighing Mode at the Touch of a Key

* **Precise Counting** ; The Automatic Reference Weight Optimisation of Reference Weight Gradually Improves the Average Piece Weight Value

Specifications & Ordering information

<L1>

Cat. No.	Weighing Capacity Max kg	Read-ability d g	Smallest Part Weight Normal g / piece	Counting Resolution Point	Dimensions & Weight (W×D×H, kg)	Weighing Plate (W×D)	Price(₩)
KE.CKE.36K0.1	36kg	0.1g	1g	360,000	350×390×120mm, 7kg	340×240mm, Stainless Steel	
KE.CKE.65K0.2	65kg	0.2g	2g	325,000			

Optional Accessories

KE.FKB.A02S05 Protective Working Cover, Weighing Plate Size 340×240mm, 5ea

KE.YKS.01 (1) Statistics Printer, "GLP" Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232, 165×140×50mm, Adapter 230V/50Hz

KE.YKB.01 (2) Standard Printer, Capable, 0.75 lines/sec, max 40 Characters/line, with 112mm Paper & RS232

Cs. Price, ₩

5ea

Ea. Price, ₩



(1)
Statistics Printer
KERN YKS.01



(2)
Thermal Printer
KERN YKB.01

Articles

Cat. No

Description

DAIHAN™ [d] 0.1 & 1g, max 610 & 6,100g Economical Compact Scale "BAS", 17.5 × 24 × h6cm

Good for Saving Money, with Pcs. Counting, Back Light LCD, Battery, without Adapter, 경제형 컴팩트 스케일/저울, 계수 기능

"BAS" Eco-Scale



CE

3-Year Warranty

■ Supplied by DAIHAN Korea

* Ergonomic Design for Intuitive User Interface

* Ambient Temperature 5°C~35°C

■ CE Certified

* Pcs. Counting & Weighing Mode

* Numeric Linearity Calibration

Specifications & Ordering information

<L1>

Spec.	Model	"BAS610"	"BAS6.1K"
Weighing Range. g		max 610g	max 6,100g
Weighing Plate. mm		155 × 155mm	
Read-Out [d] & Reproducibility. g		0.1g & 0.1g	1g & 1g
Material		Body ; ABS, Platform ; Stainless-steel	
Dimension & Weight		175 × 240 × h60mm, 0.6kg	
Packing & Gross Weight		280 × 180 × h80mm, 0.8kg	
Power Supply		DC 9V/300mA / AA × 2Ea.	
Cat. No.		DH.Bal8021	DH.Bal8022
Price, Each. ₩			

DAIHAN™ [d] 0.1 & 1g, max 610 & 6,100g Functional Compact Scale "BAP610" & "BA6.1KS", 17 × 25 × h6.5cm

With 10~100Pcs. Counting · Bargraph · LCD Display, with DC & AC Adapter, 기본형 스케일/저울 계수기능

"BAP" Functional Scale



CE

3-Year Warranty

■ Supplied by DAIHAN Korea

■ CE Certified

■ Digital Large LCD Display and Ergonomic Design

* Easy Operation, Sturdy and Balancing with Level Gauge

* 10~100 Pcs. Counting and Bar-graph Display Function

* Material ; Stainless-steel(STS) Weighing Plate and ABS Resin Body

* 4 × Touch Key ; MODE, ZERO, TARE, and ENTER

* AA Battery and AC 230V Adapter for DC Power 12V/300mA

* Fast Stabilization Time - 3 seconds

* Permissible Ambient Temp 5°C~35°C

Specifications & Ordering information

<L1>

Spec.	Model	"BAP610"	"BA6.1KS"
Weighing Range. g		max 610g	max 6,100g
Read-Out [d] & Reproducibility. g		0.1g & 0.1g	1g & 1g
Linearity. g		0.1g	1g
Weighing Mode		Gram	
STS Weighing Plate. mm		120 × 130mm	160 × 170mm
Color Light		Silver-color	
Dimension & Weight		170 × 250 × h65mm, 0.8kg	170 × 250 × h65mm, 1.0kg
Packing & Gross Weight		210 × 315 × h115mm, 1.7kg	210 × 315 × h115mm, 1.9kg
DC Power & AC Adapter		with Rechargeable Battery & AC 220V Adapter	
Cat. No.		DH.Bal8026	DH.Bal8031
Price, Each. ₩			



(2) φ90mm Reusable Al. Dish

φ75 · 90 · 102mm Reusable Aluminum Weighing Dish

Conveniently Selectable Dishes, 재사용 알루미늄 평량 디쉬

<L3>

■ Ideal for General Weighing and Moisture Drying Application.

- DH.Dis8005 (1) φ75 Alumi Weigh-Dish, φ75 × h8mm, 4.5g, 30ml
 DH.Dis8006 (2) φ90 Alumi Weigh-Dish, φ90 × h4mm, 3.05g, 25ml **New**
 KE.MLB.A01 , φ90 × h7mm, 3.7g, 40ml
 DH.Dis8007 (3) φ102 Alumi Weigh-Dish, φ102 × h8mm, 3.8g, 57ml

Pk.	Price, ₩
100ea	
50ea	
80ea	
50ea	

Articles	Cat. No	Description
----------	---------	-------------

[d] 2 & 10g, max 10 & 30kg Economical Waterproof Scale “T-2810K” & “T-2830K”, Digital LCD, IP×2
 With Rechargeable Battery, Tray, Back Light & LCD Display, DC Adapter, Pcs. Counting, 경제형 방수 스케일/저울

CE



“T-28K” Waterproof Scale
; 10kg & 30kg

230×240×h170mm

3-Year Warranty

- Supplied by DAIHAN Korea
 ■ Digital LCD Display
 * Included Rechargeable Battery (5,000mAh)
 * Waterproofing (Protection Against Falling Water)
 * Body & Platform ; Made of ABS

■ CE Certified
 ■ Pcs. Counting
 * AC 220V Adapter for DC Power 12V / 1.5A
 * Smart Auto-off 10min

■ 3-Year Warranty

Specifications & Ordering information <L1>

Spec.	Model	“T-2810K”	“T-2830K”
Weighing Range. g		max 10kg	max 30kg
Read-Out [d]. g		2g	10g
Reproducibility. g		2g	10g
Display		g, kg, PCS ; LCD, 5 Digit (100×h40mm), Back Light On/Off	
Platform Size. mm		230×190 mm, ABS	
Dimension & Weight		230×240×h170mm, 2.7kg	
Packing & Gross Weight		350×260×h265mm, 3.6kg	
Electrical Requirements		Adaptor DC 12V/1.5A, 220V/60Hz	
Cat. No.		DH.Bal3026	DH.Bal3027
Price, Each. ₩			

DAIHAN™ [d] 0.1 & 0.5g, max 32 & 60kg Precision High Capacity Digital Bench Scale “Bah32P” & “Bah60P”
 With Pcs. · % · Alarm · Print & Check Function and Rechargeable Battery/Back Light & LCD Display, DC & AC Adapter
 정밀 대용량 다기능 스케일/저울, Pcs. 계수 · % · 경보음 · 인쇄 및 체크 기능, 재충전 배터리/백라이트 내장

CE



320×340×h125mm

3-Year Warranty

(1) [d] 0.1g, max 32kg, Bench Scale, “Bah32P”

- CE Certified
 ■ Digital ; 3×6 digits LCD display, 20mm Digits, and Sturdy/Ergonomic Design
 ■ 7× Touch Key for “Bah32P” ; UNIT, PRINT, FUNC., %, LIMIT, TARE, and ZERO
 ■ 6× Touch Key for “Bah60P Series” ; UNIT, M+, G/N, TARE, ZERO and ON/OFF
 ■ Easy Operation & Multi-Function ; Alarm, Check, Counting Pcs., Percent %
 * Integrated Rechargeable Battery(100Hr-Run) and Back Light
 * Very Fast Stabilization Time ; 2 seconds
 * Material ; Stainless-steel Platform and ABS Resin Body
 * External Calibration

■ 3-Year Warranty
 * with AC 220V Adapter for DC Power 9V / 850mA

* Permissible Ambient Temp. −10℃ +40℃

Specifications & Ordering information <L1>

Spec.	Model	(1) “Bah32P”	(2) “Bah60P”
Weighing Range		max 32kg	max 60kg
Read-Out [d]. g		0.1g	0.5g
Reproducibility. g		0.1g	0.5g
Weighing Plate. mm		225×300mm	355×320mm
Dimension & Weight		320×340×h125mm, 3.8kg	355×320×h130mm, 6.3kg
Packing & Gross Weight		370×440×h210mm, 4.6kg	500×450×h230mm, 7.3kg
DC Power & AC Adapter		with Rechargeable Battery & AC 220V Adapter	
Cat. No.		DH.Bal8035	DH.Bal8036
Price, Each. ₩			

CE



355×320×h130mm

3-Year Warranty

(2) [d] 0.5g, max 60kg, Bench Scale, “Bah60P”

Articles

Cat. No

Description

DAIHAN™ [d] 0.5 & 1g, max 15 & 30kg Basic High Capacity Bench Scale "BAH15L" & "BAH30L"

With Pcs. · % · Alarm · Print & Check Function and Rechargeable Battery/Back Light & LCD Display, DC & AC Adapter, 32 × 34 × h12.5

기본형 대용량 다기능 스케일/저울, Pcs. 계수 · % · 경보음 · 인쇄 및 체크 기능, 재충전 배터리/백라이트 내장



CE

"BAH15L" & "BAH30L" Bench Scale
; 15kg & 30kg32 × 34 × h12.5cm
3-Year Warranty

- Supplied by DAIHAN Korea
- CE Certified
- 3-Year Warranty
- Digital ; 3×6 digits LCD Display, 20mm Digits, and Sturdy/Ergonomic Design
- 7× Touch Key ; UNIT, PRINT, FUNC., %, LIMIT, TARE, and ZERO
- Easy Operation & Multi-Function ; Alarm, Check, Counting Pcs., Percent %, and Print. "RS-232C"-Option

- * Integrated Rechargeable Battery(100Hr-Run) and Back Light
- * Very Fast Stabilization Time ; 2 seconds
- * Permissible Ambient Temp. -10℃ +40℃
- * Material ; Stainless-steel Platform and ABS Resin Body
- * External Calibration
- * with AC 220V Adapter for DC Power 9V / 850mA

Specifications & Ordering information

<L1>

Spec.	Model	"BAH15L"	"BAH30L"
Weighing Range. g		max 15kg	max 30kg
Read-Out [d]. g		0.5 g	1 g
Reproducibility. g		0.5 g	1 g
Weighing Plate. mm		225 × 300mm	
Dimension & Weight		320 × 340 × h125mm, 3.8kg	
Packing & Gross Weight		370 × 440 × h210mm, 4.6kg	
DC Power & AC Adapter		with Rechargeable Battery & AC 220V Adapter	
Cat. No.		DH.Bal8041	DH.Bal8042
Price, Each. ₩			

DAIHAN™ [d] 5 & 10g, max 60 & 150kg High Capacity Multi-func. Platform Scale "BAB60" & "BAB150", 42.5 × 67.5 × h75.5cm

With Pcs. · Alarm · Back Light · Rechargeable Battery & AC Adapter, 대용량 다기능 스케일/저울

"BAB" Platform Scale
; 60kg & 150kg

CE

42.5 × 67.5 × h75.5cm
3-Year Warranty

- Supplied by DAIHAN Korea
- CE Certified
- 3-Year Warranty
- Digital ; 6 Digits, LCD Display, 24mm Digits
- Sturdy Ergonomic Design
- Easy Operation and Multi Function ; Pcs. Count, Weighing Alarm, Weight Storage/Plus/Call, and Auto-TARE

- * Integrated Rechargeable Battery & Back Light
- * Very Fast Stabilization Time ; 1 sec.
- * Permissible Ambient Temp. -10℃ +40℃
- * Automatic External Calibration
- * Material ; Stainless-steel Platform, ABS Resin Indicator
- * with AC 220V Adapter for DC Power 9V/800mA

Specifications & Ordering information

<L1>

Spec.	Model	"BAB60"	"BAB150"
Weighing Range		max 60kg	max 150kg
Read-Out [d]. g		5g	10g
Reproducibility. g		5g	10g
Weighing Plate. mm		420 × 520mm	
Dimension & Weight		425 × 675 × h755mm, 14kg	
Packing & Gross Weight		580 × 710 × h820mm, 15.9kg	
DC Power & AC Adapter		with Rechargeable Battery & AC 220V Adapter	
Cat. No.		DH.Bal8045	DH.Bal8046
Price, Each. ₩			

Kern® [d] 1 & 2g, max 8 & 15kg Compact Pet Dual Scale “FOB.NL”, Stainless-steel Plate 252 × 200mm
Protection against Dust and Water Splashes, IP67, with Protective Working Cover, 소형 애완동물용 저울, 방수 플레이트

■ Supplied by Kern® in Germany

■ Stainless Steel Design of the Housing and Weighing plate, Its Smooth Surfaces Make it Simple to Clean

* Secure & Non-slip Positioning with Rubber Feet

* HACCP-Compliant Quality System.

■ CE & ISO 9001 Certified

* Protective Working Cover Included with Delivery

“FOB” Pet Scale

2-Year Warranty

CE



285 × 255 × h90mm

“FOB10K.3NL” [d] 1g/2g, max 8/15kg

Specification

- Large Backlit LCD Display, Digit Height 25mm
- Batteries Included 4 × 1.5 V AA, Running 48 hours
- Permissible Ambient Temperature 5°C~35°C
- Dual-range Balance Switches Automatically to the Next Largest Weighing Capacity and Readability
- IP67 ; Protected against Dust and Water Splashes

Features



Adjusting Program CAL



Weighing Unit



Protection Against Dust and Splashes IPxx



Stainless Steel



Battery Operation



Strain Gauges

Specifications & Ordering information

<M1>

Cat. No.	Weighing Range Max kg	Readability & Reproducibility d g	Linearity g	Overall Dimensions (w × d × h) & Weight	Weighing Plate (w × d × h)	Price(W)
KE.FOB10K.3NL	8 kg / 15 kg	1 g / 2 g	20	285 × 255 × 90mm, 3.8kg	252 × 200 × 14mm	

Kern® [d] 10 & 20g, max 35 & 60kg Compact Veterinary Scale “EOB”, Stainless-steel Plate 315 × 305mm
With Wall Mount Display, Device · Non-slip Rubber Mat · Working Cover, 수의과(동물병원)용 저울

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

“EOB” Veterinary Scale

2-Year Warranty

CE



315 × 305 × h57mm

Specification

- Large LCD Display ; Digit Height 25 mm
- Dimensions of Display Device ; 235 × 114 × 51mm
- Operating Time ; up to 60h
- Permissible Ambient Temperature ; 10°C~30°C

Features



Adjusting Program CAL



Weighing Unit



Hold function



Battery Operation



AC 230V Adapter



Strain Gauges

* Simple and Convenient 4-key Operation

* Hold Function ; When the Weighing Conditions are Unstable, a Stable Weight is Calculated Determining an Average Value

* Integrated AUTO-OFF Function to Preserve the Batteries

* Non-slip Rubber Mat & Protective Working Cover Included

* Optional Battery Operation, 6 × 1.5 V AA not Included

Specifications & Ordering information

<M1>

Cat. No.	Weighing Range Max kg	Readability d g	Reproducibility g	Linearity g	Weighing Plate (w×d×h) & Weight	Price(W)
KE.EOB.35K10	35	10	10	± 20	315×305×57mm, 4kg	
KE.EOB.60K20	60	20	20	± 40		
KE.EOB.A01N	Optional ; Stand to Elevate Display Device, Approx. h450mm					



Stand ; Option

Figure ; Optional Stand(KERN EOB.A01N) Installed

Articles

Cat. No

Description

Kern® [d] 1g, max 15kg Slim Kitchen Scale "FGE-10K", with One Button Operation, 27.2 × 25 × h2.2cm
With Safety Glass Weighing Plate & Large LED Display, 슬림 강화유리 식품/주방용 저울, 간편 원버튼 영점 기능

NEW



2-Year Warranty

"FGE-10K" Slim Kitchen Scale
27.2 × 25 × h2.2cm, 500g



Package / 포장상태

CE

- Supplied by Kern® in Germany ■ CE Certified
- 원버튼 조작으로 간편하게 사용할 수 있는 슬림형 강화유리 식품용 저울 ; 주방이나 카페에 적합
- Ideal for Kitchen, Cafeteria, Guestonomy ■ Easy to Use-1 Button Control
- Safety Glass Weighing Plate - Easy and Hygienic Cleaning
- Low Overall Height for Secure Footing and Space-Saving Storage
- * Large High-contrasted Display that is Easy to Read
- * Secure and Non Slip Positioning with Rubber Feet
- * Comfortable Tare Function with Just One Keystroke ; Allows Taring of Containers and Sequential Ingredients without Calculation
- * Ideal for Accurate Portion Control in Cooking, Baking, and Everyday Use
- * High-capacity Sensor Design Ensures Stable Weighing with Large Containers and Heavy Ingredients
- * Integrated AUTO OFF Function to Preserve Battery Life * LCD Display ; Backlit, Digit Height 22mm
- * Batteries Included ; 3 × 1.5V, AAA * Dimensions & Weight (W × D × H) ; 275 × 250 × h22mm, 500g

Features



Adjusting Program CAL



Weighing Unit



Battery Operation



Strain Gauges

Specifications & Ordering information

<M1>

Cat. No.	Weighing -Capacity Max kg	Readability d g	Weighing Plate Dimensions	Ea. Price, ₩	Cs. Price, ₩
KE.FGE10K.3S05	15	1	275 × 205mm		5ea

, Medical / Healthcare - Body Scale 체중계



2-Year Warranty

CE



"MBD" Baby Scale



Optional ; Carrying Case

Kern® [d] 10g, max 15kg Ergonomically-Optimised Baby Scale "MBD", Ideal for Private
With Auto-off · Feeding · Hold · Memory-Function, Weight Displayed in kg or lb, 인기 보급형 신생아 체중계

- Supplied by Kern® in Germany ■ CE & ISO 9001 Certified ■ Using the Scales in Medicine is not Allowed.
- Particularly Flat and Lightweight Construction, So it is Easy to Transport and Store.
- * Large, Ergonomically-Optimised Lying Surface for Safe Weighing of Infants ; 650 × 350 × h55mm
- * Robust Plastic for Easy and Hygienic Cleaning. * Integrated AUTO-OFF Function to Preserve Batteries
- * Feeding Function ; Difference Function for Calculating the Weight Gain before and after Feeding
- * Hold Function ; When Patients do not lie completely still, a Stable Weight is Calculated Using an Average Weight
- * With Rubber Feet for Secure & Non-Slip Positioning * Batteries included (2 × 1.5V AAA), Operating Time up to 20h.

Specifications & Ordering information

<M1>

Cat. No.	Weighing Range Max kg	Readability d kg	Dimensions (W × D × H) & Weight	Price, ₩
KE.MBD.20K2S05	15	0.01	650 × 350 × 55mm, 2kg	

Ea. Price, ₩

KE.YTB.01S05 Optional ; Carrying Case with Handle and Shoulder Strap



2-Year Warranty

CE



"MBC" Baby Scale



* Figure ; Optional Height Rod Installed

(2) Optional Carrying Case

Kern® [d] 2 · 5 · 10g, max 6 · 15 · 20kg Elegantly-Precision Baby Scale "MBC"
Ideal for Weighing Neonates in Paediatrics, with Feeding & Hold-Function, 정밀형 신생아 체중계

- Supplied by Kern® in Germany ■ Medical approved in Germany(독일 의료기 인증 완료)
- CE & ISO 9001 Certified ■ Verification Class III (Verification is Optional) ■ ISO 13485 Certified
- Approval for Medical Use According to 93/42/EEC ■ Ideally Suited for Weighing Neonates in Paediatrics
- * Robust ABS Plastic for Easy and Hygienic Cleaning
- * Feeding Function ; Difference Function for Calculating the Weight Gain before and after Feeding
- * Hold Function ; When Patients do not lie completely still, a Stable Weight is Calculated Using an Average Weight
- * Baby Weighing Pan ; 600 × 280 × h55mm * Rubber Feet for Secure & non-Slip Positioning
- * Protective Working Cover * Batteries Included (6 × 1.5V AA)

Specifications & Ordering information

<M1>

Cat. No.	Weighing -Range Max kg	Readability d kg	Verif. Value e kg	Dimensions (W × D × H) & Weight	Price, ₩
KE.MBC.15K2DEM	6 / 15	0.002 / 0.005	0.002 / 0.005	600 × 407 × 120mm,	
KE.MBC.20K10M	20	0.01	0.01	5.5kg	

Ea. Price, ₩

KE.MBC.A01 Optional ; (1) Mechanical Height Rod, L80cm, 1mm Division
KE.MBC.A02 ; (2) Carrying Case, Only for Model without Height Rod



2-Year Warranty

"MGD" Body Scale
[d] 0.1kg, max 180 & 250kg
Tempered Glass Weighing Plate

Kern® [d] 0.1kg, max 180 & 250kg Popular Maximum Body Scale "MGD", with Large LCD
With Large · Flat Glass Weighing Plate, Ideal for Healthcare & Homecare, 가정용 스텝은 체중계

- Supplied by Kern® in Germany ■ CE & ISO 9001 Certified
- Large LCD Display and High-Contrast Display that is Easy to Read, Digit Height 40 & 50mm ; kg or lb
- Secure and Non-slip Positioning with Height-Adjustable Rubber Feet, Particularly Flat Design Facilitates Access
- * with Large · Flat · Sturdy Glass Weighing Plate * Easy and Hygienic Cleaning * Generous Weighing Capacity
- * Easy to Use and Time Efficient ; Simply Tap on it, Stand on it, Read the Weight, Integrated AUTO-OFF Function
- * Battery Included, "MGD-100K-1" ; 2ea×CR2032, "MGD-200K-1L" ; 3ea×1.5V AAA, Operating Time up to 100h

Specifications & Ordering information

<M1>

Cat. No.	Weighing -Range Max kg	Readability d kg	Weighing Plate Dimensions (W×D×H) & Weight	Price, ₩
KE.MGD.100K.1S05	180	0.1	300×300×19mm, 1.5kg	
KE.MGD.200K.1LS05	250	0.1	380×330×27mm, 2kg	



2-Year Warranty

"MFB" Body Scale
[d] 0.1kg, max 182kg
with Electrodes

Kern® [d] 0.1kg, max 182kg Multifunctional Body Scale "MFB", with Large LCD
With Large · Flat Glass Weighing Plate, Ideal for Healthcare & Homecare, 가정용 다기능 스텝은 체중계

- Supplied by Kern® in Germany ■ CE & ISO 9001 Certified
- Large LCD Display and High-Contrast Display that is Easy to Read, Digit Height 35mm ; kg or lb
- Secure and Non-slip Positioning with Height-Adjustable Rubber Feet, Particularly Flat Design Facilitates Access
- * Display Showing the Calorie Requirements for Maintaining the Current Body Weight Depending on the Level of Activity
- * with Large · Flat · Sturdy Glass Weighing Plate * Generous Weighing Capacity * Memory for up to 10 people
- * Easy to Use and Time Efficient ; Simply Tap on it, Stand on it, Read the Weight, Integrated AUTO-OFF Function
- * Battery Included ; 2ea×CR2032, Operating Time up to 100h

Specifications & Ordering information

<M1>

Cat. No.	Weighing -Range Max kg	Readability d kg	Weighing Plate Dimensions (W×D×H) & Weight	Price, ₩
KE.MFB.150K.100S05	182	0.1	300×300×25mm, 2kg	



(1) "AC/DC검용"
with Backlight
"MPD-M"

[d] 0.1kg,
max 250kg

2-Year Warranty



(2) "DC용"
"MPD-EM"

[d] 0.2kg,
max 250kg

2-Year Warranty



Step-on Function

Kern® [d] 0.1 & 0.2kg, max 250kg Step-on Floor Body Scale "MPD"
With Large Display, Ergonomically Optimized Design, 의료용 스텝은 체중계

- Supplied by Kern® in Germany
- Medical approved in Germany(독일 의료기 인증 완료)
- CE & ISO 9001 Certified ■ ISO 13485 Certified
- Approved as Medical Device according to 93/42/EEC
- * Easy to Use and Time Efficient ; Simply Tap on it, Stand on it, Read the Weight.
- * Secure and Non-slip Positioning with Height-Adjustable Rubber Feet.
- * Large Platform with Non-Slip and Wear-Resistant Surface
- * Suitable for Easy and Hygienic Cleaning
- * Level Indicator Make to Level the Balance Precisely
- * Robust and Ergonomically Optimized Design
- * Dimensions of Display ; 200×130×h60mm



- (1) AC/DC검용 with Backlight "MPD-M"
 - Large Backlight LCD Display, Digit Height 25mm
 - Mains Adapter External, Standard
 - Battery Operation Possible, 6ea×1.5V AA, Not Included, Operating Time up to 20h
- (2) DC용 "MPD-EM"
 - Large LCD Display without Backlight, Digit Height 25mm
 - Battery Operation, 6ea×1.5V AA, Included, Operating Time up to 20h

Specifications & Ordering information

<M1>

Cat. No.	Weighing -Range Max kg	Readability d kg	Verif. Value e kg	Weighing Plate Dimensions (W×D×H) & Weight	Price, ₩
KE.MPD.250K.100M	250	0.1	0.1	365×360×80mm, 10kg	
KE.MPD.200K.1EM	250	0.2	0.2		

DAIHAN™ [d] 50g, max 200kg Mobile Chair-type Scale “BAC-M501”, Hold & BMI Func., Backlight LCD
With Foldable Armrest & Footrest · Steerable Caster & Stopper · Rechargeable Battery, 의료용 이동식 의자형 체중계, 충전식 배터리, 체결량 지수

NEW

[d] 50g, max 200kg

Backlight LCD



3-Year Warranty

Mobile Chair-type Scale “BAC-M501”

- Supplied by DAIHAN Korea ■ 3-Year Warranty
- 국민건강보험공단 체중 지표(BMI) 기준을 적용하여 국내 건강검진과 연계한 관리가 용이함
- 바퀴 고정 기능 및 접이식 발판과 팔걸이로 보다 안전하게 의자에 앉고 내리기 편리함
- 환자 및 노약자 체중 측정에 최적화된 설계 및 기능 보유
- Ideal for Weighing Patients with limited Mobility in Hospital and Care Facilities.
- Foldable Armrests and Footrests for User Convenience.
- Steerable Caster & Stopper for Safe and Easy Mobility.
- Hold Function ; Ideal for Weighing Patients with Unstable Posture
- BMI(Body Mass Index) Analysis Function ;
- Classifies Body Mass Status(Underweight · Normal · Overweight)

* Rechargeable Battery for Weighing Device ;

- Included 6×AA Size Battery & 100-240V Universal Adapter for Charging
- Max. Operation Time ; 40hrs - Full Charging Time ; 12hrs

Specifications & Ordering information

<L1>


Spec.	Model	“BAC-M501”
Weighing Range. kg		0.05kg ~ 200kg
Read-Out [d]. g		50g
Reproducibility. g		50g
Weighing Platform. mm		370 × 370mm
Material	Frame	Powder Coated Steel
	Chair	PVC
Dimension & Weight		636×987×h950mm, 24kg
Packing & Gross Weight		640×990×h958mm, 25kg
DC Power & AC Adapter		Rechargeable Battery ; 6 × Size AA, Ni-MH 1.2V 2,000mAh Adapter ; AC100 ~ 240V
Cat. No.		DH.Bal8047
Price, Each. ₩		



Backlight LCD

Steerable Caster & Stopper
for Safe and Easy Mobility

Foldable Armrests and Footrests for User Convenience



Dimension Drawing

DAIHAN™ [d] 1mg, max 110g **Premium** Moisture Analyzer “BAL-110M”, Halogen Lamp, 0.00~100.00%, 40~199℃
With 100g Cali. Weight & Sample Plates, 5" Touch Screen, Counting Func., 3 Displays of %(g) · ℃ · Time, Ext-CAL, **조정밀/다용도 수분측정기**


“BAL-110M” [d] 1mg, max 110g



NEW
3-Year Warranty



Moisture Analyzer, “BAL-110M”



Intuitive User Interface
for Effortless Usability

Realtime Dynamic Curve



(2) “WBA-100” RS-232C
Thermal Printer



(3) “WBA-420” RS-232C
Thermal Label Printer

- Supplied by DAIHAN Korea
- CE Certified ■ PL(Product Liability) Insurance
- With Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-Information, and Traceable Data Base System
- 5 inch Touch Screen ■ Highlight Halogen Lamp / HBM Load Sensor
- Includes 50Pcs Aluminum Sample Plates(Φ90×h4mm, 3.05g)
- External Calibration ; Includes 100g Calibration Weight
- Display after Drying ; Moisture(%) - Moisture Content from Wet Weight
Dry Content(%) - Dry Weight from Wet Weight
ATRO(%) - 100 ~ 999%

- * Aluminum Shell and Stainless-steel Chamber
- * Real-time Dynamic Curve / Date and Time
- * Temperature / Time Set * Store History Set Value
- * Test History & Result on Intuitive Touch Screen
- * Built in RS232 for Printer (SL.WBA.100P, SL.WBA.420P)
- * Halogen Lamp Heating System for Ensuring Fast Drying
- * Supplied Multi - Point Calibration Function
- * 3 Test Modes ; Auto Mode · Soft Mode · Timing Mode
- * 3 Stop Modes ; Auto Stop · Timing Stop · Manual Stop

Specifications & Ordering information

<L1>

Spec.	Model	“BAL-110M”
Weighing Range. g		max 110g
Moisture Range / Read ability		0.00 ~ 100% / 0.01%
Dry Range / Read ability		Range ; 100.00 ~ 0.00%, Readability ; 0.01%
Read-Out [d]. g		0.001g
Weighing Plate. mm		Φ 90 mm
Ambient Temperature		5.0℃~35℃
Display		5-inch Touch Screen
Calibration		External Calibration ; 100g Calibration Weight Included
Heat Source / Temp Sensor		Halogen Lamp / PT100 Sensor
Temp Range / Resolution		40℃~199℃ / 1℃
Timer / Auto-stop		1~99min / 10s, Auto-stop
Number of Storage		20
Interface		RS232
Dimension & Net Weight		365×235×185mm / 6.8kg
Packing Size & Gross Weight		490×350×360mm / 7.5kg
Power Consumption		450 W
Power Supply & Cord Plug		AC 220V, 50 / 60Hz
Cat. No.		DH.Bal0110M
Price, Each. ₩		

Spare Sample Plates

DH.Dis8006 (1) Sample Plates, Φ90×h4mm, 3.05g, 25ml

Pk. Price, ₩
50ea

Optional Thermal Printers ; Specify the Balance Model for Printer Installation

* 프린터 구매 시 바란스 모델을 알려주어야 합니다.

Ea. Price, ₩

SL.WBA.100P (2) Thermal Printer “WBA-100”, Paper Roll
Included Cable Module, 1 Paper Roll

DH.WBA0084 Thermal Paper Package, 5 Paper Rolls, for “WBA-100”

SL.WBA.420P (3) Thermal Printer “WBA-420”, Label Paper Roll
Included Cable Module, 1 Label Paper Roll

DH.WBA0771 Thermal Label Paper Package, 5 Label Paper Rolls, for “WBA-420”

Articles

Cat. No

Description

Kern® [d] 1mg, max 110g High-Standard Moisture Analyzer "DAB", with 5 Memories & Graphic Display, 0~100%, 40°C~199°C
With Backlit LCD, 400W Halogen Quartz Glass Heater, Display of % · °C · Date · Time · Drying Program, 정밀형 다용도 수분 측정기

■ Supplied by Kern® in Germany

■ CE & ISO 9001 Certified

"DAB" [d] 1mg, max 110g

3-Year Warranty CE



Features



External Calibration



Memory



Data Interface RS 232



GLP/ISO record keeping



Weighing Unit



Power Supply



Strain Gauges


* Halogen Quartz Glass Heater 400 W

* Observation Window above the Sample, Useful During Initial Setting

* Internal Memory for Automatic Run of 5 Drying Programs and 15 Drying Processes Already Executed

* 50 Sample Plates Included

* Application Handbook ; On the Internet, You Will Find a Practical Application Handbook Containing Many Examples, Field Reports, Settings and Tips for Each KERN Moisture Analyzer



Backlit graphic display, digit height 14 mm

- ① Drying process active
- ② Active heating profile
- ③ Active switch-off criteria
- ④ Previous drying time
- ⑤ Current temperature
- ⑥ Current moisture content in %

Technical data

Readout	d	0.001g / 0.01 %
Weighing range	max	=maximum net weighing 110 g
Reproducibility, weight of sample	2 g	0.15 %
	10 g	0.02 %

Display after drying

Moisture [%] = Moisture Content (M) from Wet Weight (W) 0 - 100 %	
Dry Content [%] = Dry Weight (D) From (W) 0 - 100 %	
Moisture Content [g] (M)	Absolute value in [g]
Temperature Range	40 - 199 °C in steps of 1 °C
Display Switchable at all Times	

Drying modes

Standard drying	Gentle drying
Rapid Drying	

Switch off criteria

Automatic switch off (2 mg loss in weight in 45 s)
Time controlled switch off (3 min - 99 min 59 sec, 10 s increments)
Manual switch off at the press of a button

Call up of Measurement

Interval can be set from 1 sec - 10 min
(Only When used with Printer KERN YKB-01N or PC)

Dimensions

Overall dimensions	240 × 365 × h180 mm
Net weight	4.82 kg

Ordering information

KE.DAB100.3 Moisture Analyzer DAB 0~100%, 1mg ~ 110g

Optional Accessories

KE.DAB.A01	(1) Temperature Calibration Set Consist of Measuring Sensor and Display Device
KE.YKB.01	(2) Standard Printer, Capable, 0.75 lines/sec max 40 Characters/line, with 112mm Paper
KE.MLB.A05	Interface Cable, PC or Printer for "DAB" Model
KE.MLB.A01	Spare Sample dish, φ 90 mm Aluminum., 80/pk.
KE.RH.A02	Round Fiberglass Filter, Box with 100 pcs

<L1>

Ea. Price, ₩

Ea. Price, ₩

Pk. Price, ₩

80ea
100ea



(1) Temperature Calibration set
KERN "DAB-A01"



(2) Thermal Printer
KERN YKB.01

Kern® [d] 0.1 & 1g, max 500 & 1,000g Slim Pocket Balance for Reference "TGC", 10×13×h1.8cm, 200g
With Large Stainless Steel Weighing Plate & Shock Proof Plastic Cover, with Practical Tare Function, 슬림 포켓 디지털 바런스, 간편 영점기능



2-Year Warranty



Slim Pocket Balance ;
10×13×h1.8cm

Shock Proof
Plastic Cover



Incl. Shock Proof Plastic Cover
Can be used as a Weighing/Tare Cup

- Supplied by Kern® in Germany
- CE & ISO 9001 Certified
- High-Resolution Pocket Balance with Particularly Rapid Reaction and Stabilisation Time
- With Particularly Flat Design
- * Shock Proof Plastic Cover as Protection against Pressure and Dust & Can be Used as a Weighing/Tare Cup
- * Stainless Steel Weighing Plate, which Makes Cleaning Easy and Hygienic
- * Weighing Unit ; ct, dwt, g, gn, oz, ozt
- * LCD Display ; Backlit, Digit Height 12mm
- * Optimal Temperature Range ; 0℃~40℃
- * Batteries Included ; 2×CR2032, Operating Time up to 33h

Specifications & Ordering information

<L1>

Cat. No.	Weighing Capacity Max g	Readability d g	Dimensions & Weight (W×D×H, g)	Weighing Plate (φ)	Price, ₩
KE.TGC500.1S05	500g	0.1g	100×130×18mm, 200g	φ81mm	
KE.TGC1K.3S05	1,000g	1g			

Kern® [d] 0.1 & 1g, max 320 & 1,000g Pocket Balance & Calculator "CM", 8.5×13×h2.5cm, 180g
With Convenient 4-Key Operation & Hard Case Cover, Well Protected, 포켓 바런스 겸 계산기



2-Year Warranty



Pocket Balance & Calculator "CM" ;
8.5×13×h2.5cm

- Supplied by Kern® in Germany
- CE & ISO 9001 Certified
- Highly Portable, Convenient Use due to Combined Balance & Calculator
- * Simple and Convenient 4-Key Operation
- * Hard Case Cover as Protection against Pressure and Dust
- * Robust, Handy, Well Protected
- * Integrated Calculator
- * LCD Display ; Backlit, Digit Height 12mm
- * Optimal Temperature Range ; 18℃~25℃
- * Batteries Included ; 2×1.5V AAA, Operating Time up to 33h
- * Auto-Off Function to Preserve the Battery

Specifications & Ordering information

<L1>

Cat. No.	Weighing Capacity Max g	Readability d g	Dimensions & Weight (W×D×H, g)	Weighing Plate (W×D)	Price, ₩
KE.CM320.1N	320g	0.1g	85×130×25mm, 180g	70×80mm	
KE.CM1K.1N	1,000g	1g			

Kern® “E1”-Class Calibration Standard Weight, “OIML”-Class, International Accurate Class, High-degree of Accuracy With Certificate, DKD/GLP/IS, Standard, 최상급 표준 분동

- Supplied by Kern® in Germany
- Kern® Cylindrical Test Weight ; Comply with the “OIML R111-2014” in all Respects.
Kern® offers you the appropriate Test Weight Package ; Consisting of the Weight, Box, and DAkKS Calibration Certificate.
DAkKS = Deutsche Akkreditierungsstelle <German Accreditation Authority>.
- E1-Class ; Test Weight for Customers Who Require a High Degree of Accuracy for the Most Demanding Applications.
For High-resolution Balances with $d > 1,000,000$
Use Recommended with DAkKS Calibration Certificate Only.

DKD Certificate



“DKD” Calibration Certificate of Traceability

- * DKD (Deutscher Kalibrier Dienst = German Calibration Service)
- * ISO 9000FF, GLP, QS9000, GLM and ISO

※ **Note!!** The Detail of Accuracy Tolerance Information. Please see the Individual Weight Section



(1) Certified Weight-set

(1) “E1”-Class Calibration Certificated Standard Weight-sets with Certi. / Wooden Box / Non-magnetic Forceps / Brush / Glove

* Milligram Weight(1~500mg) in Plastic Box

KE.303.031.610	“E1”-Class Weight-Set, with Certi., 1 mg ~ 100 g (21 pcs)
KE.303.051.610	, 1 mg ~ 500 g (24 pcs)
KE.303.061.610	, 1 mg ~ 1 kg (25 pcs)
KE.303.081.610	, 1 mg ~ 5 kg (28 pcs)
KE.303.091.610	, 1 mg ~ 10 kg (29 pcs)
KE.304.031.610	“E1”-Class Weight-Set, with Certi., 1 g ~ 100 g (9 pcs)
KE.304.051.610	, 1 g ~ 500 g (12 pcs)
KE.304.061.610	, 1 g ~ 1 kg (13 pcs)
KE.304.081.610	, 1 g ~ 5 kg (16 pcs)
KE.304.091.610	, 1 g ~ 10 kg (17 pcs)

<L1>

Set. Price, ₩

(2) “E1”-Class Calibration Certificated Individual Weight with Certi. / Wooden Box



(a) Milligram Weight



(b) & (c) Gram and Kilogram Weight

KE.308.311.600	(a) “E1”-Class Weight, with Certi., 1 mg, Tol. ± 0.002 mg
KE.308.321.600	, 2 mg, Tol. ± 0.002 mg
KE.308.331.600	, 5 mg, Tol. ± 0.002 mg
KE.308.341.600	, 10 mg, Tol. ± 0.002 mg
KE.308.351.600	, 20 mg, Tol. ± 0.003 mg
KE.308.361.600	, 50 mg, Tol. ± 0.004 mg
KE.308.371.600	, 100 mg, Tol. ± 0.005 mg
KE.308.381.600	, 200 mg, Tol. ± 0.006 mg
KE.308.391.600	, 500 mg, Tol. ± 0.008 mg
KE.307.011.610	(b) “E1”-Class Weight, with Certi., 1 g, Tol. ± 0.010 mg
KE.307.021.610	, 2 g, Tol. ± 0.012 mg
KE.307.031.610	, 5 g, Tol. ± 0.015 mg
KE.307.041.610	, 10 g, Tol. ± 0.020 mg
KE.307.051.610	, 20 g, Tol. ± 0.025 mg
KE.307.061.610	, 50 g, Tol. ± 0.030 mg
KE.307.071.610	, 100 g, Tol. ± 0.05 mg
KE.307.081.610	, 200 g, Tol. ± 0.10 mg
KE.307.091.610	, 500 g, Tol. ± 0.25 mg
KE.307.111.610	(c) “E1”-Class Weight, with Certi., 1 kg, Tol. ± 0.5 mg
KE.307.121.610	, 2 kg, Tol. ± 1.0 mg
KE.307.131.610	, 5 kg, Tol. ± 2.5 mg
KE.307.141.610	, 10 kg, Tol. ± 5.0 mg
KE.307.151.610	, 20 kg, Tol. ± 10.0 mg
KE.307.161.610	, 50 kg, Tol. ± 25.0 mg

Ea. Price, ₩

Innotem® “F1”- & “M1”-Class Standard Weight, Set & Individual, “OIML”-Class, International Accurate Class

Polished Stainless-steel, <China-Made> “OIML” 표준분동세트 & 개별분동

■ Supplied by DAIHAN Korea

■ “OIML R111-2004” [Fr ; Organisation Internationale de Metrologie Legale ; International Organization of Legal Metrology] Offers the appropriate

* F1-Class Weight ; Test Weight for Analytical Balance & Precision Balance Verification (Class I / II ≤ 100,000e)

* M1-Class Weight ; Test Weight for Industrial and Commercial Scale Verification (Class III ≤ 10,000e)

※ Note! The Calibration Certificate ;
⇒ See the Next Page!

(1) “F1”- & “M1”- Class OIML R111-2004 Standard Weight-Set

Aluminum Case, Tweezer, Cotton cloth, Brush, 분동세트

<L1>

Weight Set	Pcs.	F1 Set		M1 Set	
		Cat. No.	Set. Price, ₩	Cat. No.	Set. Price, ₩
① 1mg ~ 500mg	12 Pcs	I7.We7051		I7.We7151	
② 1mg ~ 100g	21 Pcs	I7.We7052		I7.We7152	
③ 1mg ~ 500g	24 Pcs	I7.We7054		I7.We7154	
④ 1mg ~ 1kg	25 Pcs	I7.We7055		I7.We7155	
⑤ 1mg ~ 2kg	27 Pcs	I7.We7056		I7.We7156	
⑥ 1mg ~ 5kg	28 Pcs	I7.We7057		I7.We7157	
⑦ 1mg ~ 10kg	29 Pcs	I7.We7058		I7.We7158	
⑧ 1mg ~ 20kg	30 Pcs	I7.We7059		I7.We7159	
⑨ 1 kg ~ 5kg	4 Pcs	I7.We7060		I7.We7160	

Individual weight per Set	1	2	2	5	10	20	20	50	100	200	200	500	1	2	2	5	10	20	20	50	100	200	200	500	1	2	2	5	10	20																																																																				
Weight Set	mg	mg	mg	mg	mg	mg	mg	mg	mg	mg	mg	mg	g	g	g	g	g	g	g	g	g	g	g	g	g	kg	kg	kg	kg	kg	kg																																																																			
1 mg – 500 mg	Total weight												1.11g															①																																																																						
1 mg – 100 g																									211.11g							②																																																																		
1 mg – 500 g																																					1,111.11g							③																																																						
1 mg – 1 kg																																					2,111.11g							④																																																						
1 mg – 2 kg																																											6,111.11g							⑤																																																
1 mg – 5 kg																																																	11,111.11g							⑥																																										
1 mg – 10 kg																																																							21,111.11g							⑦																																				
1 mg – 20 kg																																																													41,111.11g							⑧																														
1 k g – 5 kg																																																																																					10,000g							⑨						

(2) “F1”- & “M1”- Class OIML R111-2004 Standard Individual Weight in a Case

Plastic Case, But 10 & 20kg in Aluminum Case, 개별분동



Weight Type	“F1” Individual Weight				“M1” Individual Weight			
	Cat. No.	Weight	± Tol mg	Ea. Price, ₩	Cat. No.	Weight	± Tol mg	Ea. Price, ₩
(a) Milligram	I7.We7001	1mg	± 0.020 mg		I7.We7101	1mg	± 0.20 mg	
	I7.We7002	2mg	± 0.020 mg		I7.We7102	2mg	± 0.20 mg	
	I7.We7003	5mg	± 0.020 mg		I7.We7103	5mg	± 0.20 mg	
	I7.We7004	10mg	± 0.025 mg		I7.We7104	10mg	± 0.25 mg	
	I7.We7005	20mg	± 0.03 mg		I7.We7105	20mg	± 0.3 mg	
	I7.We7006	50mg	± 0.04 mg		I7.We7106	50mg	± 0.4 mg	
	I7.We7007	100mg	± 0.05 mg		I7.We7107	100mg	± 0.5 mg	
	I7.We7008	200mg	± 0.06 mg		I7.We7108	200mg	± 0.6 mg	
	I7.We7009	500mg	± 0.08 mg		I7.We7109	500mg	± 0.8 mg	
(b) Gram	I7.We7011	1g	± 0.10 mg		I7.We7111	1g	± 1.0 mg	
	I7.We7012	2g	± 0.12 mg		I7.We7112	2g	± 1.2 mg	
	I7.We7013	5g	± 0.16 mg		I7.We7113	5g	± 1.6 mg	
	I7.We7014	10g	± 0.20 mg		I7.We7114	10g	± 2.0 mg	
	I7.We7015	20g	± 0.25 mg		I7.We7115	20g	± 2.5 mg	
	I7.We7016	50g	± 0.3 mg		I7.We7116	50g	± 3.0 mg	
	I7.We7017	100g	± 0.5 mg		I7.We7117	100g	± 5.0 mg	
	I7.We7018	200g	± 1.0 mg		I7.We7118	200g	± 10 mg	
	I7.We7019	500g	± 2.5 mg		I7.We7119	500g	± 25 mg	
(c) Kilogram	I7.We7021	1kg	± 5.0 mg		I7.We7121	1kg	± 50 mg	
	I7.We7022	2kg	± 10 mg		I7.We7122	2kg	± 100 mg	
	I7.We7023	5kg	± 25 mg		I7.We7123	5kg	± 250 mg	
	I7.We7024	10kg	± 50 mg		I7.We7124	10kg	± 500 mg	
	I7.We7025	20kg	± 100 mg		I7.We7125	20kg	± 1000 mg	
(d) Rectangle	I7.We7033	5kg	± 25 mg		I7.We7133	5kg	± 250 mg	
	I7.We7034	10kg	± 50 mg		I7.We7134	10kg	± 500 mg	
	I7.We7035	20kg	± 100 mg		I7.We7135	20kg	± 1000 mg	

DAIHAN™ KOLAS Calibration Certificate for "F1- & M1- Class" Weight, F1- & M1-급 분동 KOLAS 교정



■ "OIML R111-2004" [Fr ; Organisation Internationale de Metrologie Legale ; International Organization of Legal Metrology] Offers the appropriate Test Weight Package ; Consisting of the Weight, Box and DAIHAN Certificate.

■ Available Calibration Certificate Service for all Brand Weight

■ "KOLAS" Calibration Certificate of Traceability

* Accredited by KOLAS (Korea Laboratory Accreditation Scheme = Korea Calibration Service)

* KOLAS which Signed the ILAC-MRA (International Laboratory Accreditation Cooperation - Mutual Recognition Arrangement)

* Significant Instability or other Adverse Factor can affect the Validity of the Calibration (Overload, Temperature, Humidity etc.)

<L1>

Weight Range	F1 Weight Calibration Certificate		M1 Weight Calibration Certificate	
	Cat. No.	Price, ₩	Cat. No.	Price, ₩
1mg ~ 500g	DH.Weiz201.01		DH.Weiz201.11	
1kg ~ 5kg	DH.Weiz201.02		DH.Weiz201.12	
10kg	DH.Weiz201.03		DH.Weiz201.13	

* Weight-Set & Other Weights are available upon Customer's Request.

, Accessories 액세서리



Glass Weighing Bottle, with Ground Glass Stopper, Short- & Tall-form, 5~85ml

Excellent Chemical · Mechanical Resistance, Boro-glass 3.3, 글라스 다용도 평량병

(1) Short form, 20~82ml

<L1>

	Ea. ₩	Pk. Price, ₩
GL.264.202.01 20ml, $\phi 40 \times h30$ mm		10ea
GL.264.202.03 35ml, $\phi 50 \times h35$ mm		10ea
GL.264.202.05A 40ml, $\phi 60 \times h40$ mm		10ea
GL.264.202.04 50ml, $\phi 50 \times h50$ mm		10ea
GL.264.402.07 82ml, $\phi 70 \times h33$ mm		10ea

(2) Tall form, 5~85ml

	Ea. ₩	Pk. Price, ₩
GL.264.202.06 5ml, $\phi 20 \times h40$ mm		10ea
GL.264.202.07 15ml, $\phi 25 \times h50$ mm		10ea
GL.264.202.08 25ml, $\phi 30 \times h60$ mm		10ea
GL.264.202.09 60ml, $\phi 40 \times h80$ mm		10ea
GL.264.407.07 85ml, $\phi 50 \times h60$ mm		10ea



Kartell® Universal Use Weighing Bottle, PP, Transparent, 20~360ml

Excellent Chemical Resistance, Autoclavable, 0°C~125/140°C, <Italy-Made> PP 다용도 평량병

* PP(Polypropylene)

* Good Chemical & Temp-Resistance

* Multi-use ; Weighing, Sampling and Storage.

(1) Short form, 20~50ml

<L1>

	Ea. Price, ₩	Pk. Price, ₩
KA.342 20ml, $\phi 40 \times h30$ mm		10ea
KA.345 30ml, $\phi 50 \times h30$ mm		10ea
KA.346 50ml, $\phi 60 \times h34$ mm		10ea

(2) Tall form, 23~360ml

	Ea. Price, ₩	Pk. Price, ₩
KA.340 23ml, $\phi 30 \times h50$ mm		10ea
KA.343 60ml, $\phi 40 \times h70$ mm		10ea
KA.347 190ml, $\phi 60 \times h90$ mm		10ea
KA.348 360ml, $\phi 70 \times h120$ mm		10ea



Micro

Wheaton® 0.5 & 2.2ml Reusable Aluminum Micro-weighing Dish, with Smooth-surface

Ideal for Analytical Balance, $\phi 11$ & 20mm, <USA-Made> 알루미늄 미량 디쉬

<L1>

	dh pk. Price, ₩	Cs. Price, ₩
WH.370790 Weighing Dish, 0.5ml, $\phi 11 \times h6$ mm, 0.15 g/pcs	10ea	1000ea
WH.370792 , 2.2ml, $\phi 20 \times h8$ mm, 0.45 g/pcs	10ea	1000ea



(1) $\phi 75 \cdot 90 \cdot 102$ mm Reusable Al. Dish

$\phi 75 \cdot 90 \cdot 102$ mm Reusable Aluminum Dish, for Sampling & Weighing, 30~45ml

알루미늄 샘플/평량 디쉬, 재사용 가능

* Good for General Weighing, Testing and Moisture Drying Application for Moisture Analyzer

(1) $\phi 75 \cdot 90 \cdot 102$ mm Reusable Aluminum Dish

■ Supplied by DAIHAN Korea

DH.Dis8005	Al.Sample Dish, 30ml, $\phi 75 \times h8$ mm, 4.5g
DH.Dis8006	, 25ml, $\phi 90 \times h4$ mm, 3.05g New
DH.Dis8007	, 57ml, $\phi 102 \times h8$ mm, 3.8g

<L3>

Pk.	Price, ₩
100ea	
50ea	
50ea	

<L1>

(2) $\phi 90$ mm Premium Reusable Aluminum Dish, <Germany-Made>

KE.MLB.A01	Al.Sample Dish, 40ml, $\phi 90 \times h7$ mm, 3.7g
------------	--

Pk.	Price, ₩
80ea	



(1)

SciLab™ Disposable Aluminum Weighing Dish, Stackable, General- · Attached Handle- · Smooth Wall-type
Ideal for Weighing, Storage, Sampling and Drying, No Coating and Oil Free, 알루미늄 샘플/평량 디쉬

■ Supplied by SciLab Korea

* Ideal for Weighing, Storage, Sampling and Drying

* with Flat Bottom

<L3>

(1) Aluminum Weighing Dish, General-type

SL.Dis8030 Weighing Dish, Aluminum, General-type, 20ml, $\phi 51 \times h12.7\text{mm}$, 1g
SL.Dis8031 , 70ml, $\phi 70 \times h30.0\text{mm}$, 1.5g
SL.Dis8032 , 75ml, $\phi 71 \times h17.4\text{mm}$, 2g
SL.Dis8034 , 200ml, $\phi 127 \times h28.6\text{mm}$, 4g

Pk.	Price, ₩
100ea	
100ea	
100ea	
50ea	



(2)

(2) Aluminum Weighing Dish, Attached Handle-type

SL.Dis8040 Weighing Dish, Aluminum, Attached Handle-type, 8ml, $\phi 30 \times h10.0\text{mm}$, 0.3g
SL.Dis8042 , 20ml, $\phi 50 \times h12.7\text{mm}$, 1g
SL.Dis8043 , 40ml, $\phi 60 \times h16.0\text{mm}$, 1.5g
SL.Dis8044 , 80ml, $\phi 76 \times h16.0\text{mm}$, 2g

Pk.	Price, ₩
100ea	
100ea	
100ea	
100ea	



(3)

(3) Aluminum Weighing Dish, Smooth Wall-type

SL.Dis8051 Weighing Dish, Aluminum, Smooth Wall-type, 85ml, $\phi 100 \times h13\text{mm}$, 5g

Pk.	Price, ₩
100ea	



Low-Nitrogen Weighing Paper, 평량지/유산지

■ Supplied by DAIHAN Korea

* Ideal for Chemical, Medical, Laboratory

* Perfect Surface Treatment ; Powder Samples were not Asked

* NO Dust / Absorptiveness

<L1>

DH.WEP000 Weighing Paper, 50 × 50mm
DH.WEP0011 , 60 × 60mm
DH.WEP001 , 80 × 80mm
DH.WEP002 , 100 × 100mm
DH.WEP0021 , 120 × 120mm
DH.WEP003 , 150 × 150mm

Pk.	Price, ₩
500sheet	
500sheet	
500sheet	
500sheet	
500sheet	
500sheet	



(1) Square-type

Simport® Anti-Static Weighing Dish, PS, with Smooth Surface, Disposable, 9~355ml

Suitable for Clean Room, $-10^{\circ}\text{C} + 70/80^{\circ}\text{C}$, PS 일회용 플라스틱 평량디쉬, 정전기방지용

■ Supplied by Simport® in Canada

* Use ; Sampling for Liquid and Solid and Powder, Easily Transformed into Funnel-form

<L1>

(1) Square-type ; Accurate Pour-Outs With Minimum Sample Loss

SP.D250.1 Dish, Square type, 10ml, $40 \times 40 \times h8\text{mm}$
SP.D250.2 , 100ml, $78 \times 78 \times h25\text{mm}$
SP.D250.3 , 330ml, $127 \times 127 \times h25\text{mm}$

Pk.	Price, ₩	Cs.	Price, ₩
500ea		8pk	
500ea			
500ea			



(2) Boat-type

(2) Boat-type, Pour Dish ; for Static-affected Materials, Easy to Spout

SP.D251.1 Dish, Boat type, 12ml, $50 \times 37 \times h8\text{mm}$
SP.D251.2 , 140ml, $128 \times 76 \times h25\text{mm}$
SP.D251.3 , 270ml, $180 \times 117 \times h25\text{mm}$

Pk.	Price, ₩	Cs.	Price, ₩
250ea		8pk	
250ea			
250ea			



(3) Hexagonal-type

(3) Hexagonal-type ; Greater Balance Protection and Safety

SP.D252.1 Dish, Hexagonal type, 9ml, Top i.d. 35mm, Base i.d.25mm
SP.D252.2 , 58ml, Top i.d. 70mm, Base i.d.47mm
SP.D252.3 , 203ml, Top i.d.115mm, Base i.d.85mm
SP.D252.4 , 355ml, Top i.d.130mm, Base i.d.95mm

Pk.	Price, ₩	Cs.	Price, ₩
500ea		8pk	
500ea			
500ea			
500ea			

Black Weighing Dish, PS, Square-type & Diamond-form, Ideal for White Sample
With Smooth Surface, $-10^{\circ}\text{C} + 70/80^{\circ}\text{C}$, 검정 4각 & 다이아몬드형 웨이딩디쉬

<L3>

(1) Square-type, 10 · 100 · 300ml, 4각형

SL.Dis3021 Dish, Square-type, Black, $43 \times 43 \times h8\text{mm}$, 10ml
SL.Dis3022 , $89 \times 89 \times h25\text{mm}$, 100ml
SL.Dis3023 , $140 \times 140 \times h25\text{mm}$, 300ml

Pk.	Price, ₩
250ea	
250ea	
250ea	



(1) Square-type Black

(2) Kartell® Diamond-form, <UK-Made> 30 & 100ml, 다이아몬드형

KA.11101.09 Dish, Diamond-type, Black, $85 \times 56 \times h15\text{mm}$, 30ml
KA.11102.09 , $135 \times 96 \times h19\text{mm}$, 100ml

Pk.	Price, ₩	Cs.	Price, ₩
250ea		2pk	
250ea			

<L1>



(2) Diamond-form Black

Articles

Cat. No

Description

NEW



Lab Bristle Brush, For Cleaning Balance, Wood Handle with Metal Ferrule, 랩 브러쉬

<B3>

SL.Bru3011 Brush, Bristle ϕ 19 $\times$ L52mm, overall 28cm
 SL.Bru3012 , Bristle ϕ 24 $\times$ L56mm, overall 30cm

Ea.	Price, ₩



Weighing Forceps, High Grade Stainless-steel, with Plastic Tip, L100mm

Ideal for Handling Precision Weight, Non-Magnetic, Finished Surface, 웨이징 포셉/핀셋, 비자성/비부식

■ Supplied by SciLab Korea


■ 특급 스텐(18/10 크롬·니켈) 제품 ; 비자성, 비부식, 내산성 / 내약품성 / 내열성-우수

■ High Grade Stainless-steel of 18/10 Cr-Ni; Non-magnetic, Rust-free, Excellent Resistance of Acid / Chemicals / Heat

<L1>

SL.For7089 Forceps, Weighing, Plastic Grip Tip, STS, L100mm

Ea.	Price, ₩



Micro Spoon-Spatula, High Grade Stainless-steel, L150~200mm

With Bowl 3 $\times$ 5~5 $\times$ 9, Non-magnetic, 마이크로 스패츨러-스푼, 비자성/비부식

■ Supplied by SciLab Korea

* 18/10 Stainless-steel; Non-magnetic, Rust-free, Anti-Acid / Chemicals / Heat

<L1>

SL.Spa7011 Micro Spoon-Spatula, L150, Bowl-3 $\times$ 5, Blade-5 $\times$ 50, Stem- ϕ 2.5 mm
 SL.Spa7012 , L150, -5 $\times$ 7, -5 $\times$ 50, - ϕ 2.5 mm
 SL.Spa7013 , L180, -5 $\times$ 9, -5 $\times$ 50, - ϕ 2.5 mm
 SL.Spa7014 , L200, -5 $\times$ 9, -5 $\times$ 50, - ϕ 2.5 mm

Ea.	Price, ₩



Hammacher® Premium Spoon-Spatula, Rustproof Stainless-steel, L150~230mm

With 1-side Blade, Finished Polish, <Germany-Made> 프리미엄 스푼-스패츨러, 독일제, 비부식

* 18/8 Stainless-steel; Rust-proof, Anti-Acid / Chemicals / Heat

<L1>

HA.HSN133.15 Spoon-Spatula, L150, Blade - 10 $\times$ 50, Spoon - 17 $\times$ 38, Stem- ϕ 3mm
 HA.HSN134.18 , L180, - 10 $\times$ 55, - 17 $\times$ 38, - ϕ 3mm
 HA.HSN135.21 , L210, - 10 $\times$ 60, - 17 $\times$ 38, - ϕ 3mm
 HA.HSN136.23 , L230, - 10 $\times$ 60, - 17 $\times$ 38, - ϕ 4mm

Ea.	Price, ₩



Micro Half-round Scoop & Spatula, High Grade Stainless-steel, L130~210mm

Non-magnetic, High Polished, 미니 스텐 스크프-스패츨러, 비자성/비부식

* 18/10 Stainless-steel; Non-magnetic, Rust-free, Anti-Acid / Chemicals / Heat

<L1>

(1) Slender Stem-type * Supplied by SciLab Korea

SL.Sco7001 Scoop-Spatula, Stainless-steel, L130, Blade-9 $\times$ 40, Stem- ϕ 3.7mm
 SL.Sco7002 , L150, -9 $\times$ 45, - ϕ 3.3mm
 SL.Sco7003 , L185, -9 $\times$ 50, - ϕ 3.3mm
 SL.Sco7004 , L210, -9 $\times$ 50, - ϕ 3.3mm

Ea.	Price, ₩

(2) Bochem® Thin-wide Stem-type * Supplied by Bochem® in Germany

BC.3360 Mini Scoop, Stainless-steel, L170, Bowl-10 $\times$ 40, Stem thick-1mm

Ea.	Price, ₩



Premium Vibro-Scoop, High Grade Stainless-steel · PTFE-Coated, with Adjusting Knob, L190 & 220mm

Ideal for Handling High Value Materials, with Plastic-handle, 프리미엄 진동 스크프

■ 진동에 의한 미세한 용량 조절이 가능하여 밸런스 측정 및 고품질 시료를 핸들링하는데 적합

* Ideal for Use with Balance or Handling High Value Materials

* Adjustable Sample Offering with the Adjusting Knob

(1) Bochem® Stainless-steel Vibro-Spatula, with Plastic-handle, L190mm

* Supplied by Bochem® in Germany

* 18/10 Stainless-steel; Non-magnetic, Rust-free, Anti-Acid · Chemicals · Heat

<L1>

BC.3540 Vibro Scoop, Stainless-steel, L190, Blade-12 $\times$ 50mm

Ea.	Price, ₩

(2) Hammacher® PTFE-Coated Vibro-Spatula, with Plastic-handle, L220mm

* Made in Germany

* Teflon-PTFE ; Anti-Acid · Alkally · Chemicals · Corrosion, MP 260 °C

HA.HSN462.22 PTFE-coated Vibro Scoop, L220, Blade-7 $\times$ 80mm

Ea.	Price, ₩




(1) Stainless-steel



(2) PTFE-coated



(1) Stainless-steel



(2) PTFE-coated

Bochem® **Premium** Weighing Scoop, High Grade Stainless-steel · PTFE-Coated, L200~250mm
Non-magnetic, Rust-free, 웨칭 스푼, 비자성/비부식

■ Supplied by Bochem® in Germany

(1) **Premium** Stainless-steel Weighing Scoop, L200~250mm

* 18/10 Stainless-steel; Non-magnetic, Rust-free, Anti-Acid · Chemicals · Heat

BC.3600	Weighing-Scoop, STS, L 200, Head-31 × 70 , Stem-φ 5mm, 30g
BC.3602	, L 235, -31 × 105, -φ 5mm, 35g
BC.3603	, L 250, -31 × 120, -φ 5mm, 40g

<L1>


Ea.	Price, ₩

(2) PTFE-Coated Weighing Scoop, L235mm

* Teflon-PTFE ; Anti-Acid · Alkally · Chemicals · Corrosion, MP 260 °C

BC.8265	PTFE-coat. Weighing-Scoop, L235, Head-31 × 105, Stem-φ 5mm, 39g
---------	---

Ea.	Price, ₩



(1) Bowl/Dish only

(2) with Handling Stem

Weighing Scoop, High Grade Stainless-steel, with or without Handling Stem, L70~250mm
Non-magnetic, Rust-free, 경재형 웨칭 스푼, 비자성/비부식

■ Supplied by SciLab Korea

* 18/10 Stainless-steel; Non-magnetic, Rust-free, Anti-Acid/Chemicals/Heat

(1) Bowl/Dish-only, without Handling Stem, L70~120mm

SL.Sco7031	Weighing Scoop, Dish, STS, Dish-32 × 70mm , 13g
SL.Sco7032	, Dish-32 × 90mm , 18g
SL.Sco7033	, Dish-32 × 105mm, 23g
SL.Sco7034	, Dish-32 × 120mm, 26g

<L1>

Ea.	Price, ₩

(2) with Handling Stem, L200~250mm

SL.Sco7036	Weighing Scoop, Handled, STS, Bowl-31 × 70 , Overall-L200mm, 29g
SL.Sco7037	, Bowl-31 × 90 , -L220mm, 33g
SL.Sco7038	, Bowl-31 × 105, -L235mm, 45g
SL.Sco7039	, Bowl-31 × 120, -L250mm, 48g

Ea.	Price, ₩



Glass Weighing Scoop, 3 · 6 · 10ml, Easy to Use, Boro-glass 3.3, Glass 평량 스푼

* With Open Tubular Arm for the Filling of Powder into Narrow Mouth Bottles / Vessels

* With Flat Bottom for Safety Location

<L1>

SL.Sco1011	Weighing Scoop, Glass, 3ml , L(approx)-65mm
SL.Sco1012	, 6ml , -80mm
SL.Sco1013	, 10ml, -95mm

Ea.	Price, ₩	Pk.	Price, ₩
		10ea	
		10ea	
		10ea	



(1)



(2)

PP White Measuring Scoop, Flat-top, Autoclavable, 2~100ml

For Liquid · Powder · Solid, 0°C~125/140°C Stable, PP 백색 계량 스푼, 액체/고체겸용

■ 다양한 용량과 정밀한 용량, 안정적인 손잡이를 가진 스푼

* With Precision Formed Suitable Filling Edge and Comfortable

* 121 °C Autoclavable PP(Polypropylene)

* Good Chemical & Temp-Resistance Multi-Use ; Free-Standing and Flat-Top Construction

* Easily Readable Volume Quantities on the Upper Side of the Stable Handle.

* Suitable as Weighing Scoops

(1) DAIHAN Wisd™ PP Mini Flat Scoop **NEW**

■ Supplied by DAIHAN Korea

<L1>

	Ea.	₩	Pk.	Price, ₩
DH.Sco3001	10ml , L100mm		12ea	
DH.Sco3002	25ml , L135mm		12ea	
DH.Sco3003	50ml , L160mm		12ea	
DH.Sco3004	100ml, L200mm		12ea	

(2) VITLAB® **Premium** PP Flat Scoop

■ Supplied by VITLAB® in Germany

	Ea.	₩	Cs.	Price, ₩
VI.391.94	2ml , L60mm		12ea	
VI.392.94	5ml , L82mm		12ea	
VI.393.94	10ml , L100mm		12ea	
VI.394.94	25ml , L135mm		12ea	
VI.395.94	50ml , L160mm		12ea	
VI.396.94	100ml, L200mm		12ea	

DAIHAN Scientific

BATHS & BATH CIRCULATORS

Bath , Basic 100°C	P. 48
, Basic 100°C Circulation(Internal)	P. 49
, Basic 250°C Oil Bath	P. 50
, Cold Trap	P. 51
, Extraction Water Bath	P. 52
, Magnetic Stirring Bath	P. 53
, Shaking Water Bath	P. 54
, Viscosity Bath	P. 55
Bath Circulator , External, Chiller -30°C/-20°C~R/T New	P. 56~57
, Heated/Refrigerated -30°C/-20°C~200°C New	P. 58~59
, External/Internal, Heated +250°C	P. 60
, Heated/Refrigerated -35°C/-25°C~150°C	P. 61~62
Accessories , Circulation Hose New	P. 63~64
, Transfer Fluid	P. 65



Basic-type Bath
"WB", 100°C



Shaking Water Bath
"SB", 100°C



**External
Chiller / Cooling Circulator**
"CHI", -20°C~R/T



**External / Internal
Heated / Refrigerated Bath Circulator**
"CL", -35°C+150°C

DAIHAN™ Digital General Purpose Water Bath “WB”, Digital Fuzzy Control System, 6·11·22 Lit, up to 100℃, $\pm 0.1^\circ\text{C}$
With Stainless-steel Flat Lid, Back Light LCD, with Certi. & Traceability, 정밀 범용 항온수조, 디지털 퍼지 제어 시스템, 스텐리스 스틸 리드 포함



3-Year Warranty



(1) Stainless-steel Flat Lid Included

(2) Stainless-steel Dome Lid

(3) “BFC106”
Lid with Concentric Ring
Cover Sets for 6Lit(4) “BFC111”
Lid with Concentric Ring
Cover Sets for 11Lit(5) “BFC122”
Lid with Concentric Ring
Cover Sets for 22Lit Bath(6) “BTH106”
Tube Rack, 36Hole
for 6Lit Bath

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function

- * User's Self-Compensation Function to Control Difference between Real Temperature and Indicated Temperature ($\pm 10.0^\circ\text{C}$).
- * 304 Stainless-steel Bath & Lid for Superior Durability & High Thermal Efficiency/Constant.
- * Digital Timer Function(99hr 59min) ; ‘Delay Time’ and ‘Operation Time’ can be Controlled.
- * Storage Function for the Set Values of Temperature and Timer.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled).

Specifications & Ordering information

<L1>

Spec.	Model	“WB-6” Max. Capa. 6Lit	“WB-11” Max. Capa. 11Lit	“WB-22” Max. Capa. 22Lit
Capacity		6 Lit	11 Lit	22 Lit
Temp Range & Accuracy		Room Temp +5℃~100℃, $\pm 0.1^\circ\text{C}$		
Temp Resolution		0.1℃-Display, 0.1℃-Control		
Temp Uniformity & Probe		$\pm 1.0^\circ\text{C}$, PT100		
Heating Power		800 W	1 kW	2 kW
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller		Digital Fuzzy Controller with Jog - Shuttle Switch (Turn + Push)		
Safety Circuit		Over-Temp / -Current protection, Sensor Error Detection		
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel		
Lid		included Stainless-steel Flat Lid, Optional Dome Lid		
Drain Valve		None		Mounted($\phi 9.5\sim 11\text{mm}$)
Dimension (w×d×h)	Effective Space	300×155×h110mm	302×240×h110mm	500×295×h110mm
	Internal	300×155×h150mm	302×240×h150mm	500×295×h150mm
	External	325×180×h260mm	325×265×h260mm	535×330×h260mm
Packing Size (mm) and Gross Weight (kg)		350×210×h310mm, 8kg	350×300×h310mm, 10kg	570×360×h310mm, 14kg
Electrical Reqt.s.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz		
Cat. No.	230V	DH.WB000106	DH.WB000111	DH.WB000122
	120V	DH.WB000206	DH.WB000211	DH.WB000222
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Optional Lid and Tube Rack, Stainless-steel

- DH.WB010206 (2) Dome Lid, “BD206”, for 6 Lit, 324×176×h100mm
 DH.WB010211 “BD211”, for 11 Lit, 324×263×h100mm
 DH.WB010222 “BD222”, for 22 Lit, 533×328×h100mm
 DH.WB030106 (3) Lid with 2× Concentric Ring Cover Sets, “BFC106”, for 6 Lit, Stainless-steel, 330×182×h20mm
 DH.WB030111 (4) Lid with 2× Concentric Ring Cover Sets, “BFC111”, for 11 Lit, Stainless-steel, 329×269×h20mm
 DH.WB030122 (5) Lid with 6× Concentric Ring Cover Sets, “BFC122”, for 22 Lit, Stainless-steel, 539×334×h20mm
 DH.WB020106 (6) Tube Rack, “BTH106”, for 6 Lit, 36Hole× $\phi 17\text{mm}$, Stainless-steel, 220×125×h78.5mm
 DH.WB020206 Tube Rack, “BTH206”, for 6 Lit, 15Hole× $\phi 32\text{mm}$, Stainless-steel, 220×125×h78.5mm
 DH.WB020111 Tube Rack, “BTH111”, for 11 Lit, 63Hole× $\phi 17\text{mm}$, Stainless-steel, 220×205×h78.5mm
 DH.WB020211 Tube Rack, “BTH211”, for 11 Lit, 25Hole× $\phi 32\text{mm}$, Stainless-steel, 220×205×h78.5mm

Ea. Price, ₩

Water Bath with Validation (IQ, OQ) Service

- DH.WB040106 Validation Service(IQ, OQ), “BV106”, with “WB-6”
 DH.WB040111 “BV111”, with “WB-11”
 DH.WB040122 “BV122”, with “WB-22”

Ea. Price, ₩

DAIHAN™ Digital Precise Internal Circulation Water Bath “WCB”, 6 · 11 · 22 Lit, with Certi. & Traceability
With Stainless-steel Flat Lid, Powerful Circulation Pump, up to 100 °C, ±0.1 °C, 5 Lit/min, 정밀 내부순환 항온수조



up to 100 °C



3-Year Warranty



Stainless-steel Flat Lid Included

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function
- * Ideal for Biotechnology, Clinical, Environmental, Medical, Petroleum, Food Industry, Pharmaceutical and Industrial Applications. * Alarm function ; Error Status and Timer-end
- * Warning and Automatic Stop Function in case of Liquid Shortage.
- * Powerful Circulation Pump Assures Temp Uniformity ; Internal and External Circulation.
- * Connection with Evaporators and Viscometers is applicable.

Specifications & Ordering information

<L1>

Spec.	Model	“WCB-6” Max. Capa. 6Lit	“WCB-11” Max. Capa. 11Lit	“WCB-22” Max. Capa. 22Lit
Capacity		6 Lit	11 Lit	22 Lit
Temp Range & Accuracy		Room Temp +5 °C~100 °C, ±0.1 °C		
Temp Resolution		0.1 °C-Display, 0.1 °C-Control		
Temp Uniformity & Probe		±0.2 °C, PT100		
Heating Power		800 W	1 kW	2 kW
Pump		5 Lit / min (Internal / External)		
Circulation Nozzle		id.7×od.11 mm (id.9/od.12mm silicone tube connected for inner circulation)		
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end		
Controller & Display		Digital Fuzzy Controller with Back Light Digital LCD		
Safety Circuit		Over-Temp/-Current protection, Sensor Error & Low Level Detection		
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel		
Drain Valve		None		Mounted(φ9.5~11mm)
Dimension (w×d×h)	Effective Space	120×155×h150 mm	170×240×h150 mm	372×295×h150 mm
	Internal	300×155×h150 mm	302×240×h150 mm	500×295×h150 mm
	External	330×185×h385 mm	335×275×h385 mm	540×335×h410 mm
Packing Size (mm) and Gross Weight (kg)		350×210×h410mm, 11kg	360×300×h420mm, 13kg	580×380×h470mm, 18kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz		
Cat. No.	230V	DH.WCB00106	DH.WCB00111	DH.WCB00122
	120V	DH.WCB00206	DH.WCB00211	DH.WCB00222
Price, Each. ₩				


* Other Specifications are available upon Customer's Request.

※ 운임 별도

DAIHAN™ Digital Precise Circulation Immersion Pump Head “WCB-H”, 11 & 22 Lit, up to 100 °C
With Digital Fuzzy Control System, Certi. & Traceability, Power Circulation Pump, 5 Lit/min, 정밀 항온 순환기



up to 100 °C



3-Year Warranty

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function
- * Ideal for Biotechnology, Clinical, Environmental, Medical, Petroleum, Food Industry, Pharmaceutical and Industrial Applications. * Alarm function ; Error Status and Timer-end
- * Powerful Circulation Pump Assures Temp Uniformity ; Internal and External Circulation.

Specifications & Ordering information

<L1>

Spec.	Model	“WCB-11H” Max. Capa. 11Lit	“WCB-22H” Max. Capa. 22Lit
Temp Range & Sensor		Room Temp +5 °C~100 °C, PT100	
Heating Power & Pump		1 kW, 5 Lit/min	2 kW, 5 Lit/min
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end	
Controller & Display		Digital Fuzzy Controller with Back Light Digital LCD	
Safety Circuit		Over-Temp/-Current protector, Sensor Error & Low Level Detection	
Dimension		130×165×h320mm	
Packing Size (mm) and Gross Weight (kg)		400×280×h200mm, 8kg	400×280×h200mm, 10kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	
Cat. No.	230V	DH.WCB10011	DH.WCB10022
	120V	DH.WCB20011	DH.WCB20022
Price, Each. ₩			



Side



Rear

DAIHAN™ Digital High Temp Oil Bath “WHB”, Digital Fuzzy Control System, 6 · 11 · 22 Lit, up to 250°C, $\pm 1.5^\circ\text{C}$
With Stainless-steel Flat Lid, Back Light LCD, with Certi. & Traceability, 정밀 고온 오일 배스, 디지털 퍼지 제어 시스템, 스테인리스 스틸 리드 포함

CE

3-Year Warranty
up to 250°C



Digital High Temp Oil Bath, 11Lit, “WHB-11”



Stainless-steel Flat Lid Included

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy.
- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function.

- * Rustless 304 Stainless-steel Bath (t 1.2mm) & Ceramic Fiber Insulator (t 50mm)
- * User's Self-Compensation Function to Control Difference between Real Temp. and Indicated Temp.
- * Storage Function for the Set Values of Temperature & Timer.
- * Alarm & Digital Timer ; 'Delay Time' and 'Operation Time' can be Controlled.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled).

Specifications & Ordering information

<L1>

Spec.	Model	“WHB-6” Max. Capa. 6Lit	“WHB-11” Max. Capa. 11Lit	“WHB-22” Max. Capa. 22Lit
Capacity		6 Lit	11 Lit	22 Lit
Temp. Range & Accuracy		Room Temp +5°C to 250°C, $\pm 1.5^\circ\text{C}$		
Temp. Uniformity & Probe		$\pm 3.0^\circ\text{C}$, PT100		
Heating Power		1.5 kW	2.0 kW	2.5 kW
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)		
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection		
Material		Internal ; Stainless-steel (#304), External ; Powder Coated steel Insulation ; 50 mm Thickness Ceramic Fiber		
Lid		Stainless-steel Flat Lid Included		
Drain Valve		Mounted		
Dimension (w×d×h)	Effective Space	270×180×h110mm	300×250×h110mm	400×300×h140mm
	Mounted	270×180×h150mm	300×250×h150mm	400×300×h180mm
	External	425×305×h340mm	425×375×h340mm	525×425×h370mm
Packing Size (mm) & Gross Weight		470×400×h500mm, 23kg	600×490×h500mm, 25kg	600×490×h500mm, 39kg
Electrical Reqt.		1 Phase, AC 120V, 60Hz or 230V, 50 / 60Hz		
Cat. No.	230V	DH.WHB00106	DH.WHB00111	DH.WHB00122
	120V	DH.WHB00206	DH.WHB00211	DH.WHB00222
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

High-Temp. Oil Bath with Validation (IQ, OQ) Service

DH.WHB40106 Validation Service(IQ, OQ), “HBV106”, with “WHB-6”
DH.WHB40111 , “HBV111”, with “WHB-11”
DH.WHB40122 , “HBV122”, with “WHB-22”

Ea. Price, ₩

※ See page 65 for the Heating Transfer Fluids ;

- ▶ SUPER THERM FG-20® -10°C +300°C Heating Transfer Fluids
- ▶ SYN THERM 22L® -5°C +320°C Heating Transfer Fluids
- ▶ SYN THERM 7008® -40°C +400°C Cooling/Heating Transfer Fluids

DAIHAN™ Digital Precision Cold Trap Bath “WCT”, -40°C & -80°C , 10 Lit, with Certi. & Traceability
With 2 Glass Cold Trap, Digital PID Control System, Back-Light LCD, CFC-free Refrigeration System, 초저온 동결 트랩

CE



3-Year Warranty

Patented Jog-Shuttle Control System

※ DAIHAN Recommendation
; Combination of the Cold Trap Bath(WCT) with the Vacuum Oven and Vacuum Pump.



(1) Glass Cold Trap
“CTG1000”

- Supplied by DAIHAN Korea
 - CE Certified
 - PL(Product Liability) Insurance
 - with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
 - Digital PID Controller
 - CFC-free (R-404A) Refrigeration System
 - Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function
- * Ideal for exceptional Cooling System like Collecting Evaporation Vapor or Hazardous Vapor from Vacuum Line.
- * Durable Stainless-steel Bath and Powder Coated steel
- * Digital Timer Function ; ‘Delay Time’ and ‘Operation Time’ ; Can be Controlled.
- * User’s Self-Compensation Function to Control Difference between Real Temperature and Indicated Temperature.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

Specifications & Ordering information

<L1>

Spec.	Model	“WCT-40” Max. Temp. -40°C	“WCT-80” Max. Temp. -80°C
Temperature Range		$-40^{\circ}\text{C} \sim -20^{\circ}\text{C}$	$-80^{\circ}\text{C} \sim -65^{\circ}\text{C}$
Bath Volume		10 Lit	
Bath Volume & Cold Trap Type		Stainless-steel 304, 2 Glass Trap	
Refrigerator		1/2 HP	1 HP
Power Consumption		900 W	1.8 kW
Refrigerator & Refrigerant		R-404A	HFC
Display		Digital LCD with Back-Light Function	
Controller		PID Controller with Jog-Shuttle Switch (Turn+Push)	
Drain Valve		Mounted	
Safety Circuit		Over Current protection, Sensor Error Detection, Leakage Breaker	
Material		Internal ; Stainless-steel (# 304) External ; Powder Coated steel, 1.2t, Doublepaint & Baked Insulation ; 20mm Thickness Polyurethane	
Dimension (w×d×h)		Internal ; 220×220×h210mm, External ; 377×453×h705mm	
Cold Trap Inlet Valve Size		od. 10 mm	
Packing Dimension (mm) and Gross Weight (kg)		560×500×h1,010, 67kg	680×640×h1,130, 90kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or 230V, 50 / 60Hz	
Cat. No.	230V	DH.WCT00140	DH.WCT00180
	120V	DH.WCT00240	DH.WCT00280
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 윤임 별도

Spare Glass Cold Trap

DH.WCT11000 (1) Glass Cold Trap “CTG1000” for “WCT”
DURAN® Borosilicate Glass α 3.3, $\phi 106 \times h249\text{mm}$,
Safety-Pressure Resistant 1.5bar at 250°C


Ea. Price, ₩

DAIHAN™ Extraction Water Bath “WEB”, with Concentric Ring Cover Sets & Holding Frame, 4 · 6 · 8 Place
Ideal for Gentle Extraction · Soxhlet Extraction · COD Determination, Digital Fuzzy Control, Back Light LCD, up to 100 °C, ± 0.2 °C
다용도 추출용 항온수조, 4 · 6 · 8 구, 컨센트릭 링 커버 세트 & 홀딩 프레임 포함



3-Year Warranty

4-Place Extraction Water Bath “WEB-4”



3-Year Warranty

6-Place Extraction Water Bath “WEB-6”



(1) Lid with Concentric Ring Cover Set (Included)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
 - with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
 - Digital Fuzzy Control System Implementing Superior Temperature Accuracy
 - Concentric Ring Cover Sets make Effective Steaming Work
 - Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function
- * 4 · 6 · 8 Extraction Place are available.
* Includes Concentric Ring Cover Set & Holding Frame
* User's Self-compensation Function to Control Difference between Real Temperature and Indicated Temperature (± 10.0 °C).
* High Quality 304 Stainless-steel Bath and Holding Frame for Superior Durability.
* Storage Function for the Set Values of Temperature and Timer.
* Alarm Function ; Error Status and Timer-end.
* Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled).

Specifications & Ordering information

<L1>

Spec.	Model	“WEB-4” 4 Place	“WEB-6” 6 Place	“WEB-8” 8 Place
Flask Capacity		250 ~ 500ml almost all Round Flasks		
No. of Extraction Holes		4 Place	6 Place	8(4×2) Place
Temp. Range & Accuracy		Room Temp +5 °C to 100 °C, ± 0.2 °C		
Temp. Uniformity & Sensor		± 1.0 °C, PT100		
Heating Power		1.4 kW	1.8 kW	1.4 kW×2
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)		
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection, Leakage Breaker		
Material		Internal ; Stainless-steel (#304), External ; Powder Coated steel Holding Frame ; Stainless-steel		
Lid		Concentric Ring Cover Set Included		
Top Opening Dia.		φ36 · φ65 · φ96 · φ110mm		
Bath Depth		122mm		
Drain Valve		Mounted		
Dimension (w×d×h)	Internal	610×160×h160mm	910×140×h160mm	610×290×h160mm
	External	680×245×h275mm	975×245×h275mm	720×375×h275mm
Packing Size (mm) and Gross Weight (kg)		1,100×350×h400, 18 kg	1,250×350×h400, 25 kg	900×450×h400, 34 kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or 230V, 50 / 60Hz		
Cat. No.	230V	DH.WEB01004	DH.WEB01006	DH.WEB01008
	120V	DH.WEB02004	DH.WEB02006	DH.WEB02008
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Lid with Concentric Ring Cover Set and Holding Frame Set

- DH.WEB11000 (1) Lid with Concentric Ring Cover Set, “EBC100”, for all “WEB” models
- DH.WEB21004 Holding Frame Set, “EBF100”, for “WEB-4”, 2ea × 800mm Pole & 2ea × 700mm Pole
- DH.WEB21006 Holding Frame Set, “EBF200”, for “WEB-6”, 2ea × 800mm Pole & 2ea × 1,020mm Pole
- DH.WEB21008 Holding Frame Set, “EBF300”, for “WEB-8”, 4ea × 800mm Pole & 4ea × 740mm Pole


Ea. Price, ₩

SciLab™ Integrated Thermostatic Magnetic Stirring Bath “MBWCL”, up to 200℃, 0~2,000rpm, 4 & 6.5Lit
With PID Controller, Digital Display, Stainless-steel Heater, can be Used for Oil-/Water-Bath, 다목적 자력교반 항온 수조



3-Year Warranty

Integrated Thermostatic Magnetic Stirring Bath
“MBWCL-3”



Drawing of “MBWCL-3”

- Supplied by SciLab Korea
- Adopt ANSI 304 Stainless-steel Heater which can be used for Oil and Water Bath.
- Powerful Magnetic Stirring of High Performance Direct Current Brushless.
- Dual Sensor, Temperature can be converted and Displayed in Flume and Reaction Vessel.
- * Strong Magnetic Stirring makes the Stirrings in the Flume Rotate Synchronistically to ensure the Even and Precision of High Viscosity Solution Temperature.
- * The Performance of Direct Current Brushless Machine is Stable and can be Used Continuously.

Specifications & Ordering information

<L1>

Spec.	Model	“MBWCL-3” Max. 3L Flask Input	“MBWCL-5” Max. 5L Flask Input
Bath Volume & Size		4 Lit, $\phi 220 \times \text{depth} 110\text{mm}$	6.5 Lit, $\phi 254 \times \text{depth} 130\text{mm}$
Temp. Range & Stability, Digital		Room Temp. $5^{\circ}\text{C} \sim 200^{\circ}\text{C}$, $\pm 1^{\circ}\text{C} : 1$	
Speed, Knob Setting		0~2000 rpm	
Max. Flask Input		3,000ml	5,000ml
Heater, ANSI 204 Stainless-steel		500W	1,050W
Temp. Controller & Display		PID Controller, Digital Display	
Dimension		$260 \times 280 \times \text{h} 230\text{mm}$	$260 \times 280 \times \text{h} 260\text{mm}$
Electrical Reqts.		220V, 50/60Hz	
Cat. No.		SL.MbWCL3	SL.MbWCL5
Price, Each. ₩			

* Included Retractable Multi-Ring Opening Lid ; This Lid can Adjust Caliber to Match the External Diameter of the Reaction Vessel and Round Flask, Decrease Volatilization of the Refrigerant and Heat Exchange with the Outside Environment, Realize Stable Low-Temperature Environment.



DAIHAN™ Digital Precise Shaking Water Bath “SB”, Reciprocating Motion, 18 · 30 · 45 Lit

With Universal Spring Rack, Digital Fuzzy Control, up to 2 Lit Flask, 20~250 rpm, up to 100℃, ±0.1℃, with Certi. & Traceability

정밀 진탕 항온수조, 왕복 진탕 운동, 디지털 피드백 제어, 고성능 셰이킹 엔진 탑재



- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Controller Provides Precise Temperature Control
- Digital Feedback Controller Provides Precise Speed Control
- Hi-Performance Shaking Engine ensures the Best Smooth Shaking Motion and Minimum Noise and Vibration
- Easy to Disassemble Spring Rack for Cleaning
- Patented Shaking Mechanism ; Minimum Noise and Vibration
- Patented Jog-Dial Control System and Digital LCD with Back-Light

- * Ideal for Direct Immersion of Samples or Shaking in the Field of Biotechnology and Laboratories.
- * Hybridization Applications, Cell Cultures, Cell Aeration, Increasing Solubility Rates, Molecular Biology Assays, Bacterial Cultures
- * Stainless-steel Bath for Superior Durability & High Thermal Efficiency
- * Universal Stainless-steel Spring Rack Included ; for various kind of Flasks and Tubes (up to 2 ℓ)
- * User's Self-Compensation Function to Control Difference between Real Temperature and Indicated Temperature (± 10.0℃)
- * Over Temp, Over Current Protection and Sensor Error Detection
- * Digital Timer Function ; 'Delay Time' and 'Operation Time' ; Can be Controlled

Specifications & Ordering information

<L1>

Spec.	Model	“SB-18” Max. Capa. 18Lit	“SB-30” Max. Capa. 30Lit	“SB-45” Max. Capa. 45Lit
Capacity		18 Lit	30 Lit	45 Lit
Max.Q'ty of Erlenmeyer Flask		250ml×6ea, 500ml×4ea, 1000ml×1ea	250ml×10ea, 500ml×8ea, 1000ml×3ea	250ml×18ea, 500ml×15ea, 1000ml×8ea
Temp Range & Accuracy		Room Temp +5℃ to 100℃, ±0.1℃		
Temp Resolution		0.1℃-Display, 0.1℃-Control		
Temp Uniformity & Probe		±0.5℃, PT100		
Heating Power		1.4 kW	1.8 kW	2.4 kW
Shaking Speed & Resolution		20 ~ 250 rpm, 1rpm-Display/1rpm-Control		
Shaking Stroke (mm)		15, 20, 30 Available (20 as factory Default)		
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller Digital		Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)		
Safety Circuit		Over-Temp/-Current protection, Sensor Error Detection		
Rack		Universal-type Stainless-steel Spring Rack Included		
Drain Valve		Mounted		
Material		Internal ; Stainless-steel (#304), External ; Powder Coated steel		
Dimension (w×d×h)	Internal (mm)	300×300×h200	500×300×h200	600×380×h200
	External (mm)	563×353×h340	763×353×h340	863×433×h340
Packing Size (mm) & Gross Weight (kg)		620×420×h450, 24 kg	820×420×h450, 26 kg	920×500×h450, 32 kg
Electric Reqts.		1 Phase, AC 120V, 60 Hz or 230V, 50 / 60Hz		
Cat. No.	230V	DH.WSB03018	DH.WSB03030	DH.WSB03045
	120V	DH.WSB03218	DH.WSB03230	DH.WSB03245
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Universal Spring Rack

- DH.WSB55018 (1) Spring Rack, “SBR518”, for 18 Lit, 245×270×h110mm
 DH.WSB55030 , “SBR530”, for 30 Lit, 445×270×h110mm
 DH.WSB55045 , “SBR545”, for 45 Lit, 545×350×h110mm

Ea. Price, ₩

Optional Stainless-steel Dome Lid

- DH.WSB02018 (2) Dome Lid, “SBD018”, for 18 Lit, 363×353×h69mm
 DH.WSB02030 , “SBD030”, for 30 Lit, 563×353×h69mm
 DH.WSB02045 , “SBD045”, for 45 Lit, 663×433×h69mm

Ea. Price, ₩

Shaking Water Bath with Validation (IQ, OQ) Service

- DH.WSB31018 Validation Service(IQ, OQ), “SBV018”, with “SB-18”
 DH.WSB31030 , “SBV030”, with “SB-30”
 DH.WSB31045 , “SBV045”, with “SB-45”

Ea. Price, ₩



(1) Universal Spring Rack
(Included)



(2) Dome Lid, Stainless-steel (Optional)

DAIHAN™ Digital Precise Viscosity Bath “WVB”, for 5 Viscometers, Max. 30Lit/min, up to 100℃, $\pm 0.1^\circ\text{C}$
 Stainless-steel Lid with 5 Holes for Viscometer Holders, Available Reverse & Routine-type Viscometer, Transparent Window
 투시형 정밀 점도용 항온수조, 5×점도계 사용가능, 디지털 퍼지 제어 시스템, 5-구 스텐리스 스틸 리드포함



Max. 30Lit/min Viscosity Bath “WVB-30”
 Includes Flat Lid for Viscometer Holder.
 (Viscometer Holder ; Option)



Fig1. Flat Lid with 5 Viscometer Holders



(1) Viscometer Holder “WVB02010” (Optional)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Used for Precise Measurement with Capillary Viscometers
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Powerful Circulation Pump ensures the Highest Uniformity of Temperature
- Viscometer Holder ; Available Reverse & Routine-type Viscometer
- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function.
- * Stainless-steel Lid Contains Five Holes for Insertion of Viscometer Holders
- * User's Self-compensation Function to Control Difference between Real Temperature and Indicated Temperature
- * Storage Function for the Set Values of Temperature and Time
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Includes Flat Lid for Viscometer(without Viscometer Holder)
- * Over Temp., Over Current Protection and Sensor Error Detection
- * Digital Timer Function ; 'Delay Time' and 'Operation Time' ; Can be Controlled

Specifications & Ordering information

<L1>

Spec.	Model	“WVB-30” Max. Flow 30Lit/min
Range & Accuracy		Room Temp +5℃ to 100℃, $\pm 0.1^\circ\text{C}$
Temp. Resolution		0.1℃-Display, 0.1℃-Control
Temp. Uniformity & Probe		$\pm 0.2^\circ\text{C}$, PT100
Heating Power		2 kW
Pump		Magnetic Drive, Max. Flow 30 Lit/min
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end
Display		Digital LCD with Back-Light Function
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection
Lid		Included 5 Holes Stainless-steel Lid for Viscometer Holder
Drain Valve		Mounted
Viewing Window Size(mm)		230×180
Material		Internal ; Stainless-steel (#304), External ; Powder Coated steel Window ; Vacuum Sealed and Tempered Safety Pair Glasses
Dimension (w×d×h)	Effective Space (mm)	330×300×h260mm
	Internal (mm)	350×300×h300mm
	External (mm)	585×375×h420mm
Packing Size (mm) & Gross Weight (kg)		640×450×h500mm, 24kg
Electric Reqts.		1 Phase, AC 120V, 60Hz or 230V, 50 / 60Hz
Cat. No.	230V	DH.WVB01030
	120V	DH.WVB02030
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Stainless-steel Flat Lid for Viscometer Holder

DH.WVB11030 Flat Lid, “WVB030”, for 30 Lit, w380×d330, 5× ϕ 50mm Holes

Ea. Price, ₩

Optional Viscometer Holder

DH.WVB02010 (1) Viscometer Holder, “WVB02010”
 for Ostwald, Ubbelohde
 DH.WVB02011 Viscometer Holder, “WVB02011”
 for “witeg” Ubbelohde, Cannon-Fenske
 DH.WVB02012 Viscometer Holder, “WVB02012”
 for Thread-type Ubbelohde

Ea. Price, ₩

Viscosity Bath Validation (IQ, OQ) Service

DH.WVB31030 Validation Service(IQ, OQ), “WVB030”, with “WVB-30”

Ea. Price, ₩

DAIHAN™ -30℃~R/T Classical Chiller “CHI”, Refrigerated External Circulator, Fill-5 · 10 · 20 · 30 · 50Lit with Insulated STS Hose
Ideal for Evaporator, Reactor &c. Cooling Line, Lift 6 · 12m, Cooling Capa 0.6~3.5kW, 냉각 써클레이터/칠러, 단열 튜빙 포함



3-Year Warranty



Internal Bath



Connection Ports



Insulated Flexible STS Hose Included

- Supplied by DAIHAN Korea ■ CE Certified ■ Max. Lift ; 12m
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information.
- Ideal for Wide-range Reaction System and Evaporator.
- Precise Temperature Control System, $\pm 1^{\circ}\text{C}$
- Microprocessor-based Control ensures Stable and Consistent Thermal Performance
- Stainless-steel Internal Bath for Prevention of Corrosion
- Includes Two Insulated Flexible Stainless-steel Hoses (0.7m & 1m)
- * Various Filling Volumes Available ; 5 · 10 · 20 · 30 · 50 Lit
- * Effective Cooling System with Rapid Temperature Stabilization
- * Useful for 1~50Lit Various Reactors and Evaporators.
- * Equipped with Caster Wheels for Safe and Easy Mobility
- * High Quality Circulation Pump ; 30-Lit/min Max Flow Rate

Specifications & Ordering information

<L1>

Spec.	Model	“CHI5-30”	“CHI10-30”	“CHI20-30”	“CHI30-30”	“CHI50-30”
Temp. Range & Accuracy		-30℃-Range, ±1℃ at Set Temperature				
Cool-Fluid. Filling Volume		5Lit max	10Lit max	20Lit max	30Lit max	50Lit max
Useful for 1~50 Lit Evaporators, etc.		for 1~5Lit Evaporator	for 5~20Lit Evaporator	for 20~30Lit Evaporator	for 50Lit Evaporator	for 100Lit Evaporator
Refrigerator Capacity		0.55kW	0.75kW	1.09kW	1.5kW	3.75kW
Material		Internal ; Stainless-steel, External ; Powder Coated steel				
Lid & Top Opening Dia.		φ180mm		φ240mm	φ320mm	
Fluid In/Outlet Port		ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel				
Circulation Hose		0.7m & 1m, Included				
Pump Power, Lift Height, m	Max. Flow Rate	20Lit/min				30Lit/min
	Max. Lift	8m				12m
Cooling Capacity	at +15℃	0.9kW	2.0kW	2.8kW	4.6kW	10.5kW
	at 0℃	0.7kW	1.7kW	2.3kW	3.5kW	9.3kW
	at -30℃	0.3kW	0.7kW	0.6kW	0.6kW	0.5kW
External Dimension (mm)		370×470×h680	450×520×h800	570×490×h860	640×540×h1000	740×640×h1190
Electrical Reqts.		1 Phase 230V, 50 / 60 Hz				
Cat. No.		DH.CHI05.30	DH.CHI10.30	DH.CHI20.30	DH.CHI30.30	DH.CHI50.30
Price, Set. ₩						

※ 운임 별도

Spare Insulated Flexible Stainless-steel Hose

SL.Tub7033 Insulated Flexible STS Hose Only, 1/2", id17×od26mm, Overall od60mm, L0.7m
SL.Tub7034 , L1m

Ea. Price, ₩

※ For More Information about Circulator Hoses ;
SEE the “Bath Circulators, Accessories” Section !

DAIHAN™ -20℃~R/T Premium Chiller “CHI”, Heavy-duty Refrigerated External Circulator, Fill-17·29·47 Lit
Ideal for Evaporator/Reactor &c. Cooling Line, Lift 27m, Cooling Capa 0.87·1.3·3.0 kW, 다용도 냉각 써클레이터/칠러



Assembly-type Safety Connection Ports/Pipe for Fluid In/Out-let ;
ISO RC 1/2" Thread, Located On Rear-Upper Panel

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Wide-range Reaction System and Evaporator. ■ Max. Lift 27m
- Assembly-type Safety Connection Ports/Pipe for Fluid In/Out-let ; ISO RC 1/2" Thread, Stainless-steel.
- Digital PID Control System Implementing Superior Temperature Accuracy.
- Embedded External Nozzle Part for the Prevention of Nozzle Damaged.
- Stainless-steel Internal Bath for Prevention of Corrosion.
- * Wide-range Filling Volume 17Lit · 29Lit · 47Lit
- * Effective Cooling Capacity
- * Useful for 1~50Lit Various Reactors and Evaporators. * Emergency Stop Button for safety.
- * High Quality Circulation Pump ; 30-Lit/min Max Flow Rate and 3.3-bar Max. Pressure.
- * Built in Pump Pressure Gauge and Compressor Pressure Gauge.

Specifications & Ordering information

<L1>

Spec.	Model	“CHI-20” Filling Volume Max. Capa. 17Lit	“CHI-30” Filling Volume Max. Capa. 29Lit	“CHI-50” Filling Volume Max. Capa. 47Lit
Temp Range & Stability		-20℃~Room Temp, ±0.2℃(at 15℃)-Stability		
Cool-Fluid. Filling Volume		17 Lit max	29 Lit max	47 Lit max
Useful for 1~50Lit Evaporators, etc.		for 1~20Lit Evaporator	for 10~30Lit Evaporator	for 20~50Lit Evaporator
Refrigerator Capacity(HP)		1 HP	1.5 HP	2 HP
Lid & Top Opening Dia.		Φ 100 mm		
Fluid In/Out-let Connector		ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel		
Pump Power, Lift Height, m	Max. Flow Rate	30 Lit/min		
	Max. Pressure & Lift	3.3 bar Max. Lift ; 27m		
Cooling Capacity	at +20℃	1.7 kW	2.8 kW	6.5 kW
	at 0℃	0.87 kW	1.3 kW	3.0 kW
	at -20℃	0.32 kW	0.45 kW	1.2 kW
External Dimension		450×630×h650 mm	550×830×h975 mm	605×975×h1265 mm
Material		Internal ; Stainless-steel, External ; Powder Coated steel		
Electrical Rqts.		1 Phase 230V, 50 / 60 Hz		3 Phase 380V
Cat. No.	60Hz	DH.Chi0020	DH.Chi0030	DH.Chi0050
	50Hz	DH.Chi0025	DH.Chi0035	DH.Chi0055
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

※ Cooling/Heating Transfer Fluids ;
See the Page 65.

- ▶ SYN THERM 7008® -40℃+400℃
- ▶ SUPER THERM FG-20® -10℃+300℃
- ▶ SYN THERM 22L® -5℃+320℃

Note ! Cooling Transfer Fluid Should be an
Antifreezing Solution.

※ Insulated Circulation Hoses ; See the
“Bath Circulators, Accessories” Section !



▶ Figure ; One of the Applications in the Large Volume Evaporator Line.

DAIHAN™ -30℃+200℃ **Premium** Precise Refrigerated Heating External Circulator “CRH”, Closed-type, Fill-10·18·25 Lit Ideal for Evaporator/Reactor &c. Heating & Cooling Line, with Pre-Cooling Sys, Used with External Direct Contact K-type Temp. Probe(Optional) With High Quality Magnetic Pump, Flow 20·40·60 Lit/min, Lift 38m, ±0.2℃, 정밀 고온/냉동 겸용 외부 써클레이터-밀폐형

Advantage of Closed-type Circulators ; 밀폐형 써클레이터의 장점

- (1) Small Power Consumption ; Less than Bath Circulators
- (2) Eco-Friendly ; Non-Evaporation of Heat Transfer Fluid and Non-Smell of Fluid
- (3) Uniformity of Temperature Control



Fig. “CRH” in the Reactor System.

※ Note ; Cooling Transfer Fluid Should be an Antifreezing Solution

- Supplied by DAIHAN Korea
- PL (Product Liability) Insurance
- Max. Lift ; 38m
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Patented Jog-Dial Control System and Digital LCD with Back-Light Function.
- * Innovated Powerful Circulation Pump ensures Temp. Uniformity ; Internal and External Circulation (Max Flow Rate ; 20·40·60 Lit/min)
- * Assembly-type Safety Connection Ports/Pipe (for Fluid Drain/In/Out-let ; ISO RC 1/2" Screw Stainless-steel.)
- * With High Quality Shock-proof Caster for Convenient and Safe Moving.
- * Embedded External Nozzle Part for the Prevention of Nozzle Damaged.
- * Stainless-steel Internal Bath for Prevention of Corrosion.
- * Independently 3-systems ; Refrigerating System, Heating System, Pre-Cooling System.
- * Urethane Foam Insulation * More Faster Cool-down and Heat-up Time
- * Warning and Auto-Stop Function in case of Liquid Shortage.
- * Connection with Evaporator or Viscometer is applicable.
- * 외부 순환을 위한 추가 커넥터가 필요합니다. (ex. PT 1/2" Hose nipple Connector)



Specifications & Ordering information

<L1>

Spec.	Model	“CRH-20” Useful for 5~20Lit Reactor	“CRH-50” Useful for 20~50Lit Reactor	“CRH-100” Useful for 50~100Lit Reactor
Temp. Range & Stability		-30℃+200℃-Range, ±0.2℃(at 15℃)-Stability		
Heat-Fluid Filling Volume		10 Lit max	18 Lit max	25 Lit max
Useful for 5~100Lit Reactors		for 5~20Lit Reactor	for 20~50Lit Reactor	for 50~100Lit Reactor
Refrigerator Capacity		4.5 kW	7.5 kW	18 kW
Heater		2.5 kW	6.0 kW	12 kW
Material		Internal ; Stainless-steel, External ; Powder Coated steel		
Lid & Top Opening Dia.		φ80 mm		
Fluid Drain/In/Outlet Port		ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel		
Overflow Port		φ12mm One-Touch Fitting		
Pump Power, Lift Height, m	Max. Flow Rate	20 Lit/min	40 Lit/min	60 Lit/min
	Max. Pressure & Lift	Max Pressure ; 4 bar, Max. Lift ; 38m		
Cooling Capacity	at +15℃	8.1 kW	13.1 kW	22.7 kW
	at 0℃	4.3 kW	7.3 kW	12.6 kW
	at -25℃	1.4 kW	3.1 kW	5.45 kW
Optional External Temp. Sensor K-type Probe		Can be used with External Direct Contact K-type Temp. Probe. (Optional)		
External Dimension		600×850×h1200mm	650×900×h1350mm	750×900×h1450mm
Electrical Reqts.		3Phase 380V		
Cat. No.	60Hz	DH.WHR0020	DH.WHR0050	DH.WHR0100
	50Hz	DH.WHR0025	DH.WHR0055	DH.WHR0150
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도



► Insulated Circulation Hose


⇒ See the “Bath Circulators, Accessories” Section !

※ Cooling/Heating Transfer Fluids ⇒ See the Page 65.

DAIHAN™ -20℃+200℃ Refrigerated Heating External Circulator "HCR", Closed-type, Fill-10·20·30·50Lit
With Insulated Stainless-steel Hose, Ideal for Evaporator/Reactor &c. Heating & Cooling Line, Flow 20·25·30Lit/min, Lift Max. 8 & 12m,
고온/냉동 겸용 외부 써클레이터-밀폐형, 단열 튜빙 포함

**Advantage of Closed-type Circulators ;
밀폐형 써클레이터의 장점**

- (1) Small Power Consumption ;
Less than Bath Circulators
- (2) Eco-Friendly ; Non-Evaporation of Heat Transfer
Fluid and Non-Smell of Fluid
- (3) Uniformity of Temperature Control



Connection Ports



Insulated Flexible STS Hose Included

- Supplied by DAIHAN Korea ■ CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information.
- Precise Temperature Control System, $\pm 1^\circ\text{C}$
- Microprocessor Based Control ensures Stable and Consistent Thermal Performance.
- Ideal for Chemical, Pharmaceutical, and Food Research from Lab to Pilot Scale
- Includes Two Insulated Flexible Stainless-steel Hoses (0.7m & 1m)
- * Integrated Heating and Cooling System with Rapid Temperature Stabilization
- * Optimized for Repeated Thermal Cycling * Sealed Structure minimizes Leakage Risk.
- * Internal Stainless-steel Circulation Path ensures Corrosion Resistance and Durability.
- * External Nozzle Protection Design prevents Physical Damage.
- * Equipped with Caster Wheels for Safe and Easy Mobility
- * Connection with Evaporator or Viscometer is Applicable.

Specifications & Ordering information

<L1>

Spec.	Model	"HCR10-20"	"HCR20-20"	"HCR30-20"	"HCR50-20"
Temp. Range & Accuracy	-20℃+200℃-Range, $\pm 1^\circ\text{C}$ at Set Temperature				
Heat-Fluid Filling Volume		Max. 10Lit	Max. 20Lit	Max. 30Lit	Max. 50Lit
Useful for 5~100Lit Reactors		for 5~20Lit Reactor	for 20~30Lit Reactor	for 50Lit Reactor	for 100Lit Reactor
Refrigerator Capacity		0.75kW	1.09kW	1.5kW	3.75kW
Heater		2kW	3kW	4.5kW	6kW
Material	Internal ; Stainless-steel, External ; Powder Coated steel				
Lid & Top Opening Dia.	$\phi 45\text{mm}$				
Fluid In/Outlet Port	ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel				
Circulation Hose	0.7m & 1m, Included				
Pump Power, Lift Height, m	Max. Flow Rate	20Lit/min	25Lit/min	30Lit/min	
	Max. Lift	8m		12m	
Cooling Capacity	at +15℃	2.1kW	2.8kW	5.6kW	10.5kW
	at 0℃	1.5kW	1.9kW	5.5kW	8.5kW
	at -30℃	0.5kW	0.6kW	0.9kW	0.9kW
External Dimension (mm)	620×540×h920			710×580×h1050	770×670×h1180
Electrical Reqts.	220V, 50/60Hz				3 Phase 380V
Cat. No.	DH.HCR10.20	DH.HCR20.20	DH.HCR30.20	DH.HCR50.20	
Price, Set. ₩					

※ 운임 별도

Spare Insulated Flexible Stainless-steel Hose

SL.Tub7033 Insulated Flexible STS Hose Only, 1/2", id17×od26mm, Overall od60mm, L0.7m
SL.Tub7034 , L1m

Ea.	Price, ₩

※ For More Information about Circulator Hoses ;
SEE the "Bath Circulators, Accessories" Section !

DAIHAN™ 250℃ Digital Precise High-temp Bath Circulator “CH”, ±0.1℃, 8·12·22·30 Lit

With Stainless-steel Flat Lid, Digital Fuzzy Control System, Certi. & Traceability, Flow 16 Lit/min, Lift 2.8m

Ideal for Heating Line of Facility, External Circulating Connect od.12.5 & id.10.5mm, 정밀 고온 배스겸용 고온 써클레이터

CE up to 250℃



3-Year Warranty




External Circulation Nozzle



Flat Lid Included



Fig. Multi-use High Temp Bath Circulator with Modular Reactor System



(1) Insulation Lid (Optional)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance.
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.

- Advanced Heating System Ensures Temperature Performance.

- Patented Jog-Dial Control System and Capacitive Touch Button.

- * Digital Fuzzy Control System Implementing Superior Temperature Accuracy.

- * Innovated Powerful Circulation Pump ensures Temp Uniformity ; External and Internal Circulation (Max Flow Rate ; 16 Lit/min).

- * Controlled 5 Steps Circulation speed(Flow Rate).

- * 304 Stainless-steel Bath and Lid for Superior Durability and High Thermal Efficiency.

- * Compensation Function to Control Difference between Real Temp and Indicated Temp

- * Digital Timer Function ; 'Delay Time' and 'Operation Time' can be Controlled.

- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled).

Specifications & Ordering information

<L1>

Spec.	Model	“CH-8” Max. Capa. 8 Lit	“CH-12” Max. Capa. 12 Lit	“CH-22” Max. Capa. 22 Lit	“CH-30” Max. Capa. 30 Lit
Bath Volume		8 Lit	12 Lit	22 Lit	30 Lit
Pump		Max Flow Rate ; 16 Lit/min, (Controlled 5 Steps Flow Rate Speed) Max Pressure ; 5 psi, Max Lift ; 2.8m			
Temp Range & Accuracy		Room Temp +5℃~250℃, ±0.1℃			
Temp Resolution		0.1℃-Display, 0.5℃-Control			
Temp Uniformity & Probe		±0.2℃, PT100			
Heating Power		1.0 kW	1.5 kW	2.2 kW	2.5 kW
Heat-up Time		50 min. to 100℃			
External Circulation Nozzle		id.10.5 / od.12.5 mm			
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end			
Display		Digital LCD with Back-Light Function			
Controller		Digital Fuzzy Controller with Jog - Dial Knob (Turn+Push)			
Safety Circuit		Over-Temp/-Current protection, Sensor Error			
Material		Internal ; Stainless-steel (#304), External ; Powder Coated steel Insulation ; 50mm Thickness Ceramic Fiber			
Lid & Drain Valve		Included Stainless-steel Flat Lid & Ball Valve			
Dimension (w×d×h)	Effective Space	180×130×h155	250×130×h155	340×250×h155	340×250×h200
	Mounted	180×320×h155	250×320×h155	340×440×h155	340×440×h200
	External	303×443×h569	373×443×h569	463×563×h569	463×563×h614
Packing Size (mm) & Gross Weight(kg)		400×470×h630 19 kg	400×470×h630 21 kg	490×600×h680 23 kg	490×600×h680 37 kg
Electrical Reqt.s.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz			
Cat. No.	230V	DH.WCH00408	DH.WCH00412	DH.WCH00422	DH.WCH00430
	120V	DH.WCH00508	DH.WCH00512	DH.WCH00522	DH.WCH00530
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

(1) Optional Microporous Ceramic Insulation Lid T15mm(10mm Microporous Ceramic + 5mm Ceramic Paper)

DH.WCH.L008 Insulation Lid, for CH-8 , 288×158×T15mm
DH.WCH.L012 for CH-12, 358×158×T15mm
DH.WCH.L022 for CH-22, CH-30, 448×316×T15mm

Ea. Price, ₩

※ Cooling/Heating Transfer Fluids ; See the Page 65.

- ▶ SYN THERM 7008® -40℃ +400℃
- ▶ SYN THERM 22L® -5℃ +320℃

- ▶ SUPER THERM FG-20® -10℃ +300℃

Note ! Cooling Transfer Fluid Should be an Antifreezing Solution.

※ Insulated Circulation Hoses ; See the “Bath Circulators, Accessories” Section !

DAIHAN™ -35℃ + 150℃ Digital Precise Refrigerated/Heating Bath Circulator “CL”, ±0.2℃

With Flat Lid, Digital Fuzzy Control, CFC-free, Certi. & Traceability, Capa. 8 · 12 · 22 · 30 Lit, Flow 25Lit/min, Lift 4m

Ideal for Cooling/Heating Line of Facility, External Circulating Connect od.12 & id.10mm, 정밀 냉동/가열 배스 겸용 써클레이터

CE -35℃ + 150℃



3-Year Warranty



External Circulation Nozzle



Stainless-steel Flat Lid Included

※ Cooling/Heating Transfer Fluids ;
See the Page 65.

- ▶ SYN THERM 7008® -40℃ + 400℃
- ▶ SUPER THERM FG-20® -10℃ + 300℃
- ▶ SYN THERM 22L® -5℃ + 320℃

Note ! Cooling Transfer Fluid Should be an Antifreezing Solution.

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Advanced Refrigeration System Ensures Temperature Performance
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Patented Jog-Dial Control System and Capacitive Touch Button

- * Innovated Powerful Circulation Pump ensures Temp Uniformity ; External and Internal Circulation (Max Flow Rate ; 25 Lit/min).
- * RS232C Interface for Remote Monitoring and Controlling with PC
- * All models Built in Caster
- * Much Faster Cool-down Time
- * Urethane Foam Insulation
- * Proper Liquid Level Line in Bath
- * Connection with Evaporator or Viscometer is applicable
- * User's Self-compensation Function to Control Difference between Real Temperature and Indicated Temperature (± 10.0℃).
- * Warning and Auto-Stop Function in case of Liquid Shortage.
- * Auto-controlled Refrigeration System for Compressor Protection and Power Saving.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled)

Specifications & Ordering information

<L1>

Spec.	Model	“CL-8” Max. Capa. 8 Lit	“CL-12” Max. Capa. 12 Lit	“CL-22” Max. Capa. 22 Lit	“CL-30” Max. Capa. 30 Lit
Temp Range & Accuracy		-35℃+150℃, ±0.1℃			
Temp Resolution		0.1℃-Display, 0.1℃-Control			
Temp Uniformity & Probe		±0.2℃ at -20℃, PT100			
Pump		Max Flow Rate ; 25Lit/min(Controlled 5 Steps Flow Rate Speed) Max Pressure ; 7.5 psi, Max Lift ; 4m			
Bath Volume		8 Lit	12 Lit	22 Lit	30 Lit
Heating & Refrigerator		1 kW, 1/2HP	1.5 kW, 1/2HP	2 kW, 7/8HP	2.2 kW, 1HP
Cooling Capacity	at +20℃	500 W	530 W	860 W	1030 W
	at 0℃	340 W	370 W	620 W	760 W
	at -20℃	200 W	210 W	400 W	500 W
Refrigerant		CFC -Free(R-404A) Refrigeration System			
Total Power Consumption		1.7 kW	2.2 kW	3 kW	3.3 kW
Compressor Control		Intelligent Control of Refrigerator for preventing Compressor Overload			
External Circulation Nozzle		id.10.5 / od.12.5mm(For Inner Circulat. ; to Connect directly the in-/outlet Nozzles)			
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end			
Display		Digital LCD with Back-Light Function			
Controller		Digital Fuzzy Controller with Jog-Dial Knob (Turn+Push)			
Safety Circuit		Over-Temp/-Current protection, Sensor Error Detection, Leakage Breaker			
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel Insulation ; 40 mm Thickness Urethane Foaming			
Lid & Drain Valve		Includes Stainless-steel Flat Lid & Ball Valve			
Dimension (w×d×h)	Effective Space (mm)	175×153×150	230×181×150	290×220×200	330×290×200
	Internal (mm)	175×320×150	230×350×150	290×390×200	330×460×200
	External (mm)	400×523×807	400×523×807	483×563×850	494×633×906
Packing Size (mm)		530×630×985	530×630×985	615×670×1025	630×740×1075
& Gross Weight (kg)		50 kg	55 kg	65 kg	78 kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz			
Cat. No.	230V	DH.WCL00408	DH.WCL00412	DH.WCL00422	DH.WCL00430
	120V	DH.WCL00508	DH.WCL00512	DH.WCL00522	DH.WCL00530
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

※ Insulated Circulation Hoses ; See the
“Bath Circulators, Accessories” Section !

Refrigerated Circulator with Validation (IQ, OQ) Service


- DH.WCL30308 Validation Service(IQ, OQ), “CLBV208”, with “CL-8”
- DH.WCL30312 , “CLBV212”, with “CL-12”
- DH.WCL30322 , “CLBV222”, with “CL-22”
- DH.WCL30330 , “CLBV230”, with “CL-30”

Ea. Price, ₩

(1) Optional Microporous Ceramic Insulation Lid
T15mm(10mm Microporous Ceramic+5mm Ceramic Paper)

- DH.WCR.L008 Insulation Lid, for CL-8 , 240×178×T15mm
- DH.WCR.L012 for CL-12, 316×206×T15mm
- DH.WCR.L022 for CL-22, 396×245×T15mm
- DH.WCR.L030 for CL-30, 436×316×T15mm

Ea. Price, ₩



(1) Insulation Lid (Optional)

DAIHAN™ -25℃ + 150℃ Digital Precise Refrigerated/Heating Bath Circulator "CR", ±0.2℃

With Flat Lid, Digital Fuzzy Control, CFC-free Refrigeration, Certi. & Traceability, 8 · 12 · 22 · 30 Lit, Flow 25Lit/min, Lift 4m

Ideal for Cooling/Heating Line of Facility, External Circulating Connect od.12 & id.10mm, 정밀 냉동/가열 배스 겸용 써클레이터

CE -25℃ + 150℃



3-Year Warranty



External Circulation Nozzle



Stainless-steel Flat Lid Included

※ Cooling/Heating Transfer Fluids ; See the Page 65.

- ▶ SYN THERM 7008® -40℃ + 400℃
- ▶ SUPER THERM FG-20® -10℃ + 300℃
- ▶ SYN THERM 22L® -5℃ + 320℃

Note ! Cooling Transfer Fluid Should be an Antifreezing Solution.

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.

- Advanced Refrigerating System Ensures Temperature Performance.
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy.
- Patented Jog-Dial Control System and Capacitive Touch Button

- * Innovated Powerful Circulation Pump ensures Temp Uniformity ; External and Internal Circulation (Max Flow Rate ; 25 Lit/min).
- * RS232C Interface for Remote Monitoring and Controlling with PC. * All models Built in Caster.
- * Much Faster Cool-down Time * Urethane Foam Insulation * Proper Liquid Level Line in Bath
- * User's Self-compensation Function to Control Difference between Real Temperature and Indicated Temperature (± 10.0℃).
- * Warning and Auto-Stop Function in case of Liquid Shortage.
- * Auto-controlled Refrigeration System for Compressor Protection and Power Saving.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled).

Specifications & Ordering information

<L1>

Spec.	Model	"CR-8" Max. Capa. 8 Lit	"CR-12" Max. Capa. 12 Lit	"CR-22" Max. Capa. 22 Lit	"CR-30" Max. Capa. 30 Lit
Temp Range & Accuracy		-25℃+150℃, ±0.1℃			
Temp Resolution		0.1℃-Display, 0.1℃-Control			
Temp Uniformity & Probe		±0.2℃ at -10℃, PT100			
Pump		Max Flow Rate ; 25Lit/min(Controlled 5 Steps Flow Rate Speed) Max Pressure ; 7.5 psi, Max Lift ; 4m			
Bath Volume		8 Lit	12 Lit	22 Lit	30 Lit
Heating & Refrigerator		1 kW, 1/3HP	1.5 kW, 1/3HP	2 kW, 1/2HP	2.2 kW, 3/4HP
Cooling Capacity	at +20℃	300 W	380 W	680 W	850 W
	at 0℃	200 W	250 W	500 W	620 W
	at -20℃	110 W	130 W	300 W	400 W
Refrigerant		CFC -Free(R-404A) Refrigeration System			
Total Power Consumption		1.6 kW	2 kW	2.7 kW	3 kW
Compressor Control		Intelligent Control of Refrigerator for preventing Compressor Overload			
External Circulation Nozzle		id.10.5 / od.12.5mm(For Inner Circulat. ; to Connect directly the in-/outlet Nozzles)			
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end			
Display		Digital LCD with Back-Light Function			
Controller		Digital Fuzzy Controller with Jog-Dial Knob (Turn+Push)			
Safety Circuit		Over-Temp/-Current protection, Sensor Error Detection, Leakage Breaker			
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel Insulation ; 40mm Thickness Urethane Foaming			
Lid & Drain Valve		Includes Stainless-steel Flat Lid & Ball Valve			
Dimension (w×d×h)	Effective Space (mm)	175×153×150	230×181×150	290×220×200	330×290×200
	Internal (mm)	175×320×150	230×350×150	290×390×200	330×460×200
	External (mm)	330×474×752	400×523×807	483×563×850	494×633×906
Packing Size (mm)		420×560×925	530×630×985	615×670×1025	630×740×1075
& Gross Weight (kg)		45 kg	54 kg	63 kg	70 kg
Electrical Rqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz			
Cat. No.	230V	DH.WCR00408	DH.WCR00412	DH.WCR00422	DH.WCR00430
	120V	DH.WCR00508	DH.WCR00512	DH.WCR00522	DH.WCR00530
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 윤임 별도

Refrigerated Circulator with Validation (IQ, OQ) Service


- DH.WCR30308 Validation Service(IQ, OQ), "CRBV208", with "CR-8"
- DH.WCR30312 "CRBV212", with "CR-12"
- DH.WCR30322 "CRBV222", with "CR-22"
- DH.WCR30330 "CRBV230", with "CR-30"

Ea. Price, ₩

(1) Optional Microporous Ceramic Insulation Lid
T15mm(10mm Microporous Ceramic+5mm Ceramic Paper)


- DH.WCR.L008 Insulation Lid, for CR-8, 240×178×T15mm
- DH.WCR.L012 for CR-12, 316×206×T15mm
- DH.WCR.L022 for CR-22, 396×245×T15mm
- DH.WCR.L030 for CR-30, 436×316×T15mm

Ea. Price, ₩



(1) Insulation Lid (Optional)

I. Insulated PTFE Hose Set



SciLab™ -40+240°C Insulated Circulation Hose Set, PTFE & Flexible STS, L1.7 & 1.5m
 Used for Heating/Cooling Circulators for Reactors or Evaporator System, with STS Connect Kit,
 가열 / 냉각 순환기 연결호스, 테프론 & 유연성 스텐레스

■ **Supplied by SciLab Korea**

■ **최상의 내열·내약품성 호스로, Adapter 선택과 길이 조절에 따라 거의 모든 가열/냉각 순환기에 사용 가능.**

■ **Made of PTFE & STS 304 ; Excellent Resistance to High Temp, Oils and Chemicals**

* Stainless-steel Drain Valved "T" Connect Kit for Outlet of Heat Transfer Fluid in Reactor Vessels.

* Can be Used with Almost Any Heating/Cooling Circulator depending on Adapter Selection and Length Adjustment.

I. Insulated PTFE Hose Set, Up to 20Lit Reactors, 단열 실편조 PTFE 호스 세트

► **Connects DIN GL18 or GL32 Screw Male Port of Reactor to OD11~13mm Nozzle Tube of Circulator**

* -200°C +260°C PTFE GL Screw(Female) Connect Kit with O-Ring(PTFE / VITON) Gasket & STS Connect Tube


* "T" Drain Valve Kit & STS PT1/2" or PT3/4" Both Male Screw Adapter

* 1.75mm Thick Stainless-steel Braid and Cotton Insulated * Good Flexibility

<L1>

- | Set. | Price, ₩ |
|---|----------|
| SL.Tub4041 Insulated PTFE Hose Set, GL18(Female), L1.7m , Good for up to 5Lit Reactors
Consists of ; Connector Kit(SL.Con4004) + Adapter(SL.Ada701.6)
+ 0.7m Hose(SL.Tub4043) + T Connector Kit(SL.Con701.5)
+ 1m Hose(SL.Tub4044) + Adapter(SL.Ada701.6) | |
| SL.Tub4042 Insulated PTFE Hose Set, GL18(Female), L1.5m , Good for up to 5Lit Reactors
Consists of ; Connector Kit(SL.Con4004) + Adapter(SL.Ada701.6)
+ 1.5m Hose(SL.Tub4045) + Adapter(SL.Ada701.6) | |
| SL.Tub4051 Insulated PTFE Hose Set, GL32(Female), L1.7m , Good for up to 20Lit Reactors
Consists of ; Connector Kit(SL.Con4005) + 0.7m Hose(SL.Tub4043)
+ T Connector Kit(SL.Con701.5) + 1m Hose(SL.Tub4044)
+ Adapter(SL.Ada701.6) | |
| SL.Tub4052 Insulated PTFE Hose Set, GL32(Female), L1.5m , Good for up to 20Lit Reactors
Consists of ; Connector Kit(SL.Con4005) + 1.5m Hose(SL.Tub4045)
+ Adapter(SL.Ada701.6) | |

II-A. Insulated STS Hose Set (GL32 Port)



II. Insulated Flexible STS Hose Set, Up to 100Lit Reactors, 단열 플렉시블 스텐호스 세트

* Made of Stainless-steel 304, Core Hose ID17×OD26mm.

* 17mm Thick Silicone Foam Insulated, -40°C +240°C (Max.)

(II-A) Stainless-steel Hose Set, for 10 & 20Lit Reactors (GL32 Port)


► **Connects DIN GL32 Screw Male Port of Reactor to OD11~13mm Nozzle Tube of Circulator**

* -200°C +260°C PTFE GL32 Screw(Female) Connect Kit with O-Ring(PTFE / VITON) Gasket & STS Connect Tube

* "T" Drain Valve Kit & STS 1/2" (7/8-UNF) LOK Fitting Adapter for Connection Tube of Circulator

- | Set. | Price, ₩ |
|---|----------|
| SL.Tub7031 Insulated Flexible STS Hose Set, GL32(Female), L1.7m
Consists of ; Connector Kit(SL.Con4005) + 0.7m Hose(SL.Tub7033)
+ T Connector Kit(SL.Con701.5) + 1m Hose(SL.Tub7034)
+ Adapter(SL.Ada701.6) | |
| SL.Tub7032 Insulated Flexible STS Hose Set, GL32(Female), L1.5m
Consists of ; Connector Kit(SL.Con4005) + 1.5m Hose(SL.Tub7035)
+ Adapter(SL.Ada701.6) | |

II-B. Insulated STS Hose Set (Flat Flanged Port)



(II-B) Stainless-steel Hose Set, for 30~100Lit Reactors (Flat Flanged Port)

► **Connects OD 33.5 / 34mm Male Port of Reactor to PT1/2" or PT3/4" Female Port of Circulator**

* "T" Drain Valve Kit & STS PT1/2" or PT3/4" Both Male Screw Adapter

- | Set. | Price, ₩ |
|--|----------|
| SL.Tub7036 Insulated Flexible STS Hose Set, for PT1/2" Female Port, L1.7m
Consists of ; Connector Kit(SL.Con701.1) + 0.7m Hose(SL.Tub7033)
+ T Connector Kit(SL.Con701.5) + 1m Hose(SL.Tub7034)
+ Adapter(SL.Con702.8) | |
| SL.Tub7037 Insulated Flexible STS Hose Set, for PT1/2" Female Port, L1.5m
Consists of ; Connector Kit(SL.Con701.1) + 1.5m Hose(SL.Tub7035)
+ Adapter(SL.Con702.8) | |
| SL.Tub7038 Insulated Flexible STS Hose Set, for PT3/4" Female Port, L1.7m
Consists of ; Connector Kit(SL.Con701.1) + 0.7m Hose(SL.Tub7033)
+ T Connector Kit(SL.Con701.5) + 1m Hose(SL.Tub7034)
+ Adapter(SL.Con702.9) | |
| SL.Tub7039 Insulated Flexible STS Hose Set, for PT3/4" Female Port, L1.5m
Consists of ; Connector Kit(SL.Con701.1) + 1.5m Hose(SL.Tub7035)
+ Adapter(SL.Con702.9) | |

<Continued On Next Page>

<Continued From...SciLab™ -40+240℃ Insulated Circulation Hose Set>



III. Accessories for Circulation Hose Sets

(1) Insulated PTFE Hose Only, 1/2"(7/8-UNF) Both Female

SL.Tub4043 Insulated PTFE Hose Only, 1/2", id12 × od14mm, Overall od17.5mm, L0.7m
 SL.Tub4044 , L1m
 SL.Tub4045 , L1.5m

<L1>

Ea. Price,₩



(2) Insulated Flexible Stainless-steel Hose Only, 1/2"(7/8-UNF) Both Female

SL.Tub7033 Insulated Flexible STS Hose Only, 1/2", id17 × od26mm, Overall od60mm, L0.7m
 SL.Tub7034 , L1m
 SL.Tub7035 , L1.5m

Ea. Price,₩

(3) PTFE GL(Female) Connect Kit

* With PTFE/VITON O-Ring Gasket & PTFE or STS Connect Tube



SL.Con4002 PTFE Connect Kit, GL18, PTFE Tube od12.7mm for id10~11mm Hose
 SL.Con4007 , GL25, PTFE Tube od16mm for id13~16mm Hose
 SL.Con4003 , GL32, PTFE Tube od20mm for id18~19mm Hose
 SL.Con4004 PTFE Connect Kit, GL18, STS Tube od12.7mm for id10~11mm Hose
 SL.Con4005 , GL32, STS 1/2"(7/8-UNF) Male Connector

Kit. Price,₩



(4) STS Flat Flange Connect Kit

SL.Con701.1 STS Connect Kit, od 33.5 / 34, PT Female Clamp & 1/2"(7/8-UNF) Male
 SL.Con701.1B , PT Female Clamp & 1/2"(7/8-UNF) Male, L-type

Kit. Price,₩

(5) STS Drain Valved "T" Connect Kit

* For Outlet of Heat Transfer Fluid in Reactor Vessel



SL.Con701.2 STS "T" Connect Kit, with od12.7mm Tube for (GL18) id10mm Hose
 SL.Con701.3 , with od20mm Tube for (GL32) id18mm Hose
 SL.Con701.4 STS "T" Reducing Connect Kit, with od20mm Tube & od12.7mm Tube
 SL.Con701.5 STS "T" Connect Kit, with 1/2"(7/8-UNF) Male Connector for id1/2"(7/8-UNF) Hose

Kit. Price,₩

(6) STS "LOK Fitting" Reducing Adaptor

SL.Ada701.6 STS "LOK Fitting" Adaptor, 1/2"(7/8-UNF) Male & Female for od12~12.7mm Tubes


Ea. Price,₩



(7) STS PT "Screw to Tube" Connector

SL.Con702.1 STS PT "Screw to Tube" Connector, PT1/2" Male & od12.7mm Tube
 SL.Con702.2 , PT1/2" Male & od20mm Tube
 SL.Con702.3 , PT3/4" Male & od12.7mm Tube
 SL.Con702.4 , PT3/4" Male & od20mm Tube

Ea. Price,₩



(8) Copper "Tube to Tube" Connector

SL.Con702.5 Copper "Tube to Tube" Connector, for od12.7 & od12.7mm Tube
 SL.Con702.6 , for od20 & od12.7mm Tube(Reducing)
 SL.Con702.7 , for od20 & od20mm Tube

Ea. Price,₩



(9) STS PT Screw Both Male Adapter

SL.Con702.8 STS PT Screw Both Male Adapter, 1/2"(7/8-UNF) & PT1/2"
 SL.Con702.9 , 1/2"(7/8-UNF) & PT3/4"

Ea. Price,₩



ShinEtsu® -30°C + 250°C Cool-/Heat-ing Transfer Silicone Fluid, 1 kg/Can

열매 · 냉매 실리콘 오일, 윤활제

- **Excellent Electrical Insulative Properties**
- **Products Available in Various Viscosities**
- * Temperature has little effect on viscosity
- * Superior Thermo-oxidative Stability
- * Outstanding Lubricity

- **Colorless & Transparent**
- **High Resistance to Shear**
- * High Flash Point
- * Low Surface tension
- * Non-corrosive

- **Low Vapor Pressure**
- **Excellent Chemical Stability**
- * Releasability
- * High Compressibility

<L2>

S8.KF96.100 Silicone Fluid, -20°C + 150°C, 1kg/Can

S8.KF54 Phenylated Silicone Fluid, -30°C + 250°C, 1kg/Can

Can.	Price, ₩
1kg	
1kg	



SK™ -40°C + 400°C Best COOL-/HEAT-ing Transfer Fluid, Synthetic Silicone Mixture, 20Lit

For Low/High-Temp Baths, Clearly Harmless Liquid, 냉매/열매 겸용 실리콘계 합성 열매체유

- **Supplied by SciLab Korea**

* Colorless and Transparent Liquid.

- **The Best Circulation Properties / Thermal Stability.**

* Do not use for Silicone Tube

(1) SYN THERM 7008®, -40°C + 400°C, 저/고온용

* Kinematic Viscosity, cSt ⇒ (1) -40°C ; 51.56cSt, (2) 0°C ; 16.08cSt, (3) 100°C ; 3.43cSt, (4) 200°C ; 1.36 cSt <L1>

* Specific Gravity ; 0.9435 at 15°C * Flash Point(PMCC) ; 168°C * Pour Point ; -72°C

SL.Syn7008 Cool/Heat Transfer Fluid, -40+400°C, 20Lit

Pk.	Price, ₩
20Lit	

(2) SUPER THERM FG-20®, -10°C + 300°C, 고온용

* Food Grade Thermal Medium Oil Satisfying FDA, NSF Standards

* Kinematic Viscosity, cSt ⇒ (1) 25°C ; 35.65cSt, (2) 100°C ; 4.23cSt, (3) 200°C ; 1.27 cSt

* Specific Gravity ; 0.8374 at 15°C * Flash Point ; 230°C * Pour Point ; -15°C

SL.Fg0020 Cool/Heat Transfer Fluid, -10+300°C, 20Lit

Pk.	Price, ₩
20Lit	



SHL™ SYN THERM 22L® -5°C + 320°C Heating Transfer Fluid

Synthetic Alkylbenzene Derivative, for Low- & High-Temp Circulation, 알킬벤젠계 합성열매

- **Supplied by SciLab Korea**

* Hazards Identification of MSDS(Material Safety Data Sheet) ;

- H315 Causes Skin Irritation

- H319 Causes Serious Eye Irritation

- H335 May Cause Respiratory Irritation

- H336 May Cause Drowsiness or Dizziness

* Kinematic Viscosity & Density ;

(1) at 0°C ; 22.7cSt, 0.885kg/ℓ (2) at 40°C ; 21.80cSt, 0.859kg/ℓ (3) at 80°C ; 6.09cSt, 0.834kg/ℓ

(4) at 140°C ; 6.09cSt, 0.795kg/ℓ (5) at 200°C ; 1.10cSt, 0.756kg/ℓ (6) at 260°C ; 0.74cSt, 0.718kg/ℓ <L1>

SL.SYN22L SYN THERM 22L® Heating Transfer Fluid, -5°C + 320°C, 20Lit

Pk.	Price, ₩
20Lit	

DAIHAN™ Sample Bag Mixer “BLE-400”, Lab Stomacher, Digital Controller, Low Noise, for 330~1,620ml Bags
Max 400ml Blending Capacity, Adjustable Pedal Distance, 10 Level Speed Adjustment, 백 믹서/무균 시료 분쇄기, 미생물 실험 전처리



3-Year Warranty



Pedal Distance
Adjust Gauge



- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Food Microbiology, Veterinary, Clinical, Pharmaceutical, Molecular Extraction, Toxin Research, Environmental Analyses
- 10 Level Motor Speed Adjusting Function
- Low Noise by Innovative Rail System
- Ergonomic Door Lever ; Easy to Open & Close
- Power Stomaching Mechanism
- * Diversified Choices of Bags ; 330 ~ 1,620ml Bags(Max 400ml Blending Capacity)
- * For Blending and Filtering Solid Samples for Microbiological Analysis.
- * Pedal Distance Adjust Function Ensure Quick and Precise Sterile Blending of Samples
- * Multifunction Digital Control Display for Monitoring Operating Condition, Speed, Time
- * Easy Cleaning Pedals
- * Wide Angle Door Open by Free Stop Torque Hinge
- * Lever-type Ergonomic Door Handle
- * Powerful Blending by 10:1 Gear and Powerful BLDC Motor

Specifications & Ordering Information

<L1>

Spec.	Model	“BLE-400”
Capacity		330ml ~ 1620ml Bag
Blending Capacity		Max. 400ml
Speed		1~10 Level (1~8 strokes/sec)
Pedal Distance Adjustment		0 ~ 45mm
Display & Controller		Digital LCD with Back Light / Digital Controller
Timer, Alarm		Timer ; Continuous & 0~99min 59sec, Alarm ; Error status & Timer-end
Motor, Transmission		100W, 1/10 Ratio
Door		Viewing Window Type
Material		Pedal ; Stainless-steel(STS304), External ; Powder Coated Steel, Window ; Acrylic
Dimensions(w×d×h mm) & Weight		250×474×h264mm, 18kg
Packing Size(w×d×h mm)		380×580×h470mm
Gross Weight		19kg
Power Consumption		160W
Power Supply		AC 100V~220V, 50/60Hz
Cat. No.	220V	DH.Ble0400
Price, Each. ₩		

Sterile Sample Bags for Sample Mixer


<L2>

Model		3001	3002	3003	3004	1523W	1523FW
Volume		Up to 1,620ml				Up to 710ml	
Sterilization		Gamma Sterilization					
Thickness		0.1mm				0.09mm	
Filter	Type	None	Lateral	Lateral	None	None	Lateral
	Porosity	5~150 μ m					
	Material	None	Non-woven	Non-woven	None	None	Non-woven
Bag	Size	190×300mm				150×230mm	
	Material	PE	PE+Nylon	PE+Nylon	PE	PE	PE+Nylon

DH.Bag3001	(1) Sterile Bag, 190×h300mm, 1620ml, Plain type
DH.Bag3002	(2) , Filter type
DH.Bag3003	(3) , Filter & Wire type
DH.Bag3004	(4) , Wire type
3M.06.1523W	(5) , 150×230mm, 710ml, Wire type
3M.06.1523FW	(6) , Filter & Wire type

Pk.	Price, ₩	Cs.	Price, ₩
50ea		10pk	
50ea		10pk	
50ea		10pk	
50ea		5pk	
		500ea	
		250ea	

How to use Sample Bag



(1) Plain type
190×300mm



(2) Filter type
190×300mm



(3) Filter & Wire type
190×300mm



(4) Wire type
190×300mm



(5) Wire type
150×230mm



(6) Filter & Wire type
150×230mm

NEW



Hanil® Multi Function Blender with Grinder, with Powerful RED Motor, 600~3,200ml With Durable Titanium Blade, KTR & KETI Certification, 블렌더

- Supplied by DAIHAN Korea
- Electronic Safety Certification By KTR(Korea Testing & Research Institute), KETI(Korea Electronics Technology Institute)
- Powerful RED(Reinforced, Efficiency, Durability) Motor & Durable Titanium Blade
- Useful Multi Function ; Mixer, Grinding, Mincing ■ Easy to Use and Simple Design

(1) Double Mix Station Hand Blender, 800ml NEW

- * 전용 거치 스테이션으로 정리와 보관이 용이함.
- * Included Chopping Bowl, Double Whisk, Blade Mix Shaft, Blender and Station
- * Poly Carbonate Beaker ; 800ml * Poly Carbonate Mixer Bottle ; 500ml with Cover
- * Poly Carbonate Chop Cup ; 600ml * Mixer / Chopping
- * 5-Step Speed Control Jog-Dial, Max 14,500 RPM From Powerful DC Motor
- * Consumption of Electrical Motor Power 270W
- * Rated Voltage 220V, 60Hz * Weight ; 2.3kg

※ 운임 별도 <B4>

Pk.	Price, ₩
1ea	



H2.HMHD800 Station Hand Blender, 800ml

(2) Small Blender with Grinder, 600ml

- * Glass Mixer Cup ; 600ml, Glass Grinding Cup ; 300ml
- * Mixer / Grinding * Stainless Blade STS304
- * 4-Step Speed Control Jog-Dial
- * Consumption of Electrical Motor Power 300W
- * Rated Voltage 220V, 60Hz * Weight ; 2.2kg

※ 운임 별도

Pk.	Price, ₩
1ea	



H2.SGMF600 Blender with Grinder, Small, 600ml

(3) Middle Volume Blender with Grinder, 1,600ml

- * Glass Mixer Cup ; 1,600ml, Glass Grinding Cup ; 100ml
- * Mixer / Grinding * Safety Design
- * 4-Step Speed Control Jog-Dial
- * Consumption of Electrical Motor Power 320W
- * Rated Voltage 220V, 60Hz * Weight ; 5kg

※ 운임 별도

Pk.	Price, ₩
1ea	

H2.HMF1600PB Blender with Grinder, Middle Volume, 1,600ml

(4) Premium Blender with Grinder, 3,200ml NEW

- * Stainless-steel Mixer Cup with Transparent window ; 3,200ml
- * Mixer / Grinding * Safety Design * Stainless Blade STS304
- * 3-Step Speed Control Button
- * Consumption of Electrical Motor Power 400W
- * Rated Voltage 220V, 60Hz * Weight ; 5.8kg

※ 운임 별도

Pk.	Price, ₩
1ea	

H2.NMF3300 Blender with Grinder, Big Size, Stainless-steel, 3,200ml

NEW



(4)

DAIHAN™ Class II Type A2 Biological Safety Cabinet “BSC”, with 2×ULPA-U15 Filter & UV Lamp, w1200·1500·1800mm
BioSafety Level 1~3, “ECM” DC Motor, Real-Time Air-Flow & Air-Speed Monitoring System, 7 inch Touch TFT LCD Controller
Multi-Safety Function, Compliance with EN12469, 생물안전작업대 “BSC”, 안전등급 BSL 1~3, 안전 저소를 설계, Alarm부




3-Year Warranty

Class II Type A2 Biological Safety Cabinet(BSC),
“BSC-15” Model

Re. Approx Size of Particulate Matter

- PM2.5 초미세먼지 ; 1.0 ~ 2.5 μ m Particles
- PM1.0 극초미세먼지 ; $\leq$ 1.0 μ m Particles
- Fume(Acids) ; $\leq$ 1.0 μ m (이하) Particles
- Bacteria ; 0.5 μ m ~ 0.5m/m Particles
- Virus ; 0.018 ~ 0.4 μ m Particles
- VOCs ; Gas (기체)



Safe Air Circulation System by
EN12469 Standard

필터형 생물안전작업대 ; 별매품인 Steel Frame Stand 와 조합가능 !

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- ▶ Filter Replacement Cycle Alarm Function ; 필터 교체 시기 알람 ALARM 내장 !
- “ECM” DC 모터 사용 시간을 카운트하여 필터 교체 알람(6개월) - 필터교체 주기는 약 6~12개월로
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Compliance with Biological Safety Cabinet European Standard EN12469 ; Exhaust of 30% Filtered Airflow-rate & Recirculation of 70% Filtered Airflow-rate
- Class II Type A2 Protection for the Operator, Environment and Sample at BSL 1~3
- 2ea ULPA-U15(Ultra Low Particulate Air) Filter for Recirculation & Exhaust ;
- $\geq$ 0.12 μ m Particles Removal 99.9995% by High Quality Dense Glass Fiber
- Recommended Replacement Cycle of ULPA-U15 Filter ; Every 6 ~ 12 months
- ▶ Movable Benchtop. Ensuring High Performance Exhaust Filtering without Duct
- ▶ Electronically Commutated “ECM” DC Motor ;
- High Energy Efficient (70% Energy Efficiency Higher than General AC motor)
- Low Noise Level (Less than 60dB) - Ensuring Stable Air-speed
- ▶ Intuitive User Interface of 7” Full Touch TFT LCD Controller ;
- Monitoring Function for Real-time Circulation Airflow & Airspeed
- Maintenance of ULPA-U15 filter Operation Time - Error Status Recording Function.



“ECM” DC Motor

- * Adjustable Height UV Protection Tempered Glass Sliding Sash ;
Safety Height-zone Alarm Function(Alarm sound over 200mm height).
- * Timeless Usage based on Built-in Comfortable Armrest * Easy to Clean Working-table
- * IP68 AC Inlet allows Desktop Equipment to be able to Positioned Inside

Specifications & Ordering information

<L1>

Spec.	Model	“BSC-12” w1200	“BSC-15” w1500	“BSC-18” w1800
Working Space		1250×450×h700mm 360Lit	1550×450×h700mm 450Lit	1850×450×h700mm 540Lit
Air Flow-type		Down Flow, Vertical Recirculation Air Flow-type (Recirculation 70%, Exhaust 30% by EN12469 Standard)		
Air Volume		460 m³/h	580 m³/h	700 m³/h
Air Velocity		Inflow ; 0.53 m/s $\pm$ 0.025, Downflow ; 0.35 m/s $\pm$ 0.08m/s		
Bactericidal UV Lamp		1×Powerful 254nm UV Decontamination, 30W		
Filter and the Replacement Cycle		* 2ea ULPA-U15(Ultra Low Penetration Air) Filter (1×for Exhaust & 1×for Recirculation), $\geq$ 0.12 μ m Particles Removal 99.9995% * Replacement Cycle of ULPA Filter ; 6 ~ 12 months		
Motor		1.25kW Electronically Commutated “ECM” DC Motor		
Window		UV Protection Tempered Glass Sliding Window		
Controller		Full Touch Screen 7” TFT LCD		
Noise Level(dB)		$\leq$ 68 dB		
Internal Lighting		2×15W LED Lamp(800Lux)		
Material		Internal ; Stainless-steel(#304), External ; Powder Coated Steel		
Utility Device		2×AC Outlet(220V, IP68), 1×Vacuum Cock(3/8” NPT)		
Dimension(mm)	Outer	1360×870×h1570	1560×870×h1570	1860×870×h1570
	Without Stand	1250×450×h700	1550×450×h700	1850×450×h700
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz		
Cat. No.	230V	DH.BSC12	DH.BSC15	DH.BSC18
	120V	DH.BSC121	DH.BSC151	DH.BSC181
Price, Each. ₩				

※ 운임 별도

<Continued On Next Page>

<Continued From...DAIHAN™ Class II Type A2 Biological Safety Cabinet "BSC">



Optional Frame Stand



Spare ULPA Filter Kit(2×Filters) ;
1ea for Recirculation & 1ea for Exhaust

Optional Frame Stand for Biological Safety Cabinet, Powder Coated Steel

DH.BSC12S Stand, For Biological Safety Cabinet w1200 Model, 1360×870×h700
DH.BSC15S w1500 Model, 1560×870×h700
DH.BSC18S w1800 Model, 1860×870×h700

Spare ULPA-U15 Filter Kit, ≥0.12μm Particles Removal 99.9995%

DH.BSC12F ULPA-U15 Filter Kit, Included 2×ULPA Filter for w1200 Model,
1× For Recirculation ; 1220×457×h65mm,
1× For Exhaust ; 657×455×h85mm
DH.BSC15F ULPA-U15 Filter Kit, Included 2×ULPA Filter for w1500 Model,
1× For Recirculation ; 1520×457×h65mm,
1× For Exhaust ; 857×455×h85mm
DH.BSC18F ULPA-U15 Filter Kit, Included 2×ULPA Filter for w1800 Model,
1× For Recirculation ; 1820×457×h65mm,
1× For Exhaust ; 1057×455×h85mm

<L1>

Ea. Price, ₩

Ea. Price, ₩



Inspection Certificate of Filters
MAI (주)마디 CERTIFIED



※ Please get in touch with DAIHAN as Filter System Replacement requires Airflow Validation

※ Filter System Replacement Cycle depends on the Environment but 6-12 Month is Recommended

▶ Caution ! Please Replace the Polluted Filters in Good Time ! 오염된 필터들을 적시에 교체하세요 !

<User Convenience of "BSC">



Safety Height-Zone Alarm Function
(Alarm sound over 200mm height).



Comfortable Armrest



7" Full Touch TFT LCD Controller
for User Convenience

- Monitoring Function for Real-time Airflow & Airspeed
- Real-time Error Alarm
- Maintenance of ULPA filter Operation Time

Installed Utility Device



1× Vacuum Cock
(3/8" NPT)



2× 220V AC Outlet
(IP68)



Pressure Gauge
for Downflow
Air Circulation

DAIHAN ComFort™ Sturdily Fabricated Benchtop Clean Bench/PCR Workstation “CBP”, Vertical Flow, ULPA-U15 & Pre-Filter Clear PC 3-Side Viewing, Safely Fabricated Body & Filter, Bactericidal UV Lamp, Touch-type Controller, w700~1,500mm Ideal for Sensitive Biological Applications, LED Lamp, 2 Utility Holes, 8-Step Vertical Air Velocity(0.3~0.6m/s), 5-Years Warranty 컴포트™ 크린벤치 / PCR 워크스테이션, 투명 PC 3면 / 조립식 벤치탑, 수직 기류형, 8단계 풍속 조절, 필터 교체주기 Alarm, 5년 보증-조립기기 안전성



Fabricated Portable Clean Bench 조립식 크린벤치





Fig. 11° Slanted 2-Step Folding Window & Optional Stand


Inspection Certificate of Filters
MDI CERTIFIED

▶ Easy & Quick Replacement of Framed Filters 쉽고 신속한 필터교체 DIY



Remove the Filters in order, then Install the New Filters in Reverse Order. Reassemble the Filter Cover and Close the Front Panel.

▶ Tool-free Fabrication-Easy & Sturdy 별도의 도구 없이 쉽고 튼튼한 조립 DIY



견고한 구조의 조립식 Portable 필터형 크린벤치 ; Optional Steel Frame Stand!

■ Supplied by DAIHAN Korea

■ PL(Product Liability) Insurance

▶ 5-Year Warranty of Safely Sturdy Fabrication and Replaced Filter System 조립기기 안전성 5년 보증.

▶ ULPA-U15 Filter ; $\geq 0.12\mu\text{m}$ Particles Removal 99.9995% by High Quality Dense Glass Fiber ; Recommended Replacement Cycle of ULPA-U15 Filter ; Every 6~12 months

▶ Polyester Pre-Filter($\geq 3\mu\text{m}$ Particle) for Main Filter Protection From Dust

▶ Filter Replacement Cycle Alarm Function ; 필터 교체 시기 알람 ALARM 내장 !

- Filter Replacement Alarm by Counting AC Motor Usage Time

- Filter Replacement Cycle is Approx. 6~12 Months.

■ with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

▶ Powerful AC Motor ; Ensuring Stable Air-speed, Low Noise Level(Less than 65dB)

- Optimized Air Flow System for Uniform Air Flow & Velocity

▶ UV Protection PC(Polycarbonate) 3-Side Viewing Window(8mm Thick) ; 11° Slanted 2-Step Folding

* $2 \times \phi 50\text{mm}$ ON-OFF Hole(PC) for Various Applications * Easy to Clean Working-table

Specifications & Ordering information

<L1>

Spec.	Model	“CBP-07” w700	“CBP-09” w900	“CBP-12” w1200	“CBP-15” w1500
Working Space w×d×h(mm), Lit.		670×570×h635mm, 240Lit	870×570×h635mm, 310Lit	1170×570×h635mm, 420Lit	1470×570×h635mm, 530Lit
Air Flow-type		Vertical Air Flow-type			
Air Volume		360m³/h	480m³/h	650m³/h	820m³/h
Air Flow Velocity(m/s)		0.3m/s ~ 0.6m/s, 8 Levels			
Bactericidal UV Lamp		1×Powerful 254nm UV Decontamination, 30W			
Filter and the Replacement Cycle		* 1 ULPA-U15 Filter(≥0.12μm Particles Removal 99.9995%) * 1 Polyester Pre-Filter(≥3μm Particles) * Replacement Cycle of ULPA & Pre-Filter ; 6~12 months			
Motor & Noise Level(dB)		1×150W AC Motor, and ≤ 65dB		2×150W AC Motor, and ≤ 65dB	
Controller		Digital PID Controller with Touch-button			
Internal Lighting		1ea 17W LED Lamp(600Lux)		1ea 29W LED Lamp(600Lux)	
Clear PC Front Window		11° Slanted 2-Step Folding Transparent Window-Optimal Visibility			
Utility Device Hole		2×φ50mm ON-OFF Hole(PC) ; In-Out of Electrical Cords &c.			
Material		Head & Rear Wall ; Powder-coated Steel			
		3-Side ; UV Protection PC (Polycarbonate) Window - 8mm Thick.			
		Working Plate ; Stainless-steel #304			
Dimension (mm)	Outer	825×645×h1060	1025×645×h1060	1325×645×h1060	1625×645×h1060
	Inner	670×570×h635	870×570×h635	1170×570×h635	1470×570×h635
Net Weight		65kg	75kg	82kg	94kg
Packing & Gross Weight	Head	1160×780×h655, 58kg	1160×780×h655, 68kg	1760×780×h655, 89kg	1760×780×h655, 146kg
	Wall & Plate	1020×925×h200, 31kg	1020×925×h200, 41kg	1620×925×h200, 51kg	1620×925×h200, 61kg
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz			
Cat. No. (w/o Stand)	230V	DH.CBP07	DH.CBP09	DH.CBP12	DH.CBP15
	120V	DH.CBP07.1	DH.CBP09.1	DH.CBP12.1	DH.CBP15.1
Price, Each. ₩					

※ 운임 별도

Optional Frame Stand for PCR Workstation, Powder Coated Steel

DH.CBP07S	Stand, For w700 , 710×625×h700mm, Packing ; 1350×805×h225mm, 20kg
DH.CBP09S	, For w900 , 910×625×h700mm, Packing ; 1350×805×h225mm, 25kg
DH.CBP12S	, For w1200, 1210×625×h700mm, Packing ; 1350×805×h225mm, 35kg
DH.CBP15S	, For w1500, 1510×625×h700mm, Packing ; 1550×805×h225mm, 45kg

Spare ULPA-U15($\geq 0.12\mu\text{m}$ Particles Removal 99.9995%) & Polyester Pre-Filter($\geq 3\mu\text{m}$ Particles) Kit

DH.CBP07.F1	Prefilter for-w700 · 900 · 1200 · 1500, 457×457×h10mm
DH.CBP07.F2	ULPA-U15 Filter for-w700 , 720×457×h70mm
DH.CBP09.F2	-w900 , 920×457×h70mm
DH.CBP12.F2	-w1200, 1220×457×h70mm
DH.CBP15.F2	-w1500, 1520×457×h70mm

Ea. Price, ₩


Ea. Price, ₩

DAIHAN™ Sturdily Steel-type Filtering Vertical Clean Bench/PCR Workstation “CBV”, HEPA-H14 Filter & Pre-Filter, w1,200 & 1,500mm With Bactericidal UV Lamp, Touch-type Controller, 8 Steps Vertical Air Velocity(0.3~0.6m/s), 2 × LED Lamp, Vacuum Cock, and 2 AC Outlets Ideal for Sensitive Biological Applications, 크린벤치 / PCR 워크스테이션, 견고한 스틸 타입, 수직 기류형, 8단계 풍속 조절, 필터 Alarm




Re. Approx Size of Particulate Matter

- PM2.5 초미세먼지 ; 1.0 ~ 2.5 μ m Particles
- PM1.0 극초미세먼지 ; $\leq$ 1.0 μ m Particles
- Fume(Acids) ; $\leq$ 1.0 μ m (이하) Particles
- Bacteria ; 0.5 μ m ~ 0.5m/m Particles
- Virus ; 0.018 ~ 0.4 μ m Particles
- VOCs ; Gas (기체)



► **Caution ! Please Replace the Polluted Filters in Good Time !**
오염된 필터들을 적시에 교체하세요 !



견고한 스틸 타입 필터형 크린벤치 ; 별매품인 Steel Frame Stand 와 조합가능 !

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- Filter Replacement Cycle Alarm Function ; 필터 교체 시기 알람 ALARM 내장 !
 - AC 모터 사용 시간을 카운트하여 필터 교체 알람 (6개월) - 필터교체 주기는 약 6~12 개월로
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- HEPA-H14 Filter ;
 - $\geq$ 0.3 μ m Particles Removal 99.995% by High Quality Dense Glass Fiber
 - Recommended Replacement Cycle of HEPA-H14 Filter ; Every 6~12 months
- Polyester Pre-Filter($\geq$ 3 μ m Particle) for Main Filter Protection From Dust
- Movable Benchtop ; Ensuring High Performance Exhaust Filtering without Duct
- Powerful AC Motor ; Ensuring Stable Air-speed, Low Noise Level(Less than 65dB)
 - Optimized Air Flow System for Uniform Air Flow & Velocity
- * Utility Devices for Various Applications ; 2 × AC Outlet(220V, IP68), 1 × Vacuum Cock(3/8" NPT)
- * Adjustable Height UV Protection Tempered Glass Sliding Sash * Easy to Clean Working-table
- * IP68 AC Inlet allows Desktop Equipment to be able to Positioned Inside

Specifications & Ordering information

<L1>

Spec.	Model	"CBV-12" w1200	"CBV-15" w1500
Working Space		1200×660×h610mm / 470Lit	1500×660×h610mm / 590Lit
Air Flow-type		Vertical Air Flow-type	
Air Volume		650m ³ /h	820m ³ /h
Air Flow Velocity		8 Steps Vertical Air Flow Velocity(0.3~0.6 m/s)	
Bactericidal UV Lamp		1×Powerful 254nm UV Decontamination, 30W	
Filter and the Replacement Cycle		* 1 HEPA-H14 Filter($\geq$ 0.3 μ m Particles Removal 99.995%) * 1 Polyester Pre-Filter($\geq$ 3 μ m Particles) * Replacement Cycle of HEPA & Pre-Filter ; 6~12 months	
Motor		2×150W AC Motor	
Noise Level		$\leq$ 65dB	
Display / Controller		Back Light LCD / Microprocessor Digital Controller with Touch UI	
Internal Lighting		1×29W LED Lamp(600Lux)	
Window		UV Protection Tempered Glass Sliding Window	
Utility Device		2×AC Outlet(220V), 1×Air Cock, 1×Gas Cock	
Material		Internal ; Stainless-steel(# 304), External ; Powder Coated steel	
Dimension(mm)	Outer	1300×720×h1210	1600×720×h1210
	Inner	1200×660×h610	1500×660×h610
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz	
Cat. No.	230V	DH.CBV12	DH.CBV15
	120V	DH.CBV121	DH.CBV151
Price, Each. ₩			

※ 운임 별도

Optional Frame Stand for Clean Bench, Powder Coated Steel

DH.CBV12S Stand, For Clean Bench w1200 Model, 1300×720×h700
DH.CBV15S Stand, For Clean Bench w1500 Model, 1600×720×h700

Ea.	Price, ₩

Spare HEPA-H14($\geq$ 0.3 μ m Particles Removal 99.995%) and 3 μ m Polyester Pre-Filter ; Easy to Replace Filters

DH.CBP07.F1 Prefilter for-w1200 · 1500, 457×457×h10mm
DH.CBV12FH HEPA-H14 Filter for-w1200, 1140×400×h80mm
DH.CBV15FH HEPA-H14 Filter for-w1500, 1440×400×h80mm

Ea.	Price, ₩



DAIHAN™ Powder-coated Steel Ductless Filtering Reagent Storage Cabinet “RSC-PS”, Internal Air Circulation System, 400 · 800 · 1000-Lit.
Ideal for Storage of Acid · Formaldehyde · VOCs, With Hybrid Carbon Filter · Powder-Coated Steel Chamber & Clear PVC Window
스틸형 내산성 필터식 밀폐형 시약장, 에어필터링 순환식, 휘발성 유기화합물 · 산 · 염기성 및 유해시약 보관용, 무덕트




“RSC-800PS”

Two Powerful
Fan Motors
for Internal
Air Circulation

Two Activated
Carbon Filters ;
Acid and VOCs
+Formaldehyde

Spill-Proof Shelf



Optimized Internal Air Circulation System



(1) Spill-Proof Shelf, Stainless-steel



Inspection Certificate of Filters

MDI (주)마디 CERTIFIED



(2), (3) Carbon Filters

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Fully Enclosed Self-Circulating System ;
 - Prevents Air Inside the Cabinet from being Exposed into the Room.
 - Reagent Gases are Filtered through Two Carbon Filters(Acid and VOCs + Formaldehyde) and Purified Air Circulates Inside.
- Robust Powder-coated Steel Chamber is Resistant to most Corrosive Chemicals.
- Filter Replacement Cycle Alarm Function ; 필터 교체 시기 알람 ALARM 내장 !
 - Filter Replacement Alarm by Counting Motor Usage Time
 - Filter Replacement Cycle is Approx. 6~12 Months.
 - Filter can be Easily Replaced by Opening the Front Top Panel.
- * Fast Air Circulation Boost Starts Automatically When The Door Opens for User Safety.
- * Fan Speed can be Set from 1 to 8 Steps to Increase Filter Life and Minimize Energy Consumption.
- * Height Adjustable Spill-Proof-type STS Shelves(h50mm, T1.2mm)
- * Easy to Observe Inside with PVC Window.
- * High Efficiency Low Noise Turbo Fan with Overheat Protection System
- * Provides Key Lock for Security. * Display Internal Air Circulation Speed

Specifications & Ordering information

<L1>

Spec.		Model	“RSC-400PS”	“RSC-800PS”	“RSC-1000PS”
Capacity			390Lit	800Lit	1,020Lit
Door			1 Door	2 Door	2 Door
Material	External	Powder-coated Steel			
	Internal	Powder-coated Stainless-steel			
	Door	Powder-coated Steel, Window ; Clear PVC(T10mm)			
	Shelf	Stainless-steel			
Spill-Proof Shelf	Type	Height Adjustable			
	Size(mm)	540×490×h50mm (1.2T)	475×535×h50mm (1.2T)	625×535×h50mm (1.2T)	
	Included	5	10	10	
	Max Load	35kg / Shelf			
Lamp			1×17W LED	1×29W LED	
Fan Motor			1×70W BLDC Motor	2×70W BLDC Motor	
Air Velocity			0.4m/s ~ 1m/s, 8 Levels		
Noise Level			less than 65dB		
Controller			Digital PID Controller with Touch-button		
Dimension (w×d×h, mm)	Internal	560×520×h1350	1080×550×h1350	1380×550×h1350	
	External	600×600×h1900	1200×600×h1900	1500×600×h1900	
Filter (w×d×h, mm)	Acid	460×410×h40	760×410×h40	1060×410×h40	
	VOCs + Formaldehyde	460×410×h40	760×410×h40	1060×410×h40	
Power Consumption			100W	200W	
Electrical Regts.			1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz		
Cat. No.			DH.RSC400	DH.RSC800	DH.RSC1000
Price, Ea. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

(1) Spare Spill-Proof Shelf, Stainless-steel

DH.RSC.S400 Spill-Proof Shelf, STS, for “RSC-400PS”, 540 × 490 × h50mm(T1.2mm)
DH.RSC.S800 , for “RSC-800PS”, 475 × 535 × h50mm(T1.2mm)
DH.RSC.S1000 , for “RSC-1000PS”, 625 × 535 × h50mm(T1.2mm)

Ea. Price, ₩

(2) Spare Acid Filter, Chemical-coated Activated Carbon

DH.FHF09.F3 Acid Filter, for “RSC-400PS”, 460 × 410 × T40mm
DH.FHF12.F3 , for “RSC-800PS”, 760 × 410 × T40mm
DH.FHF15.F3 , for “RSC-1000PS”, 1060 × 410 × T40mm

Ea. Price, ₩

(3) Spare VOCs + Formaldehyde Filter, Activated Carbon

DH.FHF09.F5 VOCs+Formaldehyde CH₂O Filter, for “RSC-400PS”, 460 × 410 × T40mm
DH.FHF12.F5 , for “RSC-800PS”, 760 × 410 × T40mm
DH.FHF15.F5 , for “RSC-1000PS”, 1060 × 410 × T40mm

Ea. Price, ₩

DAIHAN™ PP Ductless Filtering Reagent Storage Cabinet “RSC”, Internal Air Circulation System, 240 · 470 · 660 · 940-Lit.

Ideal for Storage of Acid · Chemical · VOCs, With Hybrid Composite Filter · All PP Chamber & Clear PVC Window

PP 내산성 필터식 밀폐형 시약장, 에어필터링 순환식, 휘발성 유기화합물 · 산 · 염기성 및 유해시약 보관용, 무다트



I. “RSC240” & II. “RSC470”
1-Room Models

3-Year Warranty



IV. “RSC940”
2-Room Model

3-Year Warranty



- Supplied by DAIHAN Korea
- Fully Enclosed Self-Circulating System ;
 - Prevents Exposure of the Air Inside the Cabinet to the Room.
 - Reagent Gases are Filtered through 4-Layer Hybrid Composite Filter(HEPA, Acid, Alkaline, VOCs) and Circulated Inside.
- Durable PP Chamber is Resistant to most Corrosive Chemicals.
 - * Height Adjustable PP Sliding Drawer Shelf & PP Tray for Spill Storage
 - * Perforated Drawer Shelf Available(Optional)
 - * Easy to Observe Inside with PVC Window.
 - * High Efficiency Low Noise Turbo Fan with Overheat Protection System
 - * Display Temp, Humidity, VOCs Concentration
 - * Alarm Function when VOCs, Temp and Humidity are above the Set Value
 - * Fan Speed can be Set from 1 to 5 Steps to Increase Filter Life and Minimize Energy Consumption.
 - * Provides Key Lock for Security.
 - * Front Opening for Convenient Filter Management
 - * Completed Hazardous Gas Removal Capability & Door Durability Test



Fig1. Circulation Air Flow
of 2-Room Model

Specifications & Ordering information

<L1>

Spec.		Model	I. “RSC240” 1-Room	II. “RSC470” 1-Room	III. “RSC660” 2-Room ; 660Lit		IV. “RSC940” 2-Room ; 940Lit	
					Left (260Lit)	Right (400Lit)	Left (470Lit)	Right (470Lit)
Capacity			240Lit	470Lit	660Lit		940Lit	
Door			1 Door	1 Door	2 Doors		2 Doors	
Material	External	White PP(Polypropylene) 8T, Front Panel ; Gray PVC 10T						
	Internal	White PP(Polypropylene) 8T						
	Door	White PP(Polypropylene) 8T, Window ; Clear PVC 5T						
	Drawer	White PP(Polypropylene) 8T						
Dimension (w×d×h, mm)	Internal	670×470 ×h730	670×470 ×h1,450	300×470 ×h1,450	670×470 ×h1,450	670×470 ×h1,450	670×470 ×h1,450	
	External	700×600 ×h1,000	700×600 ×h2,000	1,000×600×h2,000		1,400×600×h2,000		
Shelf	Type	Sliding Drawer		Fixed Shelves	Sliding Drawer	Sliding Drawer		
	Size(mm)	568×450×h80		300×465 ×h40	568×450 ×h80	568×450×h80		
	Included	2	4	2	4	4	4	
	Max	4	8	2	8	8	8	
	Max Load	30kg / Shelf		35kg / Shelf	30kg / Shelf	30kg / Shelf		
Circulation Fan			Air Flow ; 3.48m²/min, Glass Fiber Reinforced Plastic, Power Consumption ; 21W					
Controller & Display			7" Full Touch Screen TFT LCD					
Filter	Type	4-Layer Hybrid Composite Filter(VOCs, Acid, Base, HEPA) ※ VOCs(Volatile Organic Compounds)						
	Size(mm)	500×300×h100			700×300×h100			
Net Weight			65kg	95kg	135kg		170kg	
Electrical Reqts.			AC230V 50/60Hz, Single Phase, 22W					
Cat. No.			DH.RSC240	DH.RSC470	DH.RSC660		DH.RSC940	
Price, Ea. ₩								

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Drawer Shelf, Tray, and Filter

DH.RSC.D50	(1) PP Drawer Shelf, 568×450×h80mm
DH.RSC.D50R	(2) PP Drawer Rail
DH.RSC.T50	(3) PP Tray for Spill Storage, 560×440×h50mm
DH.RSC.CF50	(4) Composite Hybrid Filter, 500×300×h100mm, for 240Lit & 470Lit
DH.RSC.CF70	(4) Composite Hybrid Filter, 700×300×h100mm, for 660 & 940Lit

Optional Drawer

DH.RSC.D50P	(5) PP Perforated Drawer Shelf, 568×450×h80mm
-------------	---

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ Acids/Corrosives Safety Cabinet "SCC", with Double Wall Fireproof Structure with Insulation Layer, 70~1,600Lit
With Manual Door, Resistant to Fume/Gas, 내부식성 캐비닛, 산 또는 부식성 유해용액 보관용, 페뉴얼타입



3-Year Warranty



2 Doors, 1100 Lit

Spare Steel Shelf

- Supplied by DAIHAN Korea
- Store Acids/Corrosives Chemicals in Safety Cabinet to Protect Operators and Reduce Fire Risk
- Double Wall ; Durable 1.2mm Thickness of Double Robust Sheet Steel Structure Prevent Deformation
- High Quality Acids/Corrosives Resistant Blue Paint; Provides High Durability against Corrosion & Stretch, and Excellent Chemical Resistance.
- 2 Vent Holes $\phi 6\text{cm}$ (2.36in) and Flashback Arrestors on Each Side to allow Fumes to Escape and Prevent Flames and Flashback into the Cabinet
- Special Features ;
 - (a) Fully Welded Construction offer greater Protection in Fire and Explosion Situation.
 - (b) 3-Point Highly Secure Locking System allows for Improved Closure and Secure Storage.
 - (c) Corrosion-Resistant Stainless Steel Piano Hinges provide Greater Durability and Stability, and make it Open and Close Door easily and smoothly.
 - (d) Prevent the Stored Material Leakage by Anti-Leakage Bottom Design
 - (e) Including Galvanized Steel Shelves, which are sloped to Direct Spills to the Back and Bottom of Leakproof sump.
 - (f) Including 4 Adjustable Feet for Leveling on Uneven Ground.

※ Optional Ductless Filtering Vapor (Fume & Gas) Purifier Box can be Installed.

Specifications & Ordering information

<L1>

Spec.	Model	“SCC-7”	“SCC-17”	“SCC-35”	“SCC-60”	“SCC-110”	“SCC-160”
Chamber Capacity, Lit		70 Lit	170 Lit	350 Lit, Tall	600 Lit	1100 Lit	1600 Lit
Door Q'ty		1 Door, Manual Closing			2 Doors, Manual Closing		
Shelf Q'ty (Galvanized Steel)		1 Shelf	2 Shelves		2 Shelves		
Structure & Material		- Double Wall - Fire Proof Powder Coated Steel; 1.2mm Thickness					
Vent		2 Apertures (ϕ60mm) ; Bottom and Opposite Top					
Optional Vapor Purifier Box		- Possible to Install Portable Filtering Vapor Purifier through 2 Apertures ; Filter Exchangeable Type - Clean Air Circulation of Cabinet Chamber					
Dimension (w×d×h)	Outer(cm)	53×46×56	70×46×89	70×46×165	115×46×165	105×86×165	145×86×165
	Inner(cm)	45×38×41	62×38×74	62×38×150	107×38×150	97×78×150	137×78×150
Packing Size(cm) & Gross Wight		58×53×70	75×53×102	75×53×180	120×53×180	120×93×180	150×93×180
		43kg	67kg	102kg	149kg	194kg	245kg
Cat. No.		DH.SCC07	DH.SCC17	DH.SCC35	DH.SCC60	DH.SCC110	DH.SCC160
Price, Each. ₩							

※ 운임 별도

Spare Steel Shelves for Flammable & Acid/Corrosives Safety Cabinets

DH.SCCT440	Steel Shelf, for 70Lit Cabinets, 440×370×h28mm
DH.SCCT600	, for 170 & 350Lit Cabinets, 600×370×h28mm
DH.SCCT1060	, for 600Lit Cabinets, 1060×370×h28mm
DH.SCCT960	, for 1100Lit Cabinets, 960×770×h28mm
DH.SCCT1360	, for 1600Lit Cabinets, 1360×770×h28mm

Ea. Price, ₩

3-Year Warranty



Safety Cabinet

▶ Figure ;
Safety Cabinet + Vapor Purifier Box

Optional Portable Ductless Filtering Vapor (Fume & Gas) Purifier Box "SCF90-F"

유해가스 종합 정화기 박스, 무덕트 필터링형, 별도로 쉽게 설치 DIY

- 설치법 ; 사용 중인(사용할) 약품장 상부에 $\phi 60\text{mm}$ Hole이 있거나 혹은 뚫어서, 본 정화기와 PVC Flexible Hose $\phi 50$ 으로 연결
- 안전 캐비닛 내에 유해가스 발생 시 덕트없이 필터 여과 정화, 청정 공기 실내에 공급.
- Powerful Forced Airflow and Low Noise by High Performance BLDC Motor
- Complete Set with 2 Filters ; (1) Acid Filter,
(2) VOCs Gas(휘발성 유기화합물) Formaldehyde Gas

DH.SCF90F Filtering Gas Purifier Box, 480×323×h340mm, Powder-Coated Steel, Complete Set with 2 Filters

Ea. Price, ₩

Spare Framed Filters For Vapor Purifier Box.

* The Size of Framed Filters ; 300×300×h20mm

DH.SCF90F3	Chemical Coated Activated Carbon Filter, for Fume(Acids)
DH.SCF90F5	Activated Carbon Filter, for Organics(VOCs) & Formaldehyde Gas

Ea. Price, ₩

▶ Filter Details ; See Page 78!

DAIHAN™ Premium Acids/Corrosives Safety Cabinet "SCC", Compliance with EN14727 EU Standard, 124~836Lit.

Powder Coated Steel with Safety Melamine Coated Anti-Corrosive Plywood Interior, Spill-Proof PP Tray, Self Closing Door & Locking Type
Ideal for Storage of Acids & Corrosives and Flammables, 내부식성 캐비닛, 산 또는 부식성 유해용액 보관용, Resistant Fume/Gas.



3-Year Warranty

I. Under Bench-Type, 2 Doors
Self Closing Door



3-Year Warranty

II. Floor-Type, 2 Doors
Self Closing Door



(1) Spill-Proof & Anti-Corrosive PP Tray, h30mm, T5mm. (Included)



(2) Safety Melamine Coated Plywood Shelf, T15.4mm. (Included)

- Supplied by DAIHAN Korea
- 부식성 화학물질을 안전 캐비닛에 보관하여 작업자와 작업실을 보호함.
- Complies with EN14727 EU Standard
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

- * Highly Secure Door Locking System
- * Prevention of Static Electricity by Ground
- * Includes 4 Adjustable Level Feet
- * Easy to move and install Using Forklift

▶ Double Wall Safety Construction

1) Safety Melamine Resin-Coated Anti-Corrosive Plywood Interior(Metal-Free) for Excellent Chemical Resistance

2) Durable 1.2mm Thickness Robust Sheet Steel Structure Prevents Deformation.

▶ Included h30mm Spill-Proof & Anti-Corrosive PP Trays Prevent Accidental Spill(T5mm).

* Trays are Placed on Each Shelf and on the Bottom.

※ Optional Ductless Filtering Vapor (Fume & Gas) Purifier Box can be Installed.

Specifications & Ordering information

<L1>

Spec.	Model	I. Under Bench-Type		II. Floor-Type		
		"SCC-100"	"SCC-200"	"SCC-400"	"SCC-500"	"SCC-800"
Effective Capacity		124 Lit	248 Lit	418 Lit	500 Lit	836 Lit
Dimension (w×d×h)	Outer (mm)	600×535×692	1200×535×692	659×535×1912	1200×535×1222	1200×535×1912
	Inner (mm)	547×440×518	547×440×518 ×2 Inner Chamber	547×440×1737	547×440×1047 ×2 Inner Chamber	547×440×1737 ×2 Inner Chamber
Material	Internal	Safety Melamine Coated Anti-Corrosive Plywood				
	External	Powder-Coated Steel				
Door ; Self Closing		1 ea	2 ea	1 ea	2 ea	2 ea
Shelf (Melamine-Coated Plywood)		1 ea	2 ea	3 ea	4 ea	6 ea
Shelf Capacity		Max. Load ; 50kg				
Tray (PP)		2 ea	4 ea	4 ea	6 ea	8 ea
Vent (Hole Diameter 60mm)		2 ea	3 ea	3 ea	5 ea	5 ea
Weight		90 kg	100 kg	150 kg	200 kg	250 kg
Cat. No.		DH.SCC100	DH.SCC200	DH.SCC400	DH.SCC500	DH.SCC800
Price, Each. ₩						

※ 운임 별도

Spare Shelf & Tray

- DH.SCCT (1) Tray, PP, 540×430×h30mm, T5mm
DH.SCCS (2) Shelf, Melamine(MF) Coated Plywood, 545×435mm, T15.4mm

Ea.	Price, ₩



▶ Figure ;
Safety Cabinet + Vapor Purifier Box

Optional Portable Ductless Filtering Vapor (Fume & Gas) Purifier Box "SCF90-F"
유해가스 종합 정화기 박스, 무덕트 필터링형, 별도로 쉽게 설치 DIY

- 설치법 ; 사용 중인(사용할) 약품장 상부에 $\phi 60\text{mm}$ Hole이 있거나 혹은 뚫어서, 본 정화기와 PVC Flexible Hose $\phi 50$ 으로 연결
- 안전 캐비닛 내에 유해가스 발생 시 덕트없이 필터 여과 정화, 청정 공기 실내에 공급.
- Powerful Forced Airflow and Low Noise by High Performance BLDC Motor
- Complete Set with 2 Filters ; (1) Acid Filter, (2) VOCs Gas(휘발성 유기화합물) Formaldehyde Gas
- 별도로 구매하여 사용자가 직접 Safety Cabinet에 쉽게 연결가능 (※ 자세한 설치 방법은 사용자 매뉴얼 참조)
- Filter Replacement Cycle Alarm Function ; 필터교체 주기는 약 6~12개월로 사용 조건에 따라 다를 수 있음

DH.SCF90F Filtering Gas Purifier Box, 480×323×h340mm, Powder-Coated Steel, Complete Set with 2 Filters

Ea.	Price, ₩

Spare Framed Filters For Vapor Purifier Box.

* The Size of Framed Filters ; 300×300×h20mm

- DH.SCF90F3 Chemical Coated Activated Carbon Filter, for Fume(Acids)
DH.SCF90F5 Activated Carbon Filter, for Organics(VOCs) & Formaldehyde Gas

Ea.	Price, ₩

▶ Filter Details ; See Page 78!

DAIHAN™ Flammable Safety Cabinet “SCF15”, with Double Wall Fireproof Structure with Insulation Layer, 70~1,600Lit
With Manual- / Self Closing-Door, 내화 안전 캐비닛, 유해 인화성 물질(용액) 안전보관용

NEW

3-Year Warranty

“SCF15”
Safe Storage



(1) Single Door,
70 Lit



(1) Single Door,
350 Lit, Tall

- Supplied by DAIHAN Korea ■ 15min Fire-Resistance
- Store Flammable Chemicals in Safety Cabinet to Protect Operators and Reduce Fire Risk
- Double Wall ; Durable 1.2mm Thickness of Double Fire Proof Sheet Steel Structure
- High Quality Heat & Fire Resistant Yellow Paint ; Provides High Durability against Corrosion & Stretch, and Excellent Fire Resistance.
- 2 Vent Holes $\phi 6\text{cm}(2.36\text{in})$ and Flashback Arrestors on Each Side to allow Fumes to Escape and Prevent Flames and Flashback into the Cabinet
- Special Features ;
 - (a) Fully Welded Construction offer greater Protection in Fire and Explosion Situation.
 - (b) 3-Point Highly Secure Locking System allows for Improved Closure and Excellent Storage.
 - (c) Corrosion-Resistant Stainless Steel Piano Hinges provide Greater Durability and Stability, and make it Open and Close Door easily and smoothly.
 - (d) Prevent the Stored Material Leakage by Anti-Leakage Bottom Design
 - (e) Including Galvanized Steel Shelves, which are sloped to Direct Spills to the Back and Bottom of Leakproof sump.
 - (f) Including 4 Adjustable Feet for Leveling on Uneven Ground.

※ Optional Ductless Filtering Vapor (Fume & Gas) Purifier Box can be Installed.

(1) With Manual-Closing Door, 70~1600Lit

<L1>

Model		“SCF15-7”	“SCF15-17”	“SCF15-35”	“SCF15-60”	“SCF15-110”	“SCF15-160”
Spec.							
Chamber Capacity, Lit		70 Lit	170 Lit	350 Lit, Tall	600 Lit	1100 Lit	1600 Lit
Door Q'ty		1 Door, Manual Closing			2 Doors, Manual Closing		
Shelf Q'ty (Galvanized Steel)		1 Shelf	2 Shelves		2 Shelves		
Structure & Material		- Double Wall - Fire Proof Powder Coated Steel; 1.2mm Thickness					
Vent		2 Apertures (ϕ 60mm) ; Bottom and Opposite Top					
Optional Vapor Purifier Box		- Possible to Install Portable Filtering Vapor Purifier through 2 Apertures ; Filter Exchangeable Type - Clean Air Circulation of Cabinet Chamber					
Dimension (w×d×h)	Outer(cm)	53×46×56	70×46×89	70×46×165	115×46×165	105×86×165	145×86×165
	Inner(cm)	45×38×41	62×38×74	62×38×150	107×38×150	97×78×150	137×78×150
Packing Size(cm) & Gross Wight		58×53×70 43kg	75×53×102 67kg	75×53×180 102kg	120×53×180 148kg	120×93×180 193kg	150×93×180 244kg
Cat. No.		DH.SCF1507	DH.SCF1517	DH.SCF1535	DH.SCF1560	DH.SCF15110	DH.SCF15160
Price, Each. ₩							

※ 운임 별도

(2) With Self-Closing Door, 600 · 1100 · 1600 Lit

Spec. Model	“SCF15-60S”	“SCF15-110S”	“SCF15-160S”
Chamber Capacity, Lit	600 Lit	1100 Lit	1600 Lit
Door Q'ty	2 Doors, Self-Closing		
Shelf Q'ty (Galvanized Steel)	2 Shelves		
Structure & Material	- Double Wall - Fire Proof Powder Coated Steel; 1.2mm Thickness		
Vent	2 Apertures ($\phi 60\text{mm}$) ; Bottom and Opposite Top		
Optional Vapor Purifier Box	- Possible to Install Portable Filtering Vapor Purifier through 2 Apertures ; Filter Exchangeable Type - Clean Air Circulation of Cabinet Chamber		
Dimension (w×d×h)	Outer(mm)	1150×460×1650	1050×860×1650
	Inner(mm)	1070×380×1500	970×780×1500
Packing Size(cm) & Gross Weight(kg)		1200×530×1800	1200×930×1800
		153kg	200kg
Cat. No.	DH.SCF1560S	DH.SCF15110S	DH.SCF15160S
Price, Each. ₩			


※ 운임 별도

(3) Spare Steel Shelves for Flammable & Acid/Corrosives Safety Cabinets

DH.SCCT440	Steel Shelf, for 70Lit Cabinets, 440×370×h28mm
DH.SCCT600	, for 170 & 350Lit Cabinets, 600×370×h28mm
DH.SCCT1060	, for 600Lit Cabinets, 1060×370×h28mm
DH.SCCT960	, for 1100Lit Cabinets, 960×770×h28mm
DH.SCCT1360	, for 1600Lit Cabinets, 1360×770×h28mm

Ea. Price, ₩

3-Year Warranty



(2) 2 Doors, 1100 Lit



(3) Spare Steel Shelf

※ Exhaust Filtering System for
Safety Cabinets
⇒ See Page 78!

DAIHAN™ Premium 90min-Resistant Flammables Safety Cabinet "SCF90", Type 90 by EN14470-1 Standard, 103~720 Lit
With 5-Layer Fireproof Structure(Powder-coated Steel · Ceramic-Fiber Hybrid · Calcium Silicate Microporous · Gypsum · Melamine-Coated Plywood)
Automatic Self Closing Door, Patented Safe Drawer-type Shelf, 90분 내화 안전 캐비닛, 유해 인화성 물질(용액) 안전 보관용

3-Year Warranty



I. Under Bench-Type, 2 Doors
Self Closing Door

3-Year Warranty



II. Floor-Type, 2 Doors
Self Closing Door

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 인화성 화학물질을 안전 캐비닛에 보관하여 작업자를 보호하고 화재 및 폭발 위험을 줄임
- Type 90(90min.) Compliance with Safety Storage Cabinets for Flammable EU Standard EN14470-1 ; The Internal Temperature of Safety Cabinet keep under 180K(Temperature Variation 180℃) for up to 90min.
- EN14470-1 시험 기준으로 한국가스안전공사의 인화성 위험물 보관함 성능평가 시험을 통한 안전성 확보.
- 5-Layer Safety Structure ; Complete Fire Resistance & Prevention of Deformation / Explosion. 5겹의 완벽한 내화 및 항폭발력의 구조임.
 - ① Powder-Coated Steel · ② Hybrid Ceramic-Fiber · ③ Microporous Calcium Silicate · ④ Gypsum · and ⑤ Melamine-Coated Plywood
- 자동잠금 장치 고장시 수동잠금이 가능함
If the Self Closing Device doesn't work, it can be also Used As a Manual Closing type.
- h30mm Spill-Proof Roller Drawer-Type Shelves(T1.2mm)
 - * Highly Secure Door Locking System
 - * Prevention of Static Electricity by Ground
 - * Prevents the Stored Material Leakage by Anti-leakage Bottom Design
 - * Easy to move and install Bottom Design Using Forklift
 - * Installed 4 Adjustable Level Feet
- ▶ 캐비닛의 도어와 환풍구가 열려 있는 상태에서 화재 발생 시 화염의 유입을 방지하기 위해서 캐비닛 주변온도가 +70℃ ±10℃가 되면 안전을 위해 자동으로 선반이 안으로 들어가고 도어와 환풍구가 닫힘.
Even if the Doors and Vents are Open and the Shelves are Out, When the Ambient Temp reaches +70℃ ±10℃, the Shelves will Automatically Retract and the Doors and Vents will Close to Prevent Flames from Entering.
- ▶ 특허 등록된 안전 편리형 서랍식 선반 장착. Patented Safe Drawer-type Shelf(특허 ; 제 10-2476972호) ; Auto-return Mechanism Roller Sliding Drawer-Shelf
- ※ Optional Ductless Filtering Vapor (Fume & Gas) Purifier Box can be Installed.


Specifications & Ordering information

<L1>

Spec.	Model	I. Under Bench-Type		II. Floor-Type		
		"SCF90-100"	"SCF90-200"	"SCF90-300"	"SCF90-400"	"SCF90-700"
Type of Fire Resistance		Type 90; 90min Fire Resistance by EN14470-1 EU Standard				
Effective Cabinet Capacity		103 Lit	221 Lit	337 Lit	405 Lit	720 Lit
Dimension (w×d×h)	Outer (mm)	592×600×690	1102×600×690	592×600×1912	1102×600×1222	1102×600×1912
	Inner (mm)	455×470×485	970×470×485	455×470×1580	970×470×890	970×470×1580
Material	Internal	Melamine Resin-Coated Plywood				
	Insulation	Ceramic-Fiber Hybrid Board, Calcium Silicate Microporous Board, Gypsum Board				
	External	Powder-Coated Steel				
Door Types & Q'ty		Self Closing type Door				
		1 ea	2 ea	1 ea	2 ea	2 ea
Shelf, Patented Ease of Use		h30mm Spill-Proof Roller Sliding Drawer-Type(T1.2mm) ; Powder Coated Steel				
Shelf Load & Q'ty		Max. Load 80kg				
		-	-	3 ea	2 ea	2 ea
Vent, φ60mm Hole		2ea Vent Holes on the Backside		2ea Vent Holes on the Upperside		
Weight		140kg	220kg	250kg	330kg	430kg
Cat. No.		DH.SCF9100	DH.SCF9200	DH.SCF9300	DH.SCF9400	DH.SCF9700
Price, Each. ₩						

※ 운임 별도

※ 5-Layer Fireproof Structure For Complete 90min Fire Resistance & Explosion Resistance.



- (1) Powder-Coated Steel ; 파우더 코팅된 스틸을 사용하여 내화특성과 열 변형을 최소화함.
- (2) Ceramic-Fiber Hybrid Board ; 100% 무기질 성분으로 연소 시 연기와 악취가 발생하지 않고 내화 단열성능이 우수함.
- (3) Calcium Silicate Microporous Board ; 섬유와 필러로 보강된 불연성 매트릭스 광물로서 습기에 강하며 폼알데하이드가 포함되어 있지 않다.
- (4) Gypsum Board ; 불연특성의 석고에 무기질 섬유를 보강한 내화성 단열재.
- (5) Melamine-Coated Plywood ; 단단하고 열에 강한 멜라민 수지로 코팅된 합판.

<Continued On Next Page>

<Continued From...DAIHAN™ Premium 90min-Resistant Flammables Safety Cabinet "SCF90">



Patented Safe Drawer-Type Shelves ; Auto-Return Mechanism

Spare Safe Drawer Shelves, Roller Sliding Form.

DH.SCF901D1 Shelf, Powder-Coated Steel, For 300 Model, 450×460×h30mm(T1.2mm)
DH.SCF902D2 , For 400 & 700 Model, 965×460×h30mm(T1.2mm)

<L1>

Ea. Price, ₩

※ 한국가스안전공사 인화성 위험물 보관함 성능평가 완료!!!



EN14470-1 Standard 기준으로 화재 발생 상황을 구현한 한국가스안전공사의 방화시험을 통해 90분 동안 Safety Cabinet 내부 온도 변화 폭이 180℃ 이하로 유지하는 성능평가를 완료함.



3-Year Warranty

Optional Ductless Filtering Vapor Purifier Box

Optional Portable Ductless Filtering Vapor (Fume & Gas) Purifier Box "SCF90-F"
유해가스 종합 정화기 박스, 무덕트 필터링형, 별도로 쉽게 설치 DIY

- 설치법 ; 사용 중인(사용할) 약품장 상부에 $\phi 60\text{mm}$ Hole이 있거나 혹은 뚫어서, 본 정화기와 PVC Flexible Hose $\phi 50$ 으로 연결
- 안전 캐비닛 내에 유해가스 발생 시 덕트없이 필터 여과 정화, 청정 공기 실내에 공급.
- Powerful Forced Airflow and Low Noise by High Performance BLDC Motor
- Complete Set with 2 Filters ; (1) Acid Filter, (2) VOCs Gas(휘발성 유기화합물) Formaldehyde Gas
- 별도로 구매하여 사용자가 직접 Safety Cabinet에 쉽게 연결가능 (※자세한 설치 방법은 사용자 매뉴얼 참조)
- Filter Replacement Cycle Alarm Function ; 필터교체 주기는 약 6~12개월로 사용 조건에 따라 다를 수 있음

Ea. Price, ₩

DH.SCF90F Filtering Gas Purifier Box, 480×323×h340mm, Powder-Coated Steel, Complete Set with 2 Filters

Spare Framed Filters For Vapor Purifier Box.

* The Size of Framed Filters ; 300×300×h20mm

Ea. Price, ₩

DH.SCF90F3 Chemical Coated Activated Carbon Filter, for Fume(Acids)
DH.SCF90F5 Activated Carbon Filter, for Organics(VOCs) & Formaldehyde Gas

※ Filters For Vapor Purifier Box.

Filter	Target	Material	Method	Description
(1) Acids Filter	Fume Vapor	Chemical-Coated Activated Carbon	Chemical	특정한 화학물질이 코팅된 활성탄에 Fume/Vapor(Acids, NH_3 , S, etc.)나 먼지 등을 흡착 여과 정화, 2nm Molecule
(2) VOCs Filter Formaldehyde	Gas	Activated Carbon	Chemical Physical	활성탄에 VOCs(휘발성 유기화합물의 총칭) CH_2O Gas나 먼지 등을 흡착 여과 정화. 2nm Molecule

► Quick & Intuitive Fabrication of Vapor Purifier Box



1. 캐비닛 배출구 홀과 정화기 흡입구 홀 위치 정렬



2. 정화기 흡입구 덮개 제거



3. 캐비닛 배출구 덮개 제거



4. 고정 홀에 맞추어 덕트 호스 조립



5. 덕트 호스 체결 상태 확인



6. 캐비닛과 정화기 연결완료 모습



Framed Filters



Inspection Certificate of Filters
MAD CERTIFIED



Fig. Vapor Purifier Box Installed on the Safety Cabinet
안전캐비닛 위에 설치된 정화기 박스

ABS Assembling Utility Carts, Heavy-Duty, with Tray Shelf · Detachable Handle · 4 Casters

Excellent Resistance Against High Impact · Chemicals · Temp(-40+90℃), ABS 조립식 다용도 카트, 튼튼한 구조의 플라스틱

- Supplied by DAIHAN Korea ■ Sturdy Construction & Best Durability. ■ Easy to Assemble(DIY); Everybody, Anywhere! 누구나 쉽게 조립!
- Made of ABS Resin "Acrylonitrile Butadiene Styrene Copolymer" 우수한 내충격·내구성·내약품성·내열성의 공중합체 플라스틱

- ▶ Multiuse in Foodstuff, Laboratory, Medical, and Industry Field.
- ▶ Compact Cardboard Packaging For Easy Transfer; Average-67×46×21cm & 10kg

(1) Basic-type
; 2 or 3 Shelf(2) Guardrail-type
; 2 or 3 Shelf(3) 1 or 2 Drawer & Guardrail-type
; 2 or 3 ShelfFigure;
Assembling Components
조립용 부품들Figure; Compact Packaging
조립식카트 1개 간편 포장 상태**(1) Basic-type ABS Cart, with 2 & 3 Shelf**

DH.Car30562	Cart, ABS, 2 Shelf, 56×w39×h82cm, Shelf Size - 44×34cm, Shelf Distance - 56cm, Pack - 62×39×14cm, 6.5kg
DH.Car30662	, 66×w46×h82cm, - 53×40cm, - 56cm, - 67×46×15cm, 8.3kg
DH.Car30822	, 82×w49×h82cm, - 68×44cm, - 56cm, - 83×52×15cm, 10.5kg
DH.Car30563	Cart, ABS, 3 Shelf, 56×w39×h82cm, Shelf Size - 44×34cm, Shelf Distance - 27cm, - 58×39×19cm, 7.9kg
DH.Car30663	, 66×w46×h82cm, - 53×40cm, - 27cm, - 67×46×21cm, 10.6kg
DH.Car30823	, 82×w49×h82cm, - 68×44cm, - 27cm, - 83×52×20cm, 13.5kg

※ 윤임 별도 <L1>

Ea.	Price, ₩

(2) Guardrail-type ABS Cart, with 2 & 3 Shelf

DH.Car30562G	Cart, ABS, 2 Shelf, 56×w39×h82cm, Shelf Size - 44×34cm, Shelf Distance - 56cm, Pack - 62×39×14cm, 6.7kg
DH.Car30662G	, 66×w46×h82cm, - 53×40cm, - 56cm, - 67×46×15cm, 8.6kg
DH.Car30563G	Cart, ABS, 3 Shelf, 56×w39×h82cm, Shelf Size - 44×34cm, Shelf Distance - 27cm, - 58×39×19cm, 8.3kg
DH.Car30663G	, 66×w46×h82cm, - 53×40cm, - 27cm, - 67×46×21cm, 11kg
DH.Car30823G	, 82×w49×h82cm, - 68×44cm, - 27cm, - 83×52×20cm, 14kg

Ea. Price, ₩

Ea.	Price, ₩

(3) Drawer & Guardrail-type ABS Cart, with 2 & 3 Shelf, 1- or 2-Drawer

DH.Car30562GD	Cart, ABS, 2 Shelf, 1 Drawer, 56×w39×h82cm, Shelf Size - 44×34cm, Shelf Distance - 45cm, Pack - 62×39×21cm, 7.5kg
DH.Car30662GD	, 2 Drawer, 66×w46×h82cm, - 53×40cm, - 47cm, - 67×46×22cm, 9.9kg
DH.Car30563GD	Cart, ABS, 3 Shelf, 1 Drawer, 56×w39×h83cm, Shelf Size - 44×34cm, Shelf Distance - 27cm, - 58×39×26cm, 8.9kg
DH.Car30663GD	, 2 Drawer, 66×w46×h83cm, - 53×40cm, - 27cm, - 67×46×27cm, 12kg
DH.Car30823GD	, 82×w49×h83cm, - 68×44cm, - 27cm, - 83×52×30cm, 15kg

Ea. Price, ₩


Ea.	Price, ₩

Spare Handles

DH.30562.HD	Spare Handle, Thick 7cm, w36×9cm, For - w39 (width) Carts All.
DH.30662.HD	, w42×9cm, - w46 (width) Carts All.
DH.30822.HD	, w47×9cm, - w49 (width) Carts All.

<L3>

Ea.	Price, ₩



(1) 4-Drawer (2) 5-Drawer

Two-way Drawer Dressing Cart, ABS, with 4 & 5 Drawer, with Handle

Ideal for Lab · Medical · Industrial, with Stop-On Casters, 양방향 서랍식 플라스틱 카트

- Supplied by DAIHAN Korea * Accessory Drawer for Storage / Handling Convenience.
- * Multi-Purpose ; Laboratories, Hospitals, and Work Site for Transport and Storage

(1) 4-Drawer Dressing Cart

- * 2 Drawer(39×43×h8cm) & 2 Drawer(39×43×h17cm)

DH.Car3064D4	Cart, ABS, 4 Drawers, 64×48×h85cm, Shelf Size(in) - 44×43cm, Drawer Size(in) - 39×43×8/17cm
--------------	---

※ 윤임 별도 <L3>

Ea.	Price, ₩


(2) 5-Drawer Dressing Cart

- * 4 Drawer(39×43×h8cm) & 1 Drawer(39×43×h17cm)

DH.Car3064D5	Cart, ABS, 5 Drawers, 64×48×h85cm, Shelf Size(in) - 44×43cm, Drawer Size(in) - 39×43×8/17cm
--------------	---

Ea. Price, ₩

Ea.	Price, ₩



Drawer Dressing Cart, ABS, with 3~6 Drawer, with Guardrails Ideal for Lab · Medical · Industrial, Stop-On Casters, 서랍식 플라스틱 카트

■ Supplied by DAIHAN Korea

* Multi-Purpose ; Laboratories, Hospitals, and Work Site for Transport and Storage

(1) 1~3 Drawer Dressing Cart

* Shelf Size(in) - 55×43cm, Drawer Size(in) - 51×36×16cm

※ 운임 별도 <L3>

Ea. Price, ₩

DH.Car3075D1	① Cart, ABS, 1 Drawer, 74×48×h89cm
DH.Car3075D2	② , 2 Drawers, 74×48×h89cm
DH.Car3075D3	③ , 3 Drawers, 74×48×h89cm

(2) 4 Drawer Dressing Cart

* 2 Drawer(51×36×h8cm) & 2 Drawer(51×36×h16cm)

Ea. Price, ₩

DH.Car3075D4	Cart, ABS, 4 Drawers, 74×48×h89cm, Shelf Size(in) - 55×43cm, Drawer Size(in) - 51×36×8/16cm
--------------	--

(3) 5 & 6 Drawer Dressing Cart, with Detachable Hanger

* Shelf Size(in)-55×43cm, Drawer Size(in)-51×36×8/16cm

Ea. Price, ₩

DH.Car3075D5	④ Cart, ABS, 5 Drawers, 74×48×h89cm
DH.Car3075D6	⑤ , 6 Drawers, 74×48×h89cm

Optional Detachable Hanger **NEW**

DH.Car3075A Detachable Hanger, 410×540mm

Ea. Price, ₩



Stainless-steel Assembling Dressing Cart, with 2 Shelf · Guardrails For Lab · Medical · Industrial, with Stop-On Caster, 조립식 2단 카트

■ Supplied by DAIHAN Korea

■ Easy to Assemble and Disassemble, Simple to Keep after Disassembling

* Easy to Clean and Long Life-time High Quality Stainless-steel ; Composed Stainless-steel Frame/Plate Shelves and can be Fixed 4ea Casters

※ 운임 별도 <L3>

Ea. Price, ₩

C9.4212	Cart, 2 Shelf, Fabricated, 72×45×h80cm
C9.4213	, 95×45×h80cm



DAIHAN mediclin™ Stainless-steel Dressing Cart, with 2 & 3 Shelf and Guardrails For Lab · Medical · Industrial, with Stop-On Caster, 다용도 2 & 3단 카트

■ Supplied by DAIHAN Korea

■ Guard Rails help Prevent Items from Slipping off Shelves.

※ 운임 별도 <L3>

Ea. Price, ₩

(1) With Guardrails 2 Shelf

DM.Car7001	Cart, 2 Shelf, 55×45×h85cm
DM.Car7002	, 65×48×h88cm
DM.Car7003	, 78×48×h92cm
DM.Car7004	, 98×48×h92cm

(2) With Guardrails 3 Shelf

DM.Car7006	Cart, 3 Shelf, 78×48×h92cm
DM.Car7007	, 98×48×h92cm

Ea. Price, ₩



DAIHAN mediclin™ Stainless-steel Handled Dressing Cart, with 2 & 3 Shelf & Guardrails For Lab · Medical · Industrial, 핸들 타입 다용도 2 & 3단 카트

■ Supplied by DAIHAN Korea

* with Stop-On Casters * Convenient to Use Built in Handle

* Easy to Clean and Long Life-time High Quality Stainless-steel ;
Composed Stainless-steel Frame/Plate Shelves.

※ 운임 별도 <L3>

Ea. Price, ₩

(1) 2 Shelf · Vertical Handle-type

DM.Car7011	Cart, 2 Shelf, 78×48×h95cm
DM.Car7012	, 98×48×h95cm

(2) 3 Shelf · Vertical Handle-type

DM.Car7015	Cart, 3 Shelf, 78×48×h95cm
DM.Car7016	, 98×48×h95cm

Ea. Price, ₩

Articles

Cat. No

Description



(1) Horizontal Handle-type 2 Shelf (2) Semi-Vertical Handle-type, 3 Shelf

Stainless-steel Handled Dressing Cart, with 2 & 3 Shelf & Handles

For Lab · Medical · Industrial, 핸들 타입 2 & 3단 카트

■ Supplied by DAIHAN Korea

* Multipurpose Cart ; Laboratory, Hospital, Industrial etc.

* with Stop-On Casters

* Composed High Quality Stainless-steel ; Easy to Clean, Long Life-time, Chemical Resistance

※ 윤임 별도 <L3>

(1) 2 Shelf · Horizontal Handle-type

C9.4223 Cart, 2 Shelf, 65×45×h90cm

C9.4224 , 75×45×h90cm

C9.4225 , 95×45×h90cm

Ea. Price, ₩

(2) 3 Shelf · Semi-Vertical Handle-type

C9.4226 Cart, 3 Shelf, 65×45×h90cm

C9.4227 , 75×45×h90cm

C9.4228 , 95×45×h90cm

Ea. Price, ₩



(1) With 2 Tray Shelf (2) With 2 Plate Shelf, Guardrails & 2 Handle

Utility Stainless-steel Tray Shelved Cart, with Stop-On Caster

For Lab · Medical · Industrial, 트레이 선반식 다용도 카트

■ Supplied by DAIHAN Korea

* 3cm Shelf Assistance Wall ; Not Drop Equipment on Cart when Move

※ 윤임 별도 <L3>

(1) 2 Tray Shelf-type, None Handle

C9.4021 Cart, 2 Shelf, 65×45×h90cm

C9.4022 , 75×45×h90cm

C9.4023 , 85×45×h90cm

Ea. Price, ₩

(2) With 2 Plate Shelf · Guardrails · Both Handle

C9.4032 Cart, 4 Side, 2 Shelf, 75×45×h95cm

C9.4033 , 85×45×h90cm

Ea. Price, ₩

(3) With 3 Tray Shelf · Both Handle

C9.4027 Cart, 3 Shelf, 75×45×h90cm

C9.4028 , 85×45×h90cm

Ea. Price, ₩

(4) With 3 Plate Shelf · Guardrails · Both Handle

C9.4035 Cart, 4 Side, 3 Shelf, 75×45×h95cm

C9.4036 , 85×45×h95cm

C9.4037 , 95×45×h95cm

Ea. Price, ₩



(1) Compact 2 Shelf (2) Heavy Duty 3 Shelf

SciLab™ Heavy-Duty Stainless-steel Cart, 2 & 3 Shelf & Handle

For Lab · Medical · Industrial, with Stop-On Casters, 2 & 3단 카트

■ Supplied by SciLab Korea

■ With Plate Shelf · Guardrail · Handle

* Consist of Stainless-steel(STS) Pipe and Plate, and Excellent Resistance to Chemicals.

* Integrated Push Handle for Easy Maneuverability Handling

* Stainless-steel Plate at the Bottom for Extra Storage

※ 윤임 별도 <L1>

(1) Compact, 2 Plate Shelf · Guardrail · 1-Handle

SL.Car7021 Cart, 2 Shelf, 1 Handle, 54×44×h90cm

Ea. Price, ₩

(2) 3 Plate Shelf · Guardrail · 1-Handle

SL.Car7022 Cart, 3 Shelf, 1 Handle, 75×45×h95cm

SL.Car7023 , 90×45×h95cm

Ea. Price, ₩



(1) with 1 Basket & 1 Shelf (2) With 2 Basket

DAIHAN medicin™ Utility Stainless-steel Cart, with 1 or 2 Wire Basket

Ideal for Drying, Storage, and Transfer, 건조 · 보관 · 운반용 다용도 카트

■ Supplied by DAIHAN Korea

* with Stop-On Casters

※ 윤임 별도 <L3>

(1) With 1 Basket & 1 Shelf

DM.Car7021 Cart, Multi-purpose, 68×48×h95cm

Ea. Price, ₩

(2) With 2 Basket

DM.Car7022 Cart, Multi-purpose, 68×48×h95cm

Ea. Price, ₩



(1) With 2-Wire Basket & 2-Vertical Handle

(2) With 3-Wire Shelf & 1-Handle

Stainless-steel Cart, with Wire-Basket & Wire-Shelf With Stop-On Caster, 와이어 바스켓 & 선반식 카트

■ Supplied by SciLab Korea

- * Clean, No Pollution, Long Life All Stainless-steel(STS) Products ; Stainless-steel(STS) Frame, STS Plate Shelves and Position Consists of the Four Wheels as possible.
- * Purposes; (1) Laboratories, Hospitals, and Work Sites for Transport and Storage.
(2) Laboratory Instruments and Accessories for Dry Cleaning.

(1) With 2-Wire Basket & 2-Vertical Handle

SL.Car7041 Cart, 2×Shelves, 75×40×h75cm

(2) With 3-Wire Shelf & 1-Handle

SL.Car7042 Cart, 3×Shelves, 70×50×h100cm

SL.Car7043 , 90×50×h100cm

※ 운임 별도 <L1>

Ea. Price, ₩

Ea. Price, ₩



(1) With 1-Drawer

(2) With 2-Drawer

SciLab™ Drawer Cart, Stainless-steel, with 1 & 2 Drawer With Stop-On Casters, 서랍식 카트

■ Supplied by SciLab Korea

- * Long Life All Stainless-steel(STS) Frame and Plate Shelves
- * Versatile ; Laboratories, Hospitals (Dressing car), Work site, etc. for the Transport and Storage.
- * Guard Rails help Prevent Items from Slipping off Shelves.
- * Accessory Drawer for Storage / Handling Convenience.

(1) 1 Drawer & 2 Shelf

SL.Car7031 Cart, 1 Drawer, 2 Shelf, 65×45×h90cm

SL.Car7032 , 75×45×h90cm

(2) 2 Drawer & 2 Shelf

SL.Car7033 Cart, 2 Drawer, 2 Shelf, 75×45×h90cm

SL.Car7034 , 90×45×h90cm

※ 운임 별도 <L1>

Ea. Price, ₩

Ea. Price, ₩



(1)



(2) 발페달로 작동되는 뚜껑

Laundry Bag Cart, with Stainless-steel Frame, 162 & 234Lit Good for Hospital, 세탁물 운반카트, 병원용에 적합

■ Supplied by DAIHAN Korea

(1) Without Lid

C9.4411 Cart, Laundry, 60×60×h85cm, 234Lit

(2) With Lid Operated Foot Pedal (ON-OFF)

C9.4415 Cart, Laundry, 50×50×h85cm, 162Lit

※ 운임 별도 <L3>

Ea. Price, ₩

Ea. Price, ₩



DAIHAN mediclin™ Bottle Safety Rack Cart, Stainless-steel, Up to 2×25Lit Round Bottles Used in Lab · Medical · Industrial, with 2 Shelf & Stop-On Casters, 바틀 랙 카트, 25Lit 병 2개 수용

■ Supplied by DAIHAN Korea


- * Integrated Push Handle for Easy Maneuverability Handling
- * Stainless-steel Plate at the Bottom for Extra Storage

DM.Car7031 Safety Rack Cart, for 1×25ℓ Bottle, 37×34×h100cm

DM.Car7032 , for 2×25ℓ Bottle, 72×34×h100cm

※ 운임 별도 <L1>

Ea. Price, ₩

Articles	Cat. No	Description				
		Bottle Safety Rack Cart, Stainless-steel, for Square Bottles 5 & 10 Lit, 3 or 2 Places With “Stop-On” Casters, 바틀 랙 카트, 4각 바틀용에 적합				
		■ Supplied by DAIHAN Korea * Strong & Safety * Bottle is Excluded * Total Height ; 74cm, Under Bottle Bottom ; 60cm				
		※ 운임 별도 <L3>				
		<table><tr><th>Ea.</th><th>Price, ₩</th></tr><tr><td></td><td></td></tr></table>	Ea.	Price, ₩		
Ea.	Price, ₩					
		C9.7431 (1) Safety Rack Cart, For-3×5 Lit Bottle , 37×33×h60/74cm C9.7432 (2) -2×10 Lit Bottle, 47×33×h60/74cm				



SciLab™ Stainless-steel Drying Cart, Adjustable 24 Places

With 3 Shelves and Removable Pegs, Room-To-Rooms, 이동식 초자 건조대, 3단 카트와 건조대

■ Supplied by SciLab Korea						
* Structure ; (1) Cart at the Top / Middle of the Stainless-steel(STS) Wire Shelves, at the Bottom of the STS Tray. (2) Drying the Top of the Stainless-steel(STS) Frame, Modular / Removable in Position Several Heat / Plastic Pack Chemical Resistance and is Composed of Various Types can be Mounted to all Glassware.						
* Multi-purpose 24-Peg, The Peg is Easy to Move Between, and Ease of Adjustment of Spacing and Quantity.						
* Plastic(PP) 24-Peg ; ϕ 10/20 $\times$ 190mm		※ 운임 별도 <L1>				
* Cart on Wheels way "Room-To-Room" is Free to Move.						
SL.Car7051	Cart, Drying, STS, with 3 Shelves & 24 Pegs, 72 $\times$ 50 $\times$ h140cm					
SL.Rac7015	Spare Pegs, for Drying Rack/Cart, PP, ϕ 10/20 $\times$ 190mm, 24ea/pk.					
		<table><tr><th>Pk.</th><th>Price, ₩</th></tr><tr><td>24ea</td><td></td></tr></table>	Pk.	Price, ₩	24ea	
Pk.	Price, ₩					
24ea						



SciLab™ Assembling Multiuse 5-Deck Storage Cart, with Protection Curtain & Stop-On Casters

Ideal for Drying and Storage, with Stainless-steel Frame & Wire- Shelves, 조립식 다용도 5단 건조/보관대

■ Supplied by SciLab Korea ■ Adjustable Shelves Distance ; 선반 간격 조절형 ■ Multi-use ; Movable Shelves and Drying Glassware Organizations such as the Test Apparatus Made of Stainless-steel Wire and Stainless-steel Pipe * Different Tray Hole Sizes (90-/55-/43-/30-mm and Shelf) to dry different type of Glassware Conveniently * Assembly and Disassembly is Easy to Carry Mobile Drying * May Prevent the Adhesion of Dust to the Curtain. * Drying, Moving, and Storage can be Comprehensively.		※ 운임 별도 <L1> <table><tr><th>Ea.</th><th>Price, ₩</th></tr><tr><td>SL.Car7061</td><td>Cart with Curtain, Displaying, 5-deck, STS, 70×45×h160cm , Pack-54×80×21cm, 23kg</td></tr><tr><td>SL.Car7062</td><td>, 90×45×h160cm , -54×100×21cm, 26kg</td></tr><tr><td>SL.Car7063</td><td>, 120×45×h160cm, -54×130×21cm, 31kg</td></tr></table>	Ea.	Price, ₩	SL.Car7061	Cart with Curtain, Displaying, 5-deck, STS, 70×45×h160cm , Pack-54×80×21cm, 23kg	SL.Car7062	, 90×45×h160cm , -54×100×21cm, 26kg	SL.Car7063	, 120×45×h160cm, -54×130×21cm, 31kg
Ea.	Price, ₩									
SL.Car7061	Cart with Curtain, Displaying, 5-deck, STS, 70×45×h160cm , Pack-54×80×21cm, 23kg									
SL.Car7062	, 90×45×h160cm , -54×100×21cm, 26kg									
SL.Car7063	, 120×45×h160cm, -54×130×21cm, 31kg									

20 & 40kg Gas Cylinder Safety Carts, for 1 & 2 Cylinder, Made of Stainless-steel or Coated Steel, 가스실린더용 안전카트

■ **Supplied by DAIHAN Korea** * Application ; Gas Cylinder Safe Moving and Storage in Hospital, Laboratory etc.

* Structure ; Special Structure for Easy to Move, Collapsible Portable



(1) C9.4402



(2) C9.4403



(3) DH.Car7006



(4) DH.Car7007

Stainless-steel Cart, for 1×20 or 40kg Cylinder

C9.4402 (1) Cart, Stainless-steel, for 1×20kg Gas Cylinder

C9.4403 (2) , for 1×40kg Gas Cylinder

Steel Cart, for 1 or 2×40kg Cylinder **NEW**

DH.Car7006 (3) Cart, Steel, for 1×40kg Gas Cylinder, Air Caster

DH.Car7007 (4) , for 2×40kg Gas Cylinder

※ **운임 별도 <L3>**

Ea. Price, ₩

<L4>

Ea. Price, ₩

DAIHAN™ 6-Place High Performance Mini-microcentrifuge Set "CF-5", Max. 5,500 rpm, *Medicaluse*

With Circular Fixed-Angle Rotor for $6 \times 0.2/0.5/1.5/2.0\text{ml}$ Tubes and Strip Rotor for $2 \times 0.2\text{ml}$ PCR 8-tubes Strips/ $16 \times 0.2\text{ml}$ PCR Tubes
 미니-마이크로 원심분리기 세트, 앵글-/스트립-로터 포함, 0.2/0.5ml 튜브 어댑터 포함

CE

Medicaluse

■ Supplied by DAIHAN Korea

* with Certificate & Traceability

* Supplied with two Rotors ;

- Circular Fixed-Angle Rotor for $6 \times 0.5/1.5/2.0\text{-ml}$ Tubes (6 $\times$ 0.2-/0.5-ml tube Adapters included)- Strip Rotor for $2 \times 0.2\text{ml}$ PCR 8-tubes Strips or 16ea of individual 0.2ml PCR Tubes

* Smooth Opening of the Lid by Simply Pressing

* Lid-Activated On/Off Safety Switch ; Close the Clear Lid to Start and Open to Stop

■ CE Certified

* Complete Set with Circular Fixed-Angle Rotor and Fixed Strip Rotor

■ PL(Product Liability) Insurance

Circular Fixed-Angle Rotor
for $6 \times 0.2\text{-}2.0\text{ml}$ Tubes

3-Year Warranty

Operation
Display LED Lamp

Specifications & Ordering information

<M1>

Spec.	Model	"CF-5"
RPM & Max. RCF		Max. 5,500 rpm, 1,775 $\times$ g
Max. Capacity		6 $\times$ 0.5-/1.5-/2.0-ml tubes, 2 $\times$ 0.2ml PCR 8-tube Strip or 16 $\times$ 0.2ml PCR 8 Tubes
Rotor (included)	Angle Rotor	6 Holes Circular Fixed-Angle Rotor for 0.5-/1.5-/2.0-ml Tubes
	Strip Rotor	16 Holes Strip Rotor for 0.2ml PCR 8-tube Strip or individual Tubes
Adapter		Included 6 $\times$ 0.2-/ 0.5-ml Tube Adapter
Material		Lid-PC, Body-ABS, Rotor-ABS
Dimension (w $\times$ d $\times$ h)		Unit ; 178 $\times$ 150 $\times$ h109 mm, Packing ; 200 $\times$ 180 $\times$ h175mm, 2kg
Power Consumption		8 W
Electric Reqts.		1 Phase, AC 100~240V, 50/60Hz
Cat. No.	230V	DH.WCF00005
	120V	DH.WCF00025
Price, Each. ₩		

DAIHAN™ 6-Place Hi-Quality Mini-microcentrifuge Set "CF-6", Max. 5,500 rpm, Rapidly Stop Spin-Motion for Safety

With Circular Fixed-Angle Rotor for $6 \times 0.2/0.5/1.5/2.0\text{ml}$ Tubes and Strip Rotor for $2 \times 0.2\text{ml}$ PCR 8-Tubes Strips/ $16 \times 0.2\text{ml}$ PCR Tubes
 미니-마이크로 원심분리기 세트, 신속한 회전동작 정지, 앵글-/스트립-로터 포함, 0.2/0.5ml 튜브 어댑터 포함

3-Year Warranty

■ Supplied by DAIHAN Korea

■ PL(Product Liability) Insurance

■ Rapidly Stop Spin-motion for user safety

* with Certificate & Traceability

* Complete Set with two Rotors ;

- Circular Fixed-Angle Rotor for $6 \times 0.5/1.5/2.0\text{-ml}$ Tubes (6 $\times$ 0.2-/0.5-ml tube Adapters included)- Fixed Strip Rotor for $2 \times 0.2\text{ml}$ PCR 8-tubes Strips or 16ea of individual 0.2ml PCR Tubes

* Smooth Opening of the Lid by Simply Pressing

* Rubberized Base

* Lid-Activated On/Off Safety Switch ; Close the Clear Lid to Start and Open to Stop

Specifications & Ordering information

<L1>

Spec.	Model	"CF-6", Rapidly Stop Spin-motion for Safety
RPM & Max. RCF		Max. 5,500 rpm, 1,775 $\times$ g
Max. Capacity		6 $\times$ 0.5-/1.5-/2.0-ml tubes, 2 $\times$ 0.2ml PCR 8-tube Strip or 16 $\times$ 0.2ml PCR 8 Tubes
Rotor (included)	Angle Rotor	6 Holes Circular Fixed-Angle Rotor for 0.5-/1.5-/2.0-ml Tubes
	Strip Rotor	16 Holes Strip Rotor for 0.2ml PCR 8-tube Strip or individual Tubes
Adapter		Included 6 $\times$ 0.2-/ 0.5-ml Tube Adapter
Material		Lid-PC, Body-ABS, Rotor-ABS
Dimension (w $\times$ d $\times$ h)		Unit ; 144 $\times$ 192 $\times$ h108mm, 1.1kg, Packing ; 175 $\times$ 255 $\times$ h215mm, 1.6kg
Power Consumption		8 W
Electric Reqts.		1 Phase, AC 100~240V, 50/60Hz
Cat. No.	230V	DH.WCF00006
	120V	DH.WCF00026
Price, Each. ₩		

Fixed Rotor & Adapter for Spare, 0.2-/0.5-ml Tube Adapter-Set, 스페어 로터/튜브 어댑터

<L1>

Ea. Price, ₩



(1)



(2)



(3)

DH.WCF00105 (1) Circular Fixed Angle Rotor, for "CF-6", "CF-5"
for $6 \times 0.2\text{-}0.5/1.5/2.0\text{-ml}$ TubesDH.WCF00205 (2) Fixed Strip Rotor, for "CF-6", "CF-5"
for $2 \times 0.2\text{ml}$ PCR8-tube Strips or $16 \times 0.2\text{ml}$ PCR TubesDH.WCF00201 (3) 0.2ml Tube Adapter-sets, 6ea Adapter/1set,
for Circular Fixed Angle Rotor, for "CF-6", "CF-5"DH.WCF00301 (3) 0.5ml Tube Adapter-sets, 6ea Adapter/1set,
for Circular Fixed Angle Rotor, for "CF-6", "CF-5"

DAIHAN™ 12-Place High Performance Pro-Microcentrifuge Set “CF-10”, Max. 13,500rpm, *Medicaluse*
 With MC Nylon Fixed Angle Rotor for $12 \times 0.2\text{-}0.5\text{-}1.5\text{-}2.0\text{-ml}$ Tubes, Compact, Quiet, Safety Electronic Lid Lock System
 프로-마이크로 원심분리기, 전자식 안전 도어 잠금 시스템, 앵글 로터 포함, 0.2/0.5ml 튜브 어댑터 포함



Max. 13,500 rpm

Medicaluse



3-Year Warranty



(2) Tube Adapters for $12 \times 0.2\text{-}0.5\text{-ml}$ Tube (Included)

- Supplied by DAIHAN Korea
- CE Certified ■ *PL(Product Liability) Insurance*
- Unbelievably Quiet Operation - Brushless DC motor
- High Performance - Digital Feedback Speed Control, Max. 13,500rpm, $12,225 \times g$
- Useful Exclusive Spin Down Function
- Included MC Nylon Fixed Angle Rotor for $1.5\text{-}2.0\text{-ml} \times 12$ Tubes and $12 \times 0.2\text{-}0.5\text{-ml}$ PCR Tube Adapters
- *Patented Jog-Shuttle Control System* and Digital LCD with Back-Light Function
- * Plasmid, DNA and RNA work, protein work, and miniprep kit, Cell Culture, cell harvesting, subcellular separations, and many other applications.
- * Unnecessary Rotor Lid for Running and Reducing Noise ($\leq 57\text{dB}$ at Max. RPM).
- * Realize Extremely Compact Design with Exceptionally Small Footprint.
- * Acceleration and Braking Times are Lower than 15 sec (at Max. Speed).
- * Storage Function for the Set Values ; RPM & Time.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be disabled).
- * Safety System for Motor Error.
- * *Safety Electronic Lid Lock System* ; Opening Lid is not applicable during running.
- * Simple Operation for Repetitive use of 'Spin Down' Function.

Specifications & Ordering information

<M1>

Spec.	Model	"CF-10"
RCF		1,000~13,500 RPM, $12,225 \times g$
Max. Capacity		$0.2\text{-}0.5\text{-}1.5\text{-}2.0\text{-ml} \times 12$ Tubes
Rotor		MC Nylon Fixed Angle Rotor for $0.2\text{-}0.5\text{-}1.5\text{-}2.0\text{-ml} \times 12$ Tubes
Drive System		Brushless DC motor, Direct drive
Noise Level		≤ 57 dB
Acceleration & Braking Time		≤ 15 sec (at Max. speed)
Cooling Device		Air Cooling System
Timer & Alarm		99hr 59min (Continuous run), Error status & Timer-end
Display		Digital LCD with Back-Light Function
Controller		Digital Feedback Control with Jog-Shuttle Switch (Turn+Push)
Safety Circuit		Motor Error Detector, Electronic Safety Lid Lock System
Material		Lid-PC, Body-ABS, Rotor-PE
Dimension & Weight		$170 \times 226 \times 183\text{mm}$, 2.3kg
Packing Size & Gross Weight		$290 \times 240 \times 220\text{mm}$, 5kg
Power Consumption		90W
Electric Reqts.		1 Phase, AC 120V, 50/60Hz or AC 230V, 50/60Hz
Cat. No.	230V	DH.WCF00010
	120V	DH.WCF00011
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

Spare Rotor & Adapter Set

- | | |
|--|---|
| <p>DH.WCF00100</p> <p>DH.WCF00200</p> <p>DH.WCF00300</p> | <p>(1) MC Nylon Fixed Angle Rotor, "CFAR1512", for "CF-10", $12 \times 0.2\text{-}0.5\text{-}1.5\text{-}2.0\text{-ml}$ Tubes</p> <p>(2) 0.2-ml Tube Adapter-sets, "CFTA1", 12ea Adapter/1set, for Fixed Angle Rotor</p> <p>(2) 0.5-ml Tube Adapter-sets, "CFTA2", 12ea Adapter/1set, for Fixed Angle Rotor</p> |
|--|---|

Ea. Price, ₩

DAIHAN™ 12-Place **Premium** Pro-Microcentrifuge Set “CF-12”, Max. 13,500rpm, 12,225 × g
With **Aluminum Die Casting Lid** & Electronic Lid Lock System, Autoclavable PP Fixed Angle Rotor for 12 × 0.2-/0.5-/1.5-/2.0-ml Tubes
프로-마이크로 원심분리기, 전자식 안전 도어 잠금 시스템, 앵글 로터 포함, 0.2/0.5ml 튜브 어댑터 포함

Max. 13,500 rpm



3-Year Warranty

“CF-12” for 12 Tubes of 0.2 ~ 2.0ml
with PP Fixed Angle Rotor

- Supplied by DAIHAN Korea
- **Patented Jog-Shuttle Control System** and **Digital LCD with Back-Light Function**
- Supplied in Complete Set ; PP Fixed Angle Rotor for 1.5-/2.0-ml × 12 Tubes and PCR Tube Adapters (12 × 0.2ml & 12 × 0.5ml)
- **Unbelievably Quiet Operation due to Optimized Air Flow Design (Noise Level ≤ 57dB)**

- * Brushless DC Motor Provides Maintenance Free Operation.
- * Perfectly Constructed Frame Structure. Strong Layered Aluminum Die Casting Lid with Auto Lock Mechanism Assuring a Safe and Secure Lock.
- * High Performance - Digital Feedback Speed Control, Max. 13,500rpm, 12,225 × g
- * Useful Exclusive Spin Down Function
- * Easy to Use Compact Design Allows for Fast Centrifugation
- * Intuitive Controls and Easy to Read Display Simplifies Operation.
- * Ideal for Fast Isolating and Pelleting of DNA, RNA, Proteins, and Cell Viruses
- * Max. 13,500rpm with 12 Position, Autoclavable PP Fixed Angle Micro-Rotor
- * Large Backlit LCD-screen for Settings and Parameters
- * The System Displays Errors, Excess Speed, and Open Lid.
- * Load Imbalance Detection.
- * Strong Steel Layered Door with Safety, Auto-Lock Function
- * Automatic Door Release
- * Biocontainment Rotor Lid Ensuring Securely Contained Samples
- * Soft Brake Function to Protect Delicate Samples
- * PP Inner Lid for Maximize Safety (Optional)

Specifications & Ordering information

<L1>



Perfectly Constructed Frame Structure,
Strong Layered Lid

Spec.	Model	“CF-12” Aluminum Die Casting Lid & Electronic Lid Lock System
RCF		1,000~13,500 RPM, 12,225×g
Max. Capacity		0.2-/0.5-/1.5-/2.0-ml × 12 Tubes
Rotor (included)		PP Fixed Angle Rotor for 0.2-/0.5-/1.5-/2.0-ml × 12 Tubes
Adapter (included)		12×0.2ml & 12×0.5ml PCR Tubes
Drive System		Brushless DC motor, Direct drive, Vibration Insulated, Counterclockwise Rotating
Noise Level		≤ 57 dB
Acceleration & Braking Time		≤ 12 sec (at Max. speed)
Cooling Device		Air Cooling System
Timer & Alarm		99hr 59min (Continuous run), Error status & Timer-end
Display		Large Digital LCD with Back-Light Function
Controller		Digital Feedback Control with Jog-Shuttle Switch (Turn+Push)
Safety Circuit		Motor Error Detector, Electronic Safety Lid Lock System
Material		Lid-Aluminum Die Casting, Body-ABS, Rotor-PE
Dimension & Weight		187×253×h155mm, 3.8kg
Packing Size & Gross Weight		230×290×h245mm, 4.3kg
Power Consumption		90W
Electric Reqts.		1 Phase, AC 120V, 50/60Hz or AC 230V, 50/60Hz

Cat. No.	230V	DH.WCF01350
	120V	DH.WCF01351
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

Spare Rotor & Adapter Set



(2) Tube Adapter(Included)
12 × 0.2/0.5ml Tube Adapter

- DH.WCF11350 (1) PP Fixed Angle Rotor, Autoclavable, for “CF-12”,
12 × 0.2-/0.5-/1.5-/2.0-ml Tubes
- DH.WCF00200 (2) 0.2ml Tube Adapter-sets, 12ea Adapter/1set,
for Fixed Angle Rotor, “CFTA1”
- DH.WCF00300 (2) 0.5ml Tube Adapter-sets, 12ea Adapter/1set,
for Fixed Angle Rotor, “CFTA2”

Ea. Price, ₩

DAIHAN™ Premium 24-Place Hi-Speed Micro Centrifuge I. “CFH-24” & II. -20℃ Refrigerated “CFH-24R”, *Medicaluse*
With Aluminum Fixed-Angle Rotor for 0.2 · 0.5 · 1.5 · 2.0ml Tubes, Fast Cooling, 27,237 × g, Max. 17,000rpm
Optional Rotors for 0.2~50ml Tubes, 프리미엄 고속 마이크로 원심분리기, 일반형 & 냉각형

Max. 17,000 rpm *Medicaluse*



I. “CFH-24”, Standard



II. “CFH-24R”, Refrigerated

Rotor & Adapters Included



(1) Fixed Angle Rotor for 24 × 0.2 ~ 2.0ml Tubes



(a) 24 × 0.2ml Adapters
24 × 0.5ml Adapters

Optional Rotors



(2) 12 × 5.0ml Tubes



(3) 18 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes



(4) 30 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes



(5) 8 × PCR 8 Strips 0.2ml Tubes



(6), (7) 4 × 15 · 50ml Tubes



(8) 8 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes



(9) 16 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes

- Supplied by DAIHAN Korea
- Available Fixed Angle Rotors & Swing Rotors
- Multi Purpose Centrifuge Designed for Various Application ;
Separating Cell & Tissues, Purifying Proteins, Isolating DNA, Centrifuging Blood & Food, etc.
- Designed for High Versatility, Accommodating Tubes up to 50ml
- Max. Speed 17,000rpm
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * Anti-Aerosol/Anti-Overflow Rotor Design for Bio-Safe (Except for 50ml Angle Rotor)
- * Automatic Sensing Rotor-type for Autotuning Max. RPM
- * Rotor Imbalanced Detection Alarm Function.
- * Intuitive 5" Full Touch TFT LCD & Quick Menu Buttons allow Easy Controlling and Monitoring.
- * Very Quiet Operation under 56dB
- * Viewing Window for RPM Calibration
- * Setting Speed Timer Mode ; Start Timer after Reaching Setting Speed for Ensuring Correct Centrifuging Time.
- * Up to 9 Program Setting for Multiple Segments Operation
- ▶ II. Refrigerated “CFH-24R” ; Precise & Stable Powerful Cooling Performance
- Fast Cooling Performance Ensure to reach 4℃ only in 5 Minutes

Specifications & Ordering information

<L1>

Spec.	Model	I. “CFH-24” 24places, Standard	II. “CFH-24R” 24places, Refrigerated
Speed Range		500~17,000rpm (Resolution ; 100rpm)	
RCF		27,237×g (Resolution ; 10×g)	
Temp Range		-	-20℃+40℃
Rotor Capacity		Included ; 24 Places Fixed Angle Rotor for 0.2 · 0.5 · 1.5 · 2.0ml Tubes And 2×24 Adapters 0.2 & 0.5ml Optional ; Fixed Angle Rotor For-12 × 5.0ml Tubes -18 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes -30 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes -8 × PCR 8 Strip 0.2ml Tubes -4 × 15 · 50ml Tubes Swing-out Rotor For-8 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes -16 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes	
Run Time		30sec~99min / Continuous	
Motor		AC Brushless Motor(3Phase), 500W	
Noise Level		≤ 56dB	≤ 65dB
Safety Devices		Dual Door Interlock, Over-speed Detection, Over-temp Detection, Automatic Internal Diagnosis	
Acceleration & Breaking Time		≤ 20sec (at Max. Speed)	
Dimension & Weight		332×480×h283mm, 38kg	332×620×h283mm, 45kg
Packing Size & Gross Weight		455×700×h400mm, 43kg	455×925×h520mm, 50kg
Power Consumption		1.5kW	2kW
Electric Reqts.		AC220V~240V, 50Hz/60Hz	
Cat. No.		DH.CFH24	DH.CFH24R
Price, Each. ₩			

Spare Rotor and Tube Adapters

- DH.CFH24R1 (1) Aluminum Fixed Angle Rotor, 24 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes, Max 17,000rpm (RCF 27,237 × g)
DH.CFH24A1 (a) 0.2ml Tube Adapter-set, 24 Adapters/set, for Fixed Angle Rotor
DH.CFH24A2 (a) 0.5ml Tube Adapter-set, 24 Adapters/set, for Fixed Angle Rotor

Optional Rotors

- DH.CFH24R2 (2) Aluminum Fixed Angle Rotor, 12 × 5.0ml Tubes, Max 15,000rpm (RCF 22,967 × g)
DH.CFH24R3 (3) Aluminum Fixed Angle Rotor, 18 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes, Max 15,000rpm (RCF 18,891 × g)
DH.CFH24R4 (4) Aluminum Fixed Angle Rotor, 30 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes, Max 15,000rpm (RCF 24,249 × g)
DH.CFH24R5 (5) Aluminum Fixed Angle Rotor, 8 × PCR 8 Strip 0.2ml Tubes, Max 12,500rpm (RCF 11,722 × g)
DH.CFH24R8 (6) Aluminum Fixed Angle Rotor, 4 × 15ml · 50ml for Conical Tube, Max 6500 rpm (50ml ; RCF 4,411 × g / 15ml ; 4,316 × g)
DH.CFH24R9 (7) Aluminum Fixed Angle Rotor, 4 × 15ml · 50ml for Round Tube, Max 6500 rpm (50ml ; RCF 4,411 × g / 15ml ; 4,316 × g)
DH.CFH24R6 (8) Aluminum Swing Rotor, 8 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes, Max 15,000rpm (RCF 16,343 × g)
DH.CFH24R7 (9) Aluminum Swing Rotor, 16 × 0.2 · 0.5 · 1.5 · 2.0ml Tubes, Max 15,000rpm (RCF 15,682 × g)

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ MicroHematocrit Centrifuge “CEF-MHT”, 24 Places 12,000rpm, for Capillaries, **Medicaluse**
Without Rotor, 13,500×g, 미세적혈구용적 원심분리기, Rotor 별도

CE

Medicaluse

Max. 12,000 rpm

**3-Year Warranty**

“CEF-MHT” MicroHematocrit Centrifuge
* Rotor is Optional !



Optional PP Readable Rotor

- Supplied by DAIHAN Korea ■ CE Certified
- 혈액 내 적혈구 용적 측정 및 용액에서 혈액의 분리를 위해 사용.
- Ideal for Determining the Value of Blood cells in the Blood and It can Separate Blood from the Solution.
- Permanently Brushless(BLDC) Motor.
- * Microcomputer Control, Stable Operation, Low Noise.
- * LED Display, User-friendly Interface, Simple and Convenient Operation.
- * The Real-time Conversation of Speed and RCF is Convenient for Operation.

Specifications & Ordering information

<M1>

Spec.	Model	“CEF-MHT”
Max. RPM & RCF		12,000rpm, 13,500×g
Speed Accuracy		±30rpm
Rotor Capacity		24 pieces of Capillaries
Time Setting		1~99min
Motor		Permanently Brushless(BLDC) Motor
Noise Level		≤65dB
Safety Devices		Self-locking, Over-speed Safety, Automatic Alarm
Dimension & Weight		230×310×h240mm, 9kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef24MHT
Price, Each. ₩		

Optional Rotor

Ea. Price, ₩

DH.Cef24MHT.R PP(Polypropylene) Readable Rotor, 12,000rpm, for 24-capillaries



(1) (2)



Package Box

DAIHAN mediclin™ 75 µl Micro Hematocrit Capillary Tube, Blue & Heparinized-Red
For Blood Taking, Disposable Hi-Quality Glass, $\phi 1.55 \times L75\text{mm}$, 일회용 Micro 헤마토크리트 튜브

- Supplied by DAIHAN Korea ■ Made of High Grade Quality Glass in Compliance with DIN ISO 12772
- 혈액 중 적혈구의 비중을 측정하기 위해 헤마토크리트 검사에 사용, 원심분리기를 사용하여 혈액 성분 분리
- Color Coded ; (1) Blue “Non-Heparinized” (2) Red “Heparinized”
- * ID ; 1.1~1.2mm, OD ; 1.5~1.6mm * 100pcs/Vial Pk., 10Pk./Cs.

<M1>

(1) Blue Coded ; Non-Heparinized, 75 µl

DM.Tub2001 Capillary Tube, without Sodium Heparine, Blue, 75 µl

Pk.	Price, ₩	Cs.	Price, ₩
100ea		10pk	

(2) Red Coded ; Heparinized, 75 µl, 헤파린 처리

- 헤파린(항응고제) 처리된 튜브는 혈액이 응고되는 것을 방지
- 채혈과 헤마토크리트 검사 수행 사이 지연이 예상되는 경우 특히 유용
- * Useful when a Delay between Blood Sampling and Hematocrit Test is Expected.
- * Heparinized Tubes Prevent Blood from Clotting.

DM.Tub2002 Capillary Tube, with Sodium Heparine, Red, 75 µl

Pk.	Price, ₩	Cs.	Price, ₩
100ea		10pk	



(1) (2)



(3)

Heinz® 75 µl Premium MicroHematocrit Capillary Tube, Disposable Glass, Blue & Red-Hepa
For Blood Taking, $\phi 1.55 \times L75\text{mm}$, 일회용 Micro 헤마토크리트 튜브

- Supplied by Heinz Herenz® in Germany
- 혈액 중 적혈구의 비중을 확인하는 헤마토크리트 검사에 사용되며, 원심분리기와 함께 사용
- * 100pcs/Vial Pk., 10Pk./Cs.

<L1>

(1) Blue Coded ; Non-Heparinized, 75 µl

HH.1080801 Micropipet, 75 µl, Capill., Blue

Pk.	Price, ₩	Cs.	Price, ₩
100ea		10pk	

(2) Red Coded ; Heparinized, 75 µl

- * 헤파린 처리된 튜브는 혈액이 응고되는 것을 방지
- * Heparinized Tubes Prevent Blood from Clotting.

HH.1080802 Micropipet, Heparin, 75 µl, Capill., Red

Pk.	Price, ₩	Cs.	Price, ₩
100ea		10pk	

(3) Sealing Waxplate

- * For Hematocrit Tubes, 2pcs/Bag Pk.

HH.1131821 Sealing Waxplate

Ea.	Price, ₩	Pk.	Price, ₩
		2ea	

DAIHAN™ 7 · 10 · 15ml Clinical Centrifuge “CEF-15”, 4,000rpm, Max. 15ml × 8, *Medicaluse*

For 7 · 10 · 15ml Tubes, 1,980 × g, without Rotor, 임상용 원심분리기, Angle Rotor 별도



Medicaluse
Max. 4,000 rpm



3-Year Warranty

“CEF-15” for 7 · 10 · 15ml Tubes
✳ Rotor is Optional !



(1) for 7ml × 12



(2) for 10ml × 12



(3) for 15ml × 8

Optional PP Angle Rotors

- Supplied by DAIHAN Korea
- 클리닉, 연구, 대학에서 사용하기에 적합한 제품.
- 휴대용, 저소음 설계로 제작되어 이동이 용이하고 최소의 소음을 보장.
- Designed with Portable Body and Low Noise
- Ideal for Clinic, Research and University.
- Permanently Brushless(BLDC) Motor.
- * Microcomputer Control
- * Programmable and LED Display.
- * Display Speed(RPM) and RCF Mode.
- * ABS Housing make the Clean Easy.

Specifications & Ordering information

<M1>

Spec.	Model	“CEF-15”
Max. RPM & RCF		4,000rpm, 1,980 × g
Speed Accuracy		± 30rpm
Rotor Capacity		Option ; (1) 7ml × 12 Tubes, (2) 10ml × 12 Tubes, (3) 15ml × 8 Tubes
Timer		1~99min
Motor & Noise Level		Permanently Brushless(BLDC) Motor, ≤ 65dB
Safety Devices		Lid Safety, Speed Safety, Errors Detection
Dimension & Weight		310 × 260 × h220mm, 14kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef15
Price, Each. ₩		

Optional PP Angle Rotors

- DH.Cef15.R1 (1) PP(Polypropylene) Angle Rotor “R1-15” ; for 7ml × 12 Tubes
 DH.Cef15.R2 (2) “R2-15” ; for 10ml × 12 Tubes
 DH.Cef15.R4 (3) “R3-15” ; for 15ml × 8 Tubes

Ea. Price, ₩

DAIHAN™ 1.5~50ml General Purpose Clinical Centrifuge “CEF-50”, 4,000rpm, Max. 50ml × 6, **Medicaluse**
For 1.5~50ml Tubes, 2,220 × g, without Rotor, 범용 원심분리기, Angle Rotor 별도

CE

Medicaluse
Max. 4,000 rpm



3-Year Warranty

“CEF-50” for 1.5~50ml Tubes
※ Rotor is Optional !



(1) for 7ml × 24



(2) for 10ml × 24



(3) for 15ml × 12



(4) for 50ml × 6

Optional Angle Rotors



① for 1.5ml



② for 5ml



③ for 5/7ml



④ for 1.5ml



⑤ for 5ml



⑥ for 5/7ml



⑦ for 5/7ml



⑧ for 15ml



⑨ for 20ml

Optional Adapters (6ea/set)

- Supplied by DAIHAN Korea ■ CE Certified
- 임상의학, 생화학, 유전공학, 면역학 등의 분야에서 널리 사용 가능.
- LCD 디스플레이, 사용자 친화적 인터페이스, 단순하고 편리한 조작.
- Ideal for the Field of Clinical Medicine, Biochemistry, Genetic Engineering, Immunology, etc.
- Permanently Brushless(BLDC) Motor.
- * Microcomputer Control, Stable Operation, Low Noise and High Speed Accuracy.
- * Touch Panel, Programmable Operation
- * Operation Data can be Set According to Demands and Stored Automatically.
- * LCD Display, User-friendly Interface(The Real-time Display of Speed and RCF)
- * Simple and Convenient for Operation.

Specifications & Ordering information

<M1>

Spec.	Model	“CEF-50”
Max. RPM & RCF		4,000rpm, 2,220 × g
Speed Accuracy		± 30rpm
Rotor Capacity		Option ; (1) 7ml × 24 Tubes, (2) 10ml × 24 Tubes, (3) 15ml × 12 Tubes, (4) 50ml × 6 Tubes
Timer		1~99min
Motor		Permanently Brushless(BLDC) Motor
Noise Level		≤ 60dB
Safety Devices		Self-Locking, Lid Safety, Over-speed Safety, Automatic Alarm
Dimension & Weight		330 × 420 × h280mm, 14kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef50
Price, Each. ₩		

Optional Angle Rotors

- DH.Cef50.R1 (1) PP(Polypropylene) Angle Rotor “R1-50” ; for 7ml × 24 Tubes
 DH.Cef50.R2 (2) “R2-50” ; for 10ml × 24 Tubes
 DH.Cef50.R3 (3) “R3-50” ; for 15ml × 12 Tubes
 DH.Cef50.R4 (4) Aluminium Alloy Angle Rotor “R4-50” ; for 50ml × 6 Tubes

Ea. Price, ₩

Optional Adapters

- DH.Cef50.A1 MC Nylon Adapter for “R2-50” ; ① 1.5ml × 24ea/set, id ϕ 11 × h37mm
 DH.Cef50.A2 ; ② 5ml × 24ea/set, id ϕ 14 × h46mm
 DH.Cef50.A3 ; ③ 5/7ml(Vacuum Tube) × 24ea/set, id ϕ 13.5 × h83mm
 DH.Cef50.A4 for “R3-50” ; ④ 1.5ml × 12ea/set, id ϕ 11 × h37mm
 DH.Cef50.A5 ; ⑤ 5ml × 12ea/set, id ϕ 14 × h46mm
 DH.Cef50.A6 ; ⑥ 5/7ml × 12ea/set, id ϕ 13.5 × h84mm
 DH.Cef50.A7 for “R4-50” ; ⑦ 5/7ml × 6ea/set, id ϕ 17.5 × h78mm
 DH.Cef50.A8 ; ⑧ 15ml × 6ea/set, id ϕ 17.5 × h100mm
 DH.Cef50.A9 ; ⑨ 20ml × 6ea/set, id ϕ 22.3 × h100mm

Set. Price, ₩

DAIHAN™ 0.2~100ml High Speed Centrifuge “CEf-100H”, 16,500rpm, Max. 100ml × 4, *Medicaluse*
Without Rotor, 18,930 × g, 고속 원심분리기, Rotor 별도

CE


Medicaluse
Max. 16,500 rpm



3-Year Warranty

“CEf-100H” for 0.2~100ml Tubes
※ Rotor is Optional !

Optional Rotors & Adapters



- Supplied by DAIHAN Korea ■ CE Certified
- 임상의학, 생화학, 유전공학, 면역학 등의 분야에서 널리 사용 가능
- LCD 디스플레이, 사용자 친화적 인터페이스, 단순하고 편리한 조작.
- Ideal for the Field of Clinical Medicine, Biochemistry, Genetic Engineering, Immunology, etc.
- Permanently Brushless(BLDC) Motor.
- * Microcomputer Control, Stable Operation, Low Noise and High Speed Accuracy.
- * Touch Panel, Programmable Operation
- * Operation Data can be Set According to Demands and Stored Automatically.
- * LCD Display, User-friendly Interface, Simple and Convenient for Operation.
- * 9 Files of Acceleration and Deceleration Levels for Choice.

Specifications & Ordering information

<M1>

Spec.	Model	“CEf-100H”
Max. RPM & RCF		16,500rpm, 18,930 × g
Speed Accuracy		± 30rpm
Rotor Capacity		(1) 1.5/2.2ml × 12 Tubes (2) 1.5/2.2ml × 24 Tubes (3) PCR 0.2ml 8-strip × 3 (4) 5ml × 12 Tubes (5) 10ml × 12 Tubes (6) 15ml × 8 Tubes (7) 50ml Round-bottom × 6 Tubes (8) 50ml Conical-bottom × 6 Tubes (9) 100ml × 4 Tubes
Timer		1~99min
Motor		Permanently Brushless(BLDC) Motor
Noise Level		≤ 65dB
Safety Devices		Self-Locking, Over-speed Safety, Automatic Alarm, etc.
Dimension & Weight		330 × 420 × h280mm, 25kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef100H
Price, Each. ₩		

Optional Aluminium Alloy Angle Rotors

		Ea. Price, ₩
DH.Cef100H.R1	(1) Option Rotor “R1-100H” ; for 1.5/2.2ml × 12, 16500rpm, 18930xg	
DH.Cef100H.R2	(2) “R2-100H” ; for 1.5/2.2ml × 24, 13200rpm, 16110xg	
DH.Cef100H.R3	(3) “R3-100H” ; for PCR 0.2ml 8-strip × 3, 15000rpm, 13800xg	
DH.Cef100H.R4	(4) “R4-100H” ; for 5ml × 12, 13000rpm, 11400xg	
DH.Cef100H.R5	(5) “R5-100H” ; for 10ml × 12, 12000rpm, 14800xg	
DH.Cef100H.R6	(6) “R6-100H” ; for 15ml × 8, 12000rpm, 14800xg	
DH.Cef100H.R7	(7) “R7-100H” ; for 50ml Round-bottom × 6, 11000rpm, 13000xg	
DH.Cef100H.R8	(8) “R8-100H” ; for 50ml Conical-bottom × 6, 11000rpm, 13000xg	
DH.Cef100H.R9	(9) “R9-100H” ; for 100ml × 4, 10000rpm, 11570xg	

Optional Adapters

		Set. Price, ₩
DH.Cef100H.A1	MC Nylon Adapter for “R1-100H” ; Ⓐ 0.2ml × 12ea/set	
DH.Cef100H.A2	; Ⓑ 0.5ml × 12ea/set	
DH.Cef100H.A3	for “R2-100H” ; Ⓒ 0.2ml × 24ea/set	
DH.Cef100H.A4	; Ⓓ 0.5ml × 24ea/set	
DH.Cef100H.A5	for “R4-100H” ; Ⓔ 0.5ml × 12ea/set	
DH.Cef100H.A6	; Ⓕ 1.5ml × 12ea/set	
DH.Cef100H.A7	for “R5-100H” ; Ⓖ 1.5ml × 12ea/set	
DH.Cef100H.A8	; Ⓗ 5ml × 12ea/set	
DH.Cef100H.A9	for “R6-100H” ; Ⓙ 1.5ml × 8ea/set	
DH.Cef100H.A10	; Ⓚ 5ml × 8ea/set	
DH.Cef100H.A11	for “R7-100H”, “R8-100H” ; Ⓛ 15ml × 6ea/set	
DH.Cef100H.A12	; Ⓜ 20ml × 6ea/set	
DH.Cef100H.A13	for “R9-100H” ; Ⓨ 15ml × 4ea/set	
DH.Cef100H.A14	; Ⓩ 20ml × 4ea/set	

※ Adapter Size (id φ × h) ;

Ⓐ φ 6.5 × h20mm, Ⓑ φ 8 × h28mm, Ⓒ φ 6.5 × h17.5mm, Ⓓ φ 8 × h28mm, Ⓔ φ 8 × h28mm,
Ⓕ φ 11 × h37mm, Ⓖ φ 11 × h37mm, Ⓗ φ 14 × h48mm, Ⓙ φ 11 × h37mm, Ⓚ φ 17 × h47mm,
Ⓛ φ 17.6 × h97mm, Ⓜ φ 22 × h93mm, Ⓨ φ 17.6 × h97mm, Ⓩ φ 22 × h91mm

DAIHAN™ 1.5~500ml High Volume Swing-out Centrifuge “CE-500”, 4,000rpm, Max. 500ml × 4, **Medicaluse**
For 1.5~500ml Tubes/Bottles and Microplates, 3,200 × g, without Rotor, 다목적 원심분리기, Rotor 별도

CE

Medicaluse
Max. 4,000 rpm



3-Year Warranty

“CE-500” for 1.5~500ml Tubes/Bottles
※ Rotor is Optional !



(2) PP Bottle
500ml × 4

(1) Swing Rotor
(Option)

(3) Microplate Rotor
(Option)



① 1.5ml × 20 hole ② 5/7ml × 20 hole ③ 10ml × 15 hole ④ 15ml × 8 hole
⑤ 20ml × 9 hole ⑥ 50ml × 4 hole ⑦ 100ml × 2 hole ⑧ 250ml × 1 hole

Optional Adapters (4ea/set)

- Supplied by DAIHAN Korea ■ CE Certified
- 임상의학, 생화학, 유전공학, 면역학 등의 분야에서 널리 사용 가능.
- LCD 디스플레이, 사용자 친화적 인터페이스, 단순하고 편리한 조작.
- Ideal for the Field of Clinical Medicine, Biochemistry, Genetic Engineering, Immunology, etc.
- * Microcomputer Control, Variable-conversion Motor, Stable Operation, Low Noise and High Speed Accuracy.
- * Touch Panel, Programmable Operation
- * Operation Data can be Set According to Demands and Stored Automatically.
- * LCD Display, User-friendly Interface(The Real-time Conversation of Speed and RCF)
- * Simple and Convenient for Operation.

Specifications & Ordering information

<M1>

Spec.	Model	“CE-500”
Max. RPM & RCF		4,000rpm, 3,200×g
Speed Accuracy		±50rpm
Rotor Capacity		(1) 500ml × 4 Buckets (2) 8 × 96well Microplates
Timer		1~99min
Motor		Variable-conversion Motor
Noise Level		≤65dB
Safety Devices		Self-Locking, Lid Safety, Over-speed Safety, Automatic Alarm
Dimension & Weight		460×540×h340mm, 31kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef500
Price, Each. ₩		

Optional Rotors and Bottle

- DH.Cef500.R1 (1) Stainless-steel Swing Rotor “R1-500” with Bucket(500ml × 4)
DH.Cef500.B1 (2) PP Bottle with Screw Cap, 500ml for “R1-500”
DH.Cef500.R2 (3) Stainless-steel Microplate Rotor(8 × 96well) “R2-500”, 2200xg

Ea. Price, ₩

Optional Polypropylene Adapters

- DH.Cef500.A1 PP(Polypropylene) Adapter for Rotor “R1-500”
; ① (1.5ml × 20hole) × 4ea/set, id ϕ 11 × h37mm
DH.Cef500.A2 PP(Polypropylene) Adapter for Rotor “R1-500”
; ② (5/7ml × 20hole) × 4ea/set, id ϕ 13 × h57mm
DH.Cef500.A3 PP(Polypropylene) Adapter for Rotor “R1-500”
; ③ (10ml × 15hole) × 4ea/set, id ϕ 14 × h37mm
DH.Cef500.A4 PP(Polypropylene) Adapter for Rotor “R1-500”
; ④ (15ml × 8hole) × 4ea/set, id ϕ 17 × h97mm
DH.Cef500.A5 PP(Polypropylene) Adapter for Rotor “R1-500”
; ⑤ (20ml × 9hole) × 4ea/set, id ϕ 21 × h97mm
DH.Cef500.A6 PP(Polypropylene) Adapter for Rotor “R1-500”
; ⑥ (50ml × 4hole) × 4ea/set, id ϕ 30 × h92mm
DH.Cef500.A7 PP(Polypropylene) Adapter for Rotor “R1-500”
; ⑦ (100ml × 2hole) × 4ea/set, id ϕ 39 × h90mm
DH.Cef500.A8 PP(Polypropylene) Adapter for Rotor “R1-500”
; ⑧ (250ml × 1hole) × 4ea/set, id ϕ 63 × h90mm

Set. Price, ₩

DAIHAN™ 0.2~100ml Refrigerated High Speed Pro-Centrifuge “CEF-100HR”, -20℃ +40℃, 18,000rpm, Max. 100ml × 4
Without Rotor, 23,755 × g, 고속 냉장 원심분리기, Rotor 별도

CE

Max. 18,000 rpm



3-Year Warranty

“CEF-100HR” for 0.2~100ml Tubes
※ Rotor is Optional !

- Supplied by DAIHAN Korea ■ CE Certified
- 과학 연구실험과 고등교육 등의 분야에서 널리 사용 가능.
- LCD 디스플레이, 사용자 친화적 인터페이스, 단순하고 편리한 조작.
- Ideal for Centrifugal Separation used by Scientific Research Units and Institutions of Higher Learning.
- Permanently Brushless(BLDC) Motor.


- * Microcomputer Control, Stable Operation, Low Noise and High Speed Accuracy.
- * Touch Panel, Programmable Operation * Simple and Convenient for Operation.
- * Operation Data can be Set According to Demands and Stored Automatically.
- * LCD Display, User-friendly Interface(The Real-time Conversation of Speed and RCF)
- * 9 Files of Acceleration and Deceleration Levels for Choice.

Specifications & Ordering information

<L1>

Spec.	Model	“CEF-100HR”
Max. RPM & RCF		18,000rpm, 23,755 × g
Speed Accuracy		± 30rpm
Rotor Capacity	(1) 1.5/2.2ml × 12 Tubes	(2) 1.5/2.2ml × 24 Tubes
	(3) 0.5ml × 48 Tubes	(4) 5ml × 48 Tubes
	(5) 10ml × 12 Tubes	(6) 15ml × 8 Tubes
	(7) 50ml Round-bottom × 6 Tubes	(8) 50ml Conical-bottom × 6 Tubes
	(9) 100ml × 4 Tubes	
Timer		1min~9h 59min
Temperature Setting		-20℃ ~ +40℃
Temperature Accuracy		± 1℃
Refrigerant Unit		Non-CFC R134a
Motor		Permanently Brushless(BLDC) Motor
Noise Level		≤ 65dB
Safety Devices		Self-Locking, Over-speed Safety, Over-temperature Safety, Imbalance Safety, Automatic Alarm, etc.
Dimension & Weight		600×470×h330mm, 50kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef100HR
Price, Each. ₩		

Optional Rotors & Adapters



Optional Aluminium Alloy Angle Rotors

DH.Cef100HR.R1	(1) Option Rotor “R1-100HR”; for 1.5/2.2ml × 12, 18000rpm, 22530xg
DH.Cef100HR.R2	(2) “R2-100HR”; for 1.5/2.2ml × 24, 16000rpm, 23755xg
DH.Cef100HR.R3	(3) “R3-100HR”; for 0.5ml × 48, 13000rpm, 16760xg
DH.Cef100HR.R4	(4) “R4-100HR”; for 5ml × 12, 13000rpm, 11400xg
DH.Cef100HR.R5	(5) “R5-100HR”; for 10ml × 12, 13000rpm, 16760xg
DH.Cef100HR.R6	(6) “R6-100HR”; for 15ml × 8, 13000rpm, 16760xg
DH.Cef100HR.R7	(7) “R7-100HR”; for 50ml Round-bottom × 6, 12000rpm, 13000xg
DH.Cef100HR.R8	(8) “R8-100HR”; for 50ml Conical-bottom × 6, 12000rpm, 13000xg
DH.Cef100HR.R9	(9) “R9-100HR”; for 100ml × 4, 12000rpm, 13000xg

Ea. Price, ₩

Optional Adapters

DH.Cef100HR.A1	MC Nylon Adapter for “R1-100HR”; Ⓐ 0.2ml × 12ea/set
DH.Cef100HR.A2	Ⓑ 0.5ml × 12ea/set
DH.Cef100HR.A3	for “R2-100HR”; Ⓒ 0.2ml × 24ea/set
DH.Cef100HR.A4	Ⓓ 0.5ml × 24ea/set
DH.Cef100HR.A5	for “R4-100HR”; Ⓔ 0.5ml × 12ea/set
DH.Cef100HR.A6	Ⓕ 1.5ml × 12ea/set
DH.Cef100HR.A7	for “R5-100HR”; Ⓖ 1.5ml × 12ea/set
DH.Cef100HR.A8	Ⓗ 5ml × 12ea/set
DH.Cef100HR.A9	for “R6-100HR”; Ⓙ 1.5ml × 8ea/set
DH.Cef100HR.A10	Ⓚ 5ml × 8ea/set
DH.Cef100HR.A11	for “R7-100HR”, “R8-100HR”; Ⓛ 15ml × 6ea/set
DH.Cef100HR.A12	Ⓜ 20ml × 6ea/set
DH.Cef100HR.A13	for “R9-100HR”; Ⓝ 15ml × 4ea/set
DH.Cef100HR.A14	Ⓟ 20ml × 4ea/set

Set. Price, ₩

※ Adapter Size (id φ × h) ;

- Ⓐ φ 6.5 × h20mm, Ⓑ φ 8 × h28mm, Ⓒ φ 6.5 × h18mm, Ⓓ φ 8 × h28mm, Ⓔ φ 11 × h28mm,
- Ⓕ φ 11 × h37mm, Ⓖ φ 14 × h37mm, Ⓗ φ 14 × h47mm, Ⓙ φ 14 × h47mm, Ⓚ φ 14 × h47mm,
- Ⓛ φ 17 × h100mm, Ⓜ φ 22 × h92mm, Ⓝ φ 17 × h76mm, Ⓟ φ 22 × h92mm

DAIHAN™ 0.2~250ml **Premium High Speed Centrifuge I. "CEF-H250" & II. -20℃ Refrigerated "CEF-H250R", Medicaluse**
For 0.2~250ml Tubes · Microtube Plates · Bottles, Max 16,584 × g, 12,000rpm, 다목적 대용량 고속 원심분리기, 일반형 & 냉각형, Rotor 별도

NEW

Max. 12,000 rpm **Medicaluse**



3-Year Warranty

I. "CEF-H250", Standard
※ Rotors-Option



3-Year Warranty

II. "CEF-H250R", Refrigerated
※ Rotors-Option

- Supplied by DAIHAN Korea ■ **PL(Product Liability) Insurance**
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Available Various & Wide Range Fixed Angle Rotors & Swing Rotors
- Multi Purpose Centrifuge Designed for Various Application ; Separating Cell & Tissues, Clinical Medicine, Biochemistry, Purifying Proteins, Isolating DNA, Centrifuging Blood & Food, etc.
- * Anti-Aerosol/Anti-Overflow Rotor Design for Bio-Safe * Rotor Imbalanced Detection Alarm Function.
- * Automatic Sensing Rotor-type for Autotuning Max. RPM
- * Intuitive 5" Full Touch TFT LCD & Quick Menu Buttons allow Easy Controlling and Monitoring.
- * Very Quiet Operation under 60dB * Viewing Window for RPM Calibration
- * Setting Speed Timer Mode ; Start Timer after Reaching Setting Speed for Ensuring Correct Centrifuging Time.
- * Up to 100 Program Setting for Multiple Segments Operation
- ▶ **II. Refrigerated "CEF-H250R" ; Precise & Stable Powerful Cooling Performance**
- Fast Cooling Performance Ensure to reach 4℃ only in 5 Minutes

Specifications & Ordering information

<L1>

Spec.	Model	I. "CEF-H250" Standard	II. "CEF-H250R" Refrigerated
Speed Range & Max RCF		Fixed Angle Rotors 500~12,000rpm(16,584×g), Swing Rotors 500~5,000rpm(5,394×g)	
Temp Range		-	-20℃+40℃
Optional Rotors Capacity		* Fixed Angle Rotor for 6×85ml Tubes, 8×50ml Tubes, 12×15ml Tubes 24×12.5ml Tubes, 30×0.2·0.5·1.5·2.0ml Tubes * Swing-Out Rotor for 6×100ml Buckets, 4×250ml Buckets, 4×MTP Buckets	
Timer & Program Mode		0~100hr, Pulse, or Continuous, Available Max 100Programs	
Motor		AC Brushless Motor(3Phase), 500W	
Noise Level		≤ 60dB	≤ 65dB
Safety Devices		Dual Door Interlock, Over-speed Detection, Over-temp Detection, Automatic Internal Diagnosis	
Acceleration & Breaking Time		≤ 20sec (at Max. Speed)	
Dimension & Weight		475×625×h392mm, 40kg	657×625×h392mm, 82kg
Packing Size & Gross Weight		760×730×h590mm, 45kg	760×730×h590mm, 87kg
Power Consumption		1.5kW	2.2kW
Electric Reqts.		AC220V~240V, 50Hz/60Hz	
Cat. No.		DH.CefH250	DH.CefH250R
Price, Each. ₩			

A. Optional Fixed Angle Rotors and Adapters

Image	Cat. No.	Item	Capacity	Max. RPM	Max. RCF	Price, ₩
	DH.CefH250.1	Fixed Angle Rotor, Aluminum	(1) 30×0.2 / 0.5 / 1.5 / 2.0ml Conical Tubes	12,000rpm	15,520×g	
	DH.CFH24A1	Adapter for "DH.CefH250.1" Rotor	for 0.2ml Tube, 30ea/set (id ϕ3.2×h10mm)	-	-	
	DH.CFH24A2		for 0.5ml Tube, 30ea/set (id ϕ3.9×h10mm)	-	-	
	DH.CefH250.2	Fixed Angle Rotor, Aluminum	(2) 24×15ml Conical Tubes	4,000rpm	3,134×g	
	DH.CefH250.3	Fixed Angle Rotor, Aluminum	(3) 12×15ml Conical Tubes	12,000rpm	16,584×g	
	DH.CefH250.4	Fixed Angle Rotor, Aluminum	(4) 8×50ml Conical Tubes	10,000rpm	12,298×g	
	DH.CefH250.5	Fixed Angle Rotor, Aluminum	(5) 8×50ml Round Tubes	12,000rpm	16,584×g	
	DH.CefH250.6	Fixed Angle Rotor, Aluminum	(6) 6×85ml Round Tubes	10,000rpm	12,298×g	

<Continued On Next Page>

Articles

Cat. No

Description

<Continued From...DAIHAN™ 0.2~250ml Premium High Speed Centrifuge I.“CEF-H250” & II.“CEF-H250R”>

B. Optional Aluminum Swing Rotors · Buckets · Tube Adapters

<L1>

Image	Cat. No.	Item	Capacity	Max. RPM	Max. RCF	Price, ₩
 (7)	DH.CefH250.7	Swing Rotor, Aluminum	(7) 4 Buckets for 15ml & 50ml Tubes or Microtube Plates	4,000rpm	3,134×g	
 (8)	DH.CefH250.8	Square Bucket for "DH.CefH250.7" Rotor, Aluminum	(8) For 15ml & 50ml Tubes, 4ea/set (85×67×h95mm)	-	-	
 (9)	DH.CefH250.9	Adapters for "DH.CefH250.8" Square Bucket, POM (85×67×h95mm)	(9) 8×15ml Conical Tubes, 4ea/set (id 17×h80mm)	-	-	
 (10)	DH.CefH250.10		(10) 8×15ml Round Tubes, 4ea/set (id 17×h80mm)	-	-	
 (11)	DH.CefH250.11		(11) 4×50ml Conical Tubes, 4ea/set (id 30×h80mm)	-	-	
 (12)	DH.CefH250.12		(12) 4×50ml Round Tubes, 4ea/set (id 30×h80mm)	-	-	
 (13)	DH.CefH250.13	MTP Bucket for "DH.CefH250.7" Rotor, Aluminum	(13) For Microtube Plates(MTP), 4ea/set (87×124×h71mm)	-	-	
 (14)	DH.CefH250.14	Swing Rotor, Aluminum	(14) 6 Buckets for 15·50·100ml Tubes	4,000rpm	3,105×g	
 (15)	DH.CefH250.15	Round Bucket for "DH.CefH250.14" Rotor, Aluminum	(15) For 100ml Tubes, 6ea/set (id 47×h100mm)	-	-	
 (16)	DH.CefH250.16	Adapters for "DH.CefH250.15" Round Bucket, POM (id 47×h100mm)	(16) 3×15ml Conical Tubes, 6ea/set (id 17×h98mm)	-	-	
 (17)	DH.CefH250.17		(17) 3×15ml Round Tubes, 6ea/set (id 17×h98mm)	-	-	
 (18)	DH.CefH250.18		(18) 1×50ml Conical Tubes, 6ea/set (id 29.5×h82.5mm)	-	-	
 (19)	DH.CefH250.19		(19) 1×50ml Round Tubes, 6ea/set (id 29.5×h82.5mm)	-	-	
 (20)	DH.CefH250.20	Adapter-Free Bucket for "DH.CefH250.14" Rotor, Aluminum	(20) 2×50ml Tubes, 6ea/set (id 30×h88mm)	-	-	
 (21)	DH.CefH250.21	Swing Rotor, Aluminum	(21) 4 Buckets for 250ml Tubes	5,000rpm	5,394×g	
 (22)	DH.CefH250.22	Adapter-Free Round Bucket for "DH.CefH250.21" Rotor, Aluminum	(22) 1×250ml Tubes, 4ea/set (id 62×h106mm)	-	-	

DAIHAN™ 1.5~250ml Refrigerated High Volume Swing-out Centrifuge “**CEF-250R**”, -5℃ +35℃, 4,300rpm, Max. 250ml × 4
For 1.5~250ml Tubes/Bottles and Microplates, 3,582 × g, without Rotor, **대용량 냉장 원심분리기, Rotor 별도**

CE

Max. 4,300 rpm



3-Year Warranty

“**CEF-250R**” for 1.5~250ml Tubes/Bottles
※ **Rotor is Optional !**

- Supplied by DAIHAN Korea ■ CE Certified
- 과학 연구실험과 고등교육 등의 분야에서 널리 사용 가능.
- LCD 디스플레이, 사용자 친화적 인터페이스, 단순하고 편리한 조작.
- Ideal for Centrifugal Separation used by Scientific Research Units and Institutions of Higher Learning.
- Permanently Brushless(BLDC) Motor.
- * Microcomputer Control, Stable Operation, Low Noise and High Speed Accuracy.
- * Touch Panel, Programmable Operation
- * Operation Data can be Set According to Demands and Stored Automatically.
- * LED Display, User-friendly Interface(The Real-time Conversation of Speed and RCF)
- * Simple and Convenient for Operation. * 9 Files of Acceleration and Deceleration Levels for Choice.

Specifications & Ordering information

<L1>


Spec.	Model	“CEF-250R”
Max. RPM & RCF		4,300rpm, 3,582×g
Speed Accuracy		±30rpm
Rotor Capacity		(1) 250ml × 4 Buckets (2) 8 × 96well Microplate
Timer		1~99min
Temperature Setting		-5℃ ~ +35℃
Temperature Accuracy		±1℃
Refrigerant Unit		Non-CFC R134a
Motor		Permanently Brushless(BLDC) Motor
Noise Level		≤65dB
Safety Devices		Self-Locking, Over-speed Safety, Over-temperature Safety, Imbalance Safety, Automatic Alarm, etc.
Dimension & Weight		460×710×h420mm, 80kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef250R
Price, Each. ₩		



(2) PP Bottle
250ml × 4



(3) Microplate Rotor
(Option)



Optional Adapters (4ea/set)

Optional Rotors and Bottle

- DH.Cef250R.R1** (1) Stainless-steel Swing Rotor “**R1-250**” with Bucket(250ml × 4)
DH.Cef250R.B1 (2) PP(Polypropylene) Bottle with Screw Cap, 250ml for “**R1-250**”
DH.Cef250R.R2 (3) Aluminium Alloy Microplate Rotor(8 × 96well) “**R2-250**”, 4000rpm, 2200xg

Ea. Price, ₩

Optional Polypropylene Adapters

- DH.Cef250R.A1** PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ① (1.5ml × 9hole) × 4ea/set, id ϕ 11 × h37mm
DH.Cef250R.A2 PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ② (5ml × 7hole) × 4ea/set, id ϕ 14 × h47mm
DH.Cef250R.A3 PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ③ (5/7ml × 9hole) × 4ea/set, Vacu-tube, id ϕ 13 × h57mm
DH.Cef250R.A4 PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ④ (10ml × 7hole) × 4ea/set, Int. Thread, id ϕ 16 × h77mm
DH.Cef250R.A5 PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ⑤ (15ml × 5hole) × 4ea/set, id ϕ 17 × h102mm
DH.Cef250R.A6 PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ⑥ (20ml × 4hole) × 4ea/set, id ϕ 20 × h91mm
DH.Cef250R.A7 PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ⑦ (50ml × 1hole) × 4ea/set, id ϕ 30 × h91mm
DH.Cef250R.A8 PP(Polypropylene) Adapter for Rotor “**R1-250**”
; ⑧ (100ml × 1hole) × 4ea/set, id ϕ 39 × h98mm

Set. Price, ₩

DAIHAN™ 1.5~1,000ml Floor-standing Refrigerated High Volume Centrifuge “CEF-1000R”, -20℃ +40℃, 5,000rpm
Used for High Capacity Tubes & Bottle, without Rotor, 4,100×g, 대용량 냉장 원심분리기, 플라스틱 타입 Rotor 별도

CE

Max. 5,000 rpm




3-Year Warranty

“CEF-1000R” for 1.5~1,000ml Tubes/Bottles
※ Rotor is Optional !




(1) Angle Rotor “R1-1000R” 300ml × 6 (Option)
(2) Swing Rotor “R2-1000R” 500ml × 4 (Option)



PP Adapters for (2) Rotor “R2-1000R”



(3) Swing Rotor “R3-1000R” 1000ml × 4 (Option)
(4) Microplate Rotor “R4-1000R” (Option)



PP Adapters for (3) Rotor “R3-1000R”

Optional Adapters (4ea/set)

- Supplied by DAIHAN Korea ■ CE Certified
- 면역학, 수처리, 생화학, 약학, 정화 등 다양한 분야에서 널리 사용 가능.
- LCD 디스플레이, 사용자 친화적 인터페이스, 단순하고 편리한 조작.
- Ideal for Radioimmunoassay, Water Treatment, Biochemistry, Biology Pharmacy, Purification and etc.
- * Microcomputer Control, Electronic Lock and AC Variable-conversion Motor, Which delivers High Torque.
- * Touch Panel, Programmable Operation
- * Operation Data can be Set According to Demands and Stored Automatically.
- * LED Display, User-friendly Interface(The Real-time Conversation of Speed and RCF)
- * Simple and Convenient for Operation.
- * 10 Files of Acceleration and Deceleration Levels for Choice.

Specifications & Ordering information

<L1>

Spec.	Model	“CEF-1000R”
Max. RPM & RCF		5,000rpm, 4,100×g
Speed Accuracy		± 50rpm
Rotor Capacity		(1) 300ml × 6 Tubes (2) 500ml × 4 Buckets (3) 1,000ml × 4 Buckets (4) 8 × 96well Microplates
Timer		1min ~ 99h 59min
Temperature Setting		-20℃ ~ +40℃
Temperature Accuracy		± 1℃
Refrigerant Unit		Non-CFC R404a
Motor		AC Variable-conversion Motor
Noise Level		≤ 65dB
Safety Devices		Self-Locking, Over-speed Safety, Over-temperature Safety, Imbalance Safety, Automatic Alarm, etc.
Dimension & Weight		650×710×h1,100mm, 170kg
Power Supply		220V, 50/60Hz
Cat. No.		DH.Cef1000R
Price, Each. ₩		

Optional Rotors · Bottles · Adapters

Set. Price, ₩

- DH.Cef1000R.R1 (1) Aluminium Alloy Angle Rotor “R1-1000R”, for 300ml × 6, 5000rpm, 3850xg
DH.Cef1000R.A1 PP Adapter for “R1-1000R”, for 50ml × 6ea/set, φ 30 × h100mm
DH.Cef1000R.A2 , for 100ml × 6ea/set, φ 39 × h135mm
- DH.Cef1000R.R2 (2) STS Swing Rotor “R2-1000R” with Bucket(500ml × 4), 4000rpm, 4100xg
DH.Cef1000R.B1 PP Bottle for Rotor “R2-1000R” ; ① with Screw Cap, 500ml
DH.Cef1000R.A3 PP Adapter for Rotor “R2-1000R” ; ① (1.5ml × 20 hole) × 4ea/set
DH.Cef1000R.A4 ; ② (5/7ml × 20 hole) × 4ea/set, Vacu-tube
DH.Cef1000R.A5 ; ③ (10ml × 15 hole) × 4ea/set, Int. Thread
DH.Cef1000R.A6 ; ④ (15ml × 8 hole) × 4ea/set
DH.Cef1000R.A7 ; ⑤ (20ml × 9 hole) × 4ea/set
DH.Cef1000R.A8 ; ⑥ (50ml × 4 hole) × 4ea/set
DH.Cef1000R.A9 ; ⑦ (100ml × 2 hole) × 4ea/set
DH.Cef1000R.A10 ; ⑧ (250ml × 1 hole) × 4ea/set
- DH.Cef1000R.R3 (3) STS Swing Rotor “R3-1000R” with Bucket(1000ml × 4), 4000rpm, 4100xg
DH.Cef1000R.B2 PP Bottle for Rotor “R3-1000R” ; ① with Screw Cap, 1,000ml
DH.Cef1000R.A11 PP Adapter for Rotor “R3-1000R” ; ② (1.5ml × 24 hole) × 4ea/set
DH.Cef1000R.A12 ; ③ (5/7ml × 24 hole) × 4ea/set, Vacu-tube
DH.Cef1000R.A13 ; ④ (5/7ml × 37 hole) × 4ea/set, RIA
DH.Cef1000R.A14 ; ⑤ (10ml × 20 hole) × 4ea/set, Snap Cap
DH.Cef1000R.A15 ; ⑥ (15ml × 14 hole) × 4ea/set
DH.Cef1000R.A16 ; ⑦ (20ml × 12 hole) × 4ea/set
DH.Cef1000R.A17 ; ⑧ (50ml × 5 hole) × 4ea/set
DH.Cef1000R.A18 ; ⑨ (100ml × 3 hole) × 4ea/set
DH.Cef1000R.A19 ; ⑩ (250ml × 1 hole) × 4ea/set
DH.Cef1000R.A20 ; ⑪ (300ml × 1 hole) × 4ea/set
DH.Cef1000R.A21 ; ⑫ (500ml × 1 hole) × 4ea/set
- DH.Cef1000R.R4 (4) Aluminium Alloy Microplate Rotor(8 × 96well) “R4-1000R”, 4000rpm, 2200xg

※ Adapter Size (id φ × h) ;

- ③ φ 11 × h37mm, ④ φ 13 × h56mm, ⑤ φ 15.8 × h78mm, ⑥ φ 17 × h97mm, ⑦ φ 20.3 × h91mm,
⑧ φ 30 × h91mm, ⑨ φ 39 × h98mm, ⑩ φ 11 × h37mm, ⑪ φ 95 × h165mm, ⑫ φ 11 × h37mm,
⑬ φ 13 × h56mm, ⑭ φ 12.5 × h100mm, ⑮ φ 17 × h104mm, ⑯ φ 21.3 × h91mm,
⑰ φ 30 × h93mm, ⑱ φ 38 × h96mm, ⑲ φ 63 × h90mm, ⑳ φ 71 × h116mm, ㉑ φ 92 × h86mm

Common Use of Liquid Nitrogen Storage Tank & Canister-type CryoSystem™, 10~50 Lit

Without LN₂ Withdrawal Device & Canister(Vial Cane Holder), 액체질소 저장용 및 원통 캐니스터용 질소탱크 겸용, 펌프·캐니스터 별도판매



10L

20L

30L

50L


Incl. Roller Base



LN₂ Storage Tank



CryoSystem™
Equipped with
"Optional Canister"
for Vial Cane Holding



(1) Canister



(2) Dipper



(3) Roller Base ; Dolly



(4) Cane Cryovial Holder



(5) Cryo-Glove

■ Supplied by DAIHAN Korea

■ Low Daily Evaporation Rate of Liquid Nitrogen

■ Strong and Light-Weight Aluminum Construction

* Made of Antirust Aluminum Alloy

* Structured with High-vacuumed Super Adiabatic Multi-layers

* Equipped with Roller Base which is capable of achieving Easy Transportation (only for 50Lit)

* Easy to Lift and Pour

* with 2 Handles except for 10Lit

Specifications & Ordering information

<L1>

Model	"YDS-10D"	"YDS-20D"	"YDS-30D"	"YDS-50D"
Spec.	Max. Capa. 10L	Max. Capa. 20L	Max. Capa. 30L	Max. Capa. 50L
Volume	10L	20L	30L	50L
Neck Opening	φ50mm			
Static Evaporation Rate	0.15 Lit/day	0.17 Lit/day	0.17 Lit/day	0.24 Lit/day
Static Storing Period	65 days	120 days	174 days	212 days
Outer Diameter × Height	φ284 × h560 mm	φ345 × h655 mm	φ455 × h685 mm	φ455 × h850 mm
Empty Weight / Full Weight	6.2 / 13.9kg	8.2 / 23.9kg	14 / 37.7kg	18.1 / 58.0kg
Roller Base	-	-	-	Included
Optional Canister	Max Number	3ea	6ea	
	Max. No. of 1~2ml Vials	-		780ea
	Outer Diameter × Height	φ38 × h180 mm		
Optional Dipper	Overall Length	L490 mm	L535 mm	L705 mm
	Outer Diameter × Height	φ38 × h55 mm		
	Overall Length	L490 mm	L535 mm	L705 mm
Cat. No.	DH.Cry7010	DH.Cry7020	DH.Cry7030	DH.Cry7050
Price, Set. ₩				

Optional Accessories

(1) Canister, Stainless-steel

DH.Can7901	Canister, for "YDS-10D"
DH.Can7902	, for "YDS-20D"
DH.Can7903	, for "YDS-30D"
DH.Can7905	, for "YDS-50D"

Ea. Price, ₩

(2) Dipper, Stainless-steel NEW

DH.Dip7903.1	Dipper, for "YDS-10D"
DH.Dip7903.2	, for "YDS-20D"
DH.Dip7903.3	, for "YDS-30D"
DH.Dip7903.4	, for "YDS-50D"

Ea. Price, ₩

DH.Rol7922	(3) Roller Base, Dolly, OD φ360mm, w/4 Caster(2 Stop-on), for "YDS-30D"
DH.Rol7923	, for "YDS-50D"

Ea. Price, ₩

* When you order, please let us know the volume of your Tank.

SP.T313	(4) Cane, Aluminum, for 5 Vials (1.2~2ml) for "DH.Cry7050"
DH.Glo3021	(5) Glove, for Cryogenic Protection, Waterproof, Mid-arm, Free-size, L380mm

Pk. Price, ₩

12ea
1pair

LN₂ Withdrawal Device

* Ideal for Low Loss during the Liquid Discharging Process * Minimizes the Vaporization of the Liquid Nitrogen in the Container through a Rubber Ball which can be manually Compressed.

* Be fitted with Pouring Spouts or Dipper to aid in LN₂ Transfer. * Infusion Quantity ; ≥3 Lit/Min.

DH.Lwd7911 LN₂ Withdrawal Device, for 50mm Neck Opening, of 10 & 20 Lit.

Adjustable Length ; 425-565mm, Weight ; 2.4kg

DH.Lwd7912 LN₂ Withdrawal Device, for 50mm Neck Opening, of 30 & 50 Lit.

Adjustable Length ; 470-610mm, Weight ; 2.8kg

DH.Bul7920 Spare Bulb, for LN₂ Withdrawal Device



Ea. Price, ₩



Fig. Equipped with
LN₂ Withdrawal Device

Optional
LN₂ Withdrawal Device

Double Pipe
for Suction & Blow



LN₂ Withdrawal Components

- | | |
|---------------|------------------|
| ① Drainpipe | ⑦ Safety Valve |
| ② Rubber Ring | ⑧ Bulb |
| ③ Sealing Nut | ⑨ Plastic Pipe |
| ④ Locking Nut | ⑩ Air Filter |
| ⑤ Valve | ⑪ Pressure Gauge |
| ⑥ Vaporizer | |

Used for Storage of Cryovial Box with Rack, Storage Capacity ; 1050~3000vials(1.2 & 2.0ml), 랙 타입 크리오시스템/크리오 바이알 저장 액체 질소통, 시각 랙 별도판매



Rack-type CryoSystem™, 47~95Lit




Fig. Equipped with Racks(CryoBox Holder)
(Rack ; Optional)



(1) 6tier×25vials (2) 6tier×100vials (3) Dipper



(4) Roller Base ; Dolly



- Supplied by DAIHAN Korea**
■ Low Daily Evaporation Rate of Liquid Nitrogen **■ Strong and Light-Weight Aluminum Construction**
 * Made of Antirust Aluminum Alloy * Structured with High-vacuumed Super Adiabatic Multi-layers
 * Equipped with Roller Base which is capable of achieving Easy Transportation (only for 65 & 95Lit)

Specifications & Ordering information

<L1>

Spec.	Model	“YDS-47-127” Max. Capa. 47L	“YDS-65-216” Max. Capa. 65L	“YDS-95-216” Max. Capa. 95L
Volume		47L	65L	95L
Neck Opening		127 mm	216 mm	216 mm
Static Storing Period		105 days	73 days	100 days
Outer Diameter×Height		φ465×820mm	φ555×785mm	φ678×780mm
Empty Weight		18.6kg	32kg	33kg
Full Weight		56.5kg	84.52kg	110kg
Roller Base		-	Included	Included
Optional Rack	No. of Vial per Box	25	100	100
	No. of Boxes per Rack	6	6	6
	No. of Max Rack Equipped	6	3	5
	Vial Capacity (pcs)	1050	2400	2400
Optional Dipper	Outer Diameter×Height	φ65mm×h100mm		
	Overall Length	L680 mm	L635 mm	L620 mm
Cat. No.		DH.Cry7147	DH.Cry7165	DH.Cry7195
Price, Set. ₩				

Optional Accessories

DH.Rac7907	(1) Rack , Stainless-steel, Square, with 6tier×25vials/box					
DH.Rac7908	(2) , with 6tier×100vials/box, for "YDS-65"					
DH.Rac7909		, for "YDS-95"				
DH.Dip7904.1	(3) Dipper , Stainless-steel, φ65mm×h100mm, for "YDS-47" New					
DH.Dip7904.2		, for "YDS-65" New				
DH.Dip7904.3		, for "YDS-95" New				
DH.Rol7924	(4) Roller Base , Dolly, w/5 Caster(2 Stop-on), up to φ490mm, for "YDS-47"					
DH.Rol7925		, up to φ560mm, for "YDS-65"				
DH.Rol7926		, up to φ705mm, for "YDS-95"				
※ When you order, please let us know the volume of your container.						
DH.Glo3021	(5) Glove , for Cryogenic Protection, Waterproof, Mid-arm, Free-size, L380mm	<table><tr><th>Pk.</th><th>Price, ¥</th></tr><tr><td>1pair</td><td></td></tr></table>	Pk.	Price, ¥	1pair	
Pk.	Price, ¥					
1pair						
DH.Box901	(6) Cardboard Box , 81-Holes×φ13mm, with Cardboard Divider					
DH.Box902	(7) , with PP Divider					
SL.Rac3026	(8) PC Freezer Box , 25-holes×φ12.3mm, 76×76×h54 mm					
SL.Rac3031	(9) , 81-holes×φ13mm, 130×130×h54 mm, Blue					
SL.Rac3035	(10) , 100-holes×φ13mm, 132×132×h54 mm, White	<table><tr><th>Ea.</th><th>Price, ¥</th></tr><tr><td></td><td></td></tr></table>	Ea.	Price, ¥		
Ea.	Price, ¥					

3-Year Warranty

(1) Short-model



(2) Tall-model

With Digital Thermo-Hygrometer



(3) With Digital Thermo-Hygrometer



(4) With Urea Tray ; for Desiccant

DAIHAN™ 33Lit General Purpose PMMA Desiccator, Short- & Tall-form
 With Digital Thermo-Hygrometer · ABS Frame · Shelf · Desiccant Tray, 범용 PMMA 데시케이터

■ Supplied by DAIHAN Korea

■ Excellent Viewing & Air Tight Structure ; Heavy-duty ABS(Acrylonitrile Butadiene Styrene Resin) Frame, 3mm Clear PMMA(Polymethyl Methacrylate Resin) Body Panel & 5mm Shelf, and Rubber Sealing Door.

* Adjustable Shelf Quantity & Height / Distance.

* 2kg Shelf Load. * 30mm Shelf Moving Interval.

* Including 1 Urea Tray 190 × 273mm for Silica-gel / Desiccant.

※ 운임 별도 <L2>

(1) 33Lit Short-model "DeLL", 1-PMMA Shelf, with Digi. Thermo-Hygrometer**DH.DeLL**

General Desiccator, Short-model, 33Lit, "DeLL"

Outer ; 450 × 330 × h350mm, Inner ; 400 × 305 × h270mm

Ea. Price, ₩

(2) 33Lit Tall-model "DeLH", 3-PMMA Shelf, with Digi. Thermo-Hygrometer**DH.DeLH**

General Desiccator, Tall-model, 33Lit, "DeLH"

Outer ; 330 × 330 × h480mm, Inner ; 270 × 305 × h400mm

Ea. Price, ₩

DH.DeSH003

Spare, PMMA Shelf, for Short-model, 397 × 302 × thick 5mm, Load 2kg

DH.DeSH002

, for Tall-model, 267 × 302 × thick 4mm, Load 2kg

DH.DeDHY

, (3) Digital Thermo-Hygrometer, DC, 1~99℃, 11%~99%RH

DH.DeSH006

, (4) Urea Tray for Desiccant 273 × 190mm

Ea. Price, ₩

3-Year Warranty

(1) Bench-model



(2) Floor-model

With Digital Thermo-Hygrometer & Gas Valve

DAIHAN™ 164 · 207 · 320 · 640 · 604Lit Gas Exchangeable PMMA Desiccator, 1 & 4 Room
 With Digital Thermo-Hygrometer · Gas Valve · Stainless-steel Shelf, 가스치환 데시케이터

■ Supplied by DAIHAN Korea

■ Excellent Viewing & Air Tight Structure ; Heavy-duty Aluminium Frame, 5mm Clear PMMA(Polymethyl Methacrylate Resin) Body Panel, Stainless-steel Shelf, Caster, and Rubber Sealing Door.

* Ideal for Sample Storage / precise Instrument / Semiconductive Articles / Low Humidity, Characterization Test, etc.

* 20kg Shelf Load. * 50mm Shelf Moving Interval.

* Gas Exchange Valve; ϕ 6mm, One-touch fitting.

* Adjustable Shelf Quantity & Height / Distance.

* With Gas Valve / Gas Exchange / Stainless-steel Shelf / 1 Urea Tray 345 × 440mm for Desiccant / 1Digital Thermo-Hygrometer(But 4ea in 604Lit 4Room)

※ 운임 별도 <L2>

(1) 164Lit Bench Model "DeBG2", 1Room - 1Door, 2Shelf, 2Gas Valve**DH.DeBG2**

164Lit Desiccator, Gas-Ex. "DeBG2", 1Room, 2Shelf

Outer ; 574 × 535 × h765mm, Inner ; 490 × 460 × h730mm

Including 2 Shelf, 1 Silica gel Tray, and 2 Gas Valve

Ea. Price, ₩

(2) 207Lit Floor Model "DeBG1K", 1Room - 2Door, 4Shelf, 2Gas Valve, Caster**DH.DeBG1K**

207Lit Desiccator, Gas-Ex. "DeBG1K", 1Room - 2Door, 4Shelf

Outer ; 574 × 535 × h1070mm, Inner ; 490 × 460 × h920mm

Including 4 Shelf, 1 Silica gel Tray, and 2 Gas Valve

Ea. Price, ₩

(3) 320Lit Model "DeSS", 1Room - 2Door, 4Shelf, 2Gas Valve, Caster**DH.DeSS**

320Lit Desiccator, Gas-Ex. "DeSS", 1Room - 2Door

Outer ; 574 × 535 × h1770mm, Inner ; 490 × 460 × h1420mm

Including 4 Shelf, 1 Silica gel Tray, and 2 Gas Valve

Ea. Price, ₩

(4) 640Lit Model "DeWS", 1Room - 4Door, 8Shelf, 2Gas Valve, Caster**DH.DeWS**

640Lit Desiccator, Gas-Ex. "DeWS", 1Room - 4Door

Outer ; 1164 × 535 × h1770mm, Inner ; 980 × 460 × h1420mm

Including 8 Shelf, 2 Silica gel Tray, and 2 Gas Valve

Ea. Price, ₩

(5) 604Lit Model "DeW4S", 4Room - 4Door, 8Shelf, 8Gas Valve, Caster**DH.DeW4S**

604Lit Desiccator, Gas-Ex. "DeW4S", 4Room - 4Door

Outer ; 1164 × 535 × h1770mm, Inner ; 490 × 460 × h670mm, 4Room

Including 8 Shelf, 4 Silica gel Tray, and 8 Gas Valve

Ea. Price, ₩

Spare Accessory**DH.DeSH001**

Spare, Stainless-steel Shelf, 20kg Load, 490 × 460mm

DH.DeDHY

, (6) Digital Thermo-Hygrometer, DC, 1~99℃, 11%~99%RH

DH.DeSH005

, (7) Urea Tray for Desiccant 440 × 345mm

Ea. Price, ₩

3-Year Warranty

(3) 320Lit 1Room-2Door



(4) 640Lit 1Room-4Door

With Digital Thermo-Hygrometer & Gas Valve



(6) With Digital Thermo-Hygrometer



(7) With Urea Tray ; for Desiccant

3-Year Warranty



(1) Short Auto-Dry



(2) Tall Auto-Dry



(3) With Dehumidifier Unit



(4) With Digital Thermo-Hygrometer

※ Note ! ~25%RH ; If No Load Inside and Depends on Environment Conditions.

DAIHAN™ 33Lit Auto-Dry PMMA Desiccator, Short- & Tall-form, Dehumidifier ~25%RH With ABS Frame · Digital Thermo-Hygrometer, 자동 습도 조절 데시케이터

■ Supplied by DAIHAN Korea

■ Excellent Viewing & Air Tight Structure ; Heavy-duty ABS(Acrylonitrile Butadiene Styrene Resin) Frame, 3mm Clear PMMA(Polymethyl Methacrylate Resin) Body Panel & 5mm Shelf, and Rubber Sealing Door.

* Adjustable Shelf Quantity & Height / Distance.

* 2kg Shelf Load. * 30mm Shelf Moving Interval.

* Built-In Analog Dehumidifier Unit ; Repeated Regeneration of Solid Silica Gel Improves Dehumidification Capacity.

* Power Supply ; 220V, 50/60Hz, 3Watts.

※ 운임 별도 <L2>

(1) 33Lit Auto-Dry Short-model "DeADLL", 1-PMMA Shelf

DH.DeADLL Auto-Dry Desiccator, Short-model, 33Lit, "DeADLL"

Outer ; 450×330×h350mm, Inner ; 400×305×h270mm

Ea. Price, ₩

(2) 33Lit Auto-Dry Tall-model "DeADLH", 3-PMMA Shelf

DH.DeADLH Auto-Dry Desiccator, Tall-model, 33Lit, "DeADLH"

Outer ; 330×330×h480mm, Inner ; 270×305×h400mm

Ea. Price, ₩

DH.DeSH003 Spare, PMMA Shelf, for Short-model, 2kg Load, 397×302×thick 5mm
DH.DeSH002 , for Tall-model , 2kg Load, 267×302×thick 4mm

DH.DeHY , (3) Analog Dehumidifier Unit ; ~25%RH ; User Installable!

DH.DeDHY , (4) Digital Thermo-Hygrometer, DC, 1~99℃, 11%~99%RH

3-Year Warranty



(2) 320Lit
1Room-2Door



(3) 604Lit
1Room-4Door

Analog Control or Digital Display



With Analog Controlled Dehumidification Unit



(7) With Digital Thermo-Hygrometer

3-Year Warranty



(6) 604 Lit, 1Room-4Door

Fully Digital Control



With Digital Controlled Dehumidification Unit



With Digital Thermo Hygrometer

DAIHAN™ 207 · 320 · 604Lit Auto-Dry PMMA Desiccator, Analog or Digital-Control, 15~60%RH With Dehumidifier · Digital °C · %RH Meter · Al-Frame · 1Room, 자동 제습 데시케이터

■ Supplied by DAIHAN Korea

■ Excellent Auto-Dry / Viewing / Air Tight Structure ; Built-In Dehumidifier 15%~60%RH in Digital Control, Heavy-duty Aluminium Frame, 5mm Clear PMMA Body Panel, Stainless-steel Shelf 20kg Load, Caster, and Rubber Sealing Door.

* Analog or Digital Controlled Dehumidifier Unit ; Repeated Regeneration of Solid Silica Gel Improves Dehumidification Capacity. * Included 1-Dehumidifier.

* Ideal for Sample Storage / precise Instrument / Semiconductive Articles / Low Humidity, Characterization Test, etc.

* 20kg Shelf Load. * Including 4 Stainless-steel Shelf.

* Adjustable Shelf Quantity & Height / Distance.

* Display ; Digital Thermo-Hygrometer ; 1~99℃, 11%~99%RH

* Power Supply ; 220V, 50/60Hz

1) 3 Step(12~60%RH) Analog Controlled, 207 · 320 · 604 Lit

* 3 Step Analog Control ; 1 Step. 12~20%RH, 2 Step. 20~40%RH, 3 Step. 40~60%RH

※ 운임 별도 <L2>

Ea. Price, ₩

DH.DeADABG1K (1) 207Lit Analog Auto-Dry Desiccator "DeADABG1K"
1Room - 1Door, 4Shelf
Outer ; 574×535×h1070mm, Inner ; 490×460×h920mm

DH.DeADASS (2) 320Lit Analog Auto-Dry Desiccator "DeADASS"
1Room - 2Door, 4Shelf
Outer ; 574×535×h1770mm, Inner ; 490×460×h1420mm

DH.DeADAWS (3) 604Lit Analog Auto-Dry Desiccator "DeADAWS"
1Room - 4Door, 8Shelf
Outer ; 1154×535×h1770mm, Inner ; 980×460×h1420mm

2) 15~60%RH Digital Hygro-Controlled, 207 · 320 · 604 Lit

Ea. Price, ₩

DH.DeADDBG1K (4) 207Lit Digi-Control Auto-Dry Desiccator "DeADDBG1K"
1Room - 1Door, 4Shelf
Outer ; 574×535×h1070mm, Inner ; 490×460×h920mm

DH.DeADDSS (5) 320Lit Digi-Control Auto-Dry Desiccator "DeADDSS"
1Room - 2Door, 4Shelf
Outer ; 574×535×h1770mm, Inner ; 490×460×h1420m

DH.DeADDWS (6) 604Lit Digi-Control Auto-Dry Desiccator "DeADDWS"
1Room - 4Door, 8Shelf
Outer ; 1154×535×h1770mm, Inner ; 980×460×h1420mm

DH.DeSH001 Spare, Stainless-steel Shelf, 20kg Load, 490×460mm
DH.DeDHY (7) Digital Thermo-Hygrometer, DC, 1~99℃, 11%~99%RH
for (1) "DeADABG1K", (2) "DeADASS", (3) "DeADAWS"

3-Year Warranty

(1) Short Auto-Dry



(2) Tall Auto-Dry



(3) With Dehumidifier Unit



(4) With Digital Thermo-Hygrometer

※ Note ! ~25%RH ; If No Load Inside and Depends on Environment Conditions.

DAIHAN™ 33Lit UV Protected Auto-Dry PMMA Desiccator, Short- & Tall-model, ~25%RH With Dehumidifier · Digital Thermo-Hygrometer · ABS Frame, UV차단 자동 습도 조절 데시케이터

■ Supplied by DAIHAN Korea

■ Good Viewing / Excellent Auto-Dry / Air Tight Structure ; Heavy-duty ABS Frame, 3mm Transparent Amber PMMA Body Panel & 5mm Shelf, ~25% Dehumidifier Unit, and Rubber Sealing Door.

■ Made of Transparent Amber PMMA Panel ; 520nm (or Less) UV Protection. * 2kg Shelf Load.

* Analog Controlled Dehumidifier Unit ; Repeated Regeneration of Solid Silica Gel Improves Dehumidification Capacity.

* Adjustable Shelf Quantity & Height / Distance. * Power Supply ; 220V, 50/60Hz, 3Watts.

※ 운임 별도 <L2>

(1) 33Lit Short-model "DeADLL.UV", UV Protect Auto-Dry, 1-PMMA Shelf

DH.DeADLL.UV UV Protect Auto-Dry Desiccator "DeADLL.UV", Short-model, 33Lit Outer ; 450×330×h350mm, Inner ; 400×305×h270mm

(2) 33Lit Tall-model "DeADLH.UV", UV Protect Auto-Dry, 3-PMMA Shelf

DH.DeADLH.UV UV Protect Auto-Dry Desiccator "DeADLH.UV", Tall-model, 33Lit Outer ; 330×330×h480mm, Inner ; 270×305×h400mm

DH.DeSH003 Spare, PMMA Shelf, for Short-model, 2kg Load, 397×302×thick 5mm

DH.DeSH002 , for Tall-model , 2kg Load, 267×302×thick 4mm

DH.DeHY , (3) Analog Dehumidifier Unit ; ~25%RH

DH.DeDHY , (4) Digital Thermo-Hygrometer, DC, 1~99℃, 11%~99%RH

Ea. Price, ₩

Ea. Price, ₩

3-Year Warranty

(1) 7Lit



(2) 9.2Lit

1) Upper Door Model



(3) 12Lit



(4) 30Lit

2) Front Door Model



(5) 48Lit

2) Front Door Model



Option ; Gas Flowmeter



Spare Digital Thermo-Hygrometer

DAIHAN™ 7~48Lit Vacuum PMMA Desiccator, Clear, with Press-Gauge With Digital Thermo-Hygrometer, Approx - 1 Torr/133Pa, 진공 데시케이터

■ Supplied by DAIHAN Korea

■ Permissible Degree of Vacuum ; Approx -1 Torr / 133Pa

■ With 20~30mm Clear PMMA Panel, 2 Vacuum Valve, 1 Vacuum Gauge, Silicone

& Rubber Sealing Door, Shelf, and 1 Desiccant Tray(270×144×h16mm for 7 / 9.2 / 12 Lit)

* Two Vacuum Valve ϕ 8mm

* 1 Pressure Dial Gauge ; 0Bar~1Bar $\pm$ 1.5%

※ 운임 별도 <L2>

1) 7 & 9.2Lit Upper Door, Short-form, 20mm PMMA Panel

DH.DeVS (1) 7Lit Vacu. Desiccator "DeVS", 1Shelf, 1Tray Outer ; 300×300×h170mm, Inner ; 260×260×h100mm

DH.DeVM (2) 9.2Lit Vacu. Desiccator "DeVM", 1Shelf, 1Tray Outer ; 260×260×h360mm, Inner ; 210×210×h210mm

2) 12 · 30 · 48Lit Front Door, Tall-form, 20~30mm PMMA Panel

DH.DeVL (3) 12Lit Vacu. Desiccator "DeVL", 20mm Panel, 2Shelf, 1Tray Outer ; 310×225×h430mm, Inner ; 260×180×h260mm

DH.DeVLH (4) 30Lit Vacu. Desiccator "DeVLH", 25mm Panel, 4Shelf, 1Tray(190×273×h18mm) Outer ; 320×335×h600mm, Inner ; 260×280×h410mm

DH.DeVW (5) 48Lit Vacu. Desiccator "DeVW", 30mm Panel 3Shelf, 1Tray(313×230×h20mm) Outer ; 470×365×h600mm, Inner ; 400×300×h400mm

Ea. Price, ₩

Ea. Price, ₩

Optional Gas Flowmeter, for Gas-Ex. Desiccator, 10ℓ/min

DH.DeFL Gas Flowmeter, with Bracket & 2 Fitting for ϕ 6mm Tube, 70×L170mm

Ea. Price, ₩

Spare Digital Thermo-Hygrometer

DH.DeDHY Built in Digital Thermo-Hygrometer, DC, 1~99℃, 11~99%RH

Ea. Price, ₩

Desiccant Silica-Gel, Indicating-type Blue & Non-Indicating-type White, 실리카겔 건조제

(1) Indicating-type, Blue

* When Absorption of Moisture the Color Changing Blue → White

* Reusable ; Revivified to Blue By the Drying of the White Silicagel

DH.Sil1001 Silica Gel, Blue, CP, 5~10mesh, 500g/PP Bottle

DH.Sil1002 , 1kg/PP Bottle

(2) Non-Indicating-type, White

S19.320 Silica Gel, Non-woven Fabric, 70×85mm, 20g

DH.Sil1005 , White, CP, 5~10mesh, 500g/PE Bottle

Pk.	Price, ₩	Cs.	Price, ₩
500g		20pk	
1kg		10pk	

Pk.	Price, ₩	Cs.	Price, ₩
10ea		100pk	
500g			



(1) Indicating-type, Blue

1kg



(2) Non-Indicating-type, White

Articles

Cat. No

Description

Application / Chemical Compatibility List for DAIHAN Wisdose™ Dispenser Series

Code explanations

• Regular = Wisdose™ Standard • HF = Wisdose™ HF

A = Good resistance B = Acceptable with limitations C = Not recommended

1 = Acid vapours (better resistance with lower concentration). Rinse the instrument in the rinse mode otherwise do not leave instrument on bottle.

2 = Risk of damage, softening or discoloration of external parts through vapours. Rinse the instrument in the rinse mode otherwise do not leave instrument on bottle.

3 = Chemical degradation of glass parts (plunger/barrel).

Reagent	Regular	HF	Reagent	Regular	HF	Reagent	Regular	HF	Reagent	Regular	HF
Acetaldehyde (Ethanal)	A	A	Chloroethanol	B/2	B/2	Hydrobromic acid	A	A	Potassium chloride	A	A
Acetic acid 96%	A	A	Chloroform (Trichloromethane)	B/2	B/2	Hydrochloric acid 20% (HCl)	A	A	Potassium dichromate	A	A
Acetic acid 100% (glacial)	B/2	B/2	Nitro-hydrochloric acid (Aqua regia)	B/2	A	Hydrochloric acid 37% (HCl)	B/1	B/1	Potassium dihydrogen phosphate	A	A
Acetic anhydride	B/2	B/2	Chlorosulfonic acid	B/2	B/2	Hydrofluoric acid (HF)	C/3	A	Potassium hydroxide	A	C/3
Acetone (Propanone)	B/2	B/2	Chlorosulfuric acid 100%	B/1/2	B/1/2	Hydrogen peroxide	A	A	Potassium iodide	A	C/3
Acetonitrile (MECN)	A	A	Chromic acid 100%	B/1/2	B/1/2	Iodine	A	A	Potassium permanganate (persulfate)	A	A
Acetophenone	B/2	B/2	Chromosulfuric acid 100%	C/1/2	C/1/2	Iodine bromide	C/2	C/2	Potassium peroxydisulfate	A	A
Acetyl Chloride	B/2	B/2	Citric acid	A	A	Iodine chloride	C/2	C/2	Potassium sulfate	A	A
Acetylacetone	A	A	Copper fluoride	A	A	Isoamyl alcohol	A	A	Propionic acid (Propanoic acid)	A	A
Acrylic acid	A	A	Copper sulfate	A	A	Isobutanol	A	A	Propylene glycol (Propane-1,2-diol)	A	A
Acrylonitrile	B/2	B/2	Covi-Ox-T70 / Mixed Tocopherol	A	A	Isooctane	A	A	Propylene oxide	A	A
Adipic acid	A	A	Cresol	A	A	Isopropanol	A	A	Picric acid (Trinitrophenol)	B/2	B/2
Allyl alcohol	A	A	Cumene (Isopropylbenzene)	B/2	B/2	Isopropyl ether	B/2	B/2	Pyridine	B/2	B/2
Aluminum chloride	A	A	Cyanoacrylate	A	A	Iso-propylamine	B/2	B/2	Pyruvic acid	A	A
Amino acids	A	A	Cyclohexane	B/2	B/2	Lactic acid	A	A	Resorcin	A	A
Ammonia 20%	B/2	B/2	Cyclohexanone	B/2	B/2	Kerosene	A	*	Salicylaldehyde	A	A
Ammonia 20~30%	B/2	B/2	Cyclopentane	B/2	B/2	2-Methoxyethanol	A	A	Scintillation fluid	A	A
Ammonium chloride	A	A	1,2-Diethylbenzene	B/2	B/2	Methanol	A	A	Silver acetate	A	A
Ammonium fluoride	A	A	1,4-Dioxane (Diethylene dioxide)	B/2	B/2	Methoxybenzene (Anisol)	B/2	B/2	Silane (Silicon hydride)	*	B/1
Ammonium molybdate	A	A	1-Decanol	A	A	Methyl benzoate	B/2	B/2	Silver nitrate	A	A
Ammonium sulfate	A	A	Decane	A	A	Methyl chloride (Chloromethane)	B/2	B/2	Sodium acetate	A	A
Amyl alcohol (Pentanol)	A	A	Di-(2-ethylhexyl) peroxydicarbonate	B/2	B/2	Methyl ethyl ketone (MEK/Butanone)	B/2	B/2	Sodium chloride (kitchen salt)	A	A
Amyl chloride (Chloropentane)	B/2	B/2	Dibenzyl ether	B/2	B/2	Methyl formate	A	A	Sodium dichromate	A	A
Aniline	A	A	Dichloroacetic acid	A	A	Methyl iodide (Iodomethane)	B/2	B/2	Sodium fluoride	A	A
Ascorbic acid	A	A	Dichlorobenzene	A	A	Methyl methacrylate (MMA)	B/2	B/2	Sodium hydroxide 30%	A	A
n-Amyl acetate	B/2	B/2	Dichloroethane	A	A	Methyl propyl ketone (2-Pentanone)	A	A	Sodium hypochlorite	A	A
Barium chloride	A	A	Dichloroethylene	B/2	B/2	Methyl tert-butyl ether	B/2	B/2	Sodium thiosulfate	A	A
Benzaldehyde	A	A	Diesel oil (Heating oil)	A	A	Methylene chloride (Dichloromethane) (DCM)	B/2	B/2	Sulfonitric acid 100%	B/2	B/2
Benzene	B/2	B/2	Diethanolamine	A	A	Methylpentanone	A	A	Sulfur dioxide	B/2	B/2
Benzine	A	A	Diethylamine	B/2	B/2	Mineral oil (Engine oil)	A	A	Sulfuric acid 100%	B/2	B/2
Benzoyl chloride	B/2	B/2	Diethylene glycol	A	A	Monochloroacetic acid	A	A	Sulfuric acid <10%	A	A
Benzyl alcohol	A	A	Diethylether	B/2	B/2	Nitric acid 100%	A	A	Sulfuric acid (10-75%)	B/1	B/1
Benzyl chloride	B/2	B/2	Dimethylacetamide	A	A	Nitric acid 30-70%	A	A	Sulfuric acid (Cold conc.)	A	A
Bis(2-ethylhexyl) phthalate	B/2	B/2	Dimethyl sulfoxide (DMSO)	B/2	B/2	Nitric acid dil. <30%	A	A	Sulfuric acid (Hot conc.)	B/2	B/2
Boric acid 10%	A	A	Dimethylaniline	A	A	Nitrobenzene	B/2	B/2	1,1,2-Trichlorotrifluoroethane	B/2	B/2
Bromine	C/2	C/2	Dimethylformamide (DMF)	B/2	B/2	Nitromethane	B/2	B/2	Tartaric acid	A	A
Bromobenzene	B/2	B/2	Ethanol	A	A	N-methyl-2-pyrrolidone (NMP)	A	A	Tetrachlorethylene	B/2	B/2
Bromonaphthalene	A	A	Ethanolamine	B/2	B/2	Octane	A	A	Tetrahydrofuran (THF)	B/2	B/2
Butanediol	A	A	Ether	B/2	B/2	Octanol	A	A	Tetramethylammonium hydroxide	A	A
Buthanol	A	A	Ethyl acetate	B/2	B/2	Oil (vegetable, animal)	B/2	B/2	Toluene	B/2	B/2
Butyl acetate	B/2	B/2	Ethylbenzene	B/2	B/2	Oil of turpentine	B/2	B/2	Trichlorethylene	B/2	B/2
Butyl methyl ether	B/2	B/2	Ethylene chloride	B/2	B/2	Oleic acid	A	A	Trichloroacetic acid	B/2	B/2
Butylamine	B/2	B/2	Ethylene diamine	A	A	Oleum (Fuming Sulfuric acid)	A	A	Trichlorobenzene	B/2	B/2
Butyric acid	B/2	B/2	Ethylene glycol	A	A	Oxalic acid	A	A	Trichloroethane	B/2	B/2
Calcium carbonate	A	A	Fluoroacetic acid	B/2	B/2	Pentane	B/2	B/2	Tricresyl phosphate	*	C/3
Calcium chloride	A	A	Formaldehyde (Formalin)	A	A	Peracetic acid	A	A	Triethanolamine	A	A
Calcium hydroxide	A	A	Formamide	A	A	Perchloric acid 100%	B/2	B/2	Triethylamine	A	A
Calcium hypochlorite	A	A	Formic acid	A	A	Perchloric acid diluted	A	A	Triethylene glycol	A	A
Carbon disulfide	B/2	B/2	Gamma-butyrolactone	A	A	Perchloroethylene	B/2	B/2	Trifluoroacetic anhydride (TFAA)	B/2	B/2
Carbon tetrachloride	B/2	B/2	Gasoline	B/2	B/2	Petroleum	B/2	B/2	Trifluoromethane (Fluoroform)	B/2	B/2
Chlorine dioxide	B/2	B/2	Glycerin <40%	A	A	Petroleum ether / spirit	B/2	B/2	Urea	A	A
Chlorine water	B/2	B/2	Glycolic acid 50%	A	A	Phenol	A	A	Uric acid	*	C/3
Chloro naphthalene	B/2	B/2	Heating oil (Diesel oil)	A	A	Phenylethanol	B/2	B/2	Xylene	B/2	B/2
Chloroacetaldehyde 45%	A	A	Heptane	A	A	Phenylhydrazine	B/2	B/2	Zinc chloride 10%	A	C/3
Chloroacetic acid	A	A	Hexane	A	A	Phosphoric acid 100%	A	A	Zinc sulfate 10%	A	C/3
Chloroacetone	B/2	B/2	Hexanoic acid	A	A	Phosphoric acid 85%	A	A			
Chlorobenzene	B/2	B/2	Hexanol	A	A	Piperidine	B/2	B/2			
Chlorobutane	B/2	B/2	Hydriodic acid	B/2	B/2	Potassium carbonate	*	C/3			

* Store the instrument and accessories only in clean conditions in a cool and dry place.

Storage temperature ; from -20°C to +50°C (from -4 °F to 122 °F)

DAIHAN Wisdoser™ 0.25~100ml Standard & HF Bottle Top Dispenser Set, with Adjustable Intake Tube & Delivery Nozzle
 With Re-Circulation Valve and Calibration Report, Fully Autoclavable, CE · ISO · DAKKS · IAF Certified, with GL45 Glass Bottle
 “Standard” & “HF”(for Hydrofluoric Acid)-models, 편리형 디스펜서, 일반형 & 불산형, 병 포함

- Supplied by DAIHAN Korea ■ With Individual Calibration Report
- “Re-Circulation” 밸브를 통해 공기가 배기되는 동안 시약의 손실을 방지함 ■ 스프링 없는 PTFE밸브는 높은 호환성 및 밸브의 보존기간 보장
- 불산형 디스펜서는 Fluoropolymer Barrel로 제작되어 불산 사용에 적합함 ■ Specially Designed Volume Adjustment Knob.
- Springless Valve with High Purity Ceramic Components Facilitate Smooth and Jam-Free Functioning.
- Certified ISO 9001, CE, DAKKS(German Accreditation Authority) and IAF(International Accreditation Forum).
- Assured Quality with ISO 8655 Conformed Calibration.



- * Easy Calibration * Specially Designed Calibration Tool is provided for Convenient and Quick in-lab user Re-calibration.
- * PTFE Piston with ETP O-Ring Ensures Very High Chemical Compatibility& Wiper Like Smooth, Effortless Piston Movement.
- * The Error Limits are well within the Limits of DIN EN ISO 8655

- ▶ Complete Set with ; 1) Dispenser, 2) Adjustable Intake Tube, 3) Calibration Tool, 4) Operation Manual
 5) Calibration Report, 6) PP GL Thread Adapter Set / 6ea of GL-28 · 30 · 32 · 38 · 40 · 45,
 7) GL45 Glass Bottle ①1,000ml 1ea for 2.5~10ml dispensers or ②2,000ml 1ea for 25~100ml dispensers

※ GL 38 Adapter is compatible with US ASTM Thread (38-400) like 4Lit Bottle from “J.T.Baker” or “Wheaton”



Incl. Calibration Report



Complete Kit



(3) GL Adapter Set, PP



(4) Adjustable Intake Tube



(5) Delivery Nozzle

(1) Wisdoser™ Standard Dispenser Kit, 0.25~100ml, with GL45 Glass Bottle, 병 포함

DH.Dos301	Dispenser, Bottle Top, Standard, 0.25~2.5 / 0.05ml, accu. ± 0.015ml, cv.0.005ml
DH.Dos302	, 0.5~5 / 0.1ml, , ± 0.030ml, .0.010ml
DH.Dos303	, 1~10 / 0.2ml, , ± 0.060ml, .0.020ml
DH.Dos304	, 2.5~30 / 0.5ml, , ± 0.180ml, .0.060ml
DH.Dos305	, 5~60 / 1ml, , ± 0.360ml, .0.120ml
DH.Dos306	, 10~100 / 2ml, , ± 0.600ml, .0.200ml

(2) Wisdoser™ HF(Hydrofluoric) Acid Dispenser Kit, 1~30ml, with Bottle, 병 포함

* Excellent Compatibility with Hydrofluoric Acid, High Purity Media, & Acids

DH.Dos311	Dispenser, Bottle Top, HF, 1.0~10.0 / 0.2ml, accu. ± 0.060ml, cv.0.020ml
DH.Dos312	, 2.5~30.0 / 0.5ml, accu. ± 0.180ml, .0.060ml

(3) Spare GL Adapter Set, PP

* Adapter Thread ; GL-28, 32, 38, 40, and 45

DH.Dos001	Adapter Set of 5pcs, Grey, for Basis Model
-----------	--

(4) Spare Adjustable Intake Telescopic Tube for All Size

DH.Dos003	Adjustable Intake Tube for Dispenser (all sizes)
-----------	--

(5) Spare Delivery Nozzle, with PTFE Screw Adaptor and FEP Tubing

DH.Dos004	Delivery Nozzle, for 2.5 · 5 · 10ml
DH.Dos005	, for 30 · 60 · 100ml

<L1>

Ea.	Price, ₩

Set.	Price, ₩

Set.	Price, ₩
5ea	

Ea.	Price, ₩

Ea.	Price, ₩

DAIHAN™ 105 & 155 Lit Multi-purpose Drying Chamber, I.Gravity-Air "DRYG" & II.Forced-Air "DRYF", 50℃~90℃, ±1.0℃
With Digital PID Control, FND-Display, Viewing Window, and 2 Stainless-steel Wire Shelves, 다용도 드라이 챔버, 자연대류식 & 강제순환식



3-Year Warranty

I. 105Lit. Gravity-Air type



3-Year Warranty

II. 155Lit. Forced-Air Type, with Casters



Stainless-steel Wire Shelf
(Included)

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Optimized Air Flow by Forced-Air & Gravity-Air Convection Mechanism
- * Suitable for Drying, Baking, Conditioning, Curing, Pre-Heating and Aging
- * Digital PID Control System Implementing Superior Temperature Accuracy
- * Seven-segment FND-display
- * Compact Design
- * Ambient to Max 90℃ Temp Range, Accuracy of ±1.0℃
- * Corrosion Resistant 304 Stainless-steel Chamber and Shelves
- * Perfect Heat Insulation made by Heat-Resistant Rubber Packing and Air Flow Layer
- * Over Temp & Over Current Protection and Sensor Error Detection
- * 2 Stainless-steel Shelves are Included

Specifications & Ordering information

<L1>

Spec.	Model	I. Gravity-Air type		II. Forced-Air type	
		(1) "DRYG-105"	(2) "DRYG-155"	(3) "DRYF-105"	(4) "DRYF-155"
Air Flow Fan		None		Yes	
Capacity		105 Lit	155 Lit	105 Lit	155 Lit
Temp	Range & Sensor	50℃~90℃, PT100 Sensor			
	Accu. & Uniform.	Accuracy ±1.0℃ at 90℃, Uniformity ±1.5℃ at 90℃			
Heater Power		1,100 W	1,200 W	1,150 W	1,300 W
Heat-up Time		30min up to 90℃			
Timer & Alarm		99hr 59min (with Delay / Continuous Run.), Error Status & Timer-end			
Display		Seven segments FND			
Controller		Digital PID Controller			
Resolution		Control ; ±1.0℃, Display ; ±1.0℃			
Shelf (Included)		2 Stainless-steel Wire Shelf			
Load per Shelf		16 kg			
Movable Caster		none	4ea	none	4ea
Permitted Total Load		32 kg (16 kg × 2 Shelves)			
Ventilation		1 of Ventilation Hole (id. ϕ40mm) with Stainless Steel Cap			
Safety Circuit		Over Temp, and Over Current Protector, Sensor Error Detector			
Material	Internal	Stainless-steel (#304)			
	External	Powder Coated steel			
	Insulation	Glass Wool			
Dimension (w×d×h)	Interior (mm)	485×410×535	550×470×600	485×410×535	550×470×600
	Exterior (mm)	653×739×783	718×804×848	653×739×783	718×804×848
Net Weight		60 kg	71 kg	62 kg	74 kg
Packing Size (w×d×h)(mm)		753×818×934	818×878×995	753×818×934	818×878×995
Gross Weight		75 kg	89 kg	77 kg	92 kg
Power Consumption		1.1 Kw	1.2 Kw	1.2 Kw	1.3 Kw
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz			
Cat. No.	230V	DH.DryG0105	DH.DryG0155	DH.DryF0105	DH.DryF0155
	120V	DH.DryG1105	DH.DryG1155	DH.DryF1105	DH.DryF1155
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Stainless-steel Wire Shelves


DH.WOF11105 Wire Shelf, for 105 Lit, "OFS105", w455×d375mm
DH.WOF11155 , for 155 Lit, "OFS155", w520×d449mm

Ea. Price, ₩

DAIHAN™ 39 places Labware Dryer “DRYL-39”, with Heating Blower & Air Filter, 230V, 1.5kW, Max 50℃ & 75℃
Ideal for Glassware · Plasticware · Small Articles, Timer Set 1~120min, Continuous 8 hours, 실험기구 건열풍 건조기



3-Year Warranty



Assorted Drying Pegs (Included)

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Rating 1 MERV Fabric Pre-Filter ; 3~10 μ m less than 20%
- * 1 Pre-Filter Installed, 2 spare Pre-Filters Included
- * High Efficiency & Long Life!
- * Made of Galvanized Steel with Powder Coating ; All Body
- * Free-Standing ; Space Saving Design.
- * a Peg has Several Blown Air Holes.
- * Blown-over Tray on the Underside of Dryer.

Specifications & Ordering Information

<L1>

Spec.	Model	“DRYL-39”
Operation mode & Filtered Heating Air Temp		Glass Mode ; for Glassware Drying, max. 75℃ $\pm$ 5℃ Plastic Mode ; for Plasticware Drying, max. 50℃ $\pm$ 5℃ Ambient mode ; Ambient Temp
Timer		Countdown Timer Interval Setting ; 1~120min Continuous Run
39 Drying Pegs		① 7 Pegs - ϕ 12.7 $\times$ L152mm ② 6 Pegs - ϕ 12.7 $\times$ L142mm ③ 13 Pegs - ϕ 8 $\times$ L117mm ④ 13 Pegs - ϕ 6 $\times$ L79mm Stainless-steel
Safety		Over Heat Protection
Ambient Temp & Humidity		5℃~45℃, up to 90% RH
Dimension & Weight (w $\times$ d $\times$ h)		550 $\times$ 350 $\times$ h515mm, 17kg
Heater power		1,500W
Power Consumption		Max. 1,600W
Electric Reqts		230V, 60Hz
Cat. No.		DH.Dry8013
Price, Each. ₩		

※ 운임 별도

DAIHAN™ Portable Gel Documentation System “WGD-20”, 6.3M Pixel CMOS Camera, with Certi. & Traceability

With Optical Lens · Portable Darkroom Hood · Capture Software · USB Interface

젤 이미지 분석기, 고성능/최저가, 디지털 CMOS 카메라(630만 화소) 장착



3-Year Warranty

Portable Gel Documentation System “WGD-20”

※ Computer System Requirements

- Intel Chipset, USB 3.0 Support Mainboard
- Over Pentium 4 Processor
- Over 256MB Main Memory
- Over 60MB Space of HDD for Installation and Running
- Window 7, 8, 10, 11



Imaging Capture S/W

; Included S/W, Real-Time Image Viewing/Capture, Manipulate and Save Images



(1) UV Transilluminator(Optional)

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- The Most Cost-Effective System
- High Precise and Large Caliber Optical Lens with Zoom, Focus and Iris Control are Used
- * DNA, RNA gel Imaging Capture and Analysis
- * Digital CMOS Camera (Resolution 6.3 Million pixel) with Portable Darkroom is Offered
- * 3,072 × 2,048 Resolution, 14-bit CMOS Sensor
- * Available for 14-bit Gray Level Image Capturing
- * Power Consumption ; Supplied USB(USB 3.0)
- * Very Compact Portable Darkroom System Enables the Existing Transilluminators in the Lab to be used together with this system. High Flexibility in use is Offered
- * USB 3.0 Interface is available, through which Real-Time Images are transferred directly to PC
- * Software ; Easy and Fast Capture of Gel Image
- * TotalLab™ Software ; Fully Automatic, Single Button Press Complete Image Analysis in Interesting Area(Optional)

Specifications & Ordering information

<L1>

Spec.		Model	“WGD-20”	“WGD-20S” Includes DNA Analysis S/W
CCD Module	CCD & Pixel Size		1/8" CMOS, 2.4 × 2.4 um	
	Max. Resolution		3,072 × 2,048	
	Lens Mount & Optical Size		C mount, 7.4 × 4.9 mm	
	Connect to PC & Power		USB 3.0, Supplied USB	
Zoom Lens	Iris & Focal Length		F1.4 - Close, 11.5~69mm	
	Ambient Temp		-10℃ to +50℃	
Overall	Dimension & Weight		300 × 260 × h470 mm, 9 kg	
	Software Capability		Windows 7, 8, 10, 11	
Standard S/W			Imaging Capture Software Included	
DNA Analysis S/W			none	TotalLab™ CLIQS Included
Packing Size & Weight			400 × 350 × h510 mm, 18 kg	
Cat. No.			DH.WGD00020	DH.WGD00024
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Optional Transilluminator & S/W

DH.WUV00020	(1) Transilluminator, “WUV-L20”, 20 × 20cm Filter 312nm, 15W UV Lamp × 6ea
DH.WUV00050	(1) Transilluminator, “WUV-L50”, 20 × 20cm Filter 365nm + 312nm, 15W UV Lamp × 6ea
DH.WUV00120	(1) Transilluminator, “WUV-M20”, 20 × 20cm Filter 312nm, 8W UV Lamp × 6ea
DH.WGD00100	(2) DNA/Protein Analysis S/W, TotalLab™ CLIQS
DH.WGD00121	(3) , CLIQS 1D Pro
DH.WGD00200	LapTop Computer, WGD-Laptop, for WGD Gel Doc.System

Ea. Price, ₩


DNA Analysis S/W

(2) TotalLab™ CLIQS

- Automatic and Manual Lane Detection, Background Subtraction
- Fully Automatic Band Detection, Molecular Size / pI Calibration
- Quantity Calibration, Normalization

(3) CLIQS 1D Pro (additional function)

- Match Band Data Patterns across Large Groups of Lanes
- Choice of Algorithms to Generate Dendrograms ; NeighborJoining, UPGMA, Complete Linkage, Single Linkage, Ward Library, Queries Function



One Stop Analysis

DAIHAN™ Gel Imaging System “GIS-2”, 6.3M Pixel, 2.4 μ m, F1.4, 11~69mm Manual Zoom

With Capture Software, High Quality CMOS Camera & Optical Lens, Auto Control System, USB Connection, Certi. & Traceability

젤 이미지 분석기, 자외선 발산기 일체형, 디지털 CMOS 카메라(630만 화소) 장착



3-Year Warranty

“GIS-2” Gel Imaging System,
Built-in Transilluminator



Imaging Capture S/W

; Included S/W, Real-Time Image Viewing/Capture,
Manipulate and Save Images



Sturdy
Camera Bracket

F1.4, 11~69mm
Manual Len

Emission Filter

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.
- Complete System with 6.3M Pixel Digital CMOS Camera, Safety Darkroom Hood and Built-in Transilluminator(UV Type or LED Type)
- High-quality F1.4 Optical Lens with 11~69mm Manual Zoom, Focus and Iris Control is Used.
- Included Imaging Capture Software.

- * Ideal for Various Applications of DNA, RNA Protein Gel Imaging Capture and Analysis.
- * Unified Darkroom System with Built-in Sliding Transilluminator offers very safety and convenience in use
- * Built-in Single Wave Transilluminator-UV Type ; 312nm, 6 $\times$ 8W UV Lamps and ETI White LED $\times$ 3Ea -LED Type ; 470nm, Blue LED $\times$ 324Ea, 35W
- * 3072 $\times$ 2048 CMOS Resolution * Can be Used with Tablet PC(Optional)
- * USB 2.0 Interface is Available through Which Real-time Images are Transferred Directly to PC

Specifications & Ordering information

<L1>

Spec.	Model	“GIS-2”	
Transilluminator Type		Standard UV Type	Blue LED Type
Camera Resolution		6.3 Mega Pixel CMOS	
Pixel Size		2.4 μ m	
Lens		F1.4 Fixed, 11~69mm Manual Zoom	
Transilluminator		312nm, 6 $\times$ 8W UV Lamps, Filter Size ; 210 $\times$ 260mm	470nm, Blue LED $\times$ 324Ea, Filter Size ; 210 $\times$ 260mm
Software		Imaging Capture Software Included	
White Light		3 $\times$ LED Lamps	
PC Interface(OS)		Window7 or Window10 64bit or Window11	
Camera Interface		USB 2.0	
Dimension & Weight		360 $\times$ 380 $\times$ h700 mm, 24kg	
Packing Size (mm) & Gross Weight (kg)		430 $\times$ 490 $\times$ h830 mm, 28kg	
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz	
Cat. No.	230V	DH.WGD00300	DH.WGD00310
	120V	DH.WGD00301	DH.WGD00311
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Optional Software/Transilluminator/White Light Table

DH.WGD00100	(1) DNA/Protein Analysis S/W, TotalLab™ CLIQS
DH.WGD00121	(2) , CLIQS 1D Pro
DH.WUV01110	UV Lamp, 365nm, 8W
DH.WUV02110	, 312nm, 8W

Ea. Price, ₩

--	--


DNA Analysis S/W

(1) TotalLab™ CLIQS

- Automatic and Manual Lane Detection, Background Subtraction
- Fully Automatic Band Detection, Molecular Size / pI Calibration
- Quantity Calibration, Normalization

(2) CLIQS 1D Pro (additional function)

- Match Band Data Patterns across Large Groups of Lanes
- Choice of Algorithms to Generate Dendrograms ; NeighborJoining, UPGMA, Complete Linkage, Single Linkage, Ward Library, Queries Function



One Stop Analysis

DAIHAN™ UV Transilluminator “WUV”, Standard- & Compact Case-type, 48W & 90W UV Lamp
With Hinged UV-Blocking Safety Cover, Robust, Versatile, 312 · 365 · 312+365 nm, with Certi. & Traceability
자외선 발산기, 기본형 핵산(DNA, RNA), 스탠다드-컴팩트-케이스 타입



3-Year Warranty



(1) Standard Case “WUV-L” Model

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

- Low Cost, Robust, Versatile, Affordable & Gentle Excitation
- Hinged UV-Blocking Safety Cover Protects User from UV Radiation
- Stainless-steel Frame and Robust Steel Housing

- * Ideal for Standard Gel Documentation and Preparative Work, especially for Detection of very small amount of DNA & RNA Stained by EtBr
- * Long-Life Electronic (Low Heat, no Flickering) * High (100%) / Low (75%) Intensity Switch

Specifications & Ordering information

<L1>

Spec. Model	(1) Standard Case-type			(2) Compact Case-type		
	“WUV-L10” Single 365 nm	“WUV-L20” Single 312 nm	“WUV-L50” Dual 365+312 nm	“WUV-M10” Single 365 nm	“WUV-M20” Single 312 nm	
Filter Size	200×200 mm					
Wave Length	Single 365 nm	Single 312 nm	Dual 365+312 nm	Single 365 nm	Single 312 nm	
UV Lamps	15W×6ea (90W)			8W×6ea (48W)		
Power Consumption	95 W			52 W		
Dimension	495×330×h95 mm			350×275×h95 mm		
Weight	10 kg			8 kg		
Packing Size (mm) & Gross Weight (kg)	520×380×h220 mm, 12 kg			380×320×h220 mm, 10 kg		
Electric Reqts.	1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz					
Cat. No.	230V	DH.WUV00010	DH.WUV00020	DH.WUV00050	DH.WUV00110	DH.WUV00120
	120V	DH.WUV02010	DH.WUV02020	DH.WUV02050	DH.WUV00210	DH.WUV00220
Price, Each. ₩						

* Other Specifications are available upon Customer's Request.



3-Year Warranty



(2) Compact Case “WUV-M” Model

Spare UV Lamp

DH.WUV01010	UV Lamp, “UVL365”, 365nm, 15W, for “WUV-L10/-L50”
DH.WUV01020	“UVL312”, 312nm, 15W, for “WUV-L20/-L50”
DH.WUV01110	“UVM365”, 365nm, 8W, for “WUV-M10”
DH.WUV02110	“UVM312”, 312nm, 8W, for “WUV-M20”

Ea. Price, ₩



Custom® Digital UV-intensity Meter “UV340”, 1999~19,990 μW/cm²
With UV Sensor, 70 × 25 × h131mm, 335g, 디지털 자외선 강도계

- * Range ; 1,999~19,990 μW/cm²
- * Wavelength Sensitivity Characteristics ; 290-390 nm
- * Size ; Body 70 × 25 × 131mm, Sensor 46 × 27 × 168mm, Cord Length 700mm
- * Measurement of UVA, UVB is Possible
- * Separated Sensor enables you to Measure from Further.

* Accuracy ; ± (4%rdg+2 dgt)

* Sampling ; 2.5time/sec

<L2>

SL.Met3972 “UV340” Meter, UV-intensity, Digital

Ea. Price, ₩

Jayu® UV-indicator Dot Label, φ 18mm, 자외선 감시 라벨

- * Original Color → Discoloration ; Ultra-high Sensitivity / Clear(White) → Red,
Medium Sensitivity / Clear(White) → Pink
- * Labels Change Color by UV and keep Remaining Once Discolored.
- * The Back has Special Adhesives so that it can Stick Easily.
- * Quantity ; 1pack (100sheet)

<L2>

SL.Lab9043 UV-indicator Dot Label, Ultra-high Sensitivity
SL.Lab9045 , Medium Sensitivity **New**

Pk. Price, ₩

100ea
100ea

DAIHAN™ Rotary Vacuum Evaporator

For 0.1~5Lit Evaporating Flasks

“EVA”

Programmable
Electric
Automatic Lifting

Up to 195mm,
Digital Display Height


Adjustable
Immersion Angle
&
Automatic
Reverse Rotation

1,600cm² Vertical High Cooling Capacity

- 1,600cm² Cooling Surface Condenser.
- Superior Evaporation Capacity ; Max. 23ml/min(Water)

10-Step Adjustable Wide Immersion Angle

- Adjustable 10-Step Wide Immersion Angle 1° ~ 50°



Automatic Reverse Rotation

- Rotated Direction of Evaporating Flask Can be Reversed (Clockwise or Counter-Clockwise)
- Programmable Rotation Cycle Time for User Convenience



Interchangeable Detachable Heating Bath

- Separable Bath Structure from Digital Controller for User Safety and Convenience



Possible to Use a Heating Bath by Itself!

Programmable Electric Automatic Lifting



- Adjustable Electric Automatic Lifting ; Up to 195mm to Accommodate Changing Flask
- Programmable Prefer Height
- Digital Display the Height for User Convenience
- Automatic Raising when Power Error for Sample Protection

Individual Digital Control

- Evaporator and Heating Bath have the Individual Digital Control System for User Convenience.
- Heating Bath can be used in Other Applications Widely.



0.1~5.0 Lit Wide Range Capacity

- Designed for Using 0.1~5Lit Evaporating Flasks Up to 200 °C
- “EVA02” ; 2.5Lit Heating Bath for 0.1~2Lit Flasks
- “EVA05” ; 5.5Lit Heating Bath for 3 & 5Lit Flasks



up to 2Lit



up to 5Lit


DAIHAN™ Premium 0.1~2Lit Rotary Vacuum Evaporator “EVA02”, Vertical-type, DIN-29/32 or ASTM-24/40 Joint With 2.5Lit Heating Bath, Up To 200℃, 195mm Auto Lifting, Cooling Surface 1,600cm², 10-Step Immersion Angle 1°~50° and Automatic Reverse Rotation Function ; Clockwise or Counter Clockwise, 10~300rpm, 0.1~2Lit 프리미엄 회전 진공 농축기



3-Year Warranty

DIN 29/32
or
ASTM 24/40

Premium Rotary Evaporator,
with 2.5Lit Heating Bath(Included)
for 0.1~2Lit Evaporating Flask



▶ Figure ; One of Rotary Vacuum Evaporator
Systematization

Rotary Evaporator is for Research Facilities of Biotechnological, Chemical, Cosmetic, Medical Industries, Pharmaceutical, and Universities. For Crystallization, Distillation, Powder Drying, Product Concentration, and Separation of one or more Solvents.

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- Patented Jog Dial & Touch Button Control System
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Designed for Using 0.1~2Lit Evaporating Flasks, Up To 200℃
 - DURAN-glass Evaporating Flask, 1,000ml, DIN-29/32 or ASTM-24/40 Joint
 - Glass Vapor Duct, L185mm, DIN-29/32 or ASTM-24/40 Joint
 - DURAN-glass Receiving Flask, 1,000ml, Spherical Socket 35/20
- Automatic Reverse Rotation Function ; Clockwise or Counter Clockwise
- Adjustable Wide Immersion Angle ; 10 Steps, 1° ~ 50°
- No Grease is needed for the Flask Connection of Vapor Duct due to Teflon & Viton O-ring Gasket.
- * Evaporator and Heating Bath have the Individual Digital Control System.
- * Separable Bath Structure from Digital Controller for User Safety and Convenience
- * Embedded Heating Elements for Bath
- * High Cooling Surface Condenser ; 1,600cm², Max. 23ml/min
- * Adjustable Electric Automatic Lifting ; Up to 195mm to Accommodate Changing Flask
 - Programmable Prefer Height
 - Digital Display the Height for User Convenience
 - Automatic Raising when Power Error for Sample Protection
- * Rotary Evaporators can be Set with a Vacuum Pump & Circulation Chiller.

Specifications & Ordering information

<L1>

Spec.		Model	“EVA02-D29” for 0.1~2Lit Evaporating Flask DIN-29/32 Joint	“EVA02-A24” for 0.1~2Lit Evaporating Flask ASTM-24/40 Joint
Condenser Type			Vertical-type(Cooling Surface ; 1,600cm ²)	
Flask (Included)	Evaporating		1,000ml DURAN-glass Evaporating Flask with DIN-29/32 Joint	1,000ml DURAN-glass Evaporating Flask with ASTM-24/40 Joint
	Receiving		1,000ml DURAN-glass Receiving Flask, Spherical Socket 35/20	
Glass Vapor Duct			L185mm, DIN-29/32 Joint	L185mm, ASTM-24/40 Joint
Vacuum Seal			Double Sealing Rings made of Teflon + Viton Material	
Evaporating Capability			Max. 23ml/min	
Rotation Speed			10~300rpm	
Controller & Display			PID Controller, Digital Display	
Heating Bath			2.5Lit Heating Bath, Up To 200℃, Dimension ; External(215×205×h179mm), Internal(197×187×h117mm)	
Lifting Type			Adjustable Automatic Electric Lifting ; Up To 195mm Programmable Prefer Height and Digital Display Height	
Immersion Angle			Adjustable 10-step Angle 1° ~ 50°	
Power Consumption			Rotation Motor(100W), Heating Bath ; 2.5Lit Bath(1,200W/220V, 1,100W/110V)	
Material	Body		Aluminum Die Casting	
	Heating Bath		External(PBT), Internal(Ceramic-coated Stainless-steel 304)	
Dimension & Net Weight			530×440×h520(Max. 714)mm, 16kg	
Electrical Reqts.			AC220V, 50/60Hz or AC110V, 60Hz	
Cat. No.	230V		DH.Eva0229	DH.Eva0224
	120V		DH.Eva0229.1	DH.Eva0224.1
Price, Set. ₩				


※ 윤임 별도

- ▶ Optional and Spare Accessories for Evaporator ; See Page 113.
- ▶ Optional Chiller ; See the “Evaporators, Access-Chiller/Cooler” Section.
- ▶ Optional Vacuum Pumps for Evaporators ; See the “Evaporators, Access-Vacuum Pumps” Section.

DAIHAN™ Premium 3 & 5Lit Rotary Vacuum Evaporator “EVA05-D29”, Vertical-type, DIN-29/32 Joint
With 5.5Lit Heating Bath, Up To 200℃, 195mm Auto Lifting, Cooling Surface 1,600cm², 10-Step Immersion Angle 1°~50°
and Automatic Reverse Rotation Function ; Clockwise or Counter Clockwise, 10~300rpm, 3 & 5Lit 프리미엄 회전 진공 농축기



Premium Rotary Evaporator,
with 5.5Lit Heating Bath(Included)
for 3~5Lit Evaporating Flask



► Figure ; One of Rotary Vacuum Evaporator
Systematization



Optional Evapo-Flask, 29/32, 3000ml

Rotary Evaporator is for Research Facilities of Biotechnological, Chemical, Cosmetic, Medical Industries, Pharmaceutical, and Universities. For Crystallization, Distillation, Powder Drying, Product Concentration, and Separation of one or more Solvents.

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- Patented Jog Dial & Touch Button Control System
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Designed for Using 3 & 5Lit Evaporating Flasks, Up To 200℃
 - DURAN-glass Evaporating Flask, DIN-29/32, 5,000ml(Included), 3,000ml(Optional) Joint
 - Glass Vapor Duct, L185mm, DIN-29/32 Joint
 - DURAN-glass Receiving Flask, 3,000ml, Spherical Socket 35/20
- Automatic Reverse Rotation Function ; Clockwise or Counter Clockwise
- Adjustable Wide Immersion Angle ; 10 Steps, 1° ~ 50°
- No Grease is needed for the Flask Connection of Vapor Duct due to Teflon & Viton O-ring Gasket.
- * Evaporator and Heating Bath have the Individual Digital Control System.
- * Separable Bath Structure from Digital Controller for User Safety and Convenience
- * Embedded Heating Elements for Bath
- * High Cooling Surface Condenser ; 1,600cm², Max. 23ml/min
- * Adjustable Electric Automatic Lifting ; Up To 195mm to Accommodate Changing Flask
 - Programmable Prefer Height
 - Digital Display the Height for User Convenience
 - Automatic Raising when Power Error for Sample Protection
- * Rotary Evaporators can be Set with a Vacuum Pump & Circulation Chiller.

Specifications & Ordering information

<L1>

Spec.		Model	“EVA05-D29” for 3 & 5Lit Evaporating Flask, DIN-29/32 Joint
Condenser Type			Vertical-type(Cooling Surface ; 1,600cm ²)
Flask (Included)	Evaporating		5,000ml DURAN-glass Evaporating Flask with DIN-29/32 Joint
	Receiving		3,000ml DURAN-glass Receiving Flask, Spherical Socket 35/20
Glass Vapor Duct			L185mm, DIN-29/32 Joint
Vacuum Seal			Double Sealing Rings made of Teflon + Viton Material
Evaporating Capability			Max. 23ml/min
Rotation Speed			10~300rpm
Controller & Display			PID Controller, Digital Display
Heating Bath			5.5Lit Heating Bath, Up To 200℃ Dimension ; External(324×269×h209mm), Internal(305×250×h118mm)
Lifting Type			Adjustable Automatic Electric Lifting ; Up To 195mm Programmable Prefer Height and Digital Display Height
Immersion Angle			Adjustable 10-step Angle 1° ~ 50°
Power Consumption			Rotation Motor(100W), Heating Bath ; 5.5Lit Bath(1,500W/220V, 1,100W/110V)
Material	Body		Aluminum Die Casting
	Heating Bath		External(PBT), Internal(Ceramic-coated Stainless-steel 304)
Dimension & Net Weight			570×440×h520mm(Max. Height 714mm), 18kg
Electrical Reqs.			AC220V, 50/60Hz or AC110V, 60Hz
Cat. No.	230V		DH.Eva0529
	120V		DH.Eva0529.1
Price, Set. ₩			

※ 운임 별도


► Optional 3000ml DURAN-glass Evaporating Flask 29/32 Joint

SL.Fla2071 Evapo-Flask, 29/32, 3000ml

Ea. Price, ₩

- Optional and Spare Accessories for Evaporator ; See Page 113.
- Optional Chiller ; See the “Evaporators, Access-Chiller/Cooler” Section.
- Optional Vacuum Pumps for Evaporators ; See the “Evaporators, Access-Vacuum Pumps” Section.

Optional & Spare Accessories for Premium Rotary Vacuum Evaporator "EVA"



Ordering Information

<L6>

Seq. No.	Applicable Model	Cat. No.	Description	Ea. Price, ₩
(1)	EVA02-D29, EVA02-A24, EVA05-D29	DH.EvaA10050	Glass Feeding Valve Stopcock, 19/38, with od. $\phi 7 \times$ L350mm PTFE Tube	
(2)		DH.Con2059	Glass Condenser , $\phi 85 \times 471$ mm, with (3) PBT GL25 Cap & (4) PP GL10 Opentop Cap and PTFE Tube, “GC05”	
(3)		DU.2924013	Spare PBT Closetop Cap, GL25, Red, for Closure	10ea
(4)		DH.EvaA30010	Spare PP Opentop Cap, GL10, Hole $\phi 3.2$ mm, with od. $\phi 3 \times$ L400mm PTFE Tube	
(5)		SL.C17454	Stainless-steel(304) Ball Joint Clamp, w/Screw, Joint 35/20	
(6)		SL.Fla2081	DURAN-glass Receiving Flask, 250 ml , with Spherical Socket 35/20	
		SL.Fla2082	DURAN-glass Receiving Flask, 500 ml , with Spherical Socket 35/20	
		SL.Fla2083	DURAN-glass Receiving Flask, 1,000 ml , with Spherical Socket 35/20	
		SL.Fla2084	DURAN-glass Receiving Flask, 2,000 ml , with Spherical Socket 35/20	
		SL.Fla2085	DURAN-glass Receiving Flask, 3,000 ml , with Spherical Socket 35/20	
(7)		DH.EvaA30020	Condensate Silicone Trap for Glass Condenser	
(8)		DH.EvaA30030	Condenser Locknut, PBT, Hole $\phi 51$ mm	
(9)		DH.EvaA30040	Circlip Spring Ring, Stainless-steel 304, od. $\phi 7.8$ mm($\phi 1.0$ mm thickness)	
(10)		DH.EvaA30050	Vacuum Seal, Rubber & PTFE, od. $\phi 50$ mm, id. $\phi 25.4$ mm	
(11)		DH.EvaA30060	Viton O-ring, Black, id. $\phi 26.5$ mm	
(12)		DH.EvaA30070	Teflon O-ring, White, id. $\phi 26.5$ mm	
(13)		DH.EvaA30080	Fixing Stand Rod for Glass Condenser, Stainless-steel(304)	
(14)		DH.EvaA10024	Glass Vapor Duct L185mm, with ASTM-24/40 Joint	
		DH.EvaA10029	Glass Vapor Duct L185mm, with DIN-29/32 Joint	
(15)		WI.0140024	Suferfix Screwcap-type Joint Clamp, for ASTM-24/40 Joint, Gray	
	WI.0140029	Suferfix Screwcap-type Joint Clamp, for DIN-29/32 Joint, Blue		
(16)	EVA02-D29	SL.Fla2060	DURAN-glass Evaporating Flask, 100 ml , with DIN-29/32 Joint	
		SL.Fla2062	DURAN-glass Evaporating Flask, 250 ml , with DIN-29/32 Joint	
		SL.Fla2064	DURAN-glass Evaporating Flask, 500 ml , with DIN-29/32 Joint	
		SL.Fla2066	DURAN-glass Evaporating Flask, 1,000 ml , with DIN-29/32 Joint	
		SL.Fla2068	DURAN-glass Evaporating Flask, 2,000 ml , with DIN-29/32 Joint	
	EVA02-A24	SL.Fla2061	DURAN-glass Evaporating Flask, 100 ml , with ASTM-24/40 Joint	
		SL.Fla2063	DURAN-glass Evaporating Flask, 250 ml , with ASTM-24/40 Joint	
		SL.Fla2065	DURAN-glass Evaporating Flask, 500 ml , with ASTM-24/40 Joint	
		SL.Fla2067	DURAN-glass Evaporating Flask, 1,000 ml , with ASTM-24/40 Joint	
		SL.Fla2069	DURAN-glass Evaporating Flask, 2,000 ml , with ASTM-24/40 Joint	
	EVA05-D29	SL.Fla2071	DURAN-glass Evaporating Flask, 3,000 ml , with DIN-29/32 Joint	
SL.Fla2073		DURAN-glass Evaporating Flask, 5,000 ml , with DIN-29/32 Joint		
(17)	EVA02-D29, EVA02-A24, EVA05-D29	SM.28623060	Red PBT Screwcapped PP & PTFE Hose Connector Kit, Bent, GL14	
(18)		SM.28623000	Red PBT Screwcapped PP & PTFE Hose Connector Kit, Straight, GL14	
(19)		DU.2924008	PBT Closetop Cap, GL14, Red, for Closure	10ea

DAIHAN™ 0.5~2Lit Digital Rotary Evaporator “WEV-1001V” & “WEV-1001L”, Vertical & Diagonal-type
With 6.5 Lit Digital Bath 100℃, PID Control & Lifting Table, 10~180rpm, Cooling Surface 1,300cm², Certi. & Traceability, 회전식 증발 농축기



3-Year Warranty

(1) Vertical-type “WEV-1001V”



3-Year Warranty


(2) Diagonal-type “WEV-1001L”

- Supplied by DAIHAN Korea
- Cooling Surface 1,300cm²
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Distilling (Mixed) Solvent, Condensation Solution and Suspension, Concentrating Liquid, Re-crystallization Dry Powder and Granulate, Synthesizing Chemical, Soxhlet Extraction, etc.
- Included Two Evaporating Flasks ; 1×0.5Lit & 1×1Lit
- The Temperature adopts PID Control, which ensures accurate Temperature Control.
- Speed Display, Operate Visually and easily
- * Digital Speed Display, Digital Temp Display on Water Bath.
- * Double Sealing with PTFE and Vinton Materials to ensure High Vacuum Degree.
- * The Air Extraction Nozzle of Condenser is set down, and it is Safer after Connection with Vacuum Pipe
- * Building Block Design is Easily to be Extended.
- * Easy Operation of Lifting, Move up and Down Smoothly and Easy to Lock.
- * Dedicated Electric Machinery, Worm and Gear Drive Steadily.
- * The Separation of Bath and the Mainframe makes Configure Flexible.
- * These Rotary Evaporators could Set with Circulation Aspirator Pump (P.120)

Specifications & Ordering information

<L1>

Spec.	Model	(1) “WEV-1001V” Vertical-type	(2) “WEV-1001L” Diagonal-type
Condenser Type		Vertical-type	Diagonal-type
Bottle Capacity	Evaporating	Included 2 × Evaporating Flasks ; 1 × 0.5Lit & 1 × 1Lit, Available 100-/250-ml, 2Lit Evaporating Flask(Optional)	
	Receiving	Included 1Lit Flask 1ea	
Rotary Speed		10~180 rpm	
Evaporating Capability		20ml/min	
Cooling Surface		1,300 cm ²	
Achievable Vacuum		399.9Pa(below 3mmHg)	
Bath Capacity		6.5 Lit / φ 254 × 130 mm	
Temp Range		0℃~100℃	
Functional of Lifting Table		Weight Balancing Gliding Elevating & Manual Lifting	
Elevating Stroke		100~250 mm	
Vacuum Sealing		Double Sealing Rings made of Teflon + Vinton Materials	
Power Consumption		Motor ; 25 W, Heater ; 1.05kW	
Controller & Display		PID Controller, Digital display	
Dimension & Net Weight		560 × 320 × h660mm / 9.5 kg	750 × 320 × h540mm / 9.5 kg
Electrical Reqts.		AC220V, 50/60Hz or AC110V, 60Hz	
Cat. No.		DH.WEV01001V	DH.WEV01001L
Price, Each. ₩			



- ※ NOTE! ▶ Optional Vacuum Pumps for Evaporators ; See the “Evaporators, Access-Vacuum Pump” Section.
▶ For the Information of Spare Parts ; Contact Us Directly.

DAIHAN™ 5 · 10 · 20 Lit Digital Rotary Evaporator “WEV-1005” · “WEV-1010” · “WEV-1020”, Large Capacity
 With Digital Controlled Stainless-steel Bath 99℃, Electric Lift Bath, Vertical-type, 20~130rpm, 대용량 회전식 증발 농축기



3-Year Warranty

- Supplied by DAIHAN Korea
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Distilling (Mixed) Solvent, Condensation Solution and Suspension, Concentrating Liquid, Re-crystallization Dry Powder and Granulate, Synthesizing Chemical, Soxhlet Extraction, etc.
- Teflon (PTFE) and Fluorine Rubber Dual Rotating Seal, Patented Technology to ensure High Vacuum.
- A Vacuum Switching Valve ; without Affecting the System Vacuum, Distillation of the Situation and the Solution can be Continuously Collected.
- Electric Lift - Stainless-steel Water Bath, Easy to Use.
- Cooling Surface ; “WEV-1005” total 2,700cm²
 “WEV-1010” total 6,500cm²
 “WEV-1020” total 11,500cm²

* Accurate and Reliable Temp. Control due to the Introduction of Advanced Frequency Control System

* PTFE Discharge Valve, Corrosion Resistance, Pollution-free.

* Rotary Evaporators can be Set with Circulation Aspirator Pump (P.120)

Specifications & Ordering information

<L1>

Spec.	Model	(1) “WEV-1005” Evaporating Flask 5Lit	(2) “WEV-1010” Evaporating Flask 10Lit	(3) “WEV-1020” Evaporating Flask 20Lit
Evaporating Flask		5 Lit, Flange Opening ϕ 50 mm	10 Lit, Flange Opening ϕ 125 mm	20 Lit, Flange Opening ϕ 125 mm
Receiving Flask		3 Lit	5 Lit	10 Lit
Condenser		Water Refrigerating	Standing, Main & Auxiliary Cooling and Highly Efficient 3 Circulating Cold Trap	
Temp. Range		0℃ ~ 99℃		
Rotary Speed		20 ~ 140 rpm	20 ~ 130 rpm	
Evaporating Capability		Water ; 2 Lit/h, Alcohol ; 4 Lit/h	Water ; 3.2 Lit/h, Alcohol ; 6.5 Lit/h	Water ; 5 Lit/h, Alcohol ; 11 Lit/h
Cooling Surface		Condenser 2,700 cm ²	Main Condenser 4,500 cm ² Auxiliary Condenser 2,000 cm ²	Main Condenser 8,000 cm ² Auxiliary Condenser 3,500 cm ²
Achievable Vacuum		399.9Pa(below 3mmHg)		
Functional of Lifting Table		Electric		
Elevating Stroke		0 ~ 150 mm	0 ~ 160 mm	0 ~ 190 mm
Vacuum Sealing		PTFE and Fluorine Rubber Dual Rotating Seal		
Motor		250 W, Phase DC Induction Motor		
Heating Boiler		Stainless-steel Water Bath		
Display		Digital display		
Power Consumption		2.3 kW	1Phase 220V, 3.8kW or 3Phase 380V, 6.3kW	
Dimension (mm)		835×400×h1090	990×550×h1655	1195×635×h2070
Electrical Reqts.		1Phase 220V, 50/60Hz	1Phase 220V, 50/60Hz	
Cat. No.		DH.WEV01005	DH.WEV01010	DH.WEV01020
Price, Each. ₩				

※ 운임 별도



3-Year Warranty

※ **NOTE!** ▶ Optional Vacuum Pumps for Evaporators ; See the “Evaporators, Access-Vacuum Pump” Section.
 ▶ For the Information of Spare Parts ; Contact Us Directly.

DAIHAN™ 50 Lit Digital Rotary Evaporator “WEV-1050”, Large Capacity, Vertical Type, with Electric & Manual Lift Bath
With Digital Controlled Stainless-steel Bath 99℃, 20~110rpm, 19Lit/h, Cooling Surface 14,500cm², 대용량 회전식 증발 농축기



3-Year Warranty

- Supplied by DAIHAN Korea
- Cooling Surface ; 14,500cm²
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Distilling (Mixed) Solvent, Condensation Solution and Suspension, Concentrating Liquid, Re-crystallization Dry Powder and Granulate, Synthesizing Chemical, Soxhlet Extraction, etc.
- Large Capacity & Diameter Rotary Evaporator Flask - 50Lit
- Teflon (PTFE) and Fluorine Rubber Dual Rotating Seal, Patented Technology to Ensure High Vacuum.
- Stainless-steel Water Bath - Electric & Manual Lift, Easy to Use
- * The Introduction of Foreign Advanced Frequency Control System, Temperature Control System, Accurate and Reliable.
- * A Vacuum Switching Valve, without Affecting the System Vacuum Distillation of the Situation and the Solution can be Continuously Collected.
- * Water Bath Electric Lift and Manual Lift Support, Easy to Use.
- * PTFE Discharge Valve, Corrosion Resistance, Pollution-free.
- * This Rotary Evaporator could Set with Circulation Aspirator Pump (P.120)





Specifications & Ordering information

<L1>

Spec.	Model	“WEV-1050”
Evaporating Flask		50 Lit, Flange Opening ϕ 125 mm
Receiving Flask		20 Lit
Condenser		Standing, Main & Auxiliary Cooling and Highly Efficient 3 Circulating Cold Trap
Temp. Range & Rotary Speed		0℃~99℃, 20~110 rpm
Evaporating Capability		Water ; 9 Lit/h, Alcohol ; 19 Lit/h
Cooling Surface		Main Condenser 9,500 cm ² Auxiliary Condenser 5,000 cm ²
Achievable Vacuum		399.9Pa(below 3mmHg)
Functional of Lifting Table & Stroke		Electric & Manual, 0~180 mm
Vacuum Sealing		PTFE and Fluorine Rubber Dual Rotating Seal
Motor		250 W, Phase DC Induction Motor
Heating Boiler		Stainless-steel Water Bath
Display		Digital display
Dimension (mm)		1345 × 770 × h2140
Electrical Reqts.		6.3 kW, 3Phase 380V, 50Hz
Cat. No.		DH.WEV01050
Price, Each. ₩		

※ 운임 별도

The Details of Spare Parts for 50 Lit Digital Rotary Evaporator “WEV-1050”



※ NOTE! ▶ Optional Vacuum Pumps for Evaporators ; See the “Evaporators, Access-Vacuum Pump” Section.
▶ For the Information of Spare Parts ; Contact Us Directly.

DAIHAN™ -30℃~R/T Classical Chiller "CHI", Refrigerated External Circulator, Fill-5 · 10 · 20 · 30 · 50Lit with Insulated STS Hose
Ideal for Evaporator, Reactor &c. Cooling Line, Lift 6 · 12m, Cooling Capa 0.6~3.5kW, 냉각 써클레이터/칠러, 단열 튜빙 포함

NEW

CE

-30℃ ~ R/T

3-Year Warranty



Internal Bath



Connection Ports



Insulated Flexible STS Hose Included

- Supplied by DAIHAN Korea ■ CE Certified ■ Max. Lift ; 12m
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information.
- Ideal for Wide-range Reaction System and Evaporator.
- Precise Temperature Control System, $\pm 1^\circ\text{C}$
- Microprocessor-based Control ensures Stable and Consistent Thermal Performance
- Stainless-steel Internal Bath for Prevention of Corrosion
- Includes Two Insulated Flexible Stainless-steel Hoses (0.7m & 1m)
- * Various Filling Volumes Available ; 5 · 10 · 20 · 30 · 50 Lit
- * Effective Cooling System with Rapid Temperature Stabilization
- * Useful for 1~50Lit Various Reactors and Evaporators.
- * Equipped with Caster Wheels for Safe and Easy Mobility
- * High Quality Circulation Pump ; 30-Lit/min Max Flow Rate

Specifications & Ordering information

<L1>

Spec.	Model	“CHI5-30”	“CHI10-30”	“CHI20-30”	“CHI30-30”	“CHI50-30”
Temp. Range & Accuracy		-30℃-Range, ±1℃ at Set Temperature				
Cool-Fluid. Filling Volume		5Lit max	10Lit max	20Lit max	30Lit max	50Lit max
Useful for 1~50 Lit Evaporators, etc.		for 1~5Lit Evaporator	for 5~20Lit Evaporator	for 20~30Lit Evaporator	for 50Lit Evaporator	for 100Lit Evaporator
Refrigerator Capacity		0.55kW	0.75kW	1.09kW	1.5kW	3.75kW
Material		Internal ; Stainless-steel, External ; Powder Coated steel				
Lid & Top Opening Dia.		φ 180mm		φ 240mm	φ 320mm	
Fluid In/Outlet Port		ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel				
Circulation Hose		0.7m & 1m, Included				
Pump Power, Lift Height, m	Max. Flow Rate	20Lit/min				30Lit/min
	Max. Lift	8m				12m
Cooling Capacity	at +15℃	0.9kW	2.0kW	2.8kW	4.6kW	10.5kW
	at 0℃	0.7kW	1.7kW	2.3kW	3.5kW	9.3kW
	at -30℃	0.3kW	0.7kW	0.6kW	0.6kW	0.5kW
External Dimension (mm)		370×470×h680	450×520×h800	570×490×h860	640×540×h1000	740×640×h1190
Electrical Reqts.		1 Phase 230V, 50 / 60 Hz				
Cat. No.		DH.CHI05.30	DH.CHI10.30	DH.CHI20.30	DH.CHI30.30	DH.CHI50.30
Price, Set. ₩						

※ 운임 별도

Spare Insulated Flexible Stainless-steel Hose

SL.Tub7033 Insulated Flexible STS Hose Only, 1/2", id17×od26mm, Overall od60mm, L0.7m
SL.Tub7034 , L1m

Ea. Price, ₩

※ For More Information about Circulator Hoses ;
SEE the "Bath Circulators, Accessories" Section !

DAIHAN™ $-20^{\circ}\text{C} \sim \text{R/T}$ Premium Chiller “CHI”, Heavy-duty Refrigerated External Circulator, Fill-17 · 29 · 47 Lit
Ideal for Evaporator/Reactor &c. Cooling Line, Lift 27m, Cooling Capa 0.87 · 1.3 · 3.0 kW, 다용도 냉각 써클레이터/칠러



- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Wide-range Reaction System and Evaporator. ■ Max. Lift 27m
- Assembly-type Safety Connection Ports/Pipe for Fluid In/Out-let ; ISO RC 1/2" Thread, Stainless-steel.
- Digital PID Control System Implementing Superior Temperature Accuracy.
- Embedded External Nozzle Part for the Prevention of Nozzle Damaged.
- Stainless-steel Internal Bath for Prevention of Corrosion.
- * Wide-range Filling Volume 17Lit · 29Lit · 47Lit
- * Effective Cooling Capacity
- * Useful for 1~50Lit Various Reactors and Evaporators. * Emergency Stop Button for safety.
- * High Quality Circulation Pump ; 30-Lit/min Max Flow Rate and 3.3-bar Max. Pressure.
- * Built in Pump Pressure Gauge and Compressor Pressure Gauge.

Specifications & Ordering information

<L1>

Spec.	Model	“CHI-20” Filling Volume Max. Capa. 17Lit	“CHI-30” Filling Volume Max. Capa. 29Lit	“CHI-50” Filling Volume Max. Capa. 47Lit
Temp Range & Stability		$-20^{\circ}\text{C} \sim \text{Room Temp.}$ $\pm 0.2^{\circ}\text{C}$ (at 15°C) -Stability		
Cool-Fluid. Filling Volume		17 Lit max	29 Lit max	47 Lit max
Useful for 1~50Lit Evaporators, etc.		for 1~20Lit Evaporator	for 10~30Lit Evaporator	for 20~50Lit Evaporator
Refrigerator Capacity(HP)		1 HP	1.5 HP	2 HP
Lid & Top Opening Dia.		$\phi 100 \text{ mm}$		
Fluid In/Out-let Connector		ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel		
Pump Power, Lift Height, m	Max. Flow Rate	30 Lit/min		
	Max. Pressure & Lift	3.3 bar Max. Lift ; 27m		
Cooling Capacity	at $+20^{\circ}\text{C}$	1.7 kW	2.8 kW	6.5 kW
	at 0°C	0.87 kW	1.3 kW	3.0 kW
	at -20°C	0.32 kW	0.45 kW	1.2 kW
External Dimension		450×630×h650 mm	550×830×h975 mm	605×975×h1265 mm
Material		Internal ; Stainless-steel, External ; Powder Coated steel		
Electrical Reqts.		1 Phase 230V, 50 / 60 Hz		
Cat. No.	60Hz	DH.Chi0020	DH.Chi0030	DH.Chi0050
	50Hz	DH.Chi0025	DH.Chi0035	DH.Chi0055
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.


※ 운임 별도

※ Cooling/Heating Transfer Fluids ;
See the Page 65.

- ▶ SYN THERM 7008® $-40^{\circ}\text{C} + 400^{\circ}\text{C}$
- ▶ SUPER THERM FG-20® $-10^{\circ}\text{C} + 300^{\circ}\text{C}$
- ▶ SYN THERM 22L® $-5^{\circ}\text{C} + 320^{\circ}\text{C}$

Note ! Cooling Transfer Fluid Should be an Antifreezing Solution.

※ Insulated Circulation Hoses ; See the
“Bath Circulators, Accessories” Section !



▶ Figure ; One of the Applications in the Large Volume Evaporator Line.

DAIHAN™ -35℃ +150℃ Digital Precise Refrigerated/Heating Bath Circulator "CL", ±0.2℃


With Flat Lid, Digital Fuzzy Control, CFC-free, Certi. & Traceability, Capa. 8 · 12 · 22 · 30 Lit, Flow 25Lit/min, Lift 4m

Ideal for Cooling/Heating Line of Facility, External Circulating Connect od.12 & id.10mm, 정밀 냉동/가열 배스 겸용 써클레이터

CE -35℃ +150℃



3-Year Warranty

-35℃ +150℃ Refrigerated / Heating Bath Circulator
Max. 8Lit, "CL-8" & Max. 12Lit, "CL-12"External Circulation
NozzleStainless-steel Flat Lid
Included※ Note ; Cooling Transfer Fluid Should be
an Antifreezing Solution※ Validation (IQ, OQ) Service
Available upon Request► Insulated Lid
(Optional)

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- Advanced Refrigeration System Ensures Temperature Performance
- Patented Jog-Dial Control System and Capacitive Touch Button

* RS232C Interface for Remote Monitoring and Controlling with PC * All models Built in Caster
 * User's Self-compensation Function to Control Difference between Real Temperature and Indicated
 Temperature (±10.0℃). * Much Faster Cool-down Time * Proper Liquid Level Line in Bath
 * Warning and Auto-Stop Function in case of Liquid Shortage.

Specifications & Ordering information

<L1>

Spec.	Model	"CL-8" Max. Capa. 8 Lit	"CL-12" Max. Capa. 12 Lit	"CL-22" Max. Capa. 22 Lit	"CL-30" Max. Capa. 30 Lit
Temp Range & Accuracy		-35℃+150℃, ±0.1℃			
Temp Resolution		0.1℃-Display, 0.1℃-Control			
Temp Uniformity & Probe		±0.2℃ at -20℃, PT100			
Pump		Max Flow Rate ; 25Lit/min(Controlled 5 Steps Flow Rate Speed) Max Pressure ; 7.5 psi, Max Lift ; 4m			
Bath Volume		8 Lit	12 Lit	22 Lit	30 Lit
Heating & Refrigerator		1 kW, 1/2HP	1.5 kW, 1/2HP	2 kW, 7/8HP	2.2 kW, 1HP
Cooling Capacity	at +20℃	500 W	530 W	860 W	1030 W
	at 0℃	340 W	370 W	620 W	760 W
	at -20℃	200 W	210 W	400 W	500 W
Total Power Consumption		1.7 kW	2.2 kW	3 kW	3.3 kW
Compressor Control		Intelligent Control of Refrigerator for preventing Compressor Overload			
External Circulation Nozzle		id.10.5 / od.12.5mm(For Inner Circulat. ; to Connect directly the in-/outlet Nozzles)			
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end			
Display		Digital LCD with Back-Light Function			
Controller		Digital Fuzzy Controller with Jog-Dial Knob (Turn+Push)			
Safety Circuit		Over-Temp/-Current protection, Sensor Error Detection, Leakage Breaker			
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel Insulation ; 40 mm Thickness Urethane Foaming			
Lid & Drain Valve		Includes Stainless-steel Flat Lid & Ball Valve			
Dimension (w×d×h)	Effective Space (mm)	175×153×150	230×181×150	290×220×200	330×290×200
	Internal (mm)	175×320×150	230×350×150	290×390×200	330×460×200
	External (mm)	400×523×807	400×523×807	483×563×850	494×633×906
Packing Size (mm) & Gross Weight (kg)		530×630×985	530×630×985	615×670×1025	630×740×1075
		50 kg	55 kg	65 kg	78 kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz			
Cat. No.	230V	DH.WCL00408	DH.WCL00412	DH.WCL00422	DH.WCL00430
	120V	DH.WCL00508	DH.WCL00512	DH.WCL00522	DH.WCL00530
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

► Optional Microporous Ceramic Insulation Lid
T15mm(10mm Microporous Ceramic + 5mm Ceramic Paper)

DH.WCR.L008 Insulation Lid, for CL-8, 240×178×T15mm
 DH.WCR.L012 for CL-12, 316×206×T15mm
 DH.WCR.L022 for CL-22, 396×245×T15mm
 DH.WCR.L030 for CL-30, 436×316×T15mm

Ea. Price, ₩



► Insulated Circulation Hose

⇒ See the "Bath Circulators, Accessories" Section !

DAIHAN™ 105mbar Oil Free Vacuum Pump “PU-300”, with Fine Vacuum Regulator and Moisture Trap, 20Lit/min With Prote™ Overflow Protection, Built-in Thermal Protection, Protection Class IP30, 진공 펌프



3-Year Warranty


105mbar Oil Free Vacuum Pump
“PU-300”

- Supplied by DAIHAN Korea
- 진공처리가 필요한 다양한 실험에서 사용이 가능
- CE, CSA, and RoHS Certified
- Used for General Filtration, Microbiological Test, Suspended Solid Test, Solid-Phase Extraction, Microplate Vacuum Manifold, Air Sampling, etc.

- * Driven by Piston without Lubricant, Regular Oil Change and Maintenance ; No Oil Pollution
- * Prote™ Overflow Protection is Installed to Prevent a Large Amount of Liquid from being Sucked into the Pump.
- * Equipped with Filter Cartridge in Air Inlet to Filter Particle and Moisture to Prolong the Life of Pump.

<L1>

Spec.	Model	“PU-300”	
Speed of Evacuation		20 Lit/min	23 Lit/min
Max. Vacuum		105 mbar, 78.8 torr	
Motor Power & Max. Current		60 W, 0.3A	65 W, 0.7A
Motor Rotation & Noise Level		1450 rpm, 50 dB	1750 rpm, 50 dB
Hose Barb		8mm	
Size & Weight		268 × 135 × h204mm, 4.1kg	
Electric Rqts.		220V, 50Hz	110V, 60Hz
Cat. No.		DH.Pu0300	DH.Pu0301
Price, Each. ₩			

DAIHAN™ 10mbar Anticorrosive Vacuum Diaphragm Pump “PU-20”, with Fine Regulator and Moisture Trap, 18Lit/min With Built-in Thermal Protection, Quiet and Low Vibration, 내부식성 저소음/저진동 진공 다이어프램 펌프

3-Year Warranty


10mbar Anticorrosive Vacuum
Diaphragm Pump “PU-20”

- Supplied by DAIHAN Korea
- 진공처리가 필요한 다양한 실험에서 사용이 가능하며 최고 10mbar의 고진공 성능을 보여줌
- Used for Solvent purification, Gel dryer, Rotary Evaporator, Vacuum Filtration, etc.
- Resistance to Chemicals ; the Parts in Contact with the Gas are PTFE Coated.

- * Driven by Diaphragm without Lubricant, Regular Oil Change and Maintenance ; No Oil Pollution
- * Quiet and Low Vibration
- * Built-in Thermal Protection

<L1>

Spec.	Model	“PU-20”	
Speed of Evacuation		18 Lit/min	
Max. Vacuum		10 mbar, 7.5 torr	
Motor Power & Max. Current		90 W, 0.5A	
Motor Rotation & Noise Level		1450 rpm, 52dB	
Hose Barb		8mm	
Size & Weight		262 × 236 × h193mm, 7.4kg	
Electric Rqts.		220V, 50Hz	
Cat. No.		DH.Pu8020	
Price, Each. ₩			

DAIHAN™ 0.098MPa Multiuse Water Jet Vacuum Pump/Aspirator “WEV-003S” Ideal for Evaporator 1~5Lit, with Multi-Taps & Safety Valve, Bath 15Lit, 다용도 아스피레이터, 워터젯 진공펌프

3-Year Warranty



“WEV-003S” ideal for 1~5Lit Evaporator

- Supplied by DAIHAN Korea
- Resist Corrosion by Acid, Alkali and Solvents.
- Cyclically use Water
- Good for Evaporation, Distillation, Crystallization, Drying, Sublimation, Filtration &c.
- Fluorine Rubber Sealing
- Need to Replace the Water in the Flume Regularly
- * Can be used alone or in Parallel 2 Tap & 2-Display
- * Stainless-steel(#304) Body, Copper Ejector, Tee Junction, PP Back Valve

<L1>

Spec.	Model	“WEV-003S” Ideal for 1~5Lit Evaporator	
Max. Vacuum Degree		0.098 Mpa	
Pump Power & Water Bath Volume		180W & 15Lit	
Flow Rate		80 Lit/min	
Taps & Air Sucking Amount		2 Tap, 10 Lit/min (1 Tap)	
External Circulation Nozzle		id. 5/od. 10mm (id. 8/od.11 mm Silicone tube Connected for Circulation)	
Safety Function		Safety Check Valve	
Drain Pipe		id. 15mm	
Material of Water Tank / Body		PP / PP	
Dimension & Net Weight		385 × 280 × h420 mm / 11 kg	
Electric Rqts.		230V, 50/60Hz	
Cat. No.		DH.WEV0003S	
Price, Each. ₩			



(1) Without Return Holes

(2) With Return Holes

Premium DURAN Glass Rotary Evaporator Safety Trap / Splash Head Adapter

With ASTM & DIN Joints, 100~500ml, 에바포레이터 트랩(헤드 어댑터)

■ Supplied by SciLab Korea

■ Custom-Made Available / 규격외 Items 자유롭게 주문하세요!

(1) Evaporator Trap, without Return Holes <L1>

① with ASTM Joint

Ea. Price, ₩

SL.Ada2280	100ml, 24/40
SL.Ada2281	250ml, 24/40
SL.Ada2282	500ml, 24/40

② with DIN Joint

Ea. Price, ₩

SL.Ada2283	100ml, 29/32
SL.Ada2284	250ml, 29/32
SL.Ada2285	500ml, 29/32

(2) Evaporator Trap, with Return Holes

① with ASTM Joint

Ea. Price, ₩

SL.Ada2290	100ml, 24/40
SL.Ada2291	, 29/42
SL.Ada2292	250ml, 24/40
SL.Ada2293	, 29/42
SL.Ada2294	500ml, 24/40
SL.Ada2295	, 29/42

② with DIN Joint

Ea. Price, ₩

SL.Ada2296	100ml, 24/29
SL.Ada2297	, 29/32
SL.Ada2298	250ml, 24/29
SL.Ada2299	, 29/32
SL.Ada2300	500ml, 24/29
SL.Ada2301	, 29/32

(3) Evaporator Trap, with Glass Filter Disk

* Filter Porosity ; "P1" 100~160 μ m ; Protection of Condenser Contamination

① with ASTM Joint

Ea. Price, ₩

SL.Ada2303	250ml, 24/40
SL.Ada2304	, 29/42

② with DIN Joint

Ea. Price, ₩

SL.Ada2305	250ml, 24/29
SL.Ada2306	, 29/32



(3) With Filter



SciLab™ Evaporating Flask, with ASTM & DIN Joint, 100~3,000ml

For Rotary Vacuum Evaporator, Boro Glass 3.3, 에바포레이팅 플라스크

■ Supplied by SciLab Korea

* Made of Boro-glass 3.3

<L1>

Capacity(ml), Size Approx., od×h(mm)	(1) With ASTM 24/40 Joint			(2) With DIN 29/32 Joint		
	Cat. No.	Ea. Price, ₩	Pk. Price, ₩	Cat. No.	Ea. Price, ₩	Pk. Price, ₩
100ml, ϕ 61×h105mm	SL.Fla2901		12ea	SL.Fla2911		12ea
250ml, ϕ 80×h145mm	SL.Fla2902		12ea	SL.Fla2912		12ea
500ml, ϕ 100×h160mm	SL.Fla2903		8ea	SL.Fla2913		8ea
1,000ml, ϕ 125×h202mm	SL.Fla2904		6ea	SL.Fla2914		6ea
2,000ml, ϕ 155×h260mm	SL.Fla2905		2ea	SL.Fla2915		2ea
3,000ml, ϕ 190×h300mm	SL.Fla2906		2ea	SL.Fla2916		2ea



100~3,000ml

5,000ml

SciLab™ Premium DURAN-glass Evaporating Flask, with ASTM & DIN Joints, 100~5,000ml

For Rotary Vacuum Evaporator, Boro Glass 3.3, 에바포레이팅 플라스크

■ Supplied by SciLab Korea

■ Custom-Made Available / 규격외 Items 자유롭게 주문하세요!

<L1>

Capacity(ml), Size Approx., od×h(mm)	(1) With ASTM 24/40 Joint		(2) With DIN 29/32 Joint	
	Cat. No.	Ea. Price, ₩	Cat. No.	Ea. Price, ₩
100ml, ϕ 60×h118mm	SL.Fla2061		SL.Fla2060	
250ml, ϕ 81×h148mm	SL.Fla2063		SL.Fla2062	
500ml, ϕ 101×h178mm	SL.Fla2065		SL.Fla2064	
1,000ml, ϕ 126×h218mm	SL.Fla2067		SL.Fla2066	
2,000ml, ϕ 158×h248mm	SL.Fla2069		SL.Fla2068	
3,000ml, ϕ 176×h290mm	SL.Fla2070		SL.Fla2071	
5,000ml, ϕ 223×h320mm New	SL.Fla2074		SL.Fla2073	



SciLab™ 5~50Lit Big Volume Evaporating Flask, High Quality Heavy Wall Made of Boro Glass 3.3 Flask, 대용량 에바포레이팅 플라스크

- Supplied by SciLab Korea
- Wide Neck Flange BUCHI & EYELA Evaporator와 완벽 호환
- For (1) DAIHAN™ / SciLab™ / etc Evaporators, (2) BUCHI Evaporators, (3) EYELA Evaporators
- Custom-Made Available ; <ex.> DURAN® Flask로도 주문제작 가능 !

<L1>

Capacity(ml), Size Approx., od×h(mm)	(1) For DAIHAN™ / SciLab™ Evaporators		(2) For BUCHI Evaporators		(3) For EYELA Evaporators	
	Flange ; id/od 100/125mm Neck ; od 110×h155mm		Flange ; id/od 120/150mm Neck ; od 135×h115mm		Flange ; id/od 100/138mm Neck ; od 106×h115mm	
	Cat. No.	Ea. Price, ₩	Cat. No.	Ea. Price, ₩	Cat. No.	Ea. Price, ₩
5Lit, φ223×h353mm	SL.Fla2074D*		SL.Fla2074B		SL.Fla2074E	
10Lit, φ279×h409mm	SL.Fla2075D		SL.Fla2075B		SL.Fla2075E	
20Lit, φ345×h475mm	SL.Fla2076D		SL.Fla2076B		SL.Fla2076E	
30Lit, φ415×h545mm	SL.Fla2078D		SL.Fla2078B		SL.Fla2078E	
50Lit, φ480×h610mm	SL.Fla2077D		SL.Fla2077B		SL.Fla2077E	

* Only SL.Fla2074D Flange ; id/od 30/50mm, Neck ; od 40×h50mm **



SciLab™ Premium DURAN-glass Receiving Flask, with 35/20 Spheric Socket, 100~3,000ml For Rotary Vacuum Evaporator, Boro Glass 3.3, 리시빙 플라스크

- Supplied by SciLab Korea

<L1>

Cat. No.	Flask-type and Capacity (ml)	Spherical Socket	Size Approx., od φ×h (mm)	Ea. Price, ₩
SL.Fla2080	Receiving, 100ml	35/20	φ64×h110	
SL.Fla2081	, 250ml	"	φ85×h140	
SL.Fla2082	, 500ml	"	φ105×h170	
SL.Fla2083	, 1,000ml	"	φ131×h210	
SL.Fla2084	, 2,000ml	"	φ166×h240	
SL.Fla2085 New	, 3,000ml	"	φ185×h260	

DAIHAN™ 3 & 6 Places Multiuse Heating Mantle, for 250~1000ml Flasks, 450°C, with Certi. & Traceability
With Glassware Supporting Rod & Hook, Ideal for Extraction · Reflux · Distillation Applications, Lower Profile
3 & 6구 다용도 히팅맨틀, 추출 · 환류 · 증류용 등에 적합, 250~1000ml 둥근 플라스크용

- Supplied by DAIHAN Korea
- Fit to DURAN, Glassco, PYREX and Other-brand Round Flasks
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Repetitive Extraction (Kjeldahl, Soxhlet and other Extraction Procedure), Reflux and Distillation Laboratory Applications

* Top Cover ; Stainless-steel 304, Body ; Powder Coated Steel

* Analog Controller Built-in System


* Rods and Hooks for Rods Connection

* Nickel Chrome Heating Element - Coil Type

* Optional Clamp & Clamp Holder for Glassware ; Order Separately

~450°C

3-Year Warranty



6 places Heating Mantle
(WHM12293)



3 places Heating Mantle
(WHM12391)

<L1>

Working Places	Catalog No.		Round Flask 250~1000ml			Watts(W)	Inside Depth (mm)	Dimension w×d×h (mm)	Ea. Price, ₩
	230V	120V	Capa. (ml)	Dia. Max. ϕ (mm)	Bottom Style				
6 places	DH.WHM12291	DH.WHM12292	250/300	83 ~ 86	Flat	840	50	800×240×110	
6 places	DH.WHM12293	DH.WHM12294	250/300	83 ~ 86	Round	840	50	800×240×110	
6 places	DH.WHM12295	DH.WHM12296	500	102 ~ 105	Round	1440	65	950×270×130	
3 places	DH.WHM12391	DH.WHM12392	500	102 ~ 105	Round	720	65	550×270×130	
3 places	DH.WHM12393	DH.WHM12394	1,000	130 ~ 133	Round	1050	75	690×340×170	

Spare Support Rod & Hook Set

- DH.WHM30500 Support Rod & Hook Set, "MHR01", for 3Places 500ml, 2×750-/1×500-mm Rod and 2×Hook for Rod
 DH.WHM31000 "MHR02", for 3Places 1000ml, 2×750-/1×600-mm Rod and 2×Hook for Rod
 DH.WHM60250 "MHR03", for 6Places 250-/300-ml, 2×750-/1×900-mm Rod and 2×Hook for Rod
 DH.WHM60500 "MHR04", for 6Places 500ml, 2×750-/1×940-mm Rod and 2×Hook for Rod

Ea. Price, ₩




※ Optional Clamps & Clamp Holders
⇒ Visit DAIHAN-Sci.com Shopping Mall

DAIHAN™ Extraction Water Bath “WEB”, with Concentric Ring Cover Sets & Holding Frame, 4 · 6 · 8 Place
Ideal for Gentle Extraction · Soxhlet Extraction · COD Determination, Digital Fuzzy Control, Back Light LCD, Up to 100℃, ±0.2℃
다용도 추출용 항온수조, 4 · 6 · 8 구, 컨센트릭 링 커버 세트 & 홀딩 프레임 포함



3-Year Warranty

4-Place Extraction Water Bath “WEB-4”



3-Year Warranty

6-Place Extraction Water Bath “WEB-6”



(1) Lid with Concentric Ring Cover Set (Included)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
 - with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
 - Digital Fuzzy Control System Implementing Superior Temperature Accuracy
 - Concentric Ring Cover Sets make Effective Steaming Work
 - Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function
- * 4 · 6 · 8 Extraction Place are available.
* Includes Concentric Ring Cover Set & Holding Frame
* User's Self-compensation Function to Control Difference between Real Temperature and Indicated Temperature (±10.0℃).
* High Quality 304 Stainless-steel Bath and Holding Frame for Superior Durability.
* Storage Function for the Set Values of Temperature and Timer.
* Alarm Function ; Error Status and Timer-end.
* Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled).

Specifications & Ordering information

<L1>

Spec.	Model	“WEB-4” 4 Place	“WEB-6” 6 Place	“WEB-8” 8 Place
Flask Capacity		250 ~ 500ml almost all Round Flasks		
No. of Extraction Holes		4 Place	6 Place	8(4×2) Place
Temp. Range & Accuracy		Room Temp +5℃ to 100℃, ±0.2℃		
Temp. Uniformity & Sensor		±1.0℃, PT100		
Heating Power		1.4 kW	1.8 kW	1.4 kW×2
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)		
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection, Leakage Breaker		
Material		Internal ; Stainless-steel (#304), External ; Powder Coated steel Holding Frame ; Stainless-steel		
Lid		Concentric Ring Cover Set Included		
Top Opening Dia.		φ36 · φ65 · φ96 · φ110mm		
Bath Depth		122mm		
Drain Valve		Mounted		
Dimension (w×d×h)	Internal	610×160×h160mm	910×140×h160mm	610×290×h160mm
	External	680×245×h275mm	975×245×h275mm	720×375×h275mm
Packing Size (mm) and Gross Weight (kg)		1,100×350×h400, 18 kg	1,250×350×h400, 25 kg	900×450×h400, 34 kg
Electrical Reqt.		1 Phase, AC 120V, 60Hz or 230V, 50 / 60Hz		
Cat. No.	230V	DH.WEB01004	DH.WEB01006	DH.WEB01008
	120V	DH.WEB02004	DH.WEB02006	DH.WEB02008
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Lid with Concentric Ring Cover Set and Holding Frame Set

- DH.WEB11000 (1) Lid with Concentric Ring Cover Set, “EBC100”, for all “WEB” models
- DH.WEB21004 Holding Frame Set, “EBF100”, for “WEB-4”, 2ea×800mm Pole & 2ea×700mm Pole
- DH.WEB21006 Holding Frame Set, “EBF200”, for “WEB-6”, 2ea×800mm Pole & 2ea×1,020mm Pole
- DH.WEB21008 Holding Frame Set, “EBF300”, for “WEB-8”, 4ea×800mm Pole & 4ea×740mm Pole

Ea. Price, ₩



DAIHAN Scientific

FOOD ANALYZERS **NEW**

Food Fat Analyzer	P. 126
Food Fiber Analyzer	P. 127
Kjeldahl Digestion Furnace System for Food Nitrogen Analyzer	P. 128
Kjeldahl Nitrogen Analyzers of Food	P. 129~130



Food Fat Analyzer
“FFA-6”



Food Crude Fiber Analyzer
“FFI-6”



**Kjeldahl Digestion Furnace System
for Food Nitrogen Analyzer**
“FKDN”



Kjeldahl Nitrogen Analyzer of Food
“FKN-M”

DAIHAN™ Food Fat Auto-Analyzer “FFA-6”, Fully Sealed Electric Heating, 6-Sample Simultaneous Test 6-Stage Program Heating, 4-Stage Extraction System, “Non-Stopcock Design”, 식품지방 자동 분석기, 밀폐형 통합시스템



- Supplied by DAIHAN Korea
- Soxhlet 추출 방법을 기반으로 식품, 사료, 환경 제품 및 기타 광범위한 산업 분야에서 지방 함량 분석에 사용.
- 건조, 분쇄, 탈수(황산나트륨), 침전 등의 시료 전처리 필요.
- 자동 온도조절 방식의 완전 밀폐형 전기 가열 시스템을 채택하여 균일한 가열과 안전성 및 신뢰성 보장.
- A Crude Fat Analyzer Measures the Crude Fat Contents in Food · Feed · Grains, and Environmental Samples, Base on the Principles of Soxhlet Extraction.
- Sample Preparation Process such as Drying, Grinding, Dehydration and Sedimentation are Required.
- Integrates Immersion, Extraction, Condensation, and Solvent Recovery Process into a Single Unit.
- Supports Common Extraction Solvent (e.g., Diethyl Ether) and Recovers the Solvent for Reuse during the Solvent Recovery Stage.
- Fully Sealed Electric Heating System ensures Uniform Heating, Safety, and Reliability.
- Test 1-6 Samples Simultaneously without Leakage of Waste.
- Temp can be Selected Based on the Boiling Point of the Reagents and Ambient Temp during Operation, Achieving Rapid Results.
- Reagents can be Recycled, Reducing Testing Costs.
- * 2.5 inch Blue LCD Screen Simultaneously Displays Set Temp, Actual Temp, Countdown Timer for Holding Time, and Alarm for Real-Time Operation Monitoring.
- * Non-Stopcock Design ; Avoids Leakage and Glass Tube Damage Caused by Frequent Stopcock Rotation, and Improves Sealing and Service Life.
- * Automatic 6-Stage Heating ; Six Temp and Holding Time Step to Run Automatically in the Preset Order.
- * Over-Temperature Protection ; Automatically Alarms and Stops Heating when the Actual Temp Exceeds the Upper Limit.
- * The Integrated Aluminum Heating Block with Insulation Provides Uniform Heating.
- * Lift Handle & Cup Holder Enable Efficient Batch Transport and Processing of Multiple Samples.
- ▶ Complete Set with ; (1) Fat Analyzer, (2) Extraction Tube, (3) Aluminum Cup, (4) Cup Holder, (5) Filter Paper Tube, (6) Filter Paper Tube Rod, (7) Rubber Tube, (8) White Rubber Tube, (9) Power Code, (10) Fuse, (11) Instruction Manual, (12) Inspection Certificate

Fat Analysis 분석과정






Fig. Cup Holder Enables Efficient Batch Transport and Processing of Multiple Samples.

Accessories



Specifications & Ordering information

<L1>

Spec.	Model	“FFA-6”
Measuring Range		Crude Fat ≥ 0.1%
Sample Mass		Solids ; 1~10g (Particle Size ≤ 1mm)
Sample Capacity		6 Samples/Batch
Precision		Repeatability ≤ 0.8%
Temp. Range & Accuracy		R/T ~200℃, ±0.3℃
Heating		6 Stages of Heating, Drying
Heating Module		Well Aluminum Heating Block (with Insulation)
Measuring Time		Average 3 Hours
Solvent Recovery Rate		≥ 85%
Recovery Rate		≥ 99% (Cooling Water Temp > 4℃, ≤ 26℃)
Heat Characteristic		Over-Temperature Protection
Material		Liquid Spraying Coated Steel
Sealing Material		Polytetrafluoroethylene (PTFE)
Display		2.5-Inch Blue-Screen LCD
Extraction Cup		Aluminum, 80ml, id φ45, od φ52×h60mm
Dimension & Weight		560×370×h700mm, 32kg
Packing Size & Gross Weight		710×540×h786mm, 45kg
Power Consumption		500 W
Electrical Reqts.		220V ±10%, 50/60Hz
Cat. No.		DH.FFA6
Price, Each. ₩		

※ Note! Installation Charge is Extra. (설치비 별도)

Spare Accessories for “FFA-6”

- DH.FFA6.A01 (3) Aluminum Cup, 80ml, id φ45, od φ52×h60mm
- DH.FFA6.A02 (5) Filter Paper Tube, Aluminum Alloy, od φ33×h58mm
- DH.FFA6.A03 (6) Filter Paper Tube Rod, Nylon, od φ21.5×h58mm

Optional Filter Paper for “FFA-6”

- C1.F002.125 Qualitative Analysis Filter Paper, Medium Slow, φ125mm Thickness 190μm, Retention 5~8μm, Weight 97g/m²

Ea.	Price, ₩

Pk.	Price, ₩
100ea	

DAIHAN™ Food Crude Fiber Analyzer “FFI-6”, Acid-Base Washing Method, 6-Sample Simultaneous Test

High-Pressure Liquid Injection, Non-Sample Transfer Operation, with Boiling Flask, 식품조섬유 분석기, 산-염기 세척방식, 통합 밀폐형



3-Year Warranty

“FFI-6” Crude Fiber Analyzer

Fig. Crude Fiber Analyzer
with (3) Boiler & Boiling Flasks

(2) Boiler & 3 Boiling Flasks



(3) Glass Crucible




(4) Cleaning Column

■ Supplied by DAIHAN Korea

■ 곡물, 사료, 식품 및 농업 부산물의 조섬유, 리그닌, 셀룰로스, 헤미셀룰로스 함량 측정에 적합.

■ 산-염기 세척방식을 기반으로 처리 및 행균과정을 통합한 일체형 구조로 효율적으로 조섬유 함량 측정.

■ 밀폐형 시스템으로 전기 가열, 시약 주입, 분해, 여과, 행균을 시료 이동 없이 효과적으로 조섬유 추출.



■ Suitable for Determining the Content of Crude Fiber, Lignin, Cellulose, and Hemicellulose in Food · Feed · Grains, and Other Agricultural Byproducts.

■ Uses the Acid-Base Washing Principle, Integrating Acid and Alkali Treatment and Rinsing.

■ Fully Enclosed Method of Electric Heating, Liquid Addition, Digestion, Filtration, and Rinsing to Extract Crude Fiber Content Effectively without Sample Transfer.

* Supports CF(Crude Fiber), NDF(Neutral Detergent Fiber), ADF(Acid Detergent Fiber), and ADL(Lignin) Analysis.

* Heating Power can be Individually Adjusted for each Sample, Ensuring Fast and Stable Reflux for Every Sample.

* Pressure-Injection Liquid Feeding System Delivers Reagents Quickly and Smoothly, Overcoming the Poor Flow Problems in Gravity-Flow Systems.

* Automatic Electric Furnace Preheating Ensures Safety, Energy Efficiency, and Prevents Dry Burning.

* Truly Unobstructed Filtration with High-Performance, Corrosion-Resistant Components.

* Enhanced Sealing Structure with Upper and Lower Teflon Insulation Provides Excellent Leak Prevention and Corrosion Resistance.

* Over-Heating Protection(Boiler) ; When an Over-Heat Condition is Detected, the Instrument Automatically Cuts Off the Heating Power.

* Separate-Type Design Allows Boiling Flasks to be Placed on the Table for Convenient Liquid Addition and Makes the System more User-Oriented.

▶ Complete Set with ; (1) Crude Fiber Analyzer, (2) Boiler & 3 Boiling Flasks,
(3) Glass Crucible (35×60mm), (4) Cleaning Column, (5) Rubber Tube,
(6) Power Code, (7) Fuse, (8) Instruction Manual, (9) Inspection Certificate

Specifications & Ordering information

<L1>

Spec.	Model	“FFI-6”
Measuring Range		0~100%
Sample Mass		0.5~3.0g (Particle Size ≤ 1mm)
Sample Capacity		6 Samples/Batch
Precision		Repeatability ≤ 0.4%
Measuring Time		≤ 90min
Heating Method		Semi-Circular Silicon Carbide Heating Heating Power Individual Adjustment
Sealing Method		Upper and Lower Tetrafluoroethylene(PTFE) Insulation
Liquid Addition Method		Pressure Injection Type
Reagent Heating Method (Boiler)		Automatic Thermostat Adjustable, Up to 100°C Over-Heating Protection, Low-Water Cut-Off
Corrosion Resistance		Resistant to Acid, Alkali, Acetone, etc.
Dimension & Weight		560×460×h660mm, 45kg
Packing Size & Gross Weight		660×660×h666mm, 57kg
Power Consumption		Main Unit 1600 W, Anti-Dry-Burning Boiler 1350 W
Electrical Reqts.		220V ± 10%, 50/60Hz
Cat. No.		DH.FFI6
Price, Each. ₩		

※ Note! Installation Charge is Extra. (설치비 별도)

Spare Accessories for “FFI-6”

	Ea.	Price, ₩
DH.FFI6.A01 (2) Boiler & 3 Boiling Flasks Boiler ; Glass, 542×182×h55mm Boiling Flask ; Glass, 2Lit, φ 160×h163mm		
DH.FFI6.A02 (3) Glass Crucible, φ 35×h60mm		
DH.FFI6.A03 (4) Cleaning Glass Column, Boro-3.3, φ 30×h60mm		

DAIHAN™ Kjeldahl Digestion Furnace System “FKDN” For Food Nitrogen Analyzer, 8 & 14Holes, 250ml Tubes, Max. 600°C
With 6-Stage Temp Control · Vapors Hood Exhaust · Acid Vapors Scrubber “FKDN-GS”, 킬달 분해 장치-식품질소분석기용, 유해가스 중화, 6단 조절



Spare Accessories



Note !

The Tube Containing the Digested Sample is Used for Nitrogen Analyzers (“FKN-M” & “FKN-A”)



II. “FKDN-GS” Acid Vapors Scrubber

3-Year Warranty

- Supplied by DAIHAN Korea
- ▶ Fume Hood 내에서 작업 시에는 “The Acid Vapors Scrubber”는 없어도 됨.
- 프로그램 설정과 6-단계 온도 제어로 간편한 작동, 안전성 및 신뢰성을 제공하는 킬달 분해장치.
- 킬달 분해 공정에서 발생하는 산성 미스트 및 가스를 중화·여과하는 친환경 가스 중화장치 사용.
- Kjeldahl Nitrogen Method-based Sample Digestion System for Protein / Nitrogen Analysis with Programmable Settings and 6-Stage Temp Control for Simple, Safe and Reliable Digestion.
- Environmental Friendly Alkali Neutralization Filtration System for Acid Mist and Acidic Gases Generated during Kjeldahl Digestion
- * Digests Samples Using Sulfuric Acid(H₂SO₄) for Kjeldahl Protein / Nitrogen Analysis.
- * Continuously Adjustable Furnace Temperature with High Accuracy and Stability
- * Over-Temperature Protection ; Automatically Alarms and Cut Off the Heating Power
- * Program Control Provides Both Curve Heating and Linear Heating Modes
- * Heat Isolation between the Control Panel and the Furnace Body Reduces the Impact of the High-Temp Radiation from the Furnace Body on the Control System.
- * Acid Gas Scrubber and Gas Hood Exhaust Eliminates the Need to Place the Instrument in a Fume Hood.

I. Kjeldahl Digestion Furnace “FKDN-Series”

<L1>

Spec.	Model	“FKDN-8” 8-Hole	“FKDN-14” 14-Hole
Measuring Range		0.1 ~ 240.0mg of Nitrogen	
Measurement Quantity / Batch		8 Samples	14 Samples
Measurement Speed		45min/Batch	
Digestion Glass Tube		250ml, Φ42×h300mm, Boro-glass 3.3	
Temp Range & Accuracy		R/T ~600°C, ±1°C	R/T ~450°C, ±1°C
Average Heating Rate		30°C/min	
Temperature Control Methods		Program Control with Both Curve Heating and Linear Heating Modes	
Exhaust Gas Sealing Material		Teflon	
Insulation Method		Separate Air Duct Insulation	
Heating Method		Ceramic Furnace Core	Integrated Aluminum Alloy
Safety		Overheat Protection	
Display		2.5-Inch LCD Screen ; Actual Temp, Holding Time & Count Down	
Dimension & Weight		430×370×h190mm, 14kg	440×360×h170mm, 16kg
Packing Size & Gross Weight		770×500×h450mm, 30kg	770×500×h450mm, 33kg
Power Consumption		2,600 W	1,700 W
Electrical Reqt.s.		220V ±10%, 50/60Hz	
Cat. No.		DH.FKDN8	DH.FKDN14
Price, Each. ₩			

※ Note! Installation Charge is Extra. (설치비 별도)

DH.FKDN.A01	(1) Spare ; Digestion Glass Tube, 250ml, Φ42×h300mm	
DH.FKDN.A02	(2) ; Gas Hood Glass Plug, Φ42×h38mm	

Ea.	Price, ₩

II. Acid Vapors Scrubber “FKDN-GS”

- 킬달 분해 또는 기타 공정에서 발생하는 산성 미스트 및 가스를 제거하는 친환경 가스 중화장치.
- 이중 알칼리 중화 & 여과 시스템을 적용하여 산성가스를 보다 완전하고 안전하게 흡수할 수 있도록 설계됨.
- * Used for Acid Mist and Acidic Vapors Generated during Kjeldahl Digestion or Other Processes.
- * Environmental Friendly Dual Alkali Neutralization Filtration System Provides Excellent Neutralization Effect, Ensuring More Thorough Absorption of Acidic Gases
- * Neutralization Bottle Cap Features One-Button Opening and Locking, Making it Convenient to Use.

Spec.	Model	“FKDN-GS”
Neutralizing Volume		2×1Lit / Bottle
Operating Temp & Humidity		+1°C ~ 40°C, Max Relative Humidity is 80%
Minimum Vacuum Pressure		≤ -0.075Mpa, -0.75 bar
Vacuum Pump Suction		≥ 17 L/min (max.)
Maximum Number of Positions		20 Positions
Maximum Amount of Processing		80 Digestion Tubes (Theoretical Value)
Flow Rate Adjustment		Supports Stepless Speed Control
Dimension & Weight		220×215×h65mm, 3.5kg
Packing Size & Gross Weight		500×280×h300mm, 5.3kg
Power Consumption		≤ 20 W
Electrical Reqt.s.		220V ±10%, 50/60Hz
Cat. No. & Price ₩		DH.FKDN GS ₩ Ea
DH.FKDN GS.A01	Spare ; Regent Bottle & Cap, PP, 1Lit, Φ60×h200mm	₩ Ea



DAIHAN™ Kjeldahl Nitrogen Analyzer of Food “FKN-M”, Microprocessor Control, Electric Steam Generator / Over Pressure Protection With Dual Bottle Holder, Auto-Save the Last Set Parameter, 킬달 식품질소 분석기, 데이터 자동 저장, 절수형



3-Year Warranty

“FKN-M” Kjeldahl Nitrogen Analyzer

Kjeldahl Nitrogen Analyzer 분석과정



(2) Digestion Glass Tube
; Same as The Kjeldahl Digestion Furnace Tube

- Supplied by DAIHAN Korea
- 킬달 분석화법(Kjeldahl Method) 기반으로 유기 및 일부 무기 시료 내의 질소함량을 정량분석.
- 식품, 사료, 화학, 농업 등 다양한 분야에서 곡물, 육류, 유제품, 비료, 토양, 석탄, 약품 등의 질소 분석에 적합.
- 절수형 프로그램 자동 응축 시스템이 냉각수 사용을 자동으로 제어하여 운전 중 물 소비를 줄여줌.
- 마지막으로 설정한 데이터를 자동으로 기억해 서로 다른 조건의 증류에도 바로 적용할 수 있음.
- Designed for Quantitative Analysis of Nitrogen Content in Samples Based on the Kjeldahl Nitrogen Determination Method.
- Used in a Wide Range of Fields such as Food, Feed, Chemicals, Grains, Meat, Dairy Products, Fertilizers, Soil, Coal, and Pharmaceuticals.
- Microcontroller Digital Timing Control ; Alkali Addition and Distillation Times are Digitally Displayed with Countdown for a More Intuitive Operation.
- Program Controlled Water-Saving ; The Condensate Drain is Activated Simultaneously during Sample Testing and Automatically Shuts Off upon Completion for Energy Conservation and Environmental Protection.
- Memorizes the Last Set Data, Suitable for Distillation with Different Requirements.
- * Dual Receiving Bottle Holders for a Better User Experience
- * Soft-Touch Film Buttons for a Comfortable Tactile Response and a Clean, Simple Panel
- * Alarm Reminder upon Completion of Testing
- * Alkali Addition and Distillation can be Performed Independently without Affecting Each Other, Ensuring Convenient Operation.
- * Water Shortage Protection Program for the Evaporator to Prevent Dry Burning and Improve Safety.
- * A Special Anticorrosion Pump Performs Timed Alkali Injections, and the Timed Distillation Section, and Acid · Alkali · Heat-Resistant Splash-Proof Parts Ensures Stable Long-Term Use.
- * Backflow Prevention for Receiving and Digesting Solutions Ensures Test Quality and Safety.
- ▶ Complete Set with ; (1) Nitrogen Analyzer Distillation Device, (2) Digestion Tube, (3) Rubber Tube, (4) Water Tray, (5) Power Code, (6) Fuse, (7) Instruction Manual, (8) Inspection Certificate

Specifications & Ordering information

<L1>

Spec.	Model	“FKN-M”
Operation Mode		Thin-Film Touch Button Operating
Measurement Range		0.1mg~240.0mg of Nitrogen
Sample Mass		Solid 0.3~8.0g (Particle Size ≤ 0.42mm), Liquid 2.0~30.0ml
Accuracy		Relative Deviation ≤ 0.3%
Repeatability		Parallel Deviation ≤ 0.2%
Recovery Rate		≥ 99.5%
Measurement Speed		3~8min/sample
Condensation Water Control		Water-Saving Design with Automatic Program Control
Steam Generator Protection		Dual Water Shortage Protection Overpressure Relief Protection
Alkali Addition Mode		Corrosion-Resistant Pump with Timed Injection
Distillation Mode		Timed Distillation
Splash-Proof Tube Material		Strong Acid & Alkali Resistant Plastic
Display		Digital Display
Dimension & Weight		412×376×h740mm, 25kg
Packing Size & Gross Weight		545×538×h870mm, 30kg
Power Consumption		1,700 W
Electric Reqts.		220V±10%, 50/60Hz
Cat. No.		DH.FKDN
Price, Each. ₩		

※ Note! Installation Charge is Extra. (설치비 별도)

Spare Accessory for “FKN-M”

DH.FKDN.A01 (2) Digestion Tube, Glass, 250ml, $\phi 42 \times h300$ mm

Ea. Price, ₩


DAIHAN™ Automatic Kjeldahl Nitrogen Analyzer of Food “FKN-A”, Fully Automatic Control of the Kjeldahl Process Automatic Lifting of Receiving Bottle Holder, Splash-Proof Tube, 자동 킬달 식품질소 분석기, 일체형 자동공정



3-Year Warranty

“FKN-A” Automatic Kjeldahl Nitrogen Analyzer

Kjeldahl Nitrogen Analyzer 분석과정



(2) Digestion Glass Tube
; Same as The Kjeldahl Digestion
Furnace Tube

- Supplied by DAIHAN Korea
 - 킬달 분석화합법(Kjeldahl Method) 기반으로 유기 및 일부 무기 시료 내의 질소함량을 정량분석.
 - 식품, 사료, 화학, 농업 등 다양한 분야에서 곡물, 육류, 유제품, 비료, 토양, 석탄, 약품 등의 질소 분석에 적합.
 - 부식 방지형 정밀 펌프를 적용한 자동 알칼리 주입 시스템으로 알칼리 용액을 정확한 간격으로 자동 주입하여, 강산과 강알칼리에 의한 부식 위험을 줄이고 시약 주입의 정확성과 장비 내구성을 높임.
 - Designed for Quantitative Analysis of Nitrogen Content in Samples Based on the Kjeldahl Nitrogen Determination Method.
 - Used in a Wide Range of Fields such as Food, Feed, Chemicals, Grains, Meat, Dairy Products, Fertilizers, Soil, Coal, and Pharmaceuticals.
 - Fully Automatic Kjeldahl Workflow Performs Bottle Holder Lifting, Water and Boric-Acid Addition, Distillation, Delayed Alkali Dosing, Tube Rinsing and Fault Self-Checking, and even Allows Alkali to be Added during Operation.
 - 4-Level Electric Steam Generator Safety System, with Dual Water-Shortage Protection, Over-Temperature Protection, and Over-Pressure Relief, Prevents Dry Burning and Pressure-Related Risks and Supports Long-Term Stable Operation.
- * 7-Inch HD Color Touchscreen Allows for Real-Time Control of the Experimental Process.
- * Automatic Lifting System for the Receiving Bottle Holder ; Smooth Movement, Low Noise, and Rinses the Receiving Tube, Making the Human-Machine Interface More User-Friendly.
- * Program Controlled Water-Saving ; The Condensate Drain is Activated Simultaneously during Sample Testing and Automatically Shuts Off upon Completion for Energy Conservation and Environmental Protection.
- * Condensate Monitoring ; An Alarm will Sound and the Instrument will Stop Operating when the Condensate Flow Rate is Low to Protect the Instrument and Ensure Reliable Experimental Results.
- * The Splash-Proof Inlet Pipe is Made of Special High-Performance Acid · Alkali · Heat-Resistant Plastic, Providing Excellent Sealing Performance and Extending its Service Life.
- * Backflow Prevention for Receiving and Digesting Solutions Ensures Test Quality and Safety.
- ▶ Complete Set with ; (1) Nitrogen Analyzer, (2) Digestion Tube, (3) Inlet Pipe, (4) Rubber Tube, (5) Transparent Rubber Hose, (6) Water Tray, (7) Power Code, (8) Fuse, (9) Instruction Manual, (10) Inspection Certificate

Specifications & Ordering information

<L1>

Spec.	Model	“FKN-A”
Operation Mode		Automatic & Manual Operating System
Measurement Range		0.1mg~240.0mg Nitrogen
Sample Mass		Solid ≤ 8.0g (Particle Size ≤ 0.42mm) ; Liquid ≤ 30.0mL
Accuracy		Relative Deviation ≤ 0.3%
Repeatability		Parallel Deviation ≤ 0.2%
Recovery Rate		≥ 99.5%
Measurement Speed		3~8min/sample
Condensation Water Control		Water-Saving Design with Automatic Program Control
Condensate Monitoring		Condensate Flow Rate Monitoring Alarm when below Set Value
Lifting System		Automatic Lifting and Lowering of the Receiving Bottle Holder, and Cleaning of the Receiving Tube
4-Level Steam Generator Protection		Double Protection of Water Shortage, Overpressure Relief Protection, Overtemperature Protection
Alkali Adding Method		Special Anti-Corrosion Pump Regular Injection
Splash-Proof Inlet Pipe		Material ; Strong Acid and Alkali Resistant Plastic
Operation Interface		New UI Dynamic Interface, Easier Operation
Display		7-Inch LCD Ultra-Fine Color Widescreen
Dimension & Weight		440×440×h735mm, 28kg
Packing Size & Gross Weight		570×570×h870mm, 33kg
Power Consumption		1,700 W
Electric Reqts.		220V±10%, 50/60Hz
Cat. No.		DH.FKNA
Price, Each. ₩		

※ Note! Installation Charge is Extra. (설치비 별도)

Spare Accessory for “FKN-A”

DH.FKDN.A01 (2) Digestion Tube, Glass, 250mL, φ42×h300mm

Ea. Price, ₩

DAIHAN™ Freeze Dryer, Benchtop

FREEZE DRY SYSTEM “FRD”

Patented
Jog-Dial
&
Touch-Button
Control

Compact-size
Lab Scale
Bench-top
Type

Fastest
Freezing Time

Ergonomic & Compact Body

- Lab Scale, Bench-top Type Model
- Compact size for all laboratory table
- Ergonomic Design

Company	DAIHAN	Competitor
Chamber Capacity	8 Lit	8 Lit
External Dimension (mm)	460×620×h520	500×646×h976

Various Accessories

- Used with 12-port Manifold or Various Dry Chambers, Various Size Flasks, Flask Top, Vacuum Valve, Vacuum Pumps, Table Cart, etc.



12-port T-type
Manifold



Various
Dry Chambers


Intuitive User Interface

- Patented Jog-Dial Controller(Turn and Push)
- Capacitive Touch Button for user convenience
- Precise control by microprocessor
- Digital Back-Light LCD



Fastest Freezing System

- DAIHAN Freezing System realizes Quick Freezing and Quick Defrosting Function




Freezing System

Others

- High Sublimation Rate ; 3 & 6 Lit / day
- Manual Defrost Function ; Hot Gas Defrost Function for quick next cycle turn-around
- Automatic or Manual Process Control
- Control Temp, Vacuum and Timer at Automatic Mode
- LED Indication Light On in case of suitable Temp and Vacuum
- Can be Used with Various Accessories

Freezing speed comparison with Competing Product



DAIHAN™ 3 & 6Lit/day Freeze Dry System “FRD”, Lab Scale Benchtop-type, Automatic & Manual Process, -90℃
Various Optional Acc.; T-type Manifold, Round · Square · Stoppering-Acrylic Dry Chamber & Stainless-steel Dry Chamber
실험실용 동결 건조기, 컴팩트 바디 사이즈, 터치식 컨트롤러



I. 3Lit/day Freeze Dryer “FRD-8”

※ The Figure ; Installation of Optional 12-Port T-type Manifold and 12 Vacuum Handle Valve



II. 6Lit/day Freeze Dryer “FRD-16”

※ The Figure; Installation of Optional Large Volume Square-type Acrylic Dry Chamber



(1) 12-Port T-type Manifold with 12 Silicone Stoppers



(2)-㉠㉡ 8-Port Round-type Acrylic Dry Chamber (2)-㉢㉣ 3-Stage STS Shelf

- Supplied by DAIHAN Korea
- Patented Jog-Dial & Touch Button Control System
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- PL (Product Liability) Insurance
- Ideal for a Wide Range of Applications Including Pharmaceutical, Biological, Food Industry, Petrochemical and Semiconductor Industry, etc.

■ Various Optional Accessories;

- (1) T-type Manifold with 12-port
- (2) Round-type Acrylic Dry Chamber with 8-port
- (3) Large Volume Square-type Acrylic Dry Chamber
- (4) Heated Stoppering Acrylic Dry Chamber
- (5) Stainless-steel Dry Chamber with 12-port

- * Ergonomic Design and Compact Body Size
- * Allows Automatic or Manual Process Control
- * Rapid Integrated Hot Gas Defrost Function for Quick Next Cycle Turn-around
- * LCD Indicates Temp and Vacuum
- * LED Indication Light On in Case of Suitable Temp and Vacuum
- * Can be Controlled Temp and Timer at Automatic Mode
- * Homogeneous Temp Distribution across all Shelves in Dry Chamber

Specifications & Ordering information

<L1>

Spec.	Model	I. “FRD-8”	II. “FRD-16”
Sublimation Rate(Lit / day)		3Lit / day	6Lit / day
Chamber & Max Condenser Capa.		Chamber ; 8 Lit, Condenser ; 10 Lit	Chamber ; 16 Lit, Condenser ; 20 Lit
Set Temp Range		-90℃~-60℃	
Compressor & Defrost-type		1/2 HP, Manual	3/4 HP, Manual
Display		Digital LCD with Back-Light Function	
Material		Stainless-steel 304	
Dimension(w×d×h) & Weight		Exterior ; 460×620×h520 mm, Condenser ; ϕ300×h150 mm, Net Weight ; 55 kg	Exterior ; 520×680×h600 mm, Condenser ; ϕ300×h300 mm, Net Weight ; 60 kg
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	
Cat. No.	230V	DH.Frd00208	DH.Frd00216
	120V	DH.Frd00108	DH.Frd00116
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Optional Manifold & Dry Chamber

(1) 12-Port T-type Manifold

- DH.Frd01041 12-port T-type Manifold “TMF08”, for “FRD-8” with 12 Silicone Stoppers, for Connecting 3/4” Vacuum Handle Valve
- DH.Frd01042 12-port T-type Manifold “TMF16”, for “FRD-16” with 12 Silicone Stoppers, for Connecting 3/4” Vacuum Handle Valve

Ea. Price, ₩

(2) 8-Port Round-type Acrylic Dry Chamber & 3-Stage Shelf

- DH.Frd01051 ㉠ 8-port Round-type Acrylic Dry Chamber “ADC08”, for “FRD-8” with 8 Silicone Stoppers, for Connecting 3/4” Vacuum Handle Valve, Chamber Size ; Int. ϕ284×h390mm, Ext. ϕ300×h463mm
- DH.Frd01052 ㉡ 8-port Round-type Acrylic Dry Chamber “ADC16”, for “FRD-16” with 8 Silicone Stoppers, for Connecting 3/4” Vacuum Handle Valve, Chamber Size ; Int. ϕ334×h390mm, Ext. ϕ350×h453mm
- DH.Frd01061 ㉢ 3-stage Stainless-steel(STS) Shelf “SHF08”, for “FRD-8” ϕ259×11mm, Total h225mm
- DH.Frd01062 ㉣ 3-stage Stainless-steel(STS) Shelf “SHF16”, for “FRD-16” ϕ310×11mm, Total h225mm

Ea. Price, ₩

<Continued On Next Page>

<Continued From...DAIHAN™ 3 & 6Lit/day Freeze Dry System "FRD">

<L1>



(3) Large Volume Square-type Acrylic Dry Chamber With 5-Stage STS Shelf



(4) Heated Stoppering Acrylic Dry Chamber with 2-Stage STS Shelf



(5) 12-Port STS Dry Chamber



(6) Vacuum Handle Valve



Tube Adapter for Connecting Flask Top



(10) Silicone Flask Top



(11) Glass Flask



(12) Vacuum Pump



(13) 2-Stage Table Cart

(3) Large Volume Square-type Acrylic Dry Chamber with 5-Stage Shelf

- DH.Frd02100** Square-type Acrylic Dry Chamber "SADC08", for "FRD-8" with 5-stage STS shelf ; 375×375×h25mm
Chamber Size ; Int. 380×410×h390mm, Ext. 510×480×h490mm
- DH.Frd02101** Square-type Acrylic Dry Chamber "SADC16", for "FRD-16" with 5-stage STS shelf ; 425×450×h25mm
Chamber Size ; Int. 435×485×h420mm, Ext. 555×540×h475mm

(4) Heated Stoppering Acrylic Dry Chamber, with Height Adjustable 2-Stage Shelf and Analog Temp Controller Package-set

- DH.Frd07001** Heated Stoppering Acrylic Dry Chamber "HADC08", for "FRD-8" with Height Adjustable 2-stage Shelf ; $\phi 255 \times 13\text{mm}$, Total h278mm
Chamber Size ; Int. $\phi 284 \times h320\text{mm}$, Ext. $\phi 300 \times h523\text{mm}$
- DH.Frd07002** Heated Stoppering Acrylic Dry Chamber "HADC16", for "FRD-16" with Height Adjustable 2-stage Shelf ; $\phi 305 \times 13\text{mm}$, Total h278mm
Chamber Size ; Int. $\phi 334 \times h320\text{mm}$, Ext. $\phi 350 \times h523\text{mm}$,
- DH.Frd09001** ② Analog Temp Controller for Heated Stoppering Acrylic Dry Chamber, 230V, Temp Range 0°C~40°C

(5) Stainless-steel Dry Chamber & 3-Stage Shelf

- DH.Frd08001** 12-port Stainless-steel Dry Chamber "SDC08", for "FRD-8"
Chamber Size ; Int. $\phi 296 \times h350\text{mm}$, Ext. $\phi 320 \times h370\text{mm}$
- DH.Frd08002** 12-port Stainless-steel Dry Chamber "SDC16", for "FRD-16"
Chamber Size ; Int. $\phi 346 \times h350\text{mm}$, Ext. $\phi 370 \times h370\text{mm}$
- DH.Frd01061** ② 3-stage Stainless-steel(STS) Shelf "SHF08", for "FRD-8"
 $\phi 259 \times 11\text{mm}$, Total h225mm
- DH.Frd01062** ② 3-stage Stainless-steel(STS) Shelf "SHF16", for "FRD-16"
 $\phi 310 \times 11\text{mm}$, Total h225mm

Optional Vacuum Handle Valve & Tube Adapter for Connecting Flask

- DH.Frd03010** (6) 3/4" Vacuum Handle Valve, Neoprene Rubber & PP, for Connecting Flask
- DH.GL28091** (7) 3/4" Glass Tube Adapter, $\phi 19/1.8\text{t} \times \text{L}80\text{mm}$ for Connecting Middle & Large-size Flask Top
- DH.GL28092** (8) 3/4" Curved Glass Tube Adapter, $\phi 19/1.8\text{t} \times \text{L}80\text{mm}$, for Connecting Middle & Large-size Flask Top
- DH.Frd04010** (9) 3/4" PP Tube Adapter, $\phi 19/1.8\text{t} \times \text{L}80\text{mm}$, for Connecting Small-size Flask Top

Optional DURAN® Borosilicate α 3.3 Glass Flask & Silicone Flask Top

- DH.Frd05010** (10) Silicone Flask Top, 3/4" Small-size Flask Top, for 40-/80-ml Flask
- DH.Frd05020** , 3/4" Middle-size Flask Top, for 200-/300-ml Flask
- DH.Frd05030** , 3/4" Large-size Flask Top, for 500-/700-/1000-ml Flask
- DH.GL12041** (11) Glass Flask, DURAN® Glass α 3.3, 40ml , $\phi 40/2.3\text{t} \times h80\text{mm}$
- DH.GL12042** , 80ml , $\phi 40/2.3\text{t} \times h120\text{mm}$
- DH.GL12043** , 200ml , $\phi 65/2.2\text{t} \times h100\text{mm}$
- DH.GL12044** , 300ml , $\phi 65/2.2\text{t} \times h140\text{mm}$
- DH.GL12045** , 500ml , $\phi 90/2.5\text{t} \times h125\text{mm}$
- DH.GL12046** , 700ml , $\phi 90/2.5\text{t} \times h165\text{mm}$
- DH.GL12047** , 1000ml , $\phi 90/2.5\text{t} \times h225\text{mm}$

Optional Vacuum Pump & 2-Stage Table Cart for Vacuum Pump

- V3.VOP100** (12) Vacuum Pump, "VOP-100", 120(100)Lit/min, NW25, 220V
- V3.VOP150** , "VOP-150", 180(150)Lit/min, NW25, 220V
- DH.Frd06025** NW25 Vacuum Pump Connection-set, Tube · O-ring · Clamp
- DH.Frd06016** NW16 Vacuum Pump Connection-set, Tube · O-ring · Clamp · Adapter
- DH.Ca1002** (13) 2-stage Table Cart for Freeze Dryer & Vacuum Pump, with Stop-on Casters and Vacuum Pump Oil Tray, w504×d604×h780mm

DAIHAN Scientific

FUME HOODS

Ducted Fume Hoods	P. 137~138
Acids/Corrosives Resistance, PP	P. 137
Multiuse, Aluminum Frame	P. 138
Ductless Filtering Fume Hoods	P. 135~136
Filtering Purifiers	P. 139~141
Dust Collector	P. 141
Fume Multi Gas Purifiers	P. 139
Toxic Gas Multi-Purifiers	P. 140
Vapor Purifier Box	P. 141



Ductless Filtering Fume Hood
“FHF”



PP Ducted Fume Hood
“HD”



Ducted Multiuse Fume Hood
“FHD”



Filtering Fume Multi Gas Purifier
“CAFU-19N”

DAIHAN™ **Premium** Sturdily Fabricated Filtering Utility-Fume Hood “FHF”, For Acids · VOCs · Dust · Smell & Bacteria, w900 · 1200 · 1500mm Ductless, Safely Fabricated Body & Filters, 5-Filters of 3 μ m Pre · 0.3 μ m HEPA · 2 × Coated Activated Carbon · Activated Carbon With 3-Side Clear Viewing, 11° Slanted Front Sash Door of 2 Step Opening, and 2 × ϕ 220mm Holes on Lower Door for Easy Hand Work, **5-Years Warranty** 필터식 복합 흡후드/조립식 Portable 만능 흡후드-산·암모니아·휘발성 유기화합물(VOCs)·포름알데히드·미세먼지·악취·세균 등을 종합 정화, 필터 교체주기 알람, 무덕트



Ductless Filtering Benchtop ; 조립식 흡후드

- 좌·우·전면 투명 내약품 PVC Panel.
- 11° 경사의 전면 2단 투명 개·폐 Window
- 하부문에 2개의 작업 Hole(ϕ 220mm)이 있어 편안하고 안전한 실험

Re. Approx Size of Particulate Matter

- PM2.5 초미세먼지 ; 1.0 ~ 2.5 μ m Particles
- PM1.0 극초미세먼지 ; $\leq$ 1.0 μ m Particles
- Fume(Acids) ; $\leq$ 1.0 μ m (이하) Particles
- Bacteria ; 0.5 μ m ~ 0.5m/m Particles
- Virus ; 0.018 ~ 0.4 μ m Particles
- VOCs ; Gas (기체)




Fig. 11° Slanted 2-Step Folding Window & Optional Stand

견고한 구조의 조립식 Portable 필터형 흡후드, 무덕트 ; Optional Steel Frame Stand!

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- ▶ 5-Year Warranty of Safely Sturdy Fabrication and Replaced Filter System 조립기기 안전성 5년 보증.
- ▶ Filter Replacement Cycle Alarm Function ; 필터 교체 시기 알람 ALARM 내장 !
 - BLDC 모터 사용 시간을 카운트하여 필터 교체 알람(6개월) - 필터교체 주기는 약 6~12개월로
- Compliance with Laboratory Fume Hood Standard ASHRAE 110-1995
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Safe and Easy to move Workspace for Harmful Gas without Duct Facility
- ▶ Movable Benchtop. Ensuring High Performance Exhaust Filtering without Duct
- * 8 Levels Air Velocity for Various Test Application with Powerful BLDC Fan Motor
- * Operation Time Monitoring Function For Filter Replacement.
- ▶ 5×Framed Filters - Interchangeable System ; 2-Type of I. Standard Hood Unit with 3 μ m Prefilter(Polyester) Only, and II. Complete Hood-Attached 5 Filters All.

Standard Specifications

※ 운임 별도 <L1>

Spec.	Model	“FHF-09” w900	“FHF-12” w1200	“FHF-15” w1500
Working Space w×d×h (mm), Lit.		884×580×h800mm, 410Lit	1184×580×h800mm, 549Lit	1484×580×h800mm, 688Lit
Air Velocity(m/s)		0.4m/s ~ 1m/s, 8 Levels		
Motor		70W, BLDC Motor		
5×Framed Filters		- (1) 3μm Prefilter(Polyester) Included ; Protects the Main Filters From Dust. - Optional 4×Filters of (2) 0.3μm HEPA-H14, (3) Acids, (4) Ammonia+Formaldehyde, and (5) Organics(VOCs)+Formaldehyde.		
Material		Head ; Powder-coated Steel, Working Plate ; Stainless-steel #304 Chemical Resistant Clear PVC(Polyvinyl Chloride) - 8mm Thick.		
Clear PVC Front Window Working Holes, ϕmm		- 11° Slanted Sash Door of 2-Steps Opening(Clear PVC)-Optimal Visibility. - 2×ϕ220mm Holes on Lower Door for Easy Hand Work.		
Noise Level(dB)		≤ 60dB		
Controller		Digital PID Controller with Touch-button		
Internal Lighting		17W LED (600Lux)	29W LED (600Lux)	
Utility Device		Air Velocity Meter		
Utility Device Hole		2×ϕ50mm ON-OFF Hole(PC) ; 기타장비의 코드선 등의 In-Out 용		
Dimension(mm)	Outer	900×610×h1180	1200×610×h1180	1500×610×h1180
Without Stand	Inner	884×580×h800	1184×580×h800	1484×580×h800
Net Weight		82kg	103kg	122kg
Packing Size & Gross Weight	Head	1160×780×h655, 54kg	1760×780×h655, 61.5kg	1760×780×h655, 69kg
	Wall & Plate	1020×925×h200, 31kg	1620×925×h200, 51kg	1620×925×h200, 61kg
Electric Reqts.		1 Phase, AC 230V, 50/60Hz		

▶ Optional Frame Stand for Hoods, Powder Coated Steel.

- DH.FHF09S Stand, For w900 , 900×610×h730mm , Packing ; 1350×805×h225mm, 25kg
 DH.FHF12S , For w1200, 1200×610×h730mm, Packing ; 1350×805×h225mm, 35kg
 DH.FHF15S , For w1500, 1500×610×h730mm, Packing ; 1550×805×h225mm, 45kg

Ea. Price, ₩

I. Standard Hood Unit, Attached Pre-Filter Only. **ECONOMICAL OPERATION**

경제적 운영 ; 실험 여건상 필요한 필터만을 구입 사용, Stand 별도

- DH.FHF09 Fume Hood Unit, with Prefilter, “FHF-09”, w900 Model.
 DH.FHF12 , “FHF-12”, w1200 Model.
 DH.FHF15 , “FHF-15”, w1500 Model.

Ea. Price, ₩

<Continued On Next Page>

<Continued From...DAIHAN™ Premium Sturdily Fabricated Filtering Utility-Fume Hood “FHF”>

Fabricated Body & Filters 조립형 본체와 교체형 필터

► Filters & Air Flow

Max. 5 Framed Filters
Interchangeable

- (5) VOCs / CH₂O Filter
- (4) NH₃ / CH₂O Filter
- (3) Acids Filter
- (2) HEPA-H14 Filter
- (1) Polyester Pre-filter

2-Step Front Window
; Clear PVC

Outside Air
Inflow

Filtered
Clean Air

Fume/Gas




Shape of Framed Filters



Inspection Certificate of Filters
MADI (주)마디 CERTIFIED

► Tool-free Fabrication-Easy & Sturdy 별도의 도구 없이 쉽고 튼튼한 조립 DIY



II. Complete Hood, Attached 5 Filters All. The Best SAFETY WORKSTATION

가장 안전 후드 ; 1대의 Hood로 모든 유해물질의 Vapors · 먼지 · 세균 등을 정화, Stand 별도

5 Framed Filters ;

- (1) 3 μ m Pre-Filter(Polyester) ; Protects the Main Filters From Dust.
- (2) 0.3 μ m HEPA-H14 Filter ; 0.3 μ m 이상의 모든 유해성 Particles를 99.995% 여과 정화.
- (3) Chemical Coated Activated Carbon Filter ; Acids의 Fume/Vapors를 Adsorption & Filtering 흡착여과, 2nm Molecule
- (4) Chemical Coated Activated Carbon Filter ; Ammonia & Formaldehyde의 Fume/Vapors를 Adsorption & Filtering 흡착여과, 2nm Molecule
- (5) Activated Carbon Filter ; Organics(VOCs) & Formaldehyde Gas를 Adsorption & Filtering 흡착여과, 2nm Molecule

DH.FHF09C Fume Hood Complete with 5-Filters, “FHF-09C”, w900 Model.
DH.FHF12C “FHF-12C”, w1200 Model.
DH.FHF15C “FHF-15C”, w1500 Model.

<L1>
Ea. Price, ₩

◆ Selection Guide of Framed Filters

Filter	Target	Material	Method	Description
(1) Polyester Pre-Filter	Particles	Synthetic Micro Fiber	Physical	≥3 μ m Particles을 필터링 하여 메인필터를 보호하여 필터수명을 증대시킴
(2) HEPA-H14 Filter	Particles	Ultra Micro Glass Fiber	Physical	≥0.3 μ m 모든입자(세균·먼지·흙 등)를 99.995% 여과 정화.
(3) Acids Filter	Fume Vapors	Chemical-coated Activated Carbon	Chemical	특정한 화학물질이 코팅된 활성탄에 Fume/Vapor(Acids, NH ₃ , S, etc.)나 먼지 등을 흡착 여과 정화, 2nm Molecule
(4) Ammonia NH ₃ + Formaldehyde CH ₂ O Filter	Fume Vapors	Chemical-coated Activated Carbon	Chemical	특정한 화학물질이 코팅된 활성탄에 암모니아, 아민, 포름알데하이드 등의 Fume/Vapor를 흡착 여과 정화. 2nm Molecule
(5) VOCs + Formaldehyde CH ₂ O Filter	Gas	Activated Carbon	Chemical Physical	활성탄에 VOCs(휘발성 유기화합물의 총칭) CH ₂ O Gas나 먼지 등을 흡착 여과 정화. 2nm Molecule

■ Size of Framed Filters ; * Prefilters Thick-10mm, HEPA-H14 Filter Thick-80mm

For w900 ; 460×410×T40mm For w1200 ; 760×410×T40mm For w1500 ; 1060×410×T40mm

Ordering Information for Framed Filters

Framed Filter Description	For w900 Models		For w1200 Models		For w1500 Models	
	Cat. No.	Price, ₩	Cat. No.	Price, ₩	Cat. No.	Price, ₩
(1) 3 μ m Prefilter ; Particles	DH.FHF09.F1		DH.FHF12.F1		DH.FHF15.F1	
(2) HEPA-H14 Filter ; ≥0.3 Particles 99.995% 정화	DH.FHF09.F2		DH.FHF12.F2		DH.FHF15.F2	
(3) Acids Filter ; Chemical-coated Activated Carbon	DH.FHF09.F3		DH.FHF12.F3		DH.FHF15.F3	
(4) Ammonia / CH ₂ O Filter ; Chemical-coated Activated Carbon	DH.FHF09.F4		DH.FHF12.F4		DH.FHF15.F4	
(5) VOCs / CH ₂ O Filter ; Activated Carbon	DH.FHF09.F5		DH.FHF12.F5		DH.FHF15.F5	

► Easily Quick Replacement of Framed Filters 쉽고 신속한 필터교체 DIY

* Hood의 전면 상부 패널을 열어 직접 손쉽게 필터를 교체 가능



DAIHAN™ **Premium** PP Ducted Acids/Corrosives Safety Fume Hood “HD”, 4 Models of w1,200 · 1,500 · 1,800 · 2,100mm
ALL PP Strong Double Wall, I. Bypass and II. Air Curtain-Type, with Air · Gas · Water-Cock, Cup Sink, Drain, and Explosion Proof Lamp
덕트형 PP 프리미엄 내부식형 흡후드, 강산(Conc Acids) 및 강부식형 휘발성 유해 용액 시험에 적합, 일반 배기형 or 에어커튼형의 선택성, 내약품성 & 내열성 우수 125℃



3-Year Warranty

- Supplied by DAIHAN Korea ■ Durable, Hazardous Chemical Gas Remove.
- PP(Polypropylene) Material for Excellent Resistance to Acid · Alkali · Chemicals · High Temp. 125℃
- ALL PP Strong Double Wall Structure ; Exterior and Interior ; 3mm thick Plate and 10mm thick Front Part, 0~125℃ Stable.
- * With Air-/Gas-/Water-Cock, Cup Sink, Drain, Explosion Proof Lamp, Receptacle, and Multi-purpose Cabinet.
- * **By Pass Air Flow-type** ; After the Suction Forcing the Indoor Air to Discharge the Gas to the Outside
- * **Air Curtain Air Flow-type** ; The Outdoor Air Sucked Form the Air Curtain to the Front Chamber, and the Exhaust Gas Back to Suction the Air
- * Discharging the Polluted Air and Toxic Gas Generated during the Experiment in a Laboratory to the Outside through the Duct Structure
- * Process Pipe & Duct Silicone Sealing to Completely Block the Gas Leakage to the Outside
- * The Mechanical Design is Manufactured with Superior Workability


Specification

Models & Air Flow-type	I. By-Pass Air Flow Type			
	“HD-1200”	“HD-1500”	“HD-1800”	“HD-2100”
	II. Air Curtain Air Flow Type			
	“HD-C1200”	“HD-C1500”	“HD-C1800”	“HD-C2100”
Work Area (mm)	980×615×h1100	1280×615×h1100	1580×615×h1100	1880×615×h1100
Duct Size	φ 200mm	φ 250mm		
Material*	- PP Double Wall 8/10mm T. ; (1) External · (3) Internal · (4) Work Top - PVC Clear Window ; (6) Column · (2) Sash Frame			
Frame	(6) PVC Stability Column Frame			
Cabinet-Lower	(5) Built-in Multi-purpose PP Cabinet-Lower Part			
Lamp	LED Lamp(1ea×20 W)			
Cock	1ea×Air-/Gas-/Water-Cock			
Receptacle	4ea 220V Concent or 4ea 110V Concent			
Dimension (w×d×h) (mm)	1200×800×2200	1500×800×2200	1800×800×2200	2100×800×2200
Weight (kg)	100	130	160	190
Power Supply *	1 Phase, AC 220V / 60Hz			

* Other Specifications(50Hz) are available upon Customer's Request.

※ 운임 별도

Diagram of Air Flow



I. By-Pass Air Flow-type

II. Air Curtain Air Flow-type

Welcome! Custom-Made Fume Hood!!

Ordering Information

I. By-Pass Air Flow-type Ducted PP Fume Hood

DH.WHD02112	Fume Hood, 230V, 1200×800×h2200mm, “HD-1200”
DH.WHD02115	, 1500×800×h2200mm, “HD-1500”
DH.WHD02118	, 1800×800×h2200mm, “HD-1800”
DH.WHD02121	, 2100×800×h2200mm, “HD-2100” New
DH.WHD01112	Fume Hood, 120V, 1200×800×h2200mm, “HD-1200”
DH.WHD01115	, 1500×800×h2200mm, “HD-1500”
DH.WHD01118	, 1800×800×h2200mm, “HD-1800”
DH.WHD01121	, 2100×800×h2200mm, “HD-2100” New

II. Air Curtain Air Flow-type Ducted PP Fume Hood

DH.WHD02212	Fume Hood, 230V, 1200×800×h2200mm, “HD-C1200”
DH.WHD02215	, 1500×800×h2200mm, “HD-C1500”
DH.WHD02218	, 1800×800×h2200mm, “HD-C1800”
DH.WHD02221	, 2100×800×h2200mm, “HD-C2100” New
DH.WHD01212	Fume Hood, 120V, 1200×800×h2200mm, “HD-C1200”
DH.WHD01215	, 1500×800×h2200mm, “HD-C1500”
DH.WHD01218	, 1800×800×h2200mm, “HD-C1800”
DH.WHD01221	, 2100×800×h2200mm, “HD-C2100” New

- ※ Custom-Made Fume Hood ; Wide Variety of Material for User Selection.
- ※ Note! Freight and Installation Charges Separately.

<L4>

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ Ducted Multiuse Fume Hood “FHD”, made of PVC & Bakelite Plate with Aluminum Frame, w1,200 · 1,500 · 1,800 · 2,100mm Usable for Acids/Chemicals, I. By-Pass Air and II. Air-Curtain Flow Types, with Air/Gas/Water-Cock, Cup Sink, Drain, and Anti-Explo Light Lamp 덕트형 다용도 흙후드 산용 가능, 견고한 알루미늄 프레임에 내약품성 2중벽 (PVC 외벽+페놀판 내벽 & 작업대+하부 캐비닛)



- Supplied by DAIHAN Korea
- Durable, Acids / Hazardous Chemical Fume / Gas Remove.
- Safety Double Wall Structure ; (1) Exterior - 4mm thick PVC Plate and (2) Interior - 4mm thick Bakelite Plate (Phenol Formaldehyde Hard Board)
- Sash-type Tempered Safety Glass Window. 강화 내열 내충격 유리 윈도우.
- * General Chemical Gas Hood with High Thermal, Chemical Resistance, Various Material for User selection.
- * With Air-/Gas-/Water-Cock, Cup Sink, Drain, Explosion Proof Lamp, Receptacle, Multi-purpose Cabinet.
- * By Pass Air Flow-type ; After the Suction Forcing the Indoor Air to Discharge the Gas to the Outside
- * Air Curtain Air Flow-type ; The Outdoor Air Sucked Form the Air Curtain to the Front Chamber, and the Exhaust Gas Back to Suction the Air
- * Discharging the Polluted Air and Toxic Gas Generated during the Experiment in a Laboratory to the Outside through the Duct Structure
- * Process Pipe & Duct Silicone Sealing to Completely Block the Gas Leakage to the Outside


Specification

Models & 2 Air Flow Type	I. By-Pass Air Flow Type			
	“FHD-12” w1200mm	“FHD-15” w1500mm	“FHD-18” w1800mm	“FHD-21” w2100mm
	II. Air Curtain Air Flow Type			
	“FHD-12C” w1200mm	“FHD-15C” w1500mm	“FHD-18C” w1800mm	“FHD-21C” w2100mm
Work Area (mm)	990 × 650 × h1200	1290 × 650 × h1200	1590 × 650 × h1200	1890 × 650 × h1200
Duct Size	φ 200mm			
Material* Double Wall	(1) External ; PVC 4T, (2) Internal ; Bakelite(Phenol) Plate 4T and (4) Work Top ; Phenol Plate			
Frame	(6) Durable Aluminum			
Sash Frame & Window	(3) Aluminum Framed Tempered Safety Glass Sash Window			
Cabinet-Lower	(5) Powder Coated Steel Cabinet with Wood Door			
Lamp / Light	LED Lamp(1ea × 20 W)			
Cock	1ea × Air-/Gas-/Water-Cock			
Receptacle	2ea 220V Concent or 2ea 110V Concent			
Dimension (w × d × h) (mm)	1200 × 800 × 2350	1500 × 800 × 2350	1800 × 800 × 2350	2100 × 800 × 2350
Weight (kg)	150	170	200	240
Power Supply *	1 Phase, AC 220V / 60Hz			

* Other Specifications(50Hz) are available upon Customer's Request.

※ 운임 별도

Diagram of Air Flow



I. By-Pass Air Flow-type

II. Air Curtain Air Flow-type

Welcome! Custom-Made Fume Hood!!

Ordering Information

I. By-Pass Air Flow Type

DH.FHD12	Fume Hood, 230V, w1200 Model, “FHD-12”
DH.FHD15	, w1500 Model, “FHD-15”
DH.FHD18	, w1800 Model, “FHD-18”
DH.FHD21	, w2100 Model, “FHD-21”
DH.FHD121	Fume Hood, 120V, w1200 Model, “FHD-121”
DH.FHD151	, w1500 Model, “FHD-151”
DH.FHD181	, w1800 Model, “FHD-181”
DH.FHD211	, w2100 Model, “FHD-211”

II. Air Curtain Flow Type

DH.FHD12C	Fume Hood, Curtain, 230V, w1200 Model, “FHD-12C”
DH.FHD15C	, w1500 Model, “FHD-15C”
DH.FHD18C	, w1800 Model, “FHD-18C”
DH.FHD21C	, w2100 Model, “FHD-21C”
DH.FHD121C	Fume Hood, Curtain, 120V, w1200 Model, “FHD-121C”
DH.FHD151C	, w1500 Model, “FHD-151C”
DH.FHD181C	, w1800 Model, “FHD-181C”
DH.FHD211C	, w2100 Model, “FHD-211C”

※ Custom-Made Fume Hood ; Wide Variety of Material for User Selection.

※ Note! Freight and Installation Charges Separately.


<L4>

Ea. Price, ₩

Ea. Price, ₩

CAFU™ Filtering Fume Multi Gas Purifier of Acids · Dust · Organics(VOCs) · Smell · Bacteria, $\geq 0.3\mu\text{m}$, 100 & 430 CMH(m³/h)
4-Layer Filters of 1 μm Micro Bag · Ceramic Carbon · Active Carbon · HEPA-H13, 흡 복합정화기·산·초미세먼지·휘발성 유기화합물·악취·세균·필터식 이동형

4-Layer Perfect Filter



- ① - 1 μm Micro Fiber Bag Filter ;
For CAFU-19N
- 1 μm Fiber Bag(5 Room) Filter ;
For CAFU-A14
- ② Ceramic Honeycomb Carbon Filter
- ③ Impregnated Active Carbon Filter
- ④ 0.3 μm Sterilization HEPA-H13 Filter

Re. Approx Size of Particulate Matter

- PM2.5 초미세먼지 ; 1.0 ~ 2.5 μm Particles
- PM1.0 극초미세먼지 ; $\leq 1.0\mu\text{m}$ Particles
- Fume(Acids) ; $\leq 1.0\mu\text{m}$ (이하) Particles
- Bacteria ; 0.5 μm ~ 0.5m/m Particles
- Virus ; 0.018 ~ 0.4 μm Particles
- VOCs ; Gas (기체)

1-Year Warranty

1m Flexible Hose
 $\phi 65\text{mm}$

Suction Dome
 $\phi 150\text{mm}$

Controller



I. CAFU-19N, 기본형

1-Year Warranty

Suction Dome
 $\phi 250\text{mm}$

1m Flexible Hose
 $\phi 110\text{mm}$

Controller

II. CAFU-A14, 대용량

- Supplied by DAIHAN Korea ■ CE · KC · PSE · FCC · ISO-14001(환경) / 45001(안전) CERTIFIED.
- Suction of Polluted air / Toxic Gas Through Flexible Hose and Purifying by Filtering & Adsorption.
유해가스를 플렉시블 호스로 흡입 복합 필터로 정화(여과 / 흡착)하여 청정공기로 배출.
- KCL(한국건설 생활 환경 시험 연구원) 톨루엔, 자일렌, 아세트산, 폼알데하이드, 암모니아 등 95% 제거 테스트.
- 연구 / 실험실, VOCs 및 Acids 취급, 납땜 및 각종 도색 작업장 등에서 발생하는 유해가스를 효과적으로 정화하여 배출.
- 4-Layer Safety Filtering System / 4겹 필터링 구조 ; Purification and Sterilization of Contamination Source
 - ① - 1 μm Micro Fiber Bag Filters ; 1 μm 이상의 모든 유해를 제거 ; CAFU-19N용
- 1 μm Fiber Bag(5 Room) Filters ; 1 μm 이상의 모든 유해를 제거 ; CAFU-A14용
- 3-Layer Main Filter Kit ;
 - ② Ceramic Honeycomb Carbon Filter ; For the Vapors of VOCs / C₂H₄ / CH₂O / Acids / NH₃, O₃, and Bacteria.
휘발성 유기화합물 VOCs(벤젠·아세치렌, etc) 독성 가스, 산 / 암모니아 등의 흡, 악취, 세균 등을 흡착 / 여과 제거.
 - ③ Impregnated Activated Carbon Filter ; For the Fume / Vapor of Acids / Ammonia / H₂S / VOCs / C₂H₄ / CH₂O, and Bacteria.
산 / 암모니아 / 황화수소 등의 흡, VOCs 휘발성 유기화합물의 독성가스, 악취, 세균 등을 흡착 / 여과 제거.
 - ④ Sterilizing HEPA-H13 Filter ; $\geq 0.3\mu\text{m}$ Particles(세균, 흡, 미세먼지, etc)를 99.95% 여과 제거, 극세 유리섬유.
- - Filter 교체 Indicator 및 교체 주기 설정 ; CAFU-19N
- Filter 교체 차압계 Alarm Sensor 작동 ; CAFU-A14
- * Configurable Filter Replacement Cycle * Easy and Simple Filter Replacement
- * Power Consumption Reduced with a Brushless Direct Current(BLDC) Motor
- * Excellent Chemical · Heat · Weather · Wear Resistance



CAFU-A14 필터교체 알람

Specifications

Spec.	Model	I. "CAFU-19N", 기본형	II. "CAFU-A14", 대용량
Airflow Rate (Suction Capacity)		100 CMH(m ³ /h)	430 CMH(m ³ /h)
Air Speed		14m/s	15m/s
Air Suction Dome (PP)		$\phi 150\text{mm}$	$\phi 250\text{mm}$
Flexible Hose (PP)		$\phi 65 \times 1,000\text{mm}$	$\phi 110 \times 1,000\text{mm}$
Filters and the Replacement Cycle		- 1 μm Micro Fiber Bag Filter ; 분진이 1/3 찼을 때 - 3-Layer Filter Kit ; 3~6 Month	- 1 μm Fiber Bag(5 Room) Filter ; when Sensor Sounds - 3-Layer Filter Kit ; 3~6 Month
Noise Level		$\leq 45\text{db}$	$\leq 50\text{db}$
Controller		Touch PCB, Anti-Static	
Material		Steel Plate 1.2t with Epoxy Powder Coated	
Dimension (w×d×h) & Weight		340×340×h560mm, 18kg	470×425×h790mm, 30kg
Power Consumption		120W	290W
Electrical Reqts.		220V, 50/60Hz	

I. 기본형 "CAFU-19N" 흡 정밀 복합정화기, 100 CMH(m³/h)

※ 운임 별도

4-Layer Filter System

- ① 1.0 μm Micro Fiber Bag Filter ; Particles & Fume 여과 정화
- ② & ③ Ceramic Carbon and Impregnated Carbon Filter ; Fume(Acids) & VOCs Gas 흡착 & 여과 정화 <B3>
- ④ 0.3 μm HEPA-H13 Filter ; Bacteria/Particles/Fume 살균 여과 정화

Ea. Price, ₩

CC.5CAFU.19N Fume Multi-Purifier, Portable

II. 대용량 "CAFU-A14" 흡 복합정화기, 430 CMH(m³/h)

4-Layer Filter System

- ① 1.0 μm Fiber Bag(5 Room) Filter ; Particles 여과 정화
- ② & ③ Ceramic Carbon and Impregnated Carbon Filter ; Fume(Acids) & VOCs Gas 흡착 & 여과 정화
- ④ 0.3 μm HEPA-H13 Filter ; Bacteria/Particles/Fume 살균 여과 정화

※ 운임 별도

<B1>

Ea. Price, ₩

CC.7CAFU.A14 Fume Multi-Purifier, Portable

Spare Filters

- CC.5CAFU.FM2** ① 1 μm Micro Fiber Bag Filter, for "CAFU-19N"
- CC.7CAFU.FM** ① 1.0 μm Fiber Bag(5 Room) Filter, for "CAFU-A14"
- CC.5CAFU.FK** ②+③+④ 3-Layer Filter Kit ; Ceramic Carbon + Active Carbon + HEPA, for "CAFU-19N"
- CC.7CAFU.FK** , for "CAFU-A14"

▶ **Caution ! Please Replace the Polluted Filters in Good Time ! 오염된 필터들을 적시에 교체하세요 !**

◆ How to Replace Filter, 필터 교체방법



1. Before Replacing the Filter,
Power Off with the POWER
Button, and Fix the Wheel.



2. Open the Cover ① and
Replace the Bag Filter.




3. After Opening Cover ② and
Replace Main Filter,
Close Cover ①, ②

CAFU™ Portable Filtering Toxic Gas Multi-Purifier of Organics(VOCs) · Acids · Smell & Dust, 170 & 500 CMH(m³/h)

4-Layer Filters of 0.3µm HEPA · Ceramic Carbon · Active Carbon & Antibacteria, 이동형 유해가스 복합정화기-휘발성 유기화합물 · 산 · 악취 · 미세먼지-필터식

4-Layer Perfect Filter



- ① 10µm Anti-bacterial Prefilter
- ② Ceramic Honeycomb Carbon Filter
- ③ Impregnated Active Carbon Filter
- ④ 0.3µm Sterilization HEPA-H13 Filter

Re. Approx Size of Particulate Matter

- PM2.5 초미세먼지 ; 1.0 ~ 2.5µm Particles
- PM1.0 극초미세먼지 ; ≤1.0µm Particles
- Fume(Acids) ; ≤1.0µm (이하) Particles
- Bacteria ; 0.5µm ~ 0.5m/m Particles
- Virus ; 0.018 ~ 0.4µm Particles
- VOCs ; Gas (기체)

1-Year Warranty



I. 170 CMH, 기본형

- Supplied by DAIHAN Korea ■ CE · KC · PSE · FCC · ISO-14001(환경) / 45001(안전) CERTIFIED.
- Suction of Polluted air / Toxic Gas Through Flexible Hose and Purifying by Filtering & Adsorption.
유해가스를 플렉시블 호스로 흡입 복합 필터로 정화(여과 / 흡착)하여 청정공기로 배출.
- KCL(한국건설 생활 환경 시험 연구원) 톨루엔, 자일렌, 아세트산, 폼알데하이드, 암모니아 등 95% 제거 테스트.
- 연구 / 실험실, VOCs 및 Acids 취급, 납땀 및 각종 도색 작업장 등에서 발생하는 유해가스를 효과적으로 정화하여 배출.
- 4-Layer Safety Filtering System / 4겹 필터링 구조 ; Purification and Sterilization of Contamination Source
 - ① Anti-bacterial Prefilter ; Protects the Main Filter Kit from ≥10µm Dust & Bacteria 메인필터 보호 수명연장
 - ② Ceramic Honeycomb Carbon Filter ; For the Vapors of VOCs / C₂H₂ / CH₂O / Acids / NH₃, O₂, and Bacteria.
휘발성 유기화합물 VOCs(벤젠·아세치렌, etc) 독성 가스, 산 / 암모니아 등의 흡, 악취, 세균 등을 흡착 / 여과 제거.
 - ③ Impregnated Activated Carbon Filter ; For the Fume / Vapor of Acids / Ammonia / H₂S / VOCs / C₂H₂ / CH₂O, and Bacteria.
산 / 암모니아 / 황화수소 등의 흡, VOCs 휘발성 유기화합물의 독성가스, 악취, 세균 등을 흡착 / 여과 제거.
 - ④ Sterilizing HEPA-H13 Filter ; ≥0.3µm Particles(세균, 흡, 미세먼지, etc)를 99.95% 여과 제거, 극세 유리섬유.
- 필터수명/교체시기 Alarm Lamp. ■ Filter Replacement ; See the Previous Page!
- * Configurable Filter Replacement Cycle * Easy and Simple Filter Replacement
- * Power Consumption Reduced with a Brushless Direct Current(BLDC) Motor
- * Excellent Chemical · Heat · Weather · Wear Resistance



필터수명 교체 알람(점등)

I. 170 CMH(m³/h) Models CAFU-A02 ; 기본형

<B3>

Spec.	Model	"CAFU-A02"(기본형)
Airflow Rate (Suction Capacity)		170 CMH(m³/h)
Air Speed		13m/s
Air Suction Dome(PP)		φ250mm
Flexible Hose(POM)		φ77 × 1,000mm
Filters and the Replacement Cycle		- Anti-bacterial Prefilter ; 200 hours - 3-Layer Filter Kit ; 3~6 Month
Noise Level		≤50db
Controller		Push Button Controller
Material		Steel Plate 1.2t with Epoxy Powder Coated
Dimension (w×d×h) & Weight		390 × 390 × 705mm, 26kg
Power Consumption		120W
Electrical Reqts.		220V, 50/60Hz
Cat. No.		CC.2CAFU.A02
Price, Each. ₩		

※ 운임 별도

Spare Filters

- CC.2CAFU.FP ① Anti-Bacterial Prefilter, for "CAFU-A02"
CC.2CAFU.FK ②+③+④ 3-Layer Filter Kit ; Ceramic Carbon & Active Carbon & HEPA, for "CAFU-A02"

Ea. Price, ₩

II. 500 CMH(m³/h) Model CAFU-A13 ; 대용량/Heavy Duty

<B2>



II. 500 CMH, 대용량

Spec.	Model	"CAFU-A13"(대용량)
Airflow Rate (Suction Capacity)		500 CMH(m³/h)
Air Speed		17m/s
Air Suction Dome(PP)		φ250mm
Flexible Hose(PP)		φ110 × 1,000mm
Filters and the Replacement Cycle		- Anti-bacterial Prefilter ; 1 Month - 3-Layer Filter Kit ; 3~6 Month
Noise Level		≤60db
Controller		Touch PCB, Anti-Static
Material		Steel Plate 1.2t with Epoxy Powder Coated
Dimension (w×d×h) & Weight		440 × 455 × 770mm, 32kg
Power Consumption		290W
Electrical Reqts.		220V, 50/60Hz
Cat. No.		CC.3CAFU.A13
Price, Each. ₩		

※ 운임 별도

Spare Filters

- CC.3CAFU.FP ① Anti-Bacterial Prefilter for "CAFU-A13"
CC.3CAFU.FK ②+③+④ 3-Layer Filter Kit ; Ceramic Carbon & Active Carbon & HEPA, for "CAFU-A13"

Ea. Price, ₩

▶ Caution ! Please Replace the Polluted Filters in Good Time ! 오염된 필터들을 적시에 교체하세요 !

CAFU™ Potable Filtering Dust Collector of Fine Dust · Fume · Bacteria "CAFU-21", 150 CMH(m³/h)

2-Layer Filters of 1µm Micro Bag & 0.3µm HEPA-H13, 이동형 집진기-작업장의 먼지·흄·세균등의 정화, 필터형식

2-Layer Perfect Filter



- ① 1µm Fiber Bag(5 Room) Filter
② 0.3µm Sterilization HEPA-H13 Filter

Re. Approx Size of Particulate Matter

- PM2.5 초미세먼지 ; 1.0 ~ 2.5µm Particles
- PM1.0 극초미세먼지 ; ≤1.0µm Particles
- Fume(Acids) ; ≤1.0µm (이하) Particles
- Bacteria ; 0.5µm ~ 0.5m/m Particles
- Virus ; 0.018 ~ 0.4µm Particles
- VOCs ; Gas (기체)

1-Year Warranty

Suction Dome ;
φ 150mm1m Flexible Hose ;
φ 77mmController
차압계
필터 교체 알람

CAFU-21, 150 CMH

▶ 가공, 연마, 레이저 작업장 등에서 발생하는 먼지 및 흄, 박테리아를 효과적으로 정화하여 배출.

■ Supplied by DAIHAN Korea ■ CE · KC · PSE · FCC · ISO-14001(환경) / 45001(안전) CERTIFIED.

■ Suction of Polluted air / Through Flexible Hose and Purifying by Filtering.

작업장 먼지를 플렉시블 호스로 흡입 복합 필터로 정화(여과)하여 청정공기로 배출.

■ KCL(한국건설 생활 환경 시험 연구원) 부유미생물 저감 시험 세균감소율 99.9% 테스트 통과.

■ 2Layer Filtering System / 2겹 필터링 구조 ;

① 1µm Fiber Bag(5 Room) Filters ; 1µm 이상의 모든 유해물 입자 95% 제거

② Sterilizing HEPA-H13 Filter ; ≥0.3µm Particles(세균, 흄, 미세먼지, etc)를 99.95% 여과 제거, 극세 유리섬유.

* Filter 교체 차압계 알람 Sensor 작동 * Filter Replacement ; See Page 139!

* Power Consumption Reduced with a Brushless Direct Current(BLDC) Motor

Specifications

<B2>

Spec.	Model	"CAFU-21", 중형
Airflow Rate (Suction Capacity)		150 CMH(m³/h)
Air Speed		14m/s
Air Suction Dome(PP)		φ 150mm
Flexible Hose(POM)		φ 77 × 1,000mm
Filters and the Replacement Cycle		- 1µm Fiber Bag(5 Room) Filter ; ≥1µm 입자 95% 제거 ; when Sensor Sounds(차압계 센서가 울릴때 교체!) - 0.3µm Sterilizing HEPA-H13 Filter; 3~6 Month에 교체!
Noise Level		≤70db
Controller		POWER 버튼 ; Operation and Stop (Toggle Switch) VOLUME 버튼 ; Wind Speed Adjustment (Jog-controller)
Material		Plate 1.2t with Epoxy Powder Coated
Dimension (w × d × h) & Weight		390 × 390 × h600mm, 25kg
Power Consumption		140W
Electrical Reqts.		220V, 50/60Hz
Cat. No.		CC.6CAFU.21
Price, Each. ₩		

Spare Filters

CC.6CAFU.FM 1.0µm Fiber Bag(5 Room) Filter, for "CAFU-21"

CC.6CAFU.FH 0.3µm Sterilization HEPA-H13 Filter, for "CAFU-21"

Ea. Price, ₩

▶ Caution ! Please Replace the Polluted Filters in Good Time ! 오염된 필터들을 적시에 교체하세요 !

, Toxic Gas Purifier Box

Ductless Filtering Vapor(Fume & Gas) Purifier Box for Chemicals Storage Cabinets "SCF90-F", Additional Apparatus

With 4 Filters of HEPA · Acids · VOCs · NH₃/CH₂O and 2 × φ 60mm Connecting Device Hole for In- & Out-let

약품장용 유해가스 종합 정화기 박스, 부가 설치용, 무덕트 필터링형, 사용종인(사용할) 약품장에 부가적으로 설치함

3-Year Warranty

Ductless Filtering Vapor(Fume & Gas)
Purifier Box■ 설치법 ; 사용 종인(사용할) 약품장 상부에 φ 60mm Hole이 있거나 혹은 뚫어서,
본 정화기와 PVC Flexible Hose φ 50으로 연결

■ 안전 캐비닛 내에 유해가스 발생 시 덕트없이 필터 여과 정화, 청정 공기 실내에 공급.

■ Powerful Forced Airflow and Low Noise by High Performance BLDC Motor

■ Complete Set with 2 Filters ; (1) Acid Filter, (2) VOCs Gas(휘발성 유기화합물) Formaldehyde Gas

■ 별도로 구매하여 사용자가 직접 Safety Cabinet에 쉽게 연결가능 (※ 자세한 설치 방법은 사용자 매뉴얼 참조)

■ Filter Replacement Cycle Alarm Function ; 필터교체 주기는 약 6~12개월로 사용 조건에 따라 다를 수 있음

<L1>

DH.SCF90F

Filtering Gas Purifier Box, 480 × 323 × h340mm, Powder-Coated Steel,
Complete Set with 2 Filters

Ea. Price, ₩

Spare Framed Filters For Vapor Purifier Box.

* The Size of Framed Filters ; 300 × 300 × h20mm

DH.SCF90F3 Chemical Coated Activated Carbon Filter, for Fume(Acids)


DH.SCF90F5 Activated Carbon Filter, for Organics(VOCs) & Formaldehyde Gas

Ea. Price, ₩

▶ Filter Selection Guide ; See Page 78!

Spare & Optional
Filter KitsInspection
Certificate of Filters
MDI CERTIFIED

DAIHAN™ 1,100℃ 0.4Lit Dental Compact Porcelain Furnace “FUR040.P”, with Ceramic Fiber Dish, Cylindrical $\phi 90 \times h65$ mm
For Pigmentation on the Outside of the Prosthesis, Short Heat-up Time, High Quality Insulation ; Ceramic Fiber Board, 치과용 전기로, 원통형



- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 세라믹, 지르코니아 등 보철물 외면에 착색, 색상 발현 용도에 적합
- * Suitable for Coloring and Color Expression on the Outer Surface of Prosthesis
- * Door Opening can be Fixed with a Door Fixing Handle
- * Light Weight and Easy to move

※ 운임 별도 <L1>

Spec.	Model	“FUR040.P”
Temp Range		400℃ ~ 1,100℃
Capacity (Chamber)		0.4 Lit, Cylindrical $\phi 90 \times h65$ mm
Heater / Power		Kanthal A1, Inside Spiral Quartz Tube / 1.1kW
Thermocouple		K-Type Thermocouple
Timer		99hr 59min
Display		Digital LCD Display
Controller		Digital PID Programmable Controller
Program Function		5 Patterns with 9 Segments
Material / Insulation		External - Powder Coated Steel / Internal - Ceramic Fiber Insulation
Dimension (w×d×h)	Internal (mm)	$\phi 90 \times h65$
	External (mm)	w270×d290×h340
Gross Weight		16kg
Electrical Reqts.		1Phase, AC120V, 60Hz or AC 230V, 50/60Hz
Cat. No.	230V	DH.FurP040
	120V	DH.FurP041
Price, Each. ₩		

DAIHAN™ 1,650℃ 0.7 & 1.8Lit Dental Sintering Furnace “FUR070.S” & “FUR180.S”, with Ceramic Fiber Dish, 2-Side Heating
For Firmly Fixing the Prosthesis, Exposed Heating Elements, Short Heat-up Time, High Quality Insulation ; Ceramic Fiber Board, 치과용 고온 전기로



- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 세라믹, 지르코니아 등 보철물을 가열하여 단단히 고정하는데 적합
- * Suitable for Heating and Solidifying Prosthesis
- * Door Opening can be Fixed with a Door Fixing Handle
- * Light Weight and Easy to move

※ 운임 별도 <L1>

Spec.	Model	“FUR070.S”	“FUR180.S”
Temp Range		500℃ ~ 1,650℃	
Capacity		0.7 Lit	1.8 Lit
Heater / Power		MoSi ₂ Heater, 2-Side Heating / 600W	MoSi ₂ Heater, 2-Side Heating / 1.4kW
Thermocouple		B-Type Thermocouple	
Timer		99hr 59min	
Display		Digital LCD Display	
Controller		Digital PID Programmable Controller, EPC3004	
Program Function		20 Patterns with 8 Segments	
Material / Insulation		External - Powder Coated Steel / Internal - Ceramic Fiber Insulation	
Dimension (w×d×h)	Internal (mm)	90×95×90	120×130×120
	External (mm)	320×380×490	350×420×560
Gross Weight		31kg	41kg
Electrical Reqts.		1Phase, AC120V, 60Hz or AC 230V, 50/60Hz	
Cat. No.	230V	DH.FurS070	DH.FurS180
	120V	DH.FurS071	DH.FurS181
Price, Each. ₩			
Spare ; DH.FurS5003	(1) Ceramic Fiber Dish, 80×83×h16mm, for “Fur.S”	Ea	
; DH.FurS5004	, 80×83×h32mm, for “Fur.S”	Ea	
; DH.FurS5005	(2) Ceramic Fiber Stand, 90×25×h10mm, for Dish, 2ea/set	Set	

※ NOTE! Optional Heat Resistant Gloves
⇒ See the Next Page!

DAIHAN™ 1,000℃ Programmable Digital Muffle Furnace “FX”, Exposed Heating Elements-type
With Digital PID Control, Short Heat-up Time, Ceramic Fiber, 4-Side Heating, with Ceramic Fiber Plate, 3~63 Lit
디지털 전기로, 디지털 PID 컨트롤 시스템, 4면 가열 방식, 히터 노출형

CE

3-Year Warranty



Jog-Shuttle with Push Button

3~27Lit Model

CE

3-Year Warranty



63Lit Model

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function
- Digital PID Programmable Controller ; 4 Patterns with 10 Segments
- * Designed for Ash-Determinations, Enameling, Fusions, Precipitate Drying
- * User's Self-Compensation Function for Temperature ($\pm 50^\circ\text{C}$) * Resolution ; 1°C , 1min
- * Safety Switch for Door Opening * Ceramic Fiber exposed Heating Wire, 4-Side Heating
- * Chamber Aging is required before using the Furnace to Remove Impurities of the Insulation material and to Prevent Damage ; 500°C 6hr, 900°C 1hr

Specifications & Ordering information

<L1>

Spec.	Model	“FX-03” Max.Capa. 3 Lit	“FX-05” Max.Capa. 4.5 Lit	“FX-12” Max.Capa. 12 Lit	“FX-14” Max.Capa. 14 Lit	“FX-27” Max.Capa. 27 Lit	“FX-63” Max.Capa. 63 Lit
Temp. Range		300℃ ~ 1000℃					
Capacity		3 Lit	4.5 Lit	12 Lit	14 Lit	27 Lit	63 Lit
Heater		Kanthal A1, Expose Wire Heating Elements, 4-Side Heating					
Thermocouple		K-type Thermocouple					
Heating Power		1.4 kW	1.7 kW	3.3 kW	4.0 kW	4.6 kW	9 kW
Required time to be reached 800℃		16 min		20 min			23 min
Timer & Alarm		99hr 59min (Continuous run), Error status & Timer-end					
Display		Digital LCD with Back-Light Function					
Controller		Digital PID Controller with Jog-Dial Switch (Turn+Push)					
Program Function		4 Patterns with 10 Segments					
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection Safety Switch for Door Opening					
Material		Ceramic Fiber with Heating Wire-Internal, Powder Coated steel-External Ceramic Fiber-Insulation					
Dimension (w×d×h)	Internal (cm)	13×25×9	15×30×10	20×30×20	25×35×16	30×30×30	40×40×40
	External (cm)	41×48×53	44×53×55	50×55×66	55×59×62	60×59×76	100×75×95
Packing Size(w×d×h) (cm)		63×54×77	68×57×78	70×63×89	74×68×85	74×73×99	116×90×110
Gross Weight		47 kg	54 kg	59 kg	73 kg	83 kg	225 kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz					
Cat. No.	230V	DH.WF21.03	DH.WF21.05	DH.WF21.12	DH.WF21.14	DH.WF21.27	DH.WF21.63
	120V	DH.WF20.03	DH.WF20.05	DH.WF20.12	DH.WF20.14	DH.WF20.27	-
	3Phase	-	-	DH.WF23.12	DH.WF23.14	DH.WF23.27	DH.WF23.63
Price, Each. ₩							

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Ceramic Fiber Plate

- DH.WF211.03 (1) Ceramic Fiber Plate, 120×215×h10mm, “FS1103” for 3Lit
 DH.WF211.05 , 140×265×h10mm, “FS1105” for 4.5Lit
 DH.WF211.12 , 190×265×h10mm, “FS1112” for 12Lit
 DH.WF211.14 , 240×310×h10mm, “FS1114” for 14Lit
 DH.WF211.27 , 265×290×h10mm, “FS1127” for 27Lit
 DH.WF211.63 , 400×400×h10mm, “FS1163” for 63Lit

Ea. Price, ₩

Optional Glove, Heat Resistant

- SL.Glo6635 (2) Glove, Heat Resistant, 1000℃, Standard Type, L350mm
 SL.Glo6636 , Long Type, L450mm
 SL.Glo6421 (3) , 400℃, L340mm
 SL.Glo6422 , L450mm

Ea. Price, ₩



Inner Chamber of Exposed Heating Elements-type



(1) Ceramic Fiber Plate (Included)



(2) Heat-resistant Glove (1,000℃)



(3) Heat-resistant Glove (400℃)

DAIHAN™ 1,200°C Programmable Digital Muffle Furnace “FHX”, Exposed Heating Elements-type
With Digital PID Control, Short Heat-up Time, Ceramic Fiber, 4-Side Heating, with Ceramic Fiber Plate, 3~63 Lit
디지털 전기로, 디지털 PID 컨트롤 시스템, 4면 가열 방식, 히터 노출형



3-Year Warranty



Jog-Shuttle with Push Button

3~27Lit Model



3-Year Warranty



63Lit Model



Inner Chamber of Exposed Heating Elements-type

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Patented Jog-Dial Control System and Digital LCD with Back-Light Function
- Digital PID Programmable Controller ; 4 Patterns with 10 Segments

- * Digital Jog-Dial Control System
- * Included Ceramic Fiber Plate
- * Ceramic Fiber Exposed Heating Wire, 4-Side Heating
- * Safety Switch for Door Opening
- * User's Self-Compensation Function for Temperature ($\pm 50^\circ\text{C}$)
- * Resolution ; 1°C , 1min

Specifications & Ordering information

<L1>

Spec.	Model	"FHX-03" Max.Capa. 3 Lit	"FHX-05" Max.Capa. 4.5 Lit	"FHX-12" Max.Capa. 12 Lit	"FHX-14" Max.Capa. 14 Lit	"FHX-27" Max.Capa. 27 Lit	"FHX-63" Max.Capa. 63 Lit
Temp. Range		300°C ~ 1200°C					
Capacity		3 Lit	4.5 Lit	12 Lit	14 Lit	27 Lit	63 Lit
Heater		Kanthal A1, Expose Wire Heating Elements, 4-Side Heating					
Thermocouple		K-type Thermocouple					
Heating Power		1.4 kW	1.7 kW	3.3 kW	4.0 kW	4.6 kW	9 kW
Required time to be reached 800°C		16 min		20 min			23 min
Timer & Alarm		99hr 59min (Continuous run), Error status & Timer-end					
Display		Digital LCD with Back-Light Function					
Controller		Digital PID Controller with Jog-Dial Switch (Turn+Push)					
Program Function		4 Patterns with 10 Segments					
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection Safety Switch for Door Opening					
Material		Ceramic Fiber with Heating Wire-Internal, Powder Coated steel-External Ceramic Fiber-Insulation					
Dimension (w×d×h)	Internal (cm)	13×25×9	15×30×10	20×30×20	25×35×16	30×30×30	40×40×40
	External (cm)	41×48×53	44×53×55	50×55×66	55×59×62	60×59×76	100×75×95
Packing Size(w×d×h) (cm)		63×54×77	68×57×78	70×63×89	74×68×85	74×73×99	116×90×110
Gross Weight		47 kg	54 kg	59 kg	73 kg	83 kg	225 kg
Electrical Reqts.		1 Phase AC 120V, 60Hz or AC 230V, 50/60Hz or 3 Phase					
Cat. No.	230V	DH.WFH22.03	DH.WFH22.05	DH.WFH22.12	DH.WFH22.14	DH.WFH22.27	DH.WFH22.63
	120V	DH.WFH20.03	DH.WFH20.05	DH.WFH20.12	DH.WFH20.14	DH.WFH20.27	-
	3Phase	-	-	DH.WFH23.12	DH.WFH23.14	DH.WFH23.27	DH.WFH23.63
Price, Each. ₩							

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Ceramic Fiber Plate

- DH.WF211.03 (1) Ceramic Fiber Plate, 120×215mm, "FS1103" for 3Lit
- DH.WF211.05 , 140×265mm, "FS1105" for 4.5Lit
- DH.WF211.12 , 190×265mm, "FS1112" for 12Lit
- DH.WF211.14 , 240×310mm, "FS1114" for 14Lit
- DH.WF211.27 , 265×290mm, "FS1127" for 27Lit
- DH.WF211.63 , 400×400mm, "FS1163" for 63Lit

Ea. Price, ₩

Optional Glove, Heat Resistant

- SL.Glo6635 (2) Glove, Heat Resistant, 1000°C, Standard Type, L350mm
- SL.Glo6636 , Long Type, L450mm
- SL.Glo6421 (3) , 400°C, L340mm
- SL.Glo6422 , L450mm

Ea. Price, ₩



(1) Ceramic Fiber Plate (Included)



(2) Heat-resistant Glove (1,000°C)



(3) Heat-resistant Glove (400°C)

DAIHAN™ 1,450℃ Programmable High-Temp Muffle Furnace “FSC”, Exposed Heating Elements-type
With SIC Heater, Digital PID Control, Short Heat-up Time, 2-Side Heating, without Ceramic Fiber Plate, 1.9~22Lit
고온 디지털 전기로, 디지털 PID 컨트롤 시스템, 2면 가열 방식

3-Year Warranty



1.9Lit Model

- Supplied by DAIHAN Korea
- PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital PID Programmable Controller ; 2 Patterns with 10 Segments
- Digital LCD Display
- High Quality Insulation ; Ceramic Fiber Board
- SIC(Silicon Carbide) Heater

- * Safety Door Lock System
- * Excluding Ceramic Fiber Plate
- * Ceramic Fiber Exposed Heating Wire, 2-Side Heating
- * Safety Switch for Door Opening

Specifications & Ordering information

<L1>

Spec.	Model	FSC-2 Max. Capa. 1.9Lit	FSC-5 Max. Capa. 4.5Lit	FSC-11 Max. Capa. 11Lit	FSC-22 Max. Capa. 22Lit
Temp. Range	500℃~1,450℃				
Capacity		1.9 Lit	4.5 Lit	11 Lit	22 Lit
Heater	SIC Heater, 2-Side Heating				
Thermocouple	R-type Thermocouple				
Heating Power		2.8 kW	3.3 kW	7.7 kW	9.6 kW
Required time to be reached 800℃	30 min				
Timer & Alarm	99hr 59min or Continuous run				
Display	Digital LCD Display				
Controller	Digital PID Controller, UP35A, SCR Power Controller with Push Button				
Program Function	2 Patterns with 10 Segments				
Safety Circuit	Over-Temp./-Current protection, Sensor Error Detection, Safety Switch for Door Opening				
Material	Ceramic Fiber board with SIC Heater-Internal, Powder Coated steel-External, Ceramic Fiber-Insulation				
Dimension (w×d×h)	Internal (mm)	120×130×120	150×200×150	200×270×200	250×350×250
	External (mm)	440×482×700	820×555×655	870×625×715	920×705×785
Packing Size(w×d×h) (mm)		680×540×940	1350×600×890	1400×770×950	1450×900×1100
Gross Weight		57 kg	70 kg	102 kg	132 kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz, without Cord and Plug			
Cat. No.	230V	DH.FSC02002	DH.FSC02005	DH.FSC02011	DH.FSC02022
	120V	DH.FSC01002	-	-	-
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Optional Ceramic Fiber Plate

- DH.FSC2002 (1) Ceramic Fiber Plate, w110×d120mm, for 1.9Lit
 DH.FSC2005 , w140×d180mm, for 4.5Lit
 DH.FSC2011 , w190×d250mm, for 11Lit
 DH.FSC2022 , w240×d330mm, for 22Lit

Ea. Price, ₩

Optional Glove, Heat Resistant

- SL.Glo6635 (2) Glove, Heat Resistant, 1000℃, Standard Type, L350mm
 SL.Glo6636 , Long Type, L450mm
 SL.Glo6421 (3) , 400℃, L340mm
 SL.Glo6422 , L450mm

Ea. Price, ₩

Inner Chamber of
Exposed Heating Elements-type(1) Ceramic Fiber Plate
(Included)(2) Heat-resistant Glove
(1,000℃)(3) Heat-resistant Glove
(400℃)

DAIHAN™ 1,650℃ Programmable High-Temp Muffle Furnace "FSK", Exposed Heating Elements-type
With MoSi₂ Heater, Digital PID Control, Short Heat-up Time, 2-Side Heating, without Ceramic Fiber Plate, 1.9~11 Lit
고온 디지털 전기로, 디지털 PID 컨트롤 시스템, 2면 가열 방식



1.9Lit Model

3-Year Warranty

- Supplied by DAIHAN Korea
- PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital PID Programmable Controller ; 2 Patterns with 10 Segments
- Digital LCD Display
- High Quality Insulation ; Ceramic Fiber Board

- * Safety Door Lock System
- * Excluding Ceramic Fiber Plate
- * Ceramic Fiber Exposed Heating Wire, 2-Side Heating
- * Safety Switch for Door Opening

Specifications & Ordering information

<L1>

Spec.	Model	FSK-2 Max. Capa. 1.9 Lit	FSK-5 Max. Capa. 4.5 Lit	FSK-11 Max. Capa. 11 Lit
Temp. Range		600℃ ~ 1,650℃		
Capacity		1.9 Lit	4.5 Lit	11 Lit
Heater		MoSi ₂ Heater, 2-Side Heating		
Thermocouple		B-type Thermocouple		
Heating Power		1.3 kW	3.0 kW	5.2 kW
Required time to be reached 800℃		40 min		
Timer & Alarm		99hr 59min or Continuous run		
Display		Digital LCD Display		
Controller		Digital PID Controller, UP35A, SCR Power Controller with Push Button		
Program Function		2 Patterns with 10 Segments		
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection, Safety Switch for Door Opening		
Material		Ceramic Fiber board with MoSi ₂ Heater-Internal, Powder Coated steel-External, Ceramic Fiber-Insulation		
Dimension (w×d×h)	Internal (mm)	120×130×120	150×200×150	200×270×200
	External (mm)	440×482×700	900×600×695	980×670×765
Packing Size(w×d×h) (mm)		680×540×940	1430×640×930	1510×810×1000
Gross Weight		72 kg	150 kg	186 kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz, without Cord and Plug		
Cat. No.	230V	DH.FSK02002	DH.FSK02005	DH.FSK02011
	120V	DH.FSK01002	-	-
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도



4.5 & 11Lit Model

3-Year Warranty



Inner Chamber of
Exposed Heating Elements-type

Optional Ceramic Fiber Plate

- DH.FSK2002 (1) Ceramic Fiber Plate, w110×d120mm, for 1.9Lit
- DH.FSK2005 , w140×d180mm, for 4.5Lit
- DH.FSK2011 , w190×d250mm, for 11Lit

Ea. Price, ₩

Optional Glove, Heat Resistant

- SL.Glo6635 (2) Glove, Heat Resistant, 1000℃, Standard Type, L350mm
- SL.Glo6636 , Long Type, L450mm
- SL.Glo6421 (3) , 400℃, L340mm
- SL.Glo6422 , L450mm

Ea. Price, ₩



(1) Ceramic Fiber Plate
(Included)



(2) Heat-resistant Glove
(1,000℃)



(3) Heat-resistant Glove
(400℃)

DAIHAN™ Programmable Tube Furnaces “FT-series”, Embedded Heating Elements-type, 1,000℃
 With Digital PID Programmable Control, 1- or 3-Heating Zone, $\phi 6 \sim \phi 12$ cm Tubes, 1,000℃ 프로그램식 튜브 전기로

3-Year Warranty

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

- * Digital PID Programmable Controller ; 3 Patterns with 16 Segments
- * Designed for Ash-Determinations, Horizontal Application, Protective Gas Atmosphere
- * Available Order Various Size Tubes, Without Working Tube
- * Programmable Timer ; 99hr 59min * Resolution ; 1℃, 1min



(1) 1 Heating Zone Model



Figure ; Open Door



(1) Side Figure of 1 Heating Zone Model



(a) Heat-resistant Glove (1,000℃)



(b) Heat-resistant Glove (400℃)

Specifications & Ordering information

<L1>

Model	“FT-series” 1- or 3-Heating Zone						
Temp. Range & Heater	300℃~1,000℃, Kantal A1, Embedded Wire Heating Elements						
Thermocouple & Display	K-type Thermocouple, LED						
Safety Circuit	Over-Temp./-Current protection, Sensor Error Detection						
Material	Ceramic Fiber with Heating Wire-Internal, Powder Coated steel-External						

(1) 1Heating Zone Model	“FT-830”	“FT-860”	“FT-1230”	“FT-1260”	“FT-1430”	“FT-1460”	“FT-1630”	“FT-1660”
Heating Zone	1×Heating Zone							
Available Tube Size(cm)	Max. $\phi 6$ cm		Max. $\phi 8$ cm		Max. $\phi 10$ cm		Max. $\phi 12$ cm	
Heating Power	2.2kW	3.7kW	2.6kW	5.3kW	3.7kW	6.0kW	4.0kW	6.9kW
Heat-up Time (to 800℃)	17min	20min	18min	22min	19min	23min	20min	25min
Dimension (w×d×h cm)	50×37×54	80×37×54	50×40×57	80×40×57	50×43×60	80×43×60	50×46×63	80×46×63
Controller	1×Programmable PID Controller ; 5 Patterns 16 Segments							

Cat. No.	230V	DH.WFT12.0830	DH.WFT12.0860	DH.WFT12.1230	DH.WFT12.1260	DH.WFT12.1430	DH.WFT12.1460	DH.WFT12.1630	DH.WFT12.1660
	120V	DH.WFT11.0830	DH.WFT11.0860	DH.WFT11.1230	DH.WFT11.1260	DH.WFT11.1430	DH.WFT11.1460	DH.WFT11.1630	DH.WFT11.1660

Price, Each. ₩									
----------------	--	--	--	--	--	--	--	--	--

(2) 3Heating Zone Model	“FT-3860”	“FT-31260”	“FT-31460”	“FT-31660”
Heating Zone	3×Heating Zone			
Available Tube Size(cm)	Max. $\phi 6$ cm	Max. $\phi 8$ cm	Max. $\phi 10$ cm	Max. $\phi 12$ cm
Heating Power	3.7kW	5.3kW	6.0kW	6.9kW
Heat-up Time (to 800℃)	20min	22min	23min	25min
Dimension (w×d×h cm)	70×28×44	70×34×50	70×36×54	70×40×56
Controller	Independent 3×Programmable PID Controller ; 3 Patterns 16 Segments			

Cat. No.	230V	DH.WFT12.860.3	DH.WFT12.1260.3	DH.WFT12.1460.3	DH.WFT12.1660.3
	120V	DH.WFT11.860.3	DH.WFT11.1260.3	DH.WFT11.1460.3	DH.WFT11.1660.3
Price, Each. ₩					

- * Other Specifications are available upon Customer's Request.
- * Tubes Size & Other Specifications are available upon Customer's Request

※ 운임 별도

Optional Glove, Heat Resistant

- SL.Glo6635 (a) Glove, Heat Resistant, 1000℃, Standard Type, L350mm
 SL.Glo6636 , Long Type, L450mm
 SL.Glo6421 (b) , 400℃, L340mm
 SL.Glo6422 , L450mm


Ea. Price, ₩

DAIHAN™ Laboratory Assembly Furniture, High Quality Steel-Frame & Side Panel · Phenol Work Top · Stainless-steel Bolted Joint
 With Transfer Cabinet, Utility Box, 실험실용 조립식 실험대 개요, 고품질 스틸 프레임, 내열성/내충격성/내화학성 페놀 상판, 볼트식 결합

- Supplied by DAIHAN Korea
- Easy to Assemble with Stainless-steel(STS 304) Bolts and Stable
- Phenol Work Top ; Good for Heat-/Shock-Resistant, Excellent Chemical-Resistant(Acetic Acid, Alcohol, Hydrochloric Acid, Sulfuric Acid, Organic Chemicals, Base, Biological Stains, Salt, Solvent, etc.)
- High Quality Steel-Frame(60×40 & 50×30mm), Side-/Back-Steel Panel(t1.6mm)
- The Work Top are made of Melamine, Granite, Phenol(Color Choices Available ; Black or Gray)
- You can Choose the Materials of the Work Top according to Usage and Change the Height of Table.

Features of Work Top Materials


- * Phenol(Phenolic Laminate) ; Excellent for Heat-/ Chemical-/ Shock-Resistance, Strong Plate Consist of Melamine-resin and Phenolic-resin, Surface Strengthening through Chemical and Chemical Treatment, Good for All Purpose Laboratory.
- * Melamine ; Good for Chemical and Heat-Resistance. The Surface is Resistant against Scratches and Abrasion. Suitable for Analytical Laboratory.
- * Granite ; Excellent for Heat and Shock-Resistance. Suitable for Special Laboratory like Instrumental Analysis Laboratory.




Structures and Features of Assembly Furniture

- * (a) Phenol Work Top ; Good for Heat-/Shock-Resistant Phenolic Laminate(t20mm), Excellent Chemical-Resistant (Acetic Acid, Alcohol, Hydrochloric Acid, Sulfuric Acid, Organic Chemicals, Base, Biological Stains, Salt, Solvent, etc.), Color Choices Available ; Black or Gray
- * (b) Frame ; The Most Stable Steel Under Construction, Suitable for Long-term Loading of Heavy Analysis Equipment, High Quality Steel, 60×40mm & 50×30mm
- * (c) Shelves ; Shelf Fall Prevention Design, Movable Up and Down, Aluminum Frame, PVC Shelf Plate(t8mm), Epoxy Powder Coating Steel Shelf Supporter(t1.6mm)
- * (d) Utility Box ; Electrical safety design, Attachment Emergency Breaker and Switch Box for Each Line, Steel Plate(t1.6mm), Cover type Power Outlet(220V)
- * (e) Reagent Structure ; Al Frame & Steel Plate(1.6mm), Columns bend(Edge Rounding) to Secure LAN Line, Electricity, Gas, Water etc., Maintenance and Safety Enhancement
- * (f) Reagent Front / Rear Frame ; Anodizing and Electro-deposition Treatment on Environmentally Friendly Aluminum Surface to Enhance Strength, Abrasion-Resistance, Corrosion-Resistance and Electrical insulation
- * (g) Modular System Back Panel
 - Easily Removable and Easy to Follow and Check (Scalability)
 - Eco-friendly Epoxy Powder Coating to Prevent Scratches and Corrosion
 - Using Steel Plate t1.6mm
- * (h) Transfer Cabinet ; Particle Board Plate(t15/18mm) / Al-handle / Brake Caster / 450×500×h700mm, Basic 4-stage(2-/3-stage possible)
- * Attach Inner Support ; Use of Steel Pipe for Loaded Area and Flexural Reinforcement
- * Sink Tables ; We can Offer DAIHAN Custom Sink Table, for which you can Specify the Specification

DAIHAN™ Laboratory Sink Table, High Quality Steel-Frame & -Side Panel, Phenol Sink Top, PP Sink Bowl, Glass Water Guard
 With ① Cold & Hot Water Mixed Faucet, 실험실용 싱크대, 3면 강화유리, 재질 변경 가능



- * ② 3-mouth Chemical Faucet Available(Optional)
- * Sink Top ; Chemical-resistance Phenolic Laminate t20mm
- * Sink Bowl ; PP(Selectable Material)
- * Water Guard ; Composed of 3-sided Tempered Glass t8mm(Selectable Material)



Selectable Material of Custom-Made Sink Table <L4>

	Material			
(a) Sink Top	Phenol	P.P	PVC	STS
(b) Sink Bowl	-	P.P	PVC	STS
(c) Water Guard	Glass	P.C	PVC	PVC

- DH.WFE4031 Sink Table, Laboratory, 900×750×h800 / 950mm
- DH.WFE4032 , 1200×750×h800 / 950mm
- DH.WFE4033 , 1500×750×h800 / 950mm
- DH.WFE4034 , 1800×750×h800 / 950mm

Optional Eye-washer

- DH.WFE4040 (d) Option ; Eye-washer, Single Nozzle, Table Mounted, L1.4m

※ Note! Freight and Installation Charges Separately.

Ea.	Price, ₩
Ea.	Price, ₩

Welcome! Custom-Made Sink Table!!

We can Offer DAIHAN Custom Sink Table to Your Specifications.

DAIHAN™ Laboratory Assembly Center Table, High Quality Steel-Frame & -Side Panel · Phenol Work Top · Stainless-steel Bolted Joint
 With Transfer Cabinet, Utility Box, 실험실용 조립식 중앙 실험대, 고품질 스틸 프레임, 내열성/내충격성/내화학성 페놀 상판, 볼트식 결합

- Supplied by DAIHAN Korea ■ Easy to Assemble with Stainless-steel(STS 304) Bolts and Stable
- Phenol Work Top ; Good for Heat-Shock-Resistant, Excellent Chemical-Resistant(Acetic Acid, Alcohol, Hydrochloric Acid, Sulfuric Acid, Organic Chemicals, Base, Biological Stains, Salt, Solvent, etc.)
- High Quality Steel-Frame(60×40 & 50×30mm), Side-/Back-Steel Panel(t1.6mm)
- * The Most Stable Steel Under Construction, Suitable for Long-term Loading of Heavy Analysis Equipment (60×40 & 50×30mm)
- * Chemical-resistance Phenolic Laminate Work Top, t20mm (Color Choices Available ; Black or Gray)
- * Could be Change the Work Top, Height of Table on User Applications. Please Call to specify the Work Top.

(1) Assembly Basic Center Table, with Transfer Cabinet, Steel-Frame & -Side Panel, Phenol Work Top

- * Work Top ; Phenol Plate, t20mm
- * Back Panel ; Steel Plate, t1.6mm
- * Frame ; Steel, 60×40 & 50×30mm



WFE301 Series



- DH.WFE3011 Center Table, 1500×1500×h800mm, with Cabinet(2ea)
 DH.WFE3012 , 1800×1500×h800mm, with Cabinet(2ea)
 DH.WFE3013 , 2100×1500×h800mm, with Cabinet(4ea)
 DH.WFE3014 , 2400×1500×h800mm, with Cabinet(4ea)
 DH.WFE3015 , 3000×1500×h800mm, with Cabinet(4ea)
 DH.WFE3016 , 3600×1500×h800mm, with Cabinet(4ea)

※ Note! Freight and Installation Charges Separately.

<L4>


Ea. Price, ₩

(2) Assembly Shelf-type Center Table, with Transfer Cabinet & Shelves, Steel-Frame & -Side Panel, Phenol Work Top

- * Work Top ; Phenol Plate, t20mm
- * Back Panel ; Steel Plate, t1.6mm
- * Frame ; Steel, 60×40 & 50×30mm
- * Size ; d1500×h800/1900mm



WFE302 Series



- DH.WFE3021 Center Table, w1500×h800/1900mm, with Cabinet(2ea) / Shelves
 DH.WFE3022 , w1800×h800/1900mm, with Cabinet(2ea) / Shelves
 DH.WFE3023 , w2100×h800/1900mm, with Cabinet(4ea) / Shelves
 DH.WFE3024 , w2400×h800/1900mm, with Cabinet(4ea) / Shelves
 DH.WFE3025 , w3000×h800/1900mm, with Cabinet(4ea) / Shelves
 DH.WFE3026 , w3600×h800/1900mm, with Cabinet(4ea) / Shelves

※ Note! Freight and Installation Charges Separately.


Ea. Price, ₩

(3) Assembly Comprehensive Center Table, with Cabinet · Sink Table · Shelves, Steel-Frame & -Side Panel, Phenol Work Top

- * Work Top ; Phenol Plate, t20mm
- * Sink Top ; Phenol Plate, t20mm
- * Back Panel ; Steel Plate, t1.6mm
- * Water Guard ; 3-sided Tempered Glass
- * Frame ; Steel, 60×40 & 50×30mm
- * with Cabinet(4ea) / Shelves / Sink Table
- * Size ; d1500×h800/1900mm
- * Sink Bowl ; PP(Polypropylene)



WFE303 Series



- DH.WFE3031 Center Table, w2250×h1900mm, with Shelves / Sink Table
 DH.WFE3032 , w2550×h1900mm, with Shelves / Sink Table
 DH.WFE3033 , w2850×h1900mm, with Shelves / Sink Table
 DH.WFE3034 , w3150×h1900mm, with Shelves / Sink Table
 DH.WFE3035 , w3750×h1900mm, with Shelves / Sink Table
 DH.WFE3036 , w4350×h1900mm, with Shelves / Sink Table

※ Note! Freight and Installation Charges Separately.

Ea. Price, ₩

DAIHAN™ Laboratory Assembly Side Table, High Quality Steel-Frame & -Side Panel · Phenol Work Top · Stainless-steel Bolted Joint
With Transfer Cabinet, Utility Box, 실험실용 조립식 벽면 실험대, 고품질 스틸 프레임, 내열성/내충격성/내화학성 페놀 상판, 볼트식 결합


- Supplied by DAIHAN Korea
- Easy to Assemble with Stainless-steel(STS 304) Bolts and Stable
- Phenol Work Top ; Good for Heat-/Shock-Resistant, Excellent Chemical-Resistant(Acetic Acid, Alcohol, Hydrochloric Acid, Sulfuric Acid, Organic Chemicals, Base, Biological Stains, Salt, Solvent, etc.)
- High Quality Steel-Frame, Side-/Back-Steel Panel(t1.6mm)
- * The Most Stable Steel Under Construction, Suitable for Long-term Loading of Heavy Analysis Equipment (60×40mm & 50×30mm)
- * Chemical-Resistance Phenolic Laminate Work Top, t20mm (Black or Gray)
- * Could be Change the Work Top, Height of Table on User Applications.

(1) Assembly Basic Side Table, with Transfer Cabinet, Steel-Frame & -Side Panel, Phenol Work Top

- * Work Top ; Phenol Plate, t20mm
- * Back Panel ; Steel Plate, t1.6mm
- * Frame ; Steel, 60×40 & 50×30mm



WFE401 Series



- DH.WFE4011 Side Table, 1200×750×h800mm, with Cabinet(1ea)
- DH.WFE4012 , 1500×750×h800mm, with Cabinet(1ea)
- DH.WFE4013 , 1800×750×h800mm, with Cabinet(1ea)
- DH.WFE4014 , 2100×750×h800mm, with Cabinet(2ea)
- DH.WFE4015 , 2400×750×h800mm, with Cabinet(2ea)

※ Note! Freight and Installation Charges Separately.

<L4>


Ea. Price, ₩

(2) Assembly Shelf-type Side Table, with Transfer Cabinet & Shelves, Steel-Frame & -Side Panel, Phenol Work Top

- * Work Top ; Phenol Plate, t20mm
- * Back Panel ; Steel Plate, t1.6mm
- * Shelves ; PVC Shelf Plate(t8mm) & Steel / Aluminum Shelf Supporter(t1.6mm)



WFE402 Series



- DH.WFE4021 Side Table, 1200×750×h800/1900mm, with Cabinet / Shelves
- DH.WFE4022 , 1500×750×h800/1900mm, with Cabinet / Shelves
- DH.WFE4023 , 1800×750×h800/1900mm, with Cabinet / Shelves
- DH.WFE4024 , 2100×750×h800/1900mm, with Cabinet / Shelves
- DH.WFE4025 , 2400×750×h800/1900mm, with Cabinet / Shelves

※ Note! Freight and Installation Charges Separately.

Ea. Price, ₩

DAIHAN™ Laboratory Basic Corner Table, High Quality Steel-Frame & Phenol Work Top

Custom-Made available, 실험실용 코너 실험대, 고품질 스틸 프레임, 내열성/내충격성/내화학성 페놀 상판, 주문제작 가능



- Supplied by DAIHAN Korea
- Phenol Work Top ;
Good for Heat-/Shock-Resistant, Excellent Chemical-Resistant(Acetic Acid, Alcohol, Hydrochloric Acid, Sulfuric Acid, Organic Chemicals, Base, Biological Stains, Salt, Solvent, etc.)
- High Quality Steel-Frame, Side-/Back-Steel Panel(t1.6mm)
- 벽모서리 연결을 통해 넓은 실험공간 확보
- * Frame ; Steel Pipe 60×40 & 50×30mm, Steel Plate 1.6T Welding & Powder Finished
- * Chemical-Resistance Phenolic Laminate Work Top, t20mm (Black or Gray)
- * Non-standard orders can be made upon customer request

- DH.WFE3027 Furniture Corner Table, w1000×d1000×h800mm

<L4>

Ea. Price, ₩

Articles	Cat. No	Description						
Examining Stool, with Soft Cushion and Footrest, Adjustable Height, Swivel, Durable and Stable, $\phi 350 \times h450 \sim 600\text{mm}$ Ideal for Laboratory & Hospital, With or Without Caster, 실험 · 진찰 의자/결상								
		■ Supplied by DAIHAN Korea ■ 고급 인조가죽으로 제작되어 내구성이 좋고 편안하고 폭신한 쿠션으로 착석감이 좋음 * Adjustment Height Up to 600mm * Pre-Assembled Control * Light Weight * Easy to Move * without Back Cushion * Comfort & Soft Cushion						
(1) With Caster	(2) Without Caster	<L3> <table><tr><th>Ea.</th><th>Price, ₩</th></tr><tr><td>C9.4461</td><td>(1) Examining Stool, With Caster, $\phi 350 \times h450 \sim 600\text{mm}$</td></tr><tr><td>C9.4462</td><td>(2) , Without Caster, $\phi 350 \times h450 \sim 600\text{mm}$</td></tr></table>	Ea.	Price, ₩	C9.4461	(1) Examining Stool, With Caster, $\phi 350 \times h450 \sim 600\text{mm}$	C9.4462	(2) , Without Caster, $\phi 350 \times h450 \sim 600\text{mm}$
Ea.	Price, ₩							
C9.4461	(1) Examining Stool, With Caster, $\phi 350 \times h450 \sim 600\text{mm}$							
C9.4462	(2) , Without Caster, $\phi 350 \times h450 \sim 600\text{mm}$							

Ergonomic Desk Chair, with Mesh Back and Soft Cushion, Adjustable Height, Swivel, Stable, $600 \times 480 \times h1030 \sim 1110\text{mm}$
Ideal for Office, Laboratory, Home &c., with or without Armrest, with Caster, 사무용 의자/결상

A. With Armrest

B. Without Armrest

■ Supplied by DAIHAN Korea

■ 인체공학적 설계로 자세를 바르게 잡아줌으로써 피로감이 적음

■ 메쉬 소재로 된 등받이를 사용하여 내구성을 높이고 통기성을 향상시킴

■ 고급 패브릭원단으로 제작되어 편안하고 폭신한 쿠션으로 착석감이 좋음

* Adjustment Height Up to 1110mm

* Light Weight

* Easy to Move

* Mesh Back Cushion

* Soft Cushion

<L3>

A. With Armrest

DH.Cha3001 (1) Desk Chair, w/ Armrest,

600×480×h1030~1110mm, Black

DH.Cha3002 (2) Desk Chair, w/ Armrest,

600×480×h1030~1110mm, Blue

Ea.	Price, ₩

B. Without Armrest

DH.Cha3005 (3) Desk Chair, w/o Armrest,

600×480×h1030~1110mm, Black

DH.Cha3006 (4) Desk Chair, w/o Armrest,

600×480×h1030~1110mm, Blue

Ea.	Price, ₩

Folding Office Chair, with Mesh Back · Soft Cushion · Armrest · Caster, Stable, $585 \times 610 \times h650/860\text{mm}$
Ideal for Office, Laboratory, Home &c., Waterproof, 사무용 의자/결상, 접이식



(1) Black



(2) Gray

- Supplied by DAIHAN Korea
- 메쉬 소재로 된 등받이를 사용하여 내구성을 높이고 통기성을 향상시킴
- 두꺼운 스펀지에 방수기능이 있는 고급 패브릭원단으로 제작하여 착석감이 좋음
- 좌판이 위로 올라가서 접혀 수평적재가 가능

* 높이(좌석, 팔걸이, 전체) ; 450, 650, 860mm

* Mesh Back Cushion * Folding Seat

* Light Weight * Easy to Move

* Soft Cushion * Waterproof



Folding Seat

<L3>

Ea.	Price, ₩
DH.Cha3011	(1) Office Chair, w/ Armrest, 585×610×h650/860mm, Black
DH.Cha3012	(2) , Gray



Folding Seat

DAIHAN™ 10mm PMMA Glove Box Complete Set, with Natural Rubber Glove and Transfer Chamber, for 1 Operator
 With Internal Outlet, Foot S/W for Internal Pressure Regulation, Vacuum Pressure 1 Torr, 표준형 글러브박스/음압박스, 1인용

3-Year Warranty



150Lit & 179Lit



(1) Natural Rubber Gloves(Included)

- Supplied by DAIHAN Korea
- Compact, Portable, Lightweight and Self-contained
- Ergonomic Design to Fit the Average Person
- For Isolating Sensitive Research Studies from a Hostile Exterior Environment
- The Front Slopes Viewing Panel Minimizes Glare
- * Light Weight and Easy to move from Lab to Lab
- * Excellent Viewing & Air Tight Structure
- * Vacuum Gauge on Transfer Chamber
- * Main Chamber Foot Switch for Control Internal Pressure, Humidity and Oxygen

Specifications & Ordering information

<L2>

Spec.		Model	"GLOB150" 150 Lit	"GLOB180" 179 Lit
Main Chamber (w×d×h)	Outer D.		760×450×550mm	800×500×550mm
	Inner D.		740×430×530mm	780×480×530mm
	Thickness		10mm	
Transfer Chamber (w×d×h)	Outer D.		290×180×180mm	
	Inner D.		250×140×140mm	
	Thickness		20mm	
Material			PMMA(Polymethyl Methacrylate Resin)	
Vacuum Gauge			-1 to 1 Bar	
Vacuum Pressure			Pass Box : 1Torr	
Vacuum Outlet			O.D. ϕ9mm	
Gas Inlet			One touch Fitting, 1Ea.(I.D. ϕ6mm)	
Electric Outlet			2×220V and 2×110V	
Gloves			Natural Rubber Gloves, L640mm, 0.4mm thick	
Glove Port			2 Ports, 8"(ϕ203mm)	
Weight			32 kg	35 kg
Cat. No.			DH.GloB150	DH.GloB180
Price, Each. ₩				


Spare Gloves

SL.Glo6623

(1) Natural Rubber Glove, 8"(ϕ200mm), L780mm, 0.6mm thick

Pair. Price, ₩

※ Detailed Parts of Glove Box



DAIHAN™ Premium 15mm PMMA Heavy Duty Glove Box Complete Set, Natural Rubber Glove & Transfer Chamber, 1 & 2 Operators
With Internal Outlet, Foot S/W for Internal Pressure Regulation, Vacuum Pressure 1 Torr, 프리미엄 안전 글러브박스/음압박스, 1 & 2인용



182L, for 1 Operator

3-Year Warranty

- Supplied by DAIHAN Korea
- Compact, Portable, Lightweight and Self-contained
- Ergonomic Design to Fit the Average Person
- For Isolating Sensitive Research Studies from a Hostile Exterior Environment
- The Front Slopes Viewing Panel Minimizes Glare
- * Light Weight and Easy to move from Lab to Lab
- * Excellent Viewing & Air Tight Structure
- * Vacuum Gauge on Transfer Chamber
- * Main Chamber Foot Switch for Control Internal Pressure, Humidity and Oxygen

Specifications & Ordering information

<L2>

Spec.		Model	1 Operator type "GLOB190" 182 Lit	2 Operators type "GLOB420" 416 Lit
Main Chamber (w×d×h)	Outer D.		800×500×550mm	1530×530×630mm
	Inner D.		770×470×530mm	1500×500×600mm
	Thickness		15mm	
Transfer Chamber (w×d×h)	Outer D.		390×250×250mm	
	Inner D.		350×210×185mm	
	Thickness		20mm	
Material			PMMA(Polymethyl Methacrylate Resin)	
Vacuum Gauge			-1 to 1 Bar	
Vacuum Pressure			Pass Box ; 1Torr	
Vacuum Outlet			O.D. ϕ9mm	
Gas Inlet			One touch Fitting, 1Ea.(I.D. ϕ6mm)	
Electric Outlet			2×220V and 2×110V	
Gloves			Natural Rubber Gloves, L640mm, 0.4mm thick	
Glove Port			2 Ports, 8"(ϕ203mm)	4 Ports, 8"(ϕ203mm)
Weight			51 kg	86 kg
Cat. No.			DH.GloB190	DH.GloB420
Price, Each. ₩				



416L, for 2 Operators

3-Year Warranty



(1) Natural Rubber Gloves(Included)

※ Detailed parts of Glove Box
⇒ See the Previous Page.

Spare Gloves

Pair. Price, ₩

SL.Glo6623 (1) Natural Rubber Glove, 8"(ϕ200mm), L780mm,
0.6mm thick

Hypalon Gloves For Glove Box, for Port ϕ200/203, Length 800mm
Ideal for Acid, Thickness 0.4mm, 글러브 박스용 장갑

- 뛰어난 내약품성
- 재질 ; 겉면-하이퍼론, 안면-니트릴 고무
- Chemical Resistant - Hydrofluoric Acid / Aqua Regia / Strong Sulfuric Acid
- * Material ; Front / Hypalon, Back / Nitrile Rubber

<L2>

SL.Glo6601 "K-82" Gloves For Glove Box, Port ϕ200/203mm

Pair. Price, ₩



Natural Rubber Gloves for Glove Box, for Port ϕ130~200mm, Length 630~780mm
Powdered, Thickness 0.4 & 0.6mm, 천연 고무 글러브박스용 장갑

- 재질 ; 천연고무

<L2>

SL.Glo6621 Gloves For Glove Box, Thickness 0.4mm, ϕ130×L630mm
SL.Glo6622 , Thickness 0.6mm, ϕ150×L700mm
SL.Glo6623 , ϕ200×L780mm

Pair. Price, ₩



DAIHAN™ **Premium** 40mm PMMA Vacuum Glove Box Complete Set, Natural Rubber Glove & Transfer Chamber, for 1 Operator
With Internal Outlet, Foot S/W for Internal Pressure Regulation, Vacuum Pressure 10 Torr, 프리미엄 진공 클러브박스/음압박스, 1인용

3-Year Warranty



212L Vacuum Glove Box Complete Set



(1) Natural Rubber Gloves(Included)

- Supplied by DAIHAN Korea
- Ergonomic Design to Fit the Average Person
- Compact, Portable, Lightweight and Self-contained
- For Isolating Sensitive Research Studies from a Hostile Exterior Environment
- The Front Slopes Viewing Panel Minimizes Glare

- * Light Weight and Easy to move from Lab to Lab
- * Excellent Viewing & Air Tight Structure
- * Vacuum Gauge on Main Chamber and Transfer Chamber
- * Foot S/W for Internal Pressure Regulation
- * The Main Chamber Foot Switch can be Used to Effectively Control Internal Humidity and Oxygen

Specifications & Ordering information

<L2>


Spec.		Model	"GLOB220" 212 Lit
Main Chamber (w×d×h)	Outer D.		1,000×580×700mm
	Inner D.		840×500×550mm
	Thickness		40mm
Transfer Chamber (w×d×h)	Outer D.		390×250×250mm
	Inner D.		350×210×185mm
	Thickness		20mm
Material			PMMA(Polymethyl Methacrylate Resin)
Vacuum Gauge			-1 to 1 Bar
Vacuum Pressure			Main Chamber ; 1Torr, Pass Box ; 1Torr
Vacuum Outlet			O.D. ϕ9mm
Gas Inlet			One touch Fitting, 1Ea.(I.D. ϕ6mm)
Electric Outlet			2×220V and 2×110V
Gloves			Natural Rubber Gloves, L640mm, 0.4mm thick
Glove Port			2Ports, 8"(ϕ203mm)
Weight			137 kg
Cat. No.			DH.GloB220
Price, Each. ₩			

Spare Gloves

Pair. Price, ₩

SL.Glo6623 (1) High Grade Natural Rubber Glove, 8"(ϕ200mm), L780mm, 0.6mm thick

※ Detailed Parts of Vacuum Glove Box



DAIHAN™ High Temp Bath Circulator 250°C & Refrigerated/Heating Bath Circulator -35°C + 200°C

정밀 고온 배스겸용 써클레이터와 냉동/가열 써클레이터

CE



up to 250°C

External / Internal Heated

CE



-35°C + 150°C

External / Internal Heated / Refrigerated

CE



-30°C + 200°C

External Heated / Refrigerated



-20°C + 200°C

※ See the "Bath Circulators" Section for the Details.

- ▶ External/Internal Heated ;
250°C High-temp Bath Circulator, 8 · 12 · 22 · 30Lit
- ▶ External/Internal Heated/Refrigerated ;
-35°C + 150°C Refrigerated/Heating Bath Circulator, 8 · 12 · 22 · 30 Lit
-25°C + 150°C Refrigerated/Heating Bath Circulator, 8 · 12 · 22 · 30 Lit
- ▶ External Heated/Refrigerated ;
-30°C + 200°C Refrigerated/Heating Circulator, Fill-10 · 18 · 25Lit
-20°C + 200°C Refrigerated/Heating Circulator, Fill-10 · 20 · 30 · 50Lit



SK™ -40°C + 400°C Best COOL-/HEAT-ing Transfer Fluid, Synthetic Silicone Mixture, 20Lit
For Low/High-Temp Baths, Clearly Harmless Liquid, 냉매/열매 겸용 실리콘계 합성 열매체유

■ Supplied by SciLab Korea

* Colorless and Transparent Liquid.

■ The Best Circulation Properties / Thermal Stability.

* Do not use for Silicone Tube

(1) SYN THERM 7008®, -40°C + 400°C, 저/고온용

* Kinematic Viscosity, cSt ⇒ (1) -40°C ; 51.56cSt, (2) 0°C ; 16.08cSt, (3) 100°C ; 3.43cSt, (4) 200°C ; 1.36 cSt <L1>

* Specific Gravity ; 0.9435 at 15°C * Flash Point(PMCC) ; 168°C * Pour Point ; -72°C

SL.Syn7008 Cool/Heat Transfer Fluid, -40+400°C, 20Lit

Pk.	Price, ₩
20Lit	

(2) SUPER THERM FG-20®, -10°C + 300°C, 고온용

* Food Grade Thermal Medium Oil Satisfying FDA, NSF Standards

* Kinematic Viscosity, cSt ⇒ (1) 25°C ; 35.65cSt, (2) 100°C ; 4.23cSt, (3) 200°C ; 1.27 cSt

* Specific Gravity ; 0.8374 at 15°C * Flash Point ; 230°C * Pour Point ; -15°C

SL.Fg0020 Cool/Heat Transfer Fluid, -10+300°C, 20Lit

Pk.	Price, ₩
20Lit	



SHL™ SYN THERM 22L® -5°C + 320°C Heating Transfer Fluid

Synthetic Alkylbenzene Derivative, for Low- & High-Temp Circulation, 알킬벤젠계 합성열매

■ Supplied by SciLab Korea

* Hazards Identification of MSDS(Material Safety Data Sheet) ;

- H315 Causes Skin Irritation

- H319 Causes Serious Eye Irritation

- H335 May Cause Respiratory Irritation

- H336 May Cause Drowsiness or Dizziness

* Kinematic Viscosity & Density ;

(1) at 0°C ; 22.7cSt, 0.885kg/ℓ

(2) at 40°C ; 21.80cSt, 0.859kg/ℓ

(3) at 80°C ; 6.09cSt, 0.834kg/ℓ

(4) at 140°C ; 6.09cSt, 0.795kg/ℓ

(5) at 200°C ; 1.10cSt, 0.756kg/ℓ

(6) at 260°C ; 0.74cSt, 0.718kg/ℓ <L1>

SL.SYN22L SYN THERM 22L® Heating Transfer Fluid, -5°C + 320°C, 20Lit

Pk.	Price, ₩
20Lit	

DAIHAN™ High Performance Compact Heating Block/Dry Bath Incubator “HB-48”, Up to 150℃, ±0.1℃
With Large LCD Display, Fuzzy Control, Molded Heater, Modular Anodized Aluminum Blocks, Touch Button Controller
히팅 블록, 디지털 퍼지 컨트롤 시스템, 터치버튼식 조절기, 일체형 주물 히터, 블록 교체 가능, 블록 추출용 공구 포함



Block “BLC548”,
48 holes

~150℃



3-Year Warranty

Block Removal Tool

Heating Block/Dry Bath “HB-48-SET”


- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Designed for Multi-functional Dry Block Heater
- Ergonomic Design
- Patented Jog-Dial and Main Touch Button - Best Durability, Best Mobility, Best Convenient to Operate with Large Back-Light LCD
- With a Screw-type Block Removal Tool ; Convenient Handling
- * High-performance Molded Heater ensures Best Temperature Accuracy.
- * Internal Bath Material(PPS) ensures Best Heat Insulation
- * Rapid heat up time.
- * Ideal for Enzyme Reactions, Inactivation of Sera, Incubation and Other Laboratory Procedures
- * Precise Temperature Control ; ±0.1℃ at 37℃, ±0.2℃ at 90℃, ±0.4℃ at 135℃
- * User Self-Compensation Function to offset between Actual and Displayed Temperature. (±10.0℃)
- * Digital Timer Included ; up to 99hr 59min (Operation Time + Delay Time)
- * Storage Function of Temperature and Timer setting values
- * Convenient integral timer for time-sensitive incubations
- * Alarm Function ; Error Status and Timer-end
- * Easy to Read, Wide LCD Screen
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled)
- * Easily calibrated for accurate temperature control and display

Specifications & Ordering information

<L1>

Spec.	Model	“HB-48” (Only Unit)	“HB-48-SET” (includes Standard Block)
Controller		Single Digital Fuzzy Controller	
Useful Block		Interchangeable Block(90×120×h60mm) for Various Tubes	
included Standard Block		None(Only Heating Block Unit)	48 holes×for 1.5 / 2ml Centri. Tube
Heater		Molded Heater(300W)	
Dimension		Unit ; 178×226×h138mm, Block ; 90×120×h60mm	
Temp Range, & Accuracy		Ambient +5℃~150℃, ±0.1℃ at 37℃, ±0.2℃ at 90℃, ±0.4℃ at 135℃	
Heat-Up Speed & Sensor		Approx. 5℃/min, PT100	
Safety Circuit		Over Temp and Over Current Protector, Sensor Error Detector	
Timer with Alarm		99hr 59min, Alarm (Error Status and Timer - end)	
Display		Digital LCD with Back-Light Function	
Material		Unit ; PPS & BAS(HF380), Block ; Anodized Aluminum	
Included Standard Acce.		Block Extraction Tool	
Packing Size & Gross Weight		230×290×h245mm, 2.2kg	230×355×h245mm, 3.5kg
Power Consumption		300 W	
Electric Rqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	
Cat. No.	230V	DH.WHB02000	DH.WHB02100
	120V	DH.WHB02001	DH.WHB02101
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.




Controller of “HB-48”



Included Block (3) “BLC 548” 48 Holes
for 1.5/2ml & Removal Tool

Optional Interchangeable Anodized Aluminum Block (90×120×h60mm)



- DH.WHB510596 (1) for 0.2ml PCR Tube, (φ6×h16mm)×96 holes, “BLC596”
- DH.WHB510570 (2) for 0.5ml Centri. Tube, (φ7.7×h29.5mm)×63 holes, “BLC570”
- DH.WHB511548 (3) for 1.5/2ml Centri. Tube, (φ11×h35mm)×48 holes, “BLC548”
(* included in the “HB-48-SET”)
- DH.WHB521248 (4) for φ12mm Test Tube, (φ12×h35mm)×48 holes, “BLT248”
- DH.WHB521335 (5) for φ13mm Test Tube, (φ13×h47mm)×30 holes, “BLT335”
- DH.WHB521535 (6) for φ15mm Test Tube, (φ15×h47mm)×30 holes, “BLT535”
- DH.WHB521628 (7) for φ16.5mm Test Tube, (φ16.5×h47mm)×24 holes, “BLT628”
- DH.WHB521828 (8) for φ18mm Test Tube, (φ18×h47mm)×24 holes, “BLT828”
- DH.WHB531515 (9) for 15ml Conical Tube, (φ16×h47mm)×15 holes, “BLC515”
- DH.WHB535008 (10) for 50ml Conical Tube, (φ28.7×h47mm)×8 holes, “BLC008”
- DH.WHB500102 (11) Spare ; Block Removal Tool, Easy Handling

Ea. Price, ₩

DAIHAN™ High Performance Dual Heating Block/Dry Bath Incubator "HB-96D", Up to 150°C, ±0.1°C

With Fuzzy Control, Molded Heater, Modular Anodized Aluminum Blocks, Acrylic Lid, Touch Button Controller

히팅 블록, 디지털 퍼지 컨트롤 시스템, 듀얼 조절기/블록, 일체형 주물 히터, 블록 교체 가능, 아크릴 리드 & 블록용 집게 포함

CE

~150°C

2×Block "BLC548",
2×48 holes



Dual Control System

Block Cover

Tongs "BLE101"

3-Year Warranty

Heating Block/Dry Bath, Dual Control System
"HB-96D-SET"



(3) "BLC548"



(11) "BLE101"

Included Anodized Aluminum Block & Tongs

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Designed for Multi-functional Dry Block Heater
- Dual Controller & Block

■ Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function

- * Ideal for Enzyme Reactions, Inactivation of Sera, Incubation and Other Laboratory Procedures
- * Precise Temperature Control ; ±0.1°C at 37°C, ±0.2°C at 90°C, ±0.4°C at 135°C
- * User Self-Compensation Function to offset between Actual and Displayed Temperature. (±10.0°C)
- * Over Temperature & Over Current Protection and Sensor Error Detection
- * Digital Timer Included ; up to 99hr 59min (Operation Time + Delay Time)
- * Storage Function of Temperature and Timer setting values
- * Modular Block
- * Alarm Function ; Error Status and Timer-end
- * Various kind of Blocks are offered (Optional)
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * The Block Extraction Tongs and Included for Easy Handling

Specifications & Ordering information

<L1>


Spec.	Model	"HB-96D" (Only Unit)	"HB-96D-SET" (included Standard Block)
Controller		Dual Digital Fuzzy Controller	
Useful Block		Interchangeable Block(90×120×h60mm) for Various Tubes	
included Standard Block		None(Only Heating Block Unit)	2×48 holes×for 1.5/2ml Centri. Tube
Heater		Molded Heater(300W)×2ea	
Dimension(mm)		Unit ; 285×250×h155, Block ; (90×120×h60)×2ea	
Temp Range, & Accuracy		Ambient +5°C to 150°C, ±0.1°C at 37°C, ±0.2°C at 90°C, ±0.4°C at 135°C	
Heat-Up Speed & Sensor		Approx. 5°C/min, PT100	
Safety Circuit		Over Temp and Over Current Protector, Sensor Error Detector	
Timer with Alarm		99hr 59min, Alarm (Error Status and Timer - end)	
Display		Digital LCD with Back-Light Function	
Material		Unit ; Stainless-steel (#304) & Powder Coated steel, Block ; Anodized Aluminum	
Included Standard Acce.		Block Tongs "BLE101" and Block Cover	
Packing Size & Gross Weight		400×270×h200 mm, 6kg	400×270×h200 mm, 7.2kg
Power Consumption		600 W	
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	

Cat. No.	230V	DH.WHB00198	DH.WHB00199
	120V	DH.WHB02198	DH.WHB02199
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Optional Interchangeable Anodized Aluminum Block (90×120×h60mm)

Ea. Price, ₩



- DH.WHB510596 (1) for 0.2ml PCR Tube, (φ6×h16mm)×96 holes, "BLC596"
- DH.WHB510570 (2) for 0.5ml Centri. Tube, (φ7.7×h29.5mm)×63 holes, "BLC570"
- DH.WHB511548 (3) for 1.5/2ml Centri. Tube, (φ11×h35mm)×48 holes, "BLC548"
(+ included in the "HB-96D-SET")
- DH.WHB521248 (4) for φ12mm Test Tube, (φ12×h35mm)×48 holes, "BLT248"
- DH.WHB521335 (5) for φ13mm Test Tube, (φ13×h47mm)×30 holes, "BLT335"
- DH.WHB521535 (6) for φ15mm Test Tube, (φ15×h47mm)×30 holes, "BLT535"
- DH.WHB521628 (7) for φ16.5mm Test Tube, (φ16.5×h47mm)×24 holes, "BLT628"
- DH.WHB521828 (8) for φ18mm Test Tube, (φ18×h47mm)×24 holes, "BLT828"
- DH.WHB531515 (9) for 15ml Conical Tube, (φ16×h47mm)×15 holes, "BLC515"
- DH.WHB535008 (10) for 50ml Conical Tube, (φ28.7×h47mm)×8 holes, "BLC008"
- DH.WHB500101 (11) Spare ; Block Extraction Tongs, "BLE101"

DAIHAN™ Precise Heating • Mixing • Cooling Block/Dry Bath Incubator “HB-M24” & “HB-MR24”, Mixing Freq. 200~1500rpm With Standard Aluminum Block, Magnet Adhesion Technology, Lid for Heat Preservation, 15℃~100℃ (HB-M24) / 0.1℃~100℃ (HB-MR24), ±0.5℃ 히팅 • 믹싱 • 쿨링 블록, TFT 디스플레이, 마그넷 블록 간편 탈부착 기능, 기본 알루미늄 블록(1.5/2ml Tubes × 24 holes용) 포함



■ Supplied by DAIHAN Korea ■ CE Certified

I. Heating & Mixing Block/Dry Bath Incubator “HB-M24”

Temp Setting Range ; 15℃~100℃

II. Heating • Mixing & Cooling Block/Dry Bath Incubator “HB-MR24”

Temp Setting Range ; 0.1℃~100℃

- 정확하고 효율적인 가열 및 혼합은 물론 탁월한 온도 정확성과 균일성을 보장
- 임상 검사실의 시료 전처리, 효소-기질 반응, 생화학 테스트, 혈청 응고 테스트 등 다양한 용도에 적합.
- 공구의 도움 없이 자기접착기술로 쉽고 빠르게 블록 교체를 가능하게 함
- Allow for Precise and Efficient Heating and Mixing, as well as Excellent Temperature Accuracy and Uniformity
- Suitable for Clinical Laboratories, for Sample Preparation, Enzyme-substrate Reactions, Biochemistry Tests, Serum Coagulation Tests and Many Other.
- Precise Temperature Control for Heating ■ Flexible Adapter Selection.
- Quick Interchange of Blocks with Magnet Adhesion Technology without any Tools
- * Programmable * Excellent Mixing Result
- * Stable and Wide Range of Speed Adjustment * Compatible with Various Tubes
- * Block Equipped with a Lid for Heat Preservation * Three-point Temperature Calibration
- * Including Heating Block for 1.5 / 2ml Tubes (24-Holes)

3-Year Warranty

I. Heating & Mixing

II. Heating • Mixing & Cooling

Spare Block



(1) 1.5/2ml Tubes × 24 holes
(Included Standard Block)

Optional Block



Specifications & Ordering information

<L1>

Spec.	Model	I. Heating & Mixing “HB-M24”	II. Heating • Mixing & Cooling “HB-MR24”
Standard Block		1.5/2ml Tubes × 24 holes	
Mixing Frequency		200 ~ 1500rpm	
Included Standard Block		1.5 / 2ml Tubes × 24 holes, 1500rpm, hole ϕ 10.9 × 30.8mm	
Temp Range		Room Temp ~ 100℃	Heating ; Room Temp ~ 100℃ Cooling ; Room Temp ~ Below Room Temp 15℃
Temp Setting Range		15℃ ~ 100℃	0.1℃ ~ 100℃
Temp Control Accuracy		±0.5℃ at 20℃ ~ 45℃	
Temp Uniformity		±0.5℃ at 20℃ ~ 45℃	
Max. Heating Rate		5.5℃/min	
Max. Cooling Rate		-	5℃/min (100℃ ~ Room Temp) 0.5℃/min (Below Room Temp)
Mixing Orbit		ϕ 3mm	
Time Setting Range		0 ~ 99h 30min	
Program		9	
Overheating Protection		150℃	
Display		TFT (Thin Film Transistor)	
Adapter Block Material		Aluminum	
Voltage & Frequency		AC220V-240V & 50/60Hz	
Power		200W	
Dimension (W×D×H)		200 × 235 × h120mm (without Heating Block)	
Weight		7.3kg	
Cat. No.		DH.HBM24	DH.HBMR24
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Spare Block

DH.HBM24.B1 (1) 1.5 / 2ml × 24 holes Tubes, 1500rpm, hole ϕ 10.9 × 30.8mm

Optional Block

DH.HBM24.B2 (2) 0.5ml Tubes × 24 holes, 1500rpm, hole ϕ 8.1 × 26mm
 DH.HBM24.B3 (3) 5ml Tubes × 8 holes(Round Bottom), 1500rpm, hole ϕ 13.6 × 42mm
 DH.HBM24.B4 (4) 5ml Tubes × 8 holes(Cone Bottom), 1400rpm, hole ϕ 17 × 51.5mm
 DH.HBM24.B5 (5) 15ml Tubes × 8 holes, 800rpm, hole ϕ 16.9 × 80mm
 DH.HBM24.B6 (6) 50ml Tubes × 4 holes, 600rpm, hole ϕ 29.6 × 80mm
 DH.HBM24.B7 (7) PCR Plate, 0.2ml Tubes × 96 holes, 1500rpm
 DH.HBM24.B8 (8) Microplate, 1500rpm
 DH.HBM24.B9 (9) 1.5ml Chromatography Vial × 24 holes, 1500rpm, hole ϕ 12 × 20.5mm **New**
 DH.HBM24.B10 (10) Deep well Plate **New**

Ea. Price, ₩

DAIHAN™ High Performance Heating · Cooling Block/Dry Bath Incubator “HB-R48”, $-5^{\circ}\text{C} + 95^{\circ}\text{C}$, $\pm 0.1^{\circ}\text{C}$

With Digital PID Control, Acrylic Lid, Precise Temperature Control, Modular Anodized Aluminum Block

Peltier Cooling System, Digital Timer, Interchangeable Blocks, 히팅 / 쿨링 블록, 블록 교체 가능, 아크릴 리드 & 블록용 집게 포함

CE

$-5^{\circ}\text{C} + 95^{\circ}\text{C}$

Block “BLC548”,
48 holes



Patented Jog-Shuttle
Control System

Block Cover

Tongs “BLE101”

3-Year Warranty

“HB-R48-SET”

with Block “BLC548” for 1.5/2ml Centrifuge Tube $\times$ 48 holes



(3) “BLC548”



(11) “BLE101”

Included Anodized Aluminum Block & Tongs

- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

- Designed for Multi-Function ; Enzyme Reactions, Inactivation of Sera, Incubation and Other Laboratory Procedures.

- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function

* User Self-Compensation Function to offset between Actual and Displayed Temperature. ($\pm 10.0^{\circ}\text{C}$)

* Storage Function of Temperature and Timer setting values

* Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

Specifications & Ordering information


<L1>

Spec.	Model	“HB-R48-SET” (included Standard Block)
Controller		Digital Fuzzy Controller
Useful Block		Interchangeable Block(90×120×h60mm) for Various Tubes
included Standard Block		48 holes×for 1.5 / 2ml Centri. Tube
Dimension ; Unit & Block		185×230×h220mm, and 90×120×h60mm
Cooling		Two Peltier Modules
Cool-Down Speed		2.3°C/min. between 95°C to 25°C, 0.5°C/min. between 25°C to -5°C
Heater and Sensor		300W (Molded Heater), PT100
Temp Range, & Accuracy		-5°C to 95°C, and $\pm 0.1^{\circ}\text{C}$ at 37°C, $\pm 0.2^{\circ}\text{C}$ at 90°C
Heat-Up Speed		Approx. 5°C/min
Safety Device		Over Temp and Over Current Protector, Sensor Error Detector
Timer with Alarm		99hr 59min, Alarm (Error Status and Timer - end)
Display		Digital LCD with Back-Light Function
Material	Internal & External Block	Stainless steel (# 304), and Powder Coated steel Anodized Aluminum
Included Standard Acce.		Block Tongs “BLE101” and Block Cover
Packing Size, Gross Weight		220×350×h310mm, 10.2kg
Power Consumption		450 W
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz

Cat. No.	230V	DH.WHB00249
	120V	DH.WHB02249
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

Optional Interchangeable Anodized Aluminum Block (90×120×h60mm)



- DH.WHB510596 (1) for 0.2ml PCR Tube, ($\phi 6 \times h16\text{mm}$) $\times$ 96 holes, “BLC596”
- DH.WHB510570 (2) for 0.5ml Centri. Tube, ($\phi 7.7 \times h29.5\text{mm}$) $\times$ 63 holes, “BLC570”
- DH.WHB511548 (3) for 1.5/2ml Centri. Tube, ($\phi 11 \times h35\text{mm}$) $\times$ 48 holes, “BLC548”
(+ included in the “HB-R48-SET”)
- DH.WHB521248 (4) for $\phi 12\text{mm}$ Test Tube, ($\phi 12 \times h35\text{mm}$) $\times$ 48 holes, “BLT248”
- DH.WHB521335 (5) for $\phi 13\text{mm}$ Test Tube, ($\phi 13 \times h47\text{mm}$) $\times$ 30 holes, “BLT335”
- DH.WHB521535 (6) for $\phi 15\text{mm}$ Test Tube, ($\phi 15 \times h47\text{mm}$) $\times$ 30 holes, “BLT535”
- DH.WHB521628 (7) for $\phi 16.5\text{mm}$ Test Tube, ($\phi 16.5 \times h47\text{mm}$) $\times$ 24 holes, “BLT628”
- DH.WHB521828 (8) for $\phi 18\text{mm}$ Test Tube, ($\phi 18 \times h47\text{mm}$) $\times$ 24 holes, “BLT828”
- DH.WHB531515 (9) for 15ml Conical Tube, ($\phi 16 \times h47\text{mm}$) $\times$ 15 holes, “BLC515”
- DH.WHB535008 (10) for 50ml Conical Tube, ($\phi 28.7 \times h47\text{mm}$) $\times$ 8 holes, “BLC008”
- DH.WHB500101 (11) Spare ; Block Extraction Tongs, “BLE101”

Ea. Price, ₩

DAIHAN™ Aluminum-case Beaker Heating Mantle, (1) Basic & (2) Stirring-type, 450℃, 100~5,000ml
With Built-in Temp Controller, with/without Mag-stir Speed Control, with Certi. & Traceability
비커용 히팅맨틀, 온도 조절기 내장, “기본형” 및 “자석교반형”, Ni-Cr열선 내장

- Supplied by DAIHAN Korea
- Ideal for All type Beakers
- Non-magnetic Aluminum Case
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Built-in Temp-Controller. Nickel-Chrome Heating Element “Coil-type” ; Acid Resistant & Long Life Time.

~450℃

3-Year Warranty



(1) Basic-type Beaker Mantle with Controller

(1) Basic-type Beaker Mantle, Built-in Temp. Controller

* With a Device* for Holding a Optional Stainless-steel Rod ϕ 12.7mm for Article Support.

<L1>

Beakers, ml Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
100ml, ϕ 50~52	80 W	45 mm	ϕ 150 $\times$ h145	DH.WHM12141	DH.WHM12151	
250ml, ϕ 68~70	140 W	60 mm	ϕ 150 $\times$ h145	DH.WHM12142	DH.WHM12152	
400ml, ϕ 79~81	240 W	60 mm	ϕ 150 $\times$ h145	DH.WHM12143	DH.WHM12153	
600ml, ϕ 88~90	240 W	75 mm	ϕ 150 $\times$ h145	DH.WHM12144	DH.WHM12154	
800ml, ϕ 98~100	350 W	90 mm	ϕ 200 $\times$ h170	DH.WHM12145	DH.WHM12155	
1,000ml, ϕ 105~108	350 W	90 mm	ϕ 200 $\times$ h170	DH.WHM12146	DH.WHM12156	
2,000ml, ϕ 130~132	450 W	115 mm	ϕ 200 $\times$ h195	DH.WHM12147	DH.WHM12157	
3,000ml, ϕ 151~153	550 W	140 mm	ϕ 270 $\times$ h240	DH.WHM12148	DH.WHM12158	
5,000ml, ϕ 169~171	800 W	160 mm	ϕ 270 $\times$ h240	DH.WHM12149	DH.WHM12159	

50~1200rpm
~450℃

3-Year Warranty



(2) Stirring-type Beaker Mantle with a Stir-Bar

(2) Stirring-type Beaker Mantle, Built-in 50~1200rpm Mag-stir Speed & Temp. Controller

* With a Device* for Holding a Optional Stainless-steel Rod ϕ 12.7mm for Article Support.

* With 1ea $\times$ PTFE Octagonal Magnetic Stirrer Bar, L38mm

Beakers, ml Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
100ml, ϕ 50~52	80 W	55 mm	ϕ 150 $\times$ h145	DH.WHM121601	DH.WHM121611	
250ml, ϕ 68~70	140 W	55 mm	ϕ 150 $\times$ h145	DH.WHM121602	DH.WHM121612	
400ml, ϕ 79~81	240 W	85 mm	ϕ 150 $\times$ h190	DH.WHM121603	DH.WHM121613	
600ml, ϕ 88~90	240 W	100 mm	ϕ 150 $\times$ h190	DH.WHM121605	DH.WHM121615	
1,000ml, ϕ 105~108	350 W	85 mm	ϕ 200 $\times$ h170	DH.WHM121606	DH.WHM121616	
2,000ml, ϕ 130~132	450 W	105 mm	ϕ 200 $\times$ h195	DH.WHM121607	DH.WHM121617	
3,000ml, ϕ 151~153	550 W	150 mm	ϕ 270 $\times$ h240	DH.WHM121608	DH.WHM121618	
5,000ml, ϕ 169~171	800 W	180 mm	ϕ 270 $\times$ h270	DH.WHM121610	DH.WHM121620	

(3) Stand Rod, “RD300”, ϕ 12.7 $\times$ L520mm

(4) Holder / Clamp / Clip for Temp. Probe

Optional Stainless-steel Rod, Holder & Clamp

Ea. Price, ₩

- DH.WHP503037 (3) Stand Rod, ϕ 12.7 $\times$ L520mm, “RD300”, for Heating Mantles
DH.WHP503012 (4) Holder / Clamp / Clip, “CL220”, for Temp. Sensor

DAIHAN™ Aluminum-case Beaker Heating Mantles, Remote Controlled, 450℃, 100~10,000ml

For Beaker, with Nickel Chrome Heating Element, K-type Thermo-sensor Integrated, Option-Controller, with Certi. & Traceability

비커용 히팅맨틀, K-type 온도센서와 Ni-Cr 열선 내장, 자력교반기와 사용 가능, 조절기와 외장 온도 프로브 별도

- Supplied by DAIHAN Korea ■ Ideal for Use with Magnetic Stirrer ■ 100~5,000ml ; Non-magnetic Aluminum Case, 10Lit ; Powder-coated Steel Case
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Built-in Temp-Controller. Nickel-Chrome Heating Element "Coil-type" ; Acid Resistant & Long Life Time.
- Optional Remote Thermo Controllers are available (Analog ; WHM-C10A, Digital ; WHM-C10D)

<L1>



Beakers, ml Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
100ml, ϕ 50~52	80 W	40 mm	ϕ 150 $\times$ h100	DH.WHM12170	DH.WHM12180	
250ml, ϕ 68~70	140 W	60 mm	ϕ 150 $\times$ h110	DH.WHM12171	DH.WHM12181	
400ml, ϕ 79~81	240 W	70 mm	ϕ 150 $\times$ h110	DH.WHM12172	DH.WHM12182	
400ml, ϕ 76~78	240 W	70 mm	ϕ 150 $\times$ h110	DH.WHM121721	DH.WHM121821	
600ml, ϕ 88~90	240 W	90 mm	ϕ 150 $\times$ h110	DH.WHM12173	DH.WHM12183	
800ml, ϕ 98~100	350 W	100 mm	ϕ 200 $\times$ h140	DH.WHM12174	DH.WHM12184	
1,000ml, ϕ 105~108	350 W	110 mm	ϕ 200 $\times$ h140	DH.WHM12175	DH.WHM12185	
2,000ml, ϕ 130~132	450 W	130 mm	ϕ 200 $\times$ h160	DH.WHM12176	DH.WHM12186	
3,000ml, ϕ 151~153	550 W	160 mm	ϕ 270 $\times$ h200	DH.WHM12177	DH.WHM12187	
5,000ml, ϕ 169~171	800 W	180 mm	ϕ 270 $\times$ h205	DH.WHM12179	DH.WHM12189	
10Lit, ϕ 216~219	1400 W	245 mm	ϕ 320 $\times$ h270	DH.WHM121791	DH.WHM121891	

Optional Remote Controller

Cat. No	Ea. Price, ₩
230V	120V
DH.WHM13500	DH.WHM23500
DH.WHM12500	DH.WHM22500

(a) Analog Remote Controller, 5~100% of Rated Voltage, Up to 1.2kW

(b) Digital Remote Controller, Ambient +5℃~700℃, $\pm 1^\circ\text{C}$ Uniformity at Set Temp.

Optional Temp. Probe & Stand Rod, Holder, Clamp for Digital Temp. Remote Controller

Cat. No	Ea. Price, ₩
DH.WHP503040	Direct Contact Temp. Probe, K-type, "SS300", L250 $\times$ ϕ 4.8 mm
DH.WHP503050	"SS400", L380 $\times$ ϕ 4.8 mm
DH.WHP503041	"SS310", Teflon Coated, L250 $\times$ ϕ 5 mm
DH.WHP503051	"SS410", Teflon Coated, L380 $\times$ ϕ 5 mm
DH.WHP503037	(1) Stand Rod, "RD300", ϕ 12.7 $\times$ L520mm, for Heating Mantles
DH.WHP503012	(2) Holder / Clamp / Clip, "CL220", for Temp. Sensor

- (1) Stand Rod, "RD300", ϕ 12.7 $\times$ L520mm
- (2) Holder / Clamp / Clip for Temp. Probe

DAIHAN™ Fabric-Housed Beaker Heating Mantles, Remote Controlled, with Nickel Chrome Heating Element, 450℃, 50~3,000ml

For Beaker Self-Standing, with Certi. & Traceability, Option-Controller, 비커용 직물케이스 히팅맨틀, Ni-Cr 열선 내장, 조절기와 외장 온도 프로브 별도

- Supplied by DAIHAN Korea ■ Ideal for All type Round Flask and Use with Magnetic Stirrer
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Nickel-Chrome Heating Element "Coil-type" ; Acid Resistant & Long Life Time.

<L1>



Beakers, ml Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
50ml, ϕ 41~43	80 W	35 mm	ϕ 80 $\times$ h50	DH.WHM12190	DH.WHM12200	
100ml, ϕ 50~52	80 W	40 mm	ϕ 90 $\times$ h55	DH.WHM12191	DH.WHM12201	
250ml, ϕ 68~70	140 W	60 mm	ϕ 120 $\times$ h75	DH.WHM12192	DH.WHM12202	
400ml, ϕ 79~81	240 W	70 mm	ϕ 130 $\times$ h85	DH.WHM12193	DH.WHM12203	
400ml, ϕ 76~78	240 W	70 mm	ϕ 130 $\times$ h85	DH.WHM121931	DH.WHM122031	
600ml, ϕ 88~90	240 W	85 mm	ϕ 160 $\times$ h100	DH.WHM12194	DH.WHM12204	
800ml, ϕ 98~100	350 W	90 mm	ϕ 170 $\times$ h105	DH.WHM12195	DH.WHM12205	
1,000ml, ϕ 105~108	350 W	105 mm	ϕ 180 $\times$ h120	DH.WHM12196	DH.WHM12206	
2,000ml, ϕ 130~132	450 W	130 mm	ϕ 200 $\times$ h145	DH.WHM12197	DH.WHM12207	
3,000ml, ϕ 151~153	550 W	150 mm	ϕ 220 $\times$ h165	DH.WHM12198	DH.WHM12208	

Optional Remote Controller

Cat. No	Ea. Price, ₩
230V	120V
DH.WHM13500	DH.WHM23500
DH.WHM12500	DH.WHM22500

(a) Analog Remote Controller, 5~100% of Rated Voltage, Up to 1.2kW

(b) Digital Remote Controller, Ambient +5℃~700℃, $\pm 1^\circ\text{C}$ Uniformity at Set Temp.

DAIHAN™ Aluminum-case Round Flask Heating Mantle, (1) Basic & (2) Stirring-type, 450℃, 50~20,000ml
With Built-in Temp Controller; with/without Mag-stir Speed Control, with Certi. & Traceability
라운드플라스크용 히팅맨틀, 온도 조절기와 Ni-Cr 열선 내장, “기본형” 및 “자석교반형”

- Supplied by DAIHAN Korea
- Ideal for All type Round Flask
- 50~5,000ml ; Non-magnetic Aluminum Case, 10Lit & 20Lit ; Powder-coated Steel Case
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Built-in Temp-Controller. Nickel-Chrome Heating Element “Coil-type” ; Acid Resistant & Long Life Time.

~450℃

3-Year Warranty



(1) Basic-type

Rod Device*

<Ex.>
Setting Figure;
(1) Basic Mantle
+Flask+Stirrer
+Stand+Access.


50~1200rpm
~450℃

3-Year Warranty


with
1×Oval-type
Stirrer Bar

Rod Device*

(2) Stirring-type Flask Mantle with a Stir-Bar

(1) Basic-type Flask Mantle, Built-in Temp. Controller

* Heating Mantles for 50ml~5,000ml Flask are with a Device* for Holding a Optional Stainless-steel Rod
φ 12.7 mm for Article Support.

<L1>

Flasks, ml Dia. φ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside φ × h(mm)	230V	120V	
50ml, φ 49 ~ 52	80 W	35 mm	φ 150 × h150	DH.WHM12010	DH.WHM12020	
100ml, φ 63 ~ 66	80 W	45 mm	φ 150 × h150	DH.WHM12011	DH.WHM12021	
250ml, φ 83 ~ 86	160 W	60 mm	φ 150 × h150	DH.WHM12012	DH.WHM12022	
500ml, φ 102 ~ 105	250 W	75 mm	φ 200 × h175	DH.WHM12013	DH.WHM12023	
1,000ml, φ 130 ~ 133	350 W	90 mm	φ 200 × h175	DH.WHM12014	DH.WHM12024	
2,000ml, φ 165 ~ 169	450 W	110 mm	φ 270 × h200	DH.WHM12015	DH.WHM12025	
3,000ml, φ 184 ~ 187	550 W	115 mm	φ 270 × h200	DH.WHM12016	DH.WHM12026	
5,000ml, φ 221 ~ 224	800 W	135 mm	φ 325 × h230	DH.WHM12017	DH.WHM12027	
10Lit, φ 278 ~ 282	1,400 W	190 mm	φ 400 × h240	DH.WHM12019	DH.WHM12029	
20Lit, φ 344 ~ 348	1,800 W	215 mm	φ 440 × h265	DH.WHM12070	DH.WHM12071	

(2) Stirring-type Flask Mantle, Built-in 50~1200rpm Mag-stir Speed & Temp. Controller

* Heating Mantles for 50ml~5,000ml Flask are with a Device* for Holding a Optional Stainless-steel Rod
φ 12.7 mm for Article Support.

* 10-/20-Lit models Square-type Outside Shape.

* With 1ea × PTFE Oval-type Magnetic Stirrer Bar, L25mm.

Flasks, ml Dia. φ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside φ × h(mm)	230V	120V	
50ml, φ 49 ~ 52	90 W	50 mm	φ 150 × h150	DH.WHM12030	DH.WHM12040	
100ml, φ 63 ~ 66	90 W	50 mm	φ 150 × h150	DH.WHM12031	DH.WHM12041	
250ml, φ 83 ~ 86	160 W	60 mm	φ 150 × h150	DH.WHM12032	DH.WHM12042	
500ml, φ 102 ~ 105	250 W	85 mm	φ 200 × h175	DH.WHM12033	DH.WHM12043	
1,000ml, φ 130 ~ 133	350 W	85 mm	φ 200 × h175	DH.WHM12034	DH.WHM12044	
2,000ml, φ 165 ~ 169	450 W	110 mm	φ 270 × h200	DH.WHM12035	DH.WHM12045	
3,000ml, φ 184 ~ 187	550 W	110 mm	φ 270 × h200	DH.WHM12036	DH.WHM12046	
5,000ml, φ 221 ~ 224	800 W	135 mm	φ 325 × h230	DH.WHM12037	DH.WHM12047	
10Lit, φ 278 ~ 282	1,400 W	190 mm	400×500×h265	DH.WHM8060	DH.WHM8060.1	
20Lit, φ 344 ~ 348	1,800 W	215 mm	440×540×h290	DH.WHM8061	DH.WHM8061.1	

Optional Stainless-steel Rod, Holder & Clamp

Ea. Price, ₩

- DH.WHP503037 (3) Stand Rod, φ 12.7×L520mm, “RD300”, for Heating Mantles
DH.WHP503012 (4) Holder / Clamp / Clip, “CL220”, for Temp. Sensor

(3) Stand Rod, “RD300”, φ 12.7×L520mm

(4) Holder / Clamp / Clip for Temp. Probe

DAIHAN™ Aluminum-case Remote Controlled Round Flask Heating Mantle, (1) 450°C Basic & (2) 650°C Hi-Temp, 50ml~100Lit
With Nickel Chrome Heating Element, K-type Thermo-Sensor Integrated, with Certi. & Traceability, Option-Controller
라운드플라스크용 히팅맨틀, K-type 온도센서와 Ni-Cr 열선 내장, 자력교반기와 사용가능, 조절기와 외장 온도 프로브 별도

- Supplied by DAIHAN Korea ■ Ideal for All type Round Flask and Use with Magnetic Stirrer
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Optional Remote Thermo Controller, K-type Thermo-Sensor, and Nickel-Chrome Heating Element "Coil-type"; Acid Resistant & Long Life Time.
- Heating Mantles for 50ml~5,000ml Flask are with a Device* for Holding a Optional Stainless-steel Rod $\phi 12.7$ mm for Article Support.

~450°C

3-Year Warranty



(1) 450°C Basic-type

▶ Controller ; Option!

~650°C

3-Year Warranty



(2) 650°C High Temp.-type

▶ Controller ; Option!



(a)

(b)

(3) Stand Rod, "RD300", $\phi 12.7 \times L520$ mm

(4) Holder / Clamp / Clip for Temp. Probe

(1) 450°C Basic-type, 50ml~100Lit

* 50ml~5000ml ; Non-magnetic Aluminum Case, 10Lit~100Lit ; Powder-coated Steel Case

<L1>

Flasks, ml Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
50ml, $\phi 49 \sim 52$	90 W	30 mm	$\phi 150 \times h90$	DH.WHM12050	DH.WHM12060	
100ml, $\phi 63 \sim 66$	90 W	50 mm	$\phi 150 \times h90$	DH.WHM12051	DH.WHM12061	
250ml, $\phi 83 \sim 86$	160 W	50 mm	$\phi 150 \times h90$	DH.WHM12052	DH.WHM12062	
500ml, $\phi 102 \sim 105$	250 W	60 mm	$\phi 200 \times h105$	DH.WHM12053	DH.WHM12063	
1,000ml, $\phi 130 \sim 133$	350 W	85 mm	$\phi 200 \times h140$	DH.WHM12054	DH.WHM12064	
2,000ml, $\phi 165 \sim 169$	490 W	105 mm	$\phi 270 \times h135$	DH.WHM12055	DH.WHM12065	
3,000ml, $\phi 184 \sim 187$	600 W	110 mm	$\phi 270 \times h135$	DH.WHM12056	DH.WHM12066	
5,000ml, $\phi 221 \sim 224$	850 W	130 mm	$\phi 325 \times h165$	DH.WHM12057	DH.WHM12067	
10Lit, $\phi 278 \sim 282$	1,500 W	195 mm	$\phi 390 \times h240$	DH.WHM12073	DH.WHM12083	
20Lit, $\phi 344 \sim 348$	2,000 W	230 mm	$\phi 460 \times h280$	DH.WHM12074	DH.WHM12084	
30Lit, $\phi 410 \sim 420$	3,200 W	280 mm	$\phi 545 \times h330$	DH.WHM8085	DH.WHM8085.1	
50Lit, $\phi 489 \sim 494$	4,000 W	310 mm	$\phi 610 \times h390$	DH.WHM12075	DH.WHM12085	
100Lit, $\phi 609 \sim 614$	6,000 W	360 mm	$\phi 740 \times h450$	DH.WHM8084	DH.WHM8084.1	

(2) 650°C High Temp.-type, 100ml~10Lit

* 100ml~5000ml ; Non-magnetic Aluminum Case, 10Lit ; Powder-coated Steel Case

Flasks, ml Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
100ml, $\phi 63 \sim 66$	160 W	50 mm	$\phi 150 \times h90$	DH.WHM12091	DH.WHM12101	
250ml, $\phi 83 \sim 86$	250 W	50 mm	$\phi 150 \times h90$	DH.WHM12092	DH.WHM12102	
500ml, $\phi 102 \sim 105$	350 W	60 mm	$\phi 200 \times h105$	DH.WHM12093	DH.WHM12103	
1,000ml, $\phi 130 \sim 133$	490 W	85 mm	$\phi 200 \times h140$	DH.WHM12094	DH.WHM12104	
2,000ml, $\phi 165 \sim 169$	600 W	105 mm	$\phi 270 \times h135$	DH.WHM12095	DH.WHM12105	
3,000ml, $\phi 184 \sim 187$	750 W	105 mm	$\phi 270 \times h135$	DH.WHM12096	DH.WHM12106	
5,000ml, $\phi 221 \sim 224$	1,500 W	130 mm	$\phi 325 \times h165$	DH.WHM12097	DH.WHM12107	
10Lit, $\phi 278 \sim 282$	2,000 W	160 mm	$\phi 400 \times h260$	DH.WHM8100	DH.WHM8100.1	

Optional Remote Controller

- (a) Analog Remote Controller, 5~100% of Rated Voltage, Up to 1.2kW
(b) Digital Remote Controller, Ambient +5°C~700°C,
 $\pm 1^\circ\text{C}$ Uniformity at Set Temp.

Cat. No		Ea. Price, ₩
230V	120V	
DH.WHM13500	DH.WHM23500	
DH.WHM12500	DH.WHM22500	

Optional Temp. Probe & Stand Rod, Holder, Clamp for Digital Temp. Remote Controller

- DH.WHP503040 Direct Contact Temp. Probe, K-type, "SS300", L250 $\times \phi 4.8$ mm
DH.WHP503050 "SS400", L380 $\times \phi 4.8$ mm
DH.WHP503041 "SS310", Teflon Coated, L250 $\times \phi 5$ mm
DH.WHP503051 "SS410", Teflon Coated, L380 $\times \phi 5$ mm
DH.WHP503037 (3) Stand Rod, $\phi 12.7 \times L520$ mm, "RD300", for Heating Mantles
DH.WHP503012 (4) Holder / Clamp / Clip, "CL220", for Temp. Sensor

Ea. Price, ₩

DAIHAN™ Fabric-Housed Round Flask Heating Mantle, Remote Controlled, (1) Base & (2) Top Cover, 450℃, 50ml~100Lit
With Nickel Chrome Heating Element, Option-Controller, with Certi. & Traceability

라운드 플라스크용 직물케이스 히팅맨틀, Ni-Cr 열선 내장, 자력교반기와 사용가능, 조절기와 외장 온도 프로브 별도

- Supplied by DAIHAN Korea
- Ideal for All type Round Flask and Use with Magnetic Stirrer
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Nickel-Chrome Heating Element "Coil-type" ; Acid Resistant & Long Life Time.

~450℃

3-Year Warranty



(1) Self Standing-base

(1) Self Standing-base Mantle

<L1>

Flasks, ml Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside ϕ × h(mm)	230V	120V	
50ml, ϕ 49 ~ 52	90 W	30 mm	ϕ 80 × h55	DH.WHM12110	DH.WHM12120	
100/125ml, ϕ 63 ~ 66	90 W	40 mm	ϕ 100 × h65	DH.WHM12111	DH.WHM12121	
250ml, ϕ 83 ~ 86	160 W	45 mm	ϕ 125 × h75	DH.WHM12112	DH.WHM12122	
500ml, ϕ 102 ~ 105	250 W	60 mm	ϕ 155 × h85	DH.WHM12113	DH.WHM12123	
1,000ml, ϕ 130 ~ 133	350 W	70 mm	ϕ 191 × h100	DH.WHM12114	DH.WHM12124	
2,000ml, ϕ 165 ~ 169	490 W	80 mm	ϕ 226 × h120	DH.WHM12115	DH.WHM12125	
3,000ml, ϕ 184 ~ 187	600 W	90 mm	ϕ 255 × h130	DH.WHM12116	DH.WHM12126	
5,000ml, ϕ 221 ~ 224	850 W	100 mm	ϕ 293 × h155	DH.WHM12117	DH.WHM12127	

~450℃

3-Year Warranty



(2) Top Cover

(1) Self Standing base

<Ex.> (2) Top Cover
with (1) Self Standing-base
▶ Controller ; Option!

(2) Top Cover only

Flasks, ml Dia. ϕ mm	Mantle		Catalog No.		Ea. Price, ₩
	Watts (W)	Outside ϕ (mm)	230V	120V	
250ml, ϕ 83 ~ 86	80 W	ϕ 136	DH.WHM121301	DH.WHM121311	
500ml, ϕ 102 ~ 105	140 W	ϕ 155	DH.WHM121302	DH.WHM121312	
1,000ml, ϕ 130 ~ 133	240 W	ϕ 183	DH.WHM121303	DH.WHM121313	
2,000ml, ϕ 165 ~ 169	350 W	ϕ 220	DH.WHM121304	DH.WHM121314	
3,000ml, ϕ 184 ~ 187	450 W	ϕ 236	DH.WHM121305	DH.WHM121315	
5,000ml, ϕ 221 ~ 224	550 W	ϕ 275	DH.WHM121306	DH.WHM121316	
10Lit, ϕ 278 ~ 282	700 W	ϕ 332	DH.WHM121307	DH.WHM121317	
20Lit, ϕ 344 ~ 348	800 W	ϕ 400	DH.WHM121308	DH.WHM121318	
50Lit, ϕ 489 ~ 494	1,000 W	ϕ 550	DH.WHM8121	DH.WHM8121.1	
100Lit, ϕ 609 ~ 614	1,200 W	ϕ 665	DH.WHM8122	DH.WHM8122.1	

(a)



(b)



Optional Remote Controller

- (a) Analog Remote Controller, 5~100% of Rated Voltage, Up to 1.2kW
(b) Digital Remote Controller, Ambient +5℃~700℃,
±1℃ Uniformity at Set Temp.

Cat. No		Ea. Price, ₩
230V	120V	
DH.WHM13500	DH.WHM23500	
DH.WHM12500	DH.WHM22500	


DAIHAN™ Remote Controlled Reaction Round Flask Heating Mantle, Bottom Outlet-type, 450°C, 500mL~50Lit
With Nickel Chrome Heating Element, K-type Thermo-sensor Integrated, with Certi. & Traceability, Option-Controller
반응조용 히팅맨틀, 반응조 라운드 플라스크용, K-type 온도센서와 Ni-Cr 열선 내장, 조절기와 외장 온도 프로브 별도



- Supplied by DAIHAN Korea
- Ideal for All type Reaction Flasks
- 500mL~5000mL ; Non-magnetic Aluminum-Case, 10Lit~50Lit ; Powder-coated Steel Case
- Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Nickel-Chrome Heating Element "Coil-type" ; Acid Resistant & Long Life Time.
- * With a Device* for Holding a Optional Stainless-steel Rod ϕ 12.7 mm for Article Support.
- * Optional Remote Thermo Controller, K-type Thermo-Sensor

<L1>


Flasks, mL Dia. ϕ mm	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
500mL, ϕ 105	240 W	60 mm	ϕ 200 $\times$ h105	DH.WHM8220	DH.WHM8220.1	
1,000mL, ϕ 131	350 W	85 mm	ϕ 200 $\times$ h140	DH.WHM8221	DH.WHM8221.1	
2,000mL, ϕ 165 ~ 169	490 W	130 mm	ϕ 270 $\times$ h160	DH.WHM8222	DH.WHM8222.1	
3,000mL, ϕ 184 ~ 186	600 W	130 mm	ϕ 270 $\times$ h160	DH.WHM8223	DH.WHM8223.1	
5,000mL, ϕ 224 ~ 226	850 W	150 mm	ϕ 325 $\times$ h180	DH.WHM8225	DH.WHM8225.1	
10Lit, ϕ 278 ~ 282	1,400 W	195 mm	ϕ 400 $\times$ h240	DH.WHM8226	DH.WHM8226.1	
20Lit, ϕ 348 ~ 352	2,000 W	230 mm	ϕ 440 $\times$ h280	DH.WHM8227	DH.WHM8227.1	
30Lit, ϕ 410 ~ 420	3,200 W	280 mm	ϕ 545 $\times$ h330	DH.WHM8229	DH.WHM8229.1	
50Lit, ϕ 480 ~ 484	4,000 W	310 mm	ϕ 740 $\times$ h450	DH.WHM8228	DH.WHM8228.1	



Optional Remote Controller

- (a) Analog Remote Controller, 5~100% of Rated Voltage, Up to 1.2kW
(b) Digital Remote Controller, Ambient +5°C~700°C,
 $\pm 1^\circ\text{C}$ Uniformity at Set Temp.

Cat. No		Ea. Price, ₩
230V	120V	
DH.WHM13500	DH.WHM23500	
DH.WHM12500	DH.WHM22500	



Optional Temp. Probe & Stand Rod, Holder, Clamp for Digital Temp. Remote Controller

- DH.WHP503040 Direct Contact Temp. Probe, K-type, "SS300", L250 $\times$ ϕ 4.8 mm
DH.WHP503050 "SS400", L380 $\times$ ϕ 4.8 mm
DH.WHP503041 "SS310", Teflon Coated, L250 $\times$ ϕ 5 mm
DH.WHP503051 "SS410", Teflon Coated, L380 $\times$ ϕ 5 mm
DH.WHP503037 (1) Stand Rod, ϕ 12.7 $\times$ L520mm, "RD300", for Heating Mantles
DH.WHP503012 (2) Holder / Clamp / Clip, "CL220", for Temp. Sensor

Ea. Price, ₩

DAIHAN™ Aluminum-Case Cylindrical Reaction Vessel Heating Mantle, (1) Basic & (2) Stirring-type, 450℃, 100ml~20Lit
With Built-in Temp Controller; with/without Mag-stir Speed Control, with Certi. & Traceability
반응조용 히팅맨틀, 온도 조절기와 Ni-Cr 열선 내장, “기본형” 및 “자석교반형”

- Supplied by DAIHAN Korea
- Ideal for All type Reaction Vessel
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Built-in Temp-Controller. Nickel-Chrome Heating Element “Coil-type” ; Acid Resistant & Long Life Time.
- Heating Mantles for 100ml~6,000ml Flask are with a Device* for Holding an Optional Stainless-steel Rod ϕ 12.7mm for Article Support.

~450℃

3-Year Warranty



(1) Basic-type Reaction Vessel Mantle with Temp. Controller

(1) Basic-type Reaction Vessel Mantle, Built-in Temp. Controller

■ 100ml~6000ml ; Non-magnetic Aluminum Case, 10Lit~20Lit ; Powder-coated Steel Case

<L1>

Vessels, ml Dia. ϕ × h(cm)	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside ϕ × h(mm)	230V	120V	
100ml, ϕ 7 × h9	80 W	40 mm	ϕ 150 × h145	DH.WHM8130	DH.WHM8130.1	
250ml, ϕ 7 × h13	140 W	60 mm	ϕ 150 × h145	DH.WHM8131	DH.WHM8131.1	
500ml, ϕ 11 × h12	350 W	90 mm	ϕ 200 × h170	DH.WHM12146	DH.WHM12156	
1,000ml, ϕ 11 × h20	350 W	100 mm	ϕ 200 × h175	DH.WHM8133	DH.WHM8133.1	
2,000ml, ϕ 13 × h24	450 W	120 mm	ϕ 200 × h195	DH.WHM8134	DH.WHM8134.1	
2,000ml, ϕ 14 × h27	450 W	120 mm	ϕ 200 × h195	DH.WHM8135	DH.WHM8135.1	
3,000ml, ϕ 13 × h32	550 W	200 mm	ϕ 200 × h265	DH.WHM8136	DH.WHM8136.1	
4,000ml, ϕ 20 × h29	800 W	150 mm	ϕ 270 × h240	DH.WHM8138	DH.WHM8138.1	
5,000ml, ϕ 20 × h29	800 W	150 mm	ϕ 270 × h240			
6,000ml, ϕ 21.5 × h32	800 W	160 mm	ϕ 325 × h270	DH.WHM8139	DH.WHM8139.1	
10Lit, ϕ 21.5 × h45	1,400 W	220 mm	ϕ 325 × h310	DH.WHM8140	DH.WHM8140.1	
10Lit, ϕ 24 × h41	1,400 W	220 mm	ϕ 350 × h270	DH.WHM8141	DH.WHM8141.1	
20Lit, ϕ 25 × h47	2,000 W	350 mm	ϕ 350 × h400	DH.WHM8142	DH.WHM8142.1	
20Lit, ϕ 27 × h45	2,000 W	300 mm	ϕ 400 × h350	DH.WHM8143	DH.WHM8143.1	

~450℃

3-Year Warranty



(2) Stirring-type Reaction Vessel Mantle with Temp./Stirring-Controller and 1x Mag Stir bar

(2) Stirring-type Reaction Vessel Mantle, Built-in 50~1200rpm Mag-stir Speed & Temp. Controller

■ 100ml~6000ml ; Non-magnetic Aluminum Case, 10Lit~20Lit ; Powder-coated Steel Case

* With 1ea × PTFE Oval-type Magnetic Stirrer Bar, L25mm.

* 10-/20-Lit models Square-type Outside Shape.

Vessels, ml Dia. ϕ × h(cm)	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside ϕ × h(mm)	230V	120V	
100ml, ϕ 7 × h9	80 W	60 mm	ϕ 150 × h145	DH.WHM8150	DH.WHM8150.1	
250ml, ϕ 7 × h13	140 W	60 mm	ϕ 150 × h145	DH.WHM8151	DH.WHM8151.1	
500ml, ϕ 11 × h12	240 W	85 mm	ϕ 200 × h175	DH.WHM8152	DH.WHM8152.1	
1,000ml, ϕ 11 × h20	350 W	100 mm	ϕ 200 × h195	DH.WHM8153	DH.WHM8153.1	
2,000ml, ϕ 13 × h24	450 W	105 mm	ϕ 200 × h195	DH.WHM8154	DH.WHM8154.1	
2,000ml, ϕ 14 × h27	450 W	105 mm	ϕ 200 × h195	DH.WHM8155	DH.WHM8155.1	
3,000ml, ϕ 13 × h32	550 W	195 mm	ϕ 200 × h265	DH.WHM8156	DH.WHM8156.1	
4,000ml, ϕ 20 × h29	800 W	150 mm	ϕ 270 × h240	DH.WHM8158	DH.WHM8158.1	
5,000ml, ϕ 20 × h29	800 W	150 mm	ϕ 270 × h240			
6,000ml, ϕ 21.5 × h32	800 W	160 mm	ϕ 325 × h230	DH.WHM8159	DH.WHM8159.1	
10Lit, ϕ 21.5 × h45	1,400 W	230 mm	320 × 420 × h325	DH.WHM8160	DH.WHM8160.1	
10Lit, ϕ 24 × h41	1,400 W	220 mm	350 × 450 × h315	DH.WHM8161	DH.WHM8161.1	
20Lit, ϕ 25 × h47	2,000 W	350 mm	350 × 450 × h445	DH.WHM8162	DH.WHM8162.1	

Optional Stainless-steel Rod, Holder & Clamp

Ea. Price, ₩

DH.WHP503037 (3) Stand Rod, ϕ 12.7 × L520mm, “RD300”, for Heating Mantles

DH.WHP503012 (4) Holder / Clamp / Clip, “CL220”, for Temp. Sensor

(3) Stand Rod, “RD300”, ϕ 12.7 × L520mm

(4) Holder / Clamp / Clip for Temp. Probe

DAIHAN™ Remote Controlled Cylindrical Reaction Vessel Heating Mantle, (1) Plain & (2) Bottom Outlet-type, 450°C With Nickel Chrome Heating Element, K-type Thermo-Sensor Integrated, with Certi. & Traceability, Option-Controller, 100mL~100Lit 반응조용 히팅맨틀, K-type 온도센서와 Ni-Cr 열선 내장, 자력교반기와 사용가능, 조절기와 외장 온도 프로브 별도

- Supplied by DAIHAN Korea
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Optional Remote Thermo Controller, K-type Thermo-Sensor, and Nickel-Chrome Heating Element "Coil-type"; Acid Resistant & Long Life Time.
- Heating Mantles for 100mL~6,000mL Flask are with a Device* for Holding a Optional Stainless-steel Rod $\phi 12.7$ mm for Article Support.

~450°C

3-Year Warranty



(1) Plain Bottom-type
Remotecontrolled Reaction Vessel
Heating Mantle

▶ Controller ; Option!

(1) Plain Bottom-type

* 100mL~6000mL ; Non-magnetic Aluminum Case, 10Lit~100Lit ; Powder-coated Steel Case

<L1>

Vessels, mL Dia. $\phi \times h$ (cm)	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
100mL, $\phi 7 \times h9$	80 W	30 mm	$\phi 150 \times h100$	DH.WHM8170	DH.WHM8170.1	
250mL, $\phi 7 \times h13$	160 W	65 mm	$\phi 150 \times h100$	DH.WHM8171	DH.WHM8171.1	
500mL, $\phi 11 \times h12$	240 W	60 mm	$\phi 200 \times h105$	DH.WHM8172	DH.WHM8172.1	
1,000mL, $\phi 11 \times h20$	350 W	100 mm	$\phi 200 \times h140$	DH.WHM8173	DH.WHM8173.1	
2,000mL, $\phi 13 \times h24$	450 W	135 mm	$\phi 200 \times h160$	DH.WHM8174	DH.WHM8174.1	
2,000mL, $\phi 14 \times h27$	450 W	120 mm	$\phi 200 \times h160$	DH.WHM8175	DH.WHM8175.1	
3,000mL, $\phi 13 \times h32$	550 W	195 mm	$\phi 200 \times h230$	DH.WHM8176	DH.WHM8176.1	
4,000mL, $\phi 20 \times h29$	800 W	145 mm	$\phi 325 \times h220$	DH.WHM8178	DH.WHM8178.1	
5,000mL, $\phi 20 \times h29$	800 W	145 mm	$\phi 325 \times h220$	DH.WHM8178	DH.WHM8178.1	
6,000mL, $\phi 21.5 \times h32$	800 W	150 mm	$\phi 325 \times h220$	DH.WHM8179	DH.WHM8179.1	
10Lit, $\phi 21.5 \times h45$	1,400 W	230 mm	$\phi 325 \times h270$	DH.WHM8180	DH.WHM8180.1	
10Lit, $\phi 24 \times h41$	1,400 W	220 mm	$\phi 325 \times h270$	DH.WHM8181	DH.WHM8181.1	
20Lit, $\phi 25 \times h47$	2,000 W	350 mm	$\phi 350 \times h400$	DH.WHM8182	DH.WHM8182.1	
20Lit, $\phi 27 \times h45$	2,000 W	300 mm	$\phi 400 \times h350$	DH.WHM8183	DH.WHM8183.1	
50Lit, $\phi 31.5 \times h81$	4,000 W	550 mm	$\phi 440 \times h600$	DH.WHM8185	DH.WHM8185.1	
100Lit, $\phi 41.5 \times h91$	6,000 W	650 mm	$\phi 530 \times h700$	DH.WHM8186	DH.WHM8186.1	

~450°C

3-Year Warranty



(2) Bottom Outlet-type
Remotecontrolled Reaction Vessel
Heating Mantle

▶ Controller ; Option!

(2) Bottom Outlet-type

* 500mL~6000mL ; Non-magnetic Aluminum Case, 10Lit~100Lit ; Powder-coated Steel Case

* Bottom Outlet Hole Size ; L7.2 × W4.0cm(500~2,000mL) or L12.0 × W7.0cm(3,000mL~100Lit)

Vessels, mL Dia. $\phi \times h$ (cm)	Mantle			Catalog No.		Ea. Price, ₩
	Watts (W)	Inside Depth	Outside $\phi \times h$ (mm)	230V	120V	
500mL, $\phi 11 \times h12$	240 W	60 mm	$\phi 200 \times h105$	DH.WHM8190	DH.WHM8190.1	
1,000mL, $\phi 11 \times h20$	350 W	100 mm	$\phi 200 \times h140$	DH.WHM8191	DH.WHM8191.1	
2,000mL, $\phi 13 \times h24$	450 W	135 mm	$\phi 200 \times h160$	DH.WHM8192	DH.WHM8192.1	
2,000mL, $\phi 14 \times h27$	450 W	120 mm	$\phi 200 \times h160$	DH.WHM8193	DH.WHM8193.1	
3,000mL, $\phi 13 \times h32$	550 W	195 mm	$\phi 200 \times h230$	DH.WHM8194	DH.WHM8194.1	
4,000mL, $\phi 20 \times h29$	800 W	145 mm	$\phi 325 \times h220$	DH.WHM8196	DH.WHM8196.1	
5,000mL, $\phi 20 \times h29$	800 W	145 mm	$\phi 325 \times h220$	DH.WHM8196	DH.WHM8196.1	
6,000mL, $\phi 21.5 \times h32$	800 W	150 mm	$\phi 325 \times h220$	DH.WHM8197	DH.WHM8197.1	
10Lit, $\phi 21.5 \times h45$	1,400 W	230 mm	$\phi 325 \times h270$	DH.WHM8198	DH.WHM8198.1	
10Lit, $\phi 24 \times h41$	1,400 W	220 mm	$\phi 350 \times h270$	DH.WHM8199	DH.WHM8199.1	
20Lit, $\phi 25 \times h47$	2,000 W	350 mm	$\phi 350 \times h400$	DH.WHM8200	DH.WHM8200.1	
20Lit, $\phi 27 \times h45$	2,000 W	300 mm	$\phi 400 \times h350$	DH.WHM8201	DH.WHM8201.1	
50Lit, $\phi 31.5 \times h81$	4,000 W	550 mm	$\phi 440 \times h600$	DH.WHM8203	DH.WHM8203.1	
100Lit, $\phi 41.5 \times h91$	6,000 W	650 mm	$\phi 530 \times h700$	DH.WHM8204	DH.WHM8204.1	

Optional Remote Controller

- (a) Analog Remote Controller, 5~100% of Rated Voltage, Up to 1.2kW
 (b) Digital Remote Controller, Ambient +5°C~700°C,
 $\pm 1^\circ\text{C}$ Uniformity at Set Temp.

Cat. No		Ea. Price, ₩
230V	120V	
DH.WHM13500	DH.WHM23500	
DH.WHM12500	DH.WHM22500	

Optional Temp. Probe & Stand Rod, Holder, Clamp for Digital Temp. Remote Controller

Ea. Price, ₩

(3) Stand Rod, "RD300", $\phi 12.7 \times L520$ mm

(4) Holder / Clamp / Clip for Temp. Probe

- DH.WHP503040 Direct Contact Temp. Probe, K-type, "SS300", L250 × $\phi 4.8$ mm
 DH.WHP503050 "SS400", L380 × $\phi 4.8$ mm
 DH.WHP503041 "SS310", Teflon Coated, L250 × $\phi 5$ mm
 DH.WHP503051 "SS410", Teflon Coated, L380 × $\phi 5$ mm
 DH.WHP503037 (3) Stand Rod, $\phi 12.7 \times L520$ mm, "RD300", for Heating Mantles
 DH.WHP503012 (4) Holder / Clamp / Clip, "CL220", for Temp. Sensor

DAIHAN™ 3 & 6 Places Multiuse Heating Mantle, for 250~1000ml Flasks, 450℃, with Certi. & Traceability
With Glassware Supporting Rod & Hook, Ideal for Extraction · Reflux · Distillation Applications, Lower Profile
3 & 6구 다용도 히팅맨틀, 추출 · 환류 · 증류용 등에 적합, 250~1000ml 둥근 플라스크용

- Supplied by DAIHAN Korea
- Fit to DURAN, Glassco, PYREX and Other-brand Round Flasks
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Repetitive Extraction (Kjeldahl, Soxhlet and other Extraction Procedure), Reflux and Distillation Laboratory Applications

* Top Cover ; Stainless-steel 304 ; Body ; Powder Coated Steel

* Analog Controller Built-in System


* Rods and Hooks for Rods Connection

* Nickel Chrome Heating Element - Coil Type

* Optional Clamp & Clamp Holder for Glassware ; Order Separately

~450℃

3-Year Warranty



6 places Heating Mantle
(WHM12293)

3-Year Warranty



3 places Heating Mantle
(WHM12391)

<L1>

Working Places	Catalog No.		Round Flask 250~1000ml			Watts(W)	Inside Depth (mm)	Dimension w × d × h (mm)	Ea. Price, ₩
	230V	120V	Capa. (ml)	Dia. Max. Φ (mm)	Bottom Style				
6 places	DH.WHM12291	DH.WHM12292	250/300	83 ~ 86	Flat	840	50	800 × 240 × 110	
6 places	DH.WHM12293	DH.WHM12294	250/300	83 ~ 86	Round	840	50	800 × 240 × 110	
6 places	DH.WHM12295	DH.WHM12296	500	102 ~ 105	Round	1440	65	950 × 270 × 130	
3 places	DH.WHM12391	DH.WHM12392	500	102 ~ 105	Round	720	65	550 × 270 × 130	
3 places	DH.WHM12393	DH.WHM12394	1,000	130 ~ 133	Round	1050	75	690 × 340 × 170	

Spare Support Rod & Hook Set

- DH.WHM30500 Support Rod & Hook Set, "MHR01", for 3Places 500ml, 2 × 750-/1 × 500-mm Rod and 2 × Hook for Rod
 DH.WHM31000 , "MHR02", for 3Places 1000ml, 2 × 750-/1 × 600-mm Rod and 2 × Hook for Rod
 DH.WHM60250 , "MHR03", for 6Places 250-/300-ml, 2 × 750-/1 × 900-mm Rod and 2 × Hook for Rod
 DH.WHM60500 , "MHR04", for 6Places 500ml, 2 × 750-/1 × 940-mm Rod and 2 × Hook for Rod

Ea. Price, ₩



※ Optional Clamps & Clamp Holders
⇒ Visit DAIHAN-Sci.com Shopping Mall

DAIHAN™ Analog Compact Temp Controller "WHM-C10A", Up to 1.2kW, for All-Purpose Heating Instruments With Solid State, Proportional Voltage Power Control (5~100%), 아날로그 온도 조절기



3-Year Warranty

Analog Compact Temp. Controller "WHM-C10A"



<Ex.>
Heating Mantle Set with
the Analog Compact Temp.
Controller "WHM-C10A"

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Small Footprint, Compact Analog Temperature Controller
- Precise Manual Control of Heating Mantles, Tapes, Cords and Other Resistive Loads up to 1.2 kW
- * Safe and Durable Design
- * Solid state, Proportional Voltage Power Control (5~100% or Rated Voltage)

Specifications & Ordering information

<L1>

Spec.	Model	"WHM-C10A"
Control Range & Capacity		5~100% of Rated Voltage, Up to 1.2 kW
Dimension & Weight		95 × 98 × h106mm, 1kg
Packing Size & Gross Weight		190 × 160 × h130mm, 2kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz
Cat. No.	230V	DH.WHM13500
	120V	DH.WHM23500
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

DAIHAN™ Digital Feedback Temp Controller "WHM-C10D", Up to 700℃, ±1.0℃, with Back-Light LCD For All-purpose Heating Instruments, Available K-type External Direct Contact Thermocouple(Optional), 디지털 온도 조절기



3-Year Warranty

Mode Button Back Light LCD Turn+Push Jog-Shuttle

Digital Feedback Temp. Controller "WHM-C10D"



<Ex.>
Heating Mantle Set
with Temp. Probe and
the Digital Feedback Temp.
Controller "WHM-C10D"

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- User's Self-compensation Function to Control Difference between Real Temp. and Indicated Temp. (±50℃)
- Patented Jog-Shuttle Control System
- * Digital Feedback Control by High-Performance Microprocessor
- * Digital LCD with Back-Light Function * Small Footprint and User-Friendly Design.
- * Storage function, Timer, Alarm & Locking Mode
- * Heating Instruments can be controlled up to 2kW.
- * Accurate Temp. Control of the Liquid in the Vessel using Optional Direct K-type Sensor

Specifications & Ordering information

<L1>

Spec.	Model	"WHM-C10D"
Temp. Range & Accuracy		Ambient +5℃~700℃, ±1℃ at Set Temperature
Temp. Resolution		1℃-Display, 1℃-Control
Timer & Alarm		99hr 59min (Continuous run), Error status & Timer-end
Display		Digital LCD with Back Light Function (Temp., Timer)
Controller		Digital Feedback Controller with Jog Shuttle Switch (Turn+Push)
Control Capacity		up to 2kW for 230V, 1.5kW for 120V
Dimension & Weight		145 × 185 × h68mm, 2.8kg
Packing Size & Gross Weight		280 × 250 × h110mm, 4kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz
Cat. No.	230V	DH.WHM12500
	120V	DH.WHM22500
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

Optional Temp. Probe & Stand Rod, Holder, Clamp

Ea. Price, ₩

- DH.WHP503040 (1) Direct Contact Temp. Probe, K-type, "SS300", L250 × φ 4.8 mm
- DH.WHP503050 (1) Direct Contact Temp. Probe, K-type, "SS400", L380 × φ 4.8 mm
- DH.WHP503041 (2) Direct Contact Temp. Probe, K-type, "SS310", Teflon Coated, L250 × φ 5 mm
- DH.WHP503051 (2) Direct Contact Temp. Probe, K-type, "SS410", Teflon Coated, L380 × φ 5 mm
- DH.WHP503037 (3) Stand Rod, φ 12.7 × L520mm, "RD300", for Heating Mantles
- DH.WHP503012 (4) Holder / Clamp / Clip, "CL220", for Temp. Sensor

DAIHAN™ Heating Tape, Standard-type, 250℃, 13 & 25mm

<L1>

Heating Tape
250℃



3-Year Warranty

Heating Tapes

* For Glass Tubes Metal Pipet

Catalog No.	13mm Width × Length(m)	Watts & Volt	Ea. Price, ₩
DH.WHM12311	13mm × 0.6m	200W, 230V	
DH.WHM12312	× 1.0m	200W, 230V	
DH.WHM12313	× 1.5m	250W, 230V	
DH.WHM12314	× 2.0m	300W, 230V	
DH.WHM12315	× 2.5m	350W, 230V	
DH.WHM12316	× 3.0m	400W, 230V	
Catalog No.	25mm Width × Length(m)	Watts & Volt	Ea. Price, ₩
DH.WHM12321	25mm × 0.6m	200W, 230V	
DH.WHM12322	× 1.0m	200W, 230V	
DH.WHM12323	× 1.5m	250W, 230V	
DH.WHM12324	× 2.0m	300W, 230V	
DH.WHM12325	× 2.5m	350W, 230V	
DH.WHM12326	× 3.0m	400W, 230V	
DH.WHM12327	× 4.0m	550W, 230V	
DH.WHM12328	× 6.0m	700W, 230V	

► Controller ; Option!



Optional Remote Controller

(a) Analog Remote Controller, 5~100% of Rated Voltage, Up to 1.2kW
(b) Digital Remote Controller, Ambient +5℃~700℃,
± 1℃ Uniformity at Set Temp.

Cat. No		Ea. Price, ₩
230V	120V	
DH.WHM13500	DH.WHM23500	
DH.WHM12500	DH.WHM22500	

DAIHAN™ Custom-Made Heating Mantle & Insulation, Up to 450℃~650℃

Ideal for Bottom Valved Vessels or General-use, 고객 주문 사양 히팅맨틀 및 단열/보온재

We can Offer DAIHAN Custom Heating Mantles, for which you can Specify the Outside Diameter of the Vessel, Length of Mantle, Voltage, Number of Circuits and Length of Each Heating Zone, Ambient Temperature, Unusual Ambient Conditions, Size & Location of Holes for View Ports, and Tubulations.

Please Submit a Sketch for Complicated Applications and Contact our International Sales Team or our Distributor in your Country.


- Provides Uniform Heating and Insulation for Scrubbing and Fractionating Columns.
- Closes with a Zipper, with Eyelets and laces or Hooks and laces.
- Internal Operating Temperature is up to 450℃ or 650℃

Examples for Customized Heating Mantle



Welcome! Custom-Made Heating Mantle!!

DAIHAN™ Homogenizer with Direct Controller "HG-15A", Up to 27,000 rpm, 1~2,500ml
With Analog Stepless Control System, High Efficiency Dispersing Tools, High-Speed · High-Quality · Low Noise
균질기/호모지나이저, 아날로그 컨트롤러 내장, 초고속/저소음 모터



A. Homogenizer Package-set "HG-15A-SET-A"

1~50ml



(1) Dispersing Tool "HT1008"

1~100ml



(2) Dispersing Tool "HT1010"

10~1,000ml



(3) Dispersing Tool "HT1018"

50~2,500ml



(4) Dispersing Tool "HT1025"

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Useful for All Dispersing and Homogenizing Applications
- Over Heat Protection

- * Designed for Liquid and Tissue Samples
- * Useful for All Dispersing and Homogenizing Applications
- * High Speed / High Quality / Low Noise ; Max. 27,000 RPM (without load)
- * Speed & Power Control Optimized for a Sample
- * High Efficiency Dispersing Tool made of Stainless-steel and PTFE

Specifications & Ordering information

<L1>

Spec.	Model	A. "HG-15A-SET-A" with HT1018	B. "HG-15A-SET-B" with HT1025	C. "HG-15A" (Unit Only)
Speed Range & Adjustment		Max. 27,000 rpm (without load), Stepless		
Working Range		1 ~ 2,500ml		
Motor Rating (Input / Output)		Input ; 300W, Output ; 160W		
Controller		Analog Phase Controller		
Standard Component Acc.		"HT1018" Dispersing Tool, "CL200" Clamp, "ST110" Stand	"HT1025" Dispersing Tool, "CL200" Clamp, "ST110" Stand	Unit Only
Safety Device		Over Heat Protector		
Usable Dispersing Tool		φ 8mm Stator "HT1008", φ 10mm Stator "HT1010", φ 18mm Stator "HT1018", φ 25mm Stator "HT1025"		
Material		Aluminum Die-casting and ABS-Unit, Stainless-steel (# 304)-Tool		
Permissible Temp & Moisture		Ambient 5℃~50℃, 85%		
Extension Arm Size		φ 12.7 × L250mm		
Dimension & Weight		60 × 75 × h215mm, 2.5kg		
Packing Size & Gross Weight		320 × 300 × h190mm, 5kg		
Power Consumption		160W		
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz		
Cat. No.	230V	DH.WHG02016	DH.WHG02017	DH.WHG02015
	120V	DH.WHG01016	DH.WHG01017	DH.WHG01015
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

Spare Dispersing Tool & Clamp & Plate Stand


	Ea. Price, ₩
DH.WHG02108 (1) Dispersing Tool, "HT1008", Rotor φ 6, Stator φ 8mm, 1~50ml	
DH.WHG02110 (2) Dispersing Tool, "HT1010", Rotor φ 7.96, Stator φ 10mm, 1~100ml	
DH.WHG02118 (3) Dispersing Tool, "HT1018", Rotor φ 12.8, Stator φ 18mm, 10~1000ml	
DH.WHG02125 (4) Dispersing Tool, "HT1025", Rotor φ 20, Stator φ 25mm, 50~2500ml	
SL.Cla8011 Clamp, "CL200", for Stirrer or Homogenizer, Cast Aluminum, Grip Capacity φ 23mm for Stand/Rod	
DH.WOS501510 Standard Plate Stand, "ST110", for Stirrer or Homogenizer, φ 16 × h800 mm Rod, w20 × d35 cm Plate, with Silicone Pad	
DH.WOS501515 Wide Plate Stand, "ST120", for Stirrer or Homogenizer, φ 23 × h1000 mm Rod, w30 × d40 cm Plate, with Silicone Pad	

※ See page 173 for the Specifications of Dispersing Tools and Accessories.

DAIHAN™ Homogenizer with Remote Digital Controller "HG-15D", 2,000~27,000rpm

With Microprocessor Control, High-Efficiency Dispersing Tools, Back Light LCD, 1~2,500ml

디지털 균질기/호모지나이저, 분리형 리모트 컨트롤러, 초고속/저소음 모터



A. Homogenizer with Remote Digital Controller "HG-15D-SET-A"

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

■ Patented Jog-Shuttle Control System

- * Useful for All Dispersing and Homogenizing Applications
- * Designed for Liquid and Tissue Samples
- * High Speed / High Quality / Low Noise ; Max. 27,000 RPM (without load)
- * Remote Digital controller ; Fully Digital controlled by using Microprocessor
- * Speed & Power Control optimized for a Sample
- * LCD of Speed(RPM) and Power Bar Graph
- * High Efficiency Dispersing Tool made of Stainless-steel and PTFE

Specifications & Ordering information

<L1>

Spec.	Model	A. "HG-15D-SET-A" with HT1018	B. "HG-15D-SET-B" with HT1025	C. "HG-15D" (Unit Only)
Speed Range & Resolution		2,000 ~ 27,000rpm (without load), 10 rpm-Control		
Working Range		1 ~ 2,500ml		
Motor Rating (Input / Output)		Input ; 300W, Output ; 160W		
Timer & Alarm		10sec ~ 15min, Timer-end Alarm		
Display		Digital LCD with Back-Light Function		
Controller		Digital Feedback Controller (Remote Type)		
Standard Component Acc.		"HT1018" Dispersing Tool, "CL200" Clamp, "ST110" Stand	"HT1025" Dispersing Tool, "CL200" Clamp, "ST110" Stand	Unit Only
Safety Device		Over Load & Over Heat Protector		
Usable Dispersing Tool		φ 8mm Stator "HT1008", φ 10mm Stator "HT1010", φ 18mm Stator "HT1018", φ 25mm Stator "HT1025"		
Material		Aluminum Die-casting and ABS-Unit, Stainless-steel (# 304)-Tool		
Permissible Temp & Moisture		Ambient 5°C~50°C, 85%		
Extension Arm Size		φ 12.7 × L250mm		
Dimension & Weight		60 × 75 × h215mm, 2.5kg		
Packing Size & Gross Weight		580 × 340 × h230mm, 6kg		
Power Consumption		160W		
Electrical Rqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz		

Cat. No.	230V	DH.WHG02019	DH.WHG02020	DH.WHG02018
	120V	DH.WHG01019	DH.WHG01020	DH.WHG01018
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

Spare Dispersing Tool & Clamp & Plate Stand

		Ea. Price, ₩
DH.WHG02108	(1) Dispersing Tool, "HT1008", Rotor φ 6, Stator φ 8mm, 1~50ml	
DH.WHG02110	(2) Dispersing Tool, "HT1010", Rotor φ 7.96, Stator φ 10mm, 1~100ml	
DH.WHG02118	(3) Dispersing Tool, "HT1018", Rotor φ 12.8, Stator φ 18mm, 10~1000ml	
DH.WHG02125	(4) Dispersing Tool, "HT1025", Rotor φ 20, Stator φ 25mm, 50~2500ml	
SL.Cla8011	Clamp, "CL200", for Stirrer or Homogenizer, Cast Aluminum, Grip Capacity φ 23mm for Stand/Rod	
DH.WOS501510	Standard Plate Stand, "ST110", for Stirrer or Homogenizer, φ 16 × h800 mm Rod, w20 × d35 cm Plate, with Silicone Pad	
DH.WOS501515	Wide Plate Stand, "ST120", for Stirrer or Homogenizer, φ 23 × h1000 mm Rod, w30 × d40 cm Plate, with Silicone Pad	

※ See page 173 for the Specifications of Dispersing Tools and Accessories.

1~50ml



(1) Dispersing Tool "HT1008"

1~100ml



(2) Dispersing Tool "HT1010"

10~1,000ml




(3) Dispersing Tool "HT1018"

50~2,500ml



(4) Dispersing Tool "HT1025"

Articles	Cat. No	Description																																																																																								
1~50ml  (1) Dispersing Tool "HT1008"		Homogenizer Accessory, Dispersing Tool & Plate Stand 1. Dispersing Tool ■ Accurate Stator & Rotor made of Stainless-steel ■ High Efficiency Dispersing Tool made of Stainless-steel and PTFE ■ Low Cost & High Efficiency <L1> <table><tr><th>Model</th><th>(1) "HT1008"</th><th>(2) "HT1010"</th><th>(3) "HT1018"</th><th>(4) "HT1025"</th></tr><tr><td>Working Range</td><td>1 ~ 50ml</td><td>1 ~ 100ml</td><td>10 ~ 1,000ml</td><td>50 ~ 2,500ml</td></tr><tr><td>Stator Diameter</td><td>φ 8 mm</td><td>φ 10 mm</td><td>φ 18 mm</td><td>φ 25 mm</td></tr><tr><td>Rotor Diameter</td><td>φ 6 mm</td><td>φ 7.96 mm</td><td>φ 12.8mm</td><td>φ 20 mm</td></tr><tr><td>Gap between Rotor & Stator</td><td>0.25 mm</td><td>0.2 mm</td><td>0.25 mm</td><td>0.35 mm</td></tr><tr><td>Allowable Max. Speed</td><td colspan="2">27,000 rpm</td><td>20,000 rpm</td><td>15,000 rpm</td></tr><tr><td>Immersion Depth Min.</td><td colspan="2">35 mm</td><td colspan="2">60 mm</td></tr><tr><td>Immersion Depth Max</td><td colspan="2">70 mm</td><td colspan="2">130 mm</td></tr><tr><td>Shaft Length</td><td colspan="2">138 mm</td><td colspan="2">225 mm</td></tr><tr><td>Bearing Material</td><td colspan="4">PTFE</td></tr><tr><td>Suitable for Solvent</td><td colspan="4">Yes</td></tr><tr><td>Max. Temperature</td><td colspan="4">180℃</td></tr><tr><td>Sterilization Method</td><td colspan="4">All Methods</td></tr><tr><td>Suspensions & Emulsions</td><td colspan="4">Suspension ; 10 ~ 50μm, Emulsion ; 1 ~ 10μm</td></tr></table> <table><tr><th>Cat. No.</th><th>DH.WHG02108</th><th>DH.WHG02110</th><th>DH.WHG02118</th><th>DH.WHG02125</th></tr><tr><td>Price, Each. ₩</td><td></td><td></td><td></td><td></td></tr></table> • Test Condition ; 1.8 liters of water in the 2 liters capacity beaker. 2. Clamp & Plate Stand <table><tr><th></th><th>Ea. Price, ₩</th></tr><tr><td>SL.Cla8011 (5) Clamp, "CL200", for Stirrer or Homogenizer, L135mm, Cast Aluminum, Grip Capacity φ23mm for Stand/Rod</td><td></td></tr><tr><td>DH.WOS501510 (6) Standard Plate Stand, "ST110", for Stirrer or Homogenizer, φ16×h800mm Rod, w20×d35cm Plate, with Silicone Pad</td><td></td></tr><tr><td>DH.WOS501515 (7) Wide Plate Stand, "ST120", for Stirrer or Homogenizer, φ23×h1,000mm Rod, w30×d40cm Plate, with Silicone Pad</td><td></td></tr></table>	Model	(1) "HT1008"	(2) "HT1010"	(3) "HT1018"	(4) "HT1025"	Working Range	1 ~ 50ml	1 ~ 100ml	10 ~ 1,000ml	50 ~ 2,500ml	Stator Diameter	φ 8 mm	φ 10 mm	φ 18 mm	φ 25 mm	Rotor Diameter	φ 6 mm	φ 7.96 mm	φ 12.8mm	φ 20 mm	Gap between Rotor & Stator	0.25 mm	0.2 mm	0.25 mm	0.35 mm	Allowable Max. Speed	27,000 rpm		20,000 rpm	15,000 rpm	Immersion Depth Min.	35 mm		60 mm		Immersion Depth Max	70 mm		130 mm		Shaft Length	138 mm		225 mm		Bearing Material	PTFE				Suitable for Solvent	Yes				Max. Temperature	180℃				Sterilization Method	All Methods				Suspensions & Emulsions	Suspension ; 10 ~ 50μm, Emulsion ; 1 ~ 10μm				Cat. No.	DH.WHG02108	DH.WHG02110	DH.WHG02118	DH.WHG02125	Price, Each. ₩						Ea. Price, ₩	SL.Cla8011 (5) Clamp, "CL200", for Stirrer or Homogenizer, L135mm, Cast Aluminum, Grip Capacity φ23mm for Stand/Rod		DH.WOS501510 (6) Standard Plate Stand, "ST110", for Stirrer or Homogenizer, φ16×h800mm Rod, w20×d35cm Plate, with Silicone Pad		DH.WOS501515 (7) Wide Plate Stand, "ST120", for Stirrer or Homogenizer, φ23×h1,000mm Rod, w30×d40cm Plate, with Silicone Pad	
Model	(1) "HT1008"		(2) "HT1010"	(3) "HT1018"	(4) "HT1025"																																																																																					
Working Range	1 ~ 50ml		1 ~ 100ml	10 ~ 1,000ml	50 ~ 2,500ml																																																																																					
Stator Diameter	φ 8 mm		φ 10 mm	φ 18 mm	φ 25 mm																																																																																					
Rotor Diameter	φ 6 mm		φ 7.96 mm	φ 12.8mm	φ 20 mm																																																																																					
Gap between Rotor & Stator	0.25 mm	0.2 mm	0.25 mm	0.35 mm																																																																																						
Allowable Max. Speed	27,000 rpm		20,000 rpm	15,000 rpm																																																																																						
Immersion Depth Min.	35 mm		60 mm																																																																																							
Immersion Depth Max	70 mm		130 mm																																																																																							
Shaft Length	138 mm		225 mm																																																																																							
Bearing Material	PTFE																																																																																									
Suitable for Solvent	Yes																																																																																									
Max. Temperature	180℃																																																																																									
Sterilization Method	All Methods																																																																																									
Suspensions & Emulsions	Suspension ; 10 ~ 50μm, Emulsion ; 1 ~ 10μm																																																																																									
Cat. No.	DH.WHG02108	DH.WHG02110	DH.WHG02118	DH.WHG02125																																																																																						
Price, Each. ₩																																																																																										
	Ea. Price, ₩																																																																																									
SL.Cla8011 (5) Clamp, "CL200", for Stirrer or Homogenizer, L135mm, Cast Aluminum, Grip Capacity φ23mm for Stand/Rod																																																																																										
DH.WOS501510 (6) Standard Plate Stand, "ST110", for Stirrer or Homogenizer, φ16×h800mm Rod, w20×d35cm Plate, with Silicone Pad																																																																																										
DH.WOS501515 (7) Wide Plate Stand, "ST120", for Stirrer or Homogenizer, φ23×h1,000mm Rod, w30×d40cm Plate, with Silicone Pad																																																																																										
1~100ml  (2) Dispersing Tool "HT1010"																																																																																										
10 ~1,000ml  (3) Dispersing Tool "HT1018"																																																																																										
50 ~2,500ml  (4) Dispersing Tool "HT1025"																																																																																										
(5) High Quality Clamp (CL200) 																																																																																										
(6) Standard Plate Stand (ST110) 																																																																																										
(7) Wide Plate Stand (ST120) 																																																																																										




Steel Retort Base Duo Stand, for 2 × Plain Rods of Up to φ 16mm, without Rods
With Nonskid Feet Equipped, 2 × Level Screw Nut, and Base Thick-8mm, <Germany-Made>

Retort Stand Base Only, Acceptable 2 Plain Rod (below φ 16mm)

US.2040 Stand, "US2040", r=150mm
US.2041, "US2041", r=185mm

Optional Stainless-steel Rod

DH.WOS502001 Stainless-steel Rod, "RD501", φ 16 × h800mm, Non-screw
DH.WOS502002, "RD502", φ 16 × h1,000mm, Non-screw
DH.WOS502003, "RD503", φ 16 × h1,200mm, Non-screw



Stainless-steel Rod
With or Without 1.5P Screw

DAIHAN™ High-Quality Stainless-steel (#304) Rod, for Plate Stand, φ 16 & φ 23 mm

(1) Plain-type Rod, without Screw

DH.WOS502001 Stainless-steel Rod, "RD501", φ 16 × h800mm, Non-screw
DH.WOS502002, "RD502", φ 16 × h1,000mm, Non-screw
DH.WOS502003, "RD503", φ 16 × h1,200mm, Non-screw
DH.WOS502031, "RD531", φ 23 × h800mm, Non-screw
DH.WOS502033, "RD533", φ 23 × h1,200mm, Non-screw

(2) 1.5P Screw-type Rod

DH.WOS502011 Stainless-steel Rod, "RD511", φ 16 × h800mm, 1.5P Screw
DH.WOS502012, "RD512", φ 16 × h1,000mm, 1.5P Screw
DH.WOS502013, "RD513", φ 16 × h1,200mm, 1.5P Screw
DH.WOS502021, "RD521", φ 23 × h800mm, 1.5P Screw
DH.WOS502022, "RD522", φ 23 × h1,000mm, 1.5P Screw
DH.WOS502023, "RD523", φ 23 × h1,200mm, 1.5P Screw

DAIHAN™ Ultrasonic Homogenizer Full Sets, “UH-300” • “UH-500” • “UH-1000”, Up to 1,000ml, 19.5-25kHz
With Standard Probe, Automatic Data Storage, Timer(0.1sec~999sec), 디지털 초음파 호모지나이저, 사운드 프루프 박스 & 프로브 포함, 저소음



Ultrasonic Homogenizer Full Sets





Fig. Example of Use

※ Beaker is not Included



Spare & Optional Titanium Alloy Probe

- Supplied by DAIHAN Korea
- 본체, 사운드 프루프 박스, 프로브 세트 구성되어 경제적이고 편리성이 높음
- 사운드 프루프 박스 포함 - 샘플 처리 시 발생하는 소음을 최소화하여 작업 효율성을 높임
- 4.3인치 TFT 정전식 터치스크린 제어를 사용하여 쉽게 조작 가능
- Component of Set ;
 - (1) Ultrasonic Generator,
 - (2) Sonic System with Standard Horn (φ6mm for UH-300 & UH-500, φ20mm for UH-1000),
 - (3) Sound Proof Box, (4) Lab Support Jack, (5) Temp. Probe, (6) Power Cable, (7) Fuse, (8) User Manual
- Multi-Functional & Multi-Purpose Instrument that Uses Strong Ultrasound to Produce Cavitation Effect in Liquid and Ultrasonically Treat Substances
- Applications ; Cell Pulverization, Emulsification, Dispersion, Homogenization, Extraction, Defoaming, Degassing, Acceleration of Chemical Reactions
- * Large LCD Screen & Touch Operation Panel
- * Power Adjustable
- * Realize the Data Display of Man-Machine Friendly Operation Interface at a Glance
- * Intelligent Control ; Central Computer is used to Control the Total Working Time Centrally
- * Strong Sound Insulation - Made of Special Materials, which Greatly Reduces the Noise and Effectively Insulates the Noise
- * Durable Probe made by Titanium Alloy Material
- * Streamlined Appearance Design
- * Multiple Models Optional Matching Different Horn with Different Hosts
- * Data Storage Function - Can store 20sets of Data, Experimental Data and Convenient Recording
- * Liftable Operating Table Set in the Sound Proof Box - Easy to Adjust the Height of the Samples and Take Out
- * Can Realize the Voice Alarm Function caused by Malfunction, Over-Temperature or Operation Error

Specifications & Ordering information

<L1>

Spec.	Model	“UH-300”	“UH-500”	“UH-1000”
Mixing Capacity		0.1-300ml	0.1-500ml	1-1000ml
Power (Continuously Adjustable)		25-250W	10-650W	25-1200W
Frequency Range		20-25kHz		19.5-20.5kHz
Probe Included		φ6mm		φ20mm
Optional Probe		φ2mm	φ2 & φ12mm	φ6 & φ15 & φ22mm
Single Ultrasound Time		0.1-999s		
Single Gap Time		0.1-999s		
Total Time (Ultrasound + Gap)		1-99h 59m 59s		
Temperature Control Range		0-100℃		
Alarm Function		Temp. Time, Overload, No Load, Over Temperature		
Voice Alarm & Prompt Function		Yes		
Display Mode		Touch Screen		
Data Storage		20 Groups		
Power Supply		220V / 60Hz		
Dimension & Weight		410×225×h290mm, 12kg		410×225×h290mm, 15kg
Packing Size & Gross Weight		470×305×h430mm, 16kg		470×305×h430mm, 19kg
Sound Box Packing Size & Gross Weight		350×350×h550mm, 10kg		
Cat. No.		DH.Hom1001	DH.Hom1002	DH.Hom1003
Price, Set. ₩				

Spare & Optional Titanium Alloy Probe

Cat. No.	Probe Size (mm)	Power Range Suitable	Capacity Range	Power Rate	Ea. Price, ₩
DH.Pro7001	φ2	Power Less than 1000W	0.2-5ml	1-35%	
DH.Pro7002	φ6		10-100ml	1-60%	
DH.Pro7003	φ12		50-500ml	20-95%	
DH.Pro7004	φ15	Power Higher than 1200W	100-600ml	1-80%	
DH.Pro7005	φ20		100-1000ml	10-90%	
DH.Pro7006	φ22		200-1000ml	20-95%	

DAIHAN™ Homogenizing Stirrer "HS-30E", for Tissue Grinder, 200~5,000rpm
Analog Phase Control, without Tissue Grinder, with Certi. & Traceability, 호모지나이저 스테러



3-Year Warranty

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for usage with Tissue Grinders

* Designed for most Laboratory Applications such as Homogenizing or Mixing
* Compact Unit with Solid State Speed Control

Specifications & Ordering information

<L1>

Spec.	Model	A. "HS-30E" Unit Only	B. "HS-30E" Package Set
Speed Range & Max. Torque		200 ~ 5,000 rpm, 85 Ncm	
Stirring Capacity		Max. 10 Lit	
Motor Rating		Input ; 83W, Output ; 75W, 1/10 HP	
Safety Device		Over Heat Protector	
Controller		Analog Phase Controller	
Consists of		Unit Only	Unit, "CL200" Clamp, "ST110" Stand Set (Rod & Plate)
Dimension & Weight		90 × 145 × h180mm, 3.5kg	
Packing Size & Gross Weight		320 × 300 × h190mm, 5kg	
Power Consumption		75W	
Electric Reqs.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	
Cat. No.	230V	DH.WOS01010	DH.WOS01011
	120V	DH.WOS02010	DH.WOS02011
Price, Each. ₩			


* Other Specifications are available upon Customer's Request.

※ **Note! The Package does not Include a Tissue Grinder**



B. Homogenizing Stirrer Package-Set "HS-30E"

Cowie® Tissue Grinder, with PTFE Pestle Head · Borosilicate Glass Mortar · Stainless-steel Pestle Rod(23~27cm), 2~50ml
Good Chemical/Heat Resistance, Standard- & Serrated-type PTFE Head Pestles, <UK-Made> 티슈 그라인더



■ Supplied by Cowie® in U.K

* Use for the Controlled Reduction of Particle Size and Homogenisation for a Variety of Substances, especially Biological Material.

* Pestle Head ; Pure PTFE, Pestle Rod ; Stainless-steel, ϕ 6.5mm, Length 23~27cm

* STS Pestle Rod is Perfect for Motor Drive * The Mortar is precision Borosilicate Glass.

(1) Glass Mortar Only, Borosilicate, 2~50ml

CW.011.102	Mortar, 2ml, ϕ 8 × L120mm
CW.011.105	, 5ml, ϕ 12 × L135mm
CW.011.110	, 10ml, ϕ 15 × L150mm
CW.011.115	, 15ml, ϕ 19 × L155mm
CW.011.130	, 30ml, ϕ 25 × L175mm
CW.011.150	, 50ml, ϕ 32 × L195mm

(2) Standard-type Pestles Only

CW.011.202	Pestle, S/T, for 2ml
CW.011.205	, for 5ml
CW.011.210	, for 10ml
CW.011.215	, for 15ml
CW.011.230	, for 30ml
CW.011.250	, for 50ml

Ea. Price, ₩

(3) Serrated-type Pestles Only

CW.011.302	Pestle, Ser., for 2ml
CW.011.305	, for 5ml
CW.011.310	, for 10ml
CW.011.315	, for 15ml
CW.011.330	, for 30ml
CW.011.350	, for 50ml


Ea. Price, ₩

<L1>

Ea. Price, ₩

DAIHAN Wisd™ 550℃ High-temp Hi-Light Hotplate “HP-D165” & “HP-D200”, Schott NEXTREMA® Ceramic Glass Plate
With $\phi 165$ & $\phi 200$ mm Heating Zone, Digital Control, 2× Heating Control Mode, 다용도 하이라이트 고온가열판


CE



3-Year Warranty

I. “HP-D165” $\phi 165$ mm Heating Zone

※Holding Stand Acceptable!



3-Year Warranty

II. “HP-D200” $\phi 200$ mm Heating Zone

※ The Figure; Installation of
Optional Rod, Holder, and Temp Probe

- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance
- Hi-Light Heating Element & Schott NEXTREMA® Ceramic Glass Plate.
- 2× Heating Control Mode ; Max 550℃ Heating Control
 - Standard Mode ; 10-steps Heating Mode without Optional External Temp. Probe.
 - Precise Control Mode ; 5℃ Control Resolution with Optional External Temp. Probe.
- Patented Jog Dial Control System
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * High Temperature Indication LED Lamp for Prevention of Burn.
- * Can be used with Optional External Temperature Probe for Precise Control Mode

Specifications & Ordering Information

<L1>


Spec.	Model	I. “HP-D165” $\phi 165$ mm Heating Zone	II. “HP-D200” $\phi 200$ mm Heating Zone
Capacity		Max. 10Lit.	
Temp. Range & Accuracy		100℃~550℃ & $\pm 1.5^\circ\text{C}$ at 330℃	
Plate Size & Heating Zone		224×289mm, $\phi 165$ mm	259×324mm, $\phi 200$ mm
Heater & Heating Power	Hi-Light Heating Element		
		1,100W	1,800W
Timer & Alarm		99hr 59min (with Continuous Run), Error	
Controller & Display		Digital Controller & LED Display	
2× Heating Control Mode		* Standard Mode ; 10-Steps Heating Control Mode without Optional External Temp. Probe	
		* Precise Control Mode ; 100℃~500℃ Temperature Range Control (5℃ Control Resolution) with Optional External Temp. Probe	
Material		* Plate ; Schott NEXTREMA® Ceramic Glass Plate. * Housing ; PBT & PC	
Dimension(w×d×h) & Weight		232×311×h64mm, 2.13kg	267×346×h64mm, 2.56kg
Electrical Reqts.		1Phase, AC230V, 50/60Hz	
Cat. No.		DH.HtpD0165	DH.HtpD0200
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Optional Accessories

- DH.WHP503060 (1) External / Direct Contact Temp. Probe, “SS500”, Stainless-steel, $\phi 4 \times L250$ mm, Cable L190mm
- DH.WHP503061 (2) External / Direct Contact Temp. Probe, “SS510”, Teflon Coated, $\phi 4 \times L250$ mm, Cable L190mm,
- DH.WHP503036 (3) Stand Rod “RD200”, Stainless-steel, $\phi 12 \times L450$ mm, 1.5p screw
- DH.WHP503012 (4) Holder / Clamp / Clip, “CL220”, for Temp. Probe


Ea. Price, ₩



(1) Stainless-steel Temp. Probe



(2) Teflon Coated Temp. Probe



(3) Stand Rod “RD200”, Stainless-steel


(4) Holder / Clamp / Clip for Temp. Probe

DAIHAN™ Premium 500°C High-Temp Hotplate “HP500”, Solid Ceramic Plate, 200×200mm

With Large LCD, Optimum Insulation Layer, Accurate Temp Control, Touch-Button Controller, Hot-Top Indicator, Max. 500°C, Accu. $\pm 0.3^\circ\text{C}$

고온용 디지털 가열판, PID 온도 제어, 최적의 방열구조, 내열성 / 내화학성 / 내산성 솔리드 세라믹 글라스 플레이트


CE




~500°C



B. “HP500” Unit only



Controller of “HP500”



- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- High-quality Solid Ceramic Plate ; Heat-/Chemical-/Acid-Resistance and Easy to Clean
- Optimum Insulation Fourfold Layers ; BTU + Ceramic Paper + 2 STS Insulation Layer
- Maximum Temperature of 500°C
- Ergonomic Design
- Microprocessor PID control
- Patented Controller “Jog-Dial” and Main Touch Button - Best Durability, Best Mobility, Best Convenient to Operate with Back-Light Function
- 2 Internal K-type Temp. Sensor Operating Independently ; 1 Temp. Control Sensor for Temperature Accuracy and 1 Over Temp. Detection Sensor for User Safety
- Hot Top Indicator LED Lamp ; Hot Surface(Over 50°C) Warning Lamp for User Safety
- * Excellent for Sample Drying and most General Lab Heating Applications
- * User's Self -Compensation Function ; to Control Difference between Real Temp. and Indicated Temp. ($\pm 20.0^\circ\text{C}$)
- * Accurate Temp. Control of the Sample in the Vessel using PT100 Sensor (Optional)
- * Storage Function for the set Values ; Temp. and Timer
- * Minimization of Noise and Improved the Product Life
- * Use of the Special Advanced Insulation for the safety of the Body Temp.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled)

Specifications & Ordering information

<L2>

Spec.		Model	A. “HP500S” (Package-set)	B. “HP500” (Unit Only)
Material	Body		Powder Coated Aluminum	
	Plate		Solid Ceramic Plate, 200×200mm	
Stirring Capa.			Max. 20 Lit	
Temp. Range & Accuracy			Max. 500°C, $\pm 0.3^\circ\text{C}$ at Set Temperature	
Temp. Resolution			0.1 °C-Display, 0.5 °C-Control	
Temp. Uniformity			Specially Designed Heating Module ; Temp. Difference 10% less	
Heating Power / Consumption			850W	
Timer & Alarm			99hr 59min (with Continuous run.), Error status & Timer-end	
Display			Digital LCD with Back-Light Function	
Controller			Digital Feedback Controller with Jog Dial (Turn+Push) and Touch button.	
Safety Circuit			Over Temp. Protection Hot Top Indicator LED Lamp	
Included Standard Acc.			① Unit ② Temp. Probe ③ Rod ④ Holder+Clamp	① Unit Only
Permissible Temp. & RH %			Ambient +5 ~ 50°C, 85% RH	
Dimension & Weight			200×332×h114mm, 3.8kg	
Electrical Reqts.			1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	
Cat. No.	230V		DH.WHP01500	DH.WHP00500
	120V		DH.WHP01501	DH.WHP00501
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

Component of Package Set

- * Package-set ; ① Hotplate Unit ② “SS500” PT100 External / Direct Contact Temp. Probe,
③ “RD200” Stainless-steel Support Rod $\phi 12 \times L450\text{mm}$
④ “CL220” Holder+Clamp+Clip for External Temp. Probe

Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

- DH.WHP503060 (1) External / Direct Contact Temp. Probe, “SS500”, Stainless-steel, $\phi 4 \times L250\text{mm}$, Cable L190mm
DH.WHP503061 (2) External / Direct Contact Temp. Probe, “SS510”, Teflon Coated, $\phi 4 \times L250\text{mm}$, Cable L190mm
DH.WHP503036 (3) Stand Rod, “RD200”, Stainless-steel, $\phi 12 \times L450\text{mm}$, 1.5p screw, for “HP500”
DH.WHP503012 (4) Holder · Clamp · Clip, “CL220”, for Temp. Probe

Ea. Price, ₩

DAIHAN™ 380℃ Analog & Digital Hotplate “HP-A” & “HP-D”, Ceramic-Coated Plate

180×180- or 260×260-mm Plate, with High Temp Uniformity, CE Certified, with Certi. & Traceability

아날로그 & 디지털 정밀 가열판, 우수한 온도균일성 & 경제성, 세라믹 코팅 플레이트

A. 180×180mm Plate,

I. Analog Unit “HP-20A”

II. Digital Unit & III. Package Set

“HP-20D”

“HP-20D” Set



- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

- * Excellent for Sample Drying and most General Lab Heating Applications
- * Superior Temperature Uniformity ; Less than 10% on the Plate at Ambient Temperature

I. Analog Control Hotplate “HP-A”

- * Superior Linearity for Temperature Control by Specially Designed Control Circuitry (PWM control)


II. Digital Control Hotplate “HP-D”

■ Patented Jog-Shuttle Control System

- * Digital Feedback Control by High-Performance Microprocessor
- * Digital LCD Display with Back-Light Function, 0.1℃-Display, 0.5℃-Control
- * Temperature Accuracy ; $\pm 0.3^{\circ}\text{C}$ at Set Temp
- * Safety Circuit ; Over Temp Protection
- * Timer ; 99hr 59min (with Continuous run), Error status & Timer-end Alarm
- * Locking Mode for Experimental Safety Supported (Input to Jog Shuttle can be disabled)
- * User's Self-Compensation ; to Control Difference between Real Temp & Indicated Temp ($\pm 20.0^{\circ}\text{C}$)

III. Digital Control Hotplate Package-set

- * III-A. “HP-20D” Set ; Complete with ① Hotplate “HP-20D”, ② “SS500” PT100 External / Direct Contact Temp. Probe, ③ “RD200” Stainless-steel Support Rod $\phi 12 \times L450\text{mm}$, ④ “CL220” Holder+Clamp+Clip for External Temp. Probe
- * III-B. “HP-30D” Set ; Complete with ① Hotplate “HP-30D”, ② “SS500” PT100 External / Direct Contact Temp. Probe, ③ “RD100” Stainless-steel Support Rod $\phi 12.7 \times L390\text{mm}$, ④ “CL220” Holder+Clamp+Clip for External Temp. Probe



III-A. Digital Package Set
“HP-20D” Set

B. 260×260mm Plate,

I. Analog Unit “HP-30A”

II. Digital Unit & III. Package Set

“HP-30D”

“HP-30D” Set



(II-B) Digital “HP-30D”, 260×260mm Plate

Specifications & Ordering information

<L1>


Spec.		Model	A. 180×180mm Plate (I-A) Analog “HP-20A” (II-A) Digital “HP-20D”	B. 260×260mm Plate (I-B) Analog “HP-30A” (II-B) Digital “HP-30D”
Plate Size(mm)			180×180mm Plate-type	260×260mm Plate-type
Material	Body		Powder Coated Aluminum	Powder Coated steel
	Plate		Ceramic Coated Aluminum (Chemical / Acid Resistance, Easy to Clean)	
Temp. Range & Uniformity			Max. 380℃ Temp., High Uniformity ; Difference 10% less	
Heating Power/Consumption			600 W	1.2 kW
Permissible Temp. & RH %			R/T +5℃~50℃, 85% RH	
Dimension & Weight			206×307×h99mm, 2.8kg	286×387×h99mm, 3.3kg
Packing Size & Gross Weight			400×280×h200mm, 5kg	430×320×200hmm, 5kg
Electrical Reqts.			1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	
I. Analog Unit Only			(I-A) “HP-20A”	(I-B) “HP-30A”
Cat. No.	230V		DH.WHP03120	DH.WHP03121
	120V		DH.WHP02120	DH.WHP02121
Price, Each. ₩				
II. Digital Unit Only			(II-A) “HP-20D”	(II-B) “HP-30D”
Cat. No.	230V		DH.WHP03020	DH.WHP03022
	120V		DH.WHP02020	DH.WHP02022
Price, Each. ₩				
III. Digital Package-set			(III-A) “HP-20D” Set	(III-B) “HP-30D” Set
Cat. No.	230V		DH.WHP03021	DH.WHP03023
	120V		DH.WHP02021	DH.WHP02023
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.


Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

- DH.WHP503060 (1) Temp. Probe, “SS500”, “PT100”, $\phi 4 \times L250\text{mm}$
- DH.WHP503061 (2) Teflon Coated Temp. Probe, “SS510”, $\phi 4 \times L250\text{mm}$
- DH.WHP503011 Stand Rod, “RD100”, $\phi 12.7 \times L390$ for “HP-30A / 30D”
- DH.WHP503036 “RD200”, $\phi 12 \times L450$ for “HP-20A / 20D”
- DH.WHP503012 Holder · Clamp · Clip, “CL220”, for Temp. Probe


Ea. Price, ₩



DAIHAN™ 350°C Large Plate Digital Hotplate “HPLP”, 310×620mm Ceramic-Coated Plate, Superior Temp Uniformity
Built-in Digital PID Controller or Digital Remote Controller, Back Light LCD, Certificate & Traceability
디지털 대형 정밀 가열판, 디지털 피드백 온도 제어, 디지털 PID 조절기 / 디지털 원격 조절기, 백색 세라믹 코팅판



A. Large Digital Hotplate, with Built-in Digital PID Controller, Ceramic-Coated Plate-type “HPLP620”



B. Large Digital Hotplate, Digital Remote Controller, Ceramic-Coated Plate-type “HPLP620R”

- Supplied by SciLab Korea
- 3-Years Warranty
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Superior Temperature Uniformity ; Less than 10% on the Plate at Ambient Temperature
- Digital PID Controller
- Patented Jog-Shuttle Control System and Back-Light LCD

- * Excellent for Sample Drying and most General lab Heating Applications
- * Large Heating Area ; Ideal for TLC Plate Heating and Surface Heating in Electronic Industries
- * Ceramic Coated Plate -Chemical / Acid Prevention


Specifications & Ordering information

<L1>

Spec.	Models	A. “HPLP620” Built-in Controller	B. “HPLP620R” Remote Controller
Controller		Built-in Digital PID Controller	Digital Remote Controller
Plate Surface		Ceramic-coated	
Plate Size		310×620mm	
Temp. Range & Accuracy		Max. 350°C, ±0.3°C at Set Temperature	
Temp. Resolution		0.1 °C-Display, 0.5 °C-Control	
Temp. Uniformity		Specially Designed Heating Module ; Temp. Difference 10% less	
Heating Power/Consumption		2 kW	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Permissible Temp. & RH %		Ambient +5 ~ 50°C, 85% RH	
Dimension & Weight		630×360×h150mm, 16 kg	
Packing Size & Gross Weight		750×430×h220mm, 18 kg	
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	
Cat. No.	230V	DH.WHP03024	DH.WHP03025
	120V	DH.WHP02024	DH.WHP02025
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.


DAIHAN™ **Premium** 500°C High-Temp Hotplate Stirrer “**MSH500**”, Solid Ceramic Plate, 200×200mm, 80~1,500rpm With Permanently Brushless Shade Motor(BLAC), Large LCD, Optimum Insulation Layer. Accurate Temp Control, Hot-Top Indicator 고온용 디지털 가열 자력 교반기, 최적의 방열구조, 디지털 피드백 온도 제어, 내열성 · 내화학성 · 내산성 솔리드 세라믹 글라스 플레이트




A. “MSH500S” Package-set



B. “MSH500” Unit only



Controller of “MSH500” model



- Supplied by DAIHAN Korea
- **PL (Product Liability) Insurance**
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- **High-quality Solid Ceramic Plate ; Heat-/Chemical-/Acid-Resistance and Easy to Clean**
- **Optimum Insulation Fourfold Layer ; BTU + Ceramic Paper + 2 STS Insulation Layer**
- **Permanently Brushless Shade Motor & Powerful Magnet ensures Powerful Stirring and Precise Speed Control, Minimization Noise**
- * Microprocessor PID control
- * Patented Controller “Jog-Dial” and Main Touch Button - Best Durability, Best
- * Mobility, Best Convenient to Operate with Back-Light Function
- * 2 Internal K-type Temp Sensor Operating Independently ;
 - 1 Temp Control Sensor for Temperature Accuracy
 - 1 Over Temp Detection Sensor for User Safety
- * Hot Top Indicator LED Lamp ; Hot Surface(Over 50°C) Warning Lamp for User Safety
- * Excellent for Sample Drying and most General Lab Heating Applications
- * User's Self -Compensation Function ; to Control Difference between Real Temp and Indicated Temp (±20.0°C)
- * Accurate Temp Control of the Sample in the Vessel using PT100 Sensor (Optional)
- * Storage Function for the set Values ; Temp and Timer
- * Use of the Special Advanced Insulation for the safety of the Body Temp
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled)

Specifications & Ordering information

<L2>

Spec.		Model	A. “MSH500S” (Package-set)	B. “MSH500” (Unit Only)
Material	Body	Powder Coated Aluminum		
	Plate	Solid Ceramic Plate, 200×200mm		
Motor & Stirring Capa.		Permanently Brushless Shade Motor(BLAC), Max. 20 Lit		
Speed & Control Resolution		80~1,500 rpm, 5 rpm		
Temp Range & Accuracy		Max. 500℃, ±0.3℃ at Set Temperature		
Temp Resolution		0.1℃-Display, 0.5℃-Control		
Temp Uniformity		Specially Designed Heating Module ; Temp Difference 10% less		
Heating Power / Consumption		850W		
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller		Digital Feedback Controller with Jog Dial (Turn+Push) and Touch button.		
Stirring Bar		3cm bar-included, Usable ; up to 5cm bar		
Safety Circuit		Over Temp Protection, Hot Top Indicator LED Lamp		
included Standard Acc.		1ea Stir Bar, ②Temp Probe ③Rod ④Holder+Clamp		1ea Stir Bar
Permissible Temp & RH %		Ambient +5℃~50℃, 85% RH		
Dimension & Weight		200×332×h114mm, 3.8kg		
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz, with Cord/Plug		
Cat. No.	230V	DH.WMH00500		DH.WMH01500
	120V	DH.WMH00501		DH.WMH01501
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

Component of Package-set

- * Package-set ; ① Hotplate Stirrer Unit ② “**SS500**” PT100 External / Direct Contact Temp Probe,
③ “**RD200**” Stainless-steel Support Rod $\phi 12 \times L450$ mm
④ “**CL220**” Holder+Clamp+Clip for External Temp Probe

Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

- DH.WHP503060 (1) External / Direct Contact Temp Probe, “**SS500**”, Stainless-steel, $\phi 4 \times L250$ mm, Cable L190mm
DH.WHP503061 (2) External / Direct Contact Temp Probe, “**SS510**”, Teflon Coated, $\phi 4 \times L250$ mm, Cable L190mm,
DH.WHP503036 (3) Stand Rod “**RD200**”, Stainless-steel, $\phi 12 \times L450$ mm, 1.5p screw, for “**MSH500**”
DH.WHP503012 (4) Holder · Clamp · Clip, “**CL220**”, for Temp Probe

Ea. Price, ₩

DAIHAN™ 380℃ Analog & Digital Hotplate Stirrer “MSH-A” & “MSH-D”, Ceramic-Coated Plate, 80~1,500 rpm
180×180mm or 260×260mm Plate, with Accurate Temp. Control, Superior Temp Uniformity, with Certi. & Traceability
Permanently Brushless Motor(BLDC), 아날로그 & 디지털 가열 자력 교반기, 우수한 온도제어 · 온도균일성 · 세라믹 코팅 플레이트

A. 180×180mm Plate,

I. Analog Unit “MSH-20A”

II. Digital Unit & III. Package Set

“MSH-20D” “MSH-20D” Set



- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Excellent for Sample Drying and most General Lab Heating Applications
- Permanently Brushless Motor(BLDC)

I. Analog PWM Control System “MSH-A”

- * Superior Linearity PWM Control for Temperature & Stirring Speed

II. Digital Microprocessor Control System “MSH-D”

■ Patented Jog-Shuttle Control System

- * Speed & Control Resolution ; 80 ~ 1,500 rpm, 5 rpm
- * Temp Accuracy / Resolution ; $\pm 0.3^{\circ}\text{C}$ at Set Temp / 0.1°C -Display, 0.5°C -Control
- * Timer & Alarm ; 99hr 59min(with Continuous run), Error status & Timer-end
- * Safety Circuit ; Over Temp Protection
- * User's Self-Compensation ; Control Difference between Real Temp & Indicated Temp ($\pm 20.0^{\circ}\text{C}$)
- * Accurate Temp Control of the Sample in the Vessel using PT100 Sensor
- * Storage Function for the set Values ; Temp and Timer
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

III. Digital Package-set ; Complete with

- * III-A. “MSH-20D” Set ; ① Hotplate Stirrer “MSH-20D”, ② “SS500” PT100 External / Direct Contact Temp Probe, ③ “RD200” Stainless-steel Support Rod $\phi 12 \times L450\text{mm}$, ④ “CL220” Holder+Clamp+Clip for External Temp Probe
- * III-B. “MSH-30D” Set ; ① Hotplate Stirrer “MSH-30D”, ② “SS500” PT100 External / Direct Contact Temp Probe, ③ “RD100” Stainless-steel Support Rod $\phi 12.7 \times L390\text{mm}$, ④ “CL220” Holder+Clamp+Clip for External Temp Probe



① Hotplate “MSH-20D”

3-Year Warranty

III-A. Digital Package Set “MSH-20D” Set

B. 260×260mm Plate,

I. Analog Unit “MSH-30A”

II. Digital Unit & III. Package Set

“MSH-30D” “MSH-30D” Set



(II-B) Digital “MSH-30D”, 260×260mm Plate

Specifications & Ordering information

<L1>

Spec.	Model	A. 180×180mm Plate (I-A) Analog “MSH-20A” (II-A) Digital “MSH-20D”	B. 260×260mm Plate (I-B) Analog “MSH-30A” (II-B) Digital “MSH-30D”
Plate Size(mm)		180×180mm Plate-type	260×260mm Plate-type
Material	Body	Powder Coated Aluminum	Powder Coated steel
	Plate	Ceramic Coated Aluminum (Chemical / Acid Resistance, Easy to Clean)	
Motor & Stirring Capa.		Permanently Brushless Motor(BLDC, Rating ; 10w), Max. 20Lit	
Speed & Control Resolution		Digital ; 80~1500 rpm, 5 rpm / Analog ; Max. 1500 rpm	
Temp Uniformity		Specially Designed Heating Module ; Temp Difference 10% less	
Heating Power/Consumption		600W	1.2kW
Stirring Bar		3cm bar-included, Usable ; up to 5cm bar	
Permissible Temp & RH %		Ambient +5℃~50℃, 85% RH	
Dimension & Weight		206×307×h99mm, 3.3kg	286×387×h99mm, 3.5kg
Packing Size & Gross Weight		400×280×h200mm, 5kg	430×320×h200mm, 5kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	

I. Analog Unit Only		(I-A) “MSH-20A”	(I-B) “MSH-30A”
Cat. No.	230V	DH.WMH03120	DH.WMH03121
	120V	DH.WMH02120	DH.WMH02121
Price, Each. ₩			
II. Digital Unit Only		(II-A) “MSH-20D”	(II-B) “MSH-30D”
Cat. No.	230V	DH.WMH03020	DH.WMH03022
	120V	DH.WMH02020	DH.WMH02022
Price, Each. ₩			
III. Digital Package-set		(III-A) “MSH-20D” Set	(III-B) “MSH-30D” Set
Cat. No.	230V	DH.WMH03021	DH.WMH03023
	120V	DH.WMH02021	DH.WMH02023
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.


Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

- DH.WHP503060 (1) Temp Probe, “SS500”, “PT100”, $\phi 4 \times L250\text{mm}$
 DH.WHP503061 (2) Teflon Coated Temp Probe, “SS510”, $\phi 4 \times L250\text{mm}$
 DH.WHP503011 Stand Rod, “RD100”, $\phi 12.7 \times L390$ for “MSH-30A / 30D”
 DH.WHP503036 “RD200”, $\phi 12 \times L450$ for “MSH-20A / 20D”
 DH.WHP503012 Holder · Clamp · Clip, “CL220”, for Temp Probe


Ea. Price, ₩



DAIHAN™ 350℃ Systematic Multi-Hotplate Stirrer “SMHS”, 3- or 6-Places, Ceramic-Coated Plate, 80~1,500 rpm
With Digital Feedback Control, Independent Heating & Stirring Control, Permanently Brushless Motor(BLDC)
멀티 가열 자력 교반기, 우수한 온도균일성, 3- or 6- 구 개별 조절 가능




Stir Bar
A. 3 Plate-type Multi-Hotplate Stirrer “SMHS-3”



Stir Bar
B. 6(3×2) Plate-type Multi-Hotplate Stirrer “SMHS-6”

※ “SMHS-6” has 3 Ports for Connecting External Temp Probes.



- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Excellent for Multi-use and Wide Range Sample Testing
- Patented Jog-Shuttle Control System
- Permanently Brushless Motor(BLDC)

- * Digital Feedback Control by High-Performance Microprocessor
- * Ceramic Coated Plate-Chemical/Acid Resistance
- * Digital Timer Function Included ; 99 hr 59 min
- * Maximum Temperature of 350℃ with $\pm 0.3^{\circ}\text{C}$ Control Accuracy
- * Accurate Temp. Control of the Sample in the Vessel using PT100 Sensor (Optional)
- * Storage Function for the Set Values ; Temperature, RPM and Timer
- * Locking Mode for Experimental Safety Supported

Specifications & Ordering information

<L1>

Spec.	Model	A. “SMHS-3” 3 Plates	B. “SMHS-6” 6 Plates
Plate-type		3 Plate-type	6(3×2) Plate-type
Plate Size		140×140mm	
Motor & Stirring Capa.		Permanently Brushless Motor(BLDC), Max. 20Lit(H ₂ O)	
Motor Rating		9/4W×3	9/4W×6
Speed & Control Resolution		80 ~ 1,500 rpm, 5 rpm	
Temp. Range & Accuracy		Max. 350℃, $\pm 0.3^{\circ}\text{C}$ at Set Temp	
Temp. Resolution		0.1℃-Display, 0.5℃-Control	
Temp. Uniformity		Specially Designed Heating Module ; Temp Difference 10% less	
Heating Power/Consumption		400 W × 3 Plates	400 W × 6 Plates
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Digital Feedback Controller with Jog Shuttle Switch (Turn+Push)	
Stirring Bar		3cm bar-included, Usable ; up to 5cm bar	
Ports for Stand Rod		2ea	4ea
Permissible Temp. & RH %		Room Temp +5℃~50℃, 85% RH	
Material	Body	Powder Coated steel	
	Plate	Ceramic Coated Aluminum	
Dimension & Weight		510×240×h100mm, 10kg	510×410×h180mm, 18kg
Packing Size & Gross Weight		540×330×h170mm, 10kg	570×540×h230mm, 18kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	
Cat. No.	230V	DH.WMH03503	DH.WMH03506
	120V	DH.WMH02503	DH.WMH02506
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

※ “SMHS-6” has 3 Ports for Connecting External Temp Probes.

Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

- DH.WHP503060 (1) External / Direct Contact Temp Probe, “SS500”, Stainless-steel, $\phi 4 \times L250\text{mm}$, Cable L190mm
- DH.WHP503061 (2) External / Direct Contact Temp Probe, “SS510”, Teflon Coated, $\phi 4 \times L250\text{mm}$, Cable L190mm,
- DH.WHP503036 (3) Stand Rod “RD200”, Stainless-steel, 1.5p screw, $\phi 12 \times L450\text{mm}$, for “SMHS”
- DH.WHP503012 (4) Holder · Clamp · Clip, “CL220”, for Temp Probe

Ea.	Price, ₩

DAIHAN™ 50 · 80 · 120 · 250kg Automatic Ice Maker, Snow-type "IMS" & Flake-type "IMF"

With Fully Automatic System, Uniform Ice, Production & Storage, 아이스메이커



3-Year Warranty

Snow-type "IMS50" & "IMS250"

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- No Separate Storage Bin Required.
- Full Automatic System ; Continuous Ice Production
- Snow Type ; Soft Particles in the Form of Rolling Carve
- Flake Type ; A Uniform, Rigid Ice Particle in the Form of Screw Turbine

- * Ice Production and Storage Function at the Same Time
- * Useful in Medical Labs, Pharmaceutical Industry, Scientific Labs, Restaurants and Super markets in General.
- * Creates Uniform Ice with Minimal Water.
- * R404 Refrigerant and Ozone Friendly
- * Safety Water Level Sensor & Ice Level Sensor
- * Consists of Main Unit, Supply Hose, Drain Hose, Scoop, Power Cable, Operation Manual

(1) Snow Ice Type

<L3>

Spec.	Model	"IMS50"	"IMS250"
Ice Type		Snow Ice Type	
Production / Day		50 kg	250 kg
Storage Capacity		25 Lit	86 Lit
Compressor		1/2 HP, R404	1 HP, R404
Door		Front PC Door	
Safety Circuit		Power Circuit Breaker, Ice Level Sensor, Water Level Sensor	
Dimensions (w×d×h)		360×440×620mm	500×550×1195mm
Weight		40 kg	90 kg
Power Consumption		650 W	1.8 kW
Electric Reqts.		230V, 50 or 60Hz	
Cat. No.	60Hz	DH.Ims8050	DH.Ims8250
	50Hz	DH.Ims8055	DH.Ims8255
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

(2) Flake Ice Type



3-Year Warranty

Flake-type "IMF80" & "IMF120"

Spec.	Model	"IMF80"	"IMF120"
Ice Type		Flake Ice Type	
Production / Day		80 kg	120 kg
Storage Capacity		58 Lit	120 Lit
Compressor		1/3 HP, R404	1/3 HP, R404
Door		Front PC Door	
Safety Circuit		Power Circuit Breaker, Ice Level Sensor, Water Level Sensor	
Dimensions (w×d×h)		600×560×740mm	675×560×1065mm
Weight		61 kg	83 kg
Power Consumption		800 W	800 W
Electric Reqts.		230V, 50 or 60Hz	
Cat. No.	60Hz	DH.Imf8080	DH.Imf8120
	50Hz	DH.Imf8085	DH.Imf8125
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

DAIHAN Scientific INCUBATORS

CO₂	P. 185~186
General Low Temp, Mini	P. 187
General Purpose, Forced-air	P. 188~189
, Gravity-air	P. 190~192
, Multi-room, Forced-air	P. 193
Heating Block / Dry Bath Incubators	P. 154~158
Humidity / Temp Test Chambers	P. 194~196
Low Temp(B.O.D)	P. 197~198
Plant Growth Chambers	P. 199
Shaking, Bench Top-type	P. 200
, Front Door-type	P. 201
, Mini, CO ₂ & Low Temp	P. 185, 187
, Multi-stack Refrigerated	P. 203
, Top Door-type	P. 202



CO₂ Incubator Mini



CO₂ Incubator



Forced- & Gravity-air Incubator



Humidity / Temp. Chamber



Low Temp.(B.O.D) Incubator



Plant Growth Chamber



Shaking Incubator

DAIHAN™ 20Lit Mini CO₂ Incubator & CO₂ Shaking Incubator "IC-20 & ICS-20", 4℃~45℃ & 0~10% CO₂
 Programmable PID Controlled 0.1℃ & 0.1%, Compact Design for Saving Space/Money, Ideal for Cell/Tissue Culture
 With Cooling/Heating system of Peltier Element/Auto-Tuning/Safety Device, 저온 미니 CO₂ 배양기 & CO₂ 진탕 배양기

20Lit 4℃~45℃
 0~10% CO₂
 Cooling/Heating



3-Year Warranty

A. "IC-20" Mini CO₂ Incubator, 20Lit



B. "ICS-20" CO₂ Shaking Incubator Set



C. "SHK2.ODI" Separable Orbital Shaker

; It can be Used Individually and Be Accepted



Fig. Using "IC-20" as a Chest-type

A. Standard CO₂ Incubator "IC-20"

B. Assembly CO₂ Shaking Incubator Package Set "ICS-20"

= Standard CO₂ Incubator "IC-20" + Orbital Shaker "SHK2.ODI"

C. Multiuse Digital Orbital Motion Shaker "SHK2.ODI" ; Can be Used Individually.

"SHK2.ODI" is same with existing Orbital Shaker "SHK2.ODI" except Platform, and can be Accepted Optional Various Platforms ; See the Shaker Section

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- Applications ; Cell/Tissue Culture & Transfer, Embryo/Microorganism Culture
- NDIR CO₂ Sensor ; Dual Wavelength Type
Rapid Initial Response Time, not Affected by User Environment
- ABS-Resin Chamber & Housing ; Light-weight/Compact Design
Good for Money & Space Saving and Personal Use.

- * Cooling/Heating System of PELTIER ELEMENT ;
- Cooling and Heating by using Thermoelectric Semiconductor Element
- Minimizes Noise(<20dB) and Power Consumption
- * Auto-Tuning ; Rapid Response and High Stability According to the Thermal Properties and the Rate of Thermal Response of the Object
- * Safety Device ; Over Temperature & Over Current Protection, Sensor Error Detection
- * 4-Connection Device to Join a Shaker "SHK2.ODI" in the Bottom outside of the CO₂ Incubator "IC-20"

Specifications & Ordering Information

<L1>

Spec.	Model	A. CO ₂ Standard Incubator "IC-20"	B. CO ₂ Shaking Incubator Set "ICS-20"
Capacity (Lit)		20Lit Working Chamber	
Temp Range		4℃~45℃, 0.1℃ ; Digital PID Control	
Accuracy		±0.2℃ (No Load)	
Uniformity		±0.2℃ (at 37℃)	
CO ₂ Range & Control		0 ~ 10%, 0.1% ; NDIR CO ₂ Sensor	
Cooling/Heating		Peltier Element Module	
Air Flow-type		Gravity Flow Circulation	
1 Shelf and 3 Location Brackets		1 Perforated ABS Shelf(w210×d230×Thick 8mm) included Adjustable Shelf Location Height of 3-Stage	
Condensate Drain		Natural Evaporation	
Incubator Material		"ABS-Resin" Chamber & Housing, Except Shaking Base	
Incubator Interior		w213×d265×h360mm Chamber Only	
Dimension Exterior		w280×d330×h510mm	w280×d330×h670mm; Including Shaking Base
C. SHAKING-Base ; Connecting the Multiuse Shaking "SHK2.ODI"			
- Motion & Speed			- Orbital, 10 ~ 300rpm, 10mm Stroke
- Controller & Display			- Large Size Digital LCD with Back-Light Function
- Timer & Alarm			- 99hr 59min with ∞ Run, Error Status & Time End
- Motor & Max-Load			- Brushless DC Motor, Max-Load 10kg
- Safety Circuit			- Overload Protector of Motor
- Material			- Powder Coated Steel Body
- Size & Weight			- w276×d375×h138mm, 13kg
Net Weight (kg)		10kg	10kg+13kg Shaking Base ("SHK2.ODI") = 23kg
Power Consumption		56W	56W+45W (Shaking)
Electric Reqts.		Free Voltage; 1Phase, AC90-240V, 50/60Hz	
Cat. No.	230V	DH.IncC020	DH.IncCS020
	120V	DH.IncC020.1	DH.IncCS020.1
Price, Set. ₩			

C. SHAKING-Base ; Orbital Shaker

Ea. Price, ₩

DH.SHK2.ODI Shaker "SHK2.ODI", 230V
 DH.SHK2.ODI.1 , 120V

※ Note! Details of Optional Various Platforms ; See The "Shakers" Section

Optional CO₂ Regulator



Ea. Price, ₩

DH.Reg8004 (1) Regulator, Single-Stage, 30psi, CO₂

Spare Shelf

Ea. Price, ₩

DH.Inc0003 Perforated Shelf, ABS, 210×230mm

DAIHAN™ Air-Jacketed CO₂ Incubator "ICA", 101 & 150 Lit, 0~20% CO₂, Up to 50°C, ±0.1°C

With Precision CO₂ Sensor, UV Lamp 254nm, Microprocessor PID Control, Fan Forced-Convection, 2 Shelf Included, CO₂ 인큐베이터



up to 50°C

3-Year Warranty

Air-Jacketed CO₂ Incubator "ICA"



Tempered Safety Inner Glass Door & Enhanced Door Packing



(1) Perforated STS Shelf (Included)



(2) UV Lamp



(3) HEPA Filter



(4) Regulator, Single-stage

■ Supplied by DAIHAN Korea

■ Precise CO₂ Sensor in Stable & Durable Housing, CO₂ Range 0~20%

- Dual Wavelength Detection, NDIR(Non Dispersive Infra-red) Sensor
- Automatic Temperature Compensation and Auto Zero Calibration.

- * Ideal for Medicine, Biochemistry and Agricultural Science Research, Industrial Production Departments of Biological Cells, Tissues, Bacterial Culture, etc.
- * All Seamless Stainless-steel Interior for Corrosion Resistance and Decontamination
- * Tempered Safety Inner Glass Door for easy Sample View without Door Opening
- * Enhanced Door Packing to Keep Excellent Temp & CO₂ Uniformity
- * Air Filtration System(0.3μm HEPA Filter) Minimizes Contamination, Efficiency up to 99.97%.
- * Microprocessor PID Control ; Temp(-50°C), Timer, Alarm, Auto-tuning
- * High-Quality LCD Display with Back-Light Function
- * 10 Steps Programmable ; Temperature Control and 999 Cycle 10 Steps.
- * Detachable Shelf Frame ; Adjustable Shelf
- * Including Water Container, and 2 Stainless-steel Perforated Shelf
- * Temp & CO₂ Recovery Time ; 10min (Door Open 30sec)
- * Operating Temperature 18°C~32°C
- * Filter Replace Period ; Approx. 6~9 Month

Specifications & Ordering information

<L2>

Spec.	Model	"ICA-100" Max. Capa. 101Lit	"ICA-150" Max. Capa. 150Lit
Capacity		101 Lit	150 Lit
Heating Mechanism		Air-jacket	
CO ₂ Sensor		Digital Signal Processing NDIR Sensor	
Temp Range & Accuracy		Room Temp +5°C~50°C, Accuracy ±0.1°C at 20°C~25°C	
CO ₂ Range & Accuracy		0 ~ 20% CO ₂ Concentration, Accuracy ± 0.1% at 5%	
Heat up Time (to 37°C)		30 min to 37°C	
CO ₂ Concentration Time (to 37°C)		25 min to 5%	
Heating Power		300 W	570 W
CO ₂ Outlet		1 kg/cm ³	
UV Lamp		254 nm × 1ea	
Controller		Programmable Microprocessor Digital PID Controller with LCD Back Light	
Shelf (Included)		2 Stainless-steel Perforated Shelf	
Filter		2 HEPA Filter(0.3μm), CO ₂ Input & Output line	
Safety		Temp & CO ₂ Alarm, Over Temp Protector, Sensor Error Display, Leakage Breaker	
Material		Interior ; Stainless-steel Polished, Exterior ; Powder Coated Steel	
Door		Inner Door ; Tempered Safety Glass Door, External Door ; Powder Coated Steel, Magnetic Packing Door	
Dimension (w × d × h) & Weight	Interior	450 × 450 × h500mm	495 × 515 × h600mm
	Exterior	600 × 600 × h600mm, 65kg	660 × 645 × h1010mm, 77kg
	Packing	760 × 800 × h1120mm, 80kg	730 × 810 × h1130mm, 92kg
Power Consumption		340 W	610 W
Electric Reqts.		220V, 50/60Hz, 1 Phase	
Cat. No.	220V	DH.Inc8001	DH.Inc8002
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Accessory

- DH.IncA1001 (1) Perforated Shelf, STS, for "ICA-100", w440 × d400mm
 DH.IncA1002 , for "ICA-150", w485 × d445mm
 DH.IncA1005 (2) UV Lamp 254 nm, φ 30 × 205mm, for "ICA"
 DH.IncA1006 (3) HEPA Filter, for "ICA"

Ea. Price, ₩

Optional Accessory

- DH.Reg8004 (4) Regulator, Single-stage, 30psi, CO₂

Ea. Price, ₩

Air-jacketed CO₂ Incubator with Validation (IQ, OQ) Service

- DH.IncA5001 Validation Service(IQ, OQ), with "ICA-100"
 DH.IncA5002 , with "ICA-150"

Set. Price, ₩

DAIHAN™ 4℃~45℃ 20 Lit Mini-Low Temperature Incubator & Shaking Incubator “IR-20 & IRS-20”

2-Step Programmable PID Controlled 0.1℃, Compact Design for Saving Space/Money, Ideal for Culture & Storage of Microorganism/Clone
With Cooling/Heating system of Peltier Element/Auto-Tuning/Safety Device, 저온 미니-배양기 & 진탕 배양기, 20Lit, 4℃~45℃/0.1℃

4℃~45℃
Cooling / Heating



3-Year Warranty

A. “IR-20” Mini-Low Temp Incubator, 20Lit



20Lit
Orbital, Shake
10~300rpm, 8mm Str.
4℃~45℃
Cooling / Heating

3-Year Warranty

B. “IRS-20” Low Temp Shaking Incubator Set



C. “SHK2.ODI” Separable Orbital Shaker

; It can be Used Individually and Be Accepted



Fig. Using “IR-20” as a Chest-type

A. Standard Incubator “IR-20”

B. Assembly Shaking Incubator Package Set

“IRS-20” = “IR-20” Incubator + Orbital Shaker “SHK2.ODI”

C. Multiuse Digital Orbital Motion Shaker “SHK2.ODI” ; Can be Used Individually.

It is same with existing Orbital Shaker “SHK2.ODI” except Platform and The Shaker “SHK2.ODI” can be Accepted Optional Various Platforms; See “Shakers” Section

■ Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance

■ Applications Microorganism / Clone Culture & Storage, Reaction Experiment of Enzyme Immunity & Bacteriophage in Vitro Packaging, DNA Ligation, etc.

■ 2-Step Temperature Control(PID Micro Processor) ; With 2-Step Setting Temp Microorganism is Cultivated and Preserved by One Time(Setting time 99 hours)

* Cooling / Heating System of PELTIER ELEMENT ;

- Cooling and Heating by using Thermoelectric Semiconductor Element

- Minimizes Noise (<20dB) and Power Consumption

* Auto Tuning ; Rapid Response and High Stability According to the Thermal Properties and the Rate of Thermal Response of the Object

* Safety Device ; Over Temperature & Over Current Protection, Sensor Error Detection

* There are 4-Connection Device to Join a Shaker “SHK2.ODI” in the Bottom outside of the Incubator “IR-20”

Specifications & Ordering Information

<L1>

Spec.	Model	A. Standard Incubator “IR-20”	B. Shaking Incubator Set “IRS-20”
Capacity (Lit)		20Lit Working Chamber	
Temp Range		4℃~45℃, 0.1℃; Digital PID Control	
Accuracy		±0.2℃(No Load)	
Uniformity		±0.2℃(at 37℃)	
Cooling / Heating		Peltier Element Module	
2-Step Operation		Ideal for Cultivation and Preservation	
Air Flow-type		Gravity Flow Circulation	
1 Shelf and 3 Location Brackets		1 Perforated ABS Shelf(w210×d230×Thick 8mm) included Adjustable Shelf Location Height of 3-Stage	
Condensate Drain		Natural Evaporation	
Incubator Material		“ABS-Resin” Chamber & Housing, Except Shaking Base	
Incubator Interior		w213×d265×h360mm Chamber Only	
Dimension Exterior		w280×d330×h510mm	w280×d330×h670mm; Including Shaking Base
C. SHAKING-Base ; Connecting the Multiuse Shaking “SHK2.ODI”			
- Motion & Speed			- Orbital, 10 ~ 300rpm, 10mm Stroke
- Controller & Display			- Large Size Digital LCD with Back-Light Function
- Timer & Alarm			- 99hr 59min with ∞ Run, Error Status & Time End
- Motor & Max-Load			- Brushless DC Motor, Max-Load 10kg
- Safety Circuit			- Overload Protector of Motor
- Material			- Powder Coated Steel Body
- Size & Weight			- w276×d375×h138mm, 13kg
Net Weight (kg)		7kg	7kg+13kg Shaking Base (“SHK2.ODI”) = 20kg
Power Consumption		56W	56W+45W (Shaking)
Electric Reqts.		Free Voltage; 1Phase, AC90-240V, 50 / 60Hz	
Cat. No.	230V	DH.IncR020	DH.IncRS020
	120V	DH.IncR020.1	DH.IncRS020.1
Price, Set. ₩			

C. SHAKING-Base ; Orbital Shaker

DH.SHK2.ODI Shaker “SHK2.ODI”, 230V
DH.SHK2.ODI.1 , 120V

Ea. Price, ₩

※ Note! Details of Optional Various Platforms ; See The “Shakers” Section

Spare Shelf

DH.Inc0003 Perforated Shelf, ABS, 210×230mm


Ea. Price, ₩

Optional Power Cord for Vehicles

DH.Inc0004 Cigarette lighter Power Supply Cord

Ea. Price, ₩

DAIHAN™ Forced-Air Incubator “ThermoStable™ IF”, 3-Side Heating Zone, 50 · 105 · 155 Lit, **Medicaluse(230V)**
 With 2 Wire Shelf, Digital PID Control, Jog-Dial & Push Button, Digital LCD with Backlight, Certi. & Traceability, Up to 70℃, ±0.2℃
 강제 순환식 배양기/인큐베이터, 디지털 퍼지 컨트롤, 우수한 온도균일성



3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism.



WiseRemote Control Software ;

Connection of Temperature Devices with Interface RS232C to a PC

* You can Download this Software from DAIHAN-Sci.com Shopping-mall

- Supplied by DAIHAN Korea ■ GD(Good Design)-mark ■ UL, CUL, and CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Patented Jog-Dial Control System ■ PL(Product Liability) Insurance
- * Optimized Air Flow by Forced-Air Mechanism
- * 3-Side Heating Mechanism ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- * Various Application ; Microorganism Culture, Animal & Plant Cell Culture, Constant Temp, Germination, and so on
- * Digital PID Control System Implementing Superior Temperature Accuracy
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be disabled)

Specifications & Ordering information

<M1>

Spec. Model		ThermoStable™ IF-50 Max. Capa. 50Lit	ThermoStable™ IF-105 Max. Capa. 105Lit	ThermoStable™ IF-155 Max. Capa. 155Lit
Capacity		50 Lit	105 Lit	155 Lit
Dimension (w×d×h)	Interior (mm)	370×350×419	485×410×533	550×475×600
	Exterior (mm)	518×637×700	635×695×841	700×760×906
Heater Power		200W	300W	400W
Temp	Range	Room Temp +5℃~70℃		
	Fluctuation	± 0.2℃ at 37℃, ± 0.3℃ at 50℃		
	Variation	± 0.5℃ at 37℃, ± 0.9℃ at 50℃	± 0.4℃ at 37℃, ± 0.7℃ at 50℃	
	Sensor	PT 100		
Control Resolution		±0.1℃		
Heat-up Time		17 min. to 37℃, 28 min. to 50℃	17 min. to 37℃, 30 min. to 50℃	
Recovery Time (Door open 30sec)		3 min. to 37℃, 4 min. to 50℃		
Controller		Digital PID Control by Advanced Microprocessor Jog-Dial with Push Button		
RS232 Port		Available to connect with PC		
Display		Digital LCD with Back Light		
Timer		99hr 59 min (delay/continuous function)		
Storage		Temp and Time		
Material	Internal	Stainless-steel (#304), Corrosion Resistant		
	External	Powder Coated steel		
	Inner Door	Inner Tempered Glass Door		
Shelf (Included)		2 PE-coated steel Wire Shelf Load per Shelf ; 16kg		
Ventilation		1 Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm		
Safety device		Over Temp and Over Current Protector, Sensor Error Detector		
Alarm		Error Status and Timer-end		
Net Weight (kg)		46	66	81
Packing Size (w×d×h) (mm)		628×752×775	710×731×1,018	731×810×1,018
Shipping Weight (kg)		58	74	92
Power Consumption		216.3 W	302.9 W	438.1 W
Energy Consumption, at 37℃, at 50℃		28 Wh	43 Wh	58 Wh
		84 Wh	129 Wh	145 Wh
Electrical Reqs.		120V, 60Hz or 230V, 50/60Hz		
Cat. No.	230V	DH.WIF03050	DH.WIF03105	DH.WIF03155
	120V	DH.WIF04050	DH.WIF04105	DH.WIF04155
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도


Forced Convection Incubator with Validation (IQ, OQ) Service

DH.WIF31050 Validation Service(IQ, OQ), “IFV050”, with “ThermoStable IF-50”
 DH.WIF31105 “IFV105”, with “ThermoStable IF-105”
 DH.WIF31155 “IFV155”, with “ThermoStable IF-155”

Ea. Price, ₩

※ Note! Spare PE-coated Steel Wire Shelf & Optional Stainless-steel Wire Shelf
 ⇒ See The “Incubators, General Purpose, Accessories” Section

DAIHAN™ SMART Forced-Air Incubator “ThermoStable™ SIF”, 50 · 105 · 155 Lit, 3-Side Heating Zone, **Medicaluse(230V)**
 With Smart-Lab™ Controller, 4" Full Touch Screen TFT LCD Display, Fuzzy-PID Control, 2 Wire Shelf, Certi. & Traceability, Up to 70°C, $\pm 0.2^\circ\text{C}$
 스마트 강제 순환식 배양기/인큐베이터, 우수한 온도균일성, WiRe™ 서비스, 4인치 풀터치스크린 조절기



3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism.

WiRe™ App Service



4" Full Touch Screen TFT LCD (Smart-Lab™ Controller)

Controlled by Mobile-Phone

It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.



WiRe™ App Service

- Supplied by DAIHAN Korea ■ GD(Good Design)-mark ■ UL, CUL, and CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Smart-Lab™ Controller with WiRe™ App Service. ■ PL(Product Liability) Insurance

- * Self-Diagnostic Function
- * Automatic Data Recording
- * Variable Program Setting for Multiple Patterns with Multiple Segments
- * The Recorded Data Could be transferred to PC by USB memory
- * Internet Connectivity with WiFi
- * Controller SW can be upgraded through USB port.
- * Optimized Air Flow by Forced-air Mechanism
- * 3-Side Heating Mechanism ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- * Various Application ; Microorganism Culture, Animal & Plant Cell Culture, Constant Temp, Germination, and so on

Specifications & Ordering information

<M1>

Spec.		Model	ThermoStable™ SIF-50 Max. Capa. 50Lit	ThermoStable™ SIF-105 Max. Capa. 105Lit	ThermoStable™ SIF-155 Max. Capa. 155Lit
Capacity			50 Lit	105 Lit	155 Lit
Dimension (w×d×h)	Interior (mm)		370×350×419	485×410×533	550×475×600
	Exterior (mm)		518×637×700	635×695×841	700×760×906
Heater Power			200W	300W	400W
Temp	Range		Room Temp +5°C~70°C		
	Fluctuation		$\pm 0.2^\circ\text{C}$ at 37°C, $\pm 0.3^\circ\text{C}$ at 50°C		
	Variation		$\pm 0.5^\circ\text{C}$ at 37°C, $\pm 0.9^\circ\text{C}$ at 50°C	$\pm 0.4^\circ\text{C}$ at 37°C, $\pm 0.7^\circ\text{C}$ at 50°C	
	Sensor		PT 100		
Control Resolution			$\pm 0.1^\circ\text{C}$		
Heat-up Time			17 min. to 37°C, 28 min. to 50°C	17 min. to 37°C, 30 min. to 50°C	
Recovery Time (Door open 30sec)			3 min. to 37°C, 4 min. to 50°C		
Controller			Smart-Lab™ Controller		
RS232			Port Available to connect with PC		
Display			4" Full Touch Screen TFT LCD		
Timer			99hr 59 min (delay/continuous function)		
Storage			Temp and Time		
Material	Internal		Stainless-steel (#304), Corrosion Resistant		
	External		Powder Coated steel		
	Inner Door		Inner Tempered Glass Door		
Shelf (Included)			2 PE-coated steel Wire Shelf Load per Shelf ; 16 kg		
Ventilation			1 Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm		
Safety device			Over Temp and Over Current Protector, Sensor Error Detector		
Alarm			Error Status and Timer-end		
Net Weight (kg)			46	66	81
Packing Size (w×d×h) (mm)			628×752×775	710×731×1,018	731×810×1,018
Shipping Weight (kg)			58	74	92
Power Consumption			216.3 W	302.9 W	438.1 W
Energy Consumption, at 37°C, at 50°C			28 Wh	43 Wh	58 Wh
			84 Wh	129 Wh	145 Wh
Electrical Reqts.			120V, 60Hz or 230V, 50/60Hz		
Cat. No.	230V		DH.SWIF03050	DH.SWIF03105	DH.SWIF03155
	120V		DH.SWIF04050	DH.SWIF04105	DH.SWIF04155
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도


SMART Incubator with Validation (IQ, OQ) Service

DH.SWIF31050 Validation Service(IQ, OQ), "SIFV050", with "SIF-50"
 DH.SWIF31105 "SIFV105", with "SIF-105"
 DH.SWIF31155 "SIFV155", with "SIF-155"

Ea. Price, ₩

※ Note! Spare PE-coated Steel Wire Shelf & Optional Stainless-steel Wire Shelf
 ⇒ See The "Incubators, General Purpose, Accessories" Section

DAIHAN™ Gravity-Air Incubator “**ThermoStable™ IG**”, 3-Side Heating Zone, 32 · 50 · 105 · 155 Lit, **Medicaluse(230V)**
With 2 Wire Shelf, Digital PID Control, Jog-Dial & Push Button, Digital LCD with Backlight, Certi. & Traceability, Up to 70℃, ±0.2℃
자연 대류식 배양기/인큐베이터, 디지털 퍼지 컨트롤, 우수한 온도 균일도와 정확도



3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism.



WiseRemote Control Software ;
Connection of Temperature Devices
with Interface RS232C to a PC
* You can Download this Software from
DAIHAN-Sci.com Shopping-mall

- Supplied by DAIHAN Korea ■ GD(Good Design)-mark ■ UL, CUL, and CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Patented Jog-Dial Control System ■ PL(Product Liability) Insurance
- * Optimized Air Circulation by Gravity-Air preventing Cross Contamination of Samples
- * 3-Side Heating Mechanism ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- * Various Application ; Microorganism Culture, Animal & Plant Cell Culture, Constant Temp Germination, and so on
- * Digital PID Control System Implementing Superior Temperature Accuracy
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be disabled)

Specifications & Ordering information

<M1>

Spec.		Model	ThermoStable™ IG-32 Max. Capa. 32L	ThermoStable™ IG-50 Max. Capa. 50L	ThermoStable™ IG-105 Max. Capa. 105L	ThermoStable™ IG-155 Max. Capa. 155L
Capacity			32 Lit	50 Lit	105 Lit	155 Lit
Dimension (w×d×h)	Interior (mm)	312×290×360	370×350×420	485×410×535	550×475×600	
	Exterior (mm)	458×477×624	518×537×701	635×595×841	700×660×906	
Heater Power			150W	200W	300W	400W
Temp	Range	Room Temp +5℃~70℃				
	Fluctuation	± 0.2℃ at 37℃, ± 0.3℃ at 50℃				
	Variation	± 0.6℃ at 37℃, ± 1.0℃ at 50℃				
	Sensor	PT 100				
Control Resolution			±0.1℃			
Heat-up Time			25 min. to 37℃, 40 min. to 50℃			30 min. to 37℃, 40 min. to 50℃
Recovery Time (Door open 30sec)			8 min. to 37℃ 12 min. to 50℃		10 min. to 37℃ 15 min. to 50℃	
Controller			Digital PID Control by Advanced Microprocessor Jog-Dial with Push Button			
RS232			Available to connect with PC			
Display			Digital LCD with Back Light			
Timer			99hr 59 min (delay/continuous function)			
Storage			Temp and Time			
Material	Internal	Stainless-steel (#304), Corrosion Resistant				
	External	Powder Coated steel				
	Inner Door	Inner Tempered Glass Door				
Shelf (Included)			2 PE-coated Steel Wire Shelf Load per Shelf ; 16 kg			
Ventilation			1 Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm			
Safety device			Over Temp and Over Current Protector, Sensor Error Detector			
Alarm			Error Status and Timer-end			
Net Weight (kg)			32	40	59	71
Packing Size (w×d×h) (mm)			568×592×698	628×652×775	731×710×1,018	796×775×1,083
Shipping Weight (kg)			39	53	67	78
Power Consumption			159 W	169.9 W	292.9 W	420.5 W
Energy Consumption, at 37℃ 						

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Gravity Convection Incubator with Validation (IQ, OQ) Service

	Ea.	Price, ₩
DH.WIG31032 Validation Service(IQ, OQ), "IGV032", with "ThermoSatble IG-32"		
DH.WIG31050 "IGV050", with "ThermoSatble IG-50"		
DH.WIG31105 "IGV105", with "ThermoSatble IG-105"		
DH.WIG31155 "IGV155", with "ThermoSatble IG-155"		

※ Note! Spare PE-coated Steel Wire Shelf & Optional Stainless-steel Wire Shelf
⇒ See The "Incubators, General Purpose, Accessories" Section

DAIHAN™ SMART Gravity-Air Incubator “ThermoStable™ SIG”, 32 · 50 · 105 · 155 Lit, *Medicaluse*(230V)

With Smart-Lab™ Controller, 4" Full Touch Screen TFT LCD Display, Fuzzy-PID Control, 2 Wire Shelf, Certi. & Traceability, Up to 70°C, $\pm 0.2^\circ\text{C}$
 스마트 자연 대류식 배양기/인큐베이터, 우수한 온도 균일도와 정확도, WiRe™서비스, 4인치 풀터치스크린 조절기

*Medicaluse*

up to 70°C

3-Year Warranty

- Supplied by DAIHAN Korea ■ GD(Good Design)-mark ■ UL, CUL, and CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Smart-Lab™ Controller with WiRe™ App Service. ■ PL(Product Liability) Insurance

- * Self-Diagnostic Function
- * Automatic Data Recording
- * Variable Program Setting for Multiple Patterns with Multiple Segments
- * The Recorded Data could be transferred to PC by USB memory
- * Internet Connectivity with WiFi
- * Controller SW can be upgraded through USB port.
- * Optimized Air Circulation by Natural Flow preventing Cross Contamination of Samples
- * 3-Side Heating Mechanism ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- * Various Application ; Microorganism Culture, Animal & Plant Cell Culture, Constant Temp, Germination, and so on

Specifications & Ordering information

<M1>

Spec.		Model	ThermoStable™ SIG-32 Max. Capa. 32L	ThermoStable™ SIG-50 Max. Capa. 50L	ThermoStable™ SIG-105 Max. Capa. 105L	ThermoStable™ SIG-155 Max. Capa. 155L
Capacity			32 Lit	50 Lit	105 Lit	155 Lit
Dimension (w×d×h)	Interior (mm)		312×290×360	370×350×420	485×410×535	550×475×600
	Exterior (mm)		458×477×624	518×537×701	635×595×841	700×660×906
Heater Power			150W	200W	300W	400W
Temp	Range		Room Temp +5℃~70℃			
	Fluctuation		± 0.2℃ at 37℃, ± 0.3℃ at 50℃			
	Variation		± 0.6℃ at 37℃, ± 1.0℃ at 50℃			
	Sensor		PT 100			
Control Resolution			±0.1℃			
Heat-up Time			25 min. to 37℃, 40 min. to 50℃			30 min. to 37℃, 40 min. to 50℃
Recovery Time (Door open 30sec)			8 min. to 37℃ 12 min. to 50℃		10 min. to 37℃ 15 min. to 50℃	
Controller			Smart-Lab™ Controller			
RS232			Available to connect with PC			
Display			4" Full Touch Screen TFT LCD			
Timer			99hr 59 min (delay/continuous function)			
Storage			Temp and Time			
Material	Internal		Stainless-steel (#304), Corrosion Resistant			
	External		Powder Coated steel			
	Inner Door		Inner Tempered Glass Door			
Shelf (Included)			2 PE-coated Steel Wire Shelf Load per Shelf ; 16kg			
Ventilation			1 Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm			
Safety device			Over Temp and Over Current Protector, Sensor Error Detector			
Alarm			Error Status and Timer-end			
Net Weight (kg)			32	40	59	71
Packing Size (w×d×h) (mm)			568×592×698	628×652×775	731×710×1,018	796×775×1,083
Shipping Weight (kg)			39	53	67	78
Power Consumption			159 W	169.9 W	292.9 W	420.5 W
Energy Consumption, at 37℃ 						

* Other Specifications are available upon Customer's Request.


※ 운임 별도

SMART Incubator with Validation (IQ, OQ) Service

DH.SWIG31032	Validation Service(IQ, OQ), “SIGV032”, with “SIG-32”
DH.SWIG31050	“SIGV050”, with “SIG-50”
DH.SWIG31105	“SIGV105”, with “SIG-105”
DH.SWIG31155	“SIGV155”, with “SIG-155”

Ea. Price, ₩

※ Note! Spare PE-coated Steel Wire Shelf & Optional Stainless-steel Wire Shelf
 ⇒ See The “Incubators, General Purpose, Accessories” Section



3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism.

>> WiRe™ App Service



It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.

4" Full Touch Screen TFT LCD
(Smart-Lab™ Controller)

Controlled by Mobile-Phone



WiRe™ App Service



Wire Shelf, PE-Coated Steel
(Included)

Spare PE-Coated Steel Wire Shelf

(1) for Forced-Air Incubator "ThermoStable™ IF & SIF"

- DH.WIG11050** Wire Shelf, "IGS050", PE-coated Steel, with 2 Brackets, for 50 Lit, w340×d310mm
- DH.WIG11105** Wire Shelf, "IGS105", PE-coated Steel, with 2 Brackets, for 105 Lit, w455×d375mm
- DH.WIG11155** Wire Shelf, "IGS155", PE-coated Steel, with 2 Brackets, for 155 Lit, w520×d440mm

<L1>

Ea. Price, ₩

(2) for Gravity-Air Incubator "ThermoStable™ IG & SIG"

- DH.WIG11032** Wire Shelf, "IGS032", PE-coated Steel, with 2 Brackets, for 32 Lit, w280×d260mm
- DH.WIG11050** Wire Shelf, "IGS050", PE-coated Steel, with 2 Brackets, for 50 Lit, w340×d310mm
- DH.WIG11105** Wire Shelf, "IGS105", PE-coated Steel, with 2 Brackets, for 105 Lit, w455×d375mm
- DH.WIG11155** Wire Shelf, "IGS155", PE-coated Steel, with 2 Brackets, for 155 Lit, w520×d440mm

Ea. Price, ₩



Wire Shelf, Stainless-steel
(Optional)

Optional Stainless-steel Wire Shelf

(1) for Forced-Air Incubator "ThermoStable™ IF & SIF"

- DH.WON11050** Wire Shelf, "ONS050", Stainless-steel, with 2 Brackets, for 50 Lit, w340×d310mm
- DH.WON11105** Wire Shelf, "ONS105", Stainless-steel, with 2 Brackets, for 105 Lit, w455×d375mm
- DH.WON11155** Wire Shelf, "ONS155", Stainless-steel, with 2 Brackets, for 155 Lit / 305Lit, w520×d440mm

<L1>

Ea. Price, ₩

(2) for Gravity-Air Incubator "ThermoStable™ IG & SIG"

- DH.WIG12032** Wire Shelf, "IGSN032", Stainless-steel, with 2 Brackets, for 32 Lit, w280×d260mm
- DH.WON11050** Wire Shelf, "ONS050", Stainless-steel, with 2 Brackets, for 50 Lit, w340×d310mm
- DH.WON11105** Wire Shelf, "ONS105", Stainless-steel, with 2 Brackets, for 105 Lit, w455×d375mm
- DH.WON11155** Wire Shelf, "ONS155", Stainless-steel, with 2 Brackets, for 155 Lit / 305Lit, w520×d440mm

Ea. Price, ₩

DAIHAN™ Multi-Room Incubator “WIM”, Forced-Air, 125 Lit × 4 Rooms, with Certi. & Traceability, Up to 60/70°C
With 8 Wire Shelf, Digital Fuzzy Control, Independent Chamber Control, Max. Temp Uniformity & Accuracy
4룸 배양기/인큐베이터, 강제 순환식, 디지털 퍼지 제어, 개별 조절기 장착



3-Year Warranty



4 × 125Lit Rooms
Independent Fuzzy Control System
in Each Chamber

Multi-Room Incubator, Forced-Air, “WIM-R4”

- Supplied by DAIHAN Korea
- UL, CUL, and CE Certified
- PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Patented Jog-Shuttle Control System

- * Variable Application ; Microorganism Culture, Cell Culture, Constant Temp. Germination
- * Each Chamber has an Independent Controller for Different Temp and Time Control
- * Forced-air Flow by Silent Fan Motor Implements the Highest Uniformity of Temp in Each Chamber
- * Storage Function for Set Values of Temperature and Time
- * Stainless-steel Chamber and PE-coated Shelves (2ea for Each Chamber) - Corrosion Resistance
- * Inner Tempered Glass Door ; to Observe the Samples without Interruption of Incubation
- * Over Temperature & Over Current Protection, Sensor Error Detection and Leakage Breaking.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

Specifications & Ordering information

<L1>

Spec.		Model	(1) “WIM-4”	(2) “WIM-R4”	(3) “WIM-RL4” with illumination
Temp. Range			R/T +5℃~70℃	0℃ to 60℃	0℃ to 60℃
Circulation			Forced-air(Blower FAN)		
Capacity			125 Lit×4 Rooms		
Illumination			None		FRL 36W×2 (×4 Rooms)
Temp.	Accuracy	±0.2℃ at 37℃			
	Uniformity	±0.3℃ at 37℃, ±0.7℃ at 50℃	±0.5℃ at 37℃		
Heat-up Time		15min. to 37℃ 30min. to 50℃	20min. to 37℃ 30min. to 50℃		
Recovery Time (Door open 30sec)		8min. to 37℃ 10min. to 50℃	10min. to 50℃ 12min. to 50℃		
Heating / Compressor		500 W×4ea / None	650 W×4ea / 1/4 HP (×4 Rooms)		
Timer & Alarm		99hr 59min. (with delay / Continuous run.), Error status & Timer-end			
Display		Digital LCD with Back-Light Function			
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn & Push)			
Inner Door		with 4ea of Inner Tempered Glass Door			
Shelf		Included 8ea of PE-coated steel Shelf (2ea×4 Chambers)			
Permitted Total Load		Total Load 128 kg (16kg×8 Shelves)			
Safety Circuit		Over Temp. and Over Current Protector, Sensor Error Detector			
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel			
Dimension (w×d×h)	Interior (mm)	500×500×h500 mm (×4 Rooms)			
	Exterior (mm)	1,655×845×h1,710 mm			
Net Weight / Gross Weight		343kg / 441kg	350kg / 450kg		
Packing Size (w×d×h) (mm)		1,900×1,050×1,950 mm			
Power Consumption		2.6 kW	4.5 kW	4.8 kW	
Power Supply*		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz	1 Phase AC 120V, 60Hz or AC 230V,50/60Hz or 3 Phase		

Cat. No.	230V	DH.WIM04125	DH.WIM14125	DH.WIM24125
	120V	DH.WIM05125	DH.WIM15125	DH.WIM25125
	3Phase	-	DH.WIM16125	DH.WIM26125
Price, Each. ₩				

※ 운임 별도

Multi-Room Incubator with Validation (IQ, OQ) Service

DH.WIM54125 Validation Service(IQ, OQ), with “WIM-4”
DH.WIM64125 , with “WIM-R4”
DH.WIM74125 , with “WIM-RL4”

Ea. Price, ₩

Spare PE-Coated Steel Wire Shelf


DH.WIM34125 Wire Shelf, “IMS125”, PE-coated steel, with 2 Brackets,
w470×d480mm, for “WIM” 4-Room Incubator

Ea. Price, ₩

Optional Stainless-steel Wire Shelf

DH.WIM34126 Wire Shelf, “IMS126”, Stainless-steel, with 2 Brackets,
w470×d480mm, for “WIM” 4-Room Incubator

Ea. Price, ₩

Wire Shelf, PE-Coated steel
(Included)Wire Shelf, Stainless-steel
(Optional)

DAIHAN™ -20℃ + 80℃ SMART Temp/Humidity Chamber “STH-E”, with Plain Door, 155 · 305 · 420 · 800 Lit
With Smart-Lab™ System, Exterior HDPE Water Bottle, Touch-Screen LCD, CFC-Free, Up to 95% RH, 스마트 향온함습기, 경제형



3-Year Warranty

SMART Temp./Humidity Chamber “STH-E”



Installation of Included HDPE Water Bottle

WiRe™ App Service



7" Full Touch Screen TFT LCD
(Smart-Lab™ Controller)

Controlled by Mobile-Phone

It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.



WiRe™ App Service



(1) Perforated
Stainless-steel Shelf
(Included)



Optional ;
(2) Temp. & Humidity Recorder
6 Channel 2 Pens (WTH5100)

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Economical & Practical Type
- CFC Free Refrigerant (R-404A)
- Smart-Lab™ Controller ; Remote Controlled by Mobile-Phone with WiRe™ App Service
- Ergonomic 7" Full Touch Screen TFT LCD
- Internet Connectivity with WiFi
- RS232 Interface for Remote Monitoring and Controlling with PC
- Built in USB port for Data Download and Program Upgrade
- Password Protected Setting
- Variable Program Setting for Multiple Patterns with Multiple Segments
- Self-Diagnostic Function (Push Alarm to Smart Phone) ; Door Ajar and Power-/Temp.-Failure Alarm
- Min-/Max-Temperature Memory and Automatic Data Recording
- Interactive Temperature Graph
- Digital Calibration (Offset Function)
- * The Powerful and Silent Dual Fan Motor Assures Temperature Uniformity by Strong Air Circulation
- * Over Temperature and Over Current Protection for User Safety
- * Inner Tempered Glass Door with Silicone Packing
- * Low Water Level Alarm Lamp
- * Temperature & Humidity Recorder (6 Channels & 2Pens) (Optional)

Specifications & Ordering information

<L1>

Spec.	Model	"STH-E155" Max. Capa. 155L	"STH-E305" Max. Capa. 305L	"STH-E420" Max. Capa. 420L	"STH-E800" Max. Capa. 800L
Capacity		155 Lit	305 Lit	420 Lit	800 Lit
HDPE Water Bottle Capacity		10 Lit. (od. ϕ 11mm Nozzle)			
Temp.	Range	PT100 Sensor, -20℃ + 80℃			
	Accuracy	± 0.3℃			
Humidity	Range	30% ~ 95% RH			
	Accuracy	± 2%			
Dimension (w × d × h)	Interior (mm)	505 × 505 × 610	620 × 620 × 800	750 × 660 × 850	1000 × 800 × 1000
	Exterior (mm)	745 × 920 × 1645	820 × 1290 × 1720	990 × 1075 × 1885	1200 × 1100 × 1920
Heater	for Temp.	1.5 kW	2.0 kW	2.2 kW	2.5 kW
	for Humi.	2.0 kW	3.0 kW	3.0 kW	3.0 kW
Compressor		7/8HP	1.5 HP	1.5 HP	3 HP
Refrigerant		CFC-Free(R-404A) Refrigeration System			
Controller & Display		Smart-Lab Controller, 7" Full Touch Screen TFT LCD			
Program Function		Multiple Programs(Patterns by multiple Segments) Supported			
Material		Internal ; Stainless-steel(# 304), External ; Powder Coated steel			
Shelf (included)		3ea of Perforated Stainless-steel Shelf of Adjustable Position			
Safety Device		Hi-Low Temp. Protector, Over Current & Leakage Breaker, Water Level Sensor			
Others		PL(Product Liability) Insurance			
Packing Dimension (w × d × h) (cm)		99 × 112 × 185	122 × 122 × 205	128 × 128 × 216	149 × 132 × 220
Net & Gross Weight, kg		290 kg, 355 kg	309 kg, 369 kg	378 kg, 478 kg	558 kg, 678 kg
Power Consumption		4.4 kW	5.7 kW	6.3 kW	8.2 kW
Electric Reqts.		1 Phase AC 120V, 60Hz or AC 230V, 50/60Hz or 3 Phase			
Cat. No.	230V	DH.STHE0155	DH.STHE0305	DH.STHE0420	DH.STHE0800
	120V	DH.STHE4155	DH.STHE4305	DH.STHE4420	DH.STHE4800
	3Phase	DH.STHE5155	DH.STHE5305	DH.STHE5420	DH.STHE5800
Price, Each. ₩					

※ 운임 별도

Validation (IQ, OQ) Service Economic Temp./Humidity Chamber

	Ea. Price, ₩
DH.STHE3155 Validation Service(IQ, OQ), with "STH-E155", "STHEV155"	
DH.STHE3305 , with "STH-E305", "STHEV305"	
DH.STHE3420 , with "STH-E420", "STHEV420"	
DH.STHE3800 , with "STH-E800", "STHEV800"	

Spare Perforated Stainless-steel Shelf

	Ea. Price, ₩
DH.WTH11155 (1) Perforated Shelf, "THS1155", for 155 Lit, 470 × 480mm	
DH.WTH11305 , "THS1305", for 305 Lit, 590 × 580mm	
DH.WTH11420 , "THS1420", for 420 Lit, 720 × 640mm	
DH.WTH11800 , "THS1800", for 800 Lit, 970 × 980mm	

Optional Temp & Humidity Recorder and Chart Paper

	Ea. Price, ₩
DH.WTH05100 (2) Temperature & Humidity Recorder, "WTH5100"	
DH.WTH05105 Chart Paper, "WTH5105", 33 Days	

DAIHAN™ -20℃ + 100℃ SMART Temp/Humidity Chamber “STH”, with Viewing Window, 155 · 305 · 420 · 800 Lit
With Smart-Lab™, Auto Supplement Internal Water Tank, Touch-Screen LCD, CFC-Free, Up to 98% RH, 스마트 항온항습기, 투시창 도어

3-Year Warranty

Optional
Chart Recorder

SMART Precise Temp. / Humidity Chamber “STH”

- Supplied by DAIHAN Korea
- CFC Free Refrigerant (R-404A)
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Smart-Lab™ Controller ; Remote Controlled by Mobile-Phone with WiRe™ App Service.

- Ergonomic 7" Full Touch Screen TFT LCD
- Self-Diagnostic Function
- Push Alarm Service to Smart Phone ; Power-/Temperature-/Door Open-Error Status
- Built in USB port for Data Download and Program Upgrade
- Password Protected Setting
- Digital Calibration(Offset Function)
- Program Function ; Multiple Programs(Patterns by multiple Segments) Supported
- Internet Connectivity with WiFi
- Automatic Data Recording
- Interactive Temperature Graph
- Min-/Max-Temperature Memory
- * Precise Temperature Control (-20℃ to 100℃) and Humidity (30% to 98% Relative Humidity)
- * Electronic Humidity Sensor has Excellent Response Time and Durability.
- * Easy and Safety Door Lock
- * Low Water Level Alarm System
- * The Powerful and Silent Dual Fan Motor Assures Temperature Uniformity by Strong Air Circulation (2 ~ 4pcs)
- * Large Clear Viewing Window with Membrane Heating Element (Moisture-Proof)
- * Sliding Water Supply Tank ; Auto Water Supplementing System

Specifications & Ordering information

<L1>

Spec.	Model	“STH-155” Max. Capa. 155L	“STH-305” Max. Capa. 305L	“STH-420” Max. Capa. 420L	“STH-800” Max. Capa. 800L
Temp. Range & Humidity		-20℃ to 100℃, 30% to 98% RH			
Capacity		155 Lit	305 Lit	420 Lit	800 Lit
Water Tank Capacity		21 Lit			
Compressor		7/8HP	1.5HP	1.5HP	3HP
Power Consumption		4.7 kW	5.7 kW	6.3 kW	9 kW
Heater	for Temp.	1.5 kW	2.0 kW	2.2 kW	2.5 kW
	for Humi.	2.0 kW	3.0 kW	3.0 kW	3.0 kW
Temp. Accuracy		±0.3℃			
Humidity Accuracy		±0.3% RH			
Heat-up Time		35 min (-20℃ to 80℃)			
Cool-down Time		45 min (+20℃ to -20℃)			
Sensor		PT100 & Capacitive Sensor			
Material	Internal	Stainless-steel (# 304)			
	External	Powder Coated steel			
Shelf (Included)		3 Perforated Stainless-steel Shelf, Adjustable Position			
Safety Device		Over Temp. and Over Current Protector, Sensor Error and Low Level Detector, Leakage Breaker			
Dimension (w×d×h)	Interior (mm)	505×505×610	620×620×800	700×700×820	1,000×800×1,000
	Exterior (mm)	1,100×960×1,550	1,310×1,080×1,765	1,390×1,135×1,690	1,570×1,245×2,060
Packing Dimension (w×d×h) (cm)		134×116×175	152×125×207	160×130×202	178×141×237
Gross Weight		353 kg	369 kg	478 kg	678 kg
Electric Reqts.		1 Phase AC 120V, 60Hz or AC 230V, 50/60Hz or 3 Phase			

Cat. No.	230V	DH.STH01155	DH.STH01305	DH.STH01420	DH.STH01800
	120V	DH.STH05155	DH.STH05305	DH.STH05420	DH.STH05800
	3Phase	DH.STH07155	DH.STH07305	DH.STH07420	DH.STH07800
Price, Each. ₩					
Validation (IQ, OQ)	Cat. No.	DH.STH31155	DH.STH31305	DH.STH31420	DH.STH31800
	Model No.	STHV1155	STHV1305	STHV1420	STHV1800
Price, Ea. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

WiRe™ App Service

7" Full Touch Screen TFT LCD
(Smart-Lab™ Controller)

It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.



WiRe™ App Service

Controlled by Mobile-Phone



Auto Water Supplementing System

(1) Perforated
Stainless-steel Shelf
(Included)Optional ;
(2) Temp. & Humidity Recorder
6 Channel 2 Pens (WTH5100)

Spare Perforated Stainless-steel Shelf

DH.WTH11155 (1) Perforated Shelf, “THS1155”, for 155 Lit, 468.5×465mm
 DH.WTH11305 “THS1305”, for 305 Lit, 588.5×585mm
 DH.WTH11420 “THS1510”, for 420 Lit, 668.5×675mm
 DH.WTH11800 “THS1800”, for 800 Lit, 968.5×770mm

Ea. Price, ₩

Optional Temp & Humidity Recorder and Chart Paper

DH.WTH05100 (2) Temperature & Humidity Recorder, “WTH5100”
 DH.WTH05105 Chart Paper, “WTH5105”, 33 Days

Ea. Price, ₩

DAIHAN™ -40℃ + 100℃ Precise Temp/Humidity Chamber “THC”, 30~98% RH, 155 · 400 · 800 Lit
With Full Touch Screen Controller, Auto Supplement Water Tank, Viewing Window, Cable Port, 정밀 향온습기

3-Year Warranty

-40℃ + 100℃



Optional
Chart Recorder

Precise Temp / Humidity Chamber “THC-800”



Auto Water Supplementing System



(1) Perforated
Stainless-steel Shelf
(Included)



Optional ;
(2) Temp. & Humidity Recorder

- Supplied by DAIHAN Korea
- 5.7" Program Full Touch Screen Controller
- Built-in SD Card Slot for Data Download & Built-in Cable Port for Convenient Test
- * Precise Temp Control(-40℃ to 100℃) and Humidity(30% to 98% Relative Humidity)
- * Large Clear Viewing Window with Membrane Heating Element (Moisture-Proof)
- * No Temp Change in the Chamber by Separating the Humidity Heater from the Chamber
- * Automatic Defrost Function
- * Easy and Safety Door Lock
- * Sliding Water Supply Tank ; Auto Water Supplementing System
- * Water Recycling System reduce Water Consumption for Longer Operation
- * Low Water Level Alarm System
- * Built in Emergency Stop Switch
- * 6-Channel & 2-Pen Temperature & Humidity Recorder (Optional)
- * Warning Signal Lamp (Include)

- PL(Product Liability) Insurance
- Max 120 Pattern / 1200 Segment

Specifications & Ordering information

<L1>

Spec.	Model	“THC-155” Max. Capa. 155L	“THC-400” Max. Capa. 400L	“THC-800” Max. Capa. 800L
Temp. Range & Humidity		-40℃ ~ 100℃, 30 ~ 98% RH		
Capacity		155 Lit	400 Lit	800 Lit
Water Tank Capacity		22 Lit		
Compressor		1HP × 2Ea.	1.5HP × 2Ea.	2HP × 2Ea.
Power Consumption		5.5 kW	8 kW	11 kW
Heater	for Temp.	2 kW	3 kW	4 kW
	for Humi.	1.5 kW	2 kW	3 kW
Temp.	Accuracy	± 0.3℃		
Humidity	Accuracy	± 0.3% RH		
Sensor		PT100 & Capacitive Sensor		
Material	Internal	Stainless-steel (# 304)		
	External	Stainless-steel (# 430) & Powder Coated steel		
Shelf (Included)		3 Perforated Stainless-steel Shelf, Adjustable Position(THC-155 2ea)		
Safety Device		Over Temp. and Over Current Protector, Sensor Error and Low Level Detector, Leakage Breaker		
Dimension (w×d×h)	Interior (mm)	505×505×610	700×700×820	1000×800×1000
	Exterior (mm)	1000×1020×1700	1200×1220×1900	1500×1320×2000
Packing Dimension (w×d×h) (mm)		1210×1220×2000	1410×1420×2200	1710×1520×2300
Gross Weight		360 kg	480 kg	680 kg
Electric Reqts.		AC 380V 50/60Hz, 3 Phase		
Cat. No.	380V	DH.Thc0155	DH.Thc0400	DH.Thc0800
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Validation (IQ, OQ) Service Precise Temp./Humidity Chamber

DH.Thc1155 Validation Service(IQ, OQ), with “THC-155”
DH.Thc1400 , with “THC-400”
DH.Thc1800 , with “THC-800”

Ea. Price, ₩

Spare Perforated Stainless-steel Shelf, Warning Signal Lamp for “THC”

DH.ThcS0155 (1) Perforated Shelf, “THCS155”, for 155 Lit, 475×455mm
DH.ThcS0400 , “THCS400”, for 400 Lit, 670×650mm
DH.ThcS0800 , “THCS800”, for 800 Lit, 970×750mm
DH.ThcL0003 Warning Signal Lamp

Ea. Price, ₩

Optional Accessory

DH.ThcR0005 (2) Temperature & Humidity Recorder
DH.WTH05105 Chart Paper for Recorder

Ea. Price, ₩

DAIHAN™ Low Temperature(B.O.D) Incubator “ThermoStable™ IR”, 0°C ~ 60°C, $\pm 0.2^\circ\text{C}$, Medicaluse(230V)
With 3 Wire Shelf, Digital PID Control, CFC-Free Refrigeration(R-404A) System, Forced-Air Type, 150 · 250 · 420 · 700 Lit
저온(B.O.D) 배양기/인큐베이터, 강제 순환식, 우수한 온도정확성, 친환경 냉매 채용



Medicaluse



- Supplied by DAIHAN Korea
- UL, CUL, and CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- CFC-free (R-404A) Refrigeration System
- Digital PID Control System Implementing Superior Temperature Accuracy
- Patented Jog-Dial Control System

- * Variable Application ; Animal & Plant Cell Culture, Constant Temperature, Germination Test, BOD Test, and so on
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * RS232C Interface for Monitoring and Controlling with PC
- * Superior Temperature Uniformity by Forced-air
- * Over Temperature & Over Current Protection, Sensor Error Detection and Leakage Breaking.
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled)

Specifications & Ordering information

<M1>

Spec.	Model	ThermoStable™ IR-150 Max. Capa. 150L	ThermoStable™ IR-250 Max. Capa. 250L	ThermoStable™ IR-420 Max. Capa. 420L	ThermoStable™ IR-700 Max. Capa. 700L
Capacity		150 Lit	250 Lit	420 Lit	700 Lit
Temp Range & Accuracy		0°C~60°C, $\pm 0.2^\circ\text{C}$ at 37°C			
Temp Uniformity & Sensor		$\pm 0.5^\circ\text{C}$ at 37°C, $\pm 1.4^\circ\text{C}$ at 50°C, PT100 Sensor			
Heat-up Time		20 min. to 37°C 20 min. to 50°C		35 min. to 37°C 35 min. to 50°C	
Recovery Time (Door open 30sec)		10 min. to 37°C 12 min. to 50°C		12 min. to 37°C 15 min. to 50°C	
Heating Power/Consumption		700 W	1 kW	1.2 kW	2 kW
Compressor		1/4HP	1/3HP	1/2HP	1/2HP
Timer & Alarm		99hr 59min (with delay / Continuous run.), Error status & Timer-end			
Storage		Alarm & Timer			
Display		Digital LCD with Back Light			
Controller		Digital PID Controller with Jog-Shuttle Switch (Turn & Push)			
Inner Door		with 1 Inner Tempered Glass Door			
Shelf (Included)		3 PE-coated steel Shelf			
Permitted Total Load		Total Load ; 48kg (16kg x 3 Shelf)			
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel			
Dimension (w × d × h)	Interior (mm)	505 × 505 × 610	505 × 505 × 1,000	705 × 605 × 1,000	800 × 700 × 1,250
	Exterior (mm)	640 × 832 × 1,313	640 × 832 × 1,703	840 × 932 × 1,753	935 × 1,027 × 2,000
Net Weight		118 kg	157 kg	202 kg	234 kg
Packing Size (w × d × h) (mm)		900 × 1,000 × 1,650	900 × 1,000 × 2,040	1,100 × 1,100 × 2,140	1,220 × 1,220 × 2,260
& Gross Weight		191 kg	217 kg	277 kg	313 kg
Power Consumption		1.6 kW	1.8 kW	2.3 kW	2.4 kW
Electric Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz			



WiseRemote Control Software ;

Connection of Temperature Devices
with Interface RS232C to a PC

- * You can Download this Software from
DAIHAN-Sci.com Shopping-mail

Cat. No.	230V	DH.WIR03150	DH.WIR03250	DH.WIR03420	DH.WIR03700
	120V	DH.WIR04150	DH.WIR04250	DH.WIR04420	DH.WIR04700
Price, Each. ₩					
Validation (IQ, OQ)	Cat. No.	DH.WIR31150	DH.WIR31250	DH.WIR31420	DH.WIR31700
	Model No.	IRV150	IRV250	IRV420	IRV700
Price, Ea. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

(1) Spare PE-Coated Steel Wire Shelf

Ea. Price, ₩

- DH.WIR11150 Wire Shelf, PE Coated steel, “IRS150”, w470 × d480mm,
for “IR-150”, “IR-250”
DH.WIR11420 Wire Shelf, PE Coated steel, “IRS420”, w675 × d585mm, for “IR-420”
DH.WIR11700 “IRS700”, w770 × d680mm, for “IR-700”

(2) Optional Stainless-steel Wire Shelf

Ea. Price, ₩

- DH.WIR21150 Wire Shelf, Stainless-steel, “IRU150”, w470 × d480mm,
for “IR-150”, “IR-250”
DH.WIR21420 Wire Shelf, Stainless-steel, “IRU420”, w675 × d585mm, for “IR-420”
DH.WIR21700 “IRU700”, w770 × d680mm, for “IR-700”



(1) Wire Shelf,
PE Coated steel (Included)



(2) Wire Shelf,
Stainless-steel (Optional)

DAIHAN™ **SMART** Low Temperature(B.O.D) Incubator “**ThermoStable™ SIR**”, 150 · 250 · 420 · 700 Lit, **Medicaluse(230V)**
With **Smart-Lab™** Controller, 4" Full Touch Screen, Fuzzy-PID Control, CFC-Free(R-404A), 0℃ ~ 60℃, ± 0.2℃, with Certi. & Traceability
스마트 저온(B.O.D) 배양기/인큐베이터, 강제 순환식, 우수한 온도정확성, 친환경 냉매 채용, WiRe™ 서비스, 4인치 풀터치스크린 조절기



- Supplied by DAIHAN Korea
- UL, CUL, and CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Air Flow ; Forced Air-type
- CFC-free (R-404A) Refrigeration System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Smart-Lab™ Controller ; Remote Controlled by Mobile-Phone with WiRe™ App Service
 - Ergonomic 4" Full Touch Screen TFT LCD
 - Internet Connectivity with WiFi
 - RS232 Interface for Remote Monitoring and Controlling with PC
 - Automatic Data Recording
 - Built in USB port for Data Download and Program Upgrade
 - Self-Diagnostic Function
 - Variable Program Setting for Multiple Patterns with Multiple Segments
- * Variable Application ; Animal/Plant Cell Culture, Constant Temp Germination Test, BOD Test, etc.
- * Stainless-steel Chamber and PE-coated Shelves
- * Compact New Body Design
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * Over Temperature & Over Current Protection, Sensor Error Detection and Leakage Breaking.

Specifications & Ordering information

<M1>

Spec.	Model	ThermoStable™ SIR-150 Max. Capa. 150L	ThermoStable™ SIR-250 Max. Capa. 250L	ThermoStable™ SIR-420 Max. Capa. 420L	ThermoStable™ SIR-700 Max. Capa. 700L
Capacity		150 Lit	250 Lit	420 Lit	700 Lit
Temp Range & Accuracy		0℃~60℃, ± 0.2℃ at 37℃			
Temp Uniformity & Sensor		± 0.5℃ at 37℃, ± 1.4℃ at 50℃, PT100 Sensor			
Heat-up Time		20 min. to 37℃ 20 min. to 50℃		35 min. to 37℃ 35 min. to 50℃	
Recovery Time (Door open 30sec)		10 min. to 37℃ 12 min. to 50℃		12 min. to 37℃ 15 min. to 50℃	
Heating Power/Consumption		700 W	1 kW	1.2 kW	2 kW
Compressor		1/4HP	1/3HP	1/2HP	1/2HP
Compressor Control		Intelligent Control of Refrigerator for preventing Overload of Compressor			
Refrigerant		CFC-Free(R-404A) Refrigeration System			
Timer & Alarm		99hr 59min (with delay / Continuous run.), Error status & Timer-end			
Display & Controller		4" Full Touch Screen TFT LCD, Smart-Lab™ Controller			
Inner Door		with 1ea of Inner Tempered Glass Door			
Shelf		Included 3ea of PE-coated steel Shelf			
Permitted		Total Load Total Load ; 48 kg (16kg × 3 Shelf)			
Safety Circuit		Over Temp and Over Current Protector, Sensor Error Detector			
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel			
Dimension (w × d × h)	Interior (mm)	505 × 505 × 610	505 × 505 × 1,000	705 × 605 × 1,000	800 × 700 × 1,250
	Exterior (mm)	640 × 832 × 1,313	640 × 832 × 1,703	840 × 932 × 1,753	935 × 1,027 × 2,000
Net Weight / Gross Weight		118 kg / 191 kg	157 kg / 217 kg	202 kg / 277 kg	234 kg / 313 kg
Packing Size (w × d × h) (mm)		900 × 1,000 × 1,650	900 × 1,000 × 2,040	1,100 × 1,100 × 2,140	1,220 × 1,220 × 2,260
Power Consumption		1.6 kW	1.8 kW	2.3 kW	2.4 kW
Electric Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz			

Cat. No.	230V	DH.SWIR03150	DH.SWIR03250	DH.SWIR03420	DH.SWIR03700
	120V	DH.SWIR04150	DH.SWIR04250	DH.SWIR04420	DH.SWIR04700

Price, Each. ₩					
----------------	--	--	--	--	--

Validation (IQ, OQ)	Cat. No.	DH.SWIR31150	DH.SWIR31250	DH.SWIR31420	DH.SWIR31700
	Model No.	SIR-150	SIR-250	SIR-420	SIR-700

Price, Ea. ₩					
--------------	--	--	--	--	--

* Other Specifications are available upon Customer's Request.

※ 윤임 별도

(1) Spare PE-Coated Steel Wire Shelf

DH.WIR11150 Wire Shelf, "IRS150", 470 × 480mm, for "SIR-150", "SIR-250"
 DH.WIR11420 "IRS420", 675 × 585mm, for "SIR-420"
 DH.WIR11700 "IRS700", 770 × 680mm, for "SIR-700"

Ea. Price, ₩

(2) Optional Stainless-steel Wire Shelf

DH.WIR21150 Wire Shelf, "IRU150", 470 × 480mm, for "SIR-150", "SIR-250"
 DH.WIR21420 "IRU420", 675 × 585mm, for "SIR-420"
 DH.WIR21700 "IRU700", 770 × 680mm, for "SIR-700"

Ea. Price, ₩

WiRe™ App Service



4" Full Touch Screen TFT LCD
(Smart-Lab™ Controller)

Controlled by Mobile-Phone


It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.



WiRe™ App Service



(1) Wire Shelf
PE Coated steel
(Included)



(2) Wire Shelf
Stainless-steel
(Optional)

DAIHAN™ SMART Plant Growth Chamber "GC", Smart-Lab™ System, 432 · 864 Lit

With Smart-Lab™ Controller, CFC-Free(R-404A) Refrigeration System, 0~12,000 - or 15,000-Lux, 30~95% RH, 10℃~60℃

스마트 식물 성장상, 정밀한 온도/습도/조도 조절, WiRe™ 서비스, 7인치 풀터치스크린 스마트랩 조절기



3-Year Warranty

SMART Plant Growth Chamber "GC"

WiRe™ App Service

7" Full Touch Screen TFT LCD
(Smart-Lab™ Controller)

Controlled by Mobile-Phone

It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.



WiRe™ App Service



Water Tank, Silicone Tube and Lux Meter (Included)



Wire Shelf, PE-Coated Steel (Included)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Suitable for Plant Growth, Stability, Product Self-life QC Testing
- Ergonomic Screen Design 7" Full Touch Screen TFT LCD
- Smart-Lab™ Controller ; Remote Controlled by Mobile-Phone with WiRe™ App Service.

- Self-Diagnostic Function - Automatic Data Recording
- Variable Program Setting for Multiple Patterns with Multiple Segments
- The Recorded Data could be transferred to PC by USB memory
- Controller SW can be upgraded through USB port.
- Precise Temperature, Relative Humidity and Illumination Control
- Maintain Uniformity and Accuracy of Growth Circumstance

* Inner Tempered Glass Door with Silicone Packing offers easy Observation without Door Opening

* Over Temperature & Over Current Protection, Sensor Error Detection, Leakage Breaking

* Includes Water Tank(with Silicone Tube) and Lux Meter

Specifications & Ordering information

<L1>

Spec.		Model	(1) "GC-450" Max. Capa. 432 Lit	(2) "GC-1000" Max. Capa. 864 Lit
Temp.	Range & Variation		10℃~60℃, ±1℃ at 20℃	
	Fluctuation		± 0.5℃ at 10℃, ± 0.5℃ at 25℃ (0% light intensity) ± 1℃ at 10℃, ± 1℃ at 25℃ (100% light intensity)	
Humidity	Range & Variation		30% ~ 95% RH, ±3% at 60℃	
	Fluctuation		± 1%	
	Humidifier		Evaporation by Fan with Heater	
Illumination	Range		0 to 12,000 Lux	0 to 15,000 Lux
	Lamp		FPL 55W × 14	FPL 55W × 28
	Bank		3 Side Walls	
Temp. & Humidity Sensor			Temp. ; PT100, Humidity ; Capacitive Sensor	
Heat-up Time (up to 37℃)			11min. (0% illumination) / 9min. (100% illumination)	
Cool-down Time (down to 10℃)			13min. (0% illumination) / 15min. (100% illumination)	
Temp. Recovery (After Door was opened for 30sec)			3min. at 10℃, 1min. at 37℃ (0% illumination) 5min. at 10℃, 1min. at 37℃ (100% illumination)	
Heating Power /Consumption	for Temp.		1.5 kW	3.0 kW
	for Humi.		1.5 kW	3.0 kW
Compressor			1/2HP	1HP
Program Function			Temperature and Humidity ; 120 Patterns, 1,200 Segments Light Bank ; 7 steps	
Door	Internal		Tempered Safety Glass with Silicone Packing	
	External		3-Side Magnetic Packing Door with Light Bank	
Shelf (Included)			4 PE-coated Steel Shelf	8 PE-coated Steel Shelf
Permitted Total Load			64kg (16kg × 4 Shelves)	128kg (16kg × 8 Shelves)
Material			Internal ; Stainless-steel (# 304), External ; Power Coated steel	
Dimension (w×d×h)	Interior (mm)		650×700×1,050	1,200×790×1,050
	Exterior (mm)		850×1,090×1,965	1,400×1,150×1,965
Packing Size & Gross Weight			1,400×1,220×h2,210mm, 435kg	1,600×1,500×h2,200mm, 585kg
Power Consumption			4.8 kW	6.5 kW
Power Supply			AC 230V, 50/60Hz or 3 Phase	
Cat. No.	230V		DH.SWGC00450	DH.SWGC01000
	3Phase		DH.SWGC60450	DH.SWGC61000
Price, Each. ₩				
Validation (IQ, OQ)	Cat. No.		DH.SWGC30450	DH.SWGC31000
	Model No.		SGCV450	SGCV1000
Price, Ea. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Wire Shelf, PE Coated steel

Ea. Price, ₩

DH.WGC10450 Wire Shelf, PE Coated, w640×d620mm, "GCS450", for 432Lit
 DH.WGC11000 , w580×d620mm, "GCS1000", for 864Lit

DAIHAN™ Precise Shaking Incubator “IS-20”, BenchTop-type, Up to 60°C & 1000ml
With Universal Platform, Fuzzy Control, RS232C Interface for PC, Orbital Motion, 30~250 rpm, 101 Lit
진탕 배양기/인큐베이터, 벤치 탑 타입, 고정밀 디지털 퍼지 제어, Universal Platform 포함



A. “IS-20” R/T +5°C~60°C
B. “IS-20R” 10°C~60°C Refrigerated



WiseRemote Control Software ;
Download from DAIHAN-Sci.com

Connection of
Temperature
Devices with
Interface RS232C
to a PC



(1) Universal Platform “SP600” (Included) (2) Universal Spring Rack Platform “SR500” (Optional)



(3) STS Flask Holder w/Spring Ring (Optional)
(4) “FHS250” (5) “FHS100” (6) “FHS1000”

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Useful for Biological Cultures under various Temperature
- Patented Jog-Shuttle Control System ■ Stabilized Orbital Motion
- * Precise Temperature Control by Digital Fuzzy Controller, up to 60°C, PT100 Sensor
- * Advanced Shaking Mechanism Provides Quiet Shaking and Precise Speed Control
- * Automatic Stop of Shaking, Heating and Refrigeration when a Door opens
- * Over Temperature & Over Current Protection, Sensor Error Detection and Leakage Breaking
- * Accommodate Flasks up to 1,000ml * Minimum Noise and No Vibration
- * RS232C Interface for Monitoring and Controlling with PC
- * Universal Spring Rack for Various Shapes/Sizes of Flasks (Optional)

Specifications & Ordering information

<L1>

Spec.	Model	A. “IS-20”	B. “IS-20R” with Refrigerator
Temp Range		Room Temp +5°C~60°C	10°C~60°C
Temp Accuracy & Uniformity		Accuracy ; ±0.2°C at 37°C, Uniformity ; ±0.5°C at 37°C	
Capacity		101Lit	
Temp Sensor		PT100 Sensor	
Heater		500 W	700 W
Compressor		None	1/4 HP
Compressor Control		None	Intelligent Control of Refrigerator for preventing Overload of Compressor
Refrigerant		None	CFC-Free(404A) Refrigeration System
Shaking Speed & Stroke		30 ~ 250 rpm, 25mm, Orbital Motion	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)	
Room Lamp		60W Krypton Lamp	
Platform(Included) & Capacity for Optional Flask Holders		Universal Platform for Various Shapes & Sizes of Flasks Capacity ; 32ea Flask Holders for 50~100ml Flask 20ea Flask Holders for 200~300ml Flask 16ea Flask Holders for 500~1,000ml Flask	
Safety Circuit		Over Temp/ Current Protector, Leakage Breaker, Sensor Error Detector	
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel	
Dimension		Internal ; 530×530×h350mm, External ; 575×855×h850mm	
Packing Size (w×d×h) (mm)		1,000×800×h1,100mm	
& Gross Weight		147 kg	150 kg
Power Consumption		500 W	1.2 kW
Electrical Reqt.s.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	

Cat. No.	60Hz	230V	DH.WIS00220	DH.WIS01220
		120V	DH.WIS00120	DH.WIS01120
Cat. No.	50Hz	230V	DH.WIS00520	DH.WIS01520

Price, Each. ₩

* Other Specifications are available upon Customer's Request.

※ 운임 별도

C. Validation (IQ, OQ) Service Shaking Incubator

DH.WIS32010 Validation Service(IQ, OQ), “ISV2010”, with “IS-20”, 60Hz
DH.WIS32011 “ISV2011”, with “IS-20R”, 60Hz

Ea. Price, ₩

Optional Spring Rack Platform and Stainless-steel Flask Holder(with Spring Ring)

- DH.WIS504600 (1) Universal Platform, “SP600”, 478×478mm, for “IS-20(R)”, Included
- DH.WIS504500 (2) Universal Spring Rack, “SR500”, 483×483×h135mm, for “IS-20(R)”
- DH.WIS600100 (3) Stainless-steel Wire Spring Flask Holder, “FHS100”, for 100ml Flask
- DH.WIS600250 (4) “FHS250”, for 250ml Flask
- DH.WIS600500 (5) “FHS500”, for 500ml Flask
- DH.WIS601000 (6) “FHS1000”, for 1000ml Flask

Ea. Price, ₩

DAIHAN™ Precise Shaking Incubator “IS-30”, Front Door-type, Up to 60°C & 1000ml

With Universal Platform, Fuzzy Control, RS232C Interface for PC, Orbital Motion, 30~250 rpm, 127 Lit

진탕 배양기/인큐베이터, 전면도어 타입, 고정밀 디지털 피치제어, Universal Platform 포함



Up to 60°C/250rpm

Patented Jog-Shuttle Control System

Front Door

Back Light LCD



With Viewing Window and Universal Platform

3-Year Warranty

A. “IS-30” R/T +5°C~60°C

B. “IS-30R” 10°C~60°C Refrigerated

■ Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance

■ with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

■ Useful for Biological Cultures under various Temperature

■ Patented Jog-Shuttle Control System ■ Stabilized Orbital Motion

* Precise Temperature Control by Digital Fuzzy Controller, up to 60°C, PT100 Sensor

* Front Door and Sliding Platform Provides Easy Loading and Handling of Samples

* RS232C Interface for Monitoring and Controlling with PC.

* Advanced Shaking Mechanism Provides Quiet Shaking and Precise Speed Control

* Accommodate Flasks up to 1,000ml * Minimum Noise and No Vibration.

* Universal Spring Rack for Various Shapes/Size of Flasks (Optional)

* Automatic Stop of Shaking, Heating and Refrigeration when a Door opens

* Over Temperature & Over Current Protection and Leakage Breaking.

Specifications & Ordering information

<L1>

Spec.	Model	A. “IS-30”	B. “IS-30R” with Refrigerator
Temp Range		Room Temp +5°C~60°C	10°C~60°C
Temp Accuracy & Uniformity		Accuracy ; ±0.2°C at 37°C, Uniformity ; ±0.5°C at 37°C	
Capacity		127Lit	
Temp Sensor		PT100 Sensor	
Heater		600 W	800 W
Compressor		None	1/4 HP
Compressor Control		None	Intelligent Control of Refrigerator for preventing Overload of Compressor
Refrigerant		None	CFC-Free(404A) Refrigeration System
Shaking Speed, Stroke & Motion		30 ~ 250 rpm, 25mm, Orbital Motion	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)	
Platform(Included) & Capacity for Optional Flask Holders		Universal Platform for Various Shapes & Sizes of Flasks Capacity ; 35ea Flask Holders for 50~100ml Flask 16ea Flask Holders for 200~300ml Flask 12ea Flask Holders for 500~1,000ml Flask	
Safety Circuit		Over Temp/ Current Protector, Leakage Breaker	
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel	
Dimension		Internal ; 570 × 465 × h470mm, External ; 916 × 730 × h683mm	
Packing Size (w × d × h) (mm)		1,050 × 900 × h900mm	
& Gross Weight		146 kg	156 kg
Power Consumption		600 W	1.3 kW
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	

Cat. No.	60Hz	230V	DH.WIS00230	DH.WIS01230
		120V	DH.WIS00130	DH.WIS01130
	50Hz	230V	DH.WIS00530	DH.WIS01530

Price, Each. ₩

* Other Specifications (50Hz) are available upon Customer's Request.

※ 운임 별도



WiseRemote Control Software ; Download from DAIHAN-Sci.com

Connection of Temperature Devices with Interface RS232C to a PC



(1) Universal Platform “SP630” (Included)



(2) Universal Spring Rack Platform “SR530” (Optional)



Optional Spring Rack Platform and Stainless-steel Flask Holder(with Spring Ring)

- DH.WIS504630 (1) Universal Platform, “SP630”, 485 × 405mm, for “IS-30(R)”
 DH.WIS504530 (2) Universal Spring Rack, “SR530”, 498 × 408 × h135mm, for “IS-30(R)”
 DH.WIS600100 (3) Stainless-steel Wire Spring Flask Holder, “FHS100”, for 100ml Flask
 DH.WIS600250 (4) “FHS250”, for 250ml Flask
 DH.WIS600500 (5) “FHS500”, for 500ml Flask
 DH.WIS601000 (6) “FHS1000”, for 1000ml Flask

C. Validation (IQ, OQ) Service Shaking Incubator

DH.WIS33010 Validation Service(IQ, OQ), “ISV3010”, with “IS-30”, 60Hz
 DH.WIS33011 “ISV3011”, with “IS-30R”, 60Hz

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ Precise Shaking Incubator “IS-10”, Top Door-type, Orbital Motion, Up to 60°C, $\pm 0.2^\circ\text{C}$
With Universal Platform, Fuzzy Control, with or without illuminators & Recorder, 30~250 rpm, with Certi. & Traceability
진탕 배양기/인큐베이터, 탑 도어 타입, 고정밀 디지털 퍼지 제어, Universal Platform 포함



R/T +5°C~60°C/250rpm

Top Door System

Jog-Dial
with Push Button

3-Year Warranty

A. Room Temp +5°C~60°C, with Universal Platform
“IS-10”



WiseRemote Control Software ;
Download from DAIHAN-Sci.com

Connection of
Temperature
Devices with
Interface RS232C
to a PC



(1) Universal Platform
“SP610” (Included)

(6) Universal Spring Rack
Platform “SR510” (Optional)



STS Flask Holder w/Spring Ring
(Optional)

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Useful for Biological Cultures under various Temperature
- Patented Jog-Shuttle Control System

- * Precise Temperature Control by Digital Fuzzy Controller, up to 60°C
- * Digital LCD with Back-Light Function
- * Minimum Noise and No Vibration
- * Advanced Shaking Mechanism Provides Quiet Shaking and Precise Speed Control
- * Automatic Stop of Shaking, Heating and Refrigeration when a Door opens
- * Intelligent Control of Refrigerator for preventing Overload of Compressor (IS-10R / IS-10RL)
- * RS232C Interface for Monitoring and Controlling with PC

Specifications & Ordering information

<L1>

Spec.	Model	A. “IS-10”	B. “IS-10R” with Refrigerator	C. “IS-10RL” Illumination
Temp. Range		R/T +5℃~60℃	10℃~60℃	
Temp. Accuracy & Uniformity		Accuracy ; ±0.2℃ at 37℃, Uniformity ; ±0.5℃ at 37℃		
Capacity		229Lit		
Temp. Sensor		PT100 Sensor		
Heater		800 W	1 kW	1.2 kW
Compressor		None	1/3HP	
Refrigerant		CFC-Free (404A) Refrigeration System		
Shaking Speed & Stroke		30 ~ 250 rpm, 25mm Orbital Motion		
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end		
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)		
Illumination System		None		FPL 55W×4ea
Room Lamp		Krypton Lamp, 60W		None
Platform(Included)		Universal Platform for Various Shapes & Sizes of Flasks		
Safety Circuit		Over Temp./ Current Protector, Leakage Breaker, Sensor Error Detector		
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel		
Dimension		Internal ; 750×720×h400mm, External ; 1,065×900×h980mm		
Packing Size (w×d×h) (mm)		1,220×1,060×h1,200mm		
& Gross Weight		231 kg	238 kg	238 kg
Power Consumption		800 W	1.8 kW	2.1 kW
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz		

Cat. No.	60Hz	230V	DH.WIS00210	DH.WIS04210	DH.WIS05210
		120V	DH.WIS00110	DH.WIS04110	DH.WIS05110
	50Hz	230V	DH.WIS00510	DH.WIS04510	DH.WIS05510
		120V	DH.WIS00410	DH.WIS04410	DH.WIS05410

Price, Each. ₩

* Other Specifications are available upon Customer's Request.

※ 운임 별도

D. Validation (IQ, OQ) Service Shaking Incubator

DH.WIS31010 Validation Service(IQ, OQ), with “IS-10”, 60Hz
DH.WIS31011 , with “IS-10R”, 60Hz
DH.WIS31012 , with “IS-10RL”, 60Hz

Ea. Price, ₩

Optional Spring Rack Platform/Recorder/Stainless-steel Flask Holder(with Spring Ring)

- DH.WIS504610 (1) Spare Universal Platform, “SP610”, w688×h648mm, Flask Holder Capacity ; 68ea for 50~100ml Flask, 36ea for 200~300ml Flask, 21ea for 500~1000ml Flask
- DH.WIS600100 (2) Stainless-steel Wire Spring Flask Holder, “FHS100”, for 100ml Flask
- DH.WIS600250 (3) , “FHS250”, for 250ml Flask
- DH.WIS600500 (4) , “FHS500”, for 500ml Flask
- DH.WIS601000 (5) , “FHS1000”, for 1000ml Flask
- DH.WIS504510 (6) Universal Spring Rack, “SR510”, w693×d653×h135mm
- DH.WTH05100 (7) Temperature Recorder, “WTH5105”, with 6ea Ribbon Cassette & Chart Paper LED, 6ch/10sec, 144×144×h233mm
- DH.WTH05108 Ribbon Cassette, “WTH5108”, 6 Colors
- DH.WTH05105 Chart Paper, “WTH5105”, w100mm×L10m, 33 Days

Ea. Price, ₩

DAIHAN™ Multi-Stack Refrigerated Shaking Incubator “WIS-ML”, Orbital Motion, with 2 & 4 Universal Platform

With Digital Fuzzy Control, 2 or 4 Independent Shaking Speed Control, 10℃~60℃, 30~250 rpm, 2 & 4 Stack

대형 멀티 진탕 배양기, 고정밀 디지털 퍼지 제어, 2 & 4 Universal Platform 포함

CE

10℃~60℃ / 250rpm



B. Multi-Stack Refrigerated Shaking Incubator with 4 Universal Platform “WIS-ML04”

3-Year Warranty

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Excellent for Multi-use and wide range of Sample Testing under various Temperature
- Patented Jog-Shuttle Control System

- * Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * Advanced Shaking Mechanism provides Quiet Shaking and Precise Speed Control
- * Space Saving Design for Large Quantity Samples with 2 or 4 Independent Shaking Speed Controller
- * Automatic Stop of Shaking, Heating and Refrigeration when the Door is Opened
- * RS232C Interface for Monitoring and Controlling with PC
- * Including Universal Platform
- * Over Temperature & Over Current Protection and Sensor Error Detection
- * Universal Spring Rack for Various Shapes/Sizes of Flasks (Optional)

Specifications & Ordering information

<L1>

Spec.	Model	A. “WIS-ML02” 2 Shaking Platform	B. “WIS-ML04” 4 Shaking Platform
Shaking Stacks		2 Stack	4 Stack
Temp. Range & Accuracy		10℃~60℃, ±0.2℃ at 37℃	
Temp. Uniformity & Sensor		±1.0℃ at 37℃, PT100 Sensor	
Heater		1.5 kW	1.8 kW
Compressor		1/2HP	3/4HP
Compressor Control		Intelligent Control of Refrigerator for Preventing Overload of Compressor	
Refrigerant		CFC-Free(404A) Refrigeration System	
Shaking Speed & Stroke		30 ~ 250 rpm, 25mm Orbital Motion	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)	
Door		Magnetic Packing Door with Pair Glass	
Plate Size		w700×d600mm	
Platform(Included) & Capacity for Optional Flask Holders		Universal Platform for Various Shapes & Sizes of Flasks Capacity ; 68ea Flask Holder, for 50~100ml Flask 36ea Flask Holder, for 200~300ml Flask 21ea Flask Holder, for 500~1,000ml Flask	
Safety Circuit		Over Temp./ Current Protector, Sensor Error Detector	
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel	
Dimension (w×d×h)	Internal (mm)	700×700×1,100	700×700×1,790
	External (mm)	1,290×850×1,445	1,290×850×2,125
Packing Size (w×d×h) (mm)	1,500×1,200×1,750		1,500×1,200×2,350
	397 kg		555 kg
Power Consumption		2.4 kW	3.2 kW
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz	1 Phase AC 120V, 60Hz or AC 230V, 50/60Hz or 3 Phase
Cat. No.	230V	DH.WIS05002	DH.WIS05004
	120V	DH.WIS05102	DH.WIS05104
	3Phase	-	DH.WIS03004
Price, Each. ₩			

※ 운임 별도

C. Validation (IQ, OQ) Service Shaking Incubator

DH.WIS35002 Validation Service(IQ, OQ), with “WIS-ML02”
 DH.WIS35004 , with “WIS-ML04”

Ea. Price, ₩


Optional Spring Rack Platform and Stainless-steel Flask Holder(with Spring Ring)

Ea. Price, ₩



- DH.WIS504620 (1) Universal Platform, “SP620”, w645×h645mm, for “WIS-ML02/04”, Included
 DH.WIS504520 (2) Universal Spring Rack, “SR520”, w653×d653×h135mm, for “WIS-ML02/04”
 DH.WIS600100 (3) Stainless-steel Wire Spring Flask Holder, “FHS100”, for 100ml Flask
 DH.WIS600250 (4) , “FHS250”, for 250ml Flask
 DH.WIS600500 (5) , “FHS500”, for 500ml Flask
 DH.WIS601000 (6) , “FHS1000”, for 1000ml Flask

DAIHAN™ Digital Jar Tester “JT-M6C”, 6-Paddle Stirrer for Flocculation Test, 300 rpm, with Certi. & Traceability
 With Programmable Digital Feedback Control, 6 Impellers, 6 Independent Light Source, Up to 6 × 1 Lit Beakers
 프로그램/디지털 자 테스터, 정밀 직류 기어드 모터와 특화된 타이밍벨트, 디지털 피드백 컨트롤



3-Year Warranty

Jar Tester, $\phi 8 \times L400\text{mm}$ Shaft “JT-M6C”



Geared DC Motor with Special Timing Belt

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Flocculation Testing in Water Treatment Laboratories
- 6 Independent Light Source for Each Stirring Point
- * Deciding on Coagulant Chemistries * Digital Display for RPM, Power Bar-Graph and Status
- * Digital Feedback Control using Microprocessor ; Stable and Accurate Stirring Speed
- * Geared DC Motor and Specially Designed Timing Belt used
- * Stable and Accurate Stirring Speed by Feedback Control
- * Digital LCD with Back-Light Function * Programmable Operation ; 10 Steps, 6 Programs

Specifications & Ordering information

<L1>

Spec.	Model	“JT-M6C”
6 × Paddler (Stainless steel Shaft)		$\phi 8 \times L400\text{mm}$ Straight Blade w75 × h25mm
Stirring Capacity		6 × 1 Liter Beaker
Motor & Power		Geared DC Motor with Timing Belt (10:1 ratio), 40W
Speed & Accuracy		20~300 rpm, ± 5 rpm
Timer		99hr 59min 59sec (Programmable)
Display		Digital LCD with Back-Light for RPM, Power Bar-Graph and Status
Controller		Programmable Digital Feedback Controller with Jog-Shuttle Switch (Turn + Push)
Program Function		10 Steps, 6 Programs, and Storage for Set Values ; RPM and Timer
Florescent Lamp		20W × 6 (Independent Light Source for Each Sample)
Dimension (w × d × h)		1,000 × 310 × h535mm
Safety Device		Overload & Over Heat Protector for motor
Others		PL(Product Liability) Insurance and Storage Function (RPM, Timer)
Packing and Gross Weight		1,120 × 430 × h740mm, 69 kg
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz
Cat. No.	230V	DH.WJT06020
Price, Each. ₩	120V	DH.WJT06021

※ 운임 별도



Fluorescent Light Magnifier "BLM-F"

DAIHAN Wisd™ Multi-Use Fluorescent Light Magnifier "BLM-F", × 5 & × 8 Magnification
Base Stand-type, ϕ 127mm B270 High Clear White Glass Lens, Flexible Arm, 12W, 형광 조명 확대경

- Supplied by DAIHAN Korea
- Ideal for Laboratory, Hospital, Circuit Board Check-Up, etc.
- International Standard ϕ 127mm B270 White Glass Lens
- Light Source ; High Brightness Daylight Blue Fluorescent
- 20% Energy Saving Ballast Stabilizer Control System
- Adjustable Height and Lens Angle Direction ■ 2-Stage Spring Flexible Arm
- * Light Illumination ; 2,800lx * Adjustable Height ; 300 ~ 500mm
- * Base Size ; 220 × 330mm * Electric Power ; AC 220V, 12W

Ordering Information

DH.Mag3105 Magnifier, Fluorescent Light, × 5, "BLM-F105"
DH.Mag3108 , × 8, "BLM-F108"

<L2>

Ea.	Price, ₩
-----	----------



LED Light Magnifier "BLM-L"

DAIHAN Wisd™ Multi-Use LED Light Magnifier "BLM-L", × 5 & × 8 Magnification
Base Stand-type, ϕ 127mm B270 High Clear White Glass Lens, Flexible Arm, 12W, LED 조명 확대경

- Supplied by DAIHAN Korea
- Ideal for Laboratory, Hospital, Circuit Board Check-Up, etc.
- International Standard ϕ 127mm B270 White Glass Lens
- Light Source ; High Brightness White LED Bulbs (36ea)
- Adjustable Height and Lens Angle Direction
- 2-Stage Spring Flexible Arm ■ Scratch and Nick-Resistant PP Lens Cover
- * Ultra Bright LED Light ; 5,000lx * Adjustable Height 300 ~ 350mm
- * Base Size ; 220 × 330mm * Electric Power ; AC 220V, 12W

Ordering Information

DH.Mag3205 Magnifier, Fluorescent Light, × 5, "BLM-L105"
DH.Mag3208 , × 8, "BLM-L108"

<L2>

Ea.	Price, ₩
-----	----------



Fluorescent Light Magnifier Precision-type, "SK102A"

SEKI-Optical® Precision Fluorescent Light Magnifier "SK102A", × 5 & × 8
With Objective Lens × 13 Magnification, ϕ 50mm High Clear White Glass Lens, 초정밀 조명 확대경

- Supplied by DAIHAN Korea
- High Brightness Daylight Blue Fluorescent Light Source
- ϕ 50mm B270 white Glass Lens with × 13 Objective Lens
- 2-Stage Spring Flexible Arm
- Adjustable Height 300 ~ 500mm and Lens Angle Direction
- Ideal for Near-work, Jewelry, Laboratory, Hospital, Circuit Board Check-Up, etc.
- * Base Size ; 220 × 330mm * Electric Power ; AC 220V, 22W

S6.SK102A05 Magnifier, Fluorescent Light, Precision, × 5
S6.SK102A08 , × 8

<L2>

Ea.	Price, ₩
-----	----------

DAIHAN™ Pocket Anemometer of Flow · Temp · RH% · Dew-point · Compass “ANE5”, IP54 Water Resistance & Shock Proof
1.1~20m/s, -15℃+50℃, 0~100% RH, Dew-point & Compass 360°, Tripod-stand Mountable, 소형 다기능 아네모메타



Pocket Anemometer for 3× Measuring

■ Supplied by DAIHAN Korea

- * Small Pocket Knife Shape & Compact Design. * IP54 Water Resistant & Shock Proof.
- * Temp. Units °C/°F Selectable. * Auto Sleep Mode Disable Included. * Tripod Mountable Receptacle.
- * Measures 7 Parameters ; Air Velocity, Compass, Temperature, Humidity, Windchill, Wet Bulb & Dew Point
- * 6 Wind Units Selectable ; KNT(Knot), BF(Beaufort), FPM(Feet Per Minute), MPH(Mile Per Hour), MPS(Meter Per Second), KPH(Kilometer Per Hour).

Specifications & Ordering information

<L1>

Spec.	Model	“ANE5”
3× Measuring Range		(1) Wind Speed ; 0.5 ~ 44.7m/h, 60 ~ 3937ft/m, 0.4 ~ 38.8knt, 1.1 ~ 20.0m/s, 0.7 ~ 72.0km/h, 1 ~ 8bf, (2) Temp. ; -15℃+50℃, (3) RH% ; 0%RH ~ 100%RH
Accuracy		(1) Wind Speed ; ±(3% of Reading +0.2m/s), (2) Temp. ; ±1℃, (3) RH% ; ±3%RH(at 25℃, 10 ~ 90%RH, other ±5%RH)
Resolution		(1) Temp. ; 0.1℃/°F, (2) RH% ; 0.1%RH
Response time		Temp. & RH% ; 60 seconds (typical)
Display		(1) Windchill ; -20.0℃+50.0℃, (2) Dew Point ; -20.0℃+50.0℃, (3) Wet Bulb ; -5.0℃+50.0℃, (4) Compass ; 0~360 Degree
Operating Temp. & RH%		-15℃+50℃ & Humidity < 80%
Storage Temp. & RH%		-20℃+50℃ & Humidity < 90%
Battery		1×CR2032 Battery
Dimension & Weight		140×45×h25 mm (90° Folded), 235×45×h25 mm (180° Full Opened) / 90g LCD Size ; 16×28mm
Standard Package		Meter / Battery / Manual / Strap / Paper Box
Cat. No. & Price ₩		DH.Ane3005 Set
Spare ; DH.Ane3005.1	Mini Fan for Anemometer	Ea

DAIHAN™ Combo-Anemometer of Flow · Temp · RH% · Dew-point “ANE3”, PC-interface, Datalogging
0.5m/s, -20℃+60℃, 0.1~99.9% RH, -20℃+59.9℃ D-point, & Wet Bulb, with 18mm Vane Probe & Cable, 정밀 아네모메타



Combo - Anemometer Hard Carrying Case

■ Supplied by DAIHAN Korea

- * Measures Air Velocity, Volume, Temperature, Humidity, Dew Point & Wet Bulb Temperature.
- * Thermistor Temp. Sensor and Capacitive Humidity Sensor for Quick and Accurate Measurement.
- * 18mm Mini Vane Probe with Built-in Temp. & Humidity Sensor. * 99 Single Point Memories.
- * Winds Peed Metric and Imperial Units Selectable. * Up to 12000 Points Auto Datalogging Points.
- * PC Interface for Data download * Date and Time Stamp (12-Hour or 24-Hour)

Specifications & Ordering information

<L1>

Spec.	Model	“ANE3”
3× Measuring Range		(1) Wind Speed ; 0.5 ~ 20M/S, (2) Temp. ; -20℃+60℃, (3) RH% ; 0.1 ~ 99.9%RH
Accuracy		(1) Wind Speed ; ±(3% of reading + 0.2M/S), (2) Temp. ; ±0.6℃, (3) RH% ; ±3%RH(at 25℃, 10 ~ 90%RH, other ±5%RH)
Resolution		(1) Temp. ; 0.1℃/°F, (2) RH% ; 0.1%RH
Response time		Temp. & RH% ; 60 seconds (typical)
Display		(1) Wet Bulb & Dew Point Temp. ; -20℃+59.9℃, (2) Air Volume ; 0 to 99999m³/minute
Operating Temp. & RH%		0℃~50℃ & Humidity < 80%
Storage Temp. & RH%		-10℃+50℃ & Humidity < 90%
Battery		6×AAA Battery or 9V DC Adaptor
Dimension & Weight(g)		180×58×h24.5 mm (LCD Size ; 54×37 mm) / 200g
Standard Package		Meter / Probe / Battery / Manual / Hard Carrying Case / Software Kit
Cat. No. & Price ₩		DH.Ane3003 Set
Spare ; DH.Ane3003.1	Expanded Vane Probe, 273×18.5mm (1055mm Extended), w/Rubber Holder	Ea

DAIHAN™ Digital Anemometer Set "A9201", for Air Velocity & Temperature, Max/Min/Average
 With 10 Memories Data Storage · Recalling · Clearing, m/s · m/hr · ft/min · knots, and °C/°F, 디지털 아네모메타 세트



Digital Anemometer Set, for Air Velocity / Temp.

■ Supplied by DAIHAN Korea

- * Ideal for Professional in Aviation, Laboratory and Ventilation Systems. * CE Approval.
- * Instant Display of the Air Velocity and Temperature. * Maximum, Minimum and Average Calculation.
- * Selectable Air Velocity Unit m/s, km/hr, mph, ft/min or Knots. * Switchable Temperature Unit °C/°F.
- * 10 Memories Data Storage/Recalling/Clearing. * 1/4in. 20 UNC Threading for Tripod Mounting.

Specifications & Ordering information

<L1>

Spec.	Model	"A9201"
2 × Measuring Range		(1) Temp.; -50°C + 70°C or -58°F + 158°F, (2) Air Velocity; 0.5 to 45.0m/s · 1.8 to 162.0km/hr · 1.1 to 100.0mph, 98 to 880ft/min · 1.0 to 88.0knots
Accuracy		(1) Temp.; ± 0.5°C (for 0°C ~ 50°C) or ± 0.9°F (for 32°F ~ 122°F), otherwise ± 1°C or 1.8°F, (2) Air Velocity; ± 2%
Resolution		(1) Temp.; 0.1 for whole Range, (2) Air Velocity; 0.1
Ambient Temp.		0°C ~ 50°C or 32°F ~ 122°F
Probe Type		NTC-Probe; 76 × 67 × 32mm (1Meter Cable)
Color		Grey
Battery		3 × 1.5V AAA Battery
Dimension / Weight(g)		Unit; 85 × 175 × h39mm (LCD Size; 54 × 35 mm) / 200g Vane-Probe; 76 × 67 × 32mm, Length 125mm / 130g
Cat. No. & Price ₩	A1.A9201	Ea

DAIHAN™ Precise 10cm Vane Anemometer of Flow · Temp · RH% · Dew-point · Wet Bulb · Air Volume "ANE7"
 With Optional Air Flow Cones, 0.2~30m/s, -20°C + 60°C, 0.1~99.9% RH, D-point & 99,999m³/min., 정밀 아네모메타



Precise 10cm Vane Anemometer



■ Supplied by DAIHAN Korea

- * Measures Air Velocity, Volume, Temperature, Humidity, Dew Point & Wet Bulb Temperature.
- * Built-in Thermistor Temp. & Capacitive Humidity Sensor for Quick and Accurate Measurement.
- * Unlimited Points & Time Weighted Average Auto Calculation. * Display Max./Min. Value.
- * Hold Function Freezes Current Readings * Inlet and Outlet Measurement is Doable
- * Auto Power off 20mins Inactivity. * Low Battery Indicator * Blue Backlight for Dark Area Operation.

Specifications & Ordering information

<L1>

Spec.	Model	"ANE7"
3 × Measuring Range		(1) Wind Speed; 0.2 ~ 30M/S, (2) Temp.; -20°C + 60°C (3) RH%; 0.1%RH ~ 99.9%RH
Accuracy		(1) Wind Speed; ± (1.5% of Reading + 0.3m/s) for under 20m/s, ± (3% of Reading + 0.3m/s) for above 20m/s, (2) Temp.; ± 0.6°C, (3) RH%; / ± 3%RH (at 25°C, 10 ~ 90%RH, other ± 5%RH)
Resolution		(1) Temp; 0.1°C/°F, (2) RH%; 0.1%RH
Response time		Temp. & RH%; 60 seconds (typical)
Display		(1) Wet Bulb Temp.; -20°C + 59.9°C, (2) Dew Point Temp.; -5°C + 59.9°C, (3) Air Volume; 0 to 99999m³/minute
Operating Temp. & RH%		0°C ~ 50°C & Humidity < 80%
Storage Temp. & RH%		-10°C + 50°C & Humidity < 90%
Battery		4 × AAA Battery
Dimension & Weight		269 × 106 × h51 mm (LCD Size; 32.5 × 54 mm) / 200g
Standard Package		Meter / Manual / Battery / Hard Carrying Case
Cat. No. & Price ₩	DH.Ane3007	Set
Option; DH.Ane3007.1 Air Flow Cone, Round(φ210mm) + Square(346 × 346mm) Set		Set



DAIHAN™ Methane & Propane Gas Leak Detector "GAS3", with Alarm, 40~640 ppm With 5-Level Leakage Indicator & Auto-Calibration Function, 메탄 & 프로판가스 누출 검지기

- * Low Cost Handheld type Detector Checking Methane and Propane Gas Leakage.
- * Bright 5-Level Leakage LED Indicator.
- * Quality & Reliable Semi-Conductor Sensor.
- * Built in Auto-Calibration Function while Meter Powers on.
- * Long & Flexible Stainless Goose Neck Detect Probe.
- * Status Ready and Battery Low LED Indicator.
- * Acoustic Intermittent Beep.
- * Auto Power off Timeout about 10 mins.
- * Built in Earphone Jack for Mute Mode Use.

Specifications & Ordering information

<L1>

Spec.	Model	"GAS3"
Methane / Propane Range		40 ~ 640 ppm
Sensitivity		40 ppm (Methane & Propane)
Response Time & Warm Up		within 10 seconds / within 60 seconds
Battery		4 × AA Battery or 9V Adaptor
Flexible Probe Length		18" (465mm)
Dimension & Weight(g)		175 × 47 × h28 mm / 120g
Standard Package		Meter / Manual / Battery / Carrying Case
Cat. No. & Price ₩	DH.Gas3003	Set



DAIHAN™ Pocket type CO Detector "GAS5", with Dual LCD Display, Alarm, 0~999 ppm With Wrist Strap/Soft Carrying Pouch, 핸디 일산화탄소(CO) 디텍터

- Supplied by DAIHAN Korea
- * Easy to Manually Calibrate at 0 or 500ppm Carbon Monoxide Standard Gas.
- * Low Battery Indicator.
- * Dual LCD Display with Current Carbon Monoxide and Maximum CO Readings.
- * Long Life and Reliable Electrochemical Carbon Monoxide Sensor.
- * Backlight Function Works in Dark Area.
- * Warning Carbon Monoxide Audible Alert. Alarm Adjustable, Preset at 25ppm.

Specifications & Ordering information

<L1>

Spec.	Model	"GAS5"
CO Range & Resolution		0 ~ 999 ppm / 1 ppm
Accuracy (at 20 ± 5°C, 50 ± 20%RH)		±20% at 0 ~ 100ppm ; ±15% at 100 ~ 500ppm
Response Time / Battery		within 60 seconds / 3 × AAA Battery
Dimension & Weight(g)		175 × 47 × h28 mm (LCD Size ; 11 × 26 mm) / 120g
Standard Package		Meter / Manual / Battery / Plain Box / Wrist Strap / Soft Carrying Pouch
Cat. No. & Price ₩	DH.Gas3005	Set

DAIHAN™ Combo Air Quality Meter of CO · CO₂ · Temp · RH% · Dew-point · Wet Bulb "GAS7", PC Data Analysis 0~1000ppm-CO, 0~9999ppm-CO₂, -20°C+60°C, 0.1~99.9% RH, -20°C+59.9°C-DP, -5°C+59.9°C-WB, Alarm 콤보 에어테스터



- * Manual Calibration is available for CO₂ / CO / RH%
- * Designed with electrochemical CO Sensor and NDIR(Non-dispersive Infrared) Waveguide Technology CO₂ Sensor
- * Warning Carbon Monoxide Level Programmable
- * Max, Min, Average Functions Included
- * USB Output to PC for Data Analysis
- * Audible Alarm (~80dB) Threshold Setting
- * 99 Points Memory Review

Specifications & Ordering information

<L1>

Spec.	Model	"GAS7"
6 × Measuring Range		(1) CO ₂ ; 0 ~ 9999 ppm(5001 ~ 9999 ppm out of scale range), (2) CO ; 0 ~ 1000 ppm, (3) Temp ; -20°C + 60°C (-5°F + 140°F), (4) RH ; 0.1 ~ 99.9%RH, (5) DP ; -20°C + 59.9°C, (6) WB ; -5°C + 59.9°C
Accuracy		(1) CO ₂ ; ±30 ppm ±5% of reading (0 ~ 5000 ppm), Other ranges are not specified, (2) CO ; ±10 ppm for less than 100 ppm / ±10% of reading for 101 ~ 500 ppm, ±20% of reading for 501 ppm and above, (3) Temp ; ±0.6°C/±0.9°F, (4) RH ; ±3%RH(at 25°C, 10~90%RH), ±5%RH(at 25°C, other range)
Resolution		1 ppm, 0.1°C/°F, 0.1%RH
Response Time		(1) CO ₂ ; < 30 seconds (90% step change), (2) CO ; < 60 seconds (90% step change), (3) Tair ; < 2 mins (90% step change), (4) RH ; < 10 mins (90% step change)
Power Supply		4 × AA Battery or 9V adaptor (>24hours, Alkaline Battery)
Dimension & Weight(g)		205 × 70 × h56 mm(LCD Size ; 26 × 44 mm) / 200g
Standard Package		Meter / Manual / Battery / USB Cable & CD / Carrying case
Cat. No. & Price ₩	DH.Gas3007	Set

Option ; DH.Gas3007.1 ① Cal. Salt Bottle, for Calibration, 33%, 56ml
DH.Gas3007.2 ② , 75%, 56ml

Ea

Ea

DAIHAN™ Air Quality Meter of CO₂ · Temp · RH% · Dew-point · Wet-Bulb "GAS13", PC Data Analysis
 0~9999ppm-CO₂, -10℃+60℃, 0.1~99.9% RH, -20℃+59.9℃-DP, -5℃+59.9℃-WB, 80dB-Alarm, 다기능 에어 미터



■ Supplied by DAIHAN Korea

- * Big LCD display ; CO₂, Temp and Relative Humidity, Display Air / Dew Point / Wet Bulb Temp by Selection
 * NDIR(Non-Dispersive Infrared) Waveguide Technology CO₂ Sensor * Warning Carbon Dioxide Level

<L1>

Spec.	Model	"GAS13"
5× Measuring Range		(1) CO ₂ ; 0 ~ 9999 ppm, (2) Temp ; -10℃+60℃(14°F ~ 140°F), (3) RH% ; 0.1 ~ 99.9%RH, (4) DP ; -20.0℃+59.9℃, (5) WB ; -5.0℃+59.9℃
Accuracy		(1) CO ₂ ; ±50 ppm ±5% of reading (0 ~ 2000 ppm), (2) Temp ; ±0.6℃ / ±0.9°F, (3) RH ; ±3% RH(at 25℃, 10 ~ 90%RH)
Resolution		1 ppm, 0.1℃/°F, 0.1%RH
Response Time		(1) CO ₂ ; <30 sec, (2) Tair ; <2 min, (3) RH ; <10 min. (90% step change)
Warm-up Time / Battery		30 sec. / 4×AA Battery
Dimension / Weight(g)		205×70×h56mm(LCD Size ; 26×44mm) / 200g
Standard Package		Meter / Manual / Battery / Hard Carrying Case
Cat. No. & Price ₩		DH.Gas3013 Set
Option ; DH.Gas3013.1		① USB Cable & Software Ea
DH.Gas3013.2		② RS232 Cable & Software Ea
DH.Gas3007.1		③ Cal. Salt Bottle, for Calibration, 33%, 56ml Ea
DH.Gas3007.2		④ , 75%, 56ml Ea

DAIHAN™ Wallmount Desktop Air Quality Monitor of CO₂ · Temp · RH% · Dew-point · Wet-Bulb "GAS15", PC Data Analysis
 0~9999ppm-CO₂, -10℃+60℃, 0.1~99.9%RH, -20℃+59.9℃-DP, -5℃+59.9℃-WB, NDIR-Sensor, 다기능 벽걸이 에어 모니터



- * Super Large LCD display CO₂ Value * Warning Carbon Dioxide Level with 80dB Alarm
 * NDIR(Non-Dispersive Infrared) Waveguide Technology with ABC(Automatic Background Calibration)
 * Long Time Drift Compensation * Alarm Output Drive Relay to Control System

<L1>

Spec.	Model	"GAS15"
5× Measuring Range		(1) CO ₂ ; 0 ~ 9999 ppm (2001~9999 ppm out of Scale Range), (2) Temp ; -10℃+60℃(14°F ~140°F), (3) RH ; 0.1 ~ 99.9%RH, (4) DP ; -20.0℃+59.9℃, (5) WB ; -5.0℃+59.9℃
Accuracy		(1) CO ₂ ; ±50 ppm ±5% of reading(0 ~ 2000 ppm), (2) Temp ; ±0.6℃ / ±0.9°F, (3) RH ; ±3% RH(at 25℃, 10 ~ 90%RH) or ±5% RH(at 25℃, others)
Resolution		1 ppm, 0.1℃/°F, 0.1%RH
Response Time		(1) CO ₂ ; <30 sec, (2) Tair ; <2 min, (3) RH ; <10 min. (90% step change)
Warm-up Time		30 sec.
Power Supply		12V DC Adaptor
Dimension / Weight(g)		130×85×h60mm (LCD Size ; 50×65 mm) / 200g
Standard Package		Meter / Manual / 12V Plug Adaptor / Plain Box
Cat. No. & Price ₩		DH.Gas3015 Set
Option ; DH.Gas3013.1		① USB Cable & Software Ea
DH.Gas3013.2		② RS232 Cable & Software Ea
DH.Gas3007.1		③ Cal. Salt Bottle, for Calibration, 33%, 56ml Ea
DH.Gas3007.2		④ , 75%, 56ml Ea

DAIHAN™ "Green Life" Indoor Air Monitor of CO₂ · Temp · RH% · Time Clock "GAS9", Desktop-model, 15 Degree Tilt Angle
 With CO₂ Alarm, 0~9999ppm with Auto-Background Calibration, -10℃+60℃, 0.1~99.9%RH, "그린라이프" 룸에어모니터



Indoor Air Monitor

■ Supplied by DAIHAN Korea

- * Display CO₂ Value (with Good, Normal, Poor Air), Year-Month-Date, Hour-Minute, Air Temp & Relative Humidity
 * NDIR (Non-dispersive Infrared) Waveguide Technology * Warning Carbon Dioxide Level from 1000~5000 ppm

<L1>

Spec.	Model	"GAS9"
3× Measuring Range		(1) CO ₂ ; 0 ~ 9999 ppm (2001 ~ 9999 ppm out of scale range), (2) Temp ; -10℃+60℃(14°F ~ 140°F), (3) RH% ; 0.1 ~ 99.9%RH
Resolution		1 ppm, 0.1℃/°F, 0.1%RH
Accuracy		(1) CO ₂ ; ±50ppm ±5% of reading (0 ~ 2000 ppm), (2) Temp ; ±0.6℃/±0.9°F, (3) RH% ; ±5%RH(at 25℃, 20~90%RH), ±7%RH(at 25℃, other range)
Response time		(1) CO ₂ ; < 2 mins (90% step change), (2) Tair ; < 2 mins (90% step change), (3) RH ; < 10 mins (90% step change)
Power Supply		5V DC Adaptor(±10%), >= 500mA
Dimension & Weight(g)		120×80mm(LCD Size ; 60×82 mm) / 200g
Standard Package		Meter / Manual / Adaptor / Plain Box
Cat. No. & Price ₩		DH.Gas3009 Set



(1) "VAR"
With Controller

(2) "VAC"
Without Controller

V.D.Heijden® Mechanical Precise Vacuum Manometer "VAR" & "VAC", NO-Mercury, 0~1020 mbar
With or Without Controller, 메카니컬 정밀 마노미터, 진공조절식("VAR" 모델)

■ Supplied by V.D.Heijden® in Germany

■ This device is serially switched with an Electric Vacuum Pump (Oil or Diaphragm Pump)

* This vacuum is displayed in mbar (no more conversion necessary).

* 2× Barrel Clamps are Attached to the Back of this Meter so that the Device can be Quickly Clamped to a Supporting Wall.

* The Precise Vacuum is displayed by a Measurement-device without Any Hysteresis Problems. * Reading ; 0 ~ 1020 mbar

* Connection ; 8 mm Nickeled * Measuring device ; 72 × 72 mm * Metal Dimension ; w80 × d80 × h150 mm

(1) Fine Model "VAR" With Controller, 0~1020 mbar

You can Set Your own Vacuum Level.

VD.3000 Mechanical Manometer "VAR", with Controller, 0~1020 mbar

(2) Standard Model "VAC" Without Controller, 0~1020 mbar

VD.3001 Mechanical Manometer "VAC", w/o-Controller, 0~1020 mbar

<L1>

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ Digital Pressure / Vacuum Manometer "MAN3003", 10-Different Units ; psi, mbar, etc

+/- (±) 30psi / 1551mmHg / 2068 mbar, PC Data Analysable, Max-Min / Hold / Record, with 8mm Metal Lug, 디지털 마노미터



Digital Pressure / Vacuum Manometer



Hard Carrying Case



■ Supplied by DAIHAN Korea

■ Portable Digital Gauge ; 10 Differential Units Pressure / Vacuum Meter, with 8mm Metal Connection

* Over-range Indication with Error Messages

* Connection ; 8mm Metal Lug

* Low Battery Indicator

■ PC Interface via Cable and Software to Download Data

* 10 units or 5 units for Selection (Metric, Imperial)

* Backlight Function works in Dark Area

* Real Time Shows in Recording Mode

Specifications & Ordering information

<L1>

Spec.	Model	"MAN3003"
Range & Resolution	psi	± 30.00 & 0.02
	kPa and mmHg	kPa ; ±206.8 & 0.2, and mmHg ; ±1551 & 1
	kg/cm² and mbar	kg/cm² ; ±2.109 & 0.002, and mbar ; ±2068 & 2
	bar and inH ₂ O	Bar ; ±2.068 & 0.002, and inH ₂ O ; ±830.4 & 0.5
	cmH ₂ O and inHg	cmH ₂ O ; ±2101 & 2, and inHg ; ±61.08 & 0.04
	ftH ₂ O and oz/inch²	ftH ₂ O ; ±69.2 & 0.04, and oz/inch² ; ±480.0 & 0.3
Accuracy		±0.3% of Full Scale at 25°C
Response Time		0.5 sec
Compensated Temp Range		0°C~50°C
Battery		9V Battery or Adaptor
Dimension / Weight(g)		182 × 72 × 30 mm (LCD Size ; 29 × 51mm) / 150g
Standard Package		Meter, Battery, Manual, 2 × Connection Hose(50cm), Hard Carry Case
Cat. No. & Price ₩		DH.Man3003 Set
Option ; DH.Gas3013.1		① USB Cable & Software Ea
DH.Gas3013.2		② RS232 Cable & Software Ea
DH.Man3003.1		③ Tubing for 30psi Pressure / Vacuum and 8mm Lug, 50cm, Black Ea

※ Units Conversion ; 1psi × 6.8947=kpa, 1psi × 68.947=mbar(hpa), 1psi × 0.068966=bar, 1psi × 51.71433=torr
1psi × 70.038=cmH₂O, 1psi × 51.715=mmHg, 1psi × 16=oz/inch², 1psi × 6890=pascal(pa)

※ Note! ; Exceeding Maximum Pressure will cause Permanent Sensor Damage !



Salinity Pen-type Meter

DAIHAN™ Pocket Dual Temp & Salinity Pen-type Meter "SAL3", IP65 Waterproof

With Large Display LCD for Salinity and Temp, 0~70.0ppt, 0°C~50.0°C, 휴대용 염도계, IP65 방수

■ Supplied by DAIHAN Korea

■ 디지털 휴대용 염도계, 염도와 온도 동시 디스플레이

* One Touch Auto-calibration & Calibration Values are Adjustable * Auto Temperature Compensation

* Large Display LCD for Simultaneous Output of Salinity and Temperature (w18 × h30 mm)

* Quick Response & Auto Power off (20min / Disable)

* Built-in NaCl conductivity to TDS Conversion Factor

* Hold Function to Freeze Current Reading

* Temperature unit °C/°F Selectable

<L1>

Spec.	Model	"SAL3"
2× Measuring Range		(1) Salinity ; 0.00 ~ 10.0ppt(NaCl), 10.1 ~ 70.0ppt(NaCl), (2) Temp ; 0°C~50.0°C
2× Accuracy		(1) Salinity ; ±1%F.S ±1digit, for 0 ~ 10.00ppt, ±2%F.S for 10.1 ~ 70.0ppt, (2) Temp ; ±0.5°C
2× Resolution		(1) Salinity ; 0.01, for 0 ~ 10.00ppt, 0.1 for 10.1 ~ 70.0ppt, (2) Temp ; 0.1
Sensor life time / Battery		> 6 months (with Good Maintenance) / LR44 × 4pcs.
Dimension / Weight(g)		165 × 35 × h32mm (LCD Size ; 30 × 18mm) / 115g
Standard Package		Meter / Battery / Manual / Paper Box
Cat. No. & Price ₩		DH.Sal3003 Set

DAIHAN™ Compact Digital Sound Level Meter “SOU3”, 40~130dB, 0.1dB/1dB Resolution
 With Digital Display (30×35mm) & Bar Graph, Pocket Size Rigid Design, 휴대용 소음계



Pocket-type

■ **Supplied by DAIHAN Korea**

- * Used for Acoustic (sound that travels through air) Measurements, which had been widely used in the HVAC(Heating, Ventilation and Air Conditioning) or Other Industrial Related Application.
- * The Most Cost-effective Best Selling Meter
- * Fast / Slow Response is selectable.
- * A weighting enables the Meter to Respond in the Same Manner as the Human Ear.
- * C weighting is suitable for Flat Response Measurements with No Increase or Decrease of Amplitude over the Frequency Spectrum.
- * Auto / Manual Ranging, 4 Measurement Ranges in 10 Steps.
- * A/C Weighting is selectable
- * Tripod Mountable

Specifications & Ordering information

<L1>

Spec.	Model	“SOU3”
2× Measuring Range		(1) A Weight ; 40dB ~ 130dB (2) C Weight ; 45dB ~ 130dB
Accuracy at 94dB, 1KHz		± 2dB
2× Resolution		(1) Digital Display ; 0.1dB (2) Quasi-Analog Bar ; 1dB
Digital / Quasi-Analog Display (4 Range)		40 ~ 70dB, 60 ~ 90dB, 80 ~ 110dB, 100 ~ 130dB
Digital-/ Bar Updating Time		160mS / 40mS
Microphone		Dia.6mm Electric Condenser Microphone
Battery / Battery Life Time		AAA×4pcs / 60 Hours Typical
Dimension / Weight(g)		230×57×h44mm(LCD Size ; 30×35mm) / 130g
Standard Package		Meter / Battery / Manual / Paper Box
Cat. No. & Price ₩		DH.Sou3003 Set

DAIHAN™ Digital Sound Level Meter “SOU5”, with Deluxe Super Large LCD, 30dB~130dB, 0.1dB/1dB Resolution
 With Digital Display(42×62mm) & Analog Bar Graph, Max-Min Backlight, 디지털 소음계



Hard Carrying Case



①

②

■ **Supplied by DAIHAN Korea**

- * Meet IEC60651, Type2 Class Standard.
- * Using the Background Noise Absorber.
- * Backlight Function to Use in Dark Area.
- * Automatic Shut off.
- * Auto/Manual Ranging, 6 Measurement Ranged in 10 Steps.
- * Freezing the Maximum Sound Level Reading.
- * Easy to Replace Battery.
- * Max./Min. Recording.
- * A/C Weighting is Selectable.
- * Easy to Calibrate by Standard Calibrator
- * AC/DC Analog Output
- * Professional and Rigid Design.
- * RS232 interface for PC download.

Specifications & Ordering information

<L1>

Spec.	Model	“SOU5”
2× Measuring Range		(1) A Weight ; 30dB ~ 130dB (2) C Weight ; 35dB ~ 130dB
Accuracy at 94dB, 1KHz		± 1.5dB
2× Resolution		(1) Digital Display ; 0.1dB (2) Quasi-Analog Bar ; 1dB
Digital / Quasi-Analog Display (6 Range)		30 ~ 80dB, 40 ~ 90dB, 50 ~ 100dB, 60 ~ 110dB, 70 ~ 120dB, 80 ~ 130dB
Digital-/ Bar Updating Time		160mS / 40mS
Analog Output		AC 0.707Vrms(at Full Scale) ; DC 10mV/dB
Microphone		Dia.6mm Electric Condenser Microphone
Frequency Range		31.5Hz ~ 8KHz
Battery / Battery Life Time		9V Battery×1pcs. or 9VDC Adaptor / 30 Hours Typical
Dimension & Weight(g)		254×77×h37mm (LCD Size ; 42×62mm) / 200g
Standard Package		Meter / Battery / Manual / Carry Case
Cat. No. & Price ₩		DH.Sou3005 Set

Option ; DH.Gas3013.1

① USB Cable & Software

Ea

DH.Gas3013.2

② RS232 Cable & Software

Ea

DAIHAN™ Digital Pro-Tachometer “TAC3”, Contact & Non-Contact 2 in 1 Meter, 10.0~99,999rpm, Auto Power off in 20min
Non-Contact type with Red LED-Beam, Contact type with Convex Tip or Concave Tip, 프로급 타코미터, 접촉/비접촉 겸용



Digital Pro-Tachometer

- Supplied by DAIHAN Korea * Three Steps to Switch Operation Mode
* Speed Units are Selectable in Contact Mode ; rpm, m.min, ft/min, yd/min

Specifications & Ordering information

<L1>

Spec.	Model	“TAC3”
3× Measuring Range		(1) Rotation (Non-contact) ; 10.0 ~ 99999rpm (2) Rotation (Contact) ; 10.0 ~ 19999rpm (3) Surface Speed (Contact) ; 1.0000 ~ 9999.9M/Min.
Accuracy / Measuring Distance		± (0.1% of rdg. +2digit) / 5 ~ 30cm (Non-contact)
Resolution		rpm ; 0.001 / 0.01 / 0.1 / 1, Others ; 0.0001 / 0.001 / 0.01 / 0.1
Time Base & Sampling Time		4.0Mhz Quartz Crystal / 1sec.(> 60rpm), > 1sec.(10-60rpm)
Dimension / Weight(g)		165×51×h32mm (LCD Size ; 15×33 mm) / 90g
Standard Package		Meter / Battery / Manual / Reflective Tape / Convex Tip / Concave Tip
Cat. No. & Price ₩	DH.Tac3003	Set
Spare ; DH.Tac3003.1	(A) Convex Tip for contact Rotate Speed Use	Ea
DH.Tac3003.2	(B) Concave Tip for contact Rotate Speed Use	Ea
DH.Tac3003.3	(C) Reflective Tape for Non-Contact Rotate Speed, 16×600mm	Ea

DAIHAN™ Digital Pocket Tachometer “TAC5” & “TAC7”, Contact or Non-Contact 2-type, 10~99999rpm, Auto Power off in 20min
Non-Contact type with Red LED-Beam, Contact type with Convex Tip or Concave Tip, Hold Function, 접촉식 or 비접촉식 타코미터

Non-Contact type
(DH.TAC5)Contact type
(DH.TAC7)

- Supplied by DAIHAN Korea * Speed Units are Selectable in Contact Mode ; RPM, M/Min, Ft/Min, Yd/Min
* Contact Type ; Straight-line Contact with the Center of the Rotating Body
* Non-contact Type ; Reflective Tape Affixed to an Object and Measuring the Rotation. (50 to 150mm)

Specifications & Ordering information

<L1>

Spec.	Model	“TAC5”	“TAC7”
Measure Range		Rotation(Non-contact) ; 10.0 ~ 99999rpm	Rotation(Contact) ; 10.0 ~ 19999rpm, Surface Speed(Contact) ; 1.000~9999.9M/Min
Accuracy		± (0.1% of rdg. +2digit)	
Resolution		rpm ; 0.001 / 0.01 / 0.1 / 1, Others ; 0.0001 / 0.001 / 0.01 / 0.1	
Measure Distance		5~30cm	Contact
Power Consumption / Time Base		5uA(Idle), 18mA(Measurement) / 4.0Mhz Quartz Crystal	
Sampling Time / LCD Size		1Second(> 60RPM), > 1Second(10~60RPM) / 33×h15mm	
Battery		1×9V Battery	
Dimension / Weight(g)		124×50×h33mm / ~80g	149×50×h33mm / ~85g
Standard Package		Meter / Battery / Manual / Paper Box / Reflective Tape / Convex Tip / Concave Tip	
Cat. No. & Price ₩	DH.Tac3005	Set	DH.Tac3007
Spare ; DH.Tac3003.1	(A) Convex Tip for contact Rotate Speed Use		Ea
DH.Tac3003.2	(B) Concave Tip for contact Rotate Speed Use		Ea
DH.Tac3003.3	(C) Reflective Tape for Non-Contact Rotate Speed, 16×600mm		Ea

DAIHAN™ Handy/Digital Thermo-Hygrometer “H9237”, Rotary Sensor Protector, °C/°F & R.H%, Max/Min
With Large LCD, Rubber Coated & Black Carrying Case, -50°C + 70°C, 20~99.9% R.H, 0.1 Divi, 포켓/디지털 온습도계



Size ;
165 × 38 × 18 mm



■ 온도, 습도 측정이 필요한 모든 곳에서 편리하게 사용 가능

- * Ideal for Storage Rooms, Pharmaceutical Rooms, Greenhouse, Laboratory, etc.
- * Pocket Size with Large LCD.
- * Rubber Coated.
- * Low Battery Indicator.
- * Rotary Sensor Protector Protects The Sensor when not in Use.
- * Battery Included.
- * Easy to Change Battery.
- * with Black Carrying Case.
- * CE Approval.

Specifications & Ordering information

<L1>

Spec.	Model	“H9237”
Measurement Range		-50°C + 70°C or -58°F + 158°F / 20%~99.9% R.H (No Condensation)
Resolution		0.1° for Whole Range/ 0.1% R.H for Whole Range
Accuracy		± 1.0°C or ± 1.8°F In the Range of 0°C~50°C or -4°F + 122°F, Otherwise more than ± 1.0°C or ± 1.8°F
Sensor Type		NTC & Humidity
Ambient Temperature		0°C~50°C or 32°F~122°F
Weight (g) / Color		55g / White
Battery / Package		1 × 3 Volts “CR 2032” Size or Equivalent / White Box
Dimension		165 × 37 × 18mm (LCD Size ; 58 × 12mm)
Cat. No. & Price ₩	A1.H9237	Ea

DAIHAN™ Pocket Dew-Point Thermo-Hygrometer of Temp · RH% · DP°C · WB°C “THH5”, Max/Min, Hold
Optional External Temp Probe(-20°C + 70°C), -20°C + 50°C, 0~100% RH, -78.7°C + 50°C DP, 포켓 다기능 이슬점 온습도계



■ Supplied by DAIHAN Korea

■ 간편히 휴대하면서 온습도 측정이 가능한 제품

- * Designed with Capacitance-type Humidity Sensor to get Accurate and Quick Reading
- * Simultaneously Triple Display Temp / RH% / External Probe Temp
- * Maximum / Minimum and Reset function
- * Temperature unit °C/°F Selectable
- * External Probe for Checking Food Quality & HVAC/R Applications
- * Auto Power Off time frame Selectable
- * Hold function to Freeze current Reading

Specifications & Ordering information

<L1>

Spec.	Model	“THH5”
4 × Measuring Range		(1) Temp ; -20°C + 50°C, (2) Humidity ; 0%RH~100%RH, (3) Dew Point ; -78.7°C + 50.0°C, (4) Wet Bulb ; -21.6°C + 50.0°C
2 × Accuracy		(1) Temp ; ± 0.6°C, (2) Humidity ; ± 3%RH (at 25°C, 10~90%RH, other ± 5%RH)
2 × Resolution		(1) Temp ; 0.1°C/°F, (2) Humidity ; 0.1%RH
Battery		2 × AAA Battery
Dimension / Weight(g)		48.4 × 178.5 × 24.7mm(LCD Size ; 28 × 33 mm) / 100g
Standard Package		Meter / Battery / Manual / Paper Box
Cat. No. & Price ₩	DH.Thh3005	Set
Option ; DH.Gas3007.1	① Cal. Salt Bottle for Calibration, 33%, 56ml	Ea
DH.Gas3007.2	② , 75%, 56ml	Ea
DH.Thh3005.1	③ External Temp Probe, 116cm, with Handle	Ea

DAIHAN™ Handy Dew-Point IR Thermo-Hygrometer for Air & Surface Temp · RH% · DP °C “THH11”

Built in Laser-Beam (670nm), Backlight, PC Data Analysis, -40°C + 500°C IR, 0~100% RH, -68°C + 49.9°C DP, 휴대용 IR 온습도계



■ Supplied by DAIHAN Korea

- * The Most Popular IR+ Humidity Two in One Meter
- * Triple Display Measures Air/Surface Temp & Humidity
- * Auto Power off Time Adjustable
- * Gear Damper Device (Design Patent) for Sensor Protection
- * IR Emissivity Value Adjustable 0.3~1.0
- * Laser Beam & Backlight Function Selectable

■ IR온도계와 습도계 두가지를 동시에 측정 가능한 제품

- * Design with Precision Humidity Sensor
- * Temperature Unit °C/°F is Switchable
- * PC Output Data Download
- * Powered by 4pcs AAA Batteries

Specifications & Ordering information

<L1>

Spec.	Model	“THH11”		
Measure Range		-40°C + 500°C	-20°C + 50°C	0~100%RH
Resolution		0.1°C / 0.1°F	0.1°C / 0.1°F	0.1%RH
Accuracy		±2°C at 20~200°C, Others is ±3°C	±0.6°C	±3%RH (at 25°C, 10~90%RH, others ±5%RH)
3×Measuring Range		(1) IR Temp ; -40°C + 500°C, (2) Temp ; -20°C + 50°C, (3) Humidity ; 0~100%RH		
3×Accuracy		(1) IR Temp ; ±2°C at 20~200°C, Others is ±3°C, (2) Temp ; ±0.6°C, (3) ±3%RH (at 25°C, 10~90%RH, others ±5%RH)		
3×Resolution		(1) IR Temp ; 0.1°C/°F, (2) Temp ; 0.1°C/°F, (3) Humidity ; 0.1%RH		
Dew Point & Wet Bulb Temp		-68°C + 49.9°C & -21.6°C + 49.9°C		
IR Repeatability		±1°C		
Emissivity		0.3~1.0 adjustable		
Dimension / Weight(g)		70×50×h175mm(LCD Size ; 44×26mm) / 140g		
Standard Package		Meter / Battery / Manual / Carry Case		
Cat. No. & Price ₩		DH.Thh3011		
Option ; DH.Gas3013.1		① USB Cable & Software		
DH.Gas3013.2		② RS232 Cable & Software		
				Ea

SciLab™ Desktop/Wallmount Thermo-Hygrometer “THE3001”, Jumbo LCD Display with 3-Row for Humidity · Temperature · Time With 1.5m Probe Cable, Hourly Chime · Daily Alarm · Calendar Display-Function, 대형 온습도계

- * Ideal for Use in Incubators, Drying Chambers, Hoods, Cleanrooms, Cabinets and General Lab Applications
- * Function ; ① Temperature/Humidity Display ② Max/Min Temperature and Humidity Memory
- ③ 12/24 Hourly Clock ④ Hourly Chime and Daily Alarm
- ⑤ Calendar Display

Specifications & Ordering information

<L1>

Spec.	Model	“THE3001”	
Measuring Range		(1) Temperature ; -50°C + 70°C / -58°F + 158°F	(2) Humidity ; 10~90%RH
Accuracy		(1) Temperature ; ±1°C (1.8°F)	(2) Humidity ; +5%RH
Resolution		(1) Temperature ; 0.1°C (0.2°F)	(2) Humidity ; 1%
Dimension / Weight(g)		100×110×21.5 mm (Body), 80×60 mm(LCD)	
Sensor Type		NTC & Humidity External Sensor ; 1.5m Sensor Wire with Probe(L22.5×φ5.5mm)	
Standard Package		Meter / Manual / Plastic Blister Package	
Battery		1×AAA 1.5V	
Cat. No. & Price ₩		SL.The3001	
			Set

DAIHAN™ 2- & 3-Row Jumbo Display Thermo-Hygrometer “H9213” & “H9216”, for Indoor/Outdoor, °C/°F & R.H.%
With 3m Cord/Probe, Memory Recalling, -50°C + 70°C, 20~99.9% R.H, 0.1 Divi, 2 & 3열 실내·외 온습도계



Size ;
120 × 110 × 25 mm

(1) Jumbo 2-Row Display , “H9213”



Stand & Hanger



Size ;
120 × 110 × 25 mm

(2) Jumbo 3-Row Display , with Alarm Timer
“H9216”



Stand & Hanger

■ 온도 및 습도를 센서로 측정하거나 측정프로브를 통해 측정하는 방법 선택 가능

- * Ideal for Office, Home, Greenhouse, Laboratory. * (1) Dual or (2) Triple Row Display for Easy Reading.
- * Comfortable and Uncomfortable Index (Comfort/Warm/Cool). * Jumbo LCD.
- * 3m Sensor Wire with Probe for Outdoor Temperature Measurement. * Temperature Unit °C/°F Selectable.

(1) With 2-Row Jumbo Display ; On/Off, °C/°F, Max/Min (°C/°F), Temperature and Humidity, Indoor Humidity
Measurement Indoor/Outdoor Temperature Measurement, Memory Recalling.

Specifications & Ordering information

<L1>

Spec.	Model	“H9213”
Measurement Range		-50°C + 70°C or -58°F + 158°F / 20%~99.9% R.H (No Condensation)
Resolution		0.1° for Whole Range/ 0.1% R.H for Whole Range
Accuracy		± 1.0°C or ± 1.8°F In the Range of -20°C + 70°C or -4°F + 158°F Otherwise more than ± 1.0°C or ± 1.8°F
Ambient Temperature		0°C~50°C or 32°F~122°F
Weight (g) / Color		165g / White
Battery / Package		2 × 1.5 Volts “AAA” Size or Equivalent / Sliding Blister with Print Card
Sensor Type		NTC & Humidity External Sensor ; 3m Sensor Wire with Probe(L22.5 × φ 5.5mm)
Dimension		120 × 110 × 25mm (LCD Size ; 72 × 61mm)
Cat. No. & Price ₩	A1.H9213	Ea

(2) With 3-Row Jumbo Display, with Alarm Timer ; On/Off, °C/°F, Max/Min (°C/°F), Temperature and Humidity
Memory Recalling, Clock (HH/MM/SS), Alarm/Snooze Set, 12hr or 24hr Format Selectable, Indoor/Outdoor
Temperature Measurement, Indoor Humidity Measurement.

Specifications & Ordering information

<L1>

Spec.	Model	“H9216”
Measurement Range / Alarm Timer		-50°C + 70°C or -58°F + 158°F / 20%~99.9% R.H (No Condensation) Alarm Timer Function
Resolution		0.1° for Whole Range/ 0.1% R.H for Whole Range
Accuracy		± 1.0°C or ± 1.8°F In the Range of -20°C + 70°C or -4°F + 158°F, Otherwise more than ± 1.0°C or ± 1.8°F
Ambient Temperature		0°C~50°C or 32°F~122°F
Weight (g) / Color		165g / White
Battery / Package		2 × 1.5 Volts “AAA” Size or Equivalent / Sliding Blister with Print Card
Sensor Type		NTC & Humidity External Sensor ; 3m Sensor Wire with Probe(L22.5 × φ 5.5mm)
Dimension		120 × 110 × 25mm (LCD Size ; 72 × 61mm)
Cat. No. & Price ₩	A1.H9216	Ea

DAIHAN™ Desktop/Wallmount Dew-Point Thermo-Hygrometer of Temp · RH% · DP°C “THH7” & “THH9”, Max/Min, Hold
Programmable Danger Zone, Alarm 65dB, Time Clock, 0°C~50°C, 20~90% RH, -20°C + 70°C DP, 탁자/벽걸이용 이슬점 온습도계

■ SD카드를 통해 측정 데이터를 저장 및 재확인 가능한 온습도계

- * Simultaneously Displays Temperature and Humidity * Temperature Unit °C/°F Switchable
- * Check Humidity Condense Conditions Through Dew Point * Max/Min Function to Review Data
- * Real Time Display(Hour & Minute, AM/PM) * Programmable Danger Zone for Warning Purpose
- * Audible Beeper and Bright Red Led for Danger Zone Warning * Long Battery Life One year Via Alkaline Battery
- * Designed with Capacitance type humidity Sensor to Get Accurate and Quick Reading
- * Advanced SD Card Logging Design for Continuous Recording(DH.TH9 only)

Specifications & Ordering information

<L1>

Spec.	Model	A. “THH7”	B. “THH9”
3 × Measuring Range		(1) Temp ; 0°C~50°C(display ; -10.0°C + 70.0°C), (2) Humidity ; 20~90%RH(Display ; 0%~99.9%RH), (3) Dew Point ; -20°C + 70°C	
2 × Accuracy		(1) Temp ; ± 1°C, (2) Humidity ; ± 5%RH (at 25°C, 10~90%RH, other ± 7%RH)	
2 × Resolution		(1) Temp ; 0.1°C/°F, (2) Humidity ; 0.1%RH	
Temp Response Time		2 mins	
SD Card		-	Max. 2GB SD Card or 4GB SDHC Card
Battery		3 × AA Battery	3 × 1.5Volts AA Battery or 5V Adaptor
Dimension / Weight(g)		110 × 91 × h30mm(LCD Size ; 63 × 70mm) / 200g	
Standard Package		Meter / Battery / Manual / Paper Box	Meter / Adaptor / Manual / Paper Box / 256MB SD Card
Cat. No. & Price ₩	DH.Thh3007	Ea	DH.Thh3009 Ea



B. SD-Card Logging
Dew-point
Thermo-Hygrometer



(THH9 only)

DAIHAN™ Large Display Temp · RH% Monitor & Datalogger “THH3017”, with PDF Report, $-30 + 70^{\circ}\text{C}$

Auto-Generate a PDF & Excel Report, 대형 디스플레이 온도 및 습도 모니터 & 데이터로거



Large Display Temp · RH% Monitor & Datalogger “THH3017”

■ Supplied by DAIHAN Korea

- 온도 값을 측정하고 기록하여 사용자의 필요에 따라 Excel 또는 PDF 형식으로 보고서 생성이 가능함.
- 플러그 & 플레이 방식의 제품으로 별도의 소프트웨어 설치 필요 없이 PC에 연결만 하면 사용할 수 있음.

■ High Accuracy Humidity and Temperature Measuring.**■ Dew Point Temperature to Check Mould Growing Risk.**

* Monitoring and Logging Both can be Powered Up by Batteries and Power Adaptor.

* LED Flash & Audible Buzzer Warning Threshold are Programmable.

* Automatically Generate Statistic Report to Ease Your Work. * Salt Bottle Self-calibration Function is Integrated.

Specifications & Ordering information

<L1>

Spec.	Model	“THH3017”
3×Measure Range		(1) Temp ; $-30 + 70^{\circ}\text{C}$ / $-22 + 158^{\circ}\text{F}$ (2) Humidity ; 0.1%RH~99.9%RH, (3) Dew Point ; $-85.2^{\circ}\text{C} + 70.0^{\circ}\text{C}$
2×Accuracy		(1) Temp ; $\pm 0.5^{\circ}\text{C}$, (2) Humidity ; $\pm 3\%\text{RH}$ (at 25°C , 10~90%RH, Other $\pm 5\%\text{RH}$)
2×Resolution		(1) Temp ; $\pm 0.1^{\circ}\text{C}/^{\circ}\text{F}$, (2) Humidity ; $\pm 0.1\%\text{RH}$
Sampling Point		24,000 Readings
Sampling Interval		30 seconds, 1, 5, 10, 30, 60, 90, 120 minutes
Dimensions & Weight		120×93×h42 mm / 200g
Battery & Life Time		4×AA Battery or 9V DC Adaptor 3 Months Life Time While Alarm Function is OFF 2 Months Life Time While Alarm Function is ON
Standard Package		Logger, Battery, Manual, Paper Box, USB Cable
Cat. No. & Price ₩	DH.Thh3017	Set

DAIHAN™ Temp $^{\circ}\text{C}/^{\circ}\text{F}$ & RH% USB Datalogger with PDF Report “THH3”, Waterproof(IP65), $-30^{\circ}\text{C} + 70^{\circ}\text{C}$, 0.1~99.9% RH

Auto-Generate a PDF & Excel Report while Plugged into PC, Programmable by Users, 온습도 측정 및 Data 로거



USB Datalogger

■ Supplied by DAIHAN Korea

- * Software & Driver are not required for any operation
- * High Accuracy Thermistor and Capacitive RH Sensor
- * Programmable ; 6 Language, Sampling Interval, Start Delay, Alarm Delay, Unit, Alarm Range, Alarm Type &c.
- * Designed with Comprehensive Usage with Waterproof Function(IP65)
- * EN12830 and RoHS Compliance

Specifications & Ordering information

<L1>

Spec.	Model	“THH3”
		Temperature Humidity
Measure Range	NTC Thermistor, $-30^{\circ}\text{C} + 70^{\circ}\text{C}$, $-22^{\circ}\text{F} + 158^{\circ}\text{F}$	Capacitor, 0.1~99.9%RH
Accuracy	$\pm 0.5^{\circ}\text{C}$	$\pm 3\%$ (at 25°C , 10~90%), others $\pm 5\%$
Resolution	$0.1^{\circ}\text{C} / 0.1^{\circ}\text{F}$	0.1%RH
Battery		2×CR2032(3.0 V) Battery
Dimension / Weight(g)		110×45×h17mm / 80g
Standard Package		Logger / Manual / Paper Box
Cat. No. & Price ₩	DH.Thh3003	Ea

DAIHAN™ Temp · Humidity · Barometric Bluetooth Datalogger “THH15”, $-30^{\circ}\text{C} + 70^{\circ}\text{C}$, 0.1~99.9% RH, 300~1100hPa

Generate PDF & Excel Report Through APP, Real Time Check without Stopping, Waterproof(IP65), 온습도 및 기압 Data 로거



Bluetooth Datalogger

■ Supplied by DAIHAN Korea**■ 측정 Data를 어플리케이션을 통해 실시간으로 볼 수 있으며, 측정 후에도 결과 확인이 가능함**

- * Programming Logger & Download Data Wirelessly Through Advanced Bluetooth 4.0 Technology
- * Generate and Send Out PDF Report Through APP Report Contains Data Summary, Graph and Detailed Data
- * Please Search for “Nice Logger” in APP Store.
- * EN12830 and RoHS Compliance

Specifications & Ordering information

<L1>

Spec.	Model	“THH15”
		Temperature Humidity Barometric
Measure Range	NTC Thermistor, $-30^{\circ}\text{C} + 70^{\circ}\text{C}$, $-22^{\circ}\text{F} + 158^{\circ}\text{F}$	Capacitor, 0.1~99.9%RH 300~1100hPa 8.9~32.5inHg
Accuracy	$\pm 0.5^{\circ}\text{C}$	$\pm 3\%$ (at 25°C , 10~90%), others $\pm 5\%$ $\pm 3\text{hPa}(0\sim 60^{\circ}\text{C})$, $\pm 5\text{hPa}(-20\sim 0^{\circ}\text{C})$
Resolution	$0.1^{\circ}\text{C} / 0.1^{\circ}\text{F}$	0.1%RH 0.1hPa, 0.1inHg
Ambient Temp		$-30^{\circ}\text{C} + 70^{\circ}\text{C}$, (Humidity < 90%)
Battery		CR2(3.0 V) Battery × 1pc
Dimension / Weight(g)		70×50×h23mm (Body), 30×18(LCD) / 80g
Standard Package		Logger / Manual / Paper Box / Battery
Cat. No. & Price ₩	DH.Thh3015	Ea

DAIHAN™ WBGT Checker "THE3021", Temp 0°C ~ 50°C · 0.1~99.9%RH · WBGT 0°C ~ 50°C, with Alarm Function, 70×60×h22mm
Ideal for Outdoor & Indoor Heat Stress Measuring, 휴대용 WBGT 습구흑구 온열지수 측정기, 체감온도 측정, 열사병 예방용



Wearable Design

■ Supplied by DAIHAN Korea

■ WBGT (Wet Bulb Globe Temperature) = 습구흑구 온열지수 ; 건구온도 + 습구온도 + 흑구온도 ; 기온, 습도, 복사열, 기류 등을 종합적으로 분석해 열에 의해 인간이 받는 스트레스를 나타낸 수치.

■ 휴대용, 착용 가능한 제품으로, 야외활동을 주로하는 운동선수나 군인, 산업현장에서의 사용에 유용.

* Simultaneously Display WBGT, Temperature & RH.

* Programmable Buzzer and Red LED 4 level Danger Zone Warning

* Audible Alarm Volume Adjustable (75 dB / 60 dB / mute)

* Can be used the Tripod to Hold the Checker Vertically for Long Term use.



<L1>

Specifications & Ordering information

WBGT	Flag	Activity Intensity
< 80°F < 26.7°C	White	Normal activity. Caution should be taken.
80-84.9°F 26.7-29.4°C	Green	Discretion is required in planning intense physical activity.
85-87.9°F 29.4-31°C	Yellow	Intense activity for beginners should be curtailed.
88-89.9°F 31-32.2°C	Red	Intense exercise must be curtailed Be on high alert.
> 90°F > 32.2°C	Black	Cancel all outdoor exercise.

Spec.	Model	"THE3021"
Temp & Humidity Range (Resolution)		0°C~50°C & 0.1~99.9% RH (0.1°C & 0.1% RH)
Temp & Humidity Accuracy		±0.6°C & ±5%RH (at 25°C, 10~90%RH); others ; ±7%RH
Globe Temp & Black Ball Size		0°C~70°C & Dia. 40mm
WBGT Temp & Alarm Range		0°C~50°C & 20°C~50°C
Dimensions & Weight		70×60×h22mm (LCD Size ; 40×39mm) & 90g
Standard Package		Meter, 2×AAA Battery, Manual, Paper Box, Arm Strap
Cat. No. & Price ₩		DH.The3021 Set

DAIHAN™ Professional WBGT Monitor "THE3023", Temp 0°C ~ 50°C · 0.1~99.9%RH · WBGT 0°C ~ 50°C, with Alarm Function, 302×75×h36.7mm
Ideal for Outdoor & Indoor Heat Stress Measuring, with Backlight Indicator, 전문가용 WBGT 습구흑구 온열지수 측정기, 열사병 예방용



Pull up to Protect Sensor Design Patent Pending

■ Supplied by DAIHAN Korea

■ WBGT (Wet Bulb Globe Temperature) = 습구흑구 온열지수 ; 건구온도 + 습구온도 + 흑구온도 ; 기온, 습도, 복사열, 기류 등을 종합적으로 분석해 열에 의해 인간이 받는 스트레스를 나타낸 수치.

■ 야외활동에서 고온에 대한 적절한 노출 수준 결정을 도와줌 (열사병 예방)

■ 혹서의 작업환경이나 산업현장에서의 안전관리 목적으로 사용하기에 적합함.

* Back Light Indicator Function (Green, Yellow, Red, Flashing red) as WBGT Condition Reminder.

* Easy-to-read Big LCD Display, 6 Dust Proof Rubber Keys.

* One Button Switch to Check Air / Globe temperature & In/Out WBGT.

* Designed with Capacitance Type Humidity Sensor to get Accurate and Quick Reading.

* Hold Function to Freeze Reading and Max/Min Function Included

* Temperature unit °C/°F Switchable



Can be used with Tripod

<L1>

Specifications & Ordering information

WBGT indicator	87785 / 87782	WBGT indicator	87785 / 87782
Green backlight (26.7-29.4°C)		Red backlight (31-32.2°C)	
Yellow backlight (29.4-31°C)		Red flashing backlight (Above 32.2°C)	

Indicator Display

Spec.	Model	"THE3023"
Temp & Humidity Range(Resolution)		0°C~50°C & 0.1~99.9% RH (0.1°C & 0.1% RH)
Temp & Humidity Accuracy		±0.6°C & ±3%RH (at 25°C, 10~90%RH); others ±5%RH
Globe Temp(Inside Black Ball)		0°C~80°C
Globe Temp Acc.(Indoor & Outdoor)		±1°C at 15°C~40°C, others ±1.5°C & ±1.5°C at 15°C~40°C, others ±2°C
WBGT Temp & Alarm Range		0°C~50°C & 20°C~50°C
Dimensions & Weight		302×75×h36.7mm (LCD Size ; 50×32mm) & 300g, Black Ball ; Dia.75×h75mm
Standard Package		Meter, 4×AAA Battery, Manual, Paper Box
Cat. No. & Price ₩		DH.The3023 Set

DAIHAN™ Heat Stroke Prevention WBGT Monitor "THE3", Temp · RH% · WBGT · Alarm 65dB
Desktop/Wallmount, Programmable Danger Zone, 0°C ~ 50°C, 20~90% RH, WBGT 0°C ~ 60°C, WBGT 모니터



W.B.G.T Monitor

■ Supplied by DAIHAN Korea

■ 혹서의 작업환경이나 산업현장에서의 안전관리 목적으로 사용하기에 적합

* Checks Dangerous Heat Conditions through WBGT

* Temperature Unit °C / °F is Switchable

* Simultaneously Displays Temperature and Humidity

* Max/Min Function to Review Data

* Audible Beeper and Bright Red LED for Danger Zone Warning

* Real Time Display (Hour & Minute, AM/PM)

* Battery Low indicator

Specifications & Ordering information

<L1>

Spec.	Model	"THE3"
3× Measuring Range		(1) Temp ; 0°C~50°C (Display ; -10°C + 70°C), (2) Humidity ; 20 ~ 90%RH (Display ; 0 ~ 99.9%RH), (3) WBGT ; 0 ~ 60°C
2× Accuracy		(1) Temp ; ±1°C (2) Humidity ; ±5%RH (at 25°C, 10 ~ 90%RH, others ±7%RH)
2× Resolution		(1) Temp ; 0.1°C/0.1°F, (2) Humidity ; 0.1%RH
Temp Response Time		120 Seconds Typical
Dimension / Weight(g)		110×91×h30mm(LCD Size ; 63×70 mm) / ~200g
Standard Package		Meter / 3×AA Battery / Manual / Paper Box
Cat. No. & Price ₩		DH.The3003 Set

DAIHAN™ Non-Contact Precision Healthcare IR Thermometer “The01”, Precision 0.1 Divi.

With FND Display, Ergonomic Hand Grip, Clinical-mode(32.0℃~43.0℃/±0.3℃) & Object-mode(0℃~100℃/±1.5℃)

비접촉식 적외선 정밀 온도계, 체온계 겸용, 듀얼측정모드

2 Modes;

- Clinical-mode (32.0℃~43.0℃)
- Object-mode (0℃~100℃)

**3-Year Warranty**

Non-Contact Precision IR Thermometer with FND Display

- Supplied by DAIHAN Korea
- **PL(Product Liability) Insurance**
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 위생적인 비접촉식 측정 방식, 0.5초의 빠른 응답속도
- Non-contact Temperature Checker Designed for 32.0℃~43.0℃, 0℃~100.0℃
- Ideal for Checking Body Heat, Bath Water, Milk Temperature and Other Object.
- 다양한 온도 측정이 가능한 2종류의 측정 모드
- 2×Measuring Mode ; (1) Clinical-mode ; Ideal for Body Heat, 32.0℃~43.0℃/ ±0.3℃
- (2) Object-mode ; Ideal for Various Object, 0℃~100.0℃/ ±1.5℃
- * High-quality FND Display
- * Ergonomic Hand Grip
- * Fast Response Time ; 0.5sec
- * Temperature Units are ℉, ℃(Factory Default ; ℃)

Specifications & Ordering information

<M1>

Spec.		Model	“The01”
Measuring Range	Clinical-mode		32.0℃~43.0℃
	Object-mode		0℃~100.0℃
Resolution & Temp Accuracy	Clinical-mode		0.1℃ & ±0.3℃
	Object-mode		0.1℃ & ±1.5℃
Targeting Distance			1~5cm
Battery			2×1.5V Battery(AAA)
Dimension & Weight(g)			40×55×h168mm (FND Size ; 26×h28mm) & 108g
Standard Package			Meter, Battery, Manual
Cat. No. & Price ₩			DH.The0100 Set

DAIHAN™ 22℃~41.5℃ Non-Contact Healthcare IR Thermometer “THE19”, Precision 0.1 Divi.

Ideal for Surface Temperature Measuring, with Wrist Strap, 비접촉식 적외선 저온정밀 온도계



Incl. Soft Carry Case

Healthcare IR Thermometer, 22℃~41.5℃


- Supplied by DAIHAN Korea
- 위생적인 비접촉식 측정 방식, 0.5초의 빠른 응답속도
- Non-contact Temperature Checker Designed for 22℃~41.5℃
- Ideal for Checking Bath Water and Milk Temperature
- * Fast Response Time ; 0.5sec
- * Temperature Unit is Factory Default as ℃

**Specifications & Ordering information**

<L1>

Spec.		Model	“THE19”
Measuring Range			22℃~41.5℃
Resolution & Temp Accuracy			0.1℃ & ±0.3℃
Targeting Distance			3~5cm
Battery			1×9V Battery
Dimension & Weight(g)			44×40×h170mm (LCD Size ; 27×17mm) & 140g
Standard Package			Meter, Battery, Manual, Pouch, Wrist Strap, Paper Box
Cat. No. & Price ₩			DH.The3019 Set
Spare ; DH.The3019.1			Soft Carry Case(Pouch) 176×75×h35mm, Black for DH.The3019 Ea

DAIHAN™ -40°C +280°C Waterproof IP65 IR(Infrared)+Contact Pt-Probe Thermometer "THE15", with Laser-Pointer, MAX/MIN With 180° Rotary Piercing Probe 8cm & Timer, Data-Hold °C/°F, 접촉+비접촉 겸용 IR 온도계



- Supplied by DAIHAN Korea * Follow HACCP Principle to Check Food Safety * Certified HACCP
- * IP65 Waterproof Targeting, Washable Probe
 - * Quick to Inspect Danger Zone Temperature
 - * Big LCD with Backlight to Display Measured Temperature
 - * Flip over Probe for Internal Temperature Measurement
 - * Built-in Laser to Highlight the Area
 - * Temperature Unit C / F is Switchable
 - * Auto Power off in 20 mins Inactivity
 - * Quick and Accurate IR Response Time ; 500ms
 - * Count down Timer Ensure Food Safety Storage
 - * Max Reading Displays while in IR Scanning Status

Specifications & Ordering information

<L1>

Spec.	Model	"THE15"
2 × Measuring		(1) IR Temp ; -40°C + 280°C, (2) Probe Temp ; -40°C + 200°C
Accuracy		(1) IR Temp ; [-40°C (32°F)] ; ±1°C (2°F) ±0.1, [0°C~65°C (32°F~150°F)] ; ±1°C (2°F), >65°C (150°F)] ; ±1.5% of rdg, (2) Probe Temp ; [-5°C (23°F)] ; ±1°C (2°F), [-5°C+65°C (23°F~150°F)] ; ±0.5°C (1°F), >65°C (150°F)] ; ±1°C or 1.5% of rdg
Resolution		(1) IR Temp ; 0.1°C/°F, (2) Probe Temp ; 0.1°C/°F
D:S ratio / IR Repeatability / Emissivity		3:1 / ±1°C / 0.97
Battery		1 × 9V Battery or Adaptor
Dimension / Weight		61.5 × 35.5 × 167mm (LCD Size ; 21 × 28mm) / 140g
Standard Package		Meter / Battery / Manual / Paper Box
Cat. No. & Price ₩		DH.The3015 Set
Option ; DH.The3015.1		Soft Carry Case(Pouch) 176 × 75 × h35mm, Black for DH.The3015 Ea

DAIHAN™ -40°C +500°C Compact Gun type IR(Infrared) Thermometer "THE13", Comply with FDA Standard Class II With Ergonomic Hand Grip, °C/°F, 비접촉 적외선(IR) 온도계



- Supplied by DAIHAN Korea * Built in Laser Beam(670nm, <1mw) for Targeting.
- * Adjustable Emissivity from 0.3 to 1.0 for Metal and Non-Metal Surface Measurement.
 - * Product Complies with FDA Standard Class II
 - * Quick and accurate IR Response Time ; 500ms.

Specifications & Ordering information

<L1>

Spec.	Model	"THE13"
Measurement Range		-40°C + 500°C
Temp Accuracy		±2% or ±2°C at -20°C + 200°C, others ±3% or 3°C
Resolution		0.1°C/°F for -40°C + 110°C, 1°C/°F for others
Repeatability / Distance Ratio		±1°C / 10 : 1
Battery		1 × 9V Battery
Dimension & Weight(g)		80 × 45 × h157mm (LCD Size ; 21.8 × 31.8mm) / 140g
Standard Package		Meter / Battery / Manual / Paper Box
Cat. No. & Price ₩		DH.The3013 Set
Option ; DH.The3013.1		Soft Carry Case(Pouch), Black for DH.The3013 Ea

DAIHAN™ -40°C +816°C High Temp IR(Infrared) Thermometer "THE17", with Large LCD & Laser-Pointer With Wrist Strap & Hard Carry Case & Soft Carrying Pouch, °C/°F, 비접촉 적외선(IR) 온도계



- Supplied by DAIHAN Korea
- * Built in Laser Beam(670nm, <1mw) for Targeting.
 - * Quick and accurate IR Response Time ; 500ms.
 - * Rubber-Like Hand Grip Design.
 - * Large LCD w/ Backlight Function.
 - * Adjustable Emissivity from 0.3 to 1.0
 - * Easy Flip Battery Compartment Design.

Specifications & Ordering information

<L1>

Spec.	Model	"THE17"
Measurement Range		-40°C + 816°C
Temp Accuracy		<0°C ; ±2°C / 0°C~300°C ; ±2% or ±2°C 300°C~500°C ; ±2.5% / >500°C ; ±3%
Resolution		0.1°C/°F(<300°C), 1°C/°F(>300°C)
Distance Ratio / Battery		12 : 1 / 2 × AA Battery
Dimension & Weight(g)		195 × 134 × h50 mm (LCD Size ; 34 × 37mm) / 240g
Standard Package		Meter / Battery / Manual / Hard Carry Case / Wrist Strap / Soft Carrying Pouch
Cat. No. & Price ₩		DH.The3017 Set
Spare ; DH.The3017.1		Soft Carry Case(Pouch) 200 × 133 × h42mm, Black for DH.The3017 Ea

DAIHAN™ IP67 Waterproof Digital Thermocouple Thermometer "TC9227A", with K-Probe
 With 10 Sets Memories Data Storage, Recalling & Clearing, $-100^{\circ}\text{C} + 800^{\circ}\text{C}$, 방수/휴대용 1채널 디지털 온도계


- * Max/Min ($^{\circ}\text{C}/^{\circ}\text{F}$) Memory Recalling.
- * EL Backlight for Use in Dim Light Environment.
- * Alerting Alarm Beeps when Temperature Reaches/Beyond the Pre-Set Limit
- * Built in Holster with Flip Out Stand at the Back of the Unit.
- * Low Battery Indicator

Specifications & Ordering information

<L1>

Spec.	Model	"TC9227A"
Measurement Range		$-100^{\circ}\text{C} + 800^{\circ}\text{C}$ or $-148^{\circ}\text{F} + 1472^{\circ}\text{F}$
Resolution & Accuracy		0.1°, and $\pm 2^{\circ}\text{C}$
Water Proof & Weight		IP67 and 185g(only main unit) / Grey&Yellow
Battery		3 $\times$ 1.5Volts "AAA" Size or Equivalent
Sensor Type		K - Probe ; od. $\phi 3.5 \times \text{L}110\text{mm}$ Plastic Handle ; od. $\phi 17 \times \text{L}165\text{mm}$, Cable ; od. $\phi 4 \times \text{L}1050\text{mm}$
Dimension		54 $\times$ 165 $\times$ 25 mm (LCD Size ; 32 $\times$ 23 mm)
Cat. No. & Price ₩	A1.TC9227A	Ea

DAIHAN™ Portable Digital 2-Channel Thermocouple Thermometer "T9214T", K · J · E · R · T · N · S-Type Probe Acceptable
 With $\Delta T = T_1 - T_2$ · Max-Min · Avg Memory, Without Probe, $-200^{\circ}\text{C} + 1370^{\circ}\text{C}$, 0.1 Divi., 휴대용 2채널 디지털 온도계

■ Multiple Thermocouple (Type K, J, T, E, N, R, S) Compatible

- * 2 Channel Reading (T_1 & T_2).
- * Differential Temperature Measurement ($\Delta T = T_1 - T_2$).
- * 30 Sets Data Storage
- * 7" Large LCD with Function Indicators and Backlight.
- * Alarm Beeps when Temperature Reaches/ Beyond/ The Pre-Set Limit.
- * Audible Alarm.
- * Optional PVC Holster with Flip Out Stand at the Back of the Unit.
- * Low Battery Indicator.

Specifications & Ordering information

<L1>

Spec.	Model	"T9214T"
Measurement Range		K-Type ; $-200^{\circ}\text{C} + 1370^{\circ}\text{C}$ or $-328^{\circ}\text{F} + 2498^{\circ}\text{F}$, Clock-Timer, Alarm
Resolution		0.1° for Whole Range
Accuracy		$\pm 0.05\%$ of Reading $+0.3^{\circ}\text{C} (+0.5^{\circ}\text{F})$ Below $-100^{\circ}\text{C} (+148^{\circ}\text{F})$; Add 0.15%
Ambient Temperature		$0^{\circ}\text{C} \sim 50^{\circ}\text{C}$ or $32^{\circ}\text{F} \sim 122^{\circ}\text{F}$
Weight (g) / Color		310 / Yellow-Grey
Battery / Package		Battery Included ; 4 $\times$ 1.5Volts "AAA" Size or Equivalent / White Box
Sensor Type		Type K, J, E, T, N, R, S Thermocouple Acceptable ; Not Include
Dimension		81 $\times$ 172 $\times$ 33 mm (LCD Size ; 59 $\times$ 35 mm)
Cat. No. & Price ₩	A1.T9214T ; Without Probe	Ea

※ Option ; K-Probes $\Rightarrow$ See the 222 Page.
DAIHAN™ RS232/PC Data Analysis Thermometer "THE7", with Dual Input & Multiple-Type Probe, Compatible
 With Time Clock/Backlight, K · J · T · R · S · E-Probe in Use, $-200^{\circ}\text{C} + 1370^{\circ}\text{C}$ in K-Type, 2채널 대기종 센서 수용 온도계

■ Supplied by DAIHAN Korea

- * Non Sleep Mode Function available, Unit $^{\circ}\text{C}/^{\circ}\text{F}$ is Switchable, Super Big LCD for easy Reading.
- * RS232 interface for PC download, Data records stamped with real time clock, Backlight for dark area operation.

Specifications & Ordering information

* Hold Function to Freeze Current Reading

<L1>

Spec.	Model	"THE7"
6 $\times$ Measuring Range		(1) K-temp ; $-200^{\circ}\text{C} + 1370^{\circ}\text{C}$, (2) J-temp ; $-200^{\circ}\text{C} + 760^{\circ}\text{C}$, (3) T-temp ; $-200^{\circ}\text{C} + 390^{\circ}\text{C}$, (4) R-temp ; $0^{\circ}\text{C} \sim 1760^{\circ}\text{C}$, (5) S-temp ; $0^{\circ}\text{C} \sim 1760^{\circ}\text{C}$, (6) E-temp ; $-200^{\circ}\text{C} + 736^{\circ}\text{C}$
Accuracy		(1) K- / (2) J- / (3) T- / (6) E-temp ; $\pm (0.1\% \text{ rdg} + 0.7^{\circ}\text{C})$ (4) R- / (5) S-temp ; $\pm (0.3\% \text{ rdg} + 0.7^{\circ}\text{C})$
Resolution		(1) K-temp ; 0.1° for $-200^{\circ}\text{C} + 650^{\circ}\text{C}$, 1° for $640^{\circ}\text{C} \sim 1370^{\circ}\text{C}$, (2) J-temp ; 0.1° for $-200^{\circ}\text{C} + 500^{\circ}\text{C}$, 1° for $490^{\circ}\text{C} \sim 760^{\circ}\text{C}$, (3) T-temp ; 0.1° for $-200^{\circ}\text{C} + 390^{\circ}\text{C}$, etc.
Battery		9V $\times$ 1pc or Adaptor
Dimension / Weight(g)		182 $\times$ 72 $\times$ h30 mm (LCD Size ; 30 $\times$ 51 mm) / 150g
Standard Package		Meter / K Bead Probe / Battery / Manual
Cat. No. & Price ₩	DH.The3007	Set
Spare ; DH.The3007.1	① 4m K-Bead Thermocouple	Ea
Option ; DH.Gas3013.1	② USB Cable & Software	Ea
DH.Gas3013.2	③ RS232 Cable & Software	Ea
DH.The3007.2	④ Rubber Holster, Blue	Ea

DAIHAN™ Wallmount SD-Card Logging 4-Channel K-type Thermometer "THE9", -200℃ +1370℃

Ultimated Auto-Logging Memories in SD-Card, with 2 Bead K-Probe · Alarm · Time Clock, 벽걸이 4채널 SD로거/온도계



► No-Bothering Computer Setting! Auto-Log Anytime/Anywhere!

■ Supplied by DAIHAN Korea

- * Super Big LCD Portable / Wallmount SD-card Logger Measures Temperature via K-type Thermometer
- * Hold Function Freezers Current Readings, Checking Max / Min Value from Power On
- * One Button Press to Display T1-T2, Display Temperature Change from Baseline
- * Blue Backlight
- * 12/24 hour Timer

Specifications & Ordering information

<L1>

Spec.	Model	"THE9"
K temp range / Resolution		-200℃ +1370℃, -328°F + 2498°F / 0.1℃, 0.1°F
K temp accuracy		±(0.3% rdg + 1℃)
Sampling rate		Programmable from 1 second up
Battery		4 × AA Battery 9V adaptor ; Low Battery Indicator
Dimension / Weight(g)		152 × 100 × h39mm(LCD Size ; 47 × 104 mm) / 300g
Standard Package		Meter / 2 × Bead Probe / Battery / Manual / SD-card / Plain Box
Cat. No. & Price ₩	DH.The3009	Set
Option ; DH.The3009.1	① 1m K-Bead Thermocouple	Ea
DH.The3009.2	② Universal Plug 9V Adaptor	Ea

DAIHAN™ Foldable Portable Digital Multi-Thermometer-Set "T9218", with 4 K-type Probe

With 180° Moving Probe Connector, Timer(Count Up), -200℃ +1370℃, 0.1 Divi.

With 4 K-Probe for Air · Immers. · Penetrat. · Surface, and Carrying Case, 휴대용 접이식 디지털 온도계 세트



■ Ideal for Multi-use with Multi-Function in Lab & Industry

- * Unique Feature of Folding Probe
- * Probe Connector is Adjustable Through 45°, 90°, 135°, 180°
- * Comprehensive and Variety Use
- * Four Probes Included for Different Measurement (e.g. Surface, Semi-Solid, Air)
- * Extra Carrying Box for Easy and Protective Storage
- * Four Probes Included
- * Large LCD
- * Easy to Operate
- * CE Approval
- * Battery Included
- * High Degree of Accuracy

Function ; On/Off, Hold/Test, °C/°F, Max/Min (°C/°F) Memory Recalling, and Timer.

Specifications & Ordering information

<L1>

Spec.	Model	"T9218"
Measurement Range		-200℃ +1370℃ or -328°F +2498°F, Timer
Resolution		0.1° for Whole Range
Accuracy		±1.0℃ or ±1.8°F In the Range of -49℃ +740℃ or -56.2°F +1364°F, Otherwise more than ±2.0℃ or ±4°F
Ambient Temperature		0℃ ~50℃ or 32°F ~122°F
Weight (g) / Color		80 / Black
Battery / Package		2 × 3Volts "CR 2032" Size or Equivalent / Black Carrying Case
Sensor Type		Included 4 × Type K Thermocouple
Dimension		153 × 50 × 19 mm (LCD Size ; 50 × 22 mm)
Cat. No. & Price ₩	A1.T9218	Ea
Spare ; A1.TK069T	① Surface Probe, K-type	Ea
; A1.TK066A	② Airflow Probe, K-type	Ea
; A1.TK1480	③ Immersion Probe, Round, K-type	Ea
; A1.TK133S	④ Penetration Probe, K-type	Ea

Foldable/Portable Digital Multi-Thermometer-Set



Digital Temp Controller/Displayer for K-type Probes

With Female Minisocket for Miniplug, -50℃ +400℃, 디지털 온도측정/조절계

■ Ideal for Temp Control of Reactor-system

- * Ideal for Temperature Monitor / Controller for Reactor System
- * With K-type Probe Connector & Electronic Cable 2m
- * Jumbo LCD with Digit Height 18mm
- * Dimension ; 90 × 70 × h50mm
- * 100~240VAC Free
- * High Accuracy

* Note ; The Controller Can be used in other Cable & Plug System.

SL.The5011 Digital Temp Controller/Displayer-Set, -50℃ +400℃, K-type

<L1>

Ea. Price, ₩

Cowie® K-type PTFE Encapsulated Stainless-steel Thermocouple Probe, Can be Bent, $-180^{\circ}\text{C} + 280^{\circ}\text{C}$
 With Miniplug & PFA Insulated 1m & 2m Cable, Probe Body L100~900mm, PTFE 내약품/내부식 청결형 온도 프로브



■ **Made in UK.** ■ **Ideal for Reactor ; Especially Long Body Probes**

- * Outstanding Chemical Resistance, No-Contamination of Working Media, and Pressure Resistant
- * Fully Pure PTFE Encapsulated Stainless-steel Core Probe for Rigidity, But can be Bent
- * with Push-On Ferrule
- * Response ; 50% in 6~14 seconds, $-180^{\circ}\text{C} + 280^{\circ}\text{C}$

(1) With L100~L700mm Probe Body and 1m Cable

CW.600.101.1	PTFE Capsulated K-probe, Body- $\phi 6 \times L100$, Ferrule-L35mm, Total -L135mm
CW.600.201.1	- $\phi 7 \times L200$, -L40mm, -L240mm
CW.600.301.1	- $\phi 8 \times L300$, -L40mm, -L340mm
CW.600.401.1	- $\phi 8 \times L400$, -L50mm, -L450mm
CW.600.501.1	- $\phi 9 \times L500$, -L60mm, -L560mm
CW.600.601.1	- $\phi 9 \times L600$, -L60mm, -L660mm
CW.600.701.1	- $\phi 9 \times L700$, -L45mm, -L745mm

(2) With L900mm Probe Body and 2m Cable

CW.600.901.1.2	PTFE Capsulated K-probe, Body- $\phi 9 \times L900$, Ferrule-L45mm, Total-L945mm
-----------------------	---

<L1>

Ea. Price, ₩

K-type Thermocouple Probe, for Wide-Range Temp $-100^{\circ}\text{C} + 1000^{\circ}\text{C}$
 With Common Use Plug "Miniature 2 Flat" and Cable, "K-형" 온도계 & 각종 열기구 온도 프로브

- **K-Type Sensor Probes for Different Applications (Penetration Probe, Surface Probe, Pipe Clamp)** ■ **Cable ; od. $\phi 3 \times L1020\text{mm}$**
- **Miniature 2 Flat (L14 $\times$ Thick 1mm) Plug Type** ■ **Handle ; Thermo Plastic od. $\phi 17 \times L170\text{mm}$**

(1) For Liquid & Air, Stainless-steel Probe, Resin Handle with Cable/Plug, $-60^{\circ}\text{C} + 700^{\circ}\text{C}$



SL.The7193	Liquid/ "K", $-40^{\circ}\text{C} + 700^{\circ}\text{C}$, $\phi 1.0 \times L306\text{mm}$, 1m Cable/Plug
DH.WHP503040	, $-60^{\circ}\text{C} + 600^{\circ}\text{C}$, $\phi 4.8 \times L250\text{mm}$, 1.8m Cable/Plug
DH.WHP503050	, $-60^{\circ}\text{C} + 500^{\circ}\text{C}$, $\phi 4.8 \times L380\text{mm}$, 1.8m Cable/Plug

<L2>

Ea. Price, ₩

(2) For Semi-Solid Penetration-type, Stainless-steel Probe, Resin Handle with Cable/Plug, $-100^{\circ}\text{C} + 1000^{\circ}\text{C}$




A1.TK175S	Penetration Probe "K", $-100^{\circ}\text{C} + 800^{\circ}\text{C}$, $\phi 3.2 \times L150\text{mm}$, 1m Cable/Plug
A1.TK225S	, $-100^{\circ}\text{C} + 1000^{\circ}\text{C}$, $\phi 6.4 \times L200\text{mm}$, 1m Cable/Plug

<L1>

Ea. Price, ₩


(3) For Surface / Contact, Stainless-steel Probe, Resin Handle with Cable/Plug, $-100^{\circ}\text{C} + 400^{\circ}\text{C}$



A1.TK004T	(a) Surface Probe "K", Head ; $\phi 15 \times t1.5$, Stem ; $\phi 3.5 \times L180\text{mm}$, $-100^{\circ}\text{C} + 400^{\circ}\text{C}$, 1m Cable/Plug
A1.TK033T	(b) ; $\phi 16 \times t1.4$, Stem ; $\phi 6.5 \times L215\text{mm}$, $-100^{\circ}\text{C} + 400^{\circ}\text{C}$, 1m Cable/Plug

Ea. Price, ₩

(4) Pipe Clamping Probe, with Cable Plug, $-100^{\circ}\text{C} + 120^{\circ}\text{C}$



A1.TK0063	Pipe Clamp Probe, Pipe / Tubing, $\phi 6 \sim 22\text{mm}$ Grip $-100^{\circ}\text{C} + 120^{\circ}\text{C}$
------------------	--

Ea. Price, ₩

DAIHAN™ IP67 Waterproof Digital Thermocouple Thermometer "T9234", with 1 NTC-Probe

With 10 Sets Memories Data Storage, Recalling & Clearing, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$, 방수/휴대용 1채널 디지털 온도계



- * Max/Min ($^{\circ}\text{C}/^{\circ}\text{F}$) Memory Recalling.
- * Alerting Alarm Beeps when Temperature Reaches/Beyond the Pre-set Limit
- * EL Backlight for Use in Dim Light Environment.
- * Supplied with Sensor Probe and Plastic Handle Assembled Together
- * PVC Holster with Flip Out Stand at the Back of the Unit.
- * 1 Second Normal Temperature Sensing Time.
- * Optional One Hour Auto Shut Off
- * CE Approval.
- * Low Battery Indicator.

Specifications & Ordering information

<L1>

Spec.	Model	"T9234"
Measurement Range		$-50^{\circ}\text{C} + 300^{\circ}\text{C}$ or $-58^{\circ}\text{F} + 572^{\circ}\text{F}$
Resolution		0.1 $^{\circ}\text{C}$ for $-19.9^{\circ}\text{C} + 199.9^{\circ}\text{C}$, Otherwise 1 $^{\circ}\text{C}$
Accuracy		$\pm 1.0^{\circ}\text{C}$ or $\pm 1.8^{\circ}\text{F}$ In the Range of $-30^{\circ}\text{C} + 250^{\circ}\text{C}$ or $-22^{\circ}\text{F} + 482^{\circ}\text{F}$ Otherwise more than $\pm 2.0^{\circ}\text{C}$ or $\pm 4^{\circ}\text{F}$
Ambient Temperature		0 $^{\circ}\text{C}$ ~ 50 $^{\circ}\text{C}$ or 32 $^{\circ}\text{F}$ ~ 122 $^{\circ}\text{F}$
Weight(g) / Color / Water Proof		280 / Grey · Yellow / IP 67
Battery / Package		2 x 1.5Volts "AAA" Size or Equivalent / White Box
Sensor Type		NTC Probe ; od. $\phi 3.5 \times \text{L}110\text{mm}$ Plastic Handle ; od. $\phi 17 \times \text{L}120\text{mm}$, Cable ; od. $\phi 4 \times \text{L}1450\text{mm}$
Dimension		77 x 172 x 46mm (LCD Size ; 51 x 32 mm)
Cat. No. & Price ₩	A1.T9234	Ea

DAIHAN™ IP67 Waterproof Digital 2-Channel Thermocouple Thermometer "T9234T", with 2 NTC-Probe

With 10 Sets Memories Data Storage, Recalling & Clearing, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$, 방수/휴대용 2채널 디지털 온도계



- * Convenient Function for Calculating Differential Temp-Measurement ($\Delta T = T_1 - T_2$).
- * Alerting Alarm Beeps when Temperature Reaches/Beyond the Pre-Set Limit
- * 1 Second Normal Temperature Sensing Time.
- * Low Battery Indicator.
- * CE Approval.
- * EL Backlight for Use in Dim Light Environment.
- * Universal Kit Set with Two Stainless-steel Probes Together.

Specifications & Ordering information

<L1>

Spec.	Model	"T9234T"
Measurement Range		$-50^{\circ}\text{C} + 300^{\circ}\text{C}$ or $-58^{\circ}\text{F} + 572^{\circ}\text{F}$
Resolution		0.1 $^{\circ}\text{C}$ for $-19.9^{\circ}\text{C} + 199.9^{\circ}\text{C}$, Otherwise 1 $^{\circ}\text{C}$
Accuracy		$\pm 1.0^{\circ}\text{C}$ or $\pm 1.8^{\circ}\text{F}$ In the Range of $-30^{\circ}\text{C} + 250^{\circ}\text{C}$ or $-22^{\circ}\text{F} + 482^{\circ}\text{F}$ Otherwise more than $\pm 2.0^{\circ}\text{C}$ or $\pm 4^{\circ}\text{F}$
Ambient Temperature		0 $^{\circ}\text{C}$ ~ 50 $^{\circ}\text{C}$ or 32 $^{\circ}\text{F}$ ~ 122 $^{\circ}\text{F}$
Weight(g) / Color / Water Proof		280 / Grey · Yellow / IP 67
Battery / Package		2 x 1.5Volts "AAA" Size or Equivalent / White Box
Sensor Type		NTC Probe ; od. $\phi 3.5 \times \text{L}110\text{mm}$ Plastic Handle ; od. $\phi 17 \times \text{L}120\text{mm}$, Cable ; od. $\phi 4 \times \text{L}1450\text{mm}$
Dimension		77 x 172 x 46 mm (LCD Size ; 51 x 32 mm)
Cat. No. & Price ₩	A1.T9234T	Ea

DAIHAN™ NTC-type Probe, for NTC-type Thermocouple Thermometer

- (1) Penetration Probe ; Suitable Semi-Solid, Liquid, Air-Temperature.
- (2) Surface Probe ; Superior Contact, Plate-Temperature.
- (3) Robust Piercing Probe ; Stronger Penetration Force, Rapid Response Time, Suitable for Frozen Food.

<L1>

Cat. No.	Dimension	Temperature $^{\circ}\text{C} / ^{\circ}\text{F} / \text{Others}$	Ea. Price, ₩
(1) A1.TN0300	4Pin(Pole) Plug($\phi 38 \times \text{L}39$ mm), Cable ; od. $\phi 4 \times \text{L}1450$ mm with Protective Probe Cap Handle ; Thermo Plastic od. $\phi 17 \times \text{L}120$ mm od. $\phi 5.0 \times \text{L}275$ mm	Penetration, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$ / $-58^{\circ}\text{F} + 572^{\circ}\text{F}$, Very Sharp Tip	
(2) A1.TN133T	4Pin(Pole) Plug($\phi 38 \times \text{L}39$ mm), Cable ; od. $\phi 4 \times \text{L}1450$ mm Handle ; Thermo Plastic od. $\phi 17 \times \text{L}120$ mm od. $\phi 3.5 \times \text{L}110$ mm Surface Probe $\phi 15\text{mm} \times \text{thick } 1.5\text{mm}$	Surface, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$ / $-58^{\circ}\text{F} + 572^{\circ}\text{F}$, Thin/Flat Disk	
(3) A1.TN332C	4Pin(Pole) Plug($\phi 38 \times \text{L}39$ mm), Cable ; od. $\phi 4 \times \text{L}1450$ mm Handle ; Thermo Plastic od. $\phi 17 \times \text{L}120$ mm od. $\phi 5.0/3.5 \times \text{L}110$ mm Piercing-type Probe ; od. $\phi 3.5$ mm od. $\phi 5.0$ mm With Protective Probe Cap	Robust Piercing, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$ / $-58^{\circ}\text{F} + 572^{\circ}\text{F}$, Stronger Penetration Force	



DAIHAN™ Foldable Handy Digital Multi-Thermometer Set "T9306"

With 4 NTC-Probe & 180° Moving Probe, -50℃ + 300℃, 휴대용 접이식 온도계 세트

- 4가지 타입의 NTC 온도 프로브를 포함하고 있는 접이식 온도계 세트 * CE Approval.
- * Stainless-steel Probe Adjusts From 0° to 180°. * Jumbo LCD with Digit Height 18mm.
- * Low Battery Indicator. * Black Carrying Case with EVA Foam for Protective and Convenient Carriage.
- * 1 Second Normal Temperature Sensing Time.

Specifications & Ordering information

<L1>

Spec.	Model	"T9306"
Measurement Range		-50℃ + 300℃ or -58°F + 572°F
Resolution		0.1℃ for -19.9℃ + 199.9℃, Otherwise 1℃
Accuracy		±1.0℃ or ±1.8°F In the Range of -30℃ + 250℃ or -22°F + 482°F, Otherwise more than ±2.0℃ or ±4°F
Ambient Temperature for Plastic Case		0℃~50℃ or 32°F~122°F
Sensor Type		Included 4 Type NTC Thermocouple
Battery		1×1.5Volts "A76" Size (LR44, G13) or Equivalent
Color		Grey
Dimension & Weight		154×41×17mm (LCD Size ; 50×23mm) & 65g
Package		Black Carrying Case with EVA foam
Cat. No. & Price ₩	A1.T9306	Ea

DAIHAN™ Module Handy Solar Digital Thermometer "T9277B", Max/Min

With 2 of Internal Sensor & 1m External Cable/NTC-Probe, -50℃ + 300℃, 0.1/1.0℃ Divi., 태양광 충전 온도계



- 태양광 충전을 통해 작동 가능한 휴대용 온도계
- * 1m Wire with Stainless-steel Probe.
- * Strong Magnet for Easy Attachment.
- Unique Design of Using Solar Energy
- * 1 Second Normal Temperature Sensing Time.
- * CE Approval.


Specifications & Ordering information

<L1>

Spec.	Model	"T9277B"
Measurement Range		-50℃ + 300℃ or -58°F + 572°F
Resolution		0.1℃ for -19.9℃ + 199.9℃, Otherwise 1℃
Accuracy		±1.0℃ or ±1.8°F In the Range of -30℃ + 250℃ or -22°F + 482°F, Otherwise more than ±2.0℃ or ±4°F
Ambient Temperature		0℃~50℃ or 32°F~122°F
Sensor Type		External Sharp NTC-Probe ; od. φ3.2×L46mm with 1m Cable
Battery		1×1.5Volts "A76" Size (LR44, G13) or Equivalent
Package / Color		Sliding Blister with Print Card / Green
Dimension & Weight		103×43×16mm (LCD Size ; 41×15mm) & 70g
Cat. No. & Price ₩	A1.T9277B	Ea

DAIHAN™ Pocket Clip Digital Thermometer "T9307", Pen-Shape, L.125mm NTC-Probe Type

With Protective Plastic Sleeve/Pocket Clip, -50℃ + 300℃, 0.1/1.0℃ Divi., 포켓형 디지털 온도계(포켓 클립형)



- * 4mm Diameter of Stainless-steel Probe which Good for Solid Measurement.
- * Instant and Precision Temperature Display.
- * Large LCD Display with Digit Height 10mm.
- * 1 Second Normal Temperature Sensing Time.
- * CE Approval.
- * Pen-Shaped Design

Specifications & Ordering information

<L1>

Spec.	Model	"T9307"
Measurement Range		-50℃ + 300℃ or -58°F + 572°F
Resolution		0.1℃ for -19.9℃ + 199.9℃, Otherwise 1℃
Accuracy		±1.0℃ or ±1.8°F In the Range of -30℃ + 250℃ or -22°F + 482°F, Otherwise more than ±2.0℃ or ±4°F
Ambient Temperature		0℃~50℃ or 32°F~122°F
Weight (g) / Color		35 / Black
Battery / Package		1×1.5 Volts "A76" Size (LR44, G13) or Equivalent / Sliding Blister with Print Card
Sensor Type		NTC-Probe ; od. φ4×L125mm
Dimension		215×27×27mm (LCD Size ; 35×11mm)
Cat. No. & Price ₩	A1.T9307	Ea

DAIHAN™ Metal Housing Waterproof (IP67) Thermometer “ST9341C”, with Z-Clip & Protective Plastic Sleeve
 Ideal for Kitchen or Lab, Max/Min, with NTC-Probe, $-50^{\circ}\text{C} + 150^{\circ}\text{C}$, 메탈 디지털 방수 온도계(조절 가능 Z클립)



Metal Housing Waterproof(IP67) Thermometer

- * Unique Design of Metal Pulling Housing
- * 2 Points Block Design to Prevent Product Rolling
- * Main Unit Heat Resistance to 90°C
- * Waterproof IP67 Rated
- * Dishwasher Safe
- ▶ Z-Clip ; You can Adjust the Immersion Depth

Specifications & Ordering information

<L1>

Spec.	Model	“ST9341C”
Measurement Range		$-50^{\circ}\text{C} + 150^{\circ}\text{C}$ or $-58^{\circ}\text{F} + 302^{\circ}\text{F}$
Resolution		0.1°C
Accuracy		$\pm 1^{\circ}\text{C}$
Ambient Temperature		$0^{\circ}\text{C} \sim 50^{\circ}\text{C}$ or $32^{\circ}\text{F} \sim 122^{\circ}\text{F}$
Weight (g) / Color		30(g) / Silver
Battery		1 × 1.5 Volts “393” Size or Equivalent
Sensor Type		NTC-Probe ; od. $\phi 4 \times \text{L}122\text{mm}$
Dimension		$146 \times 37 \times 37 \text{ mm}$ (LCD Size ; $24 \times 12\text{mm}$)
Cat. No. & Price ₩	A1.ST9341C	Ea

DAIHAN™ Waterproof Pocket Digital Thermometer “T9231C”, Max/Min, Extension Handle, with NTC-Probe
 With Protective Sleeve/Pocket Clip, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$, $0.1/1.0^{\circ}\text{C}$ Divi., 포켓형 디지털 방수 온도계(손잡이 가변형)



“T9231C”

Waterproof / Pocket Digital Thermometer

- * Pen-Shaped Design.
- * Ingenious Design of Using Probe Sleeve Vent Hole as a Holder
- * Connecting Probe Sleeve to Head of Unit as a Handle Extension.
- * Press Button Membrane Keypad Function and Display.
- * Waterproof ; IP 67.
- * Convenient Use of Hanger.
- * CE Approval.
- * 1 Second Normal Temperature Sensing Time.

Specifications & Ordering information

<L1>

Spec.	Model	“T9231C”
Measurement Range		$-50^{\circ}\text{C} + 300^{\circ}\text{C}$ or $-58^{\circ}\text{F} + 572^{\circ}\text{F}$
Resolution		0.1°C for $-19.9^{\circ}\text{C} + 199.9^{\circ}\text{C}$, Otherwise 1°C
Accuracy		$\pm 1.0^{\circ}\text{C}$ or $\pm 1.8^{\circ}\text{F}$ In the Range of $-30^{\circ}\text{C} + 250^{\circ}\text{C}$ or $-22^{\circ}\text{F} + 482^{\circ}\text{F}$, Otherwise more than $\pm 2.0^{\circ}\text{C}$ or $\pm 4^{\circ}\text{F}$
Ambient Temperature		$0^{\circ}\text{C} \sim 50^{\circ}\text{C}$ or $32^{\circ}\text{F} \sim 122^{\circ}\text{F}$
Weight (g) / Color		40 / Black
Battery / Package		1 × 1.5 Volts “A76” Size (LR44, G13) or Equivalent / Sliding Blister with Print Card
Sensor Type		NTC-Probe ; od. $\phi 4 \times \text{L}115\text{mm}$
Waterproof		IP 67
Dimension		$206 \times 13 \times 19\text{mm}$ (LCD Size ; $25 \times 8\text{mm}$)
Cat. No. & Price ₩	A1.T9231C	Ea

DAIHAN™ Handy Waterproof Multi-Thermometer "T9304C", with Rubber Coated-Housing/-Hanger
 With Detachable $\phi 5 \times L280\text{mm}$ Long NTC-Probe, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$, $0.1/1.0^{\circ}\text{C}$ Division, 다용도 방수 온도계(분리형 온도센서)



- * Unique Design of Probe Changeable Unit.
- * 1 Second Normal Temperature Sensing Time.
- * Supplied with Longer and Broader ($280\text{mm} \times \phi 5\text{mm}$) Probe Together.
- * Waterproof Function ; IP 67.
- * Rubber Coated Housing and Hanger.
- * Flexible for Multi Usage.
- * CE Approval.

Specifications & Ordering information

<L1>

Spec.	Model	"T9304C"
Measurement Range		$-50^{\circ}\text{C} + 300^{\circ}\text{C}$ or $-58^{\circ}\text{F} + 572^{\circ}\text{F}$
Resolution		0.1°C for $-19.9^{\circ}\text{C} + 199.9^{\circ}\text{C}$, Otherwise 1°C
Accuracy		$\pm 1.0^{\circ}\text{C}$ or $\pm 1.8^{\circ}\text{F}$ In the Range of $-30^{\circ}\text{C} + 250^{\circ}\text{C}$ or $-22^{\circ}\text{F} + 482^{\circ}\text{F}$, Otherwise more than $\pm 2.0^{\circ}\text{C}$ or $\pm 4^{\circ}\text{F}$
Ambient Temperature		$0^{\circ}\text{C} \sim 50^{\circ}\text{C}$ or $32^{\circ}\text{F} \sim 122^{\circ}\text{F}$
Weight (g) / Color / Waterproof		275 / Grey · Yellow, IP 67
Battery / Package		1 × 1.5Volts "A76" Size (LR44, g13) or Equivalent / White Box
Sensor Type		NTC-Probe ; od. $\phi 6 \times L275\text{mm}$ (with Connection 325mm)
Dimension		59 × 127 × 22mm (LCD Size ; 39 × 35mm)
Cat. No. & Price ₩	A1.T9304C	Ea

DAIHAN™ Waterproof Handy Multi-Thermometer "T9314C", with Adjustable Hanger & NTC-Probe, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$
 With Jumbo LCD, Max/Min, High-/Low-Alarm, Jumbo LCD, $0.1/1.0^{\circ}\text{C}$ Divi., 포켓형 디지털 방수 온도계(조절 가능 행거)



- * Designed for Comprehensive of Use.
- * Rubber Coated Adjustable Hanger.
- * Press Button Membrane Keypad Function and Display.
- * Alarm Beeps when Temperature Reaches/ Beyond the Pre-Set Unit.
- * 1 Second Normal Temp Sensing Time.
- * Jumbo LCD Display with Digit Height 18mm.

Specifications & Ordering information

<L1>

Spec.	Model	"T9314C"
Measurement Range		$-50^{\circ}\text{C} + 300^{\circ}\text{C}$ or $-58^{\circ}\text{F} + 572^{\circ}\text{F}$
Resolution		0.1°C for $-19.9^{\circ}\text{C} + 199.9^{\circ}\text{C}$, Otherwise 1°C
Accuracy		$\pm 1.0^{\circ}\text{C}$ or $\pm 1.8^{\circ}\text{F}$ In the Range of $-30^{\circ}\text{C} + 250^{\circ}\text{C}$ or $-22^{\circ}\text{F} + 482^{\circ}\text{F}$, Otherwise more than $\pm 2.0^{\circ}\text{C}$ or $\pm 4^{\circ}\text{F}$
Ambient Temperature		$0^{\circ}\text{C} \sim 50^{\circ}\text{C}$ or $32^{\circ}\text{F} \sim 122^{\circ}\text{F}$
Battery / Package / Waterproof		1 × 1.5V "A67" Size(LR44, G13) or Equivalent / Sliding Blister with Print Card, IP 67
Sensor Type		NTC-Probe ; od. $\phi 4 \times L115\text{mm}$
Dimension & Weight		210 × 53 × 18mm (LCD Size ; 32 × 30mm), 70g
Cat. No. & Price ₩	A1.T9314C	Ea

DAIHAN™ Waterproof Handy Corkscrew Heavy-Duty Thermometer "T9311C", for Frozen Food & Solid
 With $\phi 8\text{mm}$ Corkscrew NTC-Probe, Protective Sleeve, $-50^{\circ}\text{C} + 300^{\circ}\text{C}$, $0.1/1.0^{\circ}\text{C}$ Divi., 나사형 방수 온도계



- * Rugged Heavy Duty Design.
- * Corkscrew Probe with 8mm Diameter.
- * 1 Second Normal Temperature Sensing Time.
- * CE Approval.
- * Protective Plastic Probe Sleeve with Pocket Clip
- * Waterproof ; IP 67.

Specifications & Ordering information

<L1>

Spec.	Model	"T9311C"
Measurement Range		$-50^{\circ}\text{C} + 300^{\circ}\text{C}$ or $-58^{\circ}\text{F} + 572^{\circ}\text{F}$
Resolution		0.1°C for $-19.9^{\circ}\text{C} + 199.9^{\circ}\text{C}$, Otherwise 1°C
Accuracy		$\pm 1.0^{\circ}\text{C}$ or $\pm 1.8^{\circ}\text{F}$ In the Range of $-30^{\circ}\text{C} + 250^{\circ}\text{C}$ or $-22^{\circ}\text{F} + 482^{\circ}\text{F}$, Otherwise more than $\pm 2.0^{\circ}\text{C}$ or $\pm 4^{\circ}\text{F}$
Ambient Temperature		$0^{\circ}\text{C} \sim 50^{\circ}\text{C}$ or $32^{\circ}\text{F} \sim 122^{\circ}\text{F}$
Weight (g) / Color / Waterproof		75 / Ivory White, IP 67
Battery / Package		1 × 1.5Volts "A76" Size (LR44, G13) or Equivalent / Sliding Blister with Print Card
Sensor Type		NTC-Probe ; od. $\phi 8 \times L100\text{mm}$
Dimension		161 × 100 × 30mm (LCD Size ; 35 × 11mm)
Cat. No. & Price ₩	A1.T9311C	Ea

Alla® Folding Oven Thermometer "91000.003F", with 1m Stainless Cable & NTC-Probe, Timer · Alarm · Clock, Up to 300°C/572°F
 With 1m Stainless Cable & Probe, "Cable-Stay-in Oven", 1° - Display, "23hr-59min-59sec" Timer, 접이식 오븐 온도계



■ **Supplied by Alla® in France**

- * Useable Thermometer Acts as Great Helper in Kitchen.
- * Monitoring of Internal Temperature by Probe
- * Comprehensive Multi-Functional Device for Temperature and Alarm Checking.
- * With Waterproof Probe (IP 67) and Supplied with 1m Stainless-steel Cable
- * Battery & Protective Plastic Probe Cap Included
- * Magnet at the Back for Easy Attachment
- * Maximum Count ; 23 : 59 : 59 (hr : min. : sec)
- * In Blister Packaging

Function ; On/Off, °C/°F, Timer, Clock, Alarm, High/Low Alarm

<L1>

Spec.	Model	"91000.003F"
Measurement Range		0°C~300°C(1°C) / 32°F~572°F
Resolution		1°C
Accuracy		$\pm 1^\circ\text{C} / 1.8^\circ\text{F}$ between 0°C~150°C / 32°F~302°F $\pm 3^\circ\text{C} / 5.4^\circ\text{F}$ between 150°C~300°C / 302°F~572°F
Battery		2 × 1.5 Volts "AAA" Size
Dimension (mm)		Folding Screen ; 125 × 70 × 10mm (LCD Size ; 50 × 35mm) Probe & Cable ; 150mm × ϕ 4mm Probe & 100cm Cable
Cat. No. & Price ₩	AL.91000.003F	Ea

CE



DAIHAN™ Digital Soil Thermometer "T9233C", with 0~180° Moving Head
 With Large LCD Display, Pocket-type, -50°C + 300°C, 0.1°C, 토양 온도계, 180° 무빙헤드

- * Unique Moving Head for Different Viewing Angle.
- * Large LCD Display with Digit Height 10mm.
- * Reduced Tip for Stronger Penetration Force.
- * Supplied with Protective Probe Cap.
- * 1 Second Normal Temperature Sensing Time.
- * CE Approval.

<L1>

Spec.	Model	"T9233C"
Dimension		191 × 53 × 53 mm (LCD Size ; 37 × 13 mm)
Measurement Range		-50°C + 300°C or -58°F + 572°F
Resolution		0.1° for -19.9° + 199.9°, Otherwise 1°
Accuracy		$\pm 1.0^\circ\text{C}$ or $\pm 1.8^\circ\text{F}$ In the Range of -30°C + 250°C or -22°F + 482°F, Otherwise more than $\pm 2.0^\circ\text{C}$ or $\pm 4^\circ\text{F}$
Ambient Temperature		0°C~50°C or 32°F~122°F
Weight (g) / Color		40 / Black
Battery / Package		1 × 1.5 Volts "A76" Size (LR44, G13) or Equivalent / Sliding Blister with Print Card
Sensor Type		NTC-Probe ; od. $\phi 5/3 \times L120\text{mm}$
Cat. No. & Price ₩	A1.T9233C	Ea

DAIHAN™ Jumbo 2-Row Display Indoor/Outdoor Digital Thermometer "T9239" & "T9267D", with NTC-Probe & 3m Cable, Max/Min
 With Flip Out Stand · Hanging Hole · Magnet for Attachment, -50°C + 70°C, 0.1°C Divi., 실내/외용 디지털 온도계(2열 표시)



- * Ideal for Office, Home, Laboratory, etc.
- * Dual Row Display for Easy Reading.
- * Indoor/Outdoor Temperature Display.
- * Pocket Size with Large LCD Display.
- * Flip Out Stand, Hanging Hole and Strong Magnet at the Back for Easy Attachment.
- * 3 Meters Sensor Wire with Probe for Outdoor Temperature Measurement.
- * CE Approval.

<L1>

Spec.	Model	"T9239"	"T9267D"
LCD Size		LCD Size ; 30 × 40 mm	
Measurement Range		-50°C + 70°C or -58°F + 158°F	
Resolution		0.1 for Whole Range	
Accuracy		$\pm 1.0^\circ\text{C}$ or $\pm 1.8^\circ\text{F}$ In the Range of -20°C + 70°C or -4°F + 158°F Otherwise more than $\pm 1.0^\circ\text{C}$ or $\pm 1.8^\circ\text{F}$	
Ambient Temperature		0°C~50°C or 32°F~122°F	
Weight (g) / Color		70(g) / White	82(g) / White
Battery		1 × 1.5 Volts "A76" Size (LR44, G13) or Equivalent	
Dimension		70 × 70 × 17mm	67 × 65 × 22mm
Cat. No. & Price ₩	A1.T9239	Ea	A1.T9267D Ea

DAIHAN™ $-100^{\circ}\text{C} + 300^{\circ}\text{C}$ Low-High HACCP Precise Pt-100 Thermometer with Timer "THE5", Max · Min · Avg.
Ideal for Foodstuff, with IP67 Waterproof, 99Points Data Record, Alarm & Temp-Programmable, $0.01^{\circ}\text{C}/^{\circ}\text{F}$, 정밀 백금 센서 저온-고온 온도계



■ Supplied by DAIHAN Korea

- * $^{\circ}\text{C}/^{\circ}\text{F}$ Switchable
- * Max./Min./Avg. Reading Display
- * Hold Function to Freeze Current Reading
- * Countdown Timer Function Included

- * Rigid Housing Design and the Most Accurate Pt100 Probe
- * Complete Set IP67 Waterproof and Washable in Washing Machine
- * Green/Red Dual Color Backlight to Indicate Temp. Hazard
- * Temp. High Alarm Threshold Programmable

Specifications & Ordering information

<L1>

Spec.	Model	"THE5"
Measure Range		$-100^{\circ}\text{C} + 300^{\circ}\text{C}$
Resolution		$0.01^{\circ}\text{C} / 0.01^{\circ}\text{F}$
Accuracy		$\pm 0.15^{\circ}\text{C}$
Temp. Response Time		3 mins
RTD Probe Size (mm)		Length - 150mm Diameter - 3.2mm Handle Length - 100mm, Cable Length - 1100mm
Battery		6 × AAA Battery
Dimension / Weight(g)		169 × 78.3 × h34.4 mm (LCD Size ; 54 × 32.5mm) / 200g
Standard Package		Meter / RTD Probe / Stainless Holder / Battery / Manual
Cat. No. & Price ₩		DH.The3005 Set

, Digital USB Dataloggers 온도 데이터 로거

DAIHAN™ Single Use USB Cold Chain Temperature Datalogger "THE3035", with PDF Report, Waterproof (IP65), $-30^{\circ}\text{C} + 70^{\circ}\text{C}$
Auto-Generate a PDF & Excel Report, Programmable by Users, 일회용 콜드체인 온도 Data 로거



■ Supplied by DAIHAN Korea

- 온도 값을 측정하고 기록하여 사용자의 필요에 따라 Excel 또는 PDF 형식으로 보고서 생성이 가능함.
- Temperature Values can be Measured and Recorded to Generate Reports in Excel or PDF Format by the User.
- USB-type Product ; Can be Used by Simply Connecting to a PC without Having to Install any Software.
- * IP65 waterproof bag could provide the moisture barrier properties of data logger.

Specifications & Ordering information

<L1>

Spec.	Model	"THE3035"
Measure Range (Resolution) & Accuracy		$-30^{\circ}\text{C} + 70^{\circ}\text{C}$, $-22^{\circ}\text{F} + 158^{\circ}\text{F}$ (0.1°C , 0.1°F) & $\pm 0.5^{\circ}\text{C}$
Sampling Point & Interval		8,192 readings & 30 seconds, 5, 10, 30, 60, 90, 120 minutes
Dimensions & Weight		80 × 43 × 2.5mm, 105 × 55 × h2.5 mm(with Bag) & 10g
Standard Package & Battery		Logger, Manual, Paper Box, 1 × 3.0V CR2032 Battery
Cat. No. & Price ₩		DH.The3035 Ea



DAIHAN™ LCD Temperature Datalogger "DL104", 62000 Sets of Data
With USB Cable, $-20^{\circ}\text{C} + 70^{\circ}\text{C}$, 데이터 로거

- * Instant current temperature read-out
- * Completed with hidden buttons to avoid unexpected operation
- * Hanging hole and magnet for easy mounting
- * Works on Windows XP / Vista / 7 or above

Specifications & Ordering information

<L1>

Spec.	Model	"DL104"
Measure Range (Resolution) & Accuracy		$-20^{\circ}\text{C} + 70^{\circ}\text{C}$ or $-4^{\circ}\text{F} + 158^{\circ}\text{F}$ (0.1°C) & $\pm 1^{\circ}\text{C}$
Mode of recording		One-time, Endless Loop (overwriting), End At
Data Capa. and Waterproof		62000 sets of data, and IP × 5
Battery		1 × 3 volts "CR2032" size or equivalent × 2
Log interval		1/2/5/10/15/30 seconds, 1/2/5/10/15/30/60/90/120 minutes
Dimension		41 × 60 × 23mm (LCD Size ; 28 × 18mm)
Cat. No. & Price ₩		A1.DL104 Ea

DAIHAN™ Temperature USB 2.0 Datalogger “THE3031”, with PDF Report, Waterproof (IP67), $-30^{\circ}\text{C} + 70^{\circ}\text{C}$
Auto-Generate a PDF & Excel Report, Programmable by Users, USB Type 온도 Data 로거



USB Datalogger

- Supplied by DAIHAN Korea
 - 온도 값을 측정하고 기록하여 사용자의 필요에 따라 Excel 또는 PDF 형식으로 보고서 생성이 가능함.
 - USB형식의 제품으로 별도의 소프트웨어 설치 필요없이 PC에 연결만 하면 사용할 수 있음.
 - Temperature Values can be Measured and Recorded to Generate Reports in Excel or PDF Format by the User.
 - USB-type Product ; Can be Used by Simply Connecting to a PC without Having to Install any Software.
- * Automatically Generate a PDF and Excel Report while Plugged into a Computer.

Specifications & Ordering information

<L1>

Spec.	Model	“THE3031”
Measure Range (Resolution) & Accuracy		$-30^{\circ}\text{C} + 70^{\circ}\text{C}$, $-22^{\circ}\text{F} + 158^{\circ}\text{F}$ (0.1°C, 0.1°F) & $\pm 0.5^{\circ}\text{C}$
Sampling Point & Interval		48,000 readings & 30 seconds, 5, 10, 30, 60, 90, 120 minutes
Dimensions & Weight		110×45×17mm & 50g
Standard Package & Battery		Logger, Manual, Paper Box, 2×3.0V CR2032 Battery
Cat. No. & Price ₩		DH.The3031 Ea

DAIHAN™ Bluetooth 4.0 Wireless Temperature Monitoring Datalogger “THE3033”, Waterproof (IP67), $-30^{\circ}\text{C} + 70^{\circ}\text{C}$
With Free Application for Smart Devices, Bluetooth Type 온도 Data 로거



Bluetooth Datalogger

- Supplied by DAIHAN Korea
- 온도 값을 측정하고 기록하여 Application을 통해 사용자가 실시간으로 데이터 확인 및 보고서 생성이 가능함.
- Users can Check the Temperature Change in Real-time through the Application “Nice Logger” Downloaded in the App Store.

* Report Contains Data Summary, Graph and Detailed Data.
* Big LCD Provides Real-time Data Every 5 seconds.

Specifications & Ordering information

<L1>

Spec.	Model	“THE3033”
Measure Range (Resolution) & Accuracy		$-30^{\circ}\text{C} + 70^{\circ}\text{C}$, $-22^{\circ}\text{F} + 158^{\circ}\text{F}$ (0.1°C, 0.1°F) & $\pm 0.5^{\circ}\text{C}$
Sampling Point & Interval		48,000 readings & 10 seconds to 120 minutes
Dimensions & Weight		70×50×h23mm & 70g
Standard Package & Battery		Logger, Manual, Paper Box, 1×3.0V CR2 Battery
Cat. No. & Price ₩		DH.The3033 Ea

DAIHAN™ Ultra Low Temperature Datalogger “THE3018”, with RTD Pt1000 Sensor, $-100^{\circ}\text{C} + 120^{\circ}\text{C}$
Auto-Generate a PDF & Excel Report, Programmable by Users, 초저온 측정 Data 로거



Ultra Low Temperature Datalogger

- Supplied by DAIHAN Korea
- 동결 건조 등의 초저온 측정을 위한 외부 RTD Pt1000 프로브 제공.
- * Plug & Play. PC Software and Driver are not Required.
- * One Key Operation to know Leftover Battery Life Time.
- * Optimize Lyophilization by Thermal Analysis
- * PDF/Excel Report Include the Summary, Graph, and the Details Data.
- * User Programmable ; 6 Languages, Sampling Interval, Start Delay, Alarm Delay, Unit, Alarm Range, Alarm Type.

Specifications & Ordering information

<L1>

Spec.	Model	“THE3018”	
		Temp (NTC Thermistor)	Temp (Pt1000 Sensor)
Measure Range		$-30^{\circ}\text{C} + 70^{\circ}\text{C}$ / $-22 + 158^{\circ}\text{F}$	$-100^{\circ}\text{C} + 120^{\circ}\text{C}$ / $-148 + 248^{\circ}\text{F}$
Accuracy		$\pm 0.5^{\circ}\text{C}$	$\pm 0.5^{\circ}\text{C}$ (at $10^{\circ}\text{C} \sim 40^{\circ}\text{C}$ Place)
Resolution		0.1°C / 0.1°F	
Sampling Point		24,000 Internal & 24,000 External	
Sampling Interval		30 seconds, 5, 10, 30, 60, 90, 120 minutes	
Alarm range & Type		$-30^{\circ}\text{C} + 70^{\circ}\text{C}$ / Single Event, Cumulative, Disable	
Dimensions & Weight		110×45×h17mm / 50g	
Battery		2×3.0V CR2032 Battery	
Standard Package		Logger, RTD Probe, Battery, Manual, Paper Box	
Cat. No. & Price ₩		DH.The3018 Ea	

DAIHAN™ Programmable Ball Mill “BML”, 2-/6- Pots, 50~600 rpm, ± 5 rpm, 10 Steps 6 Programs
With Digital Feedback Control System, Digital LCD, with Certi. & Traceability, 디지털/프로그램식 볼밀/자밀, 2구 & 6구용, 정밀형, 속도조절식

CE



3-Year Warranty

(1) 2-Places Programmable Ball Mill “BML-2”

CE



3-Year Warranty

(2) 6-Places Programmable Ball Mill “BML-6”



Minimum Noise and Safety

Optional Soundproof Cabinet “BMC-6”
for 6-Places Programmable Ball Mill “BML-6”

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Grinding / Mixing ; Crystalline Chemicals, Inks, Paints, Pigments, Plants, Racks and Soils
- Microprocessor, Programmable, Variable Speed ; 50~600 rpm
- Digital Feedback Control and LCD ; Speed (RPM), Power Bar-Graph and Timer Status
- * Antistatic Rubber Roller ; Keeps Jars / Pots Centered
- * Digital LCD with Back-Light Function
- * Over Load Protection * Ergonomic and Compact Design
- * NBR Roller with Pot Guide Ring
- * Programmable Operation ; 10 Steps 6 Programs
- * Storage Function (RPM, Timer)
- * Soundproof Cabinet for Safety and Minimum Noise (Optional)

Specifications & Ordering information

<L1>

Spec. / Model		(1) “BML-2” Max. 2 Pots	(2) “BML-6” Max. 6 Pots
Capacity		2-places, up to 2 pots	6-places, up to 6 pots
Motor		DC 40W Geared Motor (5 : 1 ratio)	DC 100W Geared Motor (5 : 1 ratio)
Maximum Load		30 kg	60 kg
Speed Range		50 ~ 600 rpm	
Speed Accuracy		± 5 rpm	
Timer		99 hr 59 min (Programmable)	
Display		Digital LCD with Back-Light Function	
Controller		Programmable Digital Feedback Controller with Jog-Shuttle Switch (Turn+Push)	
Program Function		10 Steps, 6 Programs	
Safety Device		Over Load Protector for Motor	
Dimension (w×d×h)	Roller (mm)	50×300	70×800
	External (mm)	620×300×h515	1,200×500×h795
Packing Size (w×d×h) (mm)		760×460×h700	1,300×620×h1,010
& Gross Weight		52 kg	115 kg
Power Consumption		40 W	100 W
Electric Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	
Cat. No.	230V	DH.WBM05002	DH.WBM05006
	120V	DH.WBM05102	DH.WBM05106
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Optional Soundproof Cabinet

* Designed for Minimizing Noise, and Safety Work

- DH.WBM050102 Soundproof Cabinet, “BMC2”, for 2-places“BML-2”,
1,020×570×h1,120mm
- DH.WBM050106 Soundproof Cabinet, “BMC6”, for 6-places“BML-6”,
1,600×770×h1,400mm

Ea.	Price, ₩



Porcelain Mill Pot, with Clamp, 0.5~5 Lit, 자기 포트

■ Supplied by DAIHAN Korea ■ 500 & 1000ml Pots ; No-Lid, No-Clamp., Have to Use with Rubber-Stoppers <L3>

* Ideal for Milling, 0.5 & 1 Lit. Pot is Plug Stopper

Ea.	Price, ₩
D12.0041	Pot, Porcelain, 0.5 Lit., ϕ 120×L150 mm
D12.0042	, 1.0 Lit., ϕ 150×L180 mm
D12.0043	, 2.0 Lit., ϕ 180×L210 mm (with clamp)
D12.0044	, 3.0 Lit., ϕ 210×L240 mm (with clamp)
D12.0045	, 5.0 Lit., ϕ 240×L270 mm (with clamp)



Kartell® Wide-Neck Storage Bottle, HDPE, with Heavy Wall & Handle, 5 & 10 Lit

Can be Used for Ball-Mill Pot, -50°C +105/120°C, <Italy-Made> HDPE 광구, 스토리지 바틀

■ Handle를 제거하고 Mill Pot로 사용함

* Very Resistant Bottles, Ideal for Liquid or Powders. Wide Mouth for Easy Filling and Emptying.

* Inner Cap & Blue Screwcap ; Perfect Seal.

* Non Autoclavable

<L1>

Ea.	Price, ₩	Pk.	Price, ₩
KA.1640	Bottle, HDPE, W-N, 5 Lit , ϕ 160×h 320 mm, Neck id. 90 mm	4ea	
KA.1642	, 10 Lit, ϕ 210×h 390 mm, Neck id. 120 mm	4ea	



Porcelain Ball

Porcelain Ball, ϕ 20~ ϕ 70mm, <Korea-Made> 자기 볼

■ Supplied by DAIHAN Korea

Pk.	Price, ₩
D12.1014	Ball, Porcelain, ϕ 20 mm
D12.1015	, ϕ 30 mm
D12.1016	, ϕ 40 mm
D12.1017	, ϕ 50 mm
D12.1018	, ϕ 60 mm
D12.1019	, ϕ 70 mm

<L3>

Alumina Ball, High-purity, Al₂O₃ 93%, ϕ 2~ ϕ 30mm, 알루미나 볼

■ Supplied by DAIHAN Korea

* Alumina Ball Push into the Millport, Samples were Crushed, Mixed, Decentralization <L3>

* Manufacture High-density ; Very Strong and Solid Materials Used Safely Without Fear of Damage

Pk.	Price, ₩
A8.1020	Ball, Alumina, ϕ 2 mm , Al ₂ O ₃ 93%
A8.1021	, ϕ 3 mm , Al ₂ O ₃ 93%
A8.1022	, ϕ 5 mm , Al ₂ O ₃ 93%
A8.1023	, ϕ 10 mm, Al ₂ O ₃ 93%
A8.1024	, ϕ 20 mm, Al ₂ O ₃ 93%
A8.1026	, ϕ 30 mm, Al ₂ O ₃ 93%



Alumina-and Zirconia-Balls

Zirconia Ball, ZrO₂ 94.8%, Non-Conductive, ϕ 2~ ϕ 20mm, 지르코니아 볼

■ Supplied by DAIHAN Korea

■ Excellent Resistant to Mechanical and Thermal Shock.

Pk.	Price, ₩
A8.1030	Ball, Zirconia, ϕ 2 mm , ZrO ₂ 94.8%
A8.1031	, ϕ 3 mm , ZrO ₂ 94.8%
A8.1032	, ϕ 5 mm , ZrO ₂ 94.8%
A8.1033	, ϕ 10 mm, ZrO ₂ 94.8%
A8.1034	, ϕ 15 mm, ZrO ₂ 94.8%
A8.1035	, ϕ 20 mm, ZrO ₂ 94.8%

<L3>



STS #304 Balls

Fine Stainless-steel Ball, #304, ϕ 1.6~ ϕ 25.4mm, 비자성 스텐레스 볼

* Non-Magnetic & Excellent Corrosion Resistance and Durability

* Can be Used in High/Low Temperature

Cs.	Price, ₩
SL.Bal7001	Balls, Stainless-steel, ϕ 1.6mm
SL.Bal7002	, ϕ 3.2mm
SL.Bal7003	, ϕ 6.3mm
SL.Bal7004	, ϕ 9.5mm
SL.Bal7005	, ϕ 12.7mm
SL.Bal7006	, ϕ 16mm
SL.Bal7007	, ϕ 19mm
SL.Bal7008	, ϕ 25.4mm

<L2>

DAIHAN Wisd™ Cutting Mill & Pin Mill Common Use “CML140” & “PML140”, Dry type, Max 4600rpm, Output<0.4~1.3mm
With Stainless-steel Body, Powerful 1Phase 3HP Motor, Maximum load 5kg, Input<10mm, 실험실용 커팅밀 & 핀밀 겸용



Pin Mill, “PML140”



(1) Grinding Disk Rotor



(2) Grinding Disk Stator



(3) Grinding Disk Knife



(4) Screen Sieve

3-Year Warranty

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance ■ Low Noise
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Chemicals, Pharmaceuticals, Synthetic Resins, Metals, Ceramics, Mining, etc.
- Extensive Crushing from Granular to Fine Grinding
- Excellent Quality of Output Material due to Low Heat Generated
- * High Strength and Hardness ; STS304 Body/Feeder/Rotor, STS430 Screen, STS400 Base Frame
- * The Raw Material is Diffused together with the Airflow by Rotation, and Crushed by the Impact of the Rotating Plate and the Fixed Plate.
- * Adjust the Desired Grinding Particle Size by Replacing the Internal Screen.
- * Pulverizing Various Raw Materials by Replacing the Grinding Disk
- * Easy to Separate and Assemble * Easy to Clean and Maintain * Simple Interface
- * Include ; 0.4~1.3mm 5 STS Screen Sieves, 1 STS Base Frame, 2 collection hopper

Specifications & Ordering information

<L3>

Spec.	Model	“CML140” Cutting Mill	“PML140” Pin Mill
Grinding Disc		Knife type Rotor 8ea×Cutting Knife	Pin type Rotor & Stator 10ea×Rotor Pin, 17ea×Stator Pin
Capacity		5kg(3~10kg/h)	
Hopper Size(w×d×h) / Input Raw Material		120×120×h200 / <10mm	
Screen / Output material Dia. (included)		60, 80, 100, 140, 180 / <1.3, 0.8, 0.6, 0.5, 0.4mm	
Speed & Motor	Feeder	25RPM, 40W, DC MOTER	
	Mill & Cutting	4600RPM, 1Phase 3HP	
Body Dia.		φ 165mm	
Dimension & Weight		450×550×h1150, Net. 85kg / Gross. 90kg	
Power Consumption		2.2kW	
Electrical Reqts.		3 Phase, 380V 50/60Hz	
Cat. No.	60Hz	DH.Cml7010	DH.Pml7010
	50Hz	DH.Cml7015	DH.Pml7015
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Accessory for Standard Chamber

		Ea. Price, ₩
DH.PmlR1011	(1) Grinding Disk Rotor, Stainless-steel, “R1011”, for “PML140”	
DH.PmlS1012	(2) Grinding Disk Stator, Stainless-steel, “S1012”, for “PML140”	
DH.CmlK1011	(3) Grinding Disk Knife, Stainless-steel, “K1011”, for “CML140”	
DH.PmlS1015	(4) Screen Sieve 60, Stainless-steel, Output Material Dia.<φ 1.3mm, “S1015”, for “PML140”, “CML140”	
DH.PmlS1016	(4) Screen Sieve 80, Stainless-steel, Output Material Dia.<φ 0.8mm, “S1016”, for “PML140”, “CML140”	
DH.PmlS1017	(4) Screen Sieve 100, Stainless-steel, Output Material Dia.<φ 0.6mm, “S1017”, for “PML140”, “CML140”	
DH.PmlS1018	(4) Screen Sieve 180, Stainless-steel, Output Material Dia.<φ 0.5mm, “S1018”, for “PML140”, “CML140”	
DH.PmlS1019	(4) Screen Sieve 180, Stainless-steel, Output Material Dia.<φ 0.4mm, “S1019”, for “PML140”, “CML140”	

DAIHAN Scientific

OVENS

General Purpose , Forced-air	P. 234~235
, Forced-air, SMART	P. 236~237
, Gravity-air	P. 238~239
, Gravity-air, SMART	P. 240~241
High Clean Air	P. 242
Industrial/Large-capacity	P. 243
Vacuum , Gas Exchangeable	P. 244
, Gas Exchangeable, SMART-type	P. 245
Vacuum Pumps	P. 250



Drying Oven
"SON-105"



Vacuum Drying Oven
"SOV-70"



WiRe™ App Service

Controlled by Mobile-Phone



High Clean Air Drying Oven
"WOC-560"



Industrial Drying Oven
"WOF-L1000"


DAIHAN™ Forced-Air Drying Oven “OF”, 3-Side Heating Zone, 50 · 105 · 155 · 305 Lit, Up to 250℃, ±0.3℃

With Wire Shelves, Digital PID Control, Jog-Dial & Push Button, Digital LCD, with Back-light, Pre-Heating Zone, Certi. & Traceability

강제 순환식 정밀 건조기/오븐, 고정밀 디지털 퍼지 제어, 우수한 온도 정확성/균일성



(1) Standard Model, Forced-Air Oven



3-Side Heating Zone Forced-Air Mechanism ;
The Best Temp Uniformity & Accuracy by High
Performance Heating Mechanism
(Exclude OF-305, OF-W305)



Perfect Door Lock by Double Latch

- Supplied by DAIHAN Korea
- GD(Good Design)-mark
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- Digital PID Control System Implementing Superior Temperature Accuracy
- Optimized Air Flow by Forced-air Mechanism
- Patented Jog-Dial Control System
- * Suitable for Drying, Baking, Conditioning, Curing, Pre-Heating and Aging
- * Compact New Body Design
- * RS232C Interface for Monitoring and Controlling with PC
- * Ambient 5℃~250℃ Range with Fluctuation of ±0.3℃ at 100℃
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be disabled)
- UL, CUL, and CE Certified(Exclude 305Lit-model)
- PL(Product Liability) Insurance
- * Alarm Function ; Error status and Timer-end

Specification

Spec. Model		(1) Standard Model			
		“OF-50” Max. Capa. 50L	“OF-105” Max. Capa. 105L	“OF-155” Max. Capa. 155L	“OF-305” Max. Capa. 305L
		(2) Window Model			
		-	“OF-W105”	“OF-W155”	“OF-W305”
Capacity		50 Lit	105 Lit	155 Lit	305 Lit
Dimension (w×d×h)	Interior (mm)	370×345×420	485×415×535	550×480×600	550×470×1,180
	Exterior (mm)	518×647×724	653×712×882	718×777×947	721×777×1,619
Heater Power		650 W	1.4 kW	1.6 kW	3.4 kW
Temp	Range	Room Temp +5℃~250℃			
	Fluctuation	±0.3℃ at 100℃, ±0.5℃ at 150℃			
	Variation	±1℃ at 100℃, ±1.8℃ at 150℃		±1℃ at 100℃, ±2℃ at 150℃	
	Sensor	PT 100			
Control Resolution		±1℃			
Heat-up Time		25 min. to 100℃, 45 min. to 150℃		25 min. to 100℃, 40 min. to 150℃	
Recovery Time (Door open 30sec)		8 min. to 100℃, 8 min. to 150℃	6 min. to 100℃, 6 min. to 150℃		
Controller		Digital PID Control by Advanced Microprocessor Jog-Dial with Push Button			
RS232 Port		Available to connect with PC			
Display		Digital LCD with Back Light			
Timer		99hr 59 min (delay/continuous function)			
Storage		Temperature and Timer setting values			
Material	Internal	Stainless-steel (#304), Corrosion Resistant			
	External	Powder Coated steel			
	Insulation	Glass Wool			
Shelf (included)		2 Wire Shelf			3 Wire Shelf
		Load per Stainless-steel Wire Shelf ; 16kg			
Ventilation		Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm			
		1ea	2ea		3ea
Circulation		Forced-air Type			
Safety device		Over Temp and Over Current Protector, Sensor Error Detector			
Alarm		Error Status and Timer-end			
Net Weight (kg)		44	69	78	159
Packing Size (w×d×h)(mm)		622×754×807	743×819×1,068	808×884×1,133	1,000×1,000×1,870
Shipping Weight (kg)		59	79	92	219
Power Consumption		674 W	1390 W	1566 W	3200 W
Energy Consumption, at 100℃ , at 150℃		117 Wh	187 Wh	204 Wh	-
		230 Wh	337 Wh	398 Wh	-
Electric Reqts.		120V, 60Hz or 230V, 50/60Hz			


* Other Specifications are available upon Customer's Request.

<Continued On Next Page>

<Continued From...DAIHAN™ Forced-Air Drying Oven "OF">



(2) Window Model, Forced-Air Oven

Stainless-steel Wire Shelf
(Included)

WiseRemote Control Software ;
Connection of Temperature Devices with
Interface RS232C to a PC

* You can Download this Software at
DAIHAN-Sci.com Shopping-mall

Ordering Information

(1) Standard Model, Forced-Air, with Wire Shelf

DH.WOF05050	Oven, Standard, 230V, 50 Lit, "OF-50"
DH.WOF05105	, 105 Lit, "OF-105"
DH.WOF05155	, 155 Lit, "OF-155"
DH.WOF05305	, 305 Lit, "OF-305"
DH.WOF06050	Oven, Standard, 120V, 50 Lit, "OF-50"
DH.WOF06105	, 105 Lit, "OF-105"
DH.WOF06155	, 155 Lit, "OF-155"
DH.WOF06305	, 305 Lit, "OF-305"

(2) Built in Viewing Window Model, Forced-Air, with Wire Shelf

DH.WOF07105	Oven, Window, 230V, 105 Lit, "OF-W105"
DH.WOF07155	, 155 Lit, "OF-W155"
DH.WOF07305	, 305 Lit, "OF-W305"
DH.WOF08105	Oven, Window, 120V, 105 Lit, "OF-W105"
DH.WOF08155	, 155 Lit, "OF-W155"
DH.WOF08305	, 305 Lit, "OF-W305"

Spare Stainless-steel Wire Shelf

DH.WOF11050	Wire Shelf, "OFS050", for 50 Lit, w340×d310mm
DH.WOF11105	, "OFS105", for 105 Lit, w455×d375mm
DH.WOF11155	, "OFS155", for 155-/305-Lit, w520×d440mm

Forced-air Oven with Validation (IQ, OQ) Service

DH.WOF31050	Validation Service(IQ, OQ), "OFV1050", with "OF-50"
DH.WOF31105	, "OFV1105", with "OF-105"
DH.WOF31155	, "OFV1155", with "OF-155"
DH.WOF32105	, "OFV2105", with "OF-W105"
DH.WOF32155	, "OFV2155", with "OF-W155"

※ 운임 별도 <L1>

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

➤ New Jog-Dial Controller

⑤ Digital LCD
with Back-light Function

③ Main Button

① Jog-Dial

② Push-Button

④ Sub Button


- Provides the user interface with an ergonomic design for users to operate in an easy and convenient way.
- A new 2-way Jog-Dial Knob (Jog Dial + Push Button) and presents excellent durability and mobility.
- High-quality LCD with Back-light function.
- A durable Main Button and a Sub Button that features superb and user-friendly operation.

DAIHAN™ SMART Forced-Air Drying Oven “SOF”, 3-Side Heating Zone, 50 · 105 · 155 · 305 Lit

With Smart-Lab™ Controller, 4" Full Touch Screen TFT LCD, Fuzzy-PID Control, WiRe™ Service, Certi. & Traceability, Up to 250℃, ±0.3℃
스마트 강제 순환식 정밀 건조기/오븐, 우수한 온도 정확성/균일성, WiRe™서비스, 4인치 풀터치스크린 조절기



(1) SMART Model, Forced-Air Oven



3-Side Heating Zone Forced-Air Mechanism ;
The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
(Exclude SOF-305, SOF-W305)



Perfect Door Lock by Double Latch

- Supplied by DAIHAN Korea
- UL, CUL, and CE Certified(Exclude 305 Lit-model)
- GD(Good Design)-mark
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Optimized Air Flow by Forced-air Mechanism
- Smart-Lab™ Controller WiRe™ App Service
 - Automatic Data Recording
 - Self-Diagnostic Function
 - Variable Program Setting for Multiple Patterns with Multiple Segments
 - Internet Connectivity with WiFi
 - Controller SW can be upgraded through USB port.
 - The Recorded Data could be transferred to PC by USB memory
 - RS232C Interface for Monitoring and Controlling with PC

* Ambient +5℃ to 250℃ of Temp Range, Fluctuation of ±0.3℃ at 100℃
* Suitable for Drying, Baking, Conditioning, Curing, Pre-Heating and Aging

Specification

Spec.		Model	(1) SMART Model			
			“SOF-50” Max. Capa. 50L	“SOF-105” Max. Capa. 105L	“SOF-155” Max. Capa. 155L	“SOF-305” Max. Capa. 305L
			(2) SMART-Window Model			
		-	“SOF-W105”	“SOF-W155”	“SOF-W305”	
Capacity		50 Lit	105 Lit	155 Lit	305 Lit	
Dimension (w×d×h)	Interior (mm)	370×345×420	485×415×535	550×480×600	550×470×1,180	
	Exterior (mm)	518×647×724	653×712×882	718×777×947	721×777×1,619	
Heater Power		650 W	1.4 kW	1.6 kW	3.4 kW	
Temp	Range	Room Temp +5℃~250℃				
	Fluctuation	±0.3℃ at 100℃, ±0.5℃ at 150℃				
	Variation	±1℃ at 100℃, ±1.8℃ at 150℃, ±1℃ at 100℃, ±2℃ at 150℃				
	Sensor	PT 100				
Control Resolution		±1℃				
Heat-up Time		25 min. to 100℃, 45 min. to 150℃		25 min. to 100℃, 40 min. to 150℃		
Recovery Time (Door open 30sec)		8 min. to 100℃, 8 min. to 150℃	6 min. to 100℃, 6 min. to 150℃			
Controller		Smart-Lab™ Controller				
RS232 Port		Available to connect with PC				
Display		4" Full Touch Screen TFT LCD				
Timer		99hr 59min (delay/continuous function)				
Storage		Temperature and Timer setting values				
Material	Internal	Stainless-steel (#304), Corrosion Resistant				
	External	Powder Coated steel				
	Insulation	Glass Wool				
Shelf (included)		2 Wire Shelf			3 Wire Shelf	
		Load per Stainless-steel Wire Shelf ; 16kg				
Ventilation		Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm				
		1ea	2ea		3ea	
Circulation		Forced-air Type				
Safety device		Over Temp and Over Current Protector, Sensor Error Detector				
Alarm		Error Status and Timer-end				
Net Weight (kg)		44	69	78	159	
Packing Size (w×d×h)(mm)		622×754×807	743×819×1,068	808×884×1,133	1,000×1,000×1,870	
Shipping Weight (kg)		59	79	92	219	
Power Consumption		674 W	1390 W	1566 W	3200 W	
Energy Consumption, at 100℃ , at 150℃		117 Wh	187 Wh	204 Wh	-	
		230 Wh	337 Wh	398 Wh	-	
Electric Reats.		120V, 60Hz or 230V, 50/60Hz				

* Other Specifications are available upon Customer's Request.


<Continued On Next Page>

<Continued From...DAIHAN™ SMART Forced-Air Drying Oven "SOF">



(2) SMART-Window Model, Forced-Air Oven

3-Year Warranty

Stainless-steel Wire Shelf
(Included)

Ordering Information

(1) SMART Model, Forced-Air, with Wire Shelf

DH.SWOF05050	Oven, Standard, 230V, 50 Lit, "SOF-50"
DH.SWOF05105	, 105 Lit, "SOF-105"
DH.SWOF05155	, 155 Lit, "SOF-155"
DH.SWOF05305	, 305 Lit, "SOF-305"
DH.SWOF06050	Oven, Standard, 120V, 50 Lit, "SOF-50"
DH.SWOF06105	, 105 Lit, "SOF-105"
DH.SWOF06155	, 155 Lit, "SOF-155"
DH.SWOF06305	, 305 Lit, "SOF-305"

(2) SMART Built-in Viewing Window Model, Forced-Air, with Wire Shelf

DH.SWOF07105	Oven, Window, 230V, 105 Lit, "SOF-W105"
DH.SWOF07155	, 155 Lit, "SOF-W155"
DH.SWOF07305	, 305 Lit, "SOF-W305"
DH.SWOF08105	Oven, Window, 120V, 105 Lit, "SOF-W105"
DH.SWOF08155	, 155 Lit, "SOF-W155"
DH.SWOF08305	, 305 Lit, "SOF-W305"

Spare Stainless-steel Wire Shelf

DH.WOF11050	Wire Shelf, for 50 Lit, "OFS050", w340×d310mm
DH.WOF11105	, for 105 Lit, "OFS105", w455×d375mm
DH.WOF11155	, for 155 / 305 Lit, "OFS155", w520×d440mm

Forced-air Oven with Validation (IQ, OQ) Service

DH.SWOF31050	Validation Service(IQ, OQ), "SOFV1050", with "SOF-50"
DH.SWOF31105	, "SOFV1105", with "SOF-105"
DH.SWOF31155	, "SOFV1155", with "SOF-155"
DH.SWOF32105	, "SOFV2105", with "SOF-W105"
DH.SWOF32155	, "SOFV2155", with "SOF-W155"

※ 운임 별도 <L1>

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

>> WiRe™ App Service

4" Full Touch Screen TFT LCD
(Smart-Lab™ Controller)

It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.




WiRe™ App Service

Controlled by Mobile-Phone

DAIHAN™ Gravity-Air Drying Oven "ON", 3-Side Heating Zone, 32 · 50 · 105 · 155 Lit, Up to 230℃, ±0.5℃
 With 2 Wire Shelf, Digital PID Control, Jog-Dial & Push Button, Digital LCD with Back-light, with Certi. & Traceability
 자연 대류식 정밀 건조기/오븐, 고정밀 디지털 퍼지 제어, 우수한 온도 정확성



(1) Standard Model, Gravity-Air Oven



3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism (Gravity-Air Mechanism)



Perfect Door Lock by Double Latch

- Supplied by DAIHAN Korea ■ UL, CUL, and CE Certified ■ GD(Good Design)-mark
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- Digital PID Control System Implementing Superior Temperature Accuracy
- Optimized Air Flow by Gravity-air Mechanism
- Patented Jog-Dial Control System ■ PL(Product Liability) Insurance
- * Suitable for Drying, Baking, Conditioning, Curing, Pre-Heating and Aging
- * New Jog-Dial with Push Button
- * RS232C Interface for Monitoring and Controlling with PC
- * Ambient +5℃ to 230℃ Range with Fluctuation of ±0.5℃ at 100℃
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be disabled)

Specification

Spec.		Model	(1) Standard Model			
			"ON-32"	"ON-50"	"ON-105"	"ON-155"
			Max. Capa. 32L	Max. Capa. 50L	Max. Capa. 105L	Max. Capa. 155L
			(2) Window Model			
		-	-	"ON-W105"	"ON-W155"	
Capacity		32 Lit	50 Lit	105 Lit	155 Lit	
Dimension (w×d×h)	Interior (mm)	310×290×360	370×350×420	485×409×535	550×474×600	
	Exterior (mm)	458×557×664	518×647×724	653×712×882	718×777×947	
Heater Power		500 W	650 W	1.4 kW	1.6 kW	
Temp	Range	Room Temp +5℃ to 230℃				
	Fluctuation	± 0.5℃ at 100℃, ± 0.6℃ at 150℃			± 0.5℃ at 100℃, ± 0.8℃ at 150℃	
	Variation	± 1.2℃ at 100℃, ± 2.2℃ at 150℃				
	Sensor	PT 100				
Control Resolution		± 1℃				
Heat-up Time		30 min. to 100℃, 43 min. to 150℃	30 min. to 100℃, 45 min. to 150℃	35 min. to 100℃, 50 min. to 150℃		
Recovery Time (Door open 30sec)		8 min. to 100℃, 10 min. to 150℃	10 min. to 100℃, 12 min. to 150℃			
Controller		Digital PID Control by Advanced Microprocessor Jog-Dial with Push Button				
RS232 Port		Available to connect with PC				
Display		Digital LCD with Back Light				
Timer		99hr 59 min (delay/continuous function)				
Storage		Temperature and Timer setting values				
Material	Internal	Stainless-steel (#304), Corrosion Resistant				
	External	Powder Coated steel				
	Insulation	Glass Wool				
Shelf (included)		2 Stainless-steel Wire Shelf, Load per Shelf ; 16 kg				
Ventilation		Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm				
		1ea		2ea		
Circulation		Gravity-air Type				
Safety device		Over Temp and Over Current Protector, Sensor Error Detector				
Alarm		Error Status and Timer-end				
Net Weight (kg)		36	44	69	78	
Packing Size (w×d×h) (mm)		562×664×747	622×754×807	743×819×1,068	808×884×1,133	
Shipping Weight (kg)		43	59	79	92	
Power Consumption		502 W	674 W	1390 W	1566 W	
Energy Consumption, at 100℃		90 Wh	117 Wh	187 Wh	204 Wh	
, at 150℃		181 Wh	230 Wh	337 Wh	398 Wh	
Electric Reqts.		120V, 60Hz or 230V, 50/60Hz				


* Other Specifications are available upon Customer's Request.

<Continued On Next Page>

<Continued From...DAIHAN™ Gravity-Air Drying Oven "ON">



(2) Window Model, Gravity-Air Oven

Stainless-steel Wire Shelf
(Included)**WiseRemote Control Software ;**Connection of Temperature Devices
with Interface RS232C to a PC* You can Download this Software at
DAIHAN-Sci.com Shopping-mall**Ordering Information****(1) Standard Model, Gravity-Air, with Wire Shelf**

DH.WON05032	Oven, Standard, 230V, 32 Lit, "ON-32"
DH.WON05050	, 50 Lit, "ON-50"
DH.WON05105	, 105 Lit, "ON-105"
DH.WON05155	, 155 Lit, "ON-155"
DH.WON06032	Oven, Standard, 120V, 32 Lit, "ON-32"
DH.WON06050	, 50 Lit, "ON-50"
DH.WON06105	, 105 Lit, "ON-105"
DH.WON06155	, 155 Lit, "ON-155"

(2) Built in Viewing Window Model, Gravity-Air, with Wire Shelf

DH.WON07105	Oven, Window, 230V, 105 Lit, "ON-W105"
DH.WON07155	, 155 Lit, "ON-W155"
DH.WON08105	Oven, Window, 120V, 105 Lit, "ON-W105"
DH.WON08155	, 155 Lit, "ON-W155"

Spare Stainless-steel Wire Shelf

DH.WON11032	Wire Shelf, "ONS032", for 32 Lit., w280 × d280mm
DH.WON11050	, "ONS050", for 50 Lit., w340 × d310mm
DH.WON11105	, "ONS105", for 105 Lit., w455 × d375mm
DH.WON11155	, "ONS155", for 155 Lit., w520 × d440mm

Gravity-air Oven with Validation (IQ, OQ) Service

DH.WON31032	Validation Service(IQ, OQ), "ONV1032", with "ON-32"
DH.WON31050	, "ONV1050", with "ON-50"
DH.WON31105	, "ONV1105", with "ON-105"
DH.WON31155	, "ONV1155", with "ON-155"
DH.WON32105	, "ONV2105", with "ON-W105"
DH.WON32155	, "ONV2155", with "ON-W155"

※ 윤임 별도 <L1>

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

➤ New Jog-Dial Controller⑤ Digital LCD
with Back-light Function

③ Main Button

① Jog-Dial

② Push-Button

④ Sub Button

- Provides the user interface with an ergonomic design for users to operate in an easy and convenient way.
- A new 2-way Jog-Dial Knob (Jog Dial + Push Button) and presents excellent durability and mobility.
- High-quality LCD with Back-light function.
- A durable Main Button and a Sub Button that features superb and user-friendly operation.


DAIHAN™ SMART Gravity-Air Drying Oven "SON", 3-Side Heating Zone, 32 · 50 · 105 · 155 Lit
 With Smart-Lab™ Controller, 4" Full Touch Screen, Fuzzy-PID Control, WiRe™ Service, with Certi. & Traceability, Up to 230℃, ±0.5℃
 스마트 자연 대류식 정밀 건조기/오븐, 우수한 온도 정확성/균일성, WiRe™ 서비스, 4인치 풀터치스크린 조절기



3-Year Warranty

WiRe™ App Service

(1) SMART Model, Gravity-Air Oven



3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism (Gravity-Air Mechanism)



Perfect Door Lock by Double Latch

- Supplied by DAIHAN Korea ■ UL, CUL, and CE Certified ■ GD(Good Design)-mark
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- 3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism

- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Optimized Air Flow by Gravity-air Mechanism ■ PL (Product Liability) Insurance
- Smart-Lab™ Controller WiRe™ App Service

- Automatic Data Recording
- Ergonomic Design 4" Full Touch Screen TFT LCD
- Self-Diagnostic Function
- Variable Program Setting for Multiple Patterns with Multiple Segments
- The Recorded Data could be transferred to PC by USB memory
- Controller SW can be upgraded through USB port.
- RS232C Interface for Monitoring and Controlling with PC
- Self-Diagnostic Function
- Internet Connectivity with WiFi
- Automatic Data Recording

- * Ambient +5℃ to 230℃ of Temp Range, Fluctuation of ±0.5℃ at 100℃
- * Suitable for Drying, Baking, Conditioning, Curing, Pre-Heating and Aging
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be disabled)

Specification

Spec. \ Model		(1) SMART Model			
		“SON-32” Max. Capa. 32L	“SON-50” Max. Capa. 50L	“SON-105” Max. Capa. 105L	“SON-155” Max. Capa. 155L
		(2) SMART-Window Model			
		-	-	“SON-W105”	“SON-W155”
Capacity		32 Lit	50 Lit	105 Lit	155 Lit
Dimension (w×d×h)	Interior (mm)	310×290×360	370×350×420	485×409×535	550×474×600
	Exterior (mm)	458×557×664	518×647×724	653×712×882	718×777×947
Heater Power		500 W	650 W	1.4 kW	1.6 kW
Temp	Range	Room Temp +5℃ to 230℃			
	Fluctuation	± 0.5℃ at 100℃, ± 0.6℃ at 150℃			± 0.5℃ at 100℃, ± 0.8℃ at 150℃
	Variation	± 1.2℃ at 100℃, ± 2.2℃ at 150℃			
	Sensor	PT 100			
Control Resolution		±1℃			
Heat-up Time		30 min. to 100℃, 43 min. to 150℃	30 min. to 100℃, 45 min. to 150℃	35 min. to 100℃, 50 min. to 150℃	
Recovery Time (Door open 30sec)		8 min. to 100℃, 10 min. to 150℃		10 min. to 100℃, 12 min. to 150℃	
Controller		Smart-Lab™ Controller			
RS232 Port		Available to connect with PC			
Display		4" Full Touch Screen TFT LCD			
Timer		99hr 59 min (delay/continuous function)			
Storage		Temperature and Timer setting values			
Material	Internal	Stainless-steel (#304), Corrosion Resistant			
	External	Powder Coated steel			
	Insulation	Glass Wool			
Shelf (included)		2 Stainless-steel Wire Shelf, Load per Shelf ; 16 kg			
Ventilation		Ventilation hole with Stainless-steel Cap, Hole Diameter 40mm			
		1ea		2ea	
Circulation		Gravity-air Type			
Safety device		Over Temp and Over Current Protector, Sensor Error Detector			
Alarm		Error Status and Timer-end			
Net Weight (kg)		36	44	69	78
Packing Size (w×d×h) (mm)		562×664×747	622×754×807	743×819×1,068	808×884×1,133
Shipping Weight (kg)		43	59	79	92
Power Consumption		502 W	674 W	1390 W	1566 W
Energy Consumption, at 100℃		90 Wh	117 Wh	187 Wh	204 Wh
, at 150℃		181 Wh	230 Wh	337 Wh	398 Wh
Electric Reqts.		120V, 60Hz or 230V, 50/60Hz			


* Other Specifications are available upon Customer's Request.

<Continued On Next Page>

<Continued From...DAIHAN™ SMART Gravity-Air Drying Oven "SON">



(2) SMART-Window Model, Gravity-Air Oven

Stainless-steel Wire Shelf
(Included)

Ordering Information

(1) SMART Model, Gravity-Air, with Wire Shelf

DH.SWON05032	Oven, Standard, 230V, 32 Lit, "SON-32"
DH.SWON05050	, 50 Lit, "SON-50"
DH.SWON05105	, 105 Lit, "SON-105"
DH.SWON05155	, 155 Lit, "SON-155"
DH.SWON06032	Oven, Standard, 120V, 32 Lit, "SON-32"
DH.SWON06050	, 50 Lit, "SON-50"
DH.SWON06105	, 105 Lit, "SON-105"
DH.SWON06155	, 155 Lit, "SON-155"

(2) SMART Built in Viewing Window Model, Gravity-Air, with Wire Shelf

DH.SWON07105	Oven, Window, 230V, 105 Lit, "SON-W105"
DH.SWON07155	, 155 Lit, "SON-W155"
DH.SWON08105	Oven, Window, 120V, 105 Lit, "SON-W105"
DH.SWON08155	, 155 Lit, "SON-W155"

Spare Stainless-steel Wire Shelf

DH.WON11032	Wire Shelf, "ONS032", for 32 Lit, w280×d280mm
DH.WON11050	, "ONS050", for 50 Lit, w340×d310mm
DH.WON11105	, "ONS105", for 105 Lit, w455×d375mm
DH.WON11155	, "ONS155", for 155 Lit, w520×d440mm

Gravity-air Oven with Validation (IQ, OQ) Service

DH.SWON31032	Validation Service(IQ, OQ), "SONV1032", with "SON-32"
DH.SWON31050	, "SONV1050", with "SON-50"
DH.SWON31105	, "SONV1105", with "SON-105"
DH.SWON31155	, "SONV1155", with "SON-155"
DH.SWON32105	, "SONV2105", with "SON-W105"
DH.SWON32155	, "SONV2155", with "SON-W155"

※ 운임 별도 <L1>

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

>> WiRe™ App Service

4" Full Touch Screen TFT LCD
(Smart-Lab™ Controller)

It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.



WiRe™ App Service

Controlled by Mobile-Phone

DAIHAN™ High Clean Air Drying Oven “WOC”, with HEPA-H13 Filter, 560 · 800 Lit, with Certi. & Traceability
With Chrome Coated Steel Wire Shelves, Digital Fuzzy Control System, Clean Air Convection Mechanism, Up to 250°C, ±1.0°C
청정 고온 건조기/오븐, 고정밀 디지털 퍼지 제어, 우수한 청정도



High Clean Air Oven, with HEPA Filter
“WOC-560”

3-Year Warranty

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Specially Designed Clean Air Convection Mechanism Provides Dynamic Air Flow for Exhaust of Contaminated Air from the Chamber and Supply of Filtered Clean Air into the Chamber
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy.
- Patented Jog-Shuttle Control System

- * Suitable for Critical Clean Application such as Electronics, Semiconductor, Pharmaceutical and Clinical Applications . HEPA-H13 Filter is used.
- * User's Self-Compensation Function to Control Difference between Real Temp and Indicated Temp (±20.0°C)
- * Storage Function ; Set Values of Temperature and Timer
- * Locking Mode Supported for Experimental Safety(Input to Jog-Shuttle can be Disabled)
- * Over Temperature & Over Current Protection, Sensor Error Detection

Specifications & Ordering information

<L1>

Spec.	Model	“WOC-560” Max. Capa. 560L	“WOC-800” Max. Capa. 800L
Capacity		560 Lit	800 Lit
Air Flow-type		High Clean Air Convection Mechanism	
Filter for Clean Air		HEPA-H13 Filter for High Temperature	
Temp	Range & Sensor	60°C~250°C, PT100 Sensor	
	Accuracy & Uniformity	±1.0°C, ±2.0°C at 200°C	
Timer & Alarm		99hr 59min (with delay / Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back Light	
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch(Turn+Push)	
Verification		Timer Operated Verification System	
Resolution		Control ; ±1.0°C, Display ; ±1.0°C	
Shelf		3ea of Perforated Shelf of Adjustable Position	
Load per Rack		16 kg, Total ; 48 kg (16kg × 3 Shelf)	
Safety Circuit		Over Temp and Over Current Protector, Sensor Error Detector	
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel, Insulation ; Glass Wool	
Others		Storage Function, Locking Mode	
Dimension (w×d×h)	Interior (mm)	800×700×1,000	1,000×800×1,000
	Exterior (mm)	1,000×1,155×1,735	1,500×1,200×1,900
Packing Size (w×d×h) (mm)		1,200×1,400×h2,120mm, 386 kg	1,700×1,400×h2,270mm, 446 kg
Power Consumption		6.5 kW	9.4 kW
Electric Reqts.		230V or 380V 3 Phase, 50/60Hz	
Cat. No.	230V	DH.WOC00560	DH.WOC00800
	380V 3 Phase	DH.WOC03560	DH.WOC03800
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도


Spare Perforated Shelf, Stainless-steel

Ea. Price, ₩

DH.WOC10560 (1) Perforated Shelf, “OCS560”, w780×d680, for “WOC-560”
DH.WOC10800 , “OCS800”, w980×d880, for “WOC-800”



High Efficiency and
High Clean Air Convection
Mechanism



High Temp HEPA Filter



(1) Perforated Shelf, Stainless-steel
(Included)

DAIHAN™ Industrial Forced-Air Drying Oven "WOF-L", 486 · 840 · 1176 Lit, with Certi. & Traceability
With Stainless-steel Perforated Shelves, Digital PID Control System, Superior Temp Accuracy, Up to 250°C, ±1.0°C
대용량 산업용 열풍 순환식 건조기 / 오븐, 우수한 온도 정확성



Industrial Forced-Air Oven "WOF-L1000"
with Stainless-steel Perforated Shelf

3-Year Warranty



(1) Perforated Shelf, Stainless-steel
(Included)

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Heavy Forced Air-Circulation by Strong and Reliable Fan System ; ideal for Semiconductor and Electronic Parts Baking Purpose

- * Voltage and Ampere Meter Equipped for Monitoring Safe Operation of the Machine
- * Digital PID Control System Enables Superior Temperature Accuracy
- * Digital Timer Function with Buzzer Alarm ; up to 99 hrs 59 min
- * Locker is Available for Operational Safety Control by an Administrator
- * Over Temperature & Over Current Protection, Sensor Error Detection and Leakage Breaking
- * Multi-Channel Temperature Sensors with Digital Display for Monitoring Temperature Uniformity in the Chamber is available(Optional)
- * Default Damper to Enhance Thermal Efficiency * Automatic Gas Exhaust Available (Optional)

Specifications & Ordering information

<L1>

Spec.	Model	“WOF-L400” Max. Capa. 486L	“WOF-L800” Max. Capa. 840L	“WOF-L1000” Max. Capa. 1,176L
Capacity		486 Lit	840 Lit	1,176 Lit
Air Flow-type		Forced-air Type		
Air Circulation Fan Motor		1/3HP	1/3 HP × 2 ea	1/3HP × 3 ea
Temp	Range	Room Temp +5℃ to 250℃		
	Accuracy	±1.0℃ at 150℃		
	Uniformity	±2.0℃ at 150℃		
	Sensor	PT100		
Heating Power/Consumption		3 kW	5 kW	7 kW
Timer & Alarm		99hr 59min (with Continuous run.), Timer-end Alarm		
Display		2 lines LED for Temperature and Timer		
Controller		Digital PID Control		
Resolution		Control ; ±1.0℃, Display ; ±1.0℃		
Door		Powder coated Steel Door with Silicone Packing		
		1 Door	2 Doors	
Shelf		Included Stainless-steel Perforated Shelf		
		3 ea	4 ea	6 ea
Load per Rack		16 kg		
Permitted Total Load		48 kg (16kg×3 Shelf)		96 kg (16kg×6 Shelf)
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel, Insulation ; Ceramic and Glass Wool Board		
Dimension (w×d×h)	Interior (mm)	900×600×900	1,250×700×960	1,400×700×1,200
	Exterior (mm)	1,200×750×1,700	1,500×840×1,800	1,650×840×2,150
Net Weight		316 kg	408 kg	478 kg
Others		Over Temp and Over Current Protector, Sensor Error Detector, Leakage Breaker Locking Mode by Mechanical Locker		
Packing Size (w×d×h) (mm) & Gross Weight		1,400×1,220×1,960	1,700×1,310×2,060	1,850×1,310×2,410
		376 kg	488 kg	568 kg
Power Consumption		3.6 kW	6.2 kW	8.8 kW
Electric Reqts.	230V	1Phase AC 230V, 50/60Hz		-
	3 Phase	-	3Phase AC 380V, 50/60Hz	
Cat. No.	230V	DH.WOFL0400	DH.WOFL0800	-
	3 Phase	-	DH.WOFL0802	DH.WOFL1002
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Perforated Shelf, Stainless-steel

	Ea.	Price, ₩
DH.WOF10400 (1) Perforated Shelf, "OFLS400", w880 × d580mm, for "WOF-L400"		
DH.WOF10800 "OFLS800", w1230 × d680mm, for "WOF-L800"		
DH.WOF11000 "OFLS1000", w660 × d680mm, for "WOF-L1000"		

Industrial Forced Convection Oven with Validation (IQ, OQ) Service

	Ea.	Price, ₩
DH.WOF30400 Validation Service(IQ, OQ), with "WOF-L400", "OFLV400"		
DH.WOF30800 , with "WOF-L800", "OFLV800"		
DH.WOF31000 , with "WOF-L1000", "OFLV1000", 3Phase		

DAIHAN™ Gas Exchangeable Precise Vacuum Drying Oven “OV”, 20 · 30 · 70 Lit, 10~750 mmHg, 200℃
With Highly Safety View Window, 2 Al-Shelf, Digital PID Control, Superior Temp & Vacuum Accuracy, Excellent Thermal Conductivity
가스 교환 정밀 진공 건조기/오븐, 고정밀 디지털 PID 제어, 알루미늄 선반 포함



WiseRemote Control Software ;
Connection of Temperature Devices with
Interface RS232C to a PC

* You can Download this Software at
DAIHAN-Sci.com Shopping-mall



(1) Standard Shelf, Al.(Included)

- Supplied by DAIHAN Korea
- UL, CUL, and CE Certified
- PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Suitable for Drying, Baking, Conditioning, Curing, Out-Gassing Solids and Liquids, Vacuum Embedding, Moisture Testing, Plating and Aging Test
- Patented Jog-Dial Control System and Back-Light LCD
- Gas Exchangeable, High Heat Conduction and New Jog-Dial Controller & Precise Gauge
- WiseRemote Control Software ; Download at DAIHAN-Sci.com.
- * Very Easy and Convenient for Vacuuming and Venting, Easy to read Vacuum Gauge
- * Highly Safety Viewing Window (Tempered Safety Glass + Polycarbonate Window)
- * Heavy-duty Safety Front Door for User Safety
- * High Quality Door Sealing ; High Temperature Grade Silicone Molded Gasket
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be disabled)
- * Compact New Body Design

Specifications & Ordering information

<L1>

Spec.	Model	“OV-20” Max. Capa. 20 Lit	“OV-30” Max. Capa. 30 Lit	“OV-70” Max. Capa. 70 Lit
Capacity		20 Lit	30 Lit	70 Lit
Vacuum Range & Level & Rate		10 ~ 750 mmHg, 0.63mbar, 0.39mbar/hr		
Temp	Range / Sensor	R/T +10℃ to 200℃, PT100		
	Accuracy	± 0.5℃ at 100℃ & 150℃		
	Uniformity	± 2.0℃ at 100℃, ±3.0℃ at 150℃		
Heating Power/Consumption		170W×2ea 130W×2ea	350W×2ea 450W×2ea	
Heat-up Time		60 min. to 100℃, 90 min. to 150℃	70 min. to 100℃, 100 min. to 150℃	80 min. to 100℃, 100 min. to 150℃
Timer & Alarm		99hr 59min (with delay / Continuous run.), Error status & Timer-end		
Display		Digital LCD with Back Light		
Controller		Digital PID Controller with Jog-Dial Switch (Turn & Push)		
Resolution		Control ; ±1.0℃, Display ; ±1.0℃		
Shelf (Included)		2 Aluminum Shelf		
Safety Circuit		Over Temp & Over Current Protector, Sensor Error Detector, Leakage Breaker		
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel, Window ; Tempered Safety Glass 12T, Poly Carbonate 5T for Safety		
Others		Storage Function (Temp and Timer), Locking Mode (Jog-Dial Input Disabled)		
Dimension (w×d×h)	Interior (mm)	265×290×265	300×330×300	400×435×400
	Exterior (mm)	480×475×547	515×515×582	618×620×682
Vacuum & Gas Inlet Valve Size		od. 10 mm		
Net Weight		48 kg	55 kg	85 kg
Packing Size (w×d×h) (mm) & Gross Weight		600×580×850	640×640×900	740×720×1000
		57 kg	80 kg	112 kg
Power Consumption		600 W	600 W	1.6 kW
Electric Rqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz		
Cat. No.	230V	DH.WOV03020	DH.WOV03030	DH.WOV03070
	120V	DH.WOV04020	DH.WOV04030	DH.WOV04070
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Aluminum Shelf

DH.WOV12020 (1) Standard Shelf, Al., “OVSNO20”, 260×245×h3mm, for 20 Lit
DH.WOV12030 “OVSNO30”, 290×285×h3mm, for 30 Lit
DH.WOV12070 “OVSNO70”, 390×390×h3mm, for 70 Lit

Ea.	Price, ₩

※ Note! Optional Vacuum Pump ; See the “Pumps, Vacuum” Section!

DAIHAN™ Gas Exchangeable SMART Vacuum Drying Oven “SOV”, 20 · 30 · 70 Lit, 10~750mmHg, 200℃
 With Smart-Lab™ Controller, 4" Full Touch Screen, Highly Safety View Window, 2 Al-Shelf, Digital PID Control, Superior Temp & Vacuum Accuracy
 Excellent Thermal Conductivity, 가스 교환 스마트 진공 건조기/오븐, 고정밀 디지털 PID 제어, 알루미늄 선반 포함



- Supplied by DAIHAN Korea
- UL, CUL, and CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Suitable for Drying, Baking, Conditioning, Curing, Out-Gassing Solids and Liquids, Vacuum Embedding, Moisture Testing, Plating and Aging Test
- Gas Exchangeable, High Heat Conduction and Precise Gauge
- Smart-Lab™ Controller
 - WiRe™ App Service.
 - Self-Diagnostic Function
 - Variable Program Setting for Multiple Patterns with Multiple Segments
 - The Recorded Data Could be transferred to PC by USB memory
 - Controller SW can be upgraded through USB port.

- * Very Easy and Convenient for Vacuuming and Venting, Easy to read Vacuum Gauge
- * Highly Safety Viewing Window (Tempered Safety Glass+Polycarbonate Window)
- * Heavy-duty Safety Front Door for User Safety
- * Compact New Body Design
- * High Quality Door Sealing ; High Temperature Grade Silicone Molded Gasket
- * Over Temperature & Over Current Protection and Sensor Error Detection

Specifications & Ordering information

<L1>

Spec.	Model	“SOV-20” Max. Capa. 20 Lit	“SOV-30” Max. Capa. 30 Lit	“SOV-70” Max. Capa. 70 Lit
Capacity		20 Lit	30 Lit	70 Lit
Vacuum Range & Level & Rate		10 ~ 750 mmHg, 0.63mbar, 0.39mbar/hr		
Temp	Range / Sensor	R/T +10 ℃ to 200 ℃, PT100		
	Accuracy	± 0.5℃ at 100℃ & 150℃		
	Uniformity	± 2.0℃ at 100℃, ± 3.0℃ at 150℃		
Heating Power/Consumption		170W×2ea 130W×2ea	350W×2ea 450W×2ea	
Heat-up Time		60 min. to 100 ℃, 90 min. to 150 ℃	70 min. to 100 ℃, 100 min. to 150 ℃	80 min. to 100 ℃, 100 min. to 150 ℃
Timer & Alarm		99hr 59min (with delay / Continuous run.), Error status & Timer-end		
Display		4" Full Touch Screen TFT LCD		
Controller		Smart-Lab™ Controller		
Resolution		Control ; ±1.0℃, Display ; ±1.0℃		
Shelf (Included)		2 Aluminum Shelf		
Safety Circuit		Over Temp & Over Current Protector, Sensor Error Detector, Leakage Breaker		
Material		Internal ; Stainless-steel (# 304), External ; Powder Coated steel, Window ; Tempered Safety Glass 12T, Poly Carbonate 5T for Safety		
Others		Program Function, WiRe™ App & Web Service.		
Dimension (w×d×h)	Interior (mm)	265×290×265	300×330×300	400×435×400
	Exterior (mm)	480×475×547	515×515×582	618×620×682
Vacuum & Gas Inlet Valve Size		od. 10 mm		
Net Weight		48 kg	55 kg	85 kg
Packing Size (w×d×h) (mm)		600×580×850	640×640×900	740×720×1000
& Gross Weight		57 kg	80 kg	112 kg
Power Consumption		600 W	600 W	1.6 kW
Electric Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz		

Cat. No.	230V	DH.SWOV03020	DH.SWOV03030	DH.SWOV03070
	120V	DH.SWOV04020	DH.SWOV04030	DH.SWOV04070
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Aluminum Shelf

DH.WOV12020 (1) Standard Shelf, Alu., “OVSN020”, 260×245×h3mm, for 20 Lit
 DH.WOV12030 “OVSN030”, 290×285×h3mm, for 30 Lit
 DH.WOV12070 “OVSN070”, 390×390×h3mm, for 70 Lit

Ea. Price, ₩

※ Note! Optional Vacuum Pump ; See the “Pumps, Vacuum” Section!

WiRe™ App Service



4" Full Touch Screen TFT LCD
 (Smart-Lab™ Controller)

Controlled by Mobile-Phone

It is possible to monitor and control the Smart-Lab™ equipment through all smart phone, tablet without restraint of time and place.



WiRe™ App Service



Pressure Gauge



(1) Standard Shelf, Al.(Included)

DAIHAN WisdPet™ Adjustable Micro & Macro Pipetters, with Ultra Low Dead Air Space, Fully Autoclavable, 0.2 ~ 1000 μ l, 5 & 10ml Fit to Universal-type Tips, Ergonomic Design, CE/ISO/IAF Certified, 정밀형 마이크로- & 매크로- 피펫터, 개별 검교정 성적서 포함



Incl. Calibration Report

- Supplied by DAIHAN Korea ■ With Individual Calibration Report
- Certified ISO 9001, ISO 13485, ISO 17025, ISO 8655 (International Standard Organization), IAF(International Accreditation Forum), CE(Conformity to European)
- Universal Tipcone; Compatible with Most of the Internationally Accepted Standard Tips like Gilson®, ABDOS®, Kartell®, Axygen®, Eppendorf®, etc.
- Highest Accuracy and Precision ■ Easy to Calibrate and Maintenance
- Specially Designed Technology to Ensure Least Dead Air Space to Achieve very High Accuracy and Precision
- * Optimized Spring Tension for Soft Effortless Plunger Movement * Contains a Calibration Certificate
- * Ergonomically Designed for Comfort Using during Long Pipetting Hours
- * Specially Designed Concentric Ejection Force to Ensure Least Tip Ejection Force Requirement
- * Soft Click Sound at Every Increment Ensures Perfect Volume Setting and Prevent any Accidental Volume Change
- * Built-in Streamlined Tip Ejector Facilitates Easy Tip Ejection and comfortable Access to Narrow Necked Bottles and Tubes.
- * Ideal for Pharmaceutical, Academic & Research Institutes and Biotechnology.

(1) Micropipetter

<L1>

Catalog No.	Volume μ l	Division μ l	Vol. tested μ l	Systematic Errors $\pm \mu$ l	Random Errors $\leq \mu$ l	Ea. Price, ₩
DH.Pip301	0.2 - 2	0.01	0.2 1 2	0.04	0.024	
DH.Pip302	0.5 - 10	0.05	2 5 10	0.1	0.05	
DH.Pip303	2 - 20	0.1	2 10 20	0.16	0.08	
DH.Pip304	10 - 100	0.5	10 50 100	0.6	0.2	
DH.Pip305	20 - 200	1	20 100 200	1.2	0.4	
DH.Pip306	100 - 1000	5	100 500 1000	6	2	



Stand for 6 Pipetters

(2) Macropipetter

* Fully Compatible with Gilson Tips

Catalog No.	Volume ml	Division ml	Vol. tested μ l	Systematic Errors $\pm \mu$ l	Random Errors $\leq \mu$ l	Ea. Price, ₩
DH.Pip307	0.5 - 5	0.02	500 2500 5000	30	10	
DH.Pip308	1 - 10	0.05	1000 5000 10000	60	20	


(3) Optional Stand, Material ; Polypropylene

Ea. Price, ₩

DH.Sta301 Table Stand, for 6 Pipetters



Minimum Air Gap between Piston and Tip Liquid



Use Tool for Instant Calibration




Adjust Volume Easily with Plug



Eject Tips Easily with Tip Ejector



Use Various Tips with a Universal Tipcone



Various Colors per Capacity for Easy Identification

Calibration Tool
(Beige Color)



Calibration Nut



Plunger
Calibration Tool
(Highlight in Red)

Easy Calibration Technique

1. Before Calibration, hold the Pipette at a 90 degree Angle, Aspirate the Solution(Distilled Water) to a depth of 3 to 4mm, and then Completely expel at least 3 times.
2. Set the Pipette to 50% of its Maximum Capacity. After Completely Aspirating the Solution, Measure on an Analytical Balance 5 times and Calculate the Average Value.
3. Engage the Calibration Tool, set the Volume of the Pipette to the Calculated average Value, and Remove the Calibration Tool.
4. Set the Volume of the Pipette back 50% of its Maximum Capacity, Repeat 5 times to ensure the Calibrated Volume is within Acceptable Limits.
5. If it is out of Tolerance, Repeat the above Procedure again.

DAIHAN™ 0.098MPa Multiuse Water Jet Vacuum Pump/Aspirator “WEV-003S”

Ideal for Evaporator 1~5Lit, with Multi-Taps & Safety Valve, Bath 15Lit, 다용도 아스피레이터, 워터젯 진공펌프



“WEV-003S” ideal for 1~5Lit Evaporator

3-Year Warranty

- Supplied by DAIHAN Korea
- Resist Corrosion by Acid, Alkali and Solvents.
- Cyclically use Water, Save the Precious Water Resources.
- Good for Evaporation, Distillation, Crystallization, Drying, Sublimation, Filtration &c.
- Fluorine Rubber Sealing
- Need to Replace the Water in the Flume Regularly

* Can be used alone or in Parallel 2 Tap & 2-Display

* Stainless-steel(#304) Body, Copper Ejector, Tee Junction, PP Back Valve

Specifications & Ordering Information

<L1>

Spec.	Model	“WEV-003S” Ideal for 1~5Lit Evaporator
Max. Vacuum Degree		0.098 Mpa
Pump Power & Water Bath Volume		180W & 15Lit
Flow Rate		80 Lit/min
Taps & Air Sucking Amount		2 Tap, 10 Lit/min (1 Tap)
External Circulation Nozzle		id. 5/od. 10mm (id. 8/od.11 mm Silicone tube Connected for Circulation)
Safety Function		Safety Check Valve
Drain Pipe		id. 15mm
Material of Water Tank / Body		PP / PP
Dimension & Net Weight		385×280×h420 mm / 11 kg
Electric Reqts.		230V, 50/60Hz
Cat. No.		DH.WEV0003S
Price, Each. ₩		

DAIHAN™ 105mbar Oil Free Vacuum Pump “PU-300”, with Fine Vacuum Regulator and Moisture Trap, 20Lit/min
With Prote™ Overflow Protection, Built-in Thermal Protection, Protection Class IP30, 진공 펌프



3-Year Warranty



105mbar Oil Free Vacuum Pump
“PU-300”

- 진공 처리가 필요한 다양한 실험에서 사용이 가능 ■ CE, CSA, and RoHS Certified
- Used for General Filtration, Microbiological Test, Suspended Solid Test, Solid-Phase Extraction, Microplate Vacuum Manifold, Air Sampling, etc.

- * Driven by Piston without Lubricant, Regular Oil Change and Maintenance ; No Oil Pollution
- * Prote™ Overflow Protection is Installed to Prevent a Large Amount of Liquid from being Sucked into the Pump.
- * Equipped with Filter Cartridge in Air Inlet to Filter Particle and Moisture to Prolong the Life of Pump.

<L1>

Spec.	Model	“PU-300”	
Speed of Evacuation		20 Lit/min	23 Lit/min
Max. Vacuum		105 mbar, 78.8 torr	
Motor Power & Max. Current		60 W, 0.3A	65 W, 0.7A
Motor Rotation & Noise Level		1450 rpm, 50 dB	1750 rpm, 50 dB
Hose Barb		8mm	
Size & Weight		268 × 135 × h204mm, 4.1kg	
Electric Rqts.		220V, 50Hz	110V, 60Hz
Cat. No.		DH.Pu0300	DH.Pu0301
Price, Each. ₩			

DAIHAN™ 10mbar Anticorrosive Vacuum Diaphragm Pump “PU-20”, with Fine Regulator and Moisture Trap, 18Lit/min
With Built-in Thermal Protection, Quiet and Low Vibration, 내부식성 저소음/저진동 진공 다이어프램 펌프

3-Year Warranty



10mbar Anticorrosive Vacuum
Diaphragm Pump “PU-20”

- 진공 처리가 필요한 다양한 실험에서 사용이 가능하며 최고 10mbar의 고진공 성능을 보여줌
- Used for Solvent purification, Gel dryer, Rotary Evaporator, Vacuum Filtration, etc.
- Resistance to Chemicals ; the Parts in Contact with the Gas are PTFE Coated.

- * Driven by Diaphragm without Lubricant, Regular Oil Change and Maintenance ; No Oil Pollution
- * Quiet and Low Vibration
- * Built-in Thermal Protection

<L1>

Spec.	Model	“PU-20”	
Speed of Evacuation		18 Lit/min	
Max. Vacuum		10 mbar, 7.5 torr	
Motor Power & Max. Current		90 W, 0.5A	
Motor Rotation & Noise Level		1450 rpm, 52dB	
Hose Barb		8mm	
Size & Weight		262 × 236 × h193mm, 7.4kg	
Electric Rqts.		220V, 50Hz	
Cat. No.		DH.Pu8020	
Price, Each. ₩			

DAIHAN™ 15mbar Anticorrosive Vacuum/Pressure Diaphragm Pump “PU-30”, 25psi, 30Lit/min
Ideal for Multi-Purpose, Oil-Free, Double Heads, 내부식성 진공/압력 다이어프램 펌프, 더블헤드

3-Year Warranty



15mbar Vacuum/Pressure
Diaphragm Pump “PU-30”

- 진공 처리가 필요한 다양한 실험에서 사용이 가능하며 부식성 시료의 오염에 강함
- Used for Centrifugal Concentration, Evaporating System, Reduced Pressure Distillation, Solid Phase Extraction, Vacuum Concentration, Vacuum Filtration & Drying System, etc.
- Resistance to Chemicals ; the Parts in Contact with the Gas are PTFE Coated.

- * Tolerance to almost all Strong Acid, Alkali, Oxidizer, Reductant and Organic Solvents
- * No Air Pollution due to the Gas Exchange Room with Filtering Material
- * Automatic Cooling Exhaust System ensures Continuous Operation 24hrs a day.
- * Automatic Power off to Protect Motor when the Pump Body Temp Reached at 130℃.

<L1>

Spec.	Model	“PU-30”	
Speed of Evacuation		30 Lit/min	
Max. Vacuum & Max. Pressure		15 mbar, 11.3 torr & 25 Psi.	
Motor Power		160 W	
Noise Level & Temp of the body		< 60 dB & < 55℃	
Pump Head		2	
Inlet & Outlet (mm)		id. ϕ 6 & 6 mm, od. ϕ 8 & 8 mm	
Diaphragm & Valves		PTFE	
Size & Weight		370 × 144 × h275mm, 13.5kg	
Electric Rqts.		220V, 50Hz	
Cat. No.		DH.Pu8030	
Price, Each. ₩			

TJV® Vacuum/Pressure Diaphragm Pump “GM-0.20”, Oil-Less, $\geq 0.075\text{Mpa}$, 250mbar, $\geq 30\text{ Psi}$, 12 Lit/min
Ideal for Multi-Purpose, 랩용 진공 펌프, 건식/다이어프램 펌프



“GM-0.20”, 12Lit/min, 250mbar

■ Vacuum Desiccators and Cabinets

■ High Efficiency

■ Vacuum Filtration

■ TLC Sprayers

■ Longer Operating Life

* Ideal for Medicinal Products Analysis, Tenuously Chemical Engineering, Biochemical Pharmacy, Food Examination, Investigating and Solving Criminal Case, High-accuracy Chromatogram Apparatus and HPLC Solvent Filter
* Essential Equipment in the Laboratory

Specifications & Ordering information

<L1>

Speed of Evacuation(L/Min) & Noise Level(dB)	Speed ; 12Lit / min, Noise Level ; < 60dB
Temp of the body(℃) & Temp Working Environment(℃)	Body ; < 55℃, Working Environ ; 7 ~ 40
Ultimate Pressure Vacuum / Pressure	$\geq 0.075\text{Mpa}$, 250mbar / $\geq 30\text{ Psi}$
Inlet / Outlet (mm)	$\phi 6 / \phi 6$
Motor Power(W) & Pump Head	Motor Power ; 75W, Pump Head ; Nylon
Dimensions(w×d×h) / Weight	Size ; 195×98×h156mm, Weight ; 4kg
Cat. No.	TJ.Pu0101
Price, Each. ₩	



(1) Single Head, “GM-0.33A”, 20Lit/min

TJV® Anticorrosive Vacuum/Pressure Diaphragm Pump “GM-series”, Oil-Free
Ideal for Multi-Purpose, 200/50mbar, 20~120Lit/min, 내부식성 다이어프램 펌프

■ 진공 처리가 필요한 다양한 실험에서 사용되며 부식성 시료의 오염에 강함

■ PTFE Coated as Pump's Chamber, Diaphragm & Valves, Suitable for Corrosive Gas or Organic Solvents.

■ Used for Centrifugal Concentration, Evaporating System, Reduced Pressure Distillation, Rotary Evaporation, Solid Phase Extraction, Vacuum Concentration, Vacuum Filtration & Drying System, etc.

■ With Quality Reliability, Unique Structural Design to Meet the Various Needs of the Laboratory.

■ Small and Light Weight, Easy to Move, Space Saving, Convenient Maintenance.

* Able to Tolerance almost all Strong Acid, Alkali, Oxidizer, Reductant and Variety of Organic Solvents.

* Ideal for Medicinal Products Analysis, Tenuously Chemical Engineering, Biochemical Pharmacy, Food Examination, Investigating and Solving Criminal Case, High-accuracy Chromatogram Apparatus and HPLC Solvent Filter

Specifications & Ordering Information

<L1>

Model	(1) “GM-0.33A” Max. 20L/min	(2) “GM-0.50A” Max. 30L/min	(3) “GM-0.50B” Max. 30L/min	(4) “GM-1.00” Max. 60L/min	(5) “GM-2” Max. 120L/min
Pump Head	1	1	2	2	2
Speed of Evacuation (L/min)	20	30	30	60	120
Vacuum	200mbar	200mbar	50mbar	200mbar	200mbar
Pressure	None	30 Psi	None	30 Psi	None
Motor Power(W)	160				300
Inlet / Outlet(mm)	O.D $\phi 6 / 6$				O.D $\phi 9 / 9$
Working Temp(℃)	7 ~ 40				
Voltage Rating	230Vac, 50Hz				
Noise Level(dB)	< 60				
Temp of the body(℃)	< 55				
Diaphragm & Valves	NBR(Nitrile Butadiene Rubber) with the PTFE Coating Outside				
Weight (kg)	7.5	8	10	10	20
Dimensions (w×d×h) (mm)	215×120×h235	210×160×h235	300×120×h235	300×160×h235	390×150×h250
Cat. No.	TJ.Pu0107	TJ.Pu0109	TJ.Pu0110	TJ.Pu0111	TJ.Pu0112
Price, Each. ₩					



(2) Single Head, “GM-0.50A”, 30Lit/min



(3) Double Heads, “GM-0.50B”, 30Lit/min



(4) Double Heads, “GM-1.00”, 60Lit/min



(5) Double Heads, “GM-2”, 120Lit/min



(1) VOP-50



(2) VOP-60



(3) VOP-100



(4) VOP-150



(5) VOP-200



(6)

(7)



(8) & (9)
Oil 1 & 4Lit

Vacuum[®] Precision Vacuum Pump “VOP-Series”, Gas Ballast Valve Installed, 42~240 Lit With Double-Stage, Direct Drive & Oil Sealed Rotary Type, 정밀 진공펌프, 이형

■ Supplied by DAIHAN Korea

■ 오일 역류 방지 장치가 장착된 진공펌프

- * Easy-to-use Large Finger Grips
- * Installed Filter Prevent Sludge Contamination
- * Ideal for Vacuum Evaporation Equipment / Vacuum Dryer / Vacuum Freeze Dryer, etc.

■ Oil Back Flow Prevention Mechanism Installed

- * Over Current Protection
- * Physical and Chemical Apparatus Industries

Specifications & Ordering information

<L4>

Model	(1) “VOP-50”	(2) “VOP-60”	(3) “VOP-100”	(4) “VOP-150”	(5) “VOP-200”
Pumping Speed (50/60Hz)	42 / 50 Lit/min	50 / 60 Lit/min	100 / 120 Lit/min	150 / 180 Lit/min	200 / 240 Lit/min
Ultimate Pressure (Pa)	6.7×10 ⁻² Pa				
Vacuum Degree (mbar)	6.7×10 ⁻⁴ mbar				
Motor Power	Single Phase 220V, 200W, 4poles		Single Phase 220V, 400W, 4poles	Single Phase 220V, 550W, 4poles	
Motor Speed(50/60Hz)	1430 / 1730 RPM				
Recommended Oil	COSMA-100				
Oil Capacity(ml)	800		900	1000	1100
Suction Pipe(o.d×i.d)	NW25(KF25)				
Temp Range	7℃~40℃				
Weight	17.5kg	17kg	25kg	28.5kg	30kg
Dimensions(w×d×h) (mm)	170×398×251	165×411×234	171×464×238	189×504×253	189×524×253
Electric Reqts.	220V, 50 / 60Hz				
Cat. No.	V3.VOP050	V3.VOP060	V3.VOP100	V3.VOP150	V3.VOP200
Price, Each. ₩					

* Per Customer's request various specifications are available.

* Using the Macleod gauge the ultimate pressure is measured. The measurements from the pirani gauge are show ln ().

Optional Accessory for VOP-series

Oil Mist Trap

* Available Less than 200L/min Model, Access Port ; NW25(KF25)

V3.OMT150 (6) Oil Mist Trap for VOP-50 / VOP-60 / VOP-100 / VOP-150

* Available Less than 560L/min Model, Access Port ; NW25(KF25)

V3.OMT300 (7) Oil Mist Trap for All VOP-series

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

Vacuum Pump Oil

V3.COS101 (8) Oil for Vacuum Pump VOP-series, 1Lit

V3.COS104 (9) , 4Lit

Nipple

V3.Nip012 (10) Nipple, for NW25 Suction Pipe, ϕ 12mm

V3.Nip015 , ϕ 15mm

V3.Nip018 , ϕ 18mm

Centering & Clamp

V3.Cen0025 (11) Centering, for NW25 Suction Pipe

V3.Cla0025 (12) Clamp, for NW25 Suction Pipe

※ Selection Guide ; We Recommend

“VOP-50” for 20Lit Vacuum Oven, “OV-20” & “SOV-20”

“VOP-60” for 20 & 30Lit Vacuum Oven, “OV-20 · 30” & “SOV-20 · 30”

“VOP-100” for 20~70Lit Vacuum Oven, “OV-20 · 30 · 70” & “SOV-20 · 30 · 70”



(10)



(11)



(12)

REACTORS

BioReactor / Fermentor <i>Coming Soon!</i>	P. 252~255
Pressure Reactors / Stainless-steel Reactors	P. 256
High Pressure Mini Autoclave Reactors of Stainless-steel(STS) ; 50~1000ml ...	P. 256~257
Hydrothermal Synthesis High Pressure Reactors of STS & PTFE ; 25~1000ml	P. 258
Vacuum or Atmospheric Reactors <i>ReNew</i>	P. 259~268
Glass Vacuum Reactors, with DN45° Flat Flanges	P. 259~266
Classic Round Flask Reactor Sets ; 0.25~50Lit	P. 259
Cylindrical Single Wall Vessel Reactor Sets ; 0.1~50Lit	P. 260~261
Double Wall Reactor Vessels - Self Standing Type ; 0.25~20Lit	P. 262
Double Wall Vessel Reactor Sets ; 1~5Lit and 10~100Lit	P. 263~266
- AutoLifting Double Wall Vessel Reactor Sets ; 20~100Lit	P. 266
- Tilting Double Wall Vessel Reactor Sets ; 20~100Lit	P. 265
Stainless-steel Double Wall Vacuum Reactor Sets with Glass Lid	P. 267~268
Reactor Circulators <i>New</i>	P. 269~271
Circulator Hoses	P. 63~64
Cooling / Heating Transfer Fluids	P. 65




DAIHAN™ BIOPRIME™ BioReactor/Fermentor Set, Modular Expansion System for Parallel Workflows Up to 8 Vessels
Master Control for pH · DO · Temp · Agitation · Foam, Single & Jacketed Vessel 1~14Lit, 0~2,000rpm, 5~95 °C
Best for Both Sensitive Animal Cells and Demanding Microbial Regimes, 바이오리액터/발효조 세트, 모듈식 병렬 확장 가능

► **BIOPRIME™ BioReactor / Fermentor Set ;**
Consists of a Controller and Complete Vessels



BIOPRIME™ Controller



BIOPRIME™ Complete Vessel

- **Supplied by DAIHAN Korea**
- 미생물 배양, 포유류 세포 배양, 줄기세포 관련 연구 및 일반적인 실험실 규모 생물공정 개발에 적합
- 단일 용기 실험부터 2-, 4-, 8-Vessel 워크플로우까지 하나의 플랫폼 제품군 내에서 확장 가능
- 규모 확장, 정밀 주입, Bio공정 지원 등 다양한 기능성 옵션 팩을 추가할 수 있는 랙 시스템을 통해 맞춤형으로 사용할 수 있음.
- **Suitable for Microbial Culture, Mammalian Cell Culture, Stem Cell-Related Studies, and General Lab-Scale Bioprocess Development.**
- **Scalable Platform Structure from Routine Single-Vessel Work to 2-, 4-, 8-Vessel Workflows within One Platform Family.**
- **Can be Customized through a Functional Slot-in Rack System that allows for the Addition of Functional Option Packs such as Scale Expansion, Precision Dosing, or Bio-Process Support.**
- * **Practical Parallel Research Capability ;**
Dual and Multi-Vessel Configurations enable Side-by-Side Condition Screening, Process Comparison, and Optimization Studies.
- * **Designed for Daily Bioprocess Research ;**
BIOPRIME Covers the Core Control Loops Most Commonly Required in Laboratory Bioprocessing(pH, DO, Temp, Agitation, and Foam/Level)
- * **Integrated Dosing for Routine and Advanced Studies ;**
Multi-Channel Dosing Supports Feed, Acid, Base, and Antifoam Operation, with Expansion Paths for More Advanced Feeding Workflows.
- * **Flexible Gas Control Strategy ;**
Multi-Gas and MFC-Ready Architecture Supports Standard Aeration Needs and Expandable Gas Control Strategies.
- * **Consistent User Experience Across Tiers ;**
A Common BIOPRIME Operating Concept Helps Users Expand from Standard to Pro and Enterprise with Minimal Retraining.
- * **Ready for Data Handling and Remote Operation ;**
USB Data Storage and Ethernet Connectivity Support Routine Data Handling, Remote Monitoring, and Network-Based Operation.
- * **Easy-to-Select Product Configuration ;**
The Lineup is Organized for Fast Model Selection to Fit Common Vessel and Controller Configurations.

Specifications information

Model	BIOPRIME™ Standard 1-Vessel Set	BIOPRIME™ Pro 4-Vessel Set	BIOPRIME™ Enterprise 8-Vessel Set
Spec.			
Vessel Q'ty	1 Vessel (Up to 2 Vessels)	4-Vessel Parallel	8-Vessel Parallel
Control	Temp ; 5~95°C, Pt100 Sensor Agitation ; 0~2,000rpm, AC Servo Drive, 200W pH ; 0.00~14.00 pH DO ; 0~200% Air Saturation, Agitation / Air~O2 / Gas-Mixing DO Cascade Foam / Level ; Conductivity Probe + Antifoam Control (ON/OFF)		
Integrated Dosing Pump	4 Channels	4 Channels + Add-On (8-Ch)	8 Channels
Pump Speed	Low-Flow ; 5~25rpm, High-Flow 20~100rpm		
Gas lines / MFC	2 Lines ; MFC-Optional	4 Lines ; MFC-Ready	4 Lines ; MFC-Ready
Display	10.1" Touchscreen	15" Touchscreen	17" Touchscreen
Data & Connectivity	USB Data Export + Ethernet		
Dimension(mm)	300×420×h610mm	440×470×h680mm	440×470×h680mm
Weight(kg)	28 kg	36 kg	38 kg
Electrical Reqts.	AC 100~240V, 50/60Hz		


<Continued On Next Page>

<Continued From...DAIHAN™ BIOPRIME™ BioReactor/Fermentor Set>

► BIOPRIME™ BioReactor / Fermentor Set ;
Consists of a Controller and Complete Vessels



BIOPRIME™ Controller



BIOPRIME™ Complete Vessel



Cold Finger

Sparger

Baffle

Ordering information of BIOPRIME™ BioReactor / Fermentor Set

<L1>

Spec.		Model		BIOPRIME™ Standard 1-Vessel Set		BIOPRIME™ Pro 4-Vessel Set		BIOPRIME™ Enterprise 8-Vessel Set	
Vessel Type	Volume	Cat. No.	Price, ₩	Cat. No.	Price, ₩	Cat. No.	Price, ₩	Cat. No.	Price, ₩
Single -Wall	1Lit	DH.BRS.S1		DH.BRP.S1		DH.BRE.S1			
	3Lit	DH.BRS.S3		DH.BRP.S3		DH.BRE.S3			
	5Lit	DH.BRS.S5		DH.BRP.S5		DH.BRE.S5			
	7Lit	DH.BRS.S7		DH.BRP.S7		DH.BRE.S7			
	10Lit	DH.BRS.S10		DH.BRP.S10		DH.BRE.S10			
	14Lit	DH.BRS.S14		DH.BRP.S14		DH.BRE.S14			
Double -Wall (Jacketed)	5Lit	DH.BRS.D5		DH.BRP.D5		DH.BRE.D5			
	7Lit	DH.BRS.D7		DH.BRP.D7		DH.BRE.D7			

※ 운임 별도

I. BIOPRIME™ Complete Vessel (Included)

- O링 기반 모듈식 밀봉은 GMP기준에 부합하는 위생적이고 유지보수가 용이한 구조
- Standard Top-Drive Mechanical Seal Configuration
- Wide Range of Agitation from 10~30rpm Ultra-Gentle Low-Shear Mixing to 1,500~2,000rpm High Intensity Agitation, Supporting both Sensitive Animal Cells and Demanding Microbial Regimes.
- * Clear Cylindrical Viewing Area enhances In-Process Observation and Documentation.
- * O-Ring Based Modular Sealing supports Disassembly, Cleaning, and Component Exchange, enabling a Hygienic, Maintainable Structure aligned with GMP-Oriented Practices.
- Complete Vessel ; Glass Vessel, Head Plate, Bottom Plate, Impeller, Baffle, Sparger, Cold Finger, and Condenser
 - Impeller ; Stirs the Culture Medium to ensure Uniform Mixing and Bubble Dispersion.
 - Baffle ; Suppresses Rotational Flow during Stirring to Improve Mixing Efficiency.
 - Sparger ; Supplies Air or Oxygen to the Culture Medium in the Form of Bubbles.
 - Cold Finger ; Condenses Vapor and Volatile Components through a Cooling Surface.
 - Condenser ; Condenses Vapor Generated during the Fermentation Process into a Liquid.
 - The Specifications of Accessories Vary Depending on the Model.

Spare Glass Vessels

Cat. No.	Vessel Type	Total Vol (L)	Work Vol Max/Min (L)	Glass ID (mm)	Height (mm)	H/D Ratio	Head Plate OD(mm)	Sensor L/N(mm)	Impeller D(mm)	Price, ₩
DH.BVSV010	Single -Wall	1	0.80 / 0.20	100	130	1.3	170	120 / 120	36	
DH.BVSV030		3	2.50 / 0.50	130	140	1.1	190	225 / 225	47	
DH.BVSV050		5	4.00 / 1.00	150	195	1.3	210	325 / 325	54	
DH.BVSV070		7	6.00 / 1.50	160	280	1.8	220	325 / 325	58	
DH.BVSV100		10	8.00 / 2.00	165	360	2.2	230	425 / 425	60	
DH.BVSV140		14	11.0 / 3.0	185	440	2.4	250	425 / 425	67	
DH.BVDJ050	Double -Wall	5	3.5 / 1.40	150	300	2.0	210	325 / 325	54	
DH.BVDJ070		7	5.0 / 0.60	160	345	2.2	220	325 / 325	58	

<Continued On Next Page>

<Continued From...DAIHAN™ BIOPRIME™ BioReactor/Fermentor Set>



BIOPRIME™ Head Plate

II. BIOPRIMER™ Head Plate (Included)

- 센서, 공급, 샘플링, 배기/응축 및 보조 유틸리티 연결을 상단 워크플로우로 구성
- 공정 안정성, 세척 용이성 및 장기적인 신뢰성 최적화
- *Designed as a Modular, Upgrade-Ready Interface that Organizes Critical Process Connections (Sensors, Feeding, Sampling, Exhaust/Condensing, and Auxiliary Utilities) into a Top-Side Workflow.*
- *Optimized for Process Stability, Cleanability, and Long-Term Reliability beyond Simple Port Mounting.*
- * **Hygienic O-Ring Sealing ;**
Non-Welded Sealing Provides Predictable Compression, Easy Inspection, and Cleanability, Making it Ideal for GMP Compliant Operation.
- * **Autoclave-Robust Construction ;**
Materials, Finishes and Joints are Selected to Maintain Seal Integrity and Dimensional Stability during Repeated Sterilization Cycles.
- * **Low-Profile, Unobstructed Top Geometry ;**
Compact Vertical Design improves Accessibility and reduces Interference between Sensors, Tubing, and Auxiliary Hardware in Dense Setups.
- * **Future-Proof Mechanical Interface ;**
Head Plate Dimensions Support Evolving Adapters and Process Modules without Redesigning the Vessel Body, Ensuring Platform Continuity.



(1) Septum Port



(2) Flame Cup Port



(3) Multi Feeding Port



(4) PG13.5 Port



(5) FT8 Port

(1) Septum Port for “Aseptic Inoculation”(Φ 19mm)

- Fast, Aseptic Access for Inoculation and Manual Additions Using a Septum Interface.
- Compatible with 19mm Septum, Blank, Manifold-style Adapters

(2) Flame Cup Port for “Flame-Guarded Aseptic Access”(Φ 19mm)

- For Lab Practices that Rely on Flame-Based Surface Sterilization around the Access Point

(3) Multi Feeding Port(Φ 19mm)

- Consolidates Multiple Addition Lines(Feed, Acid, Base, Antifoam, Trace) into a Single Large Hub via a Manifold-style Adapter.
- Clean Top-Side Routing reduces Tubing Clutter and simplifies SOP Setup.

(4) PG13.5 Port, “Industry-Standard Sensor Interface”

- Primary Port Standard for Core Sensors(pH, DO, etc.) and Organized Port Operation Using Plugs/Adapters
- Quality Assurance Point ; Specify PG13.5-Dedicated Tapping & Gauging to Prevent Mistakes

(5) FT8 Port, “Fast Service + Z-Adjustable Insert Port”(Φ 8mm)

- Collet-Based Feedthrough for Φ 8mm Tube/Rod Inserts (e.g., Sparger Riser, Foam Probe Rod), Optimized for Frequent Removal/Reinstallation and Height(Z) Tuning
- Built-in Z-Adjustment enables Rapid Experimental Tuning where Tip Height Matters(Foam Sensing, Sparger Depth)

<Continued On Next Page>

<Continued From...DAIHAN™ BIOPRIME™ BioReactor/Fermentor Set>



BIOPRIME™ Slot-in Rack System

III. Functional Slot-in Rack System (Option)

■ 다양한 기능성 옵션 팩을 추가하여 맞춤형으로 사용할 수 있는 랙 시스템

(1) Scale & Expansion Pack 규모 확장용 팩

- 단일 설정에서 4, 8개 베셀 규모의 병렬 워크플로우까지 제어시스템을 재구축하지 않고 센서/액추에이터 용량을 추가할 수 있는 모듈식 확장 프레임워크
- A Modular Expansion Framework for Adding More Sensor/Actuator Capacity
- Supports 1/4/8 Vessel Scaling for Phased Expansion without Redesigning the Controller

<L1>

Ea. Price, ₩

DH.BRF.A01 Scale & Expansion Pack

(2) Precision Dosing Pack 정밀 주입용 팩

- 미세한 유량 오차가 세포 생존력, 생산성 또는 재현성에 직접적인 영향을 미치는 경우 정확한 투여량 측정 지원
- Improve Dosing Confidence when Small Flow Errors Directly Affect Viability, Productivity, or Reproducibility.
- A Dosing-Focused Control Bundle ; Multi-Channel, Calibration Workflow, and Optional Micro-Dosing Interfaces for Very Low Flow Additions.
- Ideal for Mammalian/Stem-Cell Work, High-Value Reagents, and Any Process where "Stable Feed" Matters more than Raw Speed.

Ea. Price, ₩

DH.BRF.A02 Precision Dosing Pack

(3) Industrial Integration Pack 산업용 데이터운영 팩

- 데이터, 감사 가능성, 연결성이 제어 범위만큼 중요한 실제 조직 환경에서 운영 지원
- Suitable for Real Organizations where Data, Auditability, and Connectivity are as Important as Scope of Control.
- Ethernet Baseline with SCADA Connectivity Options for Lab-to-Pilot Data Workflows

Ea. Price, ₩

DH.BRF.A03 Industrial Integration Pack

(4) Microbial High-OTR Package 미생물용 팩

- 빠르게 성장하는 미생물 공정을 위해 높은 산소 전달 성능과 안정적인 가스 구동 제어
- Achieve High Oxygen Transfer Performance and Stable Gas-Driven Control for Fast-Growing Microbial Processes
- High-OTR-Ready Gas Delivery and Control, Sparging Strategy Readiness, and Exhaust Stability Hooks for Long High-Load Runs.
- Ideal for E. coli/Yeast Screening, High-Density Growth, or Metabolism-Driven Development

Ea. Price, ₩

DH.BRF.A04 Microbial High-OTR Package

(5) Mammalian Low-Shear + Perfusion-Ready Package 바이오 공정용 팩

- 관류 워크플로우로의 원활한 전환을 준비하는 동안 부드럽고 안정적인 배양 조건 유지
- Maintain Gentle, Stable Culture Conditions while Preparing a Smooth Path toward Perfusion Workflows
- Low-Shear Operation, Overlay-First Headspace Strategy Readiness, Foam/Antifoam Discipline, and Perfusion Loop Readiness
- Ideal for CHO/HEK and Stem-Cell Workflows where "Stability Beats" Aggressiveness

Ea. Price, ₩

DH.BRF.A05 Mammalian Low-Shear + Perfusion-Ready Pack

(6) Advanced DO Cascade Control Pack 고급 DO 캐스케이드 제어 팩

- 교반, 공기 흐름, 산소 농축 및 가스 혼합 로직을 통합적으로 사용하는 고급 DO 캐스케이드 전략 구현
- Extends the Platform to Advanced DO Cascade Strategies using Coordinated Agitation, Airflow, O₂ Enrichment, and Gas-Blending Logic.

Ea. Price, ₩

DH.BRF.A06 Advanced DO Cascade Control Pack

(7) Recipe & Setpoint Profile Control Pack 레시피 및 설정값 프로파일 제어 팩

- 레시피 관리 및 시간 기반 설정값 프로파일링을 지원하여 반복 실험, 설정 시간 단축 및 배치 간 일관성 향상
- Supports Recipe Management and Time-Based Setpoint Profiling for Repeatable Runs
- Reduces Setup Time, and improves Batch-to-Batch Consistency.

Ea. Price, ₩

DH.BRF.A07 Recipe & Setpoint Profile Control Pack

(8) Gravimetric Feed Control Pack 중량식 공급제어 팩

- 옵션 인터페이스 보드 및 외부 저울 통합을 통해 장기간 배양을 위한 중량 기반 공급 제어
- With Optional Interface Boards and External Balance Integration, BIOPRIME can Support Gravimetric Feeding, Weight-Based Addition, and Level-Directed Feed Control for Long-Duration Cultures.


Ea. Price, ₩

DH.BRF.A08 Gravimetric Feed Control Pack


Lab High Pressure Hydrothermal Synthesis Magnetic Stirring Mini Autoclave Reactor Set, 50 · 100 · 250 · 500 · 1000ml Compact Size, Stirring and Reacting Device, Max-1000rpm/220bar, 탁상용 미니 고압 수열합성 오토클레이브 반응기 세트, 자석교반형



Reactor Set



Components



Pressure Gauge
(2 Types Included)



Optional
PTFE Inner Reaction Vessel

■ Supplied by SciLab Korea

■ 최대압력 220bar(22MPa), 최대 온도 200℃와 같이 고압, 고온상에서 물리 화학적 반응을 수행할 수 있는 반응기

■ 용도 - 고온 & 고압 실험, 수소첨가 실험, 중합반응 실험 등에 쓰임 ■ 마그네틱바 포함

■ 사용분야; 실험실, 석유화학, 제약, 화학, 금속공학, 연주소 및 기타분야에서 주로 사용됨

■ Simple High Pressure Reactor with Bottom Heating and Stirring

■ Magnetic Stirrer Bar is Included(φ9×28mm).

■ Application ; High-Temperature & -Pressure Experiment, Hydrogenation, Polymerization, etc.

■ Ideal for Use in the Field of Laboratory, Petrochemical Industry, Chemical Plant, Pharmaceutical, Polymer Synthesis, Metallurgy, Research Institutes, and other Fields.

■ Optional PTFE Inner Reaction Vessel ; Used as Inner Vessel of Stainless-steel Kettle

* Made of Stainless-steel * Can be Equipped with Computer Communication Terminals

* Can be Placed in Electric Blast Drying Oven, Constant Temp Oil / Water Bath, Low Temp Tank and Other Instruments

* LCD Display Speed and Temp · Button to Adjust Parameters * Use More Convenient, Numerical

* Equipped with Original German Imported Miniature DC Motor, Pressure Gauge, Explosionproof Device,

Intake Valve, Sampling Valve and Temp Probe, which can Easily Realize the Micro-Reaction Test

* Kettle Body and the Heater Can be Quickly Separated, while the Kettle Body Can be Quickly Cooled

* The Outer Casing is Made of One-time Flame-retardant Reinforced PBT Injection Molded Casing, which is Resistant to High Temp, Corrosion and Good Insulation Performance

* 30° Bevel Control Panel is Suitable for Sitting and Standing Viewing Angles

* The Body is Designed with an Overflow Tank, and the Liquid Flows through the Overflow Tank to Both Sides of the Machine instead of Flowing to the Body or the Panel

* Steady Low Speed and Strong High Speed, No Extreme Low Speed

* Include Two Interchangeable Pressure Gauge ; 0~2.5MPa and 0~25MPa

※ Use 0~2.5MPa Pressure Gauge Installed as Standard for Low - Pressure Operation, Replace with 0~25MPa Pressure Gauge for High - Pressure Operation

Components

- Stirring Motor ; High Power Motor Imported from Germany, Strong, Stable Performance

- Gas Valve ; Gas Intake Needle Valve

- Sampling Valve ; It is Easy to Sample and Analyze the Reaction Process at any time in the Reaction Process

- Vent Valve ; It is Convenient to Vacuum the System and Replace the Gas before the Reaction, and It is Used as the Vent Valve after the Reaction.

- Pressure Gage ; Real-time Monitoring of Reaction Temperature in the Kettle

- Temperature Sensor ; Deep Inside the Reactor Body to Monitor the Reaction Temperature

- Heating Unit ; Module Heating, Fast Heating, Accurate Control, Digital Display Screen

- Two Type of Pressure Gauges (0~2.5MPa & 0~25MPa) ; Real-time Monitoring Reaction Pressure in the Kettle, Can be Switched Depending on Operating Pressure

Specifications & Ordering information

<L1>

Spec.	Model	TGFY-B-50	TGFY-B-100	TGFY-B-250	TGFY-B-500	TGFY-B-1000
Volume		50ml	100ml	250ml	500ml	1000ml
Working Temp.		200℃				
Working Pressure		220bar (22 MPa)				
Stirring Speed		1000 rpm				
Stirring Method		Internal Magnetic Stirring				
Motor		BLDC Motor				
Motor Power		40 W			70W	100W
Heating Power		600 W		800 W	1500 W	
Heating Mode		Mantle Heating				
Material		Body ; Stainless-steel 316 Valve and Connection ; Stainless-steel 316L Optional PTFE Inner Reaction Vessel				
Working Interface		LCD Display				
Window Diameter		φ40mm	φ44mm	φ58mm	φ70mm	φ82mm
Dimension (w×d×h)		310×330 ×h345mm	320×330 ×h400mm	330×330 ×h525mm	350×330 ×h570mm	370×330 ×h610mm
Net Weight		10kg	14kg	16.5kg	20kg	25kg
Cat. No.		SL.Rea7001	SL.Rea7002	SL.Rea7003	SL.Rea7004	SL.Rea7005
Price, Each. ₩						

※ Note/ Freight and Installation Charges Separately.

Optional PTFE Inner Reaction Vessel

SL.Con4301 PTFE Inner Reaction Vessel, 50ml , φ37×h58mm

SL.Con4302 , 100ml , φ43×h77mm

SL.Con4303 , 250ml , φ57×h102mm


SL.Con4304 , 500ml , φ69×h140mm

SL.Con4305 , 1000ml , φ80×h245mm

Ea. Price, ₩

Lab High Pressure Hydrothermal Synthesis Overhead Stirring Mini Autoclave Reactor Set, 50 · 100 · 250 · 500 · 1000ml

Compact Size, Stirring and Reacting Device, Max-1800rpm/220bar, 탁상용 미니 고압 수열합성 오토클레이브 반응기 세트, 기계식교반형



Components



Optional PTFE Inner Reaction Vessel

■ Supplied by SciLab Korea

■ 최대압력 220bar(22MPa), 최대 온도 350℃와 같이 고압, 고온상에서 물리 화학적 반응을 수행할 수 있는 반응기

■ 용도 - 고온 & 고압 실험, 수소첨가 실험, 중합반응 실험 등에 쓰임

■ 사용분야; 실험실, 석유화학, 제약, 화학, 금속공학, 연구소 및 기타분야에서 주로 사용됨

■ Simple High Pressure Reactor with Bottom Heating and Magnetically Coupled Mechanical Stirring

■ Application; High-Temperature & -Pressure Experiment, Hydrogenation, Polymerization, etc.

■ Ideal for Use in the Field of Laboratory, Petrochemical Industry, Chemical Plant, Pharmaceutical, Polymer Synthesis, Metallurgy, Research Institutes, and other Fields.

■ Optional PTFE Inner Reaction Vessel; Used as Inner Vessel of Stainless-steel Kettle

* Made of Stainless-steel

* Can be Placed in Electric Blast Drying Oven, Constant Temp Oil / Water Bath, Low Temp Tank and Other Instruments

* Can be Equipped with Computer Communication Terminals

* LCD Display Speed and Temp · Button to Adjust Parameters

* Use More Convenient, Numerical

* Equipped with BLDC Motor, Pressure Gauge, Explosionproof Device, Intake Valve, Sampling Valve and Temp Probe, which can Easily Realize the Micro-Reaction Test

* Kettle Body and the Heater Can be Quickly Separated, while the Kettle Body Can be Quickly Cooled

* The Outer Casing is Made of One-time Flame-retardant Reinforced PBT Injection Molded Casing, which is Resistant to High Temp, Corrosion and Good Insulation Performance

* 30° Bevel Control Panel is Suitable for Sitting and Standing Viewing Angles

* The Body is Designed with an Overflow Tank, and the Liquid Flows through the Overflow Tank to Both Sides of the Machine instead of Flowing to the Body or the Panel

* Steady Low Speed and Strong High Speed, No Extreme Low Speed

Components

- Stirring Motor ; High-Power BLDC Motor with Magnetic Drive for Strong and Stable Performance

- Gas Valve ; Gas Intake Needle Valve

- Sampling Valve ; It is Easy to Sample and Analyze the Reaction Process at any time in the Reaction Process

- Vent Valve ; It is Convenient to Vacuum the System and Replace the Gas before the Reaction, and It is Used as the Vent Valve after the Reaction.

- Pressure Gage ; Real-time Monitoring of Reaction Temperature in the Kettle

- Temperature Sensor ; Deep Inside the Reactor Body to Monitor the Reaction Temperature

- Heating Unit ; Module Heating, Fast Heating, Accurate Control, Digital Display Screen

Specifications & Ordering information

<L1>

Spec.	Model	“TGFY-C-50”	“TGFY-C-100”	“TGFY-C-250”	“TGFY-C-500”	“TGFY-C-1000”
Volume		50ml	100ml	250ml	500ml	1000ml
Working Temp.		350℃				
Working Pressure		220bar (22 MPa)				
Stirring Speed		1800 rpm				
Stirring Method		Internal Magnetic Stirring				
Motor		BLDC Motor				
Motor Power		40 W			70W	100W
Heating Power		600 W		800 W	1500 W	
Heating Mode		Mantle Heating				
Material		Body ; Stainless-steel 316 Valve and Connection ; Stainless-steel 316L Optional PTFE Inner Reaction Vessel				
Working Interface		LCD Display				
Window Diameter		φ40mm	φ44mm	φ58mm	φ70mm	φ82mm
Dimension (w×d×h)		310×330 ×h355mm	320×330 ×h410mm	330×330 ×h530mm	350×330 ×h585mm	370×330 ×h630mm
Net Weight		10kg	14kg	16.5kg	20kg	25kg
Cat. No.		SL.Rea7101	SL.Rea7102	SL.Rea7103	SL.Rea7104	SL.Rea7105
Price, Each. ₩						

※ Note! Freight and Installation Charges Separately.

Optional PTFE Inner Reaction Vessel


SL.Con4301 PTFE Inner Reaction Vessel, 50ml , φ37×h58mm
 SL.Con4302 , 100ml , φ43×h77mm
 SL.Con4303 , 250ml , φ57×h102mm
 SL.Con4304 , 500ml , φ69×h140mm
 SL.Con4305 , 1000ml , φ80×h245mm

Ea. Price, ₩

Hydrothermal Synthesis Autoclave Reactor Vessel Set with PTFE Inner Vessel, Made of Stainless-steel 304, 25~1,000ml
 Ideal for Using in Numerous Industries, High Quality PTFE Chamber, +200/220℃ Stable, 수열합성 고압 반응 베셀 세트



Vessel Set



Components

- Supplied by SciLab Korea
- 신뢰성이 높고 내구성이 뛰어난 구조의 스테인리스 스틸을 사용하여 누수 없이 장기간 안정적으로 유지
- 유해 물질의 유출과 오염 방지, 샘플 손실의 최소화
- 용기의 80%만 채우는 것을 권장함
- 용도 - 화학합성, 나노입자 합성, 중합반응, 수열분해, 촉매합성, 수열산화, 수열석출, 결정화과정, 중금속, 내화재료, 유기화학물질 등의 용해등에 쓰임
- 사용분야; 과학 연구소, R&D 연구소, 기관 조직, 기업 품질 관리 부분등 다양한 분야에서 널리 사용됨
- Application - Chemical Synthesis, Synthesizing of Nanoparticles, Polymerization Reaction, Hydrothermal Decomposition, Catalyst Synthesis, Crystallization Process, Dissolve of Heavy Metals, Refractory Material, Organic Chemicals etc.
- Ideal for Used in Field of Laboratory, R&D Labs, Institutional Organization, Quality Control Sections in Business.
- Insulation Resistance, High Lubrication, No Adhesion, Non-Hazardous, Anti-Pollution
- * Recommend Filling out the Reaction Material approx. 80% of Vessel Volume.
- * Dissolve Sample Rapidly without Any Lossless
- * Screw Sealing Type
- * Working Pressure ; ≤3MPa
- * Heating and Cooling Rate ; ≤5℃/min
- * Ultra Corrosion Resistance ; Resistant to Strong Acid, Alkali, Strong base, Aqua Regia and all kinds of Organic Solvent
- * Equipped With Polytetrafluoroethylene Bushing
- * Sturdy Reactor Chamber that Withstands Pressure and Temperature
- * Reliable, Durable Structure, the Sealing Effect Stays Stable for a Long Time without any Leakage

Reactor Vessel Set with PTFE Inner Reaction Vessel

* PTFE ; PolyTetraFluoroEthylene * Max Operating Temp. 220℃ * Safe Temperature 200℃ <L1>

Cat. No.	Capacity	Diameter	Height	Weight	Ea.	Price, ₩
SL.Ves7001	25ml	φ56mm	115mm	1kg		
SL.Ves7002	50ml	φ67mm	155mm	1.8kg		
SL.Ves7003	100ml	φ77mm	165mm	2.8kg		
SL.Ves7004	150ml	φ88mm	180mm	4.8kg		
SL.Ves7005	200ml	φ98mm	200mm	6.5kg		
SL.Ves7006	250ml	φ98mm	220mm	7kg		
SL.Ves7007	300ml	φ110mm	220mm	8kg		
SL.Ves7008	500ml	φ122mm	245mm	10kg		
SL.Ves7009	1000ml	φ158mm	275mm	16.3kg		

Spare; PTFE Inner Reaction Vessel Only (with Lid)

Cat. No.	Capacity	Diameter	Height	Thickness	Weight	Ea.	Price, ₩
SL.Con4201	25ml	φ33mm	64mm	4.5mm	80g		
SL.Con4202	50ml	φ40mm	86mm	5mm	112g		
SL.Con4203	100ml	φ50mm	98mm	5mm	150g		
SL.Con4204	150ml	φ56mm	118mm	6mm	292g		
SL.Con4205	200ml	φ63mm	135mm	7mm	477g		
SL.Con4206	250ml	φ65mm	152mm	7mm	560g		
SL.Con4207	300ml	φ70mm	156mm	7mm	644g		
SL.Con4208	500ml	φ85mm	158mm	8mm	720g		
SL.Con4209	1000ml	φ110mm	190mm	8mm	1200g		

SciLab™ 0.25~50Lit Classic Glass Vacuum Round Flask Reactor Sets, 45° DN Grooved Flange & FEP O-Ring, Boro-glass 3.3 With 450°C Heating Mantles, Digital Hi-Torque Stirrer Kit, and Full Support Frame & Access., 진공 글라스 라운드 플라스크 반응조 세트



■ Supplied by SciLab Korea

▶ 45° DN-Standard Flanged Lids, with 4 · 5 · 6 Necks ;
Perfect Compatibility with Flasks · Cylindrical Vessels · Quick Release Clamps

Lids	45° DN Flange	Necks	Center-Joint	Side-Joint
DN60	id 60 × od 100mm	4	24/40	2 × 15° 24/40, 1 × Vertical 14/23
DN100	id 100 × od 138mm	5	24/40	3 × 15° 24/40, 1 × Vertical 24/40
DN120	id 120 × od 158mm	5	34/45	3 × 15° 24/40, 1 × Vertical 24/40
DN150	id 150 × od 184mm	5	34/45	2 × 15° 24/40, 1 × 15° 34/35, 1 × Vertical 24/40
DN200	id 200 × od 242mm	6	34/45	1 × 15° 34/45, 1 × 15° 24/40, 1 × 15° 45/40, 2 × Vertical 24/40

* Side 34/35 & 45/40 Neck ; Good For Inlet-Outlet of Substances.

▶ High Torque Digital Stirring System ; Overhead Stirrer "SL.SOS01051" & "SL.SOS01071"

■ PL(Product Liability) Insurance

■ Patented Jog-Shuttle Control System

■ Controlled by Serial Number

- For 0.25~30Lit ; 60,000mPas Viscosity, 180Ncm Torque, 50~1000rpm, 230V "SL.SOS01051"

- For 50Lit ; 90,000mPas Viscosity, 210Ncm Torque, 50~1000rpm, 230V "SL.SOS01071"

- With Safety Flex Coupling between Stir-Motor and PTFE Shaft / Impeller, with PTFE Stir-Guide(Seal)

▶ Set Components ; (1) Digital Stirring System Full Set, (2) 45° DN Flanged Reaction Flask & Lid, (3) Quick Release Clamp, (4) Condenser, (5) Dropping Funnel, (6) A Mercury Bar Thermometer 250°C, (7) Flexible PTFE Joint Connector(Bellows), (8) PTFE Joint Clamp, (9) Necessity Clamps All {ex. ① 1 or 2 × STS Band Clamp for Vessel Neck & Condenser, ② 4 or 8 × Clamp Holders ϕ 25mm Grip, etc.}, (10) Fully Heavy Stand Frame & Access ; ① SL.ReF7101(Height 1200mm) for 0.25~30Lit and ② SL.ReF7102(Height 1500mm) for 50Lit,

(11) Heating Mantle ; ① Analog Dial Controlled for I. Classic Model


② Digital Remote Controlled with Direct Contact Thermo Probe for II. Premium Model

I. Classic Model, with 450°C Analog Dial Controlled Heating Mantle & Standard Temp Level by a Bar(Stem) Thermometer


기본형(경제형) 모델, 작업온도의 기준은 봉상온도계의 눈금 온도임.

II. Premium Model, with 450 ± 1°C Digital Remote Controlled Heating Mantle With Direct Contact Thermo Probe K-Type/PTFE Encapsulated.

정밀 고급형 모델, 1°C 단위로 반응 물질의 온도를 직접 조절(접촉 온도 센서 사용)



I. Classic Model 경제형



II. Premium / Precise Model 정밀형 Remote Controlled

ORDERING information

<L1>

		I. Classic-Model, 경제형			II. Premium-Model, 정밀 고급형		
Flask Reactor Set 45° DN Flanged		With Analog Dial Controlled Heating Mantle, 450°C, Temp Level by a Bar Thermometer			With Digital Remote Controlled Heating Mantle, 450°C ± 1°C by Direct Contact Thermo-Probe		
Capacity ml & Lit	Glass Lid & Neck	Cat. No.	Price, ₩	included Mantle No.	Cat. No.	Price, ₩	included Mantle No.
250ml	DN60, 4-Neck	SL.Rea01		SL.Man8032	SL.Rea01P		SL.Man8072
500ml		SL.Rea02		SL.Man8033	SL.Rea02P		SL.Man8073
1Lit	DN100, 5-Neck	SL.Rea03		SL.Man8034	SL.Rea03P		SL.Man8074
2Lit		SL.Rea04		SL.Man8035	SL.Rea04P		SL.Man8075
3Lit	DN120, 5-Neck	SL.Rea05		SL.Man8036	SL.Rea05P		SL.Man8077
5Lit		SL.Rea06		SL.Man8038	SL.Rea06P		SL.Man8079
10Lit	DN150, 5-Neck	SL.Rea07		SL.Man8040	SL.Rea07P		SL.Man8081
20Lit		SL.Rea08		SL.Man8041	SL.Rea08P		SL.Man8082
30Lit	DN200, 6-Neck	-	-	-	SL.Rea09P		SL.Man8085
50Lit		-	-	-	SL.Rea10P		SL.Man8083

▶ Optional Vacuum Pumps ⇒ See the "Pumps, Vacuum" Section !

▶ ALL of the Components & Necessity Accessories

⇒ Visit DAIHAN-Sci.com Shopping-mall !

▶ The Details of Mantles and Stirring System

⇒ See the "Heating Mantles" & "Stirrer" Section !

SciLab™ 0.1~20Lit Glass Vacuum Cylindrical Single Wall Reactor Sets, 45° DN Grooved Flange & FEP O-Ring, Boro-glass 3.3 With 450°C Heating Mantles, Digital Hi-Torque Stirrer Kit, and Full Support Frame & Access., 진공 글라스 베셀 반응조 세트, 싱글월, 평바닥



■ Supplied by SciLab Korea

▶ 45° DN-Standard Flanged Lids, with 4 · 5 · 6 Necks ;
Perfect Compatibility with Flasks · Cylindrical Vessels · Quick Release Clamps

Lids	45° DN Flange	Necks	Center-Joint	Side-Joint
DN60	id 60 × od 100mm	4	24/40	2 × 15° 24/40, 1 × Vertical 14/23
DN100	id 100 × od 138mm	5	24/40	3 × 15° 24/40, 1 × Vertical 24/40
DN120	id 120 × od 158mm	5	34/45	3 × 15° 24/40, 1 × Vertical 24/40
DN150	id 150 × od 184mm	5	34/45	2 × 15° 24/40, 1 × 15° 34/35, 1 × Vertical 24/40
DN200	id 200 × od 242mm	6	34/45	1 × 15° 34/45, 1 × 15° 24/40, 1 × 15° 45/40, 2 × Vertical 24/40

* Side 34/35 & 45/40 Neck ; Good For Inlet-Outlet of Substances.

▶ High Torque Digital Stirring System ; Overhead Stirrer (Cat. No. "SL.SOS01051")

■ PL (Product Liability) Insurance

■ Patented Jog-Shuttle Control System

■ Controlled by Serial Number

- 60,000mPas Viscosity, 180Ncm Torque, 50~1000rpm, 230V "SL.SOS01051"

- With Safety Flex Coupling between Stir-Motor and PTFE Shaft / Impeller, with PTFE Stir-Guide(Seal)


▶ Set Components ; (1) Digital Stirring System Full Set, (2) 45° DN Flanged Cylindrical Single Wall Vessel & Lid(5~20Lit Vessels ; Kettle-type), (3) Quick Release Clamp, (4) Condenser, (5) Dropping Funnel, (6) A Mercury Bar Thermometer 250°C, (7) Flexible PTFE Joint Connector(Bellows), (8) PTFE Joint Clamp, (9) Necessity Clamps All {ex. ① 1 or 2 × STS Band Clamp for Vessel Neck & Condenser, ② 4 or 8 × Clamp Holders ϕ 25mm Grip, etc.}, (10) Fully Heavy Stand Frame & Access ; ① SL.ReF7101(Height 1200mm) for Up To 10Lit and ② SL.ReF7102(Height 1500mm) for 20Lit,

(11) Heating Mantle ; ① 450°C Analog Dial Controlled for I. Classic Model


② 450±1°C Digital Remote Controlled with Direct Contact Thermo Probe for II. Premium Model

I. Classic Model, with 450°C Analog Dial Controlled Heating Mantle & Standard Temp Level by a Bar(Stem) Thermometer
기본형(경제형) 모델, 작업온도의 기준은 봉상온도계의 눈금 온도임.

II. Premium Model, with 450±1°C Digital Remote Controlled Heating Mantle With Direct Contact Thermo Probe K-Type/PTFE Encapsulated.
정밀 고급형 모델, 1°C 단위로 반응 물질의 온도를 직접 조절(접촉 온도 센서 사용)



I. Classic Model 경제형



II. Premium / Precise Model 정밀형 Remote Controlled

ORDERING information

<L1>

		I. Classic-Model, 경제형			II. Premium-Model, 정밀 고급형		
Cylindrical Single Wall Reactor Set 45° DN Flanged		With Analog Dial Controlled Heating Mantle, 450°C, Temp Level by a Bar Thermometer			With Digital Remote Controlled Heating Mantle, 450°C ± 1°C by Direct Contact Thermo-Probe		
Capacity ml & Lit	Glass Lid & Neck	Cat. No.	Price, ₩	included Mantle No.	Cat. No.	Price, ₩	included Mantle No.
100ml	DN60, 4-Neck	SL.Rea11		SL.Man8130	SL.Rea11P		SL.Man8170
250ml		SL.Rea12		SL.Man8131	SL.Rea12P		SL.Man8171
500ml	DN100, 5-Neck	SL.Rea13		SL.Man8015	SL.Rea13P		SL.Man8172
1Lit		SL.Rea14		SL.Man8133	SL.Rea14P		SL.Man8173
2Lit	DN120, 5-Neck	SL.Rea15		SL.Man8134	SL.Rea15P		SL.Man8174
3Lit		SL.Rea16		SL.Man8136	SL.Rea16P		SL.Man8176
5Lit	DN150, 5-Neck	SL.Rea17		SL.Man8138	SL.Rea17P		SL.Man8178
10Lit		SL.Rea18		SL.Man8141	SL.Rea18P		SL.Man8181
20Lit	DN200, 6-Neck	SL.Rea19		SL.Man8142	SL.Rea19P		SL.Man8182


▶ Optional Vacuum Pumps ⇒ See the "Pumps, Vacuum" Section !

▶ ALL of the Components & Necessity Accessories
⇒ Visit DAIHAN-Sci.com Shopping-mall !


▶ The Details of Mantles and Stirring System
⇒ See the "Heating Mantles" & "Stirrer" Section !

SciLab™ 0.5~50Lit Glass Vacuum Drain Valved Single Wall Reactor Sets, 45° DN Grooved Flange & FEP O-Ring, Boro-glass 3.3 With $450 \pm 1^\circ\text{C}$ Heating Mantles, Digital Hi-Torque Stirrer Kit, and Full Support Frame & Access., 배출밸브형 진공 반응조 세트, 맨틀가열, 환형 & 비커형

ReNew



**450 ± 1 °C Controlled
Mantle Heated Reactor Set,
I. Round Flask Type ; 0.5~50Lit
II. Cylindrical Vessel Type ; 0.5~20Lit**



■ Supplied by SciLab Korea

▶ **45° DN-Standard Flanged Lids of Glass & PTFE, with 4 · 5 · 6 Necks ;**

- Perfect Compatibility with Flasks · Vessels · Quick Release Clamps
- PTFE Lids ; Interchangeable Screw Joints, Chemically Inert, $-200^\circ\text{C} + 260^\circ\text{C}$, Non-Breakage

Lids	45° DN Flange	Necks	Center Joint	Side Joint
DN60	id 60 × od 100mm	4	24/40	2 × 15° 24/40, 1 × Vertical 14/23
DN100	id 100 × od 138mm	5	24/40	3 × 15° 24/40, 1 × Vertical 24/40
DN120	id 120 × od 158mm	5	34/45	3 × 15° 24/40, 1 × Vertical 24/40
DN150	id 150 × od 184mm	5	34/45	2 × 15° 24/40, 1 × 15° 34/35, 1 × Vertical 24/40
DN200	id 200 × od 242mm	6	34/45	1 × 15° 34/45, 1 × 15° 24/40, 1 × 15° 45/40, 2 × Vertical 24/40

* Side 34/35 & 45/40 Neck ; Good For Inlet-Outlet of Substances.

▶ **High Torque Digital Stirring System ;**

- **Patented Jog-Shuttle Control System** ■ **PL(Product Liability) Insurance**
- For 0.5~30Lit ; 60,000mPas Viscosity, 180Ncm Torque, 50~1000rpm, 230V "SL.SOS01051"
- For 50Lit ; Remote Controlled "SHT-120T" Stirrer (90,000mPas / 210Ncm)
- With Safety Flex Coupling between Stir-Motor and PTFE Shaft / Impeller, with PTFE Stir-Guide(Seal)

▶ **Set Components ;** (1) Digital Stirring System Full Set, (2) 45° DN Flanged Single Wall Flask / Vessel & Lid, (3) Quick Release Clamp, (4) Condenser, (5) Dropping Funnel, (6) A Mercury Bar Thermometer 250°C , (7) Flexible PTFE Joint Connector(Bellows), (8) PTFE Joint Clamp, (9) Necessity Clamps All, (10) Fabricated Support Frame(STS Pipe $\phi 27.2 \times$ Thick 2mm) & Access. Fully, (11) Digital Remote Controlled Heating Mantle, $450 \pm 1^\circ\text{C}$ with Direct Contact Thermo Probe K-type PTFE Encapsulated

■ **With Effluent Collection Safety STS Tray ; Used in Case of Vessel Failure.**
유출물 포집 트레이 ; 반응용기 파손 등 고장시 First Aid 용 !

I. Round Flask-Type, 0.5~50Lit Set ; ORDERING

<L1>

Round Reactor Set 45° DN Flanged		With Glass Lid		With PTFE Lid		included Mantle No.
		Cat. No.	Price, ₩	Cat. No.	Price, ₩	
0.5Lit	DN60, 4-Neck	SL.Rea21		SL.Rea21P		SL.Man8220
1Lit	DN100, 5-Neck	SL.Rea22		SL.Rea22P		SL.Man8221
2Lit		SL.Rea23		SL.Rea23P		SL.Man8222
3Lit	DN120, 5-Neck	SL.Rea24		SL.Rea24P		SL.Man8223
5Lit		SL.Rea25		SL.Rea25P		SL.Man8225
10Lit	DN150, 5-Neck	SL.Rea26		SL.Rea26P		SL.Man8226
20Lit	DN200, 6-Neck	SL.Rea27		SL.Rea27P		SL.Man8227
30Lit		SL.Rea28		SL.Rea28P		SL.Man8229
50Lit	DN200, 6-Neck	SL.Rea29		SL.Rea29P		SL.Man8228

II. Cylindrical-Type, 0.5~20Lit Set ; ORDERING

Cyl-Reactor Set 45° DN Flanged		With Glass Lid		With PTFE Lid		included Mantle No.
		Cat. No.	Price, ₩	Cat. No.	Price, ₩	
0.5Lit	DN100, 5-Neck	SL.Rea31		SL.Rea31P		SL.Man8190
1Lit		SL.Rea32		SL.Rea32P		SL.Man8191
2Lit	DN120, 5-Neck	SL.Rea33		SL.Rea33P		SL.Man8192
3Lit		SL.Rea34		SL.Rea34P		SL.Man8194
5Lit	DN150, 5-Neck	SL.Rea35		SL.Rea35P		SL.Man8196
10Lit		SL.Rea36		SL.Rea36P		SL.Man8199
20Lit	DN200, 6-Neck	SL.Rea37		SL.Rea37P		SL.Man8200

▶ **Choice of Frames 반응조 용량에 따라 사용되는 "프레임 종류"**

- Cyl-type 0.5~20Lit and Round-type 0.5~10Lit ; "ReF7104" $450 \times 450 \times h1500\text{mm}$
- Round-type 20 & 30Lit ; "ReF7105" $550 \times 550 \times h1700\text{mm}$
- Round-type 50Lit ; "ReF7106" $670 \times 670 \times h2000\text{mm}$


- ▶ Optional Vacuum Pumps ⇒ See the "Pumps, Vacuum" Section !
- ▶ ALL of the Components & Necessity Accessories ⇒ Visit DAIHAN-Sci.com Shopping-mall !
- ▶ The Details of Mantles and Stirring System ⇒ See the "Heating Mantles" & "Stirrer" Section !

SciLab™ 0.25~20Lit Glass Double Wall Self-Standing Vacuum Reactor Sets, Graduated, 45° DN Grooved Flange & FEP O-Ring With Digital Stirring System · Fabricated Heavy Duty Stand Frame · Tray, 0.25~20Lit 글라스 진공 이중 반응조 세트, “자립형”, 눈금부, 내열 연결호스 포함



0.25~20Lit Self-Standing Reactor Set

I. with Glass Lid
II. with PTFE Lid



■ **Supplied by SciLab Korea**

► Self-Standing Glass Double Wall Vacuum Vessels & Lids(Glass & PTFE)

- * Boro-glass 3.3, Approx Thick 4.5mm(0.25~5Lit) & 5.5mm(10 & 20Lit)
- * 45° DN-60/-100/-120/-150/-200 Flanged, Fit to 3-Segment STS Quick-Release Clamp
- * Precise Temp Control ; Direct Contact PTFE Capsulated STS Thermo-Probe K-Type
- * -200℃ +260℃ All PTFE Screw GL18 & GL32 Connect-Port for the Circulation Hose of Heat Transfer Fluid of Heating / Cooling Circulators

■ **VITON Circulation Hose Set Included.**

► High Torque Digital Stirring System ; 180Ncm, 60000mPas, 50~1000rpm

- **CE Certified** ■ **With Certificate & Traceability ; Controlled by Serial Number**
- **Patented Jog-Shuttle Control System** ■ **PL(Product Liability) Insurance**
- **With Safety Flex-Coupling between Stir-Motor and PTFE Shaft, PTFE Impeller.**
- * The Stirrer Model ; “SHT-50DX” Direct Controlled(SL.SOS01051), 85W, 230V, 50/60Hz(1 Phase)
- * Digital LCD Display Back Light 50~1000/1rpm, Timer & Alarm 99hr 59min & Error
- * Smooth Start and Change of Stirring Speed * Chuck - Corrosion Free
- * Storage Function for the Set Values ; RPM and Timer
- * Locking Mode for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Safety Device - Overload & Heat Protector for Motor

► Safety Work Fabricated Stand Frame System

- * Heavy Duty Casting Iron Base(40~69cm Adjustable Feet Distance) with 3×OD ϕ23mm STS Pipe(T. 2) Rod-Height 1200 & 1500mm and 2×Silicone Rubber Padded STS Band Clamps
- * Frame Models ; (a) SL.ReF7101(Height 1200mm Type) for 0.25~5Lit Reactors
(b) SL.ReF7102(Height 1500mm Type) for 10 & 20Lit Reactors

■ **With Effluent Collection Safety STS Tray ; Used in Case of Vessel Failure.**

유출물 포집 트레이 ; 반응용기 파손 등 고장시 First Aid 용 !

A. 45° DN Lids of Glass & PTFE(-200℃+260℃), 4 · 5 · 6-Socket Joint(Neck)

* Opaque PTFE Lids ; Interchangeable Screw Joints, Non-Breakage, Non-Scratch, Chemically Inert

Lid	Neck	Center Joint	Side Joints	Using
DN60	4	24/40	2 × 15° 24/40, 1 × Vertical 14/23	For 0.25Lit Vessels
DN100	5	24/40	3 × 15° 24/40, 1 × Vertical 24/40	For 0.5 & 1Lit Vessels
DN120	5	34/45	3 × 15° 24/40, 1 × Vertical 24/40	For 2 & 3Lit Vessels
DN150	5	34/45	2 × 15° 24/40, 1 × 15° 34/35, 1 × Vertical 24/40	For 5Lit Vessels
DN200	6	34/45	1 × 15° 34/45, 1 × 15° 24/40, 1 × 15° 45/40, 2 × Vertical 24/40	For 10 & 20Lit Vessels

* Side 34/35 & 45/40 Neck ; Good For Inlet-Outlet of Substances.

B. 0.25~20Lit Self-Standing Glass Double Wall Vessels

■ **With 45° DN Grooved Flange & Teflon FEP O-Ring, Graduated**

Capacity	Flange	Dimension	Port	Capacity	Flange	Dimension	Port
250/50ml	DN60	ϕ 105 × h150mm	GL18	3000/500ml	DN120	ϕ 160 × h340mm	GL18
500/100ml	DN100	ϕ 130 × h140mm	GL18	5000/500ml	DN150	ϕ 215 × h430mm	GL18
1000/100ml		ϕ 130 × h210mm	GL18	10/1Lit	DN200	ϕ 250 × h450mm	GL32
2000/200ml	DN120	ϕ 160 × h250mm	GL18	20/2Lit		ϕ 300 × h570mm	GL32

ORDERING for Reactor Sets, Without Heating Circulator & Vacu-Pump

<L1>

Reactor Set	I. With Glass Lid		II. With PTFE Lid	
	Cat. No.	Price, ₩	Cat. No.	Price, ₩
250 / 50ml, DN60	SL.Rea41		SL.Rea41P	
500 / 100ml, DN100	SL.Rea42		SL.Rea42P	
1 / 0.1Lit, DN100	SL.Rea43		SL.Rea43P	
2 / 0.2Lit, DN120	SL.Rea44		SL.Rea44P	
3 / 0.5Lit, DN120	SL.Rea45		SL.Rea45P	
5 / 0.5Lit, DN150	SL.Rea46		SL.Rea46P	
10 / 1Lit, DN200	SL.Rea47		SL.Rea47P	
20 / 1Lit, DN200	SL.Rea48		SL.Rea48P	

► **Lids · Spare All PTFE Screw GL18 & 32 Connector Kit · Spare Teflon FEP Sealing O-Ring**
⇒ Visit DAIHAN-Sci.com Shopping-mall !

► **Optional Heat/Cooling-Circulators & -Fluid**
⇒ See the “Bath Circulators, Accessories” Section !

► **Optional Vacuum Pumps** ⇒ See the “Pumps, Vacuum” Section !


SciLab™ 1~5Lit Standard Glass Vacuum Double Wall Reactor Set, Graduated, 45° DN Grooved Flange & FEP O-Ring, 5-Neck Lid With Digital Stirring System · PTFE Drain Valve · Tray · Fabricated Stand Frame, 1~5Lit 글라스 진공 2중 반응조 세트, 눈금부, 내열 연결호스 포함

ReNew

1~5Lit Double Wall Reactor Set

I. with Glass Lid

II. with PTFE Lid



Advantage of 45°DN Flange ;

1. Anti-Breakage
2. Excellent Vacuum Operation
3. Easy & Quick Installation / Removal



I. With Glass Lid
With 5 Joint Necks

II. With All PTFE Lid
* Interchangeable Screw
5 Joints, -200+260°C
* Top Safety ; Non-Breakage,
Chemically Inert

■ Supplied by SciLab Korea

▶ 1~5Lit Graduated Double-wall Vacuum Glass Vessel & 5-Neck Lid

■ 45°DN Flange & 5mm Thick High Boro-glass 3.3 with DN100 · DN120 · DN150 Flange

■ With Direct Contact PTFE Thermo-Probe K-Type, FEP O-Ring for Vacuum, PTFE Drain Valve, and Insulated VITON Circulation Hose Set.

▶ High-Torque Digital Agitator System, 50~1000rpm / 180Ncm

■ CE Certified ■ With Certificate & Traceability ; Controlled by Serial Number

■ Patented Jog-Shuttle Control System ■ PL(Product Liability) Insurance

■ With Safety Flex-Coupling between Stir-Motor and PTFE Shaft, PTFE Impeller.

* Digital LCD with Back-Light Function * Smooth Start and Change of Stirring Speed

* Storage Function for the Set Values ; RPM and Timer * Chuck - Corrosion Free

* Locking Mode for Experimental Safety (Input to Jog-Shuttle can be Disabled)

▶ Safety Work Fabricated Stand Frame System

* Heavy Duty Casting Iron Base(40~69cm Adjustable Feet Distance) with 3 × OD ϕ 23mm STS Pipe(T.2) Rod-Height 1200mm, 2 × Adjustable STS Band Clamp, 2(or 1) × STS O-Ring / T.12mm Supporter with Dual Stem, and Necessity Access.

* Frame Model ; (a) "SL.ReF7101" (Height 1200mm Type) for 1~3Lit.

(b) "SL.ReF7102" (Height 1500mm Type) for 5Lit.

* Vessel Support STS Ring with 4 Layer Cushion Pad(Partially) of 250°C Silicone Adhesive Glass Fiber Tape "3M 361"

■ With Effluent Collection Safety STS Tray ; Used in Case of Vessel Failure

유출물 포집 트레이 ; 반응용기 파손 등 고장시 First Aid 용 !

1~5Lit Specification

Models		“REA-1C”	“REA-2C”	“REA-3C”	“REA-5C”
Working Volume		1Lit	2Lit	3Lit	5Lit
DN Flange		DN 100	DN 120		DN 150
Jacket Volume Approx.		0.3Lit	0.6Lit	0.9Lit	1.5Lit
Vessel OD×h mm		150×250	180×280	200×300	230×400
Heat Hose Connector		All PTFE Screw GL 18			
Effluent Tray		2.7Lit “DM.Tra721”		5Lit “DM.Tra734”	
Digital Stirrer	Model	“SHT-50DX” Direct Controlled (SL.SOS01051) 60,000 mPas Viscosity, 180 Ncm Torque			
	Speed	50~1,000 rpm, Resolution ; 1 rpm-Display / 1 rpm-Control			
	Motor & Rating	Worm-Geared DC Motor (3.5:1 Ratio)			
		In · Output 83 / 75 W, 1 / 10 HP			
	Display	Digital LCD with Back-Light Function			
	Controller	Digital Feedback Controller with Jog Shuttle Switch (Turn + Push)			
Power		Consumption 85 W			
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end			
Shaft / Impeller		φ 10/8 PTFE Encapsulated Stainless-steel Shaft			
Safety Device		Overload & Heat Protector for Motor			
Dimension Overall		35×35×h120cm		40×40×h150cm	
Power Supply		1 Phase, AC 230 V, 50 / 60 Hz			

ORDERING for Reactor Set, without Circulator and Vacu-Pump

<L1>

Reactor Sets	I. With Glass Lid, 5 Necks		II. With PTFE Lid, 5 Necks	
	Cat. No.	Price, ₩	Cat. No.	Price, ₩
1 / 0.1Lit, DN100	SL.Rea51		SL.Rea51P	
2 / 0.2Lit, DN120	SL.Rea52		SL.Rea52P	
3 / 0.5Lit, DN120	SL.Rea53		SL.Rea53P	
5 / 0.5Lit, DN150	SL.Rea54		SL.Rea54P	

▶ Lid Details ; Visit DAIHAN-Sci.com Shopping-mall !

▶ Optional Heat/Cooling-Circulators & -Fluid

⇒ See the "Bath Circulators, Accessories" Section !

▶ Optional Vacuum Pumps ⇒ See the "Pumps, Vacuum" Section !

SciLab™ 10~100Lit Standard Glass Vacuum Double Wall Reactor Set, Graduated, 45° DN Grooved Flange & FEP O-Ring, 6-Neck Lid With Digital Stirring System · PTFE Drain Valve · Tray · Fabricated Heavy Work Frame, 10~100Lit 글라스 진공 2중 반응조 세트, 눈금부, 내열 연결호스 포함

ReNew

10~100Lit Double Wall Reactor Set
I. with Glass Lid
II. with PTFE Lid



Heavy-work Support Frame Dimension
φ 27.2 × Thick-1.5(10 & 20Lit) & -2(30~100Lit)mm STS Pipe

- * 10 & 20Lit ; "SL.ReF7107" 45 × 45 × h150cm
- * 30 & 50Lit ; "SL.ReF7108" 55 × 55 × h170cm
- * 100Lit ; "SL.ReF7109" 67 × 67 × h200cm

■ Supplied by SciLab Korea

► Graduated Double-wall Vacuum Glass Vessel & 6-Neck Lid

- 45° DN 200 Flange & 6/7mm Thick High Boro-glass 3.3 with Good Physical, Chemical Properties
- With Direct Contact PTFE Thermo-Probe K-Type, FEP O-Ring for Vacuum, PTFE Drain Valve, and Insulated VITON or Flex STS(Bellow) Hose Set

► High-Torque Digital Agitator System, 50~1000rpm(800rpm/100Lit)

- CE Certified ■ With Certificate & Traceability ; Controlled by Serial Number
- Patented Jog-Shuttle Control System ■ PL(Product Liability) Insurance
- With Safety Flex-Coupling between Stir-Motor and PTFE Shaft, PTFE Impeller.
- * Digital LCD with Back-Light Function * Smooth Start and Change of Stirring Speed
- * Storage Function for the Set Values ; RPM and Timer * Chuck - Corrosion Free
- * Locking Mode for Experimental Safety (Input to Jog-Shuttle can be Disabled)

► Heavier-Work Fabricated Frame System Adjustably ; 조립 조절식 견고 프레임

- Main Material ; OD φ 27.2 × Thick-1.5mm(10 & 20Lit) & -2mm(30~100Lit) STS Pipe
- * With 4 Stop-on Castors, STS Ring Supporters & Clamps, and All Access.
- * Vessel Support STS Ring with 4 Layer Cushion Pad(Partially) of 250 °C Silicone Adhesive Glass Fiber Tape "3M 361"

- With Effluent Collection Safety Tray STS ; Used in Case of Vessel Failure
- 유출물 포집 트레이 ; 반응용기 파손 등 고장시 First Aid 용 !

Specification of Reactor Set, DN 200 Flange & 6-Neck Lid, 10~100Lit

Models		“REA-10C”	“REA-20C”	“REA-30C”	“REA-50C”	“REA-100C”
Working Volume		10Lit	20Lit	30Lit	50Lit	100Lit
Jacket Volume Approx.		2.6Lit	5Lit	8Lit	12Lit	25Lit
Vessel OD×h mm		250×480	300×630	365×710	410×800	500×970
Inner Vessel OD mm		215	250	330	365	460
Heat Hose Connector		GL32		① od 35×id 18mm Flat Flange & Metal Clamp		
Drain Tube Size		Glass Tube ; od 24×id 19mm		Glass Tube ; od 34×id 20mm		
Effluent Tray		15Lit	33.9Lit		55.4Lit	120.8Lit
Digital Stirrer	Model, Control	Direct Controlled			Remote Controlled	
		“SHT-50DX”			“SHT-120T”	“SHT-140T”
	Viscosity & Torque	60,000 mPas			90,000 mPas	120,000 mPas
		180 Ncm			210 Ncm	300 Ncm
	Speed	50~1,000 rpm, Resolution ; 1 rpm-Display / 1 rpm-Control				50~800 rpm
	Motor & Rating	Worm-Geared DC Motor (3.5:1 Ratio)				
		In · Output 83 / 75 W, 1 / 10 HP			In · Output 130 / 120 W, 1 / 6 HP	
	Display	Digital LCD with Back-Light Function				
Controller	Digital Feedback Controller with Jog Shuttle Switch (Turn + Push)					
Power	Consumption 85W			Consumption 120W		
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end				
Shaft / Impeller		φ 12/10 PTFE Encapsulated STS			φ 14/12 PTFE Encapsulated STS	
Safety Device		Overload & Heat Protector for Motor				
Dimension Overall		50×50×h150cm	55×55×h170cm		60×60×h180cm	67×67×h200cm
Power Supply		1 Phase, AC 230 V, 50 / 60 Hz				

ORDERING for Reactor Set, without Circulator and Vacu-Pump

<L1>

	I. With Glass Lid, 6 Necks		II. With PTFE Lid, 6 Necks	
Reactor Set	Cat. No.	Price, ₩	Cat. No.	Price, ₩
"REA-10C", 10 / 0.2Lit	SL.Rea55		SL.Rea55P	
"REA-20C", 20 / 0.4Lit	SL.Rea56		SL.Rea56P	
"REA-30C", 30 / 0.5Lit	SL.Rea57		SL.Rea57P	
"REA-50C", 50 / 1Lit	SL.Rea58		SL.Rea58P	
"REA-100C", 100 / 1Lit	SL.Rea59		SL.Rea59P	

► Lid Details ; Visit DAIHAN-Sci.com Shopping-mall !

- Optional Heat/Cooling-Circulators & -Fluid
⇒ See the "Bath Circulators, Accessories" Section !
- Optional Vacuum Pumps ⇒ See the "Pumps, Vacuum" Section !

SciLab™ 360° Rotated 20~100Lit Vacuum Glass Double Wall Reactor Set, Graduated, 45° DN Grooved Flange & FEP O-Ring, 6-Neck Lid With Digital Stirring System · PTFE Drain Valve · Tray · Fabricated Safety Tilted Frame Set, 20~100Lit **틸팅 글라스 진공 2중 반응조 세트, 눈금부, 내열 연결호스 포함**

ReNew


360° Rotated 20~100Lit Double Wall Reactor Set

I. with Glass Lid

II. with PTFE Lid



Vessel을 안전하고 쉽게 360° 양방향 회전 시킴
; 순환호스 장착 길이에 따라 달라짐



I. With Glass Lid
With 6 Joint Necks

II. With All PTFE Lid

* Interchangeable Screw
6 Joints, -200+260°C

* Top Safety ; Non-Breakage,
Chemically Inert

■ Supplied by SciLab Korea

▶ 360° Safety Rotated / Tilted Vessel & Fabricated Sturdy Frame 안전 틸트형 조립식 프레임

■ 360° 2-way Rotation Vessel with "Bevel-Gear" Handle (when without Circulation Hose), Tilt Angles According to the Length of Installed Circulation Hose, and the Vessel Mount Housing Made of Casting Aluminum & PTFE Cushion (Lighter & Strong).

■ Fabricated Sturdy Main Frame ; OD 27.2 / Thick 3mm STS Pipe, Powder Coated Steel Base, and 4 Lockable Casters.

▶ Advantage of Safety Tilting System ; Ideal for Maintenance

■ Safety Cleaning the Vessel

■ Easy Feeding of Large Volume Raw Materials

▶ Graduated Double-wall Vacuum Glass Vessel, 20~100Lit

■ 45° DN Flange & 5mm Thick High Boro-glass 3.3 with Good Physical, Chemical Properties

■ With Direct Contact PTFE Thermo-Probe K-Type, FEP O-Ring for Vacuum, PTFE Drain Valve, and Insulated VITON or Flex STS(Bellow) Hose Set

* Cooling and Heating Solution in the Interlayer can be Completely Removed after Reaction.

▶ High-Torque Digital Agitator System, 50~1000rpm(800rpm / 100Lit)

■ CE Certified ■ With Certificate & Traceability ; Controlled by Serial Number

■ Patented Jog-Shuttle Control System

■ PL(Product Liability) Insurance

■ With Safety Flex-Coupling between Stir-Motor and PTFE Shaft, PTFE Impeller

* Storage Function for the Set Values ; RPM, Timer, and Locking * Chuck - Corrosion Free

Specification of Reactor Set

Models		“REA-20T”	“REA-30T”	“REA-50T”	“REA-100T”
Jacket Vessel Volume		20Lit	30Lit	50Lit	100Lit
Jacket Volume Approx.		5Lit	8Lit	12Lit	25Lit
Vessel OD×h mm		300×630	365×710	410×800	500×970
Inner Vessel OD mm		250	330	365	460
Heat Hose Connector		GL32	① od 35×id 18mm Flat Flange & Metal Clamp		
Drain Glass Tube Size		od 24×id 19mm	od 34×id 20mm		
Effluent Tray		33.9Lit		55.4Lit	120.8Lit
Digital Stirrer	Model, Control	Direct Controlled “SHT-50DX”		Remote Controlled	
	“SHT-120T”			“SHT-140T”	
	Viscosity & Torque	60,000 mPas		90,000 mPas	120,000 mPas
		180 Ncm		210 Ncm	300 Ncm
	Speed	50~1,000 rpm, Resolution ; 1 rpm-Display / 1 rpm-Control			50~800 rpm
	Motor & Rating	Worm-Geared DC Motor (3.5:1 Ratio)			
		In · Output 83 / 75 W, 1 / 10 HP		In · Output 130 / 120 W, 1 / 6 HP	
	Power	Consumption 85 W		Consumption 120 W	
Tilt Angle / Vessel		360° 2-way (when without Circulation Hose)			
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end			
Display		Digital LCD with Back-Light Function			
Controller		Digital Feedback Controller with Jog Shuttle Switch (Turn + Push)			
Shaft / Impeller		φ 12/10 PTFE Encapsulated STS		φ 14/12 PTFE Encapsulated STS	
Safety Device		Overload & Heat Protector for Motor			
Dimension Overall		820×870×h1750		870×950×h1900	950×1050×h2100
Frame Weight		120 kg		140 kg	160 kg
Power Supply		1 Phase, AC 230 V, 50 / 60 Hz			

ORDERING for Reactor Set, without Circulator & Vacu-Pump

<L1>

Reactor Set	I. With Glass Lid, 6 Necks		II. With PTFE Lid, 6 Necks	
	Cat. No.	Price, ₩	Cat. No.	Price, ₩
"REA-20T", 20 / 0.4Lit	SL.Rea62		SL.Rea62P	
"REA-30T", 30 / 0.5Lit	SL.Rea63		SL.Rea63P	
"REA-50T", 50 / 1Lit	SL.Rea64		SL.Rea64P	
"REA-100T", 100 / 1Lit	SL.Rea66		SL.Rea66P	

▶ Lid Details ; Visit DAIHAN-Sci.com Shopping-mall !

▶ Optional Heat/Cooling-Circulators & -Fluid

⇒ See the "Bath Circulators, Accessories" Section !

▶ Optional Vacuum Pumps ⇒ See the "Pumps, Vacuum" Section !


SciLab™ 20~100Lit AutoLifting Glass Double Wall Vacuum Reactor Set, Graduated, 45° DN Grooved Flange & FEP O-Ring, 6-Neck Lid With Digital Stirring System · PTFE Drain Valve · Tray · Tilted Auto Lifting Frame System, 자동리프팅 글라스 진공 2중 반응조 세트, 연결호스 포함

ReNew

**Safety AutoLifting & 30° Manual Tilted Reactor Set ; I. with Glass Lid
II. with PTFE Lid**



Vessel을 옆으로 30° 기울일 수 있음 !



I. With Glass Lid
With 6 Joint Necks

II. With All PTFE Lid
* Interchangeable Screw
6 Joints, -200+260 °C
* Top Safety ; Non-Breakage,
Chemically Inert

■ Supplied by SciLab Korea

► Safety AutoLifting “Up-Down” and 30° Manual Tilted Frame System

■ With Induction Gear Motor, TM-Screw, Timing Belt & Pulley, and Limit Switch

* Adjustable Lifting Speed 1~53rpm & Height Length 200~500mm * Tilt Angle 30 Degree

* Material ; OD ϕ 27.2 × Thick 2(1.5-20Lit)mm STS Pipe and AL-Profile 4040

* Fabricated Construction System ; Fully Adjustable

* Frames According to the Vessel Volumes * With 4 Castors Including 2 Lockable Castors

► Advantage of Safety Tilting System ; Ideal for Maintenance

■ Safety Cleaning the Vessel

■ Easy Feeding of Large Volume Raw Materials

► Graduated Double-wall Vacuum Glass Vessel & 6-Neck Lid

■ 45° DN 200 Flange & 6/7mm Thick High Boro-glass 3.3 with Good Physical, Chemical Properties

■ With (E) Direct Contact PTFE Thermo-Probe K-Type and PTFE Drain Valve

※ Have to Use Heating / Cooling Circulator & -Fluid.

► High-Torque Digital Agitator System, 50~1000rpm(800rpm / 100Lit)

■ CE Certified ■ With Certificate & Traceability ; Controlled by Serial Number

■ Patented Jog-Shuttle Control System ■ PL (Product Liability) Insurance

■ With Safety Flex-Coupling between Stir-Motor and PTFE Shaft, PTFE Impeller

* Storage Function for the set values ; RPM, Timer, and Locking * Chuck - Corrosion Free

※ With Effluent Collection Safety Tray STS ; Used in Case of Vessel Failure

유출물 포집 트레이 ; 반응용기 파손 등 고장시 First Aid 용 !

Specification of Reactor Set, DN200 Flange & 6-Neck Lid, 20~100Lit.

Models		“REA-20A”	“REA-30A”	“REA-50A”	“REA-100A”
Jacket Vessel Volume		20 Lit	30 Lit	50 Lit	100 Lit
Jacket Volume Approx.		5 Lit	8 Lit	12 Lit	25 Lit
Vessel OD×h mm		300×630	365×710	410×800	500×970
Inner Vessel OD mm		250	330	365	460
Heat Hose Connector		GL32	① od 35×id 18mm Flat Flange & Metal Clamp		
Drain Tube Glass		od 24×id 19mm	Glass Tube ; od 34×id 20mm		
Effluent Tray		33.9 Lit		55.4 Lit	126 Lit
Digital Stirrer	Model, Control	Direct Controlled “SHT-50DX”		Remote Controlled	
	“SHT-120T”			“SHT-140T”	
	Viscosity & Torque	60,000 mPas		90,000 mPas	120,000 mPas
		180 Ncm		210 Ncm	300 Ncm
	Speed	50~1,000 rpm, Resolution ; 1 rpm-Display / 1 rpm-Control			50~800 rpm
Motor & Rating	Worm-Geared DC Motor (3.5:1 Ratio)				
	In · Output 83 / 75 W, 1 / 10 HP		In · Output 130 / 120 W, 1 / 6 HP		
Lifting	Speed	1~53 rpm, 15 Step Control			
	Height	200~500mm Distance Control			
Tilt Angle		Manual 30 Degree			
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end			
Display		Digital LCD with Back-Light Function			
Controller		Digital Feedback Controller with Jog Shuttle Switch (Turn + Push)			
Shaft / Impeller		φ 12/10 PTFE Encapsulated STS		φ 14/12 PTFE Encapsulated STS	
Safety Device		Overload & Heat Protector for Motor			
Power Consumption		175 W		300 W	
Dimension Overall		915×670×h2100		1000×720×h2400	1200×790×h2445
Power Supply		1 Phase, AC 230 V, 50 / 60 Hz			

ORDERING for Reactor Set, without Circulator & Vacu-Pump

<L1>

Reactor Set	I. With Glass Lid, 6 Necks		II. With PTFE Lid, 6 Necks	
	Cat. No.	Price, ₩	Cat. No.	Price, ₩
“REA-20A”, 20 / 0.4Lit	SL.Rea72		SL.Rea72P	
“REA-30A”, 30 / 0.5Lit	SL.Rea73		SL.Rea73P	
“REA-50A”, 50 / 1Lit	SL.Rea74		SL.Rea74P	
“REA-100A”, 100 / 1Lit	SL.Rea76		SL.Rea76P	

► Lid Details ; Visit DAIHAN-Sci.com Shopping-mall !

► Optional Heat/Cooling-Circulators & -Fluid

⇒ See the “Bath Circulators, Accessories” Section !

► Optional Vacuum Pumps ⇒ See the “Pumps, Vacuum” Section !

SciLab™ 1~10Lit Glass Lidded STS Double Wall Self-Standing Vacuum Reactor Sets, 45° DN Grooved Flange & FEP O-Ring
With All PTFE GL Screw Connect Kit · Digital Stirring System · Fabricated Sturdy Stand Frame, 1~10Lit 이중 스텐 진공베셀 반응조 세트(유리뚜껑), “자립형”, 연결호스 포함

Advantage of 45°DN Flange ;

1. Anti-Breakage
2. Excellent Vacuum Operation
3. Easy & Quick Installation / Removal



1~10Lit STS Double Wall Reactor Set with Glass Lid, “Self-Standing”



■ Supplied by SciLab Korea

► Self-Standing Stainless-steel Double Wall Vacuum Vessels & Glass Lids

■ Stainless-steel(STS) 304 Self-Standing Double Wall Vessel ; Thick 1.5mm

45°DN Grooved Flange & Teflon FEP O-Ring(250℃)

■ Glass Lids ; Boro-glass 3.3, Approx Thick 6mm

45°DN Flange, Fit to 3-Segment STS Quick-Release Clamp

- * Precise Temp Control ; Direct Contact PTFE Coated Thermo-Probe K-Type
- * -200℃ +260℃ All PTFE Screw GL25 & GL32 Connect-Port for the Circulation Hose of Heat Transfer Fluid of Heating / Cooling Circulators
- * With Insulated VITON Circulation Hose Set.

► High Torque Digital Stirring System ; 180Ncm, 60000mPas, 50~1000rpm

■ CE Certified ■ With Certificate & Traceability ; Controlled by Serial Number

■ Patented Jog-Shuttle Control System

■ PL(Product Liability) Insurance

■ With Safety Flex-Coupling between Stir-Motor and PTFE Shaft, PTFE Impeller.

- * The Stirrer Model ; “SHT-50DX” Direct Controlled(SL.SOS01051), 85W, 230V, 50/60Hz(1 Phase)
- * Digital LCD Display Back Light 50~1000/1rpm, Timer & Alarm 99hr 59min & Error
- * Smooth Start and Change of Stirring Speed * Chuck - Corrosion Free
- * Storage Function for the Set Values ; RPM and Timer
- * Locking Mode for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Safety Device - Overload & Heat Protector for Motor

► Safety Work Fabricated Stand Frame System

- * Heavy Duty Casting Iron Base(40~69cm Adjustable Feet Distance) with 3×OD ϕ23mm STS Pipe(T. 2mm) Rod-Height 1200 & 1500mm and 2×Silicone Rubber Padded STS Band Clamps
- * Stand Frame Models ; SL.ReF7101(Height 1200 Type) for 1~5Lit
SL.ReF7102(Height 1500 Type) for 10Lit

A. Glass Lids, 45° DN Flanged, with 5 & 6-Socket Joints

Lids	45° DN Flange	Necks	Center-Joint	Side-Joint
DN100	id 100×od 138mm	5	24/40	3×15° 24/40, 1×Vertical 24/40
DN120	id 120×od 158mm	5	34/45	3×15° 24/40, 1×Vertical 24/40
DN150	id 150×od 184mm	5	34/45	2×15° 24/40, 1×15° 34/35, 1×Vertical 24/40
DN200	id 200×od 242mm	6	34/45	1×15° 34/45, 1×15° 24/40, 1×15° 45/40, 2×Vertical 24/40

* Side 34/35 & 45/40 Neck ; Good For Inlet-Outlet of Substances.

B. Self-Standing Stainless-steel Double Wall Vessels, Flat Bottom, 1~10Lit

■ With 45°DN Grooved Flange & Teflon FEP O-Ring for Vacuum, DN100~DN200

Vessel Capacity	Flange 45° DN	Dimension (mm)	FEP Vacu. O-Ring (mm)	Connect Port All PTFE
1Lit	DN100	OD 135(inner od 102)×h205	id ϕ 110×thick 4	Screw GL25
3Lit	DN120	OD 157(inner od 127)×h285	id ϕ 133×thick 4	
5Lit	DN150	OD 183(inner od 153)×h310	id ϕ 157×thick 5	
10Lit	DN200	OD 257(inner od 217)×h350	id ϕ 215×thick 5	Screw GL32

ORDERING for Glass Lidded STS Double Wall Self-Standing Vacuum Reactor Sets <L1>

Reactor Set Cat. No.	STS Vessel Capa. & 45° DN Flange	Connect Port All PTFE (260℃)	Digital Stirring High Torque	Price, ₩
SL.Rea81S	1Lit, DN100	Screw GL25	180Ncm 60,000mPas 50~1000rpm 230V, 50/60Hz	
SL.Rea82S	3Lit, DN120			
SL.Rea83S	5Lit, DN150			
SL.Rea84S	10Lit, DN200	Screw GL32		

► The Details of All Components

⇒ Visit DAIHAN-Sci.com Shopping-mall !

► Optional Heat/Cooling-Circulators & -Fluid

⇒ See the “Bath Circulators, Accessories” Section !

► Optional Vacuum Pumps ⇒ See the “Pumps, Vacuum” Section !

SciLab™ 3~100Lit Glass Lidded / Drain Valved Stainless-steel(STS) Double Wall Vacuum Reactor Sets

With Safety Construction & Digital High-Torque Stirring System, 3~100Lit 이중 스텐 진공베셀 반응조 세트(유리뚜껑), “배출밸브형”, 연결호스 포함

I. 3~10Lit Reactor Set with 45°DN Grooved Flange

ReNew

■ Supplied by SciLab Korea

► Safety Construction & Digital Precise Control System

- * Stainless-steel(STS) 304 Double Wall Vessels, Boro-glass 3.3 Lids, and Condensers
- * Precise Temp Control ; Direct Contact PTFE Coated Thermo-Probe K-Type
- * Certificated Digital Stirring System ; **Patented Jog-Shuttle Control**, Serial No-Traceability, **Product Liability Insurance**
- * With Safety Flex-Coupling between Stir-Motor and PTFE Encapsulated STS Shaft / PTFE Impeller

► Digital High Torque Stirring System ; 50~1000rpm(800rpm-100Lit), 230V, 50/60Hz(1 Phase)

- ① 3~30Lit ; 180Ncm / 60,000mPas Viscosity, ② 50Lit ; 210Ncm / 90,000mPas Viscosity,
- ③ 100Lit ; 300Ncm / 120,000mPas Viscosity, ④ 85Watt (3~30Lit), 120Watt (50 & 100Lit)
- ⑤ Timer & Alarm ; 99hr 59min(with Continuous Run), Error Status & Timer End, and
- ⑥ Safety Device ; Overload & Heat Protector for Motor, and ⑦ Stir Shaft ; PTFE Encapsulated STS

I. 3 · 5 · 10Lit 45°DN Flanged / PTFE Drain Valved STS Double Wall Reactor Sets With Glass Lid, Teflon FEP O-Ring and STS Quick-Release Clamp, 45°DN Flange Type

- * Stainless-steel(STS)304 Double Wall Vessels, Thick 1.5mm
- * STS 304 Pipe Frame 450×450×h1500mm Set, with All Clamping Devices

<L1>

Models	“REA92S”	“REA93S”	“REA94S”
Working Volume	3Lit	5Lit	10Lit
Jacket Volume Approx.	0.8Lit	1.5Lit	3Lit
DN Flange	DN 120	DN 150	DN 200
Vessel OD × h mm	157 × 320	183 × 310	243 × 350
Inner Vessel OD mm	127	153	217
Heat Hose Connector	All PTFE Screw GL25		All PTFE Screw GL32
Condenser	Double Wall Reflux φ 42 × h500mm, 24/40		2 Coiled Layer φ 85 × h560mm, 34/45
Stir Shaft, PTFE / STS	φ 12/10		
Dimension Overall	500 × 500 × h1500mm		
Glass Lids	Ports	5 Joints ; Center 34/45	6 Joints ; Center 34/45
	Side Joint	3 × 15° 24/40, 1 × Vertical 24/40	2 × 15° 24/40, 1 × 15° 34/35, 1 × Vertical 24/40
Cat. No.	SL.Rea92S	SL.Rea93S	SL.Rea94S
Price, Each. ₩			

II. 20~100Lit Reactor Set with Flat Flange & Link Clamping System



II. 20~100Lit Flat Flanged / STS Drain Ball Valved STS Double Wall Reactor Sets With 7 Port Glass Lid, PTFE Seal Gasket, and Link Clamping System, Flat Flange Type

- With an illumination Light on the Lid for Viewing the internal Vessel, 베셀 내부 투시용 등
- * Stainless-steel(STS)304 Double Wall Vessels Thick 3.4mm
- * Heavy Duty STS 304 Frame Set ; Fully Welded Rectangular Pipe, with All Clamping

Models	“REA95S”	“REA96S”	“REA97S”	“REA98S”
Working Volume	20Lit	30Lit	50Lit	100Lit
Jacket Volume Approx.	5Lit	8Lit	12Lit	35Lit
Flange id × od mm	350 × 396		400 × 496	500 × 596
Vessel OD × h mm	400 × 332	400 × 432	450 × 525	600 × 650
Inner Vessel OD mm	350		400	500
Heat Hose Connector	1/2" Male Stainless-steel Clamp			
Condenser	2 Coiled Layer φ 85 × h560mm, 34/45		2 Coiled Layer φ 100 × h595mm, 34/45	3 Coiled Layer φ 165 × h645mm, 45/50
Stir Shaft, PTFE / STS	φ 14/12		φ 16/14	
Dimension Overall	520 × 520 × h1830mm	520 × 520 × h1930mm	590 × 590 × h2020mm	740 × 740 × h2150mm
Glass Lids	Ports	6 Joints & a Light (배터리형 내부조명등 장착부) ; Center 34/45		
	Side Joint	1 × 15° 34/45, 1 × 15° 24/40, 1 × 15° 45/40, 2 × Vertical 24/40		
Cat. No.	SL.Rea95S	SL.Rea96S	SL.Rea97S	SL.Rea98S
Price, Each. ₩				

- Optional Heat/Cooling-Circulators & -Fluid
⇒ See the “Bath Circulators, Accessories” Section !
- Optional Vacuum Pumps ⇒ See the “Pumps, Vacuum” Section !

DAIHAN™ High Temp Bath Circulator 250°C & Refrigerated/Heating Bath Circulator -35°C +150°C
정밀 고온 배스겸용 써클레이터와 냉동/가열 배스 겸용 써클레이터

CE up to 250°C



3-Year Warranty

1. High Temp Bath Circulators

CE -35°C & -25°C ~+150°C



3-Year Warranty

2 & 3. Refrigerated/Heating Bath Circulators

※ See the "Bath Circulators" Section for the Details 배스 · 써클레이터 겸용.

1. 250°C High-temp Bath Circulator, 8 · 12 · 22 · 30 Lit
2. -25°C +150°C Refrigerated/Heating Bath Circulator, 8 · 12 · 22 · 30 Lit
3. -35°C +150°C Refrigerated/Heating Bath Circulator, 8 · 12 · 22 · 30 Lit



SK™ -40°C +400°C Best COOL-/HEAT-ing Transfer Fluid, Synthetic Silicone Mixture, 20Lit
For Low/High-Temp Baths, Clearly Harmless Liquid, 냉매/열매 겸용 실리콘계 합성 열매체유

■ Supplied by SciLab Korea

■ The Best Circulation Properties / Thermal Stability.

* Colorless and Transparent Liquid.

* Do not use for Silicone Tube

(1) SYN THERM 7008®, -40°C +400°C, 저/고온용

* Kinematic Viscosity, cSt ⇒ (1) -40°C ; 51.56cSt, (2) 0°C ; 16.08cSt, (3) 100°C ; 3.43cSt, (4) 200°C ; 1.36 cSt <L1>
* Specific Gravity ; 0.9435 at 15°C * Flash Point(PMCC) ; 168°C * Pour Point ; -72°C

SL.Syn7008 Cool/Heat Transfer Fluid, -40+400°C, 20Lit

Pk.	Price, ₩
20Lit	

(2) SUPER THERM FG-20®, -10°C +300°C, 고온용

* Food Grade Thermal Medium Oil Satisfying FDA, NSF Standards
* Kinematic Viscosity, cSt ⇒ (1) 25°C ; 35.65cSt, (2) 100°C ; 4.23cSt, (3) 200°C ; 1.27 cSt
* Specific Gravity ; 0.8374 at 15°C * Flash Point ; 230°C * Pour Point ; -15°C

SL.Fg0020 Cool/Heat Transfer Fluid, -10+300°C, 20Lit

Pk.	Price, ₩
20Lit	



SHL™ SYN THERM 22L® -5°C +320°C Heating Transfer Fluid
Synthetic Alkylbenzene Derivative, for Low- & High-Temp Circulation, 알킬벤젠계 합성열매

■ Supplied by SciLab Korea

* Hazards Identification of MSDS(Material Safety Data Sheet) ;

- H315 Causes Skin Irritation

- H319 Causes Serious Eye Irritation

- H335 May Cause Respiratory Irritation

- H336 May Cause Drowsiness or Dizziness

* Kinematic Viscosity & Density ;

(1) at 0°C ; 22.7cSt, 0.885kg/l (2) at 40°C ; 21.80cSt, 0.859kg/l (3) at 80°C ; 6.09cSt, 0.834kg/l
(4) at 140°C ; 6.09cSt, 0.795kg/l (5) at 200°C ; 1.10cSt, 0.756kg/l (6) at 260°C ; 0.74cSt, 0.718kg/l <L1>

SL.SYN22L SYN THERM 22L® Heating Transfer Fluid, -5°C +320°C, 20Lit

Pk.	Price, ₩
20Lit	

DAIHAN™ -30℃+200℃ **Premium** Precise Refrigerated Heating External Circulator “CRH”, Closed-type, Fill-10·18·25 Lit Ideal for Evaporator/Reactor &c. Heating & Cooling Line, with Pre-Cooling Sys, Used with External Direct Contact K-type Temp. Probe(Optional) With High Quality Magnetic Pump, Flow 20·40·60 Lit/min, Lift 38m, ±0.2℃, 정밀 고온/냉동 겸용 외부 써클레이터-밀폐형

Advantage of Closed-type Circulators ; 밀폐형 써클레이터의 장점

- (1) Small Power Consumption ; Less than Bath Circulators
- (2) Eco-Friendly ; Non-Evaporation of Heat Transfer Fluid and Non-Smell of Fluid
- (3) Uniformity of Temperature Control



3-Year Warranty



Fig. “CRH” in the Reactor System.

※ Note ; Cooling Transfer Fluid Should be an Antifreezing Solution

- Supplied by DAIHAN Korea
- PL (Product Liability) Insurance
- Max. Lift ; 38m
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Patented Jog-Dial Control System and Digital LCD with Back-Light Function.
- * Innovated Powerful Circulation Pump ensures Temp. Uniformity ; Internal and External Circulation (Max Flow Rate ; 20·40·60 Lit/min)
- * Assembly-type Safety Connection Ports/Pipe (for Fluid Drain/In/Out-let ; ISO RC 1/2" Screw Stainless-steel.)
- * With High Quality Shock-proof Caster for Convenient and Safe Moving.
- * Embedded External Nozzle Part for the Prevention of Nozzle Damaged.
- * Stainless-steel Internal Bath for Prevention of Corrosion.
- * Independently 3-systems ; Refrigerating System, Heating System, Pre-Cooling System.
- * Urethane Foam Insulation * More Faster Cool-down and Heat-up Time
- * Warning and Auto-Stop Function in case of Liquid Shortage.
- * Connection with Evaporator or Viscometer is applicable.
- * 외부 순환을 위한 추가 커넥터가 필요합니다. (ex. PT 1/2" Hose nipple Connector)



Connection Ports

Specifications & Ordering information

<L1>

Spec.	Model	“CRH-20” Useful for 5~20Lit Reactor	“CRH-50” Useful for 20~50Lit Reactor	“CRH-100” Useful for 50~100Lit Reactor
Temp. Range & Stability		-30℃+200℃-Range, ±0.2℃(at 15℃)-Stability		
Heat-Fluid Filling Volume		10 Lit max	18 Lit max	25 Lit max
Useful for 5~100Lit Reactors		for 5~20Lit Reactor	for 20~50Lit Reactor	for 50~100Lit Reactor
Refrigerator Capacity		4.5 kW	7.5 kW	18 kW
Heater		2.5 kW	6.0 kW	12 kW
Material		Internal ; Stainless-steel, External ; Powder Coated steel		
Lid & Top Opening Dia.		φ80 mm		
Fluid Drain/In/Outlet Port		ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel		
Overflow Port		φ12mm One-Touch Fitting		
Pump Power, Max. Flow Rate		20 Lit/min	40 Lit/min	60 Lit/min
Lift Height, m	Max. Pressure & Lift	Max Pressure ; 4 bar, Max. Lift ; 38m		
Cooling Capacity	at +15℃	8.1 kW	13.1 kW	22.7 kW
	at 0℃	4.3 kW	7.3 kW	12.6 kW
	at -25℃	1.4 kW	3.1 kW	5.45 kW
Optional External Temp. Sensor K-type Probe		Can be used with External Direct Contact K-type Temp. Probe. (Optional)		
External Dimension		600×850×h1200mm	650×900×h1350mm	750×900×h1450mm
Electrical Reqts.		3Phase 380V		
Cat. No.	60Hz	DH.WHR0020	DH.WHR0050	DH.WHR0100
	50Hz	DH.WHR0025	DH.WHR0055	DH.WHR0150
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도




※ Insulated Circulator Hoses ; See the “Bath Circulators, Accessories” Section !

DAIHAN™ -20℃+200℃ Refrigerated Heating External Circulator “HCR”, Closed-type, Fill-10·20·30·50Lit
With Insulated Stainless-steel Hose, Ideal for Evaporator/Reactor &c. Heating & Cooling Line, Flow 20·25·30Lit/min, Lift Max. 8 & 12m,
고온/냉동 겸용 외부 써클레이터-밀폐형, 단열 튜빙 포함

**Advantage of Closed-type Circulators ;
밀폐형 써클레이터의 장점**

- (1) Small Power Consumption ;
Less than Bath Circulators
- (2) Eco-Friendly ; Non-Evaporation of Heat Transfer
Fluid and Non-Smell of Fluid
- (3) Uniformity of Temperature Control



Connection Ports



Insulated Flexible STS Hose Included

- Supplied by DAIHAN Korea ■ CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information.
- Precise Temperature Control System, $\pm 1^\circ\text{C}$
- Microprocessor Based Control ensures Stable and Consistent Thermal Performance.
- Ideal for Chemical, Pharmaceutical, and Food Research from Lab to Pilot Scale
- Includes Two Insulated Flexible Stainless-steel Hoses (0.7m & 1m)
- * Integrated Heating and Cooling System with Rapid Temperature Stabilization
- * Optimized for Repeated Thermal Cycling * Sealed Structure minimizes Leakage Risk.
- * Internal Stainless-steel Circulation Path ensures Corrosion Resistance and Durability.
- * External Nozzle Protection Design prevents Physical Damage.
- * Equipped with Caster Wheels for Safe and Easy Mobility
- * Connection with Evaporator or Viscometer is Applicable.

Specifications & Ordering information

<L1>

Spec.	Model	“HCR10-20”	“HCR20-20”	“HCR30-20”	“HCR50-20”
Temp. Range & Accuracy	-20℃+200℃-Range, $\pm 1^\circ\text{C}$ at Set Temperature				
Heat-Fluid Filling Volume		Max. 10Lit	Max. 20Lit	Max. 30Lit	Max. 50Lit
Useful for 5~100Lit Reactors		for 5~20Lit Reactor	for 20~30Lit Reactor	for 50Lit Reactor	for 100Lit Reactor
Refrigerator Capacity		0.75kW	1.09kW	1.5kW	3.75kW
Heater		2kW	3kW	4.5kW	6kW
Material	Internal ; Stainless-steel, External ; Powder Coated steel				
Lid & Top Opening Dia.	$\phi 45\text{mm}$				
Fluid In/Outlet Port	ISO RC 1/2" (English System AG 1/2") Thread, Stainless-steel				
Circulation Hose	0.7m & 1m, Included				
Pump Power, Lift Height, m	Max. Flow Rate	20Lit/min	25Lit/min	30Lit/min	
	Max. Lift	8m		12m	
Cooling Capacity	at +15℃	2.1kW	2.8kW	5.6kW	10.5kW
	at 0℃	1.5kW	1.9kW	5.5kW	8.5kW
	at -30℃	0.5kW	0.6kW	0.9kW	0.9kW
External Dimension (mm)	620×540×h920			710×580×h1050	770×670×h1180
Electrical Rqts.	220V, 50/60Hz				3 Phase 380V
Cat. No.	DH.HCR10.20	DH.HCR20.20	DH.HCR30.20	DH.HCR50.20	
Price, Set. ₩					

※ 운임 별도

Spare Insulated Flexible Stainless-steel Hose

SL.Tub7033 Insulated Flexible STS Hose Only, 1/2", id17×od26mm, Overall od60mm, L0.7m
SL.Tub7034 , L1m

Ea.	Price, ₩

※ For More Information about Circulator Hoses ;
See the “Bath Circulators, Accessories” Section !

DAIHAN Scientific

REFRIGERATORS & ULT FREEZERS

Freezers, Bio -35 °C	P. 281
Freezers, Ultra Low Temp	P. 273~280
-95 °C, Upright-type ReNew	P. 275
-86 °C, Personal-type	P. 277
-86 °C, Upright-type ReNew	P. 276
-45 °C, Upright-type	P. 278
Accessories ; Racks, Freezer Boxes, Recoders, etc.	P. 279~280
Refrigerators	P. 282~285
0~10 °C Laboratory Refrigerators	P. 282
2~10 °C Blood Bank Refrigerators	P. 283
2~10 °C Pharm Refrigerators	P. 284~285



DAIHAN™ Energy Saving SMART Ultra Low Temp(ULT) Freezer, -95°C / -86°C **ReNew**

With Innovative INVERTER COMPRESSOR SYSTEM™ 절전형 혁신 초저온 스마트 냉동고

Innovative INVERTER COMPRESSOR SYSTEM™ **ReNew**

실제 냉각 수요에 맞춰 동작 속도를 조절하는 VSD(Variable Speed Drive) 기반 혁신 인버터 컴프레서 시스템을 개발·적용하여 전력 소모를 최소화 !

DAIHAN™ Ultra Low Temperature Freezers Feature DAIHAN's Innovative INVERTER COMPRESSOR SYSTEM™, Ensuring Higher Efficiency and Energy-saving Operation. Unlike Conventional Fixed-speed Compressors that Always Run at One Constant Speed, Innovative VSD(Variable Speed Drive) Compressor Technology Ensures Optimal Performance by Intelligently Adjusting Compressor Speed to Actual Cooling Demand. By Minimizing Unnecessary Energy Use, Noise, and Vibration, New Models Ensure Longer Compressor Life, Lower Operating Costs, and Reliable Performance for Today's Laboratories. Furthermore, the Integration of DAIHAN's Innovative Control System & Refrigeration Technology with the Inverter Compressor Makes this Outstanding Energy Efficiency Possible, Establishing a New Benchmark for Sustainable Ultra Low Temperature Freezer.



Benefit 1. Energy Saving - Approx 20%

정속형 컴프레서 대비 약 20%의 절전(에너지 절감) 효과

DAIHAN's Innovative INVERTER COMPRESSOR SYSTEM™ ensures Remarkable Energy Cost Savings Compared with Conventional Fixed-Speed Compressors, Achieved with Innovative Control System & Refrigeration Technology that Reduces Unnecessary Power Consumption. This Reduces Daily Power Consumption by Approx. 20% (2.2 kWh/day for the "Fre400-E95" Model).

- ▶ **Power Consumption Reduction Rate** ; Approx. 20%
- ▶ **Power Consumption** ; 11.1kwh/Day $\Rightarrow$ 8.9 kWh/day (based on "Fre400-E95" Model)
- ▶ **Annual Electricity Cost Savings** ; Approx. KRW250,000/year (based on "Fre400-E95" Model)

Benefit 2. Protection of the Environment

에너지 절약을 통해 전력 생산 시 발생하는 CO₂ 배출 감소 효과

DAIHAN's Innovative INVERTER COMPRESSOR SYSTEM™ achieves Significant CO₂ Emission Reductions Compared with Conventional Fixed-Speed Compressor Systems by Lowering Electricity Consumption, thereby Reducing CO₂ Emissions associated with Electricity Generation.

- ▶ **CO₂ Emission Reduction** ; Approx. 390 kg/year (based on "Fre400-E95" Model)
This Calculation uses the International Energy Agency (IEA) Average Emission Factor of approx. 0.45 kg CO₂/kWh of Electricity Savings. This is Equivalent to the CO₂ Emissions of a Gasoline Vehicle Driving approx. 2,660 km (2.31 kg CO₂/L Gasoline).
- ▶ **Tree Planting Effect** ; Approx. 65~70 Trees/year (based on "Fre400-E95" Model)
Equivalent to the Annual CO₂ Absorption of about 70 Mature Trees, Each Absorbing Apex. 5.5~6 kg CO₂ per Year.

CO₂ REDUCTION



DAIHAN™ Energy Saving ULT Freezer Models with Innovative INVERTER COMPRESSOR SYSTEM™

Model \ Capa.	393Lit.	503Lit.	714Lit.	812Lit.
-95°C Model	Fre400-E95	Fre500-E95	Fre700-E90	Fre800-E90
-86°C Model	Fre400-E86	Fre500-E86	Fre700-E80	Fre800-E80


DAIHAN™ Innovation of Energy Saving ULT Freezing System ReNew

Double Security DuoFreez™ & High Efficiency UniFreez™


I. DuoFreez™ $-95^{\circ}\text{C} \sim -65^{\circ}\text{C}$ Independent Dual Freezing System for Double Security

We offer a Unique High-Stability System Powered by an Independent Dual Freezing System that Operates Independently. Normally Both Compressors are Operational and the Dual Freezing System Maintains Temperature in the Inner Storage Chamber Down to -95°C . In the Event of an Unexpected Failure of One of the Innovative INVERTER COMPRESSOR, the Other Independent Freezing System continues to Operate as a 'Backup' and Maintains Ultra-Low Temperatures of -80°C . This Innovative Technology Provides Customers with the Most Important ULT Freezer Features ; Maximum Reliability and Robustness.

■ DuoFreez™ -95°C Independent Dual Freezing System Energy Saving - Approx 20%




■ Competitor's Typical Cascaded Dual Freezing System



II. UniFreez™ $-86^{\circ}\text{C} \sim -65^{\circ}\text{C}$ Single Freezing System for High Efficiency


In a General Way, as Shown in the Diagram Below, Most of the Ultra-Low Temperature Freezers have Two Compressors Working Together to Lower the Internal Temperature to -86°C , and the Overall Freezing Mechanism is Quite Complex. However, Our Advanced Single Freezing System uses Only a Single Innovative INVERTER COMPRESSOR. There is Only a Single Set of Compressor, Capillary Tube, Condenser, and Evaporator, But the Advanced Compressor Control Technology and Special Refrigerant Mixing Process Ensure that the Temperature is Lowered to -86°C . The Single Freezing System Overcomes the Technical Complexity of the Typical Cascaded Dual Compressor System. The Simplification of the Freezing Mechanism and Electricity Minimizes Noise, Power Consumption, and Operating Errors.

■ UniFreez™ -86°C Single Freezing System Energy Saving - Approx 20%



- Innovated Simple Structure
- Low Noise, Low Power Consumption and Simple Structure for Easy Maintenance

■ Competitor's Typical Cascaded Dual Freezing System



- Unmanageable Complex Structure

DAIHAN™ -95℃ ~ -65℃ SMART Multiuse ULT Freezer-DuoFreez™ Double Security System, CFC-Free, 393 · 503 · 714 · 812Lit
I. Energy Saving Models (Innovative INVERTER COMPRESSOR SYSTEM™) & **II. Medicaluse Models**, with **Smart WiRe™** App Monitoring System
 Ideal for Secure Storage of Vaccine · Viruses · Anatomy Components · Bio-Products · DNA · Blood · Tissue Samples etc. in Lab & Industry
 스마트 다용도 초저온 냉동고-절전형 & 의료용, 백신 · 바이러스 · 해부조직물 · Bio물질 · DNA · 혈액 · 티슈시료 등을 안전 보관 유지용

I. Energy Saving Models 절전형



II. Medicaluse 의료용(리코더 내장)



▶ You Can Choose 4 or 5 Inner Doors Model !

▶ Note !

Optional Freezer Racks, Boxes, LCO₂ Back Up Sys.
 and Remote Alarm Sys.
 ⇒ See the "Refrigerators, Freezer Access." Section.

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- With Certificate · Traceability · Delivery-Information ; Controlled by Serial Number
- DuoFreez™ Double Security ; Independent Dual Freezing System Maintains Temperature Down to -95℃. If One System Fails, The Other Maintains the Temperature at -80℃.
- ▶ **I. Energy Saving Models ; DAIHAN's Innovative INVERTER COMPRESSOR SYSTEM™ ;**
 Optimized VSD(Variable Speed Drive) Compressor Tech. Ensures Energy Saving.
 - Power Consumption Reduction Rate ;
 Approx. 20%(Compared with Conventional Fixed-speed Compressor)
 - Reducing CO₂ Emissions for Environmental Protection
 - Low Noise / Vibration, Ensuring Longer Compressor Life.
- ▶ **II. Medicaluse Models ; Built in Temp Chart Recorder**
- **Smart-Lab™ Control ; Remote Controlled by Mobile-Phone with WiRe™ App Service ;**
 - Internet Connectivity with WiFi. - Storage Management & Automatic Temp Data Recording.
 - Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.
 - Self-Diagnostic Function ; Built in Standard Warning and Alarm Function
- * High-Quality Insulation Panel and Inner Doors Ensure Defrosting Around Outer Surface
- * Automatic Vacuum Breaker ; Easy to Re-Open the Door
- * Robust & Flexible Design of Door Handle * External Temp Sensor Hole
- * Filter-Free Mechanism Lowers Heat Generated by the Compressor

Specifications & Ordering information

※ 운임 별도 <M1>

I. Energy Saving Model ; Inverter Comp. Sys.		“Fre400-E95”		“Fre500-E95”		“Fre700-E90”		“Fre800-E90”	
		Without Temp Chart Recorder							
Compressor		Two(2) INVERTER COMPRESSOR ; Double Security Variable Speed Inverter Type Compressor							
Power Consumption		8.9kWh/day		9.1kWh/day		10.8kWh/day		10.9kWh/day	
II. Medicaluse Approved ; Built-in Recorder		“Fre400-B95”		“Fre500-B95”		“Fre700-B90”		“Fre800-B90”	
		Built in Temp Chart Recorder “KT622” Record Range ; −100℃ ~ 0℃, divi 2℃ ; 7Days/1Chart, 60Charts/Cs.							
Compressor		Two(2) Standard Compressors ; Double Security Fixed Speed ON/OFF Type Compressor							
Power Consumption		11.1kWh/day		11.7kWh/day		14.2kWh/day		14.8kWh/day	
Capacity		393 Lit		503 Lit		714 Lit		812 Lit	
Temp Range		−95℃ ~ −65℃				−90℃ ~ −65℃			
		at R/T of 25℃ and Humidity of 70%							
Temp Resolution & Sensor		1℃ - Display, 1℃ - Control & PT100 Sensor							
Condenser-type		Block Condenser, No Pin, Filter Free Mechanism							
Refrigerant		Mixed Gas Refrigerant(DH.MGR) ; CFC Free							
Controller		Smart-Lab™ ; 7" Full Touch Screen TFT LCD							
Shelf(Included)		3ea(4-Inner Door), 4ea(5-Inner Door) Perforated Stainless-steel Shelves							
Max. Q'ty of Racks	4-Inner Door (Height 31.2 cm)	16pcs of 15(5×3)- Placed Rack (RFU4234)	16pcs of 20(5×4)- Placed Rack (RFU4257)	24pcs of 20(5×4)- Placed Rack (RFU4257)	28pcs of 20(5×4)- Placed Rack (RFU4257)				
	5-Inner Door (Height 25.2 cm)	20pcs of 12(4×3)- Placed Rack (RFU5234)	20pcs of 16(4×4)- Placed Rack (RFU5257)	30pcs of 16(4×4)- Placed Rack (RFU5257)	35pcs of 16(4×4)- Placed Rack (RFU5257)				
Material		Internal ; Stainless-steel(#304), External ; Powder Coated Steel Inner Door ; Powder Coated Stainless-steel(#304)							
Dimension (w×d×h)	Internal (cm)	60×50×131	60×64×131	88×62×131	100×62×131				
	External (cm)	84×75×190	84×89×190	112×87×190	124×87×190				
Weight		260 kg	300 kg	320 kg	380 kg				
Packing Size(w×d×h) (cm)		116×116×222	116×124×222	140×140×221	152×140×221				
& Gross Weight		352 kg	392 kg	442 kg	502 kg				
Electric Reqts.		1 Phase, AC120V, 60Hz or AC 230V, 50/60Hz							

I. Energy Saving Model ; Innovative INVERTER COMPRESSOR SYSTEM™, Without Recorder					
Cat. No.	230V	DH.FreE9400	DH.FreE9500	DH.FreE9700	DH.FreE9800
	(4-Inner Door)	DH.FreE9401	DH.FreE9501	DH.FreE9701	DH.FreE9801
Price, Each. ₩					
Cat. No.	230V	DH.FreE94005	DH.FreE95005	DH.FreE97005	DH.FreE98005
	(5-Inner Door)	DH.FreE94015	DH.FreE95015	DH.FreE97015	DH.FreE98015
Price, Each. ₩					

II. Medicaluse Approved Model ; Built-in Temp Chart Recorder					
Cat. No.	230V	DH.FreB9400	DH.FreB9500	DH.FreB9700	DH.FreB9800
	(4-Inner Door)	DH.FreB9401	DH.FreB9501	DH.FreB9701	DH.FreB9801
Price, Each. ₩					

DAIHAN™ -86℃ ~ -65℃ SMART Multiuse ULT Freezer-**UniFreez™** High Efficient Freezing System, CFC-Free, 393 · 503 · 714 · 812Lit
I. Energy Saving Models (Innovative INVERTER COMPRESSOR SYSTEM™) & **II. Medicaluse Models**, with **Smart WiRe™** App Monitoring System
 Ideal for Secure Storage of Vaccine · Viruses · Anatomy Components · Bio-Products · DNA · Blood · Tissue Samples etc. in Lab & Industry
 스마트 다용도 초저온 냉동고-절전형 & 의료용, 백신, 바이러스, 해부조직물, Bio물질, DNA, 혈액, 티슈시료 등을 안전 보관 유지용

I. Energy Saving Models 절전형



5-Year Warranty ; Compressor
(3-Year Warranty ; Unit)

WiRe™ App Service

II. Medicaluse 의료용(리코더 내장)



5-Year Warranty ; Compressor
(3-Year Warranty ; Unit)

WiRe™ App Service



▶ You Can Choose 4 or 5 Inner Doors Model !

▶ Note !
 Optional Freezer Racks, Boxes, LCO₂ Back Up Sys.
 and Remote Alarm Sys.
 ⇒ See the "Refrigerators, Freezer Access." Section.

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
 - With Certificate · Traceability · Delivery-Information ; Controlled by Serial Number
 - UniFreez™ High Efficient Freezing System ; Advanced System Uses a Single Compressor to Lower the Temperature to -86℃, and the Simple Mechanism Ensures Robust and Long Life-Cycle.
 - ▶ **I. Energy Saving Models ; DAIHAN's Innovative INVERTER COMPRESSOR SYSTEM™ ;**
 Optimized VSD(Variable Speed Drive) Compressor Tech. Ensures Energy Saving.
 - Power Consumption Reduction Rate ;
 Approx. 20%(Compared with Conventional Fixed-speed Compressor)
 - Reducing CO₂ Emissions for Environmental Protection
 - Low Noise / Vibration, Ensuring Longer Compressor Life.
 - ▶ **II. Medicaluse Models ; Built in Temp Chart Recorder**
Smart-Lab™ Control ; Remote Controlled by Mobile-Phone with WiRe™ App Service ;
 - Internet Connectivity with WiFi. - Storage Management & Automatic Temp Data Recording.
 - Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.
 - Self-Diagnostic Function ; Built in Standard Warning and Alarm Function
- * High-Quality Insulation Panel and Inner Doors Ensure Defrosting Around Outer Surface
 * Automatic Vacuum Breaker ; Easy to Re-Open the Door
 * Robust & Flexible Design of Door Handle * External Temp Sensor Hole
 * Filter-Free Mechanism Lowers Heat Generated by the Compressor

Specifications & Ordering information

※ 운임 별도 <M1>

I. Energy Saving Model ; Inverter Comp. Sys.		“Fre400-E86”	“Fre500-E86”	“Fre700-E80”	“Fre800-E80”
		Without Temp Chart Recorder			
Compressor		Single INVERTER COMPRESSOR ; Highly Advanced Efficiency Variable Speed Inverter Type Compressor			
Power Consumption		8.5kWh/day	8.7kWh/day	9.2kWh/day	9.4kWh/day
II. Medicaluse Approved ; Built-in Recorder		“Fre400-B86”	“Fre500-B86”	“Fre700-B80”	“Fre800-B80”
		Built in Temp Chart Recorder “KT622” Record Range ; -100℃~0℃, divi 2℃ ; 7Days/1Chart, 60Charts/Cs.			
Compressor		Single Standard Compressors ; Highly Advanced Efficiency Fixed Speed ON/OFF Type Compressor			
Power Consumption		10.2kWh/day	10.4kWh/day	13.3kWh/day	13.5kWh/day
Capacity		393 Lit	503 Lit	714 Lit	812 Lit
Temp Range		-86℃~-65℃		-80℃~-65℃	
		at R/T of 25℃ and Humidity of 70%			
Temp Resolution & Sensor		1℃ - Display, 1℃ - Control & PT100 Sensor			
Condenser-type		Block Condenser, No Pin, Filter Free Mechanism			
Refrigerant		Mixed Gas Refrigerant(DH.MGR) ; CFC Free			
Controller		Smart-Lab™ ; 7" Full Touch Screen TFT LCD			
Shelf(Included)		3ea(4-Inner Door), 4ea(5-Inner Door) Perforated Stainless-steel Shelves			
Max. Q'ty of Racks	4-Inner Door (Height 31.2cm)	16pcs of 15(5×3)- Placed Rack (RFU4234)	16pcs of 20(5×4)- Placed Rack (RFU4257)	24pcs of 20(5×4)- Placed Rack (RFU4257)	28pcs of 20(5×4)- Placed Rack (RFU4257)
	5-Inner Door (Height 25.2 cm)	20pcs of 12(4×3)- Placed Rack (RFU5234)	20pcs of 16(4×4)- Placed Rack (RFU5257)	30pcs of 16(4×4)- Placed Rack (RFU5257)	35pcs of 16(4×4)- Placed Rack (RFU5257)
Material		Internal ; Stainless-steel(#304), External ; Powder Coated Steel Inner Door ; Powder Coated Stainless-steel(#304)			
Dimension (w×d×h)	Internal (cm)	60×50×131	60×64×131	88×62×131	100×62×131
	External (cm)	84×75×190	84×89×190	112×87×190	124×87×190
Weight		260 kg	300 kg	320 kg	380 kg
Packing Size(w×d×h) (cm) & Gross Weight		116×116×222	116×124×222	140×140×221	152×140×221
		352 kg	392 kg	442 kg	502 kg
Electric Repts.		1 Phase, AC120V, 60Hz or AC 230V, 50/60Hz			

I. Energy Saving Model ; Innovative INVERTER COMPRESSOR SYSTEM™, Without Recorder

Cat. No.	230V	DH.FreE8400	DH.FreE8500	DH.FreE8700	DH.FreE8800
	(4-Inner Door)	DH.FreE8401	DH.FreE8501	DH.FreE8701	DH.FreE8801
Price, Each. ₩					
Cat. No.	230V	DH.FreE84005	DH.FreE85005	DH.FreE87005	DH.FreE88005
	(5-Inner Door)	DH.FreE84015	DH.FreE85015	DH.FreE87015	DH.FreE88015
Price, Each. ₩					

II. Medicaluse Approved Model ; Built-in Temp Chart Recorder

Cat. No.	230V	DH.FreB8400	DH.FreB8500	DH.FreB8700	DH.FreB8800
	(4-Inner Door)	DH.FreB8401	DH.FreB8501	DH.FreB8701	DH.FreB8801
Price, Each. ₩					

DAIHAN™ -86℃ ~ -65℃ Personal SMART Multiuse ULT Freezer-**UniFreez™** High Efficient Freezing System, **I. Medicaluse** & **II. Lab-use Smart-Lab™** with WiRe™ App, Programmable & Monitoring System, CFC-Free Refrigerant, Upright- & Under Bench-Type, 25 · 82Lit Ideal for Secure Storage of Vaccine · Viruses · Anatomy Components · Bio-products · DNA · Blood · Tissue Samples etc. in Lab & Industry 퍼스널(소형) 스마트 다용도 안전 초저온 백신냉동고, 백신 · 바이러스 · 해부조직물 · Bio물질 · DNA · 혈액 · 티슈시료 등을 안전 보관 유지용

5-Year Warranty ; Compressor
(3-Year Warranty ; Unit)

-86℃ ~ -65℃

(2) 25Lit Laboratory ULT Freezer

5-Year Warranty ; Compressor
(3-Year Warranty ; Unit)

-86℃ ~ -65℃

I. With or II. Without
Built in Recorder

(3) 82Lit Medicaluse Approved ULT Freezer

5-Year-Warranty ; Compressor
(3-Year Warranty ; Unit)

-86℃ ~ -65℃

(6) 82Lit Under Bench ULT Freezer



Fig1 ; Under-Bench Installation

Fig2 ; 2 Unit Stacked with Optional Stacking Panel-Set



WiRe™ App Service

Controlled by Mobile-Phone

It is Possible to Monitor and Control the Smart-Lab™ Equipment through all Mobile Phone, Tablet Without Restraint of Time and Place.



- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- UL, CUL Certified and GD mark(Only 82Lit Model)
- With Certificate · Traceability · Delivery-Information ; Controlled by Serial Number
- UniFreez™ High Efficient Freezing System ; Advanced System Uses a Single Compressor to Lower the Temperature to -86℃, and the Simple Mechanism Ensures Robust and Long Life-Cycle.
- Smart-Lab™ Control ; Remote Controlled by Mobile-Phone with WiRe™ App Service
 - Internet Connectivity with WiFi. - Storage Management & Automatic Temp. Data Recording.
 - Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.
 - Self-Diagnostic Function ; Built in Standard Warning and Alarm Function
- * Easy Maintenance and Improved Performance * Robust & Flexible Design of Door Handle
- * High-Quality Insulation Panel and Inner Doors Ensure Defrosting Around Outer Surface.
- * CFC-Free Mixed Gas Refrigerant is Used. * External Temp. Sensor Hole
- * Filter-Free Mechanism Lowers Heat Generated by the Compressor

I. Medicaluse Approved Model ; Built-in Temp Chart Recorder (See The Access. P280)

II. Laboratory / Standard Model ; Without Recorder

Specifications & Ordering information

※ 운임 별도 <M1>

Spec.	Model-type	Upright-type		Under Bench-type
I. Medicaluse Approved (Built-in Recorder)		(1) “Fre25-B86” Max. 25Lit	(3) “Fre80-B86” Max. 82Lit	(5) “Fre80-B86U” Max. 82Lit
		Built in Temp Chart Recorder Record Range ; -100℃ ~ 0℃, divi 2℃ ; 7Days/1Chart, 60Charts/Cs.		
II. Laboratory / Standard (Without Recorder)		(2) “Fre25-86” Max. 25Lit	(4) “Fre80-86” Max. 82Lit	(6) “Fre80-86U” Max. 82Lit
		Without Temp Chart Recorder		
Capacity		25 Lit	82 Lit	
Temp Range		-86℃~-65℃ at R/T of 25℃ and Humidity of 70%		
Temp Resolution & Sensor		1℃ - Display, 1℃ - Control & PT100 Sensor		
Compressor		Single Compressor ; Highly Advanced Efficiency		
Condenser-type		Block Condenser, No Pin, Filter Free Mechanism		
Refrigerant		Mixed Gas Refrigerant(DH.MGR) ; CFC Free		
Controller		Digital Feedback Controller	Smart-Lab™ ; 7" Full Touch Screen TFT LCD	
Shelf(Included)		None	1ea Stainless-steel Shelf	None
Inner Door		None	2 Door Height 252mm	None
Max. Q'ty of Racks		None, Only 10~15 Boxes are Applicable	4pcs of 12(4×3)-Placed Rack (RFU5234)	3pcs of 12(4×3)-Placed Rack or Max. 54 Box without Rack
Material	Internal	Stainless-steel(#304)		
	External	Powder Coated Steel		
	Inner Door	Powder Coated Stainless-steel(#304)		
Dimension (w×d×h)	Internal (cm)	20×42×30	34×48.5×52	52×48.5×34
	External (cm)	39×62.7×78.5	53×68.7×106	106.5×72.2×61.2
Weight		52 kg	60 kg	
Packing Size(w×d×h) (cm) & Gross Weight		88×86×130	96×76×135	135×76×96
		69 kg	102 kg	
Energy Consumption at -80℃		650VA	710VA	
Electric Reqts.		1Phase, AC 120V, 60Hz or AC 230V, 50/60Hz		

I. Medicaluse Approved Model ; Built-In Temp Chart Recorder

Cat. No.	230V	DH.FreB8025	DH.FreB8080	DH.FreB8080U
	120V	DH.FreB8026	DH.FreB8081	DH.FreB8081U
Price, Each. ₩				

II. Laboratory / Standard Model ; Without Recorder

Cat. No.	230V	DH.FreL8025	DH.FreL8080	DH.FreL8080U
	120V	DH.FreL8026	DH.FreL8081	DH.FreL8081U
Price, Each. ₩				

Optional Stacking Panel-Set for Under Bench-Type

DH.SWUF.SP02 (a) Stacking Panel-Set, 590×126mm(3mm Thickness)

Set.	Price, ₩
2ea	



▶ Note ! Optional Freezer Racks, Boxes, LCO₂ Back Up Sys. and Remote Alarm Sys.
⇒ See the "Refrigerators, Freezer Access." Section.

DAIHAN™ -45°C ~ -15°C **SMART** Multiuse Freezer-**UniFreez™** High Efficient Freezing System, *Medicaluse*
Smart-Lab™ with WiRe™ App, Programmable & Monitoring System, CFC-Free Refrigerant, Built-In Recorder, 393 · 503 · 714 · 812Lit
 Ideal for Secure Storage of Vaccine · Viruses · Anatomy Components · Bio-Products · DNA · Blood · Tissue Samples etc. in Lab & Industry
 스마트 다용도 안전 초저온 냉동고, 백신 · 바이러스 · 해부조직물 · Bio물질 · DNA · 혈액 · 티슈시료 등을 안전 보관 유지용



Fig. 5-Chamber & 5-Inner Door



WiRe™ App Service
Controlled by Mobile-Phone

It is Possible to Monitor and Control the Smart-Lab™ Equipment through all Mobile Phone, Tablet Without Restraint of Time and Place.

- Supplied by DAIHAN Korea ■ **PL (Product Liability) Insurance**
- With Certificate · Traceability · Delivery-Information ; Controlled by Serial Number
- UniFreez™ High Efficient Freezing System ; Advanced System Uses a Single Compressor to Lower the Temperature to -45°C, and the Simple Mechanism Ensures Robust and Long Life-Cycle.
- **Smart-Lab™ Control ; Remote Controlled by Mobile-Phone with WiRe™ App Service**
 - Internet Connectivity with WiFi. - Storage Management & Automatic Temp Data Recording.
 - Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.
 - Self-Diagnostic Function ; Built in Standard Warning and Alarm Function
- **Built in Temperature Chart Recorder (Details ; See the Access. P280)**
 - * Easy Maintenance and Improved Performance
 - * High-Quality Insulation Panel and Inner Doors Ensure Defrosting Around Outer Surface
 - * Automatic Vacuum Breaker ; Easy to Re-Open the Door
 - * Robust & Flexible Design of Door Handle * External Temp Sensor Hole
 - * Filter-Free Mechanism Lowers Heat Generated by the Compressor

Specifications & Ordering information

<M1>

Medicaluse Approved (Built-in Recorder)	Built in Temp Chart Recorder			
	Record Range ; -100°C~0°C, divi 2°C ; 7Days/1Chart, 60Charts/Cs.			
	"Fre400-B45" Max. 393Lit	"Fre500-B45" Max. 503Lit	"Fre700-B45" Max. 714Lit	"Fre800-B45" Max. 812Lit
Capacity	393 Lit	503 Lit	714 Lit	812 Lit
Temp Range	-45°C ~ -15°C at R/T of 25°C and Humidity of 70%			
Temp Resolution & Sensor	1°C - Display, 1°C - Control & PT100 Sensor			
Compressor	Single Compressor ; Highly Advanced Efficiency			
Condenser-type	Block Condenser, No Pin, Filter Free Mechanism			
Refrigerant	Mixed Gas Refrigerant(DH.MGR) ; CFC Free			
Controller	Smart-Lab™ ; 7" Full Touch Screen TFT LCD			
Shelf (Included)	3ea Perforated Stainless-steel Shelves			
Inner Door	5-Door(Height 25.2cm)			
Max. Q'ty of Racks	20pcs of 12(4×3)- Placed Rack (RFU5234)	20pcs of 16(4×4)- Placed Rack (RFU5257)	30pcs of 16(4×4)- Placed Rack (RFU5257)	35pcs of 16(4×4)- Placed Rack (RFU5257)
Material	Internal	Stainless-steel(#304)		
	External	Powder Coated Steel		
	Inner Door	Powder Coated Stainless-steel(#304)		
Dimension (w×d×h)	Internal (cm)	60×50×131	60×64×131	88×62×131
	External (cm)	84×75×190	84×89×190	112×87×190
Weight		210 kg	240 kg	320 kg
Packing Size(w×d×h) (cm)		116×116×222	116×124×222	140×140×221
& Gross Weight		296 kg	327 kg	442 kg
Power Consumption		720VA	720VA	720VA
Electric Reqts.		1 Phase, AC120V, 60Hz or AC 230V, 50/60Hz		
Cat. No.	230V	DH.FreB34005	DH.FreB35005	DH.FreB37005
	120V	DH.FreB34015	DH.FreB35015	DH.FreB37015
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

▶ Note !

Optional Freezer Racks, Boxes, LCO₂ Back Up Sys. and Remote Alarm Sys.
 ⇒ See the "Refrigerators, Freezer Access." Section.



(1) Sliding Drawer-Type Rack, 15 Place
"RFU4234SD"


Optional Freezer Rack, Sliding Drawer-Type & Side Access-Type For DAIHAN™ Freezer, Stainless-steel(304)

■ Supplied by DAIHAN Korea * Not Include Freezer Boxes

(1) Sliding Drawer-Type for 0.5~2ml, 3~4ml, 5ml Vial Box

<L1>


Cat. Number	Model	Rack Size(mm) (w×d×h)	Freezer Box & No. of Box	Freezer Capacity (No. of Rack)	Ea. Price, ₩
For 80Lit & 80Lit Under					
DH.RFU5234SD	“RFU5234SD”	140×420×232	0.5~2㎖, 12(4×3)	80Lit (4ea) 80Lit Under (3ea)	
DH.RFU5434SD	“RFU5434SD”	140×420×171	3~4㎖, 6(2×3)		
DH.RFU5534SD	“RFU5534SD”	140×420×211	5㎖, 6(2×3)		
For 3 Door-Type 300Lit					
DH.RFU5257SD	“RFU5257SD”	140×560×232	0.5~2㎖, 16(4×4)	300Lit(12ea)	
DH.RFU5457SD	“RFU5457SD”	140×560×171	3~4㎖, 8(2×4)		
DH.RFU5557SD	“RFU5557SD”	140×560×211	5㎖,8(2×4)		
For 4 Door-Type 400・500・700・800Lit					
DH.RFU4234SD	“RFU4234SD”	140×420×290	0.5~2㎖, 15(5×3)	400Lit(16ea)	
DH.RFU4434SD	“RFU4434SD”	140×420×256	3~4㎖, 9(3×3)		
DH.RFU4534SD	“RFU4534SD”	140×420×311	5㎖, 9(3×3)		
DH.RFU4257SD	“RFU4257SD”	140×560×290	0.5~2㎖, 20(5×4)	500Lit(16ea)	
DH.RFU4457SD	“RFU4457SD”	140×560×256	3~4㎖, 12(3×4)	700Lit(24ea)	
DH.RFU4557SD	“RFU4557SD”	140×560×311	5㎖, 12(3×4)	800Lit(28ea)	
For 5 Door-Type 400・500・700・800Lit					
DH.RFU5234SD	“RFU5234SD”	140×420×232	0.5~2㎖, 12(4×3)	400Lit(20ea)	
DH.RFU5434SD	“RFU5434SD”	140×420×171	3~4㎖, 6(2×3)		
DH.RFU5534SD	“RFU5534SD”	140×420×211	5㎖, 6(2×3)		
DH.RFU5257SD	“RFU5257SD”	140×560×232	0.5~2㎖, 16(4×4)	500Lit(20ea)	
DH.RFU5457SD	“RFU5457SD”	140×560×171	3~4㎖, 8(2×4)	700Lit(30ea)	
DH.RFU5557SD	“RFU5557SD”	140×560×211	5㎖, 8(2×4)	800Lit(35ea)	



(2) Side Access-Type Rack, 20 Place
"RFU4257SA"

(2) Side Access-Type for 0.5~2ml, 3~4ml, 5ml Vial Box

Cat. Number	Model	Rack Size(mm) (w×d×h)	Freezer Box & No. of Box	Freezer Capacity (No. of Rack)	Ea. Price, ₩
For 80Lit & 80Lit Under					
DH.RFU5234SA	“RFU5234SA”	140×420×232	0.5~2㎖, 12(4×3)	80Lit (4ea) 80Lit Under (3ea)	
DH.RFU5434SA	“RFU5434SA”	140×420×171	3~4㎖, 6(2×3)		
DH.RFU5534SA	“RFU5534SA”	140×420×211	5㎖, 6(2×3)		
For 3 Door-Type 300Lit					
DH.RFU5257SA	“RFU5257SA”	140×560×232	0.5~2㎖, 16(4×4)	300Lit(12ea)	
DH.RFU5457SA	“RFU5457SA”	140×560×171	3~4㎖, 8(2×4)		
DH.RFU5557SA	“RFU5557SA”	140×560×211	5㎖, 8(2×4)		
For 4 Door-Type 400・500・700・800Lit					
DH.RFU4234SA	“RFU4234SA”	140×420×290	0.5~2㎖, 15(5×3)	400Lit(16ea)	
DH.RFU4434SA	“RFU4434SA”	140×420×256	3~4㎖, 9(3×3)		
DH.RFU4534SA	“RFU4534SA”	140×420×311	5㎖, 9(3×3)		
DH.RFU4257SA	“RFU4257SA”	140×560×290	0.5~2㎖, 20(5×4)	500Lit(16ea)	
DH.RFU4457SA	“RFU4457SA”	140×560×256	3~4㎖, 12(3×4)	700Lit(24ea)	
DH.RFU4557SA	“RFU4557SA”	140×560×311	5㎖, 12(3×4)	800Lit(28ea)	
For 5 Door-Type 400・500・700・800Lit					
DH.RFU5234SA	“RFU5234SA”	140×420×232	0.5~2㎖, 12(4×3)	400Lit(20ea)	
DH.RFU5434SA	“RFU5434SA”	140×420×171	3~4㎖, 6(2×3)		
DH.RFU5534SA	“RFU5534SA”	140×420×211	5㎖, 6(2×3)		
DH.RFU5257SA	“RFU5257SA”	140×560×232	0.5~2㎖, 16(4×4)	500Lit(20ea)	
DH.RFU5457SA	“RFU5457SA”	140×560×171	3~4㎖, 8(2×4)	700Lit(30ea)	
DH.RFU5557SA	“RFU5557SA”	140×560×211	5㎖, 8(2×4)	800Lit(35ea)	



Spare Temperature Chart Recorder & Accessories

<L1>

- DH.KT622** (1) Built-in Temp Chart Recorder, "KT622", with "C657" Chart Paper
Record Range ; -100℃ ~ 0℃, divi 2℃ ; 7Days/1Chart, 60Charts/Cs.
DH.C657 (2) Chart Paper, "C657", 60 Charts, for 420 Days Recording
DH.P222 (3) Replacement Pens, "P222", Gea Red Color, for "KT622" Chart Recorder

Optional LCO₂ Back Up System and Remote Alarm System

- DH.LCO2B** (4) LCO₂ Back Up System, "LCO2B"
DH.RWAS (5) Remote Wire/Wireless Alarm System, "RWAS"

Ea.	Price, ₩

Optional DAIHAN Wisd™ Cardboard Freezer Box

■ Supplied by DAIHAN Korea

- (1) For 1.5/2.0ml Cryovial / Microtube, 81 Holes Cell ϕ 13mm,
Cardboard Storage Box with Cardboard / Polypropylene Divider

■ 시료의 위치 확인이 용이하도록 숫자와 알파벳이 박스 내/외 1면에 표기 되어있음

- DH.Box901** with Cardboard Divider Box, 81Holes $\times$ ϕ 13mm, 133 $\times$ 133 $\times$ h52mm
DH.Box902 with PP Divider Box, 81Holes $\times$ ϕ 13mm, 133 $\times$ 133 $\times$ h52mm

- (2) For 0.5~2ml Cryovials/Microtubes, 134 $\times$ 134 $\times$ h52mm

■ 효율적인 냉각을 위한 Vapor Exchange Hole 타입
■ Cooling Vapor Exchange Hole for Drainage and Rapid Freezing

- DH.Box903** 81-holes Box, White, h52mm, 81Holes $\times$ ϕ 13mm
DH.Box904 100-holes Box, White, h52mm, 100Holes $\times$ ϕ 12mm

- (3) For 3/4ml Cryovials, 134 $\times$ 134 $\times$ h77mm

■ 효율적인 냉각을 위한 Vapor Exchange Hole 타입
■ Cooling Vapor Exchange Hole for Drainage and Rapid Freezing

- DH.Box905** 81-holes Box, White, h77mm, 81Holes $\times$ ϕ 13mm

- (4) For 15 & 50ml Conical Tube, 150 $\times$ 150 $\times$ h122mm


■ 효율적인 냉각을 위한 Vapor Exchange Hole 타입
■ Cooling Vapor Exchange Hole for Drainage and Rapid Freezing

- DH.Box906** 36-holes Box, White, h122mm, 36Holes $\times$ ϕ 23.5mm, for 15ml
DH.Box907 16-holes Box, White, h122mm, 16Holes $\times$ ϕ 35.5mm, for 50ml

- (5) For 0.5~2ml Cryovial/Microtube, 3 Color Cardboard Freezer Box,
81 & 100 Holes, 133 $\times$ 133 $\times$ h52mm

■ Hinged Lid로 뚜껑의 여닫음이 편리함.
* Alphanumeric Identifications ; The Lid, Bottom and Four Sides
* 81-Holes Box Dimensions ; 133 $\times$ 133 $\times$ h52mm, 81Holes $\times$ ϕ 13mm
* 100-Holes Box Dimensions ; 133 $\times$ 133 $\times$ h52mm, 100Holes $\times$ ϕ 12mm

Cat. No.	81-Holes Color	Ea. Price, ₩	Cat. No.	100-Holes Color	Ea. Price, ₩
DH.Box908.B	Blue		DH.Box909.B	Blue	
DH.Box908.G	Green		DH.Box909.G	Green	
DH.Box908.W	White		DH.Box909.W	White	



Freezer with Validation (IQ, OQ) Service

(1) Laboratory Freezer

<L1>

Cat. No.	Model	Price, ₩	Cat. No.	Model	Price, ₩
DH.FreE9400.VS	"Fre400-E95"		DH.FreE8700.VS	"Fre700-E80"	
DH.FreE9500.VS	"Fre500-E95"		DH.FreE8800.VS	"Fre800-E80"	
DH.FreE9700.VS	"Fre700-E90"		DH.FreE84005.VS *	"Fre400-E86"	
DH.FreE9800.VS	"Fre800-E90"		DH.FreE85005.VS *	"Fre500-E86"	
DH.FreE94005.VS *	"Fre400-E95"		DH.FreE87005.VS *	"Fre700-E80"	
DH.FreE95005.VS *	"Fre500-E95"		DH.FreE88005.VS *	"Fre800-E80"	
DH.FreE97005.VS *	"Fre700-E90"		DH.FreL8025.VS	"Fre25-86"	
DH.FreE98005.VS *	"Fre800-E90"		DH.FreL8080.VS	"Fre80-86"	
DH.FreE8400.VS	"Fre400-E86"		DH.FreL8080U.VS	"Fre80-86U"	
DH.FreE8500.VS	"Fre500-E86"				

* 5-Inner Door

(2) Vaccine · Pharmaceutical & Blood Bank Freezer

Cat. No.	Model	Price, ₩	Cat. No.	Model	Price, ₩
DH.FreB9400.VS	"Fre400-B95"		DH.FreB8080U.VS	"Fre80-B86U"	
DH.FreB9500.VS	"Fre500-B95"		DH.FreB34005.VS *	"Fre400-B45"	
DH.FreB9700.VS	"Fre700-B90"		DH.FreB35005.VS *	"Fre500-B45"	
DH.FreB9800.VS	"Fre800-B90"		DH.FreB37005.VS *	"Fre700-B45"	
DH.FreB8400.VS	"Fre400-B86"		DH.FreB38005.VS *	"Fre800-B45"	
DH.FreB8500.VS	"Fre500-B86"				
DH.FreB8700.VS	"Fre700-B80"				
DH.FreB8800.VS	"Fre800-B80"				
DH.FreB8025.VS	"Fre25-B86"				
DH.FreB8080.VS	"Fre80-B86"				

* 5-Inner Door

DAIHAN™ Digital Bio Freezer “WLF”, -35℃ ~ -15℃, 330 · 420 · 600 Lit

With 3 Shelf, Digital Fuzzy Control System, Fast Freezing Function, CFC-Free Refrigeration, Certi. & Traceability

바이오 냉동고, 디지털 퍼지 컨트롤, 친환경 냉매(404A) 채용, 급속 냉동 기능



3-Year Warranty

Digital Bio Freezer “WLF”



with STS Inner Doors “WLF”

- Supplied by DAIHAN Korea ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Fast Freezing Function
- Patented Jog-Shuttle Control System
- CFC-Free (404A) Refrigeration System

- * Suitable for Storage of Bio Reagent, Seed, Enzyme and Various Samples
- * Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- * Storage Function ; Set Values of Temperature and Timer
- * Digital LCD with Back-Light Function
- * Locking Mode supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * Corrosion Resistant 304 Stainless-steel Chamber and Shelves
- * Over Current Protection, Sensor Error Detection, Leakage Breaking

Specifications & Ordering information

<L1>

Spec. / Model	“WLF-320” Max. Capa. 330 Lit	“WLF-420” Max. Capa. 420 Lit	“WLF-620” Max. Capa. 600 Lit
Capacity	330 Lit	420 Lit	600 Lit
Type	Upright Type		
Temp Range & Sensor	-35℃~-15℃, PT100		
Temp Resolution	1℃-Display, 1℃-Control		
Refrigerant & Compressor	CFC-Free(404A) Refrigeration System, 3/4HP	CFC-Free(404A) Refrigeration System, 1HP	
Display	Digital LCD with Back-Light Function		
Controller	Digital Controller with Jog-Shuttle Switch (Turn+Push)		
Shelf (Included)	3 Stainless-steel Shelf(# 304)		
Safety Device	Over Current Protector, Sensor Error Detector, Leakage Breaker		
Material	Stainless-steel (# 304)-Internal, Powder Coated steel-External		
Dimension ; Interior & Exterior	600×550×h1,000mm, 820×700×h1,700mm	600×700×h1,000mm, 820×850×h1,700mm	780×800×h1,000mm, 1000×950×h1,700mm
Packing Size & Gross Weight	1,160×910×h2,160mm 290 kg	1,160×1,060×h2,160mm 314 kg	1,340×1,160×h2,160mm 350 kg
Power Consumption	1.7 kW	1.8 kW	1.8 kW
Electrical Reqts.	1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz		
Cat. No.	230V	DH.WLF00320	DH.WLF00420
	120V	DH.WLF01320	DH.WLF01420
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Stainless-steel Flat Shelf

- DH.WLF10320 Flat Shelf, “LFS320”, w580×d530mm, for “WLF-320”
 DH.WLF10420 , “LFS420”, w580×d680mm, for “WLF-420”
 DH.WLF10620 , “LFS620”, w760×d780mm, for “WLF-620”

Ea. Price, ₩

DAIHAN™ 0°C~10°C SMART Laboratory Refrigerator "Ref.L", 150·280·600·1300Lit, Glass Door

With Smart-Lab™ System, Dual Eva-Defrost. CFC-Free, Forced-Air, Stainless-steel Perforated Shelf, Door Lock Device, 실험실용 다용도 냉장고



3-Year Warranty

(2) & (3) 280·600Lit with 1-Glass Door

WiRe™ App Service



3-Year Warranty

(4) 1300Lit with 2-Glass Door



(a) Stainless-steel Perforated Shelf



(b) Temp Chart Recorder, "KT622"



(c) Remote Alarm System "RWAS"

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance
- With Certificate · Traceability · Delivery-information ; Controlled by Serial Number
- Dual Evaporator for Natural Defrosting Ensures High Accuracy and Uniformity of Temp
- Smart-Lab™ Control ; Remote Controlled by Mobile-Phone with WiRe™ App Service
 - Internet Connectivity with WiFi. - Storage Management & Automatic Temp Data Recording.
 - Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.
 - Self-Diagnostic Function ; Built in Standard Warning and Alarm Function
- * Height-adjustable Perforated Stainless-steel Shelves
- * Forced-air Circulation Fan ; Superior Temp Uniformity and Quick Recovery Time.
- * Self-diagnostic Function for Air Circulation Fan
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * Double Fold Tempered Glass Window * Built in Battery for Power Failure Alarm
- * Alarm Set Points of Temperature Deviation at +5.5 and -1.5 Under Temperature
- * Automatic Condensate Removal, No Drain Line Required
- * Available Built in Chart Recorder and Remote Alarm System(Optional)

Specifications & Ordering information

<L1>

0~10°C Laboratory Refrigerator					
Spec.	Model	(1) "Ref150.L"	(2) "Ref300.L"	(3) "Ref600.L"	(4) "Ref1300.L"
Capacity		150 Lit	280 Lit	600 Lit	1300 Lit
Door		1-Door	1-Door	1-Door	2-Door
Temp Range		0°C~10°C at R/T of 30°C and Humidity of 70%			
Temp Resolution & Sensor		1°C - Display, 1°C - Control & PT100 Sensor			
Refrigerant		CFC-Free(R-404A) Refrigeration System			
Compressor & Evaporator		1/2HP, Auto Exchange Dual Evaporator System for Temperature Accuracy			
Air Circulation		Forced-air Type			
		2 EC Fan	3 EC Fan	4 EC Fan	6 EC Fan
Shelf (Included)		Stainless-steel Perforated Shelf with Front Guard			
		2	3	5	10
Material		Internal ; Stainless-steel, External ; Powder Coated-steel			
Safety		Door Lock System, Sensor Error Detection			
Dimension (w×d×h)	Internal (cm)	57.5×53.2×50	59×71.5×67	62×77×129	137×77×129
	External (cm)	72.2×80.3×143.6	75×100.8×162	75×98×195	149×98×195
Net Weight		90 kg	130 kg	150 kg	208 kg
Power Consumption		650 W	670 W	700 W	750 W
Electric Reqts.		1 Phase, AC120V, 60Hz or AC 230V, 50/60Hz			
Cat. No.	230V	DH.RefL0150	DH.RefL0300	DH.RefL0600	DH.RefL1300
	120V	DH.RefL0151	DH.RefL0301	DH.RefL0601	DH.RefL1301
Price, Each. ₩					
Validation(IQ, OQ)		DH.RefL0150.VS	DH.RefL0300.VS	DH.RefL0600.VS	DH.RefL1300.VS
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Stainless-steel Perforated Shelf

- DH.RefS8150 (a) STS Perforated Shelf, w477×d407 mm, for "Ref150.L"
- DH.RefS8300 , w492×d590 mm, for "Ref300.L"
- DH.RefS8600 , w532×d610 mm, for "Ref600.L"
- DH.RefS81300 , w520×d650 mm, for "Ref1300.L"

Ea. Price, ₩

Optional Accessory

- DH.KT622 (b) Temp Chart Recorder, "KT622", with "C657" Chart Paper
- DH.RWAS (c) Remote Wire/Wireless Alarm System, "RWAS"

Ea. Price, ₩

Articles

Cat. No

Description

DAIHAN™ 2°C~10°C **SMART** Blood Bank Refrigerator “Ref.B”, 150·280·600·1300Lit, *Medicaluse*
 With **Smart-Lab™** System, CFC-Free(R-404A), Forced-Air, Built in Temperature Chart Recorder, Dual Eva-Defrost, Glass Door
 Stainless-steel Perforated Drawer Shelf, Door Lock Device, 혈액 냉장고, 의료용



3-Year Warranty

WiRe™ App Service

(2) & (3) 280·600Lit with 1-Glass Door



3-Year Warranty

(4) 1300Lit with 2-Glass Door



(a) Stainless-steel Perforated Drawer Shelf



Temp Chart Recorder, “KT622”



(b) Remote Alarm System “RWAS”

- Supplied by DAIHAN Korea ■ **PL(Product Liability) Insurance**
- With Certificate · Traceability · Delivery-information ; Controlled by Serial Number
- Dual Control System for Safety ; Even if one of the Two Main Controller Boards Fails to Run, the Other One is still Running properly for Safe Storage of Sample.
- Dual Evaporator for Natural Defrosting Ensures High Accuracy and Uniformity of Temp
- **Smart-Lab™ Control ; Remote Controlled by Mobile-Phone with WiRe™ App Service**
 - Internet Connectivity with WiFi. - Storage Management & Automatic Temp Data Recording.
 - Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.
 - Self-Diagnostic Function ; Built in Standard Warning and Alarm Function
- * Forced-air Circulation Fan ; Superior Temp Uniformity and Quick Recovery Time.
- * Self-diagnostic Function for Air Circulation Fan
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * Double Fold Tempered Glass Window * Built in Battery for Power Failure Alarm
- * Alarm Set Points of Temperature Deviation at +5.5 and -1.5 Under Temperature
- * Automatic Condensate Removal, No Drain Line Required
- * Available Built in Remote Alarm System(Optional)

Specifications & Ordering information

<M1>

2°C~10°C Blood Bank Refrigerator					
Spec.	Model	(1) “Ref150.B”	(2) “Ref300.B”	(3) “Ref600.B”	(4) “Ref1300.B”
Capacity		150 Lit	280 Lit	600 Lit	1300 Lit
Door		1-Door			2-Door
Temp Range		2°C~10°C(4°C±2°C) at R/T of 30°C and Humidity of 70%			
Temp Resolution & Sensor		1°C - Display, 1°C - Control & PT100 Sensor			
Refrigerant		CFC-Free(R-404A) Refrigeration System			
Compressor & Evaporator		1/2HP, Auto Exchange Dual Evaporator System for Temperature Accuracy			
Air Circulation		Forced-air Type			
		2 EC Fan	3 EC Fan	4 EC Fan	6 EC Fan
Shelf (Included)		Stainless-steel Perforated Drawer Shelf			
		2	3	5	10
Built in Chart Recorder		for Temp Data, Includes 60 Chart Paper			
Material		Internal ; Stainless-steel, External ; Powder Coated-steel			
Safety		Door Lock System, Sensor Error Detection			
Dimension (w×d×h)	Internal (cm)	57.5×53.2×50	59×71.5×67	62×77×129	137×77×129
	External (cm)	72.2×80.3×143.6	75×100.8×162	75×98×195	149×98×195
Net Weight		90 kg	130 kg	150 kg	208 kg
Power Consumption		650 W	670 W	700 W	750 W
Electric Rqts.		1 Phase, AC120V, 60Hz or AC 230V, 50/60Hz			
Cat. No.	230V	DH.RefB0150	DH.RefB0300	DH.RefB0600	DH.RefB1300
	120V	DH.RefB0151	DH.RefB0301	DH.RefB0601	DH.RefB1301
Price, Each. ₩					
Validation(IQ, OQ)		DH.RefB0150.VS	DH.RefB0300.VS	DH.RefB0600.VS	DH.RefB1300.VS
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Drawer Shelf, Replacement Chart Paper and Pens

- DH.RefS8150D (a) Drawer STS Shelf, w517×d407mm, for “Ref150.B”
 DH.RefS8300D , w532×d590mm, for “Ref300.B”
 DH.RefS8600D , w572×d610mm, for “Ref600.B”
 DH.RefS81300D , w560×d650mm, for “Ref1300.B”

- DH.C657 Replacement Chart Paper, 60Charts for 420Days Recording
 DH.P222 Replacement Pens, 6Pens, Red Color for Chart Recorder

Optional Accessory

- DH.RWAS (b) Remote Wire/Wireless Alarm System, “RWAS”

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ 2°C~10°C Vaccine & Pharm. Refrigerator "Ref.P", 105 & 295Lit, *Medicaluse*

With Microprocessor PID Controller, 60mm(t) Polyurethane Insulation, Vacuum Sealed Glass Door, CFC/HCFC-Free, Forced Air Circulation Wire Shelf, Superior Temp Uniformity, **시약 약품 보관용 냉장고, 의료용**



- Supplied by DAIHAN Korea ■ PL (Product Liability) Insurance
- With Certificate · Traceability · Delivery-information ; Controlled by Serial Number
- * Energy Saving LED Lamp Installation
- * High Density Polyurethane Insulation(60mm thick) Provides Excellent Low Temperature
- * Invisible Heating Film Protects from water Condensation on the Glass Door
- * Self-Diagnostic and Visual & Audible Alarm System ;
 - High/Low Working Room Temp Alarm
 - Electric Power Failure Alarm
 - Status of Door Opening Warning Alarm
 - Low Battery for Controller Warning Alarm
 - Password Locking System for Set Temp

Specifications & Ordering information

<M2>

Spec.	Model	"Ref100.P" 105 Lit, Bench Top	"Ref290.P" 295 Lit, Floor Standing
Temp Range		2°C~10°C	
Shelf (Included)		1 Basket & 2 Wire Shelf	3 Wire Shelf
Compressor		1 / 8 HP	1 / 4 HP
Recorder		Option	
Defrost Heater		300 W, Automatic Defrost and Dehumidify System	
Forced Air Circulation		Turbo Fan with Automatic Stop During Door Open	
Controller & Display		Microprocessor Digital Controller / 4 Digit LED Display	
Door		Outer Door with Glass Window with Safety Door Lock System	
Safety		High / Low Temp Alarm, Electric Current Leakage Breaker	
Material		Internal ; Stainless-steel, External ; Epoxy Powder Coated-steel	
Dimension (w×d×h), (cm)	Interior	46×42×h63	56×66×h80
	Exterior	56×63×h91.5	68.3×79.6×h149
	Packing	72×63×h105	85×72×h186
Net & Gross Weight		65 & 80 kg	170 & 185 kg
Power Consumption		100 W	195 W
Electric Rqts.		220V, 50/60Hz, 1 Phase	
Cat. No.		DH.RefP0100	DH.RefP0290
Price, Each. ₩			
Validation(IQ, OQ)		DH.RefP5100	DH.RefP5290
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Spare Stainless-steel Teflon-Coated Steel Wire Shelf & Basket

- DH.RefA0005 (1) Teflon-coated Wire Shelf, w440×d350 mm, for "Ref100.P"
- DH.RefA0006 , w590×d470 mm, for "Ref290.P"
- DH.RefA0010 (2) Teflon-coated Wire Basket, w440×d340×h175 mm, for "Ref100.P"

Ea. Price, ₩

Optional Accessory

- DH.SteA1011 (3) Temp Chart Recorder, "STE A1011",
with "STE A1012" Chart Paper, for "Ref290.P"
- DH.SteA1012 Chart Paper, "STE A1012", 33 Days
- DH.SteA1005 (4) Thermal Printer, Paper Roll, "STE A1005", for "Ref100.P"
- DH.SteA1006 (5) Thermal Paper Package, 4 Paper Roll, for "STE A1005"
- DH.SteA1010 Mobile Alarm System, "STE A-Mobile"

Ea. Price, ₩



Articles

Cat. No

Description

DAIHAN™ Premium 2°C~10°C SMART Vaccine & Pharm. Refrigerator "Ref.P", *Medicaluse*
 With Smart-Lab™ System, Dual Eva-Defrost. CFC-Free(R-404A), Forced-Air, 150 · 280 · 600 · 1300Lit, Glass Door
 Stainless-steel Perforated Shelf & Perforated Drawer Shelf, Door Lock Device, 백신 & 약품 냉장고, 의료용



3-Year Warranty

(2) & (3) 280 · 600Lit with 1-Glass Door

WiRe™ App Service



3-Year Warranty

(4) 1300Lit with 2-Glass Door



(a) Stainless-steel Perforated Shelf



(b) Stainless-steel Perforated Drawer Shelf



(c) Temp Chart Recorder, "KT622"



(d) Remote Alarm System "RWAS"

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- With Certificate · Traceability · Delivery-information ; Controlled by Serial Number
- Dual Control System for Safety ; Even if one of the Two Main Controller Boards Fails to Run, the Other One is still Running properly for Safe Storage of Sample.
- Dual Evaporator for Natural Defrosting Ensures High Accuracy and Uniformity of Temp
- Smart-Lab™ Control ; Remote Controlled by Mobile-Phone with WiRe™ App Service
 - Internet Connectivity with WiFi. - Storage Management & Automatic Temp Data Recording.
 - Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.
 - Self-Diagnostic Function ; Built in Standard Warning and Alarm Function
- * Forced-air Circulation Fan ; Superior Temp Uniformity and Quick Recovery Time.
- * Self-diagnostic Function for Air Circulation Fan
- * Double Fold Tempered Glass Window
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * Alarm Set Points of Temperature Deviation at +5.5 and -1.5 Under Temperature
- * Automatic Condensate Removal, No Drain Line Required
- * Built in Battery for Power Failure Alarm
- * Available Built in Chart Recorder and Remote Alarm System(Optional)

Specifications & Ordering information

<M1>

2°C~10°C Vaccine & Pharm. Refrigerator					
Spec.	Model	(1) "Ref150.P"	(2) "Ref300.P"	(3) "Ref600.P"	(4) "Ref1300.P"
Capacity		150 Lit	280 Lit	600 Lit	1300 Lit
Door		1-Door	1-Door	1-Door	2-Door
Temp Range		2°C~10°C at R/T of 30°C and Humidity of 70%			
Temp Resolution & Sensor		1°C - Display, 1°C - Control & PT100 Sensor			
Refrigerant		CFC-Free(R-404A) Refrigeration System			
Compressor & Evaporator		1/2HP, Auto Exchange Dual Evaporator System for Temperature Accuracy			
Air Circulation		Forced-air Type			
		2 EC Fan	3 EC Fan	4 EC Fan	6 EC Fan
Shelf (Included)		Stainless-steel Perforated Shelf with Front Guard			
		2	3	3	6
		Stainless-steel Perforated Drawer Shelf			
		0	0	2	4
Material		Internal ; Stainless-steel, External ; Powder Coated-steel			
Safety		Door Lock System, Sensor Error Detection			
Dimension (w×d×h)	Internal (cm)	57.5×53.2×50	59×71.5×67	62×77×129	137×77×129
	External (cm)	72.2×80.3×143.6	75×100.8×162	75×98×195	149×98×195
Net Weight		90 kg	130 kg	150 kg	208 kg
Power Consumption		650 W	670 W	700 W	750 W
Electric Reqs.		1 Phase, AC120V, 60Hz or AC 230V, 50/60Hz			
Cat. No.	230V	DH.RefP0150	DH.RefP0300	DH.RefP0600	DH.RefP1300
	120V	DH.RefP0151	DH.RefP0301	DH.RefP0601	DH.RefP1301
Price, Each. ₩					
Validation(IQ, OQ)		DH.RefP0150.VS	DH.RefP0300.VS	DH.RefP0600.VS	DH.RefP1300.VS
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Stainless-steel Shelf

- DH.RefS8150 (a) Perforated STS Shelf, w477×d407mm, for "Ref150.P"
- DH.RefS8300 , w492×d590mm, for "Ref300.P"
- DH.RefS8600 , w532×d610mm, for "Ref600.P"
- DH.RefS81300 , w520×d650mm, for "Ref1300.P"
- DH.RefS8600D (b) Drawer STS Shelf, w572×d610, for "Ref600.P"
- DH.RefS81300D , w560×d650, for "Ref1300.P"

Optional Accessory

- DH.RefS8150D (b) Drawer STS Shelf, w517×d407mm, for "Ref150.P"
- DH.RefS8300D , w532×d590mm, for "Ref300.P"
- DH.KT622 (c) Temp Chart Recorder, "KT622", with "C657" Chart Paper
- DH.RWAS (d) Remote Wire/Wireless Alarm System, "RWAS"

DAIHAN Scientific

SHAKERS

Mixers , Rocking & Rolling Tube Mixers	P. 287, 293
, Vortex Mixers	P. 288
Shakers , Orbital Motion	P. 289~292
, Reciprocating Motion	P. 291~292
, Rockers	P. 293
, Rotators	P. 294
, Sieve Shaker	P. 295

• Mixer, Tube Mixer



• Mixer, Vortex



• Shaker, Orbital/Reciprocating



• Shaker, Rocker



• Shaker, Rotator



• Shaker, Sieve Shaker



DAIHAN™ Hematology Tube Mixer "MIXR-40", with 5 Mixing Rollers, Rolling-Motion, Fixed 40rpm

Ideal for 15~50ml Test Tube, with Minimum Noise Motor and 0~60min Timer, 튜브 롤러 믹서, 헤마톨로지용



3-Year Warranty

- Supplied by DAIHAN Korea ■ with 5 Mixing Rollers($\phi 32 \times L325\text{mm}$)
- Ideal for Mixing Blood Samples, Viscous Substances and Liquid-Solid Suspensions, Use in a Variety of Applications, Prevention of Blood Coagulation and Immune Precipitation, etc.
- Rolling Mixing-motion for Complete Sample Mixing

- * Stability Mixing Motion and Built in Minimize Noise AC Synchronous Motor
- * Can Use with 8ea $\times$ 15ml Test Tubes or 8ea $\times$ 50ml Test Tubes. * Stability Fixed 40rpm Mixing Speed
- * Built in 99hr 59min Timer for User Convenience * Can be Used in Cold Lab Chambers or Incubators

Specifications & Ordering information

<L1>

Spec.	Model	"MIXR-40"
No. of Mixing Roller		5 Mixing Rollers
Mixing Speed		Fixed 40rpm
Mixing Roller Size & Material		$\phi 32 \times L325\text{ mm}$, PVC
Motor		AC Synchronous Motor, 12W
Timer		0~60min Analog Timer/Continuous Run
Material(Body)		Powder-coated Steel
Dimension & Weight		490 $\times$ 220 $\times$ h100mm, 3 kg
Electric Reqts.		AC 230V, 50/60Hz
Cat. No.		DH.MixR1040
Price, Each. ₩		

DAIHAN™ Digital Tube Roller Mixer "MIXR-T6", with 6 Mixing Rollers, Rolling Motion, Speed Range 10~70rpm

Ideal for Mixing Blood Sample · Viscous Substance · Liquid-Solid Suspension, 디지털 튜브 롤러믹서



3-Year Warranty

Digital Tube Roller Mixer, "MIXR-T6"

- Supplied by DAIHAN Korea ■ CE Certified

- 혈액 샘플, 점성 물질 및 액체-고체 서스펜션의 혼합과 응고/침전 방지 등에 매우 효율적인 롤링 동작을 제공.
- * Used in a Variety of Applications, such as the Prevention of Blood Coagulation, Immune Precipitation, etc.
- * LCD Displays the Speed and Time
- * Six Roller Design for Gentle Rolling and Easy to Clean.

Specifications & Ordering information

<L1>

Spec.	Model	"MIXR-T6"
No. of Mixing Roller		6 Mixing Rollers
Speed Range		10~70rpm
Mixing Roller Length		L280mm
Max. Load Capacity		4 kg
Ambient Temp & Humidity		5°C~40°C & 80%RH
Motor Type		DC Motor
Operation Mode		Continuous / Timer
Timer Setting Range		1~1199 min
Dimension & Weight		450 $\times$ 260 $\times$ h120mm, 5.1kg
Power Consumption		30W
Electric Reqts.		110V-220V, 50Hz/60Hz
Cat. No.		DH.MixR1070
Price, Each. ₩		

High Quality Rotating Electric Mini Vortex Mixer "Vortex-M", Up to 3,000rpm, Economical & Compact
With Brushless DC Motor, Resistant to Chemical Corrosion, Stable Performance, Protection Class IP43, 소형 볼텍스 믹서



NEW

1-Year Warranty

- Supplied by DAIHAN Korea
- For Mixing Solutions in Test Tubes with Oscillating Movement
- Widely Used in Biochemistry, Genetic Engineering, Medicine and Other Experimental Needs
- * Exquisite Zinc-Aluminum Base, Stable Operation
- * Can be Used for Oscillating Test Tube, Centrifugal Tube, Eppendorf Tube
- * Simple and Reliable Structure, Small Size, Low Power Consumption and Low Noise.

Specifications & Ordering information

<L1>

Spec.	Model	Vortex-M
Speed		3,000rpm
Shaking mode		Orbital Shaking, Orbital Diameter 4.5mm
Operation mode		Touch
Tube Diameter & Volume		Max. ϕ 30mm, Max. 80ml
Operating Environment		5-40°C, 80%RH
Dimension & Weight		92×92×h81mm, 0.55 kg
Power Consumption		5W
Electric Reqts.		AC100-240V, 50/60Hz
Cat. No.	Blue Silicone Head	Orange Silicone Head
	DH.Vor701	DH.Vor702
Price, Each. ₩		

DAIHAN™ Hi-Performance Multifunction Vortex Mixer-Set "VM30", Circular Orbit, Constant Torque, 3,300rpm
With ϕ 76mm(3inch) Platform Head & Pop-Off Cup Head, Variable Speed Control, Continuous or Touch-On Mixing
다용도 볼텍스 믹서, 다양한 속도 조절, 인체공학적 디자인 적용, 강력하고 안정된 동작과 저소음



3-Year Warranty

Multifunction Vortex Mixer Set "VM30"

- Supplied by DAIHAN Korea
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Mixing Solutions in Test Tubes or Small Flasks with Circular, Oscillating Movement
- Variable Speed Control Allows Slow Speed Shaking Action up to High Speed Vortexing
- * Ergonomic Design with Stable Low Profile Body.
- * A Robust and Sturdy Diecast Construction.
- * Optimized Counter Balance System Minimizes Vibration.
- * Electronic Speed Control for Constant Speed even during Changes in Load
- * One Mixing Head for a Variety of Tube Sizes(Conical Tube can be Mixed Smoothly by Hand Holding)
- * NBR(Nitrile Butadiene Rubber) Feet for Ultra Stability and Vibration Damping
- * Powerful Motor for Efficient Mixing
- * Eccentric with Oil-less Ball Bearings

Specifications & Ordering information

<L1>

Spec.	Model	“VM30”	
Speed Range		0 ~ 3,300 rpm	
Shaking Motion		Orbital, Circuit Orbit ; 4.0mm	
Control		Analog Phase Controller	
Motor		Brushless AC Motor	
Head (Included)		Pop-off Cup Head and ϕ 76mm Platform Head	
Max. Load		0.5 kg, Accepting Vessels up to 76mm Diameter	
Material		Aluminum Die-casting with Powder Coating	
Others		3-position Power Switch Operates Aluminum Die-casting with Powder Coating	
Dimensions & Weight		140×190×h123 mm, 3 kg	
Packing Dimension & Weight		175×255×h215 mm, 3 kg	
Power Consumption		30 W	
Electrical Reqts.		1 Phase, AC 120V, 50 / 60Hz or AC 230V, 50 / 60Hz	
Cat. No.	230V	60Hz	DH.WVM00030
		50Hz	DH.WVM00031
	120V	60Hz	DH.WVM00032
		50Hz	DH.WVM00033
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Spare Platform Head

- DH.WVM00110 (1) Platform Head, "PM110", ϕ 76mm(3inch)
DH.WVM00210 (2) Pop-off Cup Head, "PM210"

Ea. Price, ₩

Ea. Price, ₩



(3) Microtube Platform Head,
0.5~2.0ml Tube "PM310"
(Optional)



Remove Applicator
(Included)

(3) Optional Platform Head

- DH.WVM00310 Microtube Platform Head(0.5~2.0ml), "PM310", 12-places

DAIHAN™ Compact Functional Orbital Shaker “MIXM-1500”, with the Microplate Clamp, Max. 1,500rpm
For Holding Single or Double Microplates, Continuous Operation, 소형 기능성 궤도형 웨이커, Microplates & Tubes용



3-Year Warranty

Microplate Mixer, “MIXM-1500”
with (1) Single Microplate Clamp (Included)

※ Microplate is not Included

Optional Accessory



Fig 1. with (2) Double
Microplate Clamp



Fig 2. with
(3) Universal Top Plate
& (8) Tube Adapter



(3) Universal
Top Plate



(4) Tube Adapter,
φ6mm, 48-Hole



(5) Tube Adapter,
φ10mm, 15-Hole



(6) Tube Adapter,
φ12mm, 16-Hole



(7) Tube Adapter,
φ16mm, 8-Hole



(8) Tube Adapter,
φ20mm, 8-Hole

■ Supplied by DAIHAN Korea

■ CE Certified

■ 마이크로 플레이트용 소형믹서로 생명과학분야에서 널리 사용됨.

* Widely Used in Life Sciences Including Microbiology, Cell and Molecular Biology, Immunology, and Biotechnology.

* Compact Shape of Microplate Mixer.

* Various Mixing Applications with Optional Adapters.

* Maintenance-free Brushless DC Motor

Specifications & Ordering information

<L1>

Spec.	Model	“MIXM-1500”
Speed Range		Single Microplate ; 0~1,500rpm Double Microplate ; 0~1,000rpm
Shaking Motion		Orbital (Dia. 4.5mm)
Max. Shaking Weight		0.5kg (with Accessory)
Motor Type		Brushless DC Motor
Motor Rating (Input/Output)		18W / 10W
Ambient Temp & Humidity		5℃~40℃ & 80%RH
Dimension & Weight		260×150×h80mm & 3kg
Power Consumption		20W
Electric Reqts.		110V-220V, 50Hz/60Hz
Cat. No.		DH.MixM1500
Price, Each. ₩		

Spare Microplate Clamp

DH.Cla7001 (1) Single Microplate Clamp “PS1.2”

Optional Accessory

DH.Cla7002 (2) Double Microplate Clamp “MXM-116”

DH.Cla7003 (3) Universal Top Plate “PS1.1”, Used with Tube Adapter

DH.Ada6001 (4) Tube Adapter, for Test Tube “VT1.3.1”, φ6mm, 48-Hole

DH.Ada6002 (5) “VT1.3.2”, φ10mm, 15-Hole

DH.Ada6003 (6) “VT1.3.3”, φ12mm, 16-Hole

DH.Ada6004 (7) “VT1.3.4”, φ16mm, 8-Hole

DH.Ada6005 (8) “VT1.3.5”, φ20mm, 8-Hole

Ea. Price, ₩

DAIHAN™ Analog Orbital Shaker “SHKO-300”, 285 × 285 Platform, 10 Levels Speed Control, 40~300rpm
With Standard Rubber Mat Platform, 0~60min Timer, 궤도형 셰이커, 플랫폼 포함

40~300rpm

Included
Rubber Mat Platform


3-Year Warranty

Orbital Shaker “SHKO-300”

- Supplied by DAIHAN Korea
- Stepless Speed Controller(10Levels Speed/40~300rpm)
- Ideal for Culturing Cells to Hybridization, Staining and Destaining Gels, and Combinatorial Chemistry, etc.
- Included Rubber Mat Platform
- Optional PCR Plate- / Universal Perforated-Platform
- * Brushless DC Motor ; Low Noise and Long Life time
- * Optimal Platform Rotation for Intensive Mixing
- * Built in 0~60min Timer for User Convenience

Specifications & Ordering information

<L1>

Spec.	Model	“SHKO-300”
Platform Size		285 × 285 × h10mm
Shaking Motion		Orbital Motion
Speed Range		40~300 rpm, 10 Levels Speed Control
Motor		20W DC 12V Motor
Timer		0~60min Analog Timer/Continuous Run
Material		Powder-coated Steel
Dimension (w × d × h)		250 × 280 × h180mm
Net Weight		7 kg
Packing Size (w × d × h) & Gross Weight		320 × 350 × h250mm 7.5 kg
Electric Reqts.		AC 230V, 50/60Hz
Cat. No.		DH.ShkO1300
Price, Each. ₩		

Spare Platforms

DH.ShkO3003 Rubber Mat Platform, Multi-use, w285 × d285mm


Ea. Price, ₩

Optional Platforms and Flask Holders

DH.ShkO3001 (1) PCR Plate Platform, for PCR Plate, 285 × 285mm
DH.WSO504100 (2) Universal Platform, for Flask Holder, 230 × 230mm
DH.WIS600100 (2) STS Wire Spring Flask Holder, “FHS100”, for 100ml Flask
DH.WIS600250 , “FHS250”, for 250ml Flask
DH.WIS600500 , “FHS500”, for 500ml Flask
DH.WIS601000 , “FHS1000”, for 1000ml Flask

Ea. Price, ₩

Spec.	Model	Universal Platform
Max. Capacity of Flask Holder	For 50~100ml	9 ea
	For 200~300ml	7 ea
	For 500~1000ml	4 ea



(1) PCR Plate Platform


(2) Universal Perforated
Platform and STS Wire Spring
Flask Holder

DAIHAN™ Premium Digital Orbital and Reciprocating Shaker “SHK2.O” & “SHK2.R”, 255 × 255 Platform, 10~300 rpm
Powerful Compact Design, with Programmable Function and Large LCD Display, Quiet & Low-Vibration, without Platform
디지털 케도형 & 왕복형 셰이커, 플랫폼 별도, 저소음/저진동 긴 수명, 부드러운 시동 동작 및 속도 조절



10~300rpm

(2) Optional Rubber Mat Platform



3-Year Warranty

I. Orbital Shaker “SHK2.O”



10~300rpm

(1) Optional Universal Platform & Flask Holder



3-Year Warranty

II. Reciprocating Shaker “SHK2.R”

- Supplied by DAIHAN Korea
- CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Culturing Cells to Hybridization, Staining and Destaining Gels, Chemistry
- PL(Product Liability) Insurance
- Patented Jog-Shuttle Control System
- Available Rubber Mat-/Universal-Platform(Optional)
- * Ideal for Hybridization, Washing Blots, Culturing Microbes
- * Compact Design Easily Fits in Incubators and on the Laboratory Bench.
- * Soft Start and Noise-free Operation.
- * Digital Speed Control easily adjusts from 10 to 300 rpm
- * Built-in Digital Timer. Use in Continuous Shaking or Timer Mode up to 99 hr 59 min.
- * Gentle Speed up and Slow down by PID Speed Control. Slow Start Motion decreases Risk of Spills.
- * Large LCD for Speed, Timer Function and Operating Mode.
- * Low Vibration, Low Noise
- * Clearly laid out Control Panel for easy Operation.
- * Speed Regulation by Microprocessor.
- * Excellently Applicable for Continuous Operation.
- * Constant Speed Independent of Load.
- * Comprehensive Range of Accessories.
- * Audible Buzzer Alarms on Time Expiration and when an Error Occurs.
- * Over-speed Detection and Protection.
- * Automatic Recovery from Power Break.
- * Brushless Motor Provides Maintenance Free Operation.

Specifications & Ordering information

<L1>

Spec.	Model	I. Orbital Shaker “SHK2.O”	II. Reciprocating Shaker “SHK2.R”
Shaking Platform Size		255 × 255 mm	
Shaking Motion		Orbital Motion	Reciprocating Motion
Speed Range & Max. Load		10 ~ 300 rpm, 10kg	
Amplitude of Rotation (mm)		10 mm	
Motor		Brushless DC Motor	
Permissible Temp & Moisture		Ambient 5℃ ~ 50℃, 80%	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Large Size Digital LCD with Back-Light Function	
Controller		Digital Fuzzy Controller with Touch Button, Jog-dial	
Program Function		10 Steps, 6 Programs	
Platform		Universal Platform for various Shapes & Size of Flask	
Dimension & Weight		276 × 375 × h138 mm(Base Platform ; 255 × 255 mm), 13kg	
Packing Dimension & Weight		345 × 420 × h205 mm, 13.5kg	
Safety Circuit		Overload Protector for Motor	
Power Consumption		45 W	
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	
Cat. No.	230V	DH.WSO00300	DH.WSR00300
	120V	DH.WSO00301	DH.WSR00301
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

<Optional Accessory>



(1) Universal Perforated Platform & Flask Holder



(3) Universal Spring Rack



(4) Platform for Flasks



(5) Platform for Separatory Funnels

Optional Platforms · Flask Holders · Spring Rack

- DH.WSO604100 (1) Universal Platform, 251 × 251mm
DH.WSO604200 , 346 × 346mm
- DH.WIS600100 (1) STS Wire Spring Flask Holder, “FHS100”, for 100ml Flask
DH.WIS600250 , “FHS250”, for 250ml Flask
DH.WIS600500 , “FHS500”, for 500ml Flask
DH.WIS601000 , “FHS1000”, for 1000ml Flask

Spec.		Model	251 × 251mm Universal Platform	346 × 346mm Universal Platform
Max. Capacity of Flask Holder	For 100ml		9 ea	20 ea
	For 250ml		7 ea	15 ea
	For 500ml		5 ea	9 ea
	For 1000ml		4 ea	6 ea

- DH.WSO604110 (2) Rubber Mat Platform, 295 × 295mm
DH.WSO604210 , 370 × 370mm
- DH.WSO604140 (3) Universal Spring Rack, 295 × 295 × h100mm
DH.WSO604240 , 370 × 370 × h100mm
- DH.WSO604130 (4) Platform for Round Flask, 350 × 260 × h153mm
DH.WSO604230 (5) Platform for Separatory Funnels, 350 × 260 × h153mm

Ea. Price, ₩

DAIHAN™ Large Digital **Orbital** and **Reciprocating** Shaker “SHO-2D” & “SHR-2D”, 350×350 Platform, 10~300 rpm
With Programmable Digital Feedback Control and Large LCD Display with Back Light, Built-in Timer and Alarm
Without Platform, 대형 디지털 궤도형 & 왕복형 셰이커, 플랫폼 별도, 저소음 및 긴 수명, 부드러운 시동 동작 및 속도조절

10~300rpm
3-Year Warranty



I. Orbital Shaker “SHO-2D”

10~300rpm
3-Year Warranty



II. Reciprocating Shaker “SHR-2D”

<Optional Accessory>



(1) Universal Perforated Platform & Flask Holder



(3) Universal Spring Rack “SP240”



(4) Platform for Flasks “SP220”



(5) Platform for Separatory Funnels “SP230”

- Supplied by DAIHAN Korea
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Culturing Cells to Hybridization, Staining and Destaining Gels, and Combinatorial Chemistry.
- PL(Product Liability) Insurance
- Low Noise & Long Life Cycle

■ Patented Jog-Shuttle Control System

- * Optimal Platform Rotation for Intensive Mixing
- * High Precision Digital Rotation Speed Control
- * Digital LCD with Back-Light Function
- * Smooth Start and Change of Shaking Speed for Stable Running
- * Storage Function(Program, RPM, Timer), Alarm(Error Status and Timer-end) & Overload Protection
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Brushless DC Motor(Low Noise & Long Life Cycle)
- * Programmable Operation Mode Supported

Specifications & Ordering information

<L1>

Spec.	Model	I. Orbital Shaker “SHO-2D”	II. Reciprocating Shaker “SHR-2D”
Shaking Platform Size		350×350 mm	
Shaking Motion		Orbital Motion	Reciprocating
Speed Range		10 ~ 300 rpm	
Max. Load		10 kg	
Amplitude of Rotation (mm)		15 mm	
Motor		Brushless DC Motor	
Permissible Ambient Temp		Ambient 5℃~50℃	
Permissible Relative Moisture		85%	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Programmable Digital Feedback Controller with Jog-Shuttle Switch(Turn+Push)	
Program Function		10 Steps, 6 Programs	
Safety Circuit		Overload Protector for Motor	
Dimension (w×d×h)		400×480×h160mm	
Net Weight		20 kg	
Packing Size (w×d×h)		580×500×h250mm	
& Gross Weight		23 kg	
Power Consumption		40 W	
Electrical Reqs.		Free Voltage (1 Phase, AC 80V~240V, 50/60Hz)	
Cat. No.	230V	DH.WSO04012	DH.WSR04022
	120V	DH.WSO05012	DH.WSR05022
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Optional Platforms · Flask Holders · Spring Rack

- DH.WSO504200 (1) Universal Platform, “SP200”, 350×350mm
 DH.WIS600100 (1) STS Wire Spring Flask Holder, “FHS100”, for 100ml Flask
 DH.WIS600250 , “FHS250”, for 250ml Flask
 DH.WIS600500 , “FHS500”, for 500ml Flask
 DH.WIS601000 , “FHS1000”, for 1000ml Flask

Spec.	Model	Universal Platform
Max. Capacity of Flask Holder	For 50~100ml	20 ea
	For 200~300ml	16 ea
	For 500~1000ml	9 ea

- DH.WSO504210 (2) Rubber Mat Platform, “SP210”, 370×370mm
 DH.WSO504240 (3) Universal Spring Rack, “SP240”, 370×370×h100mm
 DH.WSO504220 (4) Platform for Flasks, “SP220”, 350×350×h153mm
 DH.WSO504230 (5) Platform for Separatory Funnels, SP230, 351×355×h174mm

Ea. Price, ₩

DAIHAN™ Compact Tube Rocker/Mixer "MIXM", Smoothly Rocking Motion, 24rpm

With Rubber Cushion Platform, Continuous Operation, Ideal for Dishes ($\geq \phi 100\text{mm}$) or 1.5~50ml Tubes, 소형 다기능 믹서/락커

3-Year Warranty

- Supplied by DAIHAN Korea
- Quality Certification ; ISO 9001.
- Ideal for Dishes ($\geq \phi 100\text{mm}$) and 1.5 ~ 50ml Tubes.
- Capacity ; 3ea - $\phi 100\text{mm}$ Dishes, 30ea - 1.5ml Microtubes, 15ea - 15ml Tubes, or 8ea - 50ml Tubes.
- Easy to load and unload the Test Tubes without Turning the Mixer Off.
- Rubber Cushion preventing the Test Tubes from Slipping.

* Widely used in Hematology, Chemistry, Research and Bacteriology Laboratories.

* Smooth, Steady, Quiet, and Continuous Operation

* Rocking Angle ; 23°

* Speed ; 24rpm at 60Hz / 20rpm at 50Hz

* Platform Dimension ; w300×d115 mm <L1>

* Overall Dimension & Weight(N.W./G.W.) ; 350×140×h110 mm, 1.5 / 2.0 kg

DH.MixM0100 Compact Mixer/Rocker, "MIXM-100", for 230V

DH.MixM0101 , "MIXM-101", for 120V

Ea. Price, ₩

DAIHAN™ Programmable Platform Rocker "RK", Digital Feedback Control, 5~50 rpm, with Certi. & Traceability

With 290×200mm or 300×300mm Platform with Non-Slip Rubber Mat, Smooth Start & Wave Motion, Tilt Angle $\pm 6^\circ / \pm 9^\circ$

디지털 락커, 디지털 피드백 컨트롤 시스템, 부드러운 시동 동작 및 속도 조절

CE 5~50rpm

3-Year Warranty



Back-light LCD

Patented Jog-Shuttle Control System



(1) 290×200 mm Rocker "RK-1D" with Included Platform

- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance

- Stainless-steel Platform Comes with Gray Nonslip Rubber Mat

- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

- High Precision Digital Speed Control ; 1rpm
- Patented Jog-Shuttle Control System

* Gel Staining, Destaining, Blotting

* Cell Culture, Hybridization, Washing

* Programmable Operation Mode Supported ; 10 Steps and 6 Programs

* Digital Timer ; Max. 99hr 59min.

* Alarm Function ; Error Status and Timer-end

* Digital LCD with Back-Light Function for Timer, RPM and Power Bar-Graph

* Smooth Start and Change of Rocking Speed

* Overload Protection Function for Motor

* Storage Function for Set Values ; RPM and Timer

* Stackable Platform with Pad to Double Mixing Area (Optional)

* Locking Mode supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

Specifications & Ordering information

<L1>

Spec.		Model	(1) "RK-1D" Platform Size 290×200 mm	(2) "RK-2D" Platform Size 300×300 mm
Dimensions (w×d×h)	Platform (mm)		290×200	300×300
	Overall (mm)		290×240×h150	300×330×h150
Speed Range (rpm)			5 ~ 50	
Angle (degree)			± 6 or ± 9 (± 9 as Factory Default)	
Motor			14W Geared DC Motor	
Program Function			10 Steps, 6 Programs	
Timer & Alarm			99hr 59min (with Continuous run.), Error status & Timer-end	
Display			Digital LCD with Back-Light Function	
Controller			Programmable Digital Feedback Controller with Jog-Shuttle Switch(Turn+Push)	
Safety Circuit			Overload Protector for motor	
Platform (Included)			290×200mm Platform	300×300mm Platform
Others			PL(Product Liability) Insurance	
			Smooth Start and Change of the Stirring Speed for Stable Running	
			Storage Function (RPM, Timer)	
			Locking Mode (Jog-Shuttle Input Disabled)	
			Alarm (Error Status and Timer-end)	
Packing Dimension (w×d×h) (mm)			450×350×250	
Gross Weight			7 kg	8 kg
Power Consumption			35 W	
Electrical Reqts.			Free Voltage (1 Phase, AC 80V~240V, 50/60Hz)	
Cat. No.	230V		DH.WRK01010	DH.WRK01020
	120V		DH.WRK02010	DH.WRK02020
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

Optional Second Platform

Ea. Price, ₩

DH.WRK05011 Second Platform, "RK5011", for "RK-1D", w290×d200, with 4 Pillars(L120mm)

DH.WRK05021 Second Platform, "RK5021", for "RK-2D", w300×d300, with 4 Pillars(L120mm)

CE

Second Platform
"RK5021"(Optional)

Pillars

3-Year Warranty

(2) 300×300 mm Rocker "RK-2D"
with Second Platform(Optional)

DAIHAN™ Programmable Digital Rotator "RT-10", 0~90° Mixing Angle, 5~60rpm, with Certi. & Traceability
With Digital Feedback Control, Adjustable Speed & Angle of Rotation, Continuous or Timed Operation
디지털 로테이터, 디지털 피드백 컨트롤 시스템, 회전각 및 회전속도 조절 가능

CE

5~60rpm



Digital Rotator "RT-10"
(Optional Disk "WRT100" Installed)



(1) Disk "WRT100"
for 8 × 50ml Tubes



(2) Disk "WRT110"
for 16 × 15ml Tubes



(3) Disk "WRT120"
for 32 × 1.5ml Tubes



(4) Disk "WRT130"
for 16 × 1.5ml Tubes



(5) Disk "WRT140"
for 4 × 250ml Flasks



(6) Knob & Post "WRT300"
(Included)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Adjustable Speed and Angle of Rotation ; 0 to 90°
- Patented Jog-Shuttle Control System
- * General Mixing, Diffusing, Dialyzing and Hybridization, Aerating Culture
- * Keeping Biological Samples in Suspension
- * Digital Feedback Control ; Steady and Constant Rotating
- * Continuous or Timed Operation
- * Digital LCD with Back-Light Function
- * Programmable Operation Mode Supported ; 10 Steps and 6 Programs
- * Overload Protection for motor
- * Smooth Start and Change of Stirring Speed
- * Storage Function for Set values of RPM and Time
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Durable Construction for Long-Term Operation
- * 5 Different Disks and 3 Roller Drums(Optional)
- * Disk for 250ml Flask is available

Specifications & Ordering information

<L1>

Spec.	Model	"RT-10"
Speed range		5~ 60 rpm
Angle range		0° ~ 90°
Motor		35W
Program Function		10 Steps, 6 Programs
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end
Display		Digital LCD with Back-Light Function
Controller		Programmable Digital Feedback Controller with Jog-Shuttle
Temp Resistance		Lower than 65℃ (Dry Air)
Safety Circuit		Overload Protector for motor
Dimensions (w × d × h) (mm)		200 × 250 × h260 mm
Weight		6 kg
Others		PL(Product Liability) Insurance Storage Function (RPM, Timer) Locking Mode (Jog-Shuttle Input Disabled) User Self-compensation Function, Alarm (Error Status and Timer-end)
Packing Size (w × d × h) (mm)		370 × 240 × h430 mm
& Gross Weight		8.5 kg
Power Consumption		40 W
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz
Cat. No.	230V	DH.WRT00010
	120V	DH.WRT00020
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

Optional Disk & Roller Drum

DH.WRT00100	(1) Disk, "WRT100", for 50ml conical tube × 8ea
DH.WRT00110	(2) , "WRT110", for 15ml conical tube × 16ea
DH.WRT00120	(3) , "WRT120", for 1.5ml conical/centrifuge Tube × 32ea
DH.WRT00130	(4) , "WRT130", for 1.5ml conical/centrifuge Tube × 60ea
DH.WRT00140	(5) , "WRT140", 250ml Flask-capacity × 4ea

Ea. Price, ₩

Spare Knob & Post

DH.WRT00300	(6) 4 × knobs & 4 × Post, "WRT300", for Installing Disk
-------------	---

Ea. Price, ₩

DAIHAN™ Digital Vibrating Sieve Shaker “SKS-208”, Acceptable ϕ 100~203mm SievesCapa-8 Sieves ϕ 200/203 $\times$ h50mm, Ideal for Sieving · Classifying · Filtering, CE Certified, 시브셰이커, 최대 8단, ϕ 200/203

CE

8-Layer
Test Sieve

3-Year Warranty

8 Sieves ϕ 200/203 $\times$ h50mm

Spare; Accessory Kit

■ Supplied by DAIHAN Korea ■ CE Certified

■ 강력한 진동으로 시료의 체 분리에 이상적이며 화학 산업 및 연구소 체 시험에 적합

■ Ideal for Effective Sieving, Classifying and Filtering Materials in the Test Sieve in Chemical Industry and Lab.

■ ϕ 100 · ϕ 200 · ϕ 203mm 시브를 수용할 수 있으며, ϕ 200(203 $\times$ h50mm) 시브는 8개까지 팬과 커버와 함께 수용됨.Max. 8-Layer Test Sieve ϕ 200 (203 $\times$ h50mm Including the Sieve Pan & Cover) can be Placed.

■ Not Including Sieves

* High Efficiency, Elaborate Design and Durable

* Any Powder or Liquid can be Sieved

* Can be Shut Down Automatically at Time.

* Can Experiment Repeatedly.

* Full Functionality and Simpler Operation.

* Timer Range ; within 0~99min. and 0~59sec

* Wooden Packing

※ Must be Tightened the Locking Nuts on both Sides to Prevent Sieve's Departure.

Specifications & Ordering information

<L1>

Spec.	Model	“SKS-208”
Acceptable Sieves		ϕ 100~203mm
Capacity		8 Sieves of ϕ 200/203 $\times$ h50mm
Motor Power		0.12(0.15 Common Motor)/kw
Noise		Less than 50dB
Amplitude		0~3mm
Vibration Frequency		1400 Vibrations/min (Approx. 23Hz)
Dimensions (w $\times$ d $\times$ h) & Weight		350 $\times$ 380 $\times$ h300mm, 35kg
Power Supply		230V, 50Hz & 60Hz
Cat. No.	60Hz	DH.SKS208
	50Hz	DH.SKS209
Price, Each. ₩		

Spare Accessory Kit

Ea. Price, ₩

DH.SKS2081 Accessory kit (2 $\times$ Rod, 1 $\times$ Rod Jig, 4 $\times$ Clamps) for SKS-208

CE



(1)



(2)



(3)

CISA® Premium Electromagnetic Sieve Shaker “BA200N” & “RP200N”, for Sieve ϕ 60~203mm
Digital 3D Movement, Dry & Wet Sieving, Noiseless, 3차원 운동의 전자기 시브 셰이커

■ Supplied by CISA® in Spain

■ 강력한 회전운동 및 수직운동의 3D 모션으로 체분석에 이상적

■ CE Certified ■ Traceability and Certificate ; Controlled by Serial No.

■ 3-D ; Vertical · Amplitude · Roating Powerful Movement and Noiseless, Easy Sieving Process.

* The Latest Technology with Intermittent Sieving

* 3-Dimensional Movement & Digital Timer with Countdown

* Includes the System “CISA-PRESS” Closure with a Closing System of the Column of Sieves with Fixtures Bars, Stainless-steel

* Suitable for Both Dry and Wet Sieving, by Adapting the Special Kit. (Optional Wet Sieving Kit)

* No Maintenance ; Strong, Reliable and Polyvalent * Four Programs, For Different Kinds of Sieving

* Power Supply ; 230V/60Hz (230V/50Hz, 120V/50Hz or 60Hz on Request)

* Wooden Packing

* Only Sieves of the Same Diameter are Stackable

(1) Digital Sieve Shaker “BA200N”, Max. Amplitude 2.0mm

* Recommended for Materials with Density > 0.8gr/cm³

* Capacity ; 1) 8 sieves of ϕ 203 $\times$ h50mm, 2) 15 sieves of ϕ 203 $\times$ h25mm

<L1>

Item No.	Diameters stackable Sieves	Maximum Amplitude	Maximum Sample	Dimensions	Net Weight	Ea. Price, ₩
CE.BA200N	ϕ 60~203mm	2.0mm	3kg	400 $\times$ 310 $\times$ h130mm	27kg	

(2) Digital Sieve Shaker “RP200N”, Max. Amplitude 2.5mm

* Recommended for Chemical & Pharmaceutical Laboratories and all the Fields that Require the Granulometric Test of Precision

* Noiseless, It doesn't Transmit Vibration to the Environment

* Capacity ; 1) 8 sieves of ϕ 203 $\times$ h50mm, 2) 15 sieves of ϕ 203 $\times$ h25mm

Item No.	Diameters stackable Sieves	Maximum Amplitude	Maximum Sample	Dimensions	Net Weight	Ea. Price, ₩
CE.RP200N	ϕ 60~203mm	2.5mm	3kg	455 $\times$ 355 $\times$ h155mm	45kg	

Accessory for Sieve Shaker

Ea. Price, ₩

CE.FIT.KIT (3) Spare ; Top Sieving Clamp, Stainless-steel

CE.1000003 Option ; Wet sieving Kit(cover/pan) for ϕ 200/203mm

DAIHAN™ ϕ 200 × h50mm All Stainless-steel Standard Test Sieve, with Stainless-steel Mesh-Holes(■)

표준망체, 스테인레스 메쉬홀(■), ASTM/ISO 규격에 준함.



Dia ϕ 200 × h50mm Sieve



Cover & Pan Set

■ Supplied by DAIHAN Korea

■ Made of Stainless-steel #304

■ 다양한 크기의 구멍을 가진 시브를 조합하여 시브셰이커를 활용하여 다양한 크기의 입자를 분류 및 분석에 이상적

* According to ASTM and ISO Specification

* Standard-type for Laboratory and Industry

* High Degree of Corrosion Resistance and Easy Cleaning

* ϕ 200 Sieve & ϕ 203 Sieve are not Stackable

(1) With 6.30 ~ 125.00mm Large Opening

<L1>

Cat. No	ASTM	Sieve Opening	Price, ₩
DH.Sie7011	5"	125.00 mm	
DH.Sie7012	4.24"	106.00 mm	
DH.Sie7014	4"	100.00 mm	
DH.Sie7015	3"	75.00 mm	
DH.Sie7016	2.5"	63.00 mm	
DH.Sie7017	2.12"	53.00 mm	
DH.Sie7019	1.75"	45.00 mm	
DH.Sie7020	1.5"	37.50 mm	
DH.Sie7021	1.25"	31.50 mm	
DH.Sie7022	1.06"	26.50 mm	

Cat. No	ASTM	Sieve Opening	Price, ₩
DH.Sie7023	1.00"	25.00 mm	
DH.Sie7025	0.75"	19.00 mm	
DH.Sie7026	0.625"	16.00 mm	
DH.Sie7027	0.530"	13.20 mm	
DH.Sie7028	0.5"	12.50 mm	
DH.Sie7029	0.438"	11.20 mm	
DH.Sie7030	0.375"	9.50 mm	
DH.Sie7031	0.313"	8.00 mm	
DH.Sie7032	0.265"	6.70 mm	
DH.Sie7033	0.25"	6.30 mm	

(2) With 20 μ m ~ 5.60mm Micro-Small Opening

Cat. No	ASTM	Sieve Opening	Price, ₩
DH.Sie7034	No. 3.5	5.60 mm	
DH.Sie7035		5.00 mm	
DH.Sie7036	No. 4	4.75 mm	
DH.Sie7038	No. 5	4.00 mm	
DH.Sie7040	No. 6	3.35 mm	
DH.Sie7041		3.15 mm	
DH.Sie7042	No. 7	2.80 mm	
DH.Sie7043		2.50 mm	
DH.Sie7044	No. 8	2.36 mm	
DH.Sie7045		2.24 mm	
DH.Sie7046	No.10	2.00 mm	
DH.Sie7047		1.80 mm	
DH.Sie7048	No.12	1.70 mm	
DH.Sie7049		1.60 mm	
DH.Sie7050	No.14	1.40 mm	
DH.Sie7051		1.25 mm	
DH.Sie7052	No. 16	1.18 mm	
DH.Sie7054	No.18	1.00 mm	
DH.Sie7055		900 μ m	
DH.Sie7056	No.20	850 μ m	
DH.Sie7057		800 μ m	
DH.Sie7058	No.25	710 μ m	
DH.Sie7059		630 μ m	
DH.Sie7060	No.30	600 μ m	
DH.Sie7061		560 μ m	
DH.Sie7062	No.35	500 μ m	
DH.Sie7063		450 μ m	
DH.Sie7064	No.40	425 μ m	
DH.Sie7065		400 μ m	
DH.Sie7066	No.45	355 μ m	

Cat. No	ASTM	Sieve Opening	Price, ₩
DH.Sie7067		315 μ m	
DH.Sie7068	No.50	300 μ m	
DH.Sie7069		280 μ m	
DH.Sie7070	No.60	250 μ m	
DH.Sie7071		224 μ m	
DH.Sie7072	No.70	212 μ m	
DH.Sie7073		200 μ m	
DH.Sie7074	No.80	180 μ m	
DH.Sie7075		160 μ m	
DH.Sie7076	No.100	150 μ m	
DH.Sie7077		140 μ m	
DH.Sie7078	No.120	125 μ m	
DH.Sie7079		112 μ m	
DH.Sie7080	No.140	106 μ m	
DH.Sie7081		100 μ m	
DH.Sie7082	No.170	90 μ m	
DH.Sie7083		80 μ m	
DH.Sie7084	No.200	75 μ m	
DH.Sie7085		71 μ m	
DH.Sie7086	No.230	63 μ m	
DH.Sie7087		56 μ m	
DH.Sie7088	No.270	53 μ m	
DH.Sie7089		50 μ m	
DH.Sie7090	No.325	45 μ m	
DH.Sie7091		40 μ m	
DH.Sie7092	No.400	38 μ m	
DH.Sie7093		36 μ m	
DH.Sie7094	No.450	32 μ m	
DH.Sie7095	No.500	25 μ m	
DH.Sie7096	No.635	20 μ m	

Optional Cover & Pan Set

DH.Sie7100

Sieve-Cover & Pan Set

Ea. Price, ₩

DAIHAN™ $\phi 203 \times h41$ mm Standard Test Sieve, with Stainless-steel Mesh-Holes(■) & Brass-Frame

표준망체, 스테인레스 메쉬홀(■) 황동틀, KS/ASTM/ISO 규격에 준함.



Dia $\phi 203 \times h41$ mm Sieve



Cover & Pan Set

■ Supplied by DAIHAN Korea

- * Consist of Brass-frame with Stainless-steel mesh
- * Dia $\phi 203 \times h41$ mm for Domestic Business KS.A
- * $\phi 200$ Sieve & $\phi 203$ Sieve are not Stackable

* According to KS.A, ASTM and ISO 3310 Specification

* Standard-type for Laboratory and Industry

* Other sizes available ; On request

※ Note ; Only Sieving mesh is Fulfilled all of KS, ASTM and ISO Specification

$\phi 203 \times h41$ mm Standard Test Sieves

* Same Specifications ASTM & ISO standard, but Sieve-height is different.

<L2>

Cat. No.	ASTM No.	Sieve Opening Size	Ea. Price, ₩
DH.Si8001	5"	125 mm	
DH.Si8002	4"	106 mm	
DH.Si8003	3.5"	90 mm	
DH.Si8004	3"	75 mm	
DH.Si8005	2.5"	63 mm	
DH.Si8006	2"	53 mm	
DH.Si8007	1.75"	45 mm	
DH.Si8008	1.5"	37.5 mm	
DH.Si8009	1.25"	31.5 mm	
DH.Si8010	1.06"	26.5 mm	
DH.Si8011	0.875"	22.4 mm	
DH.Si8012	0.75"	19 mm	
DH.Si8013	0.625"	16 mm	
DH.Si8014	0.5"	13.2 mm	
DH.Si8015	0.438"	11.2 mm	
DH.Si8016	0.375"	9.5 mm	
DH.Si8017	0.313"	8 mm	
DH.Si8018	0.265"	6.7 mm	
DH.Si8019	0.25"	6.35 mm	
DH.Si8020	No.3.5	5.60 mm	
DH.Si8021	-	5.00 mm	
DH.Si8022	No. 4	4.75 mm	
DH.Si8023	No. 5	4 mm	
DH.Si8024	No. 6	3.35 mm	
DH.Si8025	No. 7	2.80 mm	
DH.Si8026	No. 8	2.36 mm	
DH.Si8027	No. 10	2.00 mm	
DH.Si8028	No. 12	1.70 mm	
DH.Si8029	No. 14	1.40 mm	
DH.Si8030	No. 16	1.18 mm	
DH.Si8031	No. 18	1.00 mm	
DH.Si8032	No. 20	850 μ m	
DH.Si8033	No. 25	710 μ m	
DH.Si8034	No. 30	600 μ m	
DH.Si8035	No. 35	500 μ m	
DH.Si8036	No. 40	425 μ m	
DH.Si8037	No. 45	355 μ m	
DH.Si8038	No. 50	300 μ m	
DH.Si8039	No. 60	250 μ m	
DH.Si8040	No. 70	212 μ m	
DH.Si8041	No. 80	200 μ m	
DH.Si8042	-	180 μ m	
DH.Si8043	No.100	150 μ m	
DH.Si8044	No.120	125 μ m	
DH.Si8045	No.140	106 μ m	
DH.Si8046	-	100 μ m	
DH.Si8047	No.170	90 μ m	
DH.Si8048	No.200	75 μ m	
DH.Si8049	No.230	63 μ m	
DH.Si8050	No.270	53 μ m	
DH.Si8051	No.325	45 μ m	
DH.Si8052	No.400	38 μ m	
DH.Si8053	No.450	32 μ m	
DH.Si8054	No.500	25 μ m	
DH.Si8055	No.635	20 μ m	

Optional Cover & Pan Set

DH.Si8100

Sieve-Cover & Pan Set

Ea. Price, ₩

DAIHAN Scientific **STERILIZERS**

Autoclaves , Bench-top Front Door.....	P. 299
, Front Door	P. 304
, Top-loading Vertical New	P. 300~303
Dry Heat Sterilizers	P. 305
Plasma Sterilizers	P. 306
Accessories	P. 303, 307



Autoclave, Bench-top Front Door



Autoclave, Top-loading / Vertical



Dry Heat Sterilizer



Plasma Sterilizer

DAIHAN™ Premium 25 & 40 Lit Programmable Bench-top Front Door Steam Sterilizer "Ste.B", Medicaluse

With Auto-Front Door, Pre- & Post-Vacuum System, Class-B Sterilization, Auto-Microprocessor Control, and Dual Door Lock, 110°C ~ 135°C SQUARE Chamber Autoclave, 2 Perforated Stainless-steel Tray, HEPA Filter, Max 2.30 bar, 사각챔버 벤치탑 고압 스팀 멸균기/오토클레이브

110°C ~ 135°C
0.4 ~ 2.3 bar

Medicaluse



3-Year Warranty

Thermal Printer
(Optional)25 & 40 Lit Bench-top Steam Sterilizer
"Ste.B"

MODE	Temperature (°C)	Pressure (bar)	Sterilizing Time (min)	Drying Time (min)
Unwrapped 1	121	1.1	20	30
Unwrapped 2	134	2.1	5	30
Wrapped 1	121	1.1	30	30
Wrapped 2	134	2.1	15	30
Liquid	121	1.1	20	0
Bowie Dick	134	2.1	4	5
User Program	110 - 135	0.4 - 2.3	4 - 30	0 - 60

Program Cycle



Square Inner Chamber



Sliding Perforated Shelf



(1) STS Wire Tray



(2) STS Removable Rack



(3) HEPA Filter-kit, 0.3µm



(6) Thermal Printer



(5) Movable Stand

- Supplied by DAIHAN Korea
- 치료, 병원, 생명공학용에 적합한 멸균기
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.
- Class-B Sterilization Cycle by EN13060-2 Standard ; Pre Vacuum Air Removal → Water Filling → Heating → Sterilizing → Ventilation → Post Vacuum Drying → End
- STS Square Chamber for Corrosion Resistance and Maximum Working Space
- * Ideal for Dental Clinic, Hospital, Biotechnology, etc.
- * Full Automatic Process Cycle by High-performance Microprocessor Controller
- * Built-in Separate Water Reservoir Supplies Water Automatically
- * Self-diagnostic and Safety System ;
 - Display Error Codes and Audible Alarm
 - Over Pressure Release Valve
 - Automatic and Manual Door Lock System
 - Over Temp Protector with Thermostat
- * Intuitive User Interface(LCD) & Indicator of Working(LED Lamp)
- ※ Note! When Autoclaving, Loosen Off the PP Screwcap!!

Specifications & Ordering information

<M2>

Spec.	Model	"Ste25.B"	"Ste40.B"
Capacity		25 Lit	40 Lit
Water Filling Tank		5 Lit (Reservoir)	
Sterilizing Cycle-type		Class-B Cycle by EN13060-2 Standard	
Temp Range		110°C ~ 135°C	
Pressure & Gauge		0.4 ~ 2.3 bar, Built-in Analog & Digital Pressure Gauge	
Timer & Alarm		Error Status & Timer-end	
Controller		Microprocessor PID Control System(Fully Automatic) with LCD Monitor	
Programmable Function		Pre-set 6 Programs (Customizable 5- & Fixed Bowie Dick 1- Program)	
Door Lock		Automatic Wheel Inter Locking with One Touch Sliding Handle	
Air Filtration System		Built-in HEPA Filter, 0.3µm Particle Removal 99.97%	
Vacuum Pump		Built-in Vacuum Pump(33LPM, 580mmHg) for Pre-vacuum & Post-vacuum Drying	
Air Removal-type		Vacuum	
Tray (Included)		2 Stainless-steel Perforated Tray	
Safety Circuit		Over-Temp / -Pressure / -Current protection, Sensor Error Detection	
Printer (Optional)		Thermal printer with 1 Record Paper Roll; w58 × φ40mm 1000 Sterilization Cycle Data / roll (Date, Temperature, Pressure, Time, Error Status, etc.)	
Material		Interior ; Stainless-steel 304 Polished, Exterior ; Epoxy Powder Coated Steel	
Dimension (w×d×h) & Weight	Interior	270 × 420 × h230mm	300 × 480 × h300mm
	Exterior	495 × 689 × h460mm, 60kg	527 × 787 × h530mm, 70kg
	Packing	580 × 750 × h580mm, 75kg	650 × 820 × h740mm, 85kg
Power Consumption		2.5kW	
Electric Reqs.		220V, 50/60 Hz, 1 Phase	
Cat. No.	220V	DH.SteB8001	DH.SteB8002
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Accessory

- DH.SteBA1013 (1) STS Wire Tray, 260 × 370mm, for "Ste25.B"
 DH.SteBA1017 , 273 × 410mm, for "Ste40.B"
 DH.SteBA1011 (2) STS Removable Rack, 260 × 385 × h155mm, for "Ste25.B"
 DH.SteBA1016 , 290 × 425 × h185mm, for "Ste40.B"
 DH.SteA1015 (3) Replacement HEPA Filter-kit, 0.3µm


Optional Accessory

- DH.SteBA1010 (4) STS Perforated Tray, 255 × 365 × h25mm, for "Ste25.B"
 DH.SteBA1015 , 285 × 405 × h30mm, for "Ste40.B"
 DH.SteBA1026 (5) Movable Stand with 1 Drawer, 580 × 700 × h850mm
 DH.SteA1005 (6) Thermal Printer, Paper Roll, "STE1005"
 DH.SteA1006 Thermal Paper Package, 4 Paper Roll, for "STE1005"

Bench-top Autoclave with Validation (IQ, OQ) Service

- DH.SteB5025 Validation Service(IQ, OQ), with "Ste25.B"
 DH.SteB5040 , with "Ste40.B"

DAIHAN™ **Premium** Top-Door Vertical Autoclave, **Medicaluse** “STE-AM”, With Drying Function & Recorder, Max 135 °C With 7” Touch Board Controller, ϕ 3.0mm Thick Tank, Max.2 kgf/cm², Electronic Door Lock System, Fast Cooling, 47 · 59 · 79 · 100Lit 6 Sterilizing Modes(Liquid, Solid, Liquid Warm, Unsealed Package, Agar Melting and Custom), Compact Footprint
의료용 탑도어 버티칼 고압 스팀 멸균기 / 오토클레이브, 선 · 후 진공 & 건조기능, 전자식 도어 잠금 시스템, 7인치 터치보드 조절기



Medicaluse, “STE-AM59”
with Recorder & Drying Vacuum Pump



※ For the Baskets & Recorder
⇒ See 3 Pages Later !

- Supplied by DAIHAN Korea ■ CE and PED Certified ■ PL(Product Liability) Insurance
- ▶ **Medicaluse Class** Sterilizing Process ;
Pre-Vacuum → Heating & Pressurizing → Sterilization → Post-Vacuum & Drying.
with Built in Powerful Vacuum Pump and Dot-type Recorder
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ergonomic 7” Full Touch Screen TFT LCD ; Easy Monitoring and Control.
- Optimized Inner Chamber Space for Maximum Loading Capacity
- Adjustable Exhaust Speed and Cooling Speed for Various Application
- * Ideal for Medical, Clinical, Biotechnology, Environmental, Pharmacy, Food Industry, etc.
- * Water Low Level Detection * 4 Caters for Easy Mobility
- * Automatic Over Pressure Release Valve ensures User Safety
- * Compact Design for Efficient Space Utilization

Specifications & Ordering information

<M1>

Spec.	Model	“STE-AM47”	“STE-AM59”	“STE-AM79”	“STE-AM100”
Capacity		47Lit	59Lit	79Lit	100Lit
Temp Range		45 °C ~ 135 °C			
Pressure Gauge Range		1.1bar at 121 °C, 2.1bar at 134 °C (Pressure Gauge 0~4 bar)			
Recorder		Built-in Dot-type Recorder with 1 Record Paper Roll ; w57 × ϕ 35mm 1000 Sterilization Cycle / Roll (Date, Temp, Pressure, Time, Error Status)			
Vacuum Pump for Drying		Max Vacuum 240mbar, 35L/min			
6 Sterilization Modes		5 Pre-Set & 1 Custom Mode ; Optimized for Every Application ① Liquid Mode ; 121 °C, 20min, Exhaust Speed 0, Cooling Fan OFF ② Solid Mode ; 134 °C, 10min, Exhaust Speed 3, Cooling Fan ON ③ Liquid Warm Mode ; 121 °C, 20min, Warming(50 °C/60min), Exhaust Speed 0, Cooling Fan OFF ④ Unsealed Package Mode ; 121 °C, 30min, Exhaust Speed 2, Cooling Fan ON ⑤ Agar Melting Mode ; 100 °C, 10min, Warming(50 °C/60min), Exhaust Speed 3, Cooling Fan OFF ⑥ Custom Mode ; Fully Programmable Parameters			
Timer & Alarm		9hr 59min(Sterilizing & Heating), 99hr 59min(Warming), Error Status & Timer End Alarm			
Controller		7” Full Touch Screen TFT LCD			
Door Lock System		Electronic Door Lock System			
Safety Circuit		Over-Temp/-Pressure/-Current Protection, Sensor Error Detection, Water Level Sensor			
Material		Pressure Vessel ; Stainless-steel 304, External ; Powder Coated steel, Door Cover ; ABS			
Basket (Included)		2 Stainless-steel Wire Basket			
Basket Dimension (mm)		ϕ 330 × h195	ϕ 330 × h285	ϕ 430 × h215	ϕ 430 × h285
Dimension (w × d × h)	Inner (mm)	ϕ 350 × h424	ϕ 350 × h622	ϕ 450 × h470	ϕ 450 × h621
	Outer (mm)	513 × 692 × h1,065		613 × 792 × h1,065	
	Packing (mm)	650 × 730 × h1,250		750 × 940 × h1,250	
Net Weight (kg)		87kg	90kg	125kg	135kg
Gross Weight (kg)		97kg	100kg	140kg	150kg
Electric Reqs.		1Phase - AC 120V, 60Hz or - AC 230V, 50/60Hz		1 Phase - AC 120V, 60Hz or - AC 230V, 50/60Hz 3 Phase with Terminal Block - AC 380V, 60Hz	

Cat. No.	230V	DH.SteAM47	DH.SteAM59	DH.SteAM79	DH.SteAM100
	120V	DH.SteAM47.1	DH.SteAM59.1	DH.SteAM79.1	DH.SteAM100.1
	3Phase	-	-	DH.SteAM79.2	DH.SteAM100.2
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

※ 운임 별도

DAIHAN™ Premium Top-Door Vertical Autoclave, **I. With Drying Function “STE-ALD” & II. Without Drying Function “STE-ALN”**, 135 °C 7" Touch Board Controller, ϕ 3.0mm Thick Tank, Max.2 kgf/cm², Electronic Door Lock System, Fast Cooling, 47 · 59 · 79 · 100Lit 6 Sterilizing Modes(Liquid, Solid, Liquid Warm, Unsealed Package, Agar Melting and Custom), Compact Footprint
 탑도어 버티칼 고압 스팀 멸균기 / 오토클레이브, 선 · 후 진공 & 건조기능(STE-ALD Models), 전자식 도어 잠금 시스템, 7인치 터치보드 조절기

NEW

CE



**I. With Drying Function, “STE-AL59D”
with Vacuum Pump**

CE



Figure ; Closed Door 3-Year Warranty

**II. Without Drying Function, “STE-AL59N”
without Vacuum Pump**

Supplied by DAIHAN Korea ■ CE and PED Certified ■ PL (Product Liability) Insurance

I. With Drying Function, with Vacuum Pump

; Pre-Vacuum → Heating & Pressurizing → Sterilization → Post-Vacuum & Drying.

II. Without Drying Function

; Heating & Pressurizing → Sterilization → Exhaust

■ with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

■ Ergonomic 7" Full Touch Screen TFT LCD ; Easy Monitoring and Control.

■ Optimized Inner Chamber Space for Maximum Loading Capacity

■ Adjustable Exhaust Speed and Cooling Speed for Various Application

* Ideal for Clinical, Biotechnology, Environmental, Pharmacy, Food Industry, etc.

* Water Low Level Detection

* 4 Caters for Easy Mobility

* Automatic Over Pressure Release Valve ensures User Safety

* Compact Design for Efficient Space Utilization

Specifications & Ordering information

※ 운임 별도 <M1>

I. With Drying Function	Model	“STE-AL47D”	“STE-AL59D”	“STE-AL79D”	“STE-AL100D”
	Drying Vacuum Pump	Max Vacuum 240mbar, 35L/min			
II. Without Drying Function	Model	“STE-AL47N”	“STE-AL59N”	“STE-AL79N”	“STE-AL100N”
Capacity		47Lit	59Lit	79Lit	100Lit
Temp Range		45°C ~ 135°C			
Pressure Gauge Range		1.1bar at 121°C, 2.1bar at 134°C (Pressure Gauge 0~4 bar)			
6 Sterilization Modes		5 Pre-Set & 1 Custom Mode ; Optimized for Every Application ① Liquid Mode ; 121°C, 20min, Exhaust Speed 0, Cooling Fan OFF ② Solid Mode ; 134°C, 10min, Exhaust Speed 3, Cooling Fan ON ③ Liquid Warm Mode ; 121°C, 20min, Warming(50°C/60min), Exhaust Speed 0, Cooling Fan OFF ④ Unsealed Package Mode ; 121°C, 30min, Exhaust Speed 2, Cooling Fan ON ⑤ Agar Melting Mode ; 100°C, 10min, Warming(50°C/60min), Exhaust Speed 3, Cooling Fan OFF ⑥ Custom Mode ; Fully Programmable Parameters			
Timer & Alarm		9hr 59min(Sterilizing & Heating), 99hr 59min(Warming), Error Status & Timer End Alarm			
Controller		7" Full Touch Screen TFT LCD			
Door Lock System		Electronic Door Lock System			
Safety Circuit		Over-Temp/-Pressure/-Current Protection, Sensor Error Detection, Water Level Sensor			
Material		Pressure Vessel ; Stainless-steel 304, External ; Powder Coated steel, Door Cover ; ABS			
Basket (Included)		2 Stainless-steel Wire Basket			
Basket Dimension (mm)		ϕ 330 × h195	ϕ 330 × h285	ϕ 430 × h215	ϕ 430 × h285
Dimension (w×d×h)	Inner (mm)	ϕ 350 × h424	ϕ 350 × h622	ϕ 450 × h470	ϕ 450 × h621
	Outer (mm)	513 × 692 × h1,065		613 × 792 × h1,065	
	Packing (mm)	650 × 730 × h1,250		750 × 940 × h1,250	
Net Weight (kg)		87kg	90kg	125kg	135kg
Gross Weight (kg)		97kg	100kg	140kg	150kg
Electric Reqts.		1Phase - AC 120V, 60Hz or - AC 230V, 50/60Hz		1 Phase - AC 120V, 60Hz or - AC 230V, 50/60Hz 3 Phase with Terminal Block - AC 380V, 60Hz	

I. With Drying Function Cat. No.	230V	DH.SteAL47D	DH.SteAL59D	DH.SteAL79D	DH.SteAL100D
	120V	DH.SteAL47D.1	DH.SteAL59D.1	DH.SteAL79D.1	DH.SteAL100D.1
	3Phase	-	-	DH.SteAL79D.2	DH.SteAL100D.2

Price, Each. ₩

II. Without Drying Function Cat. No.	230V	DH.SteAL47N	DH.SteAL59N	DH.SteAL79N	DH.SteAL100N
	120V	DH.SteAL47N.1	DH.SteAL59N.1	DH.SteAL79N.1	DH.SteAL100N.1
	3Phase	-	-	DH.SteAL79N.2	DH.SteAL100N.2

Price, Each. ₩

※ For the Baskets & Recorder
⇒ See 2 Pages Later !

DAIHAN™ Digital Fuzzy-Control Autoclave “MaXterile™”, **I. Recorder-type(Built-in Recorder)** & **II. Standard-type(w/o Recorder)**
With PED Certified ϕ 3.0mm Thick-Tank, Max.2 kgf/cm², Electronic Door Lock System, 47 · 60 · 80 · 100 Lit, Up to 132 °C
Lever-type Sliding Door, Steam Condensing, Solid-/Liquid-Models, 전자동 고압 스팀 멸균기/오토클레이브, 전자식 도어 잠금 시스템



I. Recorder-type, “MaXterile 80R”
With Built-in Chart Recorder



II. Standard-type, “MaXterile 80”
Without Recorder

※ For the Baskets ⇒ See Next Page !

- Supplied by DAIHAN Korea ■ CE and PED Certified ■ PL(Product Liability) Insurance
- With Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-Information, and Traceable Data Base System
- Patented Special Steam Condensing Mechanism, Lever-Type Sliding Door
- Built in New Generation Controller ; Patented Jog-Dial Control System
- Electronic Door Lock Sensing System for User Safety ; Door Open Button Automatically Lock Over 80 °C for User Safety
- **I. Recorder-Type Model (Only MaXterile 47R · 60R · 80R · 100R Models) ; Print Sterilization Cycle Data for Validation and Analyze Sterilization (Date, Temperature, Pressure, Time, Error Status, etc.)**
- * Ideal for Biotechnology, Clinical, Environmental, Medical and Food Industry.
- * Digital Timer Function - 'Delay Time' and 'Operation Time' can be Controlled.
- * Water Low Level Detection
- * Automatic Over Pressure Release Valve Ensures User Safety.
- * Sterilization Error Indication ; Error Alarm and Message in Case of Power Failure and Temperature Error
- ※ **Note/ When Autoclaving, Loosen Off the PP Screwcap!!**

Specifications & Ordering information

<M1>

	“MaXterile 47R” Max. 47Lit	“MaXterile 60R” Max. 60Lit	“MaXterile 80R” Max. 80Lit	“MaXterile 100R” Max. 100Lit
I. Recorder-type (Built-in Recorder)	Thermal Print Recorder with 1 Record Paper Roll ; w58 × ϕ 40mm 1000 Sterilization Cycle Data / Roll (Date, Temperature, Pressure, Time, Error Status, etc.)			
II. Standard-type (Without Recorder)	“MaXterile 47”	“MaXterile 60”	“MaXterile 80”	“MaXterile 100”
Capacity	47Lit	60Lit	80Lit	100Lit
Temp Range & Accuracy	Room Temp +5 °C ~ 132 °C, $\pm$ 0.5 °C at 121 °C			
Pressure & Gauge	1.2 kgf/cm ² at 121 °C / 2.0 kgf/cm ² at 132 °C, 0~3 kgf/cm ²			
Sterilization Mode	Solid- or Liquid-mode, Chosen by a Customer			
Heating Power/Consumption	2kW	3kW	2kW × 2	2kW × 2
Timer & Alarm	60 min, Error status & Timer-end			
Display	Digital LCD with Back-Light Function			
Controller	Digital Fuzzy-Control System (Fully Automatic)			
Door & Closing System	Lever-Type Sliding Door / Safety Electronic Door Closing System			
Safety Circuit	Over-Temp/-Pressure/-Current Protection, Sensor Error Detection			
Material	Pressure Vessel-Stainless-steel(# 304), External-Powder Coated steel			
Basket (Included)	2 Stainless-steel Wire Basket			
Basket Dimension (mm)	ϕ 270 × h240	ϕ 320 × h240	ϕ 370 × h240	ϕ 420 × h240
Dimension	Inner (mm) Outer (mm)	ϕ 300 × h670 700 × 470 × 1080	ϕ 350 × h650 700 × 470 × 1080	ϕ 400 × h650 810 × 580 × 1140
Net Weight (kg)	107kg	110kg	155kg	155kg
Packing Size (mm)	840 × 610 × h1,252mm		950 × 730 × h1,352mm	
Gross Weight (kg)	117kg	120kg	170kg	170kg
Electric Reqts.	1Phase - AC 120V, 60Hz or - AC 230V, 50/60Hz		1 Phase - AC 120V, 60Hz or - AC 230V, 50/60Hz 3 Phase with Terminal Block	
I.Recorder-type	230V	DH.WACR2047P	DH.WACR2060P	DH.WACR2080P
Cat. No.	120V	DH.WACR1047P	DH.WACR1060P	DH.WACR1080P
	3Phase	-	-	DH.WACR3080P
Price, Each. W				
II. Standard-type	230V	DH.WAC05047P	DH.WAC05060P	DH.WAC05080P
Cat. No.	120V	DH.WAC06047P	DH.WAC06060P	DH.WAC06080P
	3Phase	-	-	DH.WAC07080P
Price, Each. W				

* Other Specifications are available upon Customer's Request.

※ 운임 별도

(1) Spare Thermal Paper Roll for Recorder

Ea. Price, W

DH.WACR5004 Thermal Paper for Recording Printer, 4 × Paper Roll(w58 × ϕ 40mm) |
for 1000 Sterilization Cycle Data / Roll

Articles		Cat. No	Description	
I. Accessories for “STE-AM” · “STE-ALD” · “STE-ALN” Autoclaves				
				<L1>
(1) Spare Wire Basket, Stainless-steel				Ea. Price, ₩
 (1) Wire Basket Stainless-steel (Included)	DH.Bas71A	Wire Basket, $\phi 330 \times h195\text{mm}$, “Bas71A”, for 47 Lit		
	DH.Bas72A	, $\phi 330 \times h285\text{mm}$, “Bas72A”, for 59 Lit		
	DH.Bas73A	, $\phi 430 \times h215\text{mm}$, “Bas73A”, for 79 Lit		
	DH.Bas74A	, $\phi 430 \times h280\text{mm}$, “Bas74A”, for 100 Lit		
				<L1>
(2) Spare Recorder & Accessories				Ea. Price, ₩
 (2) Dot-type Recorder	DH.AR1	Dot-type Recorder 1 Paper Roll for 1000 Sterilization Cycle Data		
	DH.AR2	Ribbon Cartridge for Dot-type Recorder		
	DH.AR3	Paper Roll for Dot-type Recorder, w57mm $\times$ ϕ 35mm		
				Ea. Price, ₩
(3) Optional Perforated Bottom Basket, Stainless-steel				Ea. Price, ₩
 (3) Perforated Bottom Basket,Stainless-steel (Optional)	DH.Bas71AP	Basket, $\phi 330 \times h195\text{mm}$, “Bas71AP”, for 47 Lit		
	DH.Bas72AP	, $\phi 330 \times h285\text{mm}$, “Bas72AP”, for 59 Lit		
	DH.Bas73AP	, $\phi 430 \times h215\text{mm}$, “Bas73AP”, for 79 Lit		
	DH.Bas74AP	, $\phi 430 \times h280\text{mm}$, “Bas74AP”, for 100 Lit		
II. Accessories for “Maxterile™” Autoclaves				
				<L1>
(4) Spare Wire Basket, Stainless-steel				Ea. Price, ₩
 (4) Wire Basket Stainless-steel (Included)	SL.Bas7011	Wire Basket, $\phi 270 \times h240\text{mm}$, “Bas7011”, for 47 Lit , “MaXterile 47 & 47R”		
	SL.Bas7012	, $\phi 320 \times h240\text{mm}$, “Bas7012”, for 60 Lit , “MaXterile 60 & 60R”		
	SL.Bas7013	, $\phi 370 \times h240\text{mm}$, “Bas7013”, for 80 Lit , “MaXterile 80 & 80R”		
	SL.Bas7014	, $\phi 420 \times h240\text{mm}$, “Bas7014”, for 100 Lit, “MaXterile 100 & 100R”		
				<L1>
(5) Optional Perforated Bottom Basket, Stainless-steel				Ea. Price, ₩
 (5) Perforated Bottom Basket,Stainless-steel (Optional)	DH.WAC20047	Basket, $\phi 270 \times h250\text{mm}$, “WAC20047”, for 47 Lit , “MaXterile 47 & 47R”		
	DH.WAC20060	, $\phi 320 \times h250\text{mm}$, “WAC20060”, for 60 Lit , “MaXterile 60 & 60R”		
	DH.WAC20080	, $\phi 370 \times h250\text{mm}$, “WAC20080”, for 80 Lit , “MaXterile 80 & 80R”		
	DH.WAC20100	, $\phi 420 \times h250\text{mm}$, “WAC20100”, for 100 Lit, “MaXterile 100 & 100R”		
				Ea. Price, ₩
“Maxterile™” Autoclave with Validation (IQ, OQ) Service				Ea. Price, ₩
 (4) Wire Basket Stainless-steel (Included)	DH.WAC30047	Validation Service(IQ, OQ), with “MaXterile 47”		
	DH.WAC30060	, with “MaXterile 60”		
	DH.WAC30080	, with “MaXterile 80”		
	DH.WAC30100	, with “MaXterile 100”		
 (5) Perforated Bottom Basket, Stainless-steel (Optional)	DH.WACR0047	Validation Service(IQ, OQ), with “MaXterile 47R”		
	DH.WACR0060	, with “MaXterile 60R”		
	DH.WACR0080	, with “MaXterile 80R”		
	DH.WACR0100	, with “MaXterile 100R”		
Anabac® Deodorant “Classic” for Autoclaving, 100 Tablets				
Ideal for Deodorization, Deodorant Capacity ; Max 200Lit per One-Tablet, 사과향 오토클레이브용 탈취제				
	* French Bracing Fragrance	* 100 Tablets / Packed in PP/Paper Container	<L2>	
	* Deodorant Volume ; Max 200Lit per 1-Tablet			
SL.Deo1011	Deodorant for Autoclave Apple, 100 Tablets/a Container		Pk.	Price, ₩
			100ea	

DAIHAN™ 90Lit Steam Sterilizer “Ste90.F”, Floor-type Front Door, 110~135℃, **Medicaluse**

With Built-in Steam Generator, Water Ejector Vacuum Pump, Thermal Printer, 2 Perforated STS Tray, HEPA Filter

SQUARE Chamber Autoclave, Pre- & Post-Vacuum Drying, Auto-Door Lock & Alarm System, Max 2.30 bar, **사각챔버/전면문 고압 스팀 멸균기**

110~135℃
0.4 ~ 2.3 bar

Medicaluse



3-Year Warranty

90Lit Front Door Steam Sterilizer “Ste90.F”



Fig1. Vertical Sliding Door & Square Chamber

MODE	Temperature (°C)	Pressure (bar)	Sterilizing Time (min)	Drying Time (min)
Unwrapped 1	121	1.1	20	30
Unwrapped 2	134	2.1	5	30
Wrapped 1	121	1.1	30	30
Wrapped 2	134	2.1	15	30
Liquid	121	1.1	20	0
Bowie Dick	134	2.1	4	5
User Program	110 - 135	0.4 - 2.3	4 - 30	0 - 60

Program Cycle



(1) STS Perforated Tray



(2) Thermal Paper Roll

- Supplied by DAIHAN Korea ■ **PL(Product Liability)**
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.
- Class-B Sterilization Cycle by EN285 Standard ;
Pre-vacuum → Steam Generation → Steam Injection → Sterilizing → Exhaust → Post-vacuum Drying → End of Cycle

- * Ideal for Dental Clinical, Medical, Biotechnology, etc.
- * Pre-vacuum Air Removal & Post Vacuum Drying for Perfect Sterilizing
- * Built-in Steam Generator Supplies Water Automatically and Sprays Steam into the Chamber
- * Water Ejector Vacuum Pump ; Effective for Fluid Control(1800L/hr)
- * Air Filtration System ; 0.3μm HEPA Filter of Efficiency up to 99.97%
- * Automatic Dual Door Locking System for User Safety.
- * Programmable Function ; Provides 6-Program & User Customizable 1 Program
- * Intuitive User Interface(LCD) & Indicator of Working(LED Lamp)
- * Self-Diagnostic and Display Error Codes
- * Visual and Audible Alarm System
- * Sturdy and Convenient Locking Wheel helps to Move Easily
- * Built-in Thermal Printer

※ Note/ When Autoclaving, Loosen Off the PP Screwcap!!

Specifications & Ordering information

<M3>

Spec.	Model	“Ste90.F”
Capacity		90 Lit
Sterilizing Cycle-type		Class-B Cycle by EN-285 Standard
Temp Range		110℃~135℃
Pressure & Gauge		0.4~2.3 bar, Built-in Analog & Digital Pressure Gauge
Timer & Alarm		60 min, Error Status & Timer-end Alarm
Display		Digital LCD with Back-Light Function & LED Lamp (Working Process)
Controller		Microprocessor Digital PID Control System (Full Automatic)
Programmable Function		Pre-set 6 Programs (Customizable 5- & Fixed Bowie Dick 1- Program)
Door-type		Manual Vertical Sliding Door, Automatic Door Locking System
Chamber		Reinforced and Electro Polished Square Chamber with Double Jacket
Tray (Included)		2 Stainless-steel Perforated Tray
Safety Circuit		Over-Temp / -Pressure / -Current protection, Sensor Error Detection
Printer		Thermal printer with 1 Record Paper Roll ; w58 × ϕ40mm 1000 Sterilization Cycle Data / roll (Date, Temperature, Pressure, Time, Error Status, etc.)
Material		Interior ; Stainless-steel 304 Pressure Vessel, Exterior ; Stainless-steel
Dimension (w×d×h) & Weight	Interior	338 × 712 × h378mm
	Exterior	625 × 990 × h1300mm, 280kg
	Packing	850 × 1100 × h1400mm, 295kg
Power Consumption		6.0 kW
Electric Reqts.		220V, 50/60 Hz, 1 Phase
Cat. No.		DH.SteF8001
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Accessories

- DH.SteFA1001 (1) STS Perforated Tray, 323×667×h54 mm for “Ste90.F”
- DH.SteA1006 (2) Thermal Paper Package, 4 Paper Roll for Thermal Recorder
- DH.SteA1016 Replacement HEPA Filter-kit

Front Loading Sterilizer with Validation (IQ, OQ) Service

- DH.SteF5001 Validation Service(IQ, OQ), with “Ste90.F”

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ 17 · 32 · 50Lit Dry Heat Sterilizer “SteD”, with Fast Cooling System, Thermal Printer, Forced Air, Up to 230°C, *Medicaluse*
Ideal for Dental Clinic, Hospital, Biotechnology, Fast Sterilization Time(3min at 190°C & 1hr at 160°C), 건열 멸균기

Medicaluse


3-Year Warranty

A. “SteD-17” Model


Medicaluse

3-Year Warranty

B. “SteD-30” & C. “SteD-50” Model



3-Side Heating Zone Forced-Air Mechanism



(1) STS Wire Shelf (Included)

- Supplied by DAIHAN Korea ■ *PL(Product Liability) Insurance*
- 치과, 병원, 생명공학용에 적합한 멸균기
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.
- Optimized Air Flow by Forced-air Mechanism
- 3-Side Heating Zone ; The Best Temp Uniformity & Accuracy by High Performance Heating Mechanism
- Fast Cooling System ; Powerful Cooling Fan ensures only 50min Cooling Time (Down to 50°C)
- * Digital PID Control System Implementing Superior Temperature Accuracy
- * Patented Jog-Dial & Push Button Control System
- * Space Saving Compact Size
- * Efficient Square Inner Chamber
- * Alarm Function ; Error Status and Timer-end

Specifications & Ordering information

<M1>

Spec.	Model	A. “SteD-17”	B. “SteD-30”	C. “SteD-50”
Capacity		17 Lit	32 Lit	50 Lit
Heater Power		1,400 W		
Sterilization Modes		① Custom Setting Mode		
		② 3min at 190°C		
		③ 1hr at 160°C		
Temp	Range	Room Temp +5°C to 230°C		
	Fluctuation	±1.0°C at 37°C, ±1.5°C at 160°C & 190°C		
	Variation	±1.0°C at 37°C, ±3.0°C at 160°C, ±5.0°C at 190°C		
	Sensor	PT100		
Cooling Time		Approx. 50min Down to 50°C		
Temp Resolution		± 0.1°C		
Circulation		Forced-air Type		
Controller		Digital PID Control by Advanced Microprocessor Jog-Dial with Push Button		
Display		Digital LCD with Back Light		
Timer		99hr 59min (Delay / Continuous Function)		
Printer		Thermal printer with 1 Record Paper Roll ; w58×φ40mm 1000 Sterilization Cycle Data / roll (Date, Temperature, Time, Error Status, etc.)		
Material	Internal	Stainless-steel (#304)		
	External	Powder Coated Steel		
	Insulation	Glass Wool		
Shelf (Included)		2 Stainless-steel Wire Shelf Load per Shelf ; 16kg		
Safety Device		Over Temp and Over Current Protector Sensor Error Detector		
Dimension (w×d×h)	Interior	225×215×h340mm	310×290×h360mm	370×345×h420mm
	Exterior	546×485×h517mm	458×595×h664mm	518×647×h724mm
	Packing	650×592×h600mm	562×702×h747mm	622×754×h807mm
Net Weight		30kg	36kg	44kg
Shipping Weight		38kg	43kg	59kg
Power Consumption		1.5 kW		
Electric Reqts.		230V, 50/60 Hz		

Cat. No.	DH.SteD0002	DH.SteD0003	DH.SteD0005
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Stainless-steel Wire Shelf

- DH.SteD17S (1) Wire Shelf, 240×200mm, for 17 Lit
 DH.WIG12032 , 280×260mm, for 32 Lit, “IGSN032”
 DH.WON11050 , 340×310mm, for 50 Lit, “ONS050”

Dry Heat Sterilizer with Validation (IQ,OQ) Service

- DH.SteD5001 Validation Service(IQ, OQ), with “SteD-17”
 DH.SteD5003 , with “SteD-30”
 DH.SteD5005 , with “SteD-50”

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ 18Lit Compact Low Temp Plasma Sterilizer "Ste020.P", Max. 60℃

Ideal for Medical Instrument, 50% Hydrogen Peroxide Injection, Microprocessor 7" Touch Screen Controller

With Automatic Inter Lock Door & SQUARE Chamber, 2 Stainless-steel Wire Shelves, 소형 저온 플라즈마 멸균기, 빠른멸균 (Min. 30min)



Fig. 18Lit Plasma Sterilizer with Optional Perforated Trays (2 Wire Trays are Included)



- Supplied by DAIHAN Korea ■ PL(Product Liability)
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.
- 치과, 동물병원, 소방서 등 소형 병원 사용에 적합
- Ideal for Low Temperature Plasma Sterilizer to Inactivate Microorganisms for a Broad range of Metal and Nonmetal Medical Devices and Surgical Instruments at Low Temperature.
- Improved Safety Sterilization with Special Vacuum Process.
- * It does not Damage Compatible Instruments which are Sensitive to "Heat" and "Moisture"
- * Cartridge Type H₂O₂ Injection System
- * Efficient in Using Space with Compact Size Body & Squared Inner Chamber
- * Various Operating Modes ; Quick-, Standard-, Intense-Mode
- * Processing for Sterilizing ;
 - ▶ "Ste020.P" ; Vacuum → H₂O₂ Injection → Diffusion → Vacuum → Venting

Specifications & Ordering information

<M1>

Spec.	Model	"Ste020.P"
Capacity		18 Lit
Sterilization Temp		50℃
Sterilant		50% Hydrogen Peroxide(H ₂ O ₂)
Mode		Quick ; 30min Standard ; 40min Intense ; 50min
Vacuum Pump		Built-in type, 70L/min
Sterilant Removal System		Special Vacuum Process
Controller & Display		Microprocessor 7" Touch Screen Controller
Tray (Included)		2 Stainless-steel Wire Trays, 177×300×h23mm
Safety Circuit		Over-Temp Protection, Current Leakage Breaker
Printer		Thermal printer with 1 Record Paper Roll ; w58 × φ40mm 1000 Sterilization Cycle Data / roll (Date, Temperature, Time, Error Status, etc.)
Material		Chamber ; Stainless-steel(#316), Exterior ; Stainless-steel(#304)
Dimension (w×d×h) & Weight	Interior	225×320×h265mm
	Exterior	540×725×h430mm, 65kg
	Packing	600×780×h550mm, 70kg
Power Consumption		1.3 kW
Electric Reqts.		AC 220V, 50/60Hz, 1 Phase
Cat. No.		DH.SteP020
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Accessories

- DH.StePA2001 (1) STS Wire Tray, 177×300×h23mm for "Ste020.P"
- DH.StePA2006 (2) Cartridge with Hydrogen Peroxide(50%), 100ml
- DH.SteA1006 (3) Thermal Paper Package, 4 Paper Roll for Thermal Recoder

Ea. Price, ₩

Optional Accessories

- DH.StePA2002 (4) STS Perforated Tray, 177×300×h23mm for "Ste020.P"
- DH.StePA2021 (5) Movable Stand with 1 Drawer, 240×400×610mm
- DH.StePA1030 (6) Biological Indicator, 100ea.
- DH.StePA1031 (7) Chemical Indicator, 250sheet.

Ea. Price, ₩

Articles		Cat. No	Description				
 (1) for Dry-Heat  (2) for EO-Gas  (3) for Autoclave			Premium Sterilization Indicator Tape, w19mm×L50m/Roll For Dry-Heat, EO-Gas, and Steam Sterile(Autoclaving), 멸균 감지 테이프 * Changing Color after Expose to Sterilization Process(Dry-Heat, EO-gas, Steam Sterile) (1) For Dry-Heat Sterile, <Italy-Made> KA.1051 Dry-Heat Sterile Indicator, w19mm×L50m/Roll, White→Brown (2) For EO-Gas Sterile, <Germany-Made> HH.1080551 EO-Gas Sterile Indicator, w19mm×L50m/Roll, Red→Green (3) For Steam Sterile(Autoclaving), <Germany-Made> HH.1080550 Autoclave/Steam Sterile Indicator, w19mm×L50m/Roll, Green→Brown				
			<L1> <table><tr><th>Ea.</th><th>Price, ₩</th></tr><tr><td>1roll</td><td></td></tr></table>	Ea.	Price, ₩	1roll	
Ea.	Price, ₩						
1roll							
			<table><tr><th>Ea.</th><th>Price, ₩</th></tr><tr><td>1roll</td><td></td></tr></table>	Ea.	Price, ₩	1roll	
Ea.	Price, ₩						
1roll							
			<table><tr><th>Ea.</th><th>Price, ₩</th></tr><tr><td>1roll</td><td></td></tr></table>	Ea.	Price, ₩	1roll	
Ea.	Price, ₩						
1roll							

SOOSUNG® Sterilization Indicator Tape, w12 & 25mm, L50m Roll
For Steam(Autoclaving) and EO-Gas, 멸균 감지 테이프

- 색의 변화를 통해 멸균 여부를 확인하는 데 사용
- * Changing Color after Expose to Sterilization Process(EO-gas, Steam Sterile)
- * Discoloration ; Autoclave / Yellow → Black, EO-Gas / Orange → Green

(1) For Steam Sterile(Autoclaving)

DM.Tap2001 Autoclave/Steam Sterile Indicator, w12mm×L50m/Roll, Yellow → Black
DM.Tap2002 , w25mm×L50m/Roll, Yellow → Black

(2) For EO-Gas Sterile

DM.Tap2005 EO-Gas Sterile Indicator, w25mm×L50m/Roll, Orange → Green

<M1>

Ea.	Price, ₩
1roll	
1roll	

Ea.	Price, ₩
1roll	

Sterilization Indicator Roll Bags, with Printed Graduations, L200m/Roll
For Autoclave & EO-Gas Sterilization, 의료용 멸균 감지 백

- 오토클레이브 또는 EO가스 멸균 시 포장하여 포장지의 색변화를 통해 멸균 여부를 감지할 수 있음.
- When Autoclave or EO Gas Sterilization, It can be Indicated Sterilization through a Color Change of Packaging.
- It has Printed Graduations for Convenience to Cut, It is not Necessary to Worry about Object Length.

* Bag can be Sealed with Heat Sealer * Material ; PET+CPP · Paper
* Discoloration ; Autoclave / Blue → Black. EO-Gas / Red → Brown

SL.Bag9092	Sterilization Indicator Bag, w75×L200m/Roll
SL.Bag9093	, w100×L200m/Roll
SL.Bag9094	, w150×L200m/Roll
SL.Bag9095	, w200×L200m/Roll
SL.Bag9096	, w250×L200m/Roll

Roll.	Price, ₩



Dry-Heat Sterilizer Indicator Label on Sheet, 160℃/3h, 180℃/30min, Pink ⇒ Brown
For Indicating Sterilizer, Color-Changing ; from Pink to Brown, 가열 멸균 감지 라벨

* Change of Color ; Pink ► Brown * Changing Time ; 30 minutes at 180℃, 3 Hours at 160℃

SL.Lab9051 Heat-Sterile Label "D-L", w19×L32mm

<L2>

Cs.	Price, ₩
500ea	

100% Cotton Sterilization Drape, for Sterilization of Medical Tools, 2 Layer, with Strap, 30×30~150×150cm
Ideal for Dental Clinic · Hospital · etc., Various Sizes, **기계포/기구포**

NEW 	■ 가구 소독 시 겉에 감싸는 천		■ Wrapping Cloth for Sterilizing Instruments		<M1>				
	* High Quality 100% Cotton Material								
	DM.App2001	Sterilization Drape, 2 Layer,	30×30cm	, Strap - 50cm					
	DM.App2002		, 40×40cm	, - 70cm					
	DM.App2003		, 50×50cm	, - 70cm					
	DM.App2004		, 60×60cm	, - 85cm					
	DM.App2005		, 70×70cm	, - 100cm					
	DM.App2006		, 80×80cm	, - 115cm					
	DM.App2007		, 90×90cm	, - 130cm					
	DM.App2008		, 100×100cm	, - 145cm					
	DM.App2009		, 120×120cm	, - 170cm					
DM.App2010		, 150×150cm	, - 220cm						
					<table><tr><th>Ea.</th><th>Price, ₩</th></tr><tr><td colspan="2"></td></tr></table>	Ea.	Price, ₩		
Ea.	Price, ₩								

DAIHAN Scientific

STIRRERS

Magnetic Stirrers	P. 309~312
Magnetic Stirring with Heating / Hotplate Stirrer	P. 313~315
Overhead Stirrers , High Speed-type	P. 317~319
, High Torque-type	P. 320~323
, 6-Paddle Stirrer(Jar Tester)	P. 324
, Accessories (Impeller / Rod / Stand)	P. 325~330

• Magnetic Stirrer



“MS-12D”
Mini




“MS500”
Standard



“MS-MP8”
Multi-point

• Hotplate Stirrer



“SMHS-6”



“MSH500S”

• Jar Tester / 6-Paddle Stirrer



“JT-M6C”

• Overhead Stirrer




High Speed & High Torque



DAIHAN Wisd™ 120×120 Mini Hi-Quality Digital Magnetic Stirrer “MS-12D”, Max. 2,500rpm, Max. 2Lit(H₂O)
With Tempered Glass Plate & PP Housing, Permanently Brushless Motor(BLDC), Digital LED Display, Waterproof, 소형 디지털 자력교반기

CE



120×120mm Mini Magnetic Stirrer “MS-12D”

100~2,500rpm
120×120mm

- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance
- Hi-Capacity Brushless Motor(BLDC) ; High Efficiency Stirring Performance and Precise Speed Control(Minimize RPM Fluctuation), Low Maintenance Cost.
- Outstanding Chemical Resistance due to Tempered Glass Plate(120×120mm) and PP Housing
- Higher Speed Range from 100~2,500 rpm
- Waterproof Structure
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * Digital Speed LED Display.
- * Electronically Controlled Motor for more Capacity
- * Includes 1ea × 3cm Stirring Bar
- * Can be used with up to 5cm bar


Specifications & Ordering information

<L1>

Spec.	Models	“MS-12D” 120×120mm Plate
Motor & Stirring Capacity	Permanently Brushless Motor(BLDC), Max. 2Lit (H ₂ O)	
Motor Rating(Input/Output)	6 / 3W	
Speed Range & Motor	100~2,500rpm, High Performance Brushless DC-motor	
Controller & Display	Digital PID Control with Jog-knob & Digital LED display	
Stirring Bar	3cm bar-included, Usable ; up to 5cm bar	
Material	Tempered Glass Plate(5mm thick) & PP Housing	
Dimension	Plate ; 120×120mm, Body ; 136×136×h50.5mm	
Weight	0.5kg	
Electrical Reqts.	DC 12V, with AC 230V Adapter, 50/60Hz	
Cat. No.	230V	DH.StiM012
	120V	DH.StiM012.1
Price, Set. ₩		

DAIHAN Wisd™ 180×180 Midi Hi-Quality Digital Magnetic Stirrer “MS-18D”, Max. 2,500rpm, Max. 5Lit(H₂O)
With Tempered Glass Plate & PP Housing, Permanently Brushless Motor(BLDC), Digital LED Display, Waterproof, 미디 디지털 자력교반기

CE



180×180mm Midi Magnetic Stirrer “MS-18D”

100~2,500rpm
180×180mm

- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance
- Hi-Capacity Brushless Motor(BLDC) ; High Efficiency Stirring Performance and Precise Speed Control(Minimize RPM Fluctuation), Low Maintenance Cost.
- Outstanding Chemical Resistance due to Tempered Glass Plate(180×180mm) and PP Housing
- Higher Speed Range from 100~2,500 rpm
- Waterproof Structure
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * Digital Speed LED Display.
- * Electronically Controlled Motor for more Capacity
- * Includes 1ea × 3cm Stirring Bar
- * Can be used with up to 5cm bar

Specifications & Ordering information

<L1>

Spec.	Models	“MS-18D” 180×180mm Plate
Motor & Stirring Capacity	Permanently Brushless Motor(BLDC), Max. 5Lit (H ₂ O)	
Motor Rating(Input/Output)	12 / 6W	
Speed Range & Motor	100~2,500rpm, High Performance Brushless DC-motor	
Controller & Display	Digital PID Control with Jog-knob & Digital LED display	
Stirring Bar	3cm bar-included, Usable ; up to 5cm bar	
Material	Tempered Glass Plate(5mm thick) & PP Housing	
Dimension	Plate ; 180×180mm, Body ; 193×193×h80.5mm	
Weight	1.6kg	
Electrical Reqts.	DC 12V, with AC 230V Adapter, 50/60Hz	
Cat. No.	230V	DH.StiM018
	120V	DH.StiM018.1
Price, Set. ₩		

DAIHAN™ Digital & Analog Magnetic Stirrer “MS-20”, 180×180mm Ceramic-Coated Plate, Max. 1,500rpm, Max. 20 Lit
With Precise Speed Control, Permanently Brushless Motor(BLDC), with Certi. & Traceability, 디지털/아날로그 자력교반기, 세라믹 코팅 플레이트



A. Digital Magnetic Stirrer “MS-20D”



B. Analog Magnetic Stirrer “MS-20A”

Long Life Power & Mechanism!

SMPS (Switching Mode Power Supply)
Use with Safety & Regular DC !!

Stable DC !!

SSR (Solid State Relay)-Non Contact
No Electronic Harm!
Realize **Long Life Cycle!**


High Stability in Very Low & High-speed

General Phase Control
Use with AC !!
(Very sensitive Point)

PWM Control
Use with DC !!
(Super linearity)

DAIHAN Realize **Stable & Easy to Use!**

Specialized Features of “MS-20A”



- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ceramic Coated Plate-Chemical/Acid Resistance
- Permanently Brushless Motor(BLDC)
- * Ideal for Heat-sensitive Procedures such as Microbiology and Biochemistry Applications
- * Mixing for Chemical, Laboratory, Pharmaceutical Materials
- * Smart Design of Body made by Aluminum Die-Casting
- * Powerful Motor for Stirring Quantities of up to 20Lit (Water)
- * Free Voltage (AC 80V ~ 240V, 50 / 60Hz) by SMPS

A. Digital Microprocessor Control System “MS-20D”

- Digital Feedback Control by High-Performance Microprocessor
- Patented Jog-Shuttle Control System
- * Digital LCD Display with Back-Light Function
- * Storage function for the set values ; RPM and Time
- * Digital Timer Function ; 99hr 59min(with Continuous run.)
- * Alarm Function ; Error Status and Timer-End
- * Locking Mode Supported for Experimental Safety (Jog-shuttle input blocked)

B. Analog PWM Control System “MS-20A”

- * Superior Linearity for the Control of Stirring Speed by Specially Designed Control Circuitry (PWM Control)
- * Long Life Power & Mechanism
- * High Stability in Very Low & High-speed

Specifications & Ordering information

<L1>

Spec.	Model	A. “MS-20D” Digital-type	B. “MS-20A” Analog-type
Plate Size		180×180mm	
Motor & Stirring Capa. (H ₂ O)		Permanently Brushless Motor(BLDC, Rating;10W), Max. 20Lit	
Speed Range		80 ~ 1,500 rpm	Max. 1,500 rpm
Speed Accuracy (%)		± 5	-
Speed Control Resolution		5 rpm	-
Timer & Alarm		99hr 59min, Error status & Timer-end	-
Display		Digital LCD with Back-Light	-
Controller		Digital Feedback Controller with Jog Shuttle Switch(Turn+Push)	Superior Linearity PWM Controller
Stirring Bar		3cm bar-included, Usable ; up to 5cm bar	
Permissible Ambient Temp		Room Temp +5℃~50℃	
Permissible Relative Moisture		85%	
Material	Plate	Ceramic Coated Aluminum (Chemical / Acid Resistance)	
	Body	Powder Coated Aluminum	
Dimension (w×d×h)		206×307×99mm	
Weight (kg)		3.1 kg	
Packing Size (w×d×h) & Gross Weight		280×400×200mm 5 kg	
Power Consumption		10 W	
Electrical Reqts.		Free Voltage (1 Phase, AC 80 ~ 240V, 50 / 60 Hz)	
Cat. No.	230V	DH.WMS03022	DH.WMS03122
	120V	DH.WMS02022	DH.WMS02122
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Optional Stand Rod and Holder · Clamp · Clip

- DH.WHP503036 (1) Stand Rod “RD200”, Stainless-steel, 1.5p screw,
φ12×L450mm, for “MS-20A/20D”
- DH.WHP503012 (2) Holder · Clamp · Clip, “CL220”, for Temp Probe


Ea. Price, ₩

DAIHAN™ Premium Magnetic Stirrer “MS500”, Shade Motor, Max. 80~1,500 rpm, Max. 20Lit, Solid Ceramic Plate, 200×200mm With Large LCD, Permanently Brushless Shade Motor(BLAC), Touch-Button Controller, Digital Feedback Control, Precise Speed Control
고성능 디지털 자력 교반기, 디지털 피드백 제어, 초강력 쉐이드 모터

CE

80~1,500rpm

3-Year Warranty



- Supplied by DAIHAN Korea
- CE Certified ■ Ergonomic Design ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Feedback Control by High-Performance Microprocessor
- Permanently Brushless Shade Motor & Powerful Magnet ensures Powerful Stirring and Precise Speed Control, Minimization Noise
- High-quality Solid Ceramic Plate ; Heat · Chemical · Acid-Resistance and Easy to Clean
- * Large LCD Display for Current and Speed, and Timer
- * Patented Controller “Jog-Dial” and Main Touch Button - Best Durability, Best Mobility, Best Convenient to Operate.
- * Mixing for Chemical, Laboratory, Pharmaceutical Materials
- * Storage Function for the set Values ; RPM and Timer
- * Minimization of Noise and Improved the Product Life
- * Use of the Special Advanced Insulation for the safety of the Body Temp
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled)
- * Digital Timer Function(99hr 59min)
- * Ergonomic Design of Body made by Aluminum Die-casting

Specifications & Ordering information

<L2>

Spec.		Model	“MS500”
Material	Body	Powder Coated Aluminum	
	Plate	Solid Ceramic Plate, 200×200mm	
Motor & Stirring Capa.			Permanently Brushless Shade Motor(BLAC), Max. 20 Lit
Speed & Control Resolution			80 ~ 1,500 rpm, 5 rpm
Timer & Alarm			99hr 59min (with Continuous run.), Error status & Timer-end
Display			Large Digital LCD with Back-Light Function
Controller			Digital Feedback Controller with Jog Dial (Turn+Push) and Touch button.
Stirring Bar			3cm bar-included, Usable ; up to 5cm bar
Permissible Temp & RH %			Room Temp +5℃ ~ 50℃, 85% RH
Dimension & Weight			200×332×h114mm, 3.8kg
Electrical Reqts.			1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz
Cat. No.	230V	DH.WMS00500	
	120V	DH.WMS00501	
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Optional Stand Rod and Holder · Clamp · Clip

- DH.WHP503036 (1) Stand Rod, “RD200”, Stainless-steel, 1.5p screw, $\phi 12 \times L450\text{mm}$
- DH.WHP503012 (2) Holder · Clamp · Clip, “CL220”, for Temp Probe

Ea. Price, ₩


DAIHAN™ Digital Multi-point Magnetic Stirrer “MS-MP”, 4 & 8 Points, Max. 500ml per Point, 80~1,200 rpm
With Permanently Brushless Motor(BLDC), Digital Feedback Control, Digital LCD, Synchronous Operation, with Certi. & Traceability
4 & 8 멀티포인트 디지털 자력교반기, 디지털 피드백 컨트롤 시스템

CE 80~1,200rpm 3-Year Warranty



A. Digital 4-Point Magnetic Stirrer “MS-MP4”

CE 80~1,200rpm 3-Year Warranty



B. Digital 8-Point Magnetic Stirrer “MS-MP8”

- Supplied by DAIHAN Korea
- CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideally Suitable for Use in Incubators and Climatic Chambers
- Excellent for Multi-use and wide range of Sample Testing
- Patented Jog-Shuttle Control System ■ Permanently Brushless Motor(BLDC)
- * Ideal for Heat-sensitive Procedures such as Microbiology and Biochemistry Applications
- * Synchronous Operation of 4 / 8 Stirring Points * Up to 0.5 Lit. Stirring Capacity per Point
- * Increased Effectiveness by Microprocessor Control
- * Digital LCD with Back-Light Function * Digital Timer Function (99 hr 59 min)
- * Smooth Start and Change of Stirring Speed
- * Storage Function for the Set Values ; RPM and Timer
- * Extraordinary Mechanical, Physical and Chemical Resistance
- * Locking Mode supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Free Voltage (AC 80V ~ 240V, 50 / 60Hz) by SMPS


Specifications & Ordering information

<L1>

Spec.	Model	A. “MS-MP4” 4 Stirring Points	B. “MS-MP8” 8 Stirring Points
Number of Stirring Points		4 (2×2) Stirring Points	8 (4×2) Stirring Points
Plate Size		208×215mm	208×415mm
Stirring Capacity (H ₂ O)		Max. 0.5 Lit per Point	
Distance between Stirring Points		100 mm	
Speed Range		80 ~ 1,200 rpm	
Motor		Permanently Brushless Motor(BLDC, Rating ; 16/9 W)	
Speed Control Resolution		5 rpm	
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Digital Feedback Controller with Jog Shuttle Switch(Turn+Push)	
Stirring Bar		4ea 3cm bar-included, Usable ; up to 3cm bar	8ea 3cm bar-included, Usable ; up to 3cm bar
Permissible Ambient Temperature		Room Temp +5℃~50℃	
Permissible Relative Moisture		85%	
Material	Plate	Stainless-steel (# 304)	
	Body	Power Coated Steel	
Dimension (w×d×h)		208×280×75mm	208×480×75mm
Weight		2.5 kg	3.5 kg
Others		PL(Product Liability) Insurance Storage Function (RPM, Timer) Locking Mode (Jog-Shuttle Input Disabled) Alarm (Error Status and Timer-end)	
Packing Size (w×d×h) & Gross Weight		300×330×180mm 6 kg	330×540×180mm 6 kg
Power Consumption		20 W	
Electrical Reqts.		Free Voltage (1 Phase, AC 80~240V, 50/60 Hz)	
Cat. No.	230V	DH.WMS03036	DH.WMS03038
	120V	DH.WMS02036	DH.WMS02038
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.


DAIHAN™ Premium 500°C High-Temp Hotplate Stirrer “MSH500”, Solid Ceramic Plate, 200×200mm, 80~1,500rpm With Permanently Brushless Shade Motor(BLAC), Large LCD, Optimum Insulation Layer. Accurate Temp Control, Hot-Top Indicator 고온용 디지털 가열 자력 교반기, 최적의 방열구조, 디지털 피드백 온도 제어, 내열성 · 내화학성 · 내산성 솔리드 세라믹 글라스 플레이트




A. “MSH500S” Package-Set



B. “MSH500” Unit Only



Controller of “MSH500” Model



- Supplied by DAIHAN Korea
- UL, CUL, and CE Certified
- PL(Product Liability) Insurance
- Ergonomic Design
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- High-quality Solid Ceramic Plate ; Heat-/Chemical-/Acid-Resistance and Easy to Clean
- Optimum Insulation Fourfold Layer ; BTU + Ceramic Paper + 2 STS Insulation Layer
- Permanently Brushless Shade Motor & Powerful Magnet ensures Powerful Stirring and Precise Speed Control, Minimization Noise
- * Microprocessor PID control
- * Patented Controller “Jog-Dial” and Main Touch Button - Best Durability, Best
- * Mobility, Best Convenient to Operate with Back-Light Function
- * 2 Internal K-type Temp Sensor Operating Independently ;
 - 1 Temp Control Sensor for Temperature Accuracy
 - 1 Over Temp Detection Sensor for User Safety
- * Hot Top Indicator LED Lamp ; Hot Surface(Over 50℃) Warning Lamp for User Safety
- * Excellent for Sample Drying and most General Lab Heating Applications
- * User's Self -Compensation Function ; to Control Difference between Real Temp and Indicated Temp (±20.0℃)
- * Accurate Temp Control of the Sample in the Vessel using PT100 Sensor (Optional)
- * Storage Function for the set Values ; Temp and Timer
- * Use of the Special Advanced Insulation for the safety of the Body Temp
- * Locking Mode Supported for Experimental Safety (Input to Jog-Dial can be Disabled)

Specifications & Ordering information

<L2>

Spec.		Model	A. “MSH500S” (Package-Set)	B. “MSH500” (Unit Only)
Material	Body	Powder Coated Aluminum		
	Plate	Solid Ceramic Plate, 200×200mm		
Motor & Stirring Capa.		Permanently Brushless Shade Motor(BLAC), Max. 20 Lit		
Speed & Control Resolution		80~1,500 rpm, 5 rpm		
Temp Range & Accuracy		Max. 500℃, ±0.3℃ at Set Temperature		
Temp Resolution		0.1℃-Display, 0.5℃-Control		
Temp Uniformity		Specially Designed Heating Module ; Temp Difference 10% less		
Heating Power / Consumption		850W		
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller		Digital Feedback Controller with Jog Dial (Turn+Push) and Touch button.		
Stirring Bar		3cm bar-included, Usable ; up to 5cm bar		
Safety Circuit		Over Temp Protection, Hot Top Indicator LED Lamp		
included Standard Acc.		1ea Stir Bar, ②Temp Probe ③Rod ④Holder+Clamp		1ea Stir Bar
Permissible Temp & RH %		Ambient +5℃~50℃, 85% RH		
Dimension & Weight		200×332×h114mm, 3.8kg		
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz, with Cord/Plug		
Cat. No.	230V	DH.WMH00500		DH.WMH01500
	120V	DH.WMH00501		DH.WMH01501
Price, Each. ₩				

* Other Specifications are available upon Customer's Request.

Component of Package-Set

- * Package-set ; ① Hotplate Stirrer Unit ② “SS500” PT100 External / Direct Contact Temp Probe,
③ “RD200” Stainless-steel Support Rod $\phi 12 \times L450\text{mm}$
④ “CL220” Holder+Clamp+Clip for External Temp Probe

Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

Cat. No.		Ea. Price, ₩	
DH.WHP503060	(1) External / Direct Contact Temp Probe, “SS500”, Stainless-steel, $\phi 4 \times L250\text{mm}$, Cable L190mm		
DH.WHP503061	(2) External / Direct Contact Temp Probe, “SS510”, Teflon Coated, $\phi 4 \times L250\text{mm}$, Cable L190mm,		
DH.WHP503036	(3) Stand Rod “RD200”, Stainless-steel, $\phi 12 \times L450\text{mm}$, 1.5p screw, for “MSH500”		
DH.WHP503012	(4) Holder · Clamp · Clip, “CL220”, for Temp Probe		

DAIHAN™ 380°C Analog & Digital Hotplate Stirrer “MSH-A” & “MSH-D”, Ceramic-Coated Plate, 80~1,500 rpm
180×180mm or 260×260mm Plate, with Accurate Temp. Control, Superior Temp Uniformity, with Certi. & Traceability
Permanently Brushless Motor(BLDC), 아날로그 & 디지털 가열 자력 교반기, 우수한 온도제어 · 온도균일성 · 세라믹 코팅 플레이트

A. 180×180mm Plate,

I. Analog Unit “MSH-20A”

II. Digital Unit & III. Package Set

“MSH-20D” “MSH-20D” Set



- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Excellent for Sample Drying and most General Lab Heating Applications
- Permanently Brushless Motor(BLDC)

I. Analog PWM Control System “MSH-A”

- * Superior Linearity PWM Control for Temperature & Stirring Speed

II. Digital Microprocessor Control System “MSH-D”

- **Patented Jog-Shuttle Control System**
- * Speed & Control Resolution ; 80 ~ 1,500 rpm, 5 rpm
- * Temp Accuracy / Resolution ; $\pm 0.3^{\circ}\text{C}$ at Set Temp / 0.1°C -Display, 0.5°C -Control
- * Timer & Alarm ; 99hr 59min(with Continuous run), Error status & Timer-end
- * Safety Circuit ; Over Temp Protection
- * User's Self-Compensation ; Control Difference between Real Temp & Indicated Temp ($\pm 20.0^{\circ}\text{C}$)
- * Accurate Temp Control of the Sample in the Vessel using PT100 Sensor
- * Storage Function for the set Values ; Temp and Timer
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)

III. Digital Package-set ; Complete with

- * **III-A. “MSH-20D” Set** ; ① Hotplate Stirrer “MSH-20D”, ② “SS500” PT100 External / Direct Contact Temp Probe, ③ “RD200” Stainless-steel Support Rod $\phi 12 \times L450\text{mm}$, ④ “CL220” Holder+Clamp+Clip for External Temp Probe
- * **III-B. “MSH-30D” Set** ; ① Hotplate Stirrer “MSH-30D”, ② “SS500” PT100 External / Direct Contact Temp Probe, ③ “RD100” Stainless-steel Support Rod $\phi 12.7 \times L390\text{mm}$, ④ “CL220” Holder+Clamp+Clip for External Temp Probe

Specifications & Ordering information

<L1>

Spec.	Model	A. 180×180mm Plate (I-A) Analog “MSH-20A” (II-A) Digital “MSH-20D”	B. 260×260mm Plate (I-B) Analog “MSH-30A” (II-B) Digital “MSH-30D”
Plate Size(mm)		180×180mm Plate-type	260×260mm Plate-type
Material	Body	Powder Coated Aluminum	Powder Coated steel
	Plate	Ceramic Coated Aluminum (Chemical / Acid Resistance, Easy to Clean)	
Motor & Stirring Capa.		Permanently Brushless Motor(BLDC, Rating ; 10w), Max. 20Lit	
Speed & Control Resolution		Digital ; 80~1500 rpm, 5 rpm / Analog ; Max. 1500 rpm	
Temp Uniformity		Specially Designed Heating Module ; Temp Difference 10% less	
Heating Power/Consumption		600W	1.2kW
Stirring Bar		3cm bar-included, Usable ; up to 5cm bar	
Permissible Temp & RH %		Ambient +5°C~50°C, 85% RH	
Dimension & Weight		206×307×h99mm, 3.3kg	286×387×h99mm, 3.5kg
Packing Size & Gross Weight		400×280×h200mm, 5kg	430×320×h200mm, 5kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	

I. Analog Unit Only		(I-A) “MSH-20A”	(I-B) “MSH-30A”
Cat. No.	230V	DH.WMH03120	DH.WMH03121
	120V	DH.WMH02120	DH.WMH02121
Price, Each. ₩			
II. Digital Unit Only		(II-A) “MSH-20D”	(II-B) “MSH-30D”
Cat. No.	230V	DH.WMH03020	DH.WMH03022
	120V	DH.WMH02020	DH.WMH02022
Price, Each. ₩			
III. Digital Package-set		(III-A) “MSH-20D” Set	(III-B) “MSH-30D” Set
Cat. No.	230V	DH.WMH03021	DH.WMH03023
	120V	DH.WMH02021	DH.WMH02023
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

- DH.WHP503060 (1) Temp Probe, “SS500”, “PT100”, $\phi 4 \times L250\text{mm}$
 DH.WHP503061 (2) Teflon Coated Temp Probe, “SS510”, $\phi 4 \times L250\text{mm}$
 DH.WHP503011 Stand Rod, “RD100”, $\phi 12.7 \times L390$ for “MSH-30A / 30D”
 DH.WHP503036 “RD200”, $\phi 12 \times L450$ for “MSH-20A / 20D”
 DH.WHP503012 Holder · Clamp · Clip, “CL220”, for Temp Probe

Ea. Price, ₩



III-A. Digital Package Set
“MSH-20D” Set

B. 260×260mm Plate,


I. Analog Unit “MSH-30A”

II. Digital Unit & III. Package Set


“MSH-30D” “MSH-30D” Set




(II-B) Digital “MSH-30D”, 260×260mm Plate



DAIHAN™ 350°C Systematic Multi-Hotplate Stirrer “SMHS”, 3 or 6 Places, Ceramic-Coated Plate, 80~1,500 rpm
 With Digital Feedback Control, Independent Heating & Stirring Control, Permanently Brushless Motor(BLDC)
 멀티 가열 자력 교반기, 우수한 온도균일성, 3 or 6 구 개별 조절 가능




A. 3 Plate-type Multi-Hotplate Stirrer “SMHS-3”



B. 6(3×2) Plate-type Multi-Hotplate Stirrer “SMHS-6”

※ “SMHS-6” has 3 Ports for Connecting External Temp Probes.



- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Excellent for Multi-use and Wide Range Sample Testing
- Patented Jog-Shuttle Control System
- Permanently Brushless Motor(BLDC)

- * Digital Feedback Control by High-Performance Microprocessor
- * Ceramic Coated Plate-Chemical/Acid Resistance
- * Digital Timer Function Included ; 99 hr 59 min
- * Maximum Temperature of 350°C with $\pm 0.3^\circ\text{C}$ Control Accuracy
- * Accurate Temp. Control of the Sample in the Vessel using PT100 Sensor (Optional)
- * Storage Function for the Set Values ; Temperature, RPM and Timer
- * Locking Mode for Experimental Safety Supported

Specifications & Ordering information

<L1>

Spec.	Model	A. “SMHS-3” 3 Plates	B. “SMHS-6” 6 Plates
Plate-type		3 Plate-type	6(3×2) Plate-type
Plate Size		140×140mm	
Motor & Stirring Capa.		Permanently Brushless Motor(BLDC), Max. 20Lit(H ₂ O)	
Motor Rating		9/4W×3	9/4W×6
Speed & Control Resolution		80 ~ 1,500 rpm, 5 rpm	
Temp. Range & Accuracy		Max. 350°C, $\pm 0.3^\circ\text{C}$ at Set Temp	
Temp. Resolution		0.1°C-Display, 0.5°C-Control	
Temp. Uniformity		Specially Designed Heating Module ; Temp Difference 10% less	
Heating Power/Consumption		400 W × 3 Plates	400 W × 6 Plates
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Display		Digital LCD with Back-Light Function	
Controller		Digital Feedback Controller with Jog Shuttle Switch (Turn+Push)	
Stirring Bar		3cm bar-included, Usable ; up to 5cm bar	
Ports for Stand Rod		2ea	4ea
Permissible Temp. & RH %		Room Temp + 5°C~50°C, 85% RH	
Material	Body	Powder Coated steel	
	Plate	Ceramic Coated Aluminum	
Dimension & Weight		510×240×h100mm, 10kg	510×410×h180mm, 18kg
Packing Size & Gross Weight		540×330×h170mm, 10kg	570×540×h230mm, 18kg
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz	
Cat. No.	230V	DH.WMH03503	DH.WMH03506
	120V	DH.WMH02503	DH.WMH02506
Price, Each. ₩			

• Other Specifications are available upon Customer's Request.

※ 운임 별도

※ “SMHS-6” has 3 Ports for Connecting External Temp Probes.

Optional External Temp Probe, Stand Rod, and Holder · Clamp · Clip

- DH.WHP503060 (1) External / Direct Contact Temp Probe, “SS500”, Stainless-steel, $\phi 4 \times L250\text{mm}$, Cable L190mm
- DH.WHP503061 (2) External / Direct Contact Temp Probe, “SS510”, Teflon Coated, $\phi 4 \times L250\text{mm}$, Cable L190mm,
- DH.WHP503036 (3) Stand Rod “RD200”, Stainless-steel, 1.5p screw, $\phi 12 \times L450\text{mm}$, for “SMHS”
- DH.WHP503012 (4) Holder · Clamp · Clip, “CL220”, for Temp Probe

Ea.	Price, ₩

DAIHAN™ Overhead Stirrer Selection Guide of Product Line

Key Specification of Product Line

Speed or High Torque	High Speed							High Torque							
Controller Type	Analog Model		Digital Model					Digital Model							
	Integrated Controller		Integrated Controller		Remote Controller		Premium Digital Model	Integrated Controller		Remote Controller			Premium Digital Model		
Model Name	HS-50A	HS-120A	HS-30D	HS-100D	HS-30T	HS-100T	OVS-S05	HT-50DX	HT-120DX	HT-50T	HT-120T	HT-140T	OVS-T16	OVS-G30	
Max. Stirring Capacity(H ₂ O)	10Lit	20Lit	10Lit	20Lit	10Lit	20Lit	100Lit	40Lit	60Lit	40Lit	60Lit	100Lit	100Lit	100Lit	
Max. Viscosity(mPas) (Viscosity Standard Oil 1.5Lit)	5,000	10,000	5,000	10,000	5,000	10,000	50,000	60,000	90,000	60,000	90,000	120,000	200,000	350,000	
Motor Rating (Input / Output)	83W /75W	130W /120W	83W /75W	130W /120W	83W /75W	130W /120W	200W /180W	83W /75W	130W /120W	83W /75W	130W /120W	130W /120W	200W /180W	200W /180W	
Max. Torque (N·cm)	85	145	70	80	70	80	200	180	210	180	210	300	700	800	
Speed Range (RPM)	200 ~ 3,000							50 ~ 1,000						50 ~ 2,000	10 ~ 550
Max. Push Through Shafts ϕ (mm)	8						10							-	
Max. Chcuk Grip Capa. ϕ (mm)	10							12						14	
Controller	Analog Phase Controller		Digital Feedback Controller with Jog Shuttle Switch (Turn+Push)		Digital Feedback Controller with Jog Shuttle Switch (Turn+Push), Remote type		Digital Feedback Controller with Jog Shuttle Switch (Turn + Push)	Digital Feedback Controller with Jog Shuttle Switch (Turn+Push)		Digital Feedback Controller with Jog Shuttle Switch (Turn+Push), Remote type			Digital Feedback Controller with Jog Shuttle Switch (Turn+Push)		
Display	None		Back Light LCD					Back Light LCD							
Permissible Temp. & Humidity	0 ~ 40℃ / 80%														
Power Supply	1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz														

Selection Guide for Standard Viscosity

(Unit ; mPas)

Material	Viscosity	Recommended Models									
Water	1	HS-50A HS-30D HS-30T	HS-120A HS-100D HS-100T	OVS-S05	HT-50DX HT-50T	HT-120DX HT-120T	HT-140T	OVS-T16	OVS-G30		
Milk	3										
Soybean oil	160										
Varnish	500										
Lacquer	1,000										
Honey	3,000	-	-	-	-	-	-	-	-	-	-
Viscosity Standard Oil	5,000										
	10,000										
	15,000										
Ketchup	30,000										
Viscosity Standard Oil	50,000										
	60,000										
Mustard	70,000										
Viscosity Standard Oil	90,000										
Corn Syrup	100,000										
	120,000										
Viscosity Standard Oil	150,000										
	200,000										
Peanut Butter	250,000										
Viscosity Standard Oil	350,000										

DAIHAN™ **Premium** Hi-Speed Overhead Stirrer “OVS-S05”, with Permanently Brushless Motor(BLDC), 2:1, Max. 3,000rpm
 With Torque(Ncm) · Viscosity(mPas) · Temperature(°C) · Real Time Display, Optional Remote Control, 50,000mPas, Max. 100Lit
 “Push-Through” Shaft(ϕ 10mm or Less), 프리미엄 고속 교반기, 저/중점동용, 고사양 BLDC-Motor, 토크 · 점도 · 온도 표시기능



- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- Patented Jog Dial & Touch Button
- Hi-Capacity Brushless Motor(BLDC) ; High Efficiency Stirring Performance and Precise Speed Control(Minimize RPM Fluctuation), Low Maintenance Cost
- Torque Value(Ncm)-/Viscosity(mPas)-/Temperature(°C) Real-time Display
- Can be used with K-type Temp Probe for Temperature Display.
- Can be used with Optional Digital Remote Controller.
- “Push-Through” Shafts(ϕ 10mm or Less) ; Ideal for Raising or Lowering the Stir-Shafts without Moving the Motor Head - “Adjustable Height of Shaft”
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

- * Ideal for Low-/Middle-Viscosity
- * Smooth Start and Change of the Stirring Speed.
- * Large Size Digital LCD Display

Specifications & Ordering information

<L1>

Spec.	Model	“OVS-S05” Max. Speed 3,000rpm
Stirring Capacity		Max. 100Lit (H ₂ O)
Maximum Viscosity / Torque		50,000 mPas / 200 Ncm
Speed Range & Motor		200~3,000RPM, Hi-Capacity Brushless Motor(BLDC), 2:1
Motor Rating (Input / Output)		200 / 180W
Controller / Safety Device		Digital PID Control / Over Heat & Current Protector for Motor
Chuck Grip Capa. ϕ mm		dia ϕ 3~10mm Shafts, Flex-Coupling(ϕ 6~12mm)
“Push-Through” Shafts		dia ϕ 10mm or Less
Permissible Temp & Humidity		0°C~40°C, 80%
Dimension & Weight & Material		84×205×h275mm, 3.5kg, Aluminum Die Casting
Power Consumption		250W
Electrical Reqts.		1 Phase, AC 120V, 60 Hz or AC 230 V, 50 / 60 Hz

* Other Specifications are available upon Customer's Request



A. Unit-Only

DH.StiO0005 Hi-Speed Overhead Stirrer Unit “OVS-S05”, 230V
 DH.StiO0005.1 , 120V

Ea. Price, ₩

B. Package-Set

- * Complete with ① “OVS-S05” StirrerUnit
- ② “S723” Stainless-steel Impeller (Rod ϕ 8×L650mm, 3-Blade ϕ 80mm)
- ③ “CL200” Clamp
- ④ “KA1177” Safety Stand (Rod ϕ 23×L800mm, Base 34×29cm)

DH.StiO0105 Overhead Stirrer “OVS-S05” Package-set, 230V
 DH.StiO0105.1 , 120V

Ea. Price, ₩

C. Optional Digital Remote Controller, with Connecting Cable 1.5m

- * Assembly-type, Torque / Viscosity / Temp-Control & Display
- * Connecting Cable 1.5m(Included) and 2m(Optional)
- * Size ; w75×d150×h30mm, Weight ; 150g

DH.WOS50010 Digital Remote Controller with 1.5m Connecting Cable & Plug
 DH.Cabl002 Optional Connecting Cable & Plug, 2m

Ea. Price, ₩

Optional Temperature Probe


(1) Direct Contact Temp Probe, K-type

DH.WHP503040 Temp Probe, K-type, “SS300”, L250× ϕ 4.8 mm
 DH.WHP503050 , “SS400”, L380× ϕ 4.8 mm
 DH.WHP503041 , “SS310”, Teflon Coated, L250× ϕ 5 mm
 DH.WHP503051 , “SS410”, Teflon Coated, L380× ϕ 5 mm

Ea. Price, ₩

(2) K-type PTFE Encapsulated Stainless-steel Thermocouple Probe

※ Note! Optional Probes ⇒ See the “Measuring, Thermometers, K-type” Section.
 Optional Stirrer Shafts ⇒ See the “Stirrer, Accessories” Section.



DAIHAN™ High-Speed Analog Overhead Stirrer “HS-A”, for Low-/Middle-Viscosity, Max. 10 & 20Lit, Max. 3,000rpm
With “Push-Through” Shaft(φ 8mm or Less) and Chuck grip φ 3~10mm, Flex-Coupling φ 6~12mm, Up to 5,000 · 10,000mPas
아날로그 고속 교반기, 저/중점도용, 사용교반봉; ①기본은 φ 3~10mm, ②관통은 φ 3~8mm, & ③별도의 Flex-Coupling 활용시 φ 6~12mm



Max. 10Lit (H₂O) High-Speed Analog Overhead Stirrer
“HS-50A” Package-Set

- Supplied by DAIHAN Korea
- CE Certified
- PL (Product Liability) Insurance
- Ensure Minimize RPM Fluctuation ; Very Stable Stirring Speed(RPM) Performance
- * with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * “Push-Through” Shafts(φ 8mm or Less) ; Ideal for Raising or Lowering of the Stir-Shafts without Moving the Motor Head - “Adjustable Height of Shaft”
- * Easy-to-use Dial for Speed Setting
- * Smooth Start and Change of the Stirring Speed

Specifications & Ordering information

<L1>


Spec.	Models	“HS-50A” Max. Stirring Capa. 10Lit	“HS-120A” Max. Stirring Capa. 20Lit
Stirring Capacity		Max. 10Lit (H ₂ O)	Max. 20Lit (H ₂ O)
Maximum Viscosity / Torque		5,000 mPas / 85 Ncm	10,000 mPas / 145 Ncm
Speed Range & Motor		200 ~ 3,000 rpm, Direct Driven DC Motor	
Motor Rating (Input / Output)		83 / 75W	130 / 120W
Controller / Safety Device		Analog Phase Controller / Over Heat Protector for motor	
Chuck Grip Capa. φ mm		dia φ 3~10mm Shafts	
“Push-Through” Shafts		dia φ 8mm or Less	
Permissible Temp. & Humidity		0℃~40℃, 80%	
Material		ABS and Aluminum	
Dimension & Weight		90 × 145 × h180mm, 4.5kg	90 × 145 × h180mm, 5.3kg
Packing Size & Gross Weight		320 × 300 × h180mm, 5kg	320 × 300 × h180mm, 8kg
Power Consumption		85W	130W
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	

Unit-only

Cat. No.	230V	DH.WOS01210	DH.WOS01215
	120V	DH.WOS02210	DH.WOS02215
Price, Each. ₩			
Package-Set (with ① Stirrer Unit, ② Impeller, ③ Clamp, ④ Stand)			
Cat. No.	230V	DH.WOS01211	DH.WOS01216
	120V	DH.WOS02211	DH.WOS02216
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

The Figure ; Use of Coupling



Hi-flex Metal for φ 6~16mm Shafts

- * “HS-50A” Package-Set ; Complete with ① “HS-50A” Stirrer Unit, ② “S721” Stainless-steel Impeller (Rod φ 8 × L500mm, 3-Blade φ 50mm), ③ “CL200” Clamp, and ④ “KA1177” Safety Stand (Rod φ 23 × L800mm, Base w34 × d29cm)
- * “HS-120A” Package-Set ; Complete with ① “HS-120A” Stirrer Unit, ② “S723” Stainless-steel Impeller (Rod φ 8 × L650mm, 3-Blade φ 80mm), ③ “CL200” Clamp, and ④ “KA1177” Safety Stand (Rod φ 23 × L800mm, Base w34 × d29cm)

Flexible Coupling, for Safety Stirring of φ 6~16mm Shafts of Overhead Stirrers For Stirrers between Stirring Motor and Shaft, 안전교반용 유연성 커플링

- Supplied by DAIHAN Korea
- Used to Compensate Any Vibrations Between The Motor and Stirring Shafts
- * Made of DURALUMIN
- * Shaft Fastening by Screw Clamp Band
- * with Exchangeable STS Chucking Rod φ 8 / 10 / 12mm, L40mm

<L1>


	Ea.	Price, ₩
SL.StiC011	Stir-Coupling, Hi-flex, φ 8mm Chuck-rod , for-φ 6mm Stir-Shaft , 56g , φ 25 × h44mm	
SL.StiC012	-φ 8mm Stir-Shaft , 53g , φ 25 × h44mm	
SL.StiC013	-φ 10mm Stir-Shaft, 49g , φ 25 × h44mm	
SL.StiC014	, φ 10mm Chuck-rod, for-φ 8mm Stir-Shaft , 49g , φ 25 × h44mm	
SL.StiC015	-φ 10mm Stir-Shaft, 43g , φ 25 × h44mm	
SL.StiC016	-φ 12mm Stir-Shaft, 33g , φ 25 × h44mm	
SL.StiC017	, φ 12mm Chuck-rod, for-φ 10mm Stir-Shaft, 33g , φ 25 × h44mm	
SL.StiC018	-φ 12mm Stir-Shaft, 123g, φ 32 × h60mm	
SL.StiC019	-φ 14mm Stir-Shaft, 120g, φ 32 × h60mm	
SL.StiC020	-φ 16mm Stir-Shaft, 118g, φ 32 × h60mm	

※ Note! Details of Wide-Range Coupling ; See the “Stirrer, Accessories” Section.

DAIHAN™ High-Speed Digital Overhead Stirrers “HS-D” & “HS-T”, for Low-/Middle-Viscosity, Max. 10 & 20Lit, Max. 3,000rpm With “Push-Through” Shaft(ϕ 8mm or Less) and Chuck grip ϕ 3~10mm, Flex-Coupling ϕ 6~12mm, up to 5,000 · 10,000mPas
 디지털 고속 교반기, 저/중점도용, 사용교반봉 : ①기본은 ϕ 3~10mm, ②관통은 ϕ 3~8mm, & ③별도의 Flex-Coupling 활용시 ϕ 6~12mm



(I-A) Max. 10Lit (H₂O) High-Speed Digital Overhead Stirrer “HS-30D” Package-Set



(II-A) Max. 10Lit (H₂O) High-Speed Digital Overhead Stirrer, with Remote Controller “HS-30T” Package-Set

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- Ensure Minimize RPM Fluctuation ; Very Stable Stirring Speed(RPM) Performance
- Patented Jog-Shuttle Control System
- * with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * “Push-Through” Shafts(ϕ 8mm or Less) ; Ideal for Raising or Lowering of the Stir-Shafts without Moving the Motor Head - “Adjustable Height of Shaft”
- * Smooth Start and Change of the Stirring Speed

I. Digital Standard-type “HS-D”

<L1>

Spec.	Models	A. “HS-30D” Max. Stirring Capa. 10Lit	B. “HS-100D” Max. Stirring Capa. 20Lit
Maximum Viscosity / Torque		5,000 mPas / 70 Ncm	10,000 mPas / 80 Ncm
Speed Range & Resolution		200~3,000 rpm, 1 rpm-Display / 1 rpm-Control	
Motor / Safety		Device Direct Driven DC Motor / Overload & Heat Protector for Motor	
Motor Rating (Input / Output)		83 / 75W	130 / 120W
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end	
Controller / Display		Digital Feedback Controller / Digital LCD with Back-Light Function	
Chuck Grip Capa. ϕ mm		dia ϕ 3~10mm Shafts	
“Push-Through” Shafts		dia ϕ 8mm or Less	
Permissible Temp. & Humidity		0℃~40℃, 80%	
Material		ABS and Aluminum	
Dimension & Weight		90×145×h180mm, 4.5kg	90×145×h180mm, 5.3kg
Packing Size & Gross Weight		320×300×h180mm, 5kg	320×300×h180mm, 8kg
Power Consumption		85W	130W
Electrical Reqts.		1 Phase, AC 120V, 60 Hz or AC 230V, 50 / 60 Hz	

Unit-only

Cat. No.	230V	DH.WOS01015	DH.WOS01025
	120V	DH.WOS02015	DH.WOS02025

Price, Each. ₩

Digital Package-Set (with ① Stirrer Unit, ② Impeller, ③ Clamp, ④ Stand)

Cat. No.	230V	DH.WOS01016	DH.WOS01026
	120V	DH.WOS02016	DH.WOS02026

Price, Each. ₩

II. Digital Remote Controlled-type “HS-T”

Spec.	Model	A. “HS-30T” Max. Stirring Capa. 10Lit	B. “HS-100T” Max. Stirring Capa. 20Lit
Maximum Viscosity / Torque		5,000 mPas / 70 Ncm	10,000 mPas / 80 Ncm
Speed Range & Motor		200~3,000 rpm, Direct Driven DC Motor	
Motor Rating (Input / Output)		83 / 75W	130 / 120W
Timer & Alarm		99hr 59min (with Continuous run.), Timer-end Alarm	
Controller / Display		Digital Feedback Controller / Digital LCD with Back-Light Function	
Chuck Grip Capa. ϕ mm		dia ϕ 3~10mm Shafts	
“Push-Through” Shafts		dia ϕ 8mm or Less	
Permissible Temp. & Humidity		0℃~40℃, 80%	
Safety Device / Material		Overload & Heat Protector for motor / ABS and Aluminum	
Dimension & Weight*		90×100×h180mm, 3kg	90×100×h180mm, 3.8kg
Packing Size & Gross Weight		580×340×h230mm, 5kg	580×340×h230mm, 8kg
Power Consumption		75W	120W
Electrical Reqts.		1 Phase, AC 120V, 60 Hz or AC 230V, 50 / 60 Hz	

Unit & Remote Controller

Cat. No.	230V	DH.WOS01035	DH.WOS01045
	120V	DH.WOS02035	DH.WOS02045

Price, Each. ₩

Remote Package-Set (with ① Stirrer Unit + Remote Controller, ② Impeller, ③ Clamp, ④ Stand)

Cat. No.	230V	DH.WOS01036	DH.WOS01046
	120V	DH.WOS02036	DH.WOS02046

Price, Each. ₩

* Weight of Remote Controller = 2.5kg

** Other Specifications are available upon Customer's Request.

* “HS-30D(30T)” Package-Set ; Complete with ① “HS-30D(30T)” Stirrer Unit (with Remote Controller), ② “S721” Stainless-steel Impeller(Rod ϕ 8×L500mm, 3-Blade ϕ 50mm), ③ “CL200” Clamp, and ④ “ST110” Standard Plate Stand(Rod ϕ 16×L800mm, Base w20×d35cm)

* “HS-100D(100T)” Package-Set ; Complete with ① “HS-100D(100T)” Stirrer Unit (with Remote Controller), ② “S723” Stainless-steel Impeller(Rod ϕ 8×L650mm, 3-Blade ϕ 80mm), ③ “CL200” Clamp, and ④ “ST120” Wide Plate Stand(Rod ϕ 23×L1,000mm, Base w30×d40cm)

DAIHAN™ **Premium** Hi-Torque Overhead Stirrer “OVS-T16”, with Permanently Brushless Motor(BLDC), Max. 14:1, 200,000mPas With Torque(Ncm) · Viscosity(mPas) · Temperature(°C) · Real Time Display, Optional Remote Control, Max. 2,000rpm “Push-Through” Shaft(φ 10mm or Less), Chuck Grip Capa. φ 6~14mm Shafts, 프리미엄 중/고점도용 교반기, 고사양 BLDC-Motor, 토크 · 점도 · 온도 표시기능



- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- 다양한 점도 교반을 위해 φ 6~14mm 의 Impeller Shaft가 장착 가능한 Chuck Grip.
- Patented Jog Dial & Touch Button
- Hi-Capacity Brushless Motor(BLDC) ; High Efficiency Stirring Performance and Precise Speed Control(Minimize RPM Fluctuation), Low Maintenance Cost
- Torque Value(Ncm)-/Viscosity(mPas)-/Temperature(°C) Real-time Display Can be used with K-type Temp Probe for Temperature Display.
- Can be used with Optional Digital Remote Controller.
- “Push-Through” Shafts(φ 10mm or Less) ; Ideal for Raising or Lowering the Stir-Shafts without Moving the Motor Head - “Adjustable Height of Shaft”
- Wide Range Chuck Grip Capa. φ 6~14mm for High Viscosity.
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * Ideal for Low-/Middle-Viscosity
- * Smooth Start and Change of the Stirring Speed.
- * Large Size Digital LCD Display

Specifications & Ordering information

<L1>

Spec.	Model	“OVS-T16” Max. Torque 200,000mPas
Stirring Capacity		Max. 100Lit (H ₂ O)
Maximum Viscosity / Torque		200,000 mPas / 700 Ncm
Speed Range & Motor		50~2,000RPM, Hi-Capacity Brushless Motor(BLDC), Gear(2.6:1 / 14:1)
Motor Rating (Input / Output)		200 / 180W
Controller / Safety Device		Digital PID Control / Over Heat & Current Protector for Motor
Chuck Grip Capa. φ mm		dia φ 6~14mm Shafts, Flex-Coupling(φ 6~16mm)
“Push-Through” Shafts		dia φ 10mm or Less
Permissible Temp & Humidity		0°C~40°C, 80%
Dimension & Weight & Material		84 × 205 × h275mm, 3.5kg, Aluminum Die Casting
Power Consumption		250W
Electrical Reqts.		1 Phase, AC 120V, 60 Hz or AC 230V, 50 / 60 Hz

* Other Specifications are available upon Customer's Request

A. Unit-Only

DH.Sti00016 Hi-Torque Overhead Stirrer Unit “OVS-T16”, 230V
DH.Sti00016.1 , 120V

Ea. Price, ₩

B. Package-Set

- * Complete with ① “OVS-T16” Stirrer Unit
- ② “S732” Stainless-steel Impeller (Rod φ 8 × L500mm, 4-Blade φ 70mm)
- ③ “CL200” Clamp
- ④ “KA1177” Safety Stand (Rod φ 23 × L800mm, Base 34 × 29cm)

DH.Sti00116 Overhead Stirrer “OVS-T16” Package-set, 230V
DH.Sti00116.1 , 120V

Ea. Price, ₩

C. Optional Digital Remote Controller, with Connecting Cable 1.5m

- * Assembly-type, Torque / Viscosity / Temp-Control & Display
- * Connecting Cable 1.5m(Included) and 2m(Optional)
- * Size ; w75 × d150 × h30mm, Weight ; 150g

DH.WOS50010 Digital Remote Controller with 1.5m Connecting Cable & Plug
DH.CabI002 Optional Connecting Cable & Plug, 2m

Ea. Price, ₩

Optional Temperature Probe



(1) Direct Contact Temp Probe, K-type

DH.WHP503040 Temp Probe, K-type, “SS300”, L250 × φ 4.8 mm
DH.WHP503050 , “SS400”, L380 × φ 4.8 mm
DH.WHP503041 , “SS310”, Teflon Coated, L250 × φ 5 mm
DH.WHP503051 , “SS410”, Teflon Coated, L380 × φ 5 mm

Ea. Price, ₩

(2) K-type PTFE Encapsulated Stainless-steel Thermocouple Probe

※ Note! Optional Probes ⇒ See the “Measuring, Thermometers, K-type” Section.
Optional Stirrer Shafts ⇒ See the “Stirrer, Accessories” Section.



(1) Direct Contact Temp Probe, K-type (Optional)


(2) K-type PTFE Encapsulated Stainless-steel Thermocouple Probe (Optional)

DAIHAN™ **Premium Ultra** Hi-Torque Overhead Stirrer “OVS-G30”, with Permanently Brushless Motor(BLDC), Max. 60:1, 350,000mPas With Planetary-gear, Torque(Ncm) · Viscosity(mPas) · Temperature(°C) · Real Time Display, Optional Remote Control, Max. 550rpm Chuck Grip Capa. $\phi 6\sim 14$ mm Shafts, 프리미엄 초고점도용 교반기, 고사양 BLDC-Motor, 토크 · 점도 · 온도 표시기능




C. Remote Controller (Optional)

A. Max. 100Lit (H₂O) **Premium Ultra** “OVS-G30” Unit Only



3-Year Warranty

B. Max. 100Lit (H₂O) **Premium Ultra** “OVS-G30” Package-Set



(1) Direct Contact Temp Probe, K-type (Optional)

(2) K-type PTFE Encapsulated Stainless-steel Thermocouple Probe (Optional)

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- 다양한 점도 교반을 위해 $\phi 6\sim 14$ mm 의 Impeller Shaft가 장착 가능한 Chuck Grip.
- Patented Jog Dial & Touch Button
- Hi-Capacity Brushless Motor(BLDC) ; High Efficiency Stirring Performance and Precise Speed Control(Minimize RPM Fluctuation), Low Maintenance Cost
- Torque Value(Ncm)-/Viscosity(mPas)-/Temperature(°C) Real-time Display
- Can be used with K-type Temp. Probe for Temperature Display.
- Excellent Planetary Gear Design for High Torque(60:1 Ratio)
- Wide Range Chuck Grip Capa. $\phi 6\sim 14$ mm for High Viscosity.
- Can be used with Optional Wire Digital Remote Controller.
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System

* Ideal for Low-/Middle-Viscosity

* Smooth Start and Change of the Stirring Speed.

* Large Size Digital LCD Display

Specifications & Ordering information

<L1>

Spec.	Model	“OVS-G30” Max. Torque 350,000mPas
Stirring Capacity		Max. 100Lit (H ₂ O)
Maximum Viscosity / Torque		350,000 mPas / 800 Ncm
Speed Range & Motor		10~550RPM, Hi-Capacity Brushless Motor(BLDC), Gear(11:1 / 60:1)
Motor Rating (Input / Output)		200 / 180W
Controller / Safety Device		Digital PID Control / Over Heat & Current Protector for Motor
Chuck Grip Capa. ϕ mm		dia $\phi 6\sim 14$ mm Shafts, Flex-Coupling($\phi 6\sim 16$ mm)
Permissible Temp. & Humidity		0°C~40°C, 80%
Dimension & Weight & Material		84×205×h355mm, 3.5kg, Aluminum Die Casting
Power Consumption		250W
Electrical Reqts.		1 Phase, AC 120V, 60 Hz or AC 230V, 50 / 60 Hz

* Other Specifications are available upon Customer's Request

A. Unit-Only

DH.StiO0030 Ultra Hi-Torque Overhead Stirrer Unit “OVS-G30”, 230V
DH.StiO0030.1 , 120V

Ea. Price, ₩

B. Package-Set

- * Complete with ① “OVS-G30” Stirrer Unit
- ② “S732” Stainless-steel Impeller (Rod $\phi 8 \times L500$ mm, 4-Blade $\phi 70$ mm)
- ③ “CL200” Clamp
- ④ “ST200” Safety Stand (Rod $\phi 23 \times L1,200$ mm, Base 34×29 cm)

DH.StiO0130 Overhead Stirrer “OVS-G30” Package-set, 230V
DH.StiO0130.1 , 120V

Ea. Price, ₩

C. Optional Digital Remote Controller, with Connecting Cable 1.5m

- * Assembly-type, Torque / Viscosity / Temp-Control & Display
- * Connecting Cable 1.5m(Included) and 2m(Optional)
- * Size ; w75×d150×h30mm, Weight ; 150g

DH.WOS50010 Digital Remote Controller with 1.5m Connecting Cable & Plug
DH.Cabl002 Optional Connecting Cable & Plug, 2m

Ea. Price, ₩

Optional Temperature Probe

(1) Direct Contact Temp Probe, K-type

DH.WHP503040 Temp Probe, K-type, “SS300”, L250× $\phi 4.8$ mm
DH.WHP503050 , “SS400”, L380× $\phi 4.8$ mm
DH.WHP503041 , “SS310”, Teflon Coated, L250× $\phi 5$ mm
DH.WHP503051 , “SS410”, Teflon Coated, L380× $\phi 5$ mm

Ea. Price, ₩

(2) K-type PTFE Encapsulated Stainless-steel Thermocouple Probe

※ Note! Optional Probes ⇒ See the “Measuring, Thermometers, K-type” Section.
Optional Stirrer Shafts ⇒ See the “Stirrer, Accessories” Section.

DAIHAN™ High-Torque Digital Overhead Stirrer “HT-DX”, for Middle-/High-Viscosity, Max.40 & 60Lit, 60,000 & 90,000mPas With “Push-Through” Shaft(ϕ 10mm or Less) and Chuck grip ϕ 6~12mm, Flex-Coupling ϕ 6~16mm, Max.1,000rpm
중/고점도용 강력교반기, 사용교반용: ①기본은 ϕ 6~12mm, ②관통은 ϕ 6~10mm, & ③별도의 Flex-Coupling 활용시 ϕ 6~16mm



B. Max. 60Lit Digital “HT-120DX” Package-Set



HT-120DX

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- Ensure Minimize RPM Fluctuation ; Very Stable Stirring Speed(RPM) Performance
- Patented Jog-Shuttle Control System
- * with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- * “Push-Through” Shafts(ϕ 10mm or Less) ; Ideal for Raising or Lowering of the Stir-Shafts without Moving the Motor Head - “Adjustable Height of Shaft”
- * Digital LCD with Back-Light Function
- * Speed Range ; 50 ~ 1,000 rpm, Control and Display by 1 rpm
- * Timer ; 99hr 59min (with Continuous run), Error status & Timer-end Alarm

Digital Stirrer Package-Set

- * A. “HT-50DX” Set ; Complete with ① “HT-50DX” Stirrer Unit, ② “S732” Stainless-steel Impeller (Rod ϕ 8×L500mm, 4-Blade ϕ 70mm), ③ “CL200” Clamp, and ④ “ST120” Wide Plate Stand (Rod ϕ 23×L1000mm, Base w30×d40cm)
- * B. “HT-120DX” Set ; Complete with ① “HT-120DX” Stirrer Unit, ② “S736” Stainless-steel Impeller (Rod ϕ 10×L650mm, 4-Blade ϕ 90mm), ③ “CL200” Clamp, and ④ “ST200” Deluxe Safety Stand (Rod ϕ 23×L1200mm, Base w39×up to 75cm)

Specifications & Ordering information

<L1>

Spec.	Model	A. “HT-50DX” Max. Stirring Capa. 40Lit	B. “HT-120DX” Max. Stirring Capa. 60Lit
Stirring Capacity		Max. 40Lit (H ₂ O)	Max. 60Lit (H ₂ O)
Maximum Viscosity / Torque		60,000 mPas / 180 Ncm	90,000 mPas / 210 Ncm
Motor		Worm-Geared DC Motor (3.5:1 Ratio)	
Speed & Control Resolution		50~1000 rpm, 1 rpm	
Motor Rating (Input / Output)		83 / 75W	130 / 120W
Chuck Grip Capa. ϕ mm		dia ϕ 6~12mm Shafts	
“Push-Through” Shafts		dia ϕ 10mm or Less	
Permissible Temp. & Humidity		0℃~40℃, 80%	
Safety Device / Material		Overload & Heat Protector for Motor / ABS and Aluminum	
Dimension & Weight		90×180×h130mm, 3.5kg	90×210×h130mm, 4kg
Packing Size & Gross Weight		320×300×h180mm, 5kg	320×300×h180mm, 8kg
Power Consumption		85W	130W
Electrical Reqt.		1 Phase, AC 120V, 60 Hz or AC 230V, 50 / 60 Hz	

Digital Stirrer Only


Cat. No.	230V	DH.WOS01051	DH.WOS01071
	120V	DH.WOS02051	DH.WOS02071
Price, Each. ₩			

Digital Package Set (with ① Stirrer Unit, ② Impeller, ③ Clamp, ④ Stand)

Cat. No.	230V	DH.WOS01052	DH.WOS01072
	120V	DH.WOS02052	DH.WOS02072
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

DAIHAN™ Remote Control High-Torque Digital Overhead Stirrer “HT-T”, for Middle-/High-Viscosity, Max. 40 · 60 · 100Lit
 With “Push-Through” Shaft(ϕ 10mm or Less) and Chuck grip ϕ 6~12mm, Flex-Coupling ϕ 6~16mm, 120,000mPas, Max. 1,000rpm
 원격조절식 디지털 중/고점도용 강력교반기, 사용교반용 ; ①기본은 ϕ 6~12mm, ②관통은 ϕ 6~10mm, & ③별도의 Flex-Coupling 활용시 ϕ 8~16mm



Digital High-Torque Overhead Stirrer
 with Stirring Speed Remote Controller
 “HT-120T” Package-Set



“HT-140T” Digital High-Torque Overhead Stirrer
 (included Stirring Speed Remote Controller and Temp. Display)
 in Modular Reactor System

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- Ensure Minimize RPM Fluctuation ; Very Stable Stirring Speed(RPM) Performance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Patented Jog-Shuttle Control System

- * “Push-Through” Shafts(ϕ 10mm or Less) ; Ideal for Raising or Lowering of the Stir-Shafts without Moving the Motor Head - “Adjustable Height of Shaft”
- * Included Remote-type Digital Feedback Controller ;
- “HT-50T”, “HT-120T” Model ; Included Stirring-speed Remote Controller
- “HT-140T” Model ; Ideal for Reaction System.
- Included Stirring-speed Control and Temperature Display Remote Controller for Reactor Module.

Specifications & Ordering information

<L1>

Spec. / Model	“HT-50T” Max. Stirring Capa. 40L	“HT-120T” Max. Stirring Capa. 60L	“HT-140T” Max. Stirring Capa. 100L
Stirring Capacity(H ₂ O)	Max. 40Lit	Max. 60Lit	Max. 100Lit
Maximum Viscosity / Torque	60,000 mPas / 180 Ncm	90,000 mPas / 210 Ncm	120,000 mPas / 300 Ncm
Speed Range & Motor	50~1,000 rpm		50~800 rpm
	Worm-Geared Motor(3.5:1 Ratio)		
Motor Rating (Input / Output)	83 / 75W	130 / 120W	
Timer & Alarm	99hr 59min (with Continuous run.), Timer-end Alarm		
Controller / Display	Digital Feedback Remote Controller / Digital LCD with Back-Light Function		
	Stirring Speed Remote Controller		Stirring Speed & Temp. Display Remote Controller
Max. Impeller Shaft Diameter	ϕ 10mm		
Permissible Temp. & Humidity	0℃~40℃, 80%		
Safety Device	Overload & Heat Protector for motor		
Material	ABS and Aluminum		
Dimension & Weight*	155×180×h110mm, 3.5kg	155×230×h110mm, 4kg	
Packing Size & Gross Weight	580×340×h230mm, 5kg	580×340×h230mm, 8kg	
Power Consumption	75W	120W	
Electrical Reqts.	1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz		

Unit & Remote Controller

Cat. No.	230V	DH.WOS01065	DH.WOS01075	DH.WOS01085
	120V	DH.WOS02065	DH.WOS02075	DH.WOS02085
Price, Each. ₩				

Package-Set (with ① Stirrer+Remote Controller, ② Impeller, ③ Clamp, ④ Stand)

Cat. No.	230V	DH.WOS01066	DH.WOS01076	-
	120V	DH.WOS02066	DH.WOS02076	-
Price, Each. ₩				-


* Weight of Remote Controller = 2.5kg

** Other Specifications are available upon Customer's Request.

- * “HT-50T” Package-Set ; Complete with ① “HT-50T” Stirrer Unit with Remote Controller ② “S732” Stainless-steel Impeller (Rod ϕ 8×L500mm, 4-Blade ϕ 70mm) ③ “CL200” Clamp, and ④ “ST120” Wide Plate Stand(Rod ϕ 23×L1000mm, Base w30×d40cm)
- * “HT-120T” Package-Set ; Complete with ① “HT-120T” Stirrer Unit with Remote Controller ② “S736” Stainless-steel Impeller (Rod ϕ 10×L650mm, 4-Blade ϕ 90mm) ③ “CL200” Clamp, and ④ “ST200” Deluxe Safety Stand (Rod ϕ 23×L1,200mm, Base 39×up to 75cm)

※ Note! Details of Wide-Range Coupling ; See the “Stirrer, Accessories” Section.

DAIHAN™ Digital Jar Tester “JT-M6C”, 6-Paddle Stirrer for Flocculation Test, 300 rpm, with Certi. & Traceability
With Programmable Digital Feedback Control, 6 Impellers, 6 Independent Light Source, Up to 6 × 1Lit Beakers
프로그램/디지털 자 테스터, 정밀 직류 기어드 모터와 특화된 타이밍벨트, 디지털 피드백 컨트롤



3-Year Warranty

Jar Tester, φ8×L400mm Shaft “JT-M6C”



Geared DC Motor with Special Timing Belt

- Supplied by DAIHAN Korea
- CE Certified
- PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Flocculation Testing in Water Treatment Laboratories
- 6 Independent Light Source for Each Stirring Point

- * Deciding on Coagulant Chemistries
- * Digital Display for RPM, Power Bar-Graph and Status
- * Digital Feedback Control using Microprocessor ; Stable and Accurate Stirring Speed
- * Geared DC Motor and Specially Designed Timing Belt used
- * Stable and Accurate Stirring Speed by Feedback Control
- * Digital LCD with Back-Light Function
- * Programmable Operation ; 10 Steps, 6 Programs

Specifications & Ordering information

<L1>

Spec.	Model	“JT-M6C”
6 × Paddler (Stainless steel Shaft)		φ8 × L400mm
		Straight Blade w75 × h25mm
Stirring Capacity		6 × 1 Liter Beaker
Motor & Power		Geared DC Motor with Timing Belt (10:1 ratio), 40W
Speed & Accuracy		20~300 rpm, ±5 rpm
Timer		99hr 59min 59sec (Programmable)
Display		Digital LCD with Back-Light for RPM, Power Bar-Graph and Status
Controller		Programmable Digital Feedback Controller with Jog-Shuttle Switch (Turn + Push)
Program Function		10 Steps, 6 Programs, and Storage for Set Values ; RPM and Timer
Florescent Lamp		20W × 6 (Independent Light Source for Each Sample)
Dimension (w × d × h)		1,000 × 310 × h535mm
Safety Device		Overload & Over Heat Protector for motor
Others		PL(Product Liability) Insurance and Storage Function (RPM, Timer)
Packing and Gross Weight		1,120 × 430 × h740mm, 69 kg
Electrical Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz
Cat. No.	230V	DH.WJT06020
Price, Each. ₩	120V	DH.WJT06021

※ 운임 별도

Vacuum Stirrer-Guide, with DURAN® PBT Opentop Cap & PTFE/ Silicone O-Ring Seal and 24/- ~ 45/- Glass Joints
For ϕ 8~12mm Shafts, with Red PBT Screw Grip System, Chemically Inert, 200°C Stable, 스크류 스티러 쉘



■ Supplied by SciLab Korea ■ Custom-Made Available / 규격외 Items 자유롭게 주문하세요!

(1) With ASTM 24/40~45/50 Joint

SL.StiS241	for ϕ 8 mm Rod, GL25, 24/40
SL.StiS242	, 29/42
SL.StiS243	, 34/45
SL.StiS245	for ϕ 10mm Rod, GL25, 24/40
SL.StiS246	, 34/45
SL.StiS247	, GL32, 34/45
SL.StiS251	for ϕ 12mm Rod, GL32, 34/45
SL.StiS252	, 45/50

<L1>

Ea. Price, ₩

(2) With DIN 24/29~45/40 Joint

SL.StiS241D	for ϕ 8 mm Rod, GL25, 24/29
SL.StiS242D	, 29/32
SL.StiS243D	, 34/35
SL.StiS245D	for ϕ 10mm Rod, GL25, 24/29
SL.StiS246D	, 34/35
SL.StiS247D	, GL32, 34/35
SL.StiS251D	for ϕ 12mm Rod, GL32, 34/35
SL.StiS252D	, 45/40

Ea. Price, ₩



SciLab™ Stirrer-Guide/Seals, PTFE, for ϕ 8~ ϕ 14mm Shafts, with 24/- & 34/-Joints

With 2 Bearings / Screw-Retainer / Viton O-ring, -180°C +260°C Stable, <Korea-Made> PTFE 스티러 쉘

(1) for ϕ 8mm Stir-Shafts

SL.Ada405	for ϕ 8mm Rod, 24/40
SL.Ada407	, 34/45

<L1>

Ea. Price, ₩

(2) for ϕ 10mm Stir-Shafts

SL.Ada409	for ϕ 10mm Rod, 24/40
SL.Ada411	, 34/45


Ea. Price, ₩

(3) for ϕ 12 & 14mm Stir-Shafts

SL.Ada413	PTFE Stir-Guide, for ϕ 12mm Rod, 34/45
SL.Ada415	, for ϕ 14mm Rod, 34/45

Ea. Price, ₩

Cowie® Premium Taper Joint-PTFE/PEEK Vacuum Stirrer Guide, with ASTM/DIN Glass Cone Joint, for ϕ 8~16mm Shafts
Ideal for Middle/Low-Vacuum, Max-800rpm, Chemically Inert, -200°C +280°C Stable, <UK-Made> 고품질 진공 조인트 교반 쉘



High Performance Shaft Guides

■ Suitable for Glass Same Size Socket of ASTM or DIN

* with an Adjustable / Replaceable Compression Seal Sizes ; for ϕ 8~ ϕ 16mm Shafts

* PTFE in the only Media Contact Material

(1) for ASTM Socket Joint, for ϕ 8~16mm Shaft

CW.005.0.0824	Stirrer Seal, USA-type, for- ϕ 8mm Shafts, 24/40 Joint
CW.005.0.1024	- ϕ 10mm Shafts, 24/40 Joint
CW.005.0.1034	- ϕ 10mm Shafts, 34/45 Joint
CW.005.0.1434	- ϕ 14mm Shafts, 34/45 Joint
CW.005.0.1634	- ϕ 16mm Shafts, 34/45 Joint

<L1>

Ea. Price, ₩

(2) for DIN/ISO Socket Joint, for ϕ 8 & 10mm Shaft

CW.005.00.0824	Stirrer Seal, EUROPE-type, for- ϕ 8mm Shafts, 24/29 Joint
CW.005.00.1029	- ϕ 10mm Shafts, 29/32 Joint
CW.005.00.1045	- ϕ 10mm Shafts, 45/40 Joint

Ea. Price, ₩

Cowie® High-Vacuum/Pressure Stirrer Guide, PTFE/PEEK & Glass-Bearing Seal

For ϕ 8~16mm Shafts, Useable up to 800rpm, 280°C Stable, <UK-Made> 고진공/압력용 교반 쉘

■ High Vacuum 5mmHg & Pressure 3~5psi Performance

* Anti-Whip, Reduced Vibration, No-Shedding, Self-Releasing Joint Ring, Totally Inert.

* Specially Fabricated Glass Ball-Bearing & PTFE/PEEK Seal ; Useable up to 800rpm

(1) for ASTM Socket Joint Sizes ; for ϕ 8~16 mm Shafts

CW.005.1047	Stirrer Seal, USA-type, for- ϕ 8mm Shafts, 24/40 Joint
CW.005.1057	- ϕ 10mm Shafts, 24/40 Joint
CW.005.1067	- ϕ 10mm Shafts, 29/42 Joint
CW.005.1097	- ϕ 14mm Shafts, 34/45 Joint
CW.005.1117	- ϕ 16mm Shafts, 34/45 Joint


<L1>

Ea. Price, ₩



(2) for DIN/ISO Socket Joint Sizes ; for ϕ 10~16 mm Shafts

CW.005.2067	Stirrer Seal, EUROPE-type, for- ϕ 10mm Shafts, 29/32 Joint
CW.005.2077	- ϕ 10mm Shafts, 34/35 Joint
CW.005.2097	- ϕ 14mm Shafts, 34/35 Joint
CW.005.2127	- ϕ 16mm Shafts, 45/40 Joint

Ea. Price, ₩



for High Vacuum / Pressure



Flexible Coupling, for Safety Stirring of $\phi 6 \sim 16$ mm Shafts of Overhead Stirrers

For Stirrers between Stirring Motor and Shaft, 안전교반용 유연성 커플링, 교반축의 불균형등을 완전 보정

■ Supplied by DAIHAN Korea

■ Used to Compensate Any Vibrations Between The Motor and Stirring Shafts

* Made of DURALUMIN * Shaft Fastening by Screw Clamp Band

* with Exchangeable STS Chucking Rod $\phi 8 / 10 / 12$ mm, L40mm

(1) UTILITY High-Flex. Metal Coupling, with 2 $\times$ Flex-band, for $\phi 6 \sim 16$ mm Stir-Shafts Ideal for Low/High Viscosity and High/Low Speed

<L1>

Ea.	Price, ₩
SL.StiC011	Stir-Coupling, Hi-flex, $\phi 8$ mm Chuck-rod, for- $\phi 6$ mm Stir-Shaft, 56g, $\phi 25 \times h44$ mm
SL.StiC012	- $\phi 8$ mm Stir-Shaft, 53g, $\phi 25 \times h44$ mm
SL.StiC013	- $\phi 10$ mm Stir-Shaft, 49g, $\phi 25 \times h44$ mm
SL.StiC014	, $\phi 10$ mm Chuck-rod, for- $\phi 8$ mm Stir-Shaft, 49g, $\phi 25 \times h44$ mm
SL.StiC015	- $\phi 10$ mm Stir-Shaft, 43g, $\phi 25 \times h44$ mm
SL.StiC016	- $\phi 12$ mm Stir-Shaft, 33g, $\phi 25 \times h44$ mm
SL.StiC017	, $\phi 12$ mm Chuck-rod, for- $\phi 10$ mm Stir-Shaft, 33g, $\phi 25 \times h44$ mm
SL.StiC018	- $\phi 12$ mm Stir-Shaft, 123g, $\phi 32 \times h60$ mm
SL.StiC019	- $\phi 14$ mm Stir-Shaft, 120g, $\phi 32 \times h60$ mm
SL.StiC020	- $\phi 16$ mm Stir-Shaft, 118g, $\phi 32 \times h60$ mm


(2) Heavy-duty Metal Coupling, with Short Flex-band, for $\phi 10 \sim 16$ mm Stir-Shafts Low-Flexible type, Ideal for High Viscosity and Middle/Low Speed

Ea. Price, ₩

Ea.	Price, ₩
SL.StiC023	Stir-Coupling, Low-flex, $\phi 10$ mm Chuck rod, for- $\phi 10$ mm Stir-Shafts., 45g, $\phi 25 \times h30$
SL.StiC024	- $\phi 12$ mm Stir-Shafts., $\phi 25 \times h30$
SL.StiC025	, $\phi 12$ mm Chuck rod, for- $\phi 12$ mm Stir-Shafts., 105g, $\phi 32 \times h38$
SL.StiC026	- $\phi 14$ mm Stir-Shafts., $\phi 32 \times h38$
SL.StiC027	- $\phi 16$ mm Stir-Shafts., $\phi 32 \times h38$

Bochem® Stainless-steel Centrifugal-type Stirring Shaft, for Lab & Industrial Overhead Stirrers, Rod $\phi 8$ & $10 \times L400 \sim 650$ mm

Ideal for Low/Middle Viscosity, with 2 Flexible Blades, Non-magnetic/Rust-free, <Germany-Made> 스텐레스 교반봉 / 임펠러



■ Good Finished ; Surface & Stirring Center

■ Recommended to use with Flexible Coupling for Gentle Stirring & Safety Working

* 18/10 STS ; Non-magnetic, Rust-free, Anti-Acid / Chemicals / Heat

<L1>

Ea.	Price, ₩
BC.9210	Stir-Shaft, STS, Rod- $\phi 8 \times L400$ mm, Swingrotor- $\phi 90 \times h10$ mm
BC.9211	, Rod- $\phi 8 \times L500$ mm, Swingrotor- $\phi 90 \times h10$ mm
BC.9213	, Rod- $\phi 10 \times L400$ mm, Swingrotor- $\phi 90 \times h10$ mm
BC.9214	, Rod- $\phi 10 \times L500$ mm, Swingrotor- $\phi 90 \times h10$ mm
BC.9215	, Rod- $\phi 10 \times L650$ mm, Swingrotor- $\phi 90 \times h10$ mm

SciLab™ Stainless-steel Stirring Shaft, for Lab & Industrial Overhead Stirrers, Rod $\phi 8 \sim 12 \times L500 \sim 1000$ mm

Ideal for Universal-use of Wide-range Speed & Viscosity, <Korea-Made> 스텐레스 교반봉/임펠러



(1) 3-bladed Propeller-type
Universal-use



(2) 4-bladed Propeller-type
Standard-model

■ Supplied by SciLab Korea

■ Good Finished ; Surface & Stirring Center

■ Recommended to use

with Flexible Coupling for Gentle Stirring & Safety Working

■ The Rod may be Cut to Any Length with a Hacksaw

(1) 3-bladed Propeller-type, ideal for Universal-use in Middle/High Speed

with Rod $\phi 8 \sim 12 \times L500 \sim 1000$, Rotor $\phi 50 \sim 140$ mm, and Good for Low/Middle Viscosity

<L1>

* Designed for Efficient Flow ; mixed from Top to Bottom of Vessel

* Minimum Shearing

Ea.	Price, ₩
SL.StiS721	Stir-Shaft, "S721", #304, Rod- $\phi 8 \times L500$, Rotor- $\phi 50 \times h15$ mm
SL.StiS722	, "S722", - $\phi 8 \times L500$, - $\phi 60 \times h15$ mm
SL.StiS723	, "S723", - $\phi 8 \times L650$, - $\phi 80 \times h20$ mm
SL.StiS724	, "S724", - $\phi 8 \times L650$, - $\phi 100 \times h24$ mm
SL.StiS726	Stir-Shaft, "S726", #304, Rod- $\phi 10 \times L500$, Rotor- $\phi 60 \times h15$ mm
SL.StiS727	, "S727", - $\phi 10 \times L650$, - $\phi 80 \times h20$ mm
SL.StiS728	, "S728", - $\phi 10 \times L800$, - $\phi 100 \times h24$ mm
SL.StiS729	Stir-Shaft, "S729", #304, Rod- $\phi 12 \times L1000$, Rotor- $\phi 140 \times h30$ mm

(2) Standard 4-bladed Propeller-type, ideal for Low/Middle Speed & Viscosity

Popular-model, with Rod $\phi 8 \sim 12 \times L500 \sim 1000$, Rotor $\phi 50 \sim 140$ mm, 45° Angled

Ea. Price, ₩

Ea.	Price, ₩
SL.StiS731	Stir-Shaft, "S731", #304, Rod- $\phi 8 \times L500$, Rotor- $\phi 50 \times h20$ mm
SL.StiS732	, "S732", - $\phi 8 \times L500$, - $\phi 70 \times h25$ mm
SL.StiS733	, "S733", - $\phi 8 \times L650$, - $\phi 90 \times h35$ mm
SL.StiS734	, "S734", - $\phi 8 \times L650$, - $\phi 110 \times h35$ mm
SL.StiS736	Stir-Shaft, "S736", #304, Rod- $\phi 10 \times L650$, Rotor- $\phi 90 \times h35$ mm
SL.StiS737	, "S737", - $\phi 10 \times L800$, - $\phi 110 \times h40$ mm
SL.StiS738	, "S738", - $\phi 10 \times L1000$, - $\phi 110 \times h40$ mm
SL.StiS739	Stir-Shaft, "S739", #304, Rod- $\phi 12 \times L1000$, Rotor- $\phi 140 \times h40$ mm

<Continued On Next Page>

<Continued From...SciLab™ Stainless-steel Stirring Shaft>



(3) 4-bladed Dissolver Diamond-type

(3) 4-bladed Dissolver Diamond-type, Good for Particle Reduction in Middle/High Speed with Rod $\phi 8$ & $10 \times L500 \sim 650$, Rotor $\phi 70 \sim 110$ mm, and Good for Low/Middle Viscosity

<L1>

SL.StiS741	Stir-Shaft, "S741", #304, Rod- $\phi 8 \times L500$, Rotor- $\phi 70 \times h20$ mm
SL.StiS742	, "S742", - $\phi 8 \times L500$, - $\phi 90 \times h25$ mm
SL.StiS743	, "S743", - $\phi 8 \times L650$, - $\phi 110 \times h25$ mm
SL.StiS744	Stir-Shaft, "S744", #304, Rod- $\phi 10 \times L650$, Rotor- $\phi 110 \times h25$ mm

Ea. Price, ₩



(4) Heavy-duty Anchor-type

(4) Heavy-duty Anchor-type, ideal for Middle/High Viscosity, Rod $\phi 8 \sim 12$ mm Good for Low/Middle Speed, Rotor $\phi 70 \sim 180$ mm

SL.StiS751	Stir-Shaft, "S751", #304, Rod- $\phi 8 \times L500$, Rotor- $\phi 70 \times h14$ mm
SL.StiS752	, "S752", - $\phi 8 \times L650$, - $\phi 90 \times h16$ mm
SL.StiS753	Stir-Shaft, "S753", #304, Rod- $\phi 10 \times L650$, Rotor- $\phi 90 \times h16$ mm
SL.StiS754	, "S754", - $\phi 10 \times L800$, - $\phi 140 \times h25$ mm
SL.StiS755	, "S755", - $\phi 10 \times L1000$, - $\phi 160 \times h30$ mm
SL.StiS757	Stir-Shaft, "S757", #304, Rod- $\phi 12 \times L1000$, Rotor- $\phi 180 \times h30$ mm

Ea. Price, ₩

Bochem® Stainless-steel Square-type Stirring Shaft, with $6 \times \phi 12$ mm Hole, for Lab & Industrial Overhead Stirrers, Rod $\phi 8$ & $10 \times L400$ & 500 mm
Ideal for Low/Middle Viscosity, Non-magnetic/Rust-free, <Germany-Made> 스텐레스 교반봉/임펠러



Square-type

■ Good Finished ; Surface & Stirring Center

■ Recommended to Use with a Flexible Coupling for Gentle Stirring & Safety Working


* 18/10 Stainless-steel; Non-magnetic, Rust-free, Anti-Acid / Chemicals / Heat

BC.9504	Stir-Shaft, STS, Rod- $\phi 8 \times L400$ mm, Rotor- $\phi 70 \times h65$ mm $\times t1.7$ mm, $6 \times \phi 12$ mm Hole
BC.9505	, Rod- $\phi 8 \times L500$ mm, Rotor- $\phi 70 \times h65$ mm $\times t1.7$ mm, $6 \times \phi 12$ mm Hole
BC.9507	, Rod- $\phi 10 \times L400$ mm, Rotor- $\phi 70 \times h65$ mm $\times t1.7$ mm, $6 \times \phi 12$ mm Hole
BC.9508	, Rod- $\phi 10 \times L500$ mm, Rotor- $\phi 70 \times h65$ mm $\times t1.7$ mm, $6 \times \phi 12$ mm Hole

<L1>

Ea. Price, ₩

Stainless-steel Heavy-duty Helical Ribbon-type, for Middle/High Viscosity and Lab~Industry
With Rod $\phi 8 \sim 12$ mm, $L500 \sim 1100$, Rotor $\phi 60 \sim 190$ mm, Good for Reaction-system



Helical Ribbon-type

■ Supplied by DAIHAN Korea

* Patented Helical Ribbon High Efficient Mixing Performance * Generates Axial Flow.

* Designed for Middle / High Viscosity. High Flow Rate and Efficient Mixing.

* For Homogeneous Mixing and Medium is Conveyed from the Bottom upwards.

SL.Sti761	Stir-Shaft, "S761", #304, Rod- $\phi 8 \times L500$, Rotor- $\phi 60 \times h100$ mm
SL.Sti762	, "S762", - $\phi 8 \times L650$, - $\phi 80 \times h130$ mm
SL.Sti762.1	Stir-Shaft, "S762.1", #304, Rod- $\phi 10 \times L650$, Rotor- $\phi 100 \times h160$ mm
SL.Sti763	, "S763", - $\phi 10 \times L650$, - $\phi 115 \times h175$ mm
SL.Sti764	, "S764", - $\phi 10 \times L800$, - $\phi 125 \times h190$ mm
SL.Sti765	Stir-Shaft, "S765", #304, Rod- $\phi 12 \times L1000$, Rotor- $\phi 160 \times h240$ mm
SL.Sti766	, "S766", - $\phi 12 \times L1100$, - $\phi 190 \times h290$ mm

<L1>

Ea. Price, ₩



(1) (2) (3) (4)

Kartell® PP & PVC Plastic Stirring Shaft, PP-Blade, PVC-Rod with Steel Core, Rod $\phi 8 \times L350$ mm
Good for Acid & Chemicals resistance, <Italy-Made> PP & PVC 플라스틱 교반봉/임펠러

■ Ideal for Acids, Chemicals and Anti-corrosion

■ Material ; PP-Blade, PVC-Rod with Steel Core ; $\phi 8 \times L350$ mm* Polypropylene(PP) ; $-10^{\circ}\text{C} \sim 125/140^{\circ}\text{C}$ * Polyvinylchloride(PVC) ; $-20^{\circ}\text{C} \sim 80^{\circ}\text{C}$



KA.428	(1) Stirring Paddle, with 2 fixed blade (w39 $\times$ h15mm), Rod $\phi 8 \times L350$ mm
KA.430	(2) , "U" shaped blade (w65 $\times$ h78mm), Rod $\phi 8 \times L350$ mm
KA.431	(3) , with 2 flexible paddle (Swing w57mm), Rod $\phi 8 \times L350$ mm
KA.432	(4) , with 2 flexible paddle (Swing w98mm), Rod $\phi 8 \times L350$ mm

<L1>

Ea. Price, ₩

PTFE Stirring Shaft, with Stainless-steel Insert, for Lab & Industrial Overhead Stirrers, Rod $\phi 8$ & $10 \times L300 \sim 650\text{mm}$
With Centrifugal & Square-type, High Temperature Resistant and Corrosion-Proof, $-200^{\circ}\text{C} + 260^{\circ}\text{C}$, Normal-grade, PTFE 교반봉/임펠러

- Supplied by SciLab Korea ■ Recommended to Use with a Flexible Coupling for Gentle Stirring & Safety Working
- The PTFE Rod with Stainless-steel Insert ; ① OD $\phi 8\text{mm}$ with $\phi 6\text{mm}$ Steel Insert-" $\phi 8/6$ "
② OD $\phi 10\text{mm}$ with $\phi 8\text{mm}$ Steel Insert-" $\phi 10/8$ "
- Ideal for PURE-Working / Chemically Inert / Non-Contaminating / Thermally Stable
- Exposed Top Rod 4~5cm (Stainless-steel Core) ; For Strong Clamping in the Stir-Chuck or Flexible Coupling. ⇒



(1) 2-bladed Centrifugal-type

SL.Sti4091	PTFE Stirring Rod, $\phi 8/6 \times L300\text{mm}$, Swing rotor- $\phi 63\text{mm} \times h24\text{mm}$
SL.Sti4092	, $\phi 8/6 \times L400\text{mm}$, - $\phi 80\text{mm} \times h25\text{mm}$
SL.Sti4093	, $\phi 8/6 \times L500\text{mm}$, - $\phi 100\text{mm} \times h30\text{mm}$
SL.Sti4094	, $\phi 8/6 \times L600\text{mm}$, - $\phi 110\text{mm} \times h34\text{mm}$
SL.Sti4095	, $\phi 10/8 \times L650\text{mm}$, - $\phi 110\text{mm} \times h34\text{mm}$

Ea. Price, ₩


(2) Square-type, with $\phi 12\text{mm}$ Three Holes

SL.Sti4098	PTFE Stirring Rod, $\phi 10/8 \times L350\text{mm}$, Rotor $70 \times 70\text{mm} \times t4\text{mm}$, $3 \times \phi 12\text{mm}$ Hole
SL.Sti4099	, $\phi 10/8 \times L450\text{mm}$, Rotor $70 \times 70\text{mm} \times t4\text{mm}$, $3 \times \phi 12\text{mm}$ Hole

Ea. Price, ₩

SciLab™ Premium PTFE Stirring Shaft, with Stainless-steel Insert, for Overhead Stirrers, Rod $\phi 8 \sim 14 \times L500 \sim 1100\text{mm}$
Ideal for PURE Working & Universal-use of Wide-range Speed & Viscosity, Chemically Inert, $-200^{\circ}\text{C} + 280^{\circ}\text{C}$ -stable, PTFE 교반봉/임펠러

- Recommended to Use with a Flexible Coupling for Gentle Stirring & Safety Working.
- The PTFE Rod with Stainless-steel Insert ; ① OD $\phi 8\text{mm}$ with $\phi 6\text{mm}$ Steel Insert-" $\phi 8/6$ ",
② OD $\phi 10\text{mm}$ with $\phi 8\text{mm}$ Steel Insert-" $\phi 10/8$ ", ③ OD $\phi 12\text{mm}$ with $\phi 10\text{mm}$ Steel Insert-" $\phi 12/10$ ",
and ④ OD $\phi 14\text{mm}$ with $\phi 12\text{mm}$ Steel Insert-" $\phi 14/12$ ". ■ Rod may be Cut to Any Length with a Hacksaw.
- Ideal for PURE-Working / Chemically Inert / Non-Contaminating / Thermally Stable, $-200^{\circ}\text{C} + 280^{\circ}\text{C}$ -stable.
- Exposed Top Rod 4~5cm (Stainless-steel Core) ; For Strong Clamping in the Stir-Chuck or Flexible Coupling. ⇒



(1) 2-bladed Centrifugal-type
Good for Narrow-Neck
<ex> 24/- & 34/-Joint neck

(1) 2-bladed Centrifugal-type, ideal for 24/- & 34/-Joints(Narrow) Necks & Low Viscosity with Rod $\phi 8$ & $10 \times L500 \sim 650$, Rotor $\phi 50 \sim 100\text{mm}$, and Good for Middle/High Speed

SL.StiS403	Stir-Shaft, PTFE, Rod- $\phi 8/6 \times L500$, Swingrotor- $\phi 50 \times h14\text{mm}$, fit to 24/-joints
SL.StiS404	- $\phi 8/6 \times L500$, - $\phi 70 \times h17\text{mm}$, fit to 24/-joints
SL.StiS405	- $\phi 8/6 \times L650$, - $\phi 100 \times h24\text{mm}$, fit to 34/-joints
SL.StiS406	Stir-Shaft, PTFE, Rod- $\phi 10/8 \times L650$, Swingrotor- $\phi 100 \times h24\text{mm}$, fit to 34/-joints

Ea. Price, ₩

(2) Standard 4-bladed Propeller-type, ideal for Low/Middle Speed and High Viscosity Good for Reactors, with Rod $\phi 8 \sim 14 \times L500 \sim 1100$ and Rotor $\phi 50 \sim 140\text{mm}$, 45° Angled

* Standard Stirring Shafts. * Designed for Efficient Mixing ; from Top to Bottom of Vessel

SL.StiS410	Stir-Shaft, PTFE, Rod- $\phi 8/6 \times L500$, Rotor- $\phi 50 \times h20\text{mm}$
SL.StiS411	- $\phi 8/6 \times L500$, - $\phi 70 \times h25\text{mm}$
SL.StiS412	- $\phi 8/6 \times L650$, - $\phi 70 \times h25\text{mm}$
SL.StiS413	Stir-Shaft, PTFE, Rod- $\phi 10/8 \times L650$, Rotor- $\phi 90 \times h35\text{mm}$, Good for 5~10Lit Vessels
SL.StiS414	- $\phi 10/8 \times L800$, - $\phi 90 \times h35\text{mm}$, Good for 10~20Lit Vessels
SL.StiS415	- $\phi 10/8 \times L1000$, - $\phi 90 \times h35\text{mm}$, Good for 20~50Lit Vessels
SL.StiS416.1	Stir-Shaft, PTFE, Rod- $\phi 12/10 \times L800$, Rotor- $\phi 110 \times h40\text{mm}$, Good for 10~20Lit Vessels
SL.StiS416.2	- $\phi 12/10 \times L1000$, - $\phi 110 \times h40\text{mm}$, Good for 30~50Lit Vessels
SL.StiS416.3	- $\phi 12/10 \times L1000$, - $\phi 140 \times h40\text{mm}$, Good for 30~80Lit Vessels
SL.StiS417	Stir-Shaft, PTFE, Rod- $\phi 14/12 \times L1100$, Exchangeable Rotor- $\phi 140 \times h40\text{mm}$, Good for 50~100Lit Reactors

Ea. Price, ₩



(2) Standard 4-bladed
Propeller-type

(3) with 2-Exchangeable Rotors of 4-bladed Propeller-type, for 20~100Lit Reaction Vessels Ideal for High Viscosity & Large Volume Stirring, with Rod $\phi 12$ & $14 \times L 800 \sim 1100\text{mm}$

* Translucent PFA Shaft with Stainless-steel Insert, Totally Chemical Inert

SL.StiS416.4	Stir-Shaft, PTFE, for 20~30Lit, Rod- $\phi 12/10 \times L800$, 2-Rotors of (1) Exchangeable 1st Rotor $\phi 110 \times h40$ and (2) Bottom Rotor $\phi 110 \times h40\text{mm}$
SL.StiS418	Stir-Shaft, PTFE, for 30~80Lit, Rod- $\phi 12/10 \times L1000$, 2-Rotors of (1) Exchangeable 1st Rotor $\phi 110 \times h40$ and (2) Bottom Rotor $\phi 110 \times h40\text{mm}$
SL.StiS419	Stir-Shaft, PTFE, for 50~100Lit, Rod- $\phi 14/12 \times L1100$, 2-Rotors of (1) Exchangeable 1st Rotor $\phi 140 \times h40$ and (2) Bottom Rotor $\phi 140 \times h40\text{mm}$

Ea. Price, ₩



(3) with 2-Rotors, Propeller-type

(4) Heavy-duty Anchor-type, ideal for Middle/High Viscosity, Rod $\phi 8 \sim 12\text{mm}$ Good for Low/Middle Speed, Rotor $\phi 80 \sim 140\text{mm}$, <UK-Made>

CW.005.340.82	Shaft, $\phi 8 \times L400\text{mm}$ Rod, $\phi 80\text{mm}$ Anchor
CW.005.350.82	, $\phi 8 \times L500\text{mm}$ Rod, $\phi 80\text{mm}$ Anchor
CW.005.355.2	, $\phi 10 \times L550\text{mm}$ Rod, $\phi 100\text{mm}$ Anchor
CW.005.0365.102	, Ex-type, $\phi 10 \times L650\text{mm}$ Rod, $\phi 140\text{mm}$ Anchor
CW.005.0375.102	, $\phi 10 \times L750\text{mm}$ Rod, $\phi 140\text{mm}$ Anchor
CW.005.0375.122	, $\phi 12 \times L750\text{mm}$ Rod, $\phi 140\text{mm}$ Anchor
CW.005.03100.122	, $\phi 12 \times L1000\text{mm}$ Rod, $\phi 140\text{mm}$ Anchor

Ea. Price, ₩



(4) Heavy-duty Anchor-type



(1) Glass-Shaft



(2) PTFE-Shaft

(3) PTFE Blade,
with $\phi 6.5$ Connect-hole

Exchangeable GLASS and PTFE/Steel Cored Shaft and PTFE Blade, 교반 Blade 교체형

* Glass Shafts ; DURAN Borosilicate Glass 3.3 Rod $\phi 8$ & 10mm ; Made in Korea* PTFE Encapsulated Stainless-steel Shafts $\phi 8 \sim \phi 12$ mm, $-200^{\circ}\text{C} + 280^{\circ}\text{C}$ Stable ; Made in U.K.(1) GLASS-Shaft, $\phi 8/10$, L400~600mm

SL.StiS260	Glass-Shaft, $\phi 8$, L400mm
SL.StiS261	, L500mm
SL.StiS262	, L600mm

Ea. Price, ₩

SL.StiS265	Glass-Shaft, $\phi 10$, L400mm
SL.StiS266	, L500mm
SL.StiS267	, L600mm

Ea. Price, ₩

<L1>

(2) Cowie® PTFE/Steel Cored Shaft, $\phi 8 \sim \phi 12$ mm, L400~750mm, without Blade, <UK-Made>

CW.005.540.81	Shaft, PTFE, $\phi 8 \times$ L400mm
CW.005.550.81	, $\phi 8 \times$ L500mm
CW.005.555.1	, $\phi 10 \times$ L550mm

Ea. Price, ₩

CW.005.565.1	Shaft, PTFE, $\phi 10 \times$ L650mm
CW.005.0575.101	, $\phi 10 \times$ L750mm
CW.005.0575.121	, $\phi 12 \times$ L750mm

Ea. Price, ₩

(3) Cowie® PTFE Stir-Blade, Thick-3.2mm, Ideal for 24/- & 34/-Joint Necks, <UK-Made>

① with Square End, Fit to 24/-Joints

CW.002.052.1	Blade, Squ. $52 \times h14$ mm
CW.002.076.1	, $76 \times h19$ mm
CW.002.090.1	, $90 \times h28$ mm

Ea. Price, ₩

② with Plain End, Fit to 34/-Joints

CW.002.1065.1	Blade, Pla. $65 \times h25$ mm
CW.002.1075.1	, $75 \times h25$ mm
CW.002.1105.1	, $105 \times h25$ mm
CW.002.1125.1	, $125 \times h25$ mm
CW.002.1150.1	, $150 \times h25$ mm

Ea. Price, ₩

All PTFE Adjustably Exchangeable Rotor, Screw Propeller- & Turbine-type, Shaft-Option
For Shaft $\phi 8 \sim 16$ mm, $-200^{\circ}\text{C} + 280^{\circ}\text{C}$, PTFE 교체형/조절식 교반 로터, 교반봉은 별매■ The Exchangeable Rotors fit to $\phi 8 \sim 16$ mm Rods and Lock in place with Screw Nut(1) SciLab™ Exchangeable "Screw Propeller" Rotor only, 45° Angled 4 Blade

SL.StiS461	Rotor, PTFE, for Shaft- $\phi 8$, Rotor- $\phi 50 \times h20$ mm
SL.StiS462	, $\phi 70 \times h25$ mm
SL.StiS463	, for Shaft- $\phi 10$, Rotor- $\phi 70 \times h25$ mm
SL.StiS464	, $\phi 90 \times h35$ mm
SL.StiS465	, for Shaft- $\phi 12$, Rotor- $\phi 110 \times h40$ mm
SL.StiS466	, for Shaft- $\phi 14$, Rotor- $\phi 140 \times h40$ mm

Ea. Price, ₩

<L1>

(2) Cowie® Exchangeable "Turbine" Rotor only, $\phi 50 \sim 150$ mm, <UK-Made>

CW.005.90850	Rotor, PTFE, for Shaft- $\phi 8$ mm, Rotor- $\phi 50$ mm
CW.005.90875	, $\phi 75$ mm
CW.005.91075	, for Shaft- $\phi 10$ mm, Rotor- $\phi 75$ mm
CW.005.910100	, $\phi 100$ mm
CW.005.912100	, for Shaft- $\phi 12$ mm, Rotor- $\phi 100$ mm
CW.005.916150	, for Shaft- $\phi 16$ mm, Rotor- $\phi 150$ mm

Ea. Price, ₩

Glass and PFA (Stainless-steel Insert) PLAIN Shaft, For EXCHANGEABLE ROTOR
글라스 및 스텐심 PFA봉(1) DURAN® Glass Shaft, Round Ends, Boro-glass 3.3, $\phi 8$ & 10, L400~650mm

* Good for Eco-Pure Work

SL.StiS206	Glass Shaft, $\phi 8 \times$ L400mm
SL.StiS207	, $\phi 8 \times$ L500mm

Ea. Price, ₩

SL.StiS208	Glass Shaft, $\phi 10 \times$ L450mm
SL.StiS209	, $\phi 10 \times$ L550mm
SL.StiS210	, $\phi 10 \times$ L650mm

Ea. Price, ₩

<L1>

(2) PFA with Stainless-steel Insert Shaft, $\phi 8 \sim 14 \times$ L500~1100mm, $-200^{\circ}\text{C} + 260^{\circ}\text{C}$

* The Rod may be Cut to Any Length with a Hacksaw


* Translucent PFA, Totally Chemical Inert, $-200^{\circ}\text{C} + 260^{\circ}\text{C}$

SL.StiS451	Shaft, PFA, $\phi 8/6 \times$ L500mm
SL.StiS452	, $\phi 8/6 \times$ L650mm
SL.StiS453	, $\phi 8/6 \times$ L800mm
SL.StiS454	, $\phi 10/8 \times$ L650mm

Ea. Price, ₩

SL.StiS455	Shaft, PFA, $\phi 10/8 \times$ L800mm
SL.StiS456	, $\phi 12/10 \times$ L800mm
SL.StiS457	, $\phi 12/10 \times$ L1000mm
SL.StiS458	, $\phi 14/12 \times$ L1100mm

Ea. Price, ₩



DAIHAN™ Plate Stand, High Quality, with Stainless-steel Rod ϕ 16/23mm Ideal for the Stand of Homogenizer & Overhead Stirrers 대리석 뷰렛 스탠드

- High-Quality Plate-type Stand with Support Rod and Silicone Rubber Pad adhered
- with 1.5P Screwed Support Rod made of Stainless-steel specially designed


<L1>

Ea. Price, ₩

DH.WOS501510 (1) Standard Steel Plate Stand "ST110", Powder Coated, with Base Plate/Silicone Pad for Stirrer and Homogenizer(w20×35cm), and 1.5P Screw Rod "RD 511", ϕ 16×h800mm

Ea. Price, ₩

DH.WOS501515 (2) Wide Steel Plate Stand "ST120", Powder Coated, with Base Plate/Silicone Pad for Stirrer and Homogenizer(w30×40cm), and 1.5P Screw Rod "RD 522", ϕ 23×h1000mm



DAIHAN™ Heavy Duty Safety Stand, High Safety, with Stainless-steel Rod ϕ 23mm Ideal for the Stand of Overhead Stirrers

- Never Fall Down by Safety Mechanism.
- High-Safety H-U-type Stand with Stainless-steel Support Rod

<L1>

Ea. Price, ₩

(1) Adjustable Foot Distance Deluxe Safety Stand

DH.WOS501520 Deluxe Safety Stand "ST200", H-type Base w39~75cm adjust, Rod ϕ 23×h1200mm

(2) U-type Safety Stand

DH.WOS1177 Safety Stand "KA1177", U-type Base w34/29×27/21cm, Rod ϕ 23×h800mm

Ea. Price, ₩

(3) Universal Jumbo Clamp Holder, for Stand Rod

SL.Cla8011 Clamp Holder "CL200", Cast Aluminum, Grip Capacity, ϕ 23mm for Stand Rod

Ea. Price, ₩



Steel Retort Base Duo Stand

Steel Retort Base Duo Stand, for 2×Plain Rods of up to ϕ 16mm, without Rods With Nonskid Feet equipped, 2×Level Screw Nut, and Base Thick-8mm, <Germany-Made>

<L1>

Ea. Price, ₩


Retort Stand Base only, acceptable 2 Plain Rod(below ϕ 16mm)

US.2040 Stand Base only, "US2040", r=150mm
US.2041 "US2041", r=185mm

Optional Stainless-steel Rod of ϕ 16mm, without Screw

Ea. Price, ₩

DH.WOS502001 Stainless-steel Rod, "RD501", ϕ 16×h800mm, non-screw
DH.WOS502002 "RD502", ϕ 16×h1000mm, non-screw
DH.WOS502003 "RD503", ϕ 16×h1200mm, non-screw



Stainless-steel Rod
With or Without 1.5P Screw

DAIHAN™ ϕ 16 & 23mm High-Quality Stainless-steel (#304) Rod, for Stand Bases

<L1>

Ea. Price, ₩

(1) Plain-type ϕ 16 & 23mm Rod, without Screw

DH.WOS502001 Stainless-steel Rod, "RD501", ϕ 16×h800mm, non-screw
DH.WOS502002 "RD502", ϕ 16×h1000mm, non-screw
DH.WOS502003 "RD503", ϕ 16×h1200mm, non-screw
DH.WOS502031 "RD531", ϕ 23×h800mm, non-screw
DH.WOS502033 "RD533", ϕ 23×h1200mm, non-screw

Ea. Price, ₩

(2) 1.5P Screw-type ϕ 16 & 23mm Rod

DH.WOS502011 Stainless-steel Rod, "RD511", ϕ 16×h800mm, 1.5P Screw
DH.WOS502012 "RD512", ϕ 16×h1000mm, 1.5P Screw
DH.WOS502013 "RD513", ϕ 16×h1200mm, 1.5P Screw
DH.WOS502021 "RD521", ϕ 23×h800mm, 1.5P Screw
DH.WOS502022 "RD522", ϕ 23×h1000mm, 1.5P Screw
DH.WOS502023 "RD523", ϕ 23×h1200mm, 1.5P Screw

PTFE Micro Magnetic Stirrer Bar, for Lab & Industry, L2~12mm, PTFE 마이크로 마그네틱바



(1), (2) Micro-type

* for Universal Application

* Excellent for Chemical and Corrosion Resistance

(1) Micro-type

* -200℃+260℃

SL.Sti4012	Plain, Micro, $\phi 1.5 \times L9mm$
SL.Sti4013	, $\phi 3 \times L5mm$
SL.Sti4014	, $\phi 3 \times L10mm$
SL.Sti4015	, $\phi 3 \times L12mm$
SL.Sti4016	, $\phi 4 \times L10mm$

<L1>

Ea. Price, ₩

(2) Cowie® Premium Micro-type, <UK-Made>

* -200℃+280℃

CW.001.808	Micro, $\phi 1 \times L8mm$
CW.001.802	, $\phi 2 \times L2mm$
CW.001.805	, $\phi 2 \times L5mm$
CW.001.807	, $\phi 2 \times L7mm$
CW.001.803	, $\phi 3 \times L3mm$
CW.001.806	, $\phi 3 \times L6mm$
CW.001.808.3	, $\phi 3 \times L8mm$
CW.001.810	, $\phi 3 \times L10mm$
CW.001.2912	, $\phi 3 \times L12mm$

Ea. Price, ₩

PTFE Pivot Ring-type Magnetic Stirrer Bar, Efficient Spinning even on Curved or Uneven Bases, L30~159mm
Excellent for Chemical and Corrosion Resistance, for Lab & Industry, PTFE Pivot Ring-type 마그네틱바



(1), (2) Pivot Ring-type

(1) Standard-type

* -200℃+260℃

SL.Sti4041	Pivot Ring, $\phi 6 \times L30mm$
SL.Sti4042	, $\phi 7 \times L30mm$
SL.Sti4043	, $\phi 8 \times L45mm$
SL.Sti4044	, $\phi 9 \times L50mm$
SL.Sti4045	, $\phi 10 \times L60mm$
SL.Sti4046	, $\phi 12 \times L55mm$

<L1>

Ea. Price, ₩

(2) Cowie® Premium-type, <UK-Made>

* -200℃+280℃

CW.001.340	Pivot Ring, $\phi 8 \times L40mm$
CW.001.360	, $\phi 10 \times L60mm$
CW.001.370	, $\phi 10 \times L70mm$
CW.001.357	, $\phi 24 \times L57mm$
CW.001.3108	, $\phi 24 \times L108mm$
CW.001.3159	, $\phi 24 \times L159mm$

Ea. Price, ₩

PTFE Octagonal-type Magnetic Stirrer Bar, for Lab & Industry, L13~75mm, PTFE 팔각/옥타고날형 마그네틱바



(1), (2) Octagonal-type

* Excellent for Chemical and Corrosion Resistance

(1) Standard-type

* -200℃+260℃

SL.Sti4031	Octagonal, $\phi 8 \times L15mm$
SL.Sti4032	, $\phi 8 \times L20mm$
SL.Sti4033	, $\phi 8 \times L25mm$
SL.Sti4034	, $\phi 8 \times L30mm$
SL.Sti4035	, $\phi 8 \times L38mm$
SL.Sti4036	, $\phi 8 \times L40mm$
SL.Sti4037	, $\phi 8 \times L50mm$
SL.Sti4038	, $\phi 8 \times L75mm$

<L1>

Ea. Price, ₩

(2) Cowie® Premium-type, <UK-Made>

* -200℃+280℃

CW.001.513	Octagonal, $\phi 8 \times L13mm$
CW.001.515	, $\phi 8 \times L15mm$
CW.001.522	, $\phi 8 \times L22mm$
CW.001.525	, $\phi 8 \times L25mm$
CW.001.528	, $\phi 8 \times L28mm$
CW.001.538	, $\phi 8 \times L38mm$
CW.001.541	, $\phi 8 \times L41mm$
CW.001.551	, $\phi 8 \times L51mm$
CW.001.551.10	, $\phi 10 \times L51mm$
CW.001.564	, $\phi 10 \times L64mm$
CW.001.575	, $\phi 13 \times L75mm$

Ea. Price, ₩

PTFE Octagonal-type Stirrer Bar-Set, for Lab & Industry, L13~75mm, 16pcs/set

Excellent for Chemical and Corrosion Resistance, For Universal Application, PTFE 팔각/옥타고날형 마그네틱바 세트



Octagonal-type

■ Useful Assortment of the Most Popular Size of PTFE Magnetic Stirrer Bars

* Sectional Utility Box, with Assembly Partition

* Including Hanging Loop & Hinged Locking Lid

(1) Standard-type

■ -200℃+260℃

■ Component ; ($\phi 8 \times L15mm$) · ($\phi 8 \times L20mm$) · ($\phi 8 \times L25mm$) · ($\phi 8 \times L30mm$) · ($\phi 8 \times L38mm$) · ($\phi 8 \times L40mm$) · ($\phi 10 \times L50mm$) · ($\phi 10 \times L75mm$) × 2pcs

SL.Sti4101 Set of 16 PTFE Octagonal-type Stirrer Bars

(2) Cowie® Premium-type <UK-Made>

■ -200℃+280℃

■ Component ; ($\phi 8 \times L13mm$) · ($\phi 8 \times L15mm$) · ($\phi 8 \times L25mm$) · ($\phi 8 \times L28mm$) · ($\phi 8 \times L38mm$) · ($\phi 10 \times L51mm$) · ($\phi 10 \times L64mm$) · ($\phi 13 \times L75mm$) × 2pcs

SL.Sti4001 Set of 16 PTFE Octagonal-type Stirrer Bars

B3.802 Spare Box, Divider & Hinged Locking Lid, 10-divid. type ; 200×130×h35mm



Spare Box

<L1>

Set.	Price, ₩
16pcs	

Set.	Price, ₩
16pcs	
1ea	

PTFE Oval-/Egg-type Magnetic Stirrer Bar, for Round Bottom Flask, L10~150mm, PTFE 타원/달걀형 마그네틱바



(1) Pointed edge sharp



(2) Blunt edge

* Excellent for Chemical and Corrosion Resistance

(1) Standard-type

* -200℃+260℃

SL.Sti4058	Oval, $\phi 5 \times L10\text{mm}$
SL.Sti4059	, $\phi 7 \times L15\text{mm}$
SL.Sti4060	, $\phi 9 \times L20\text{mm}$
SL.Sti4061	, $\phi 10 \times L25\text{mm}$
SL.Sti4062	, $\phi 12 \times L25\text{mm}$
SL.Sti4063	, $\phi 15 \times L35\text{mm}$
SL.Sti4064	, $\phi 16 \times L48\text{mm}$
SL.Sti4065	for Production, $\phi 20 \times L40\text{mm}$
SL.Sti4066	, $\phi 20 \times L50\text{mm}$
SL.Sti4067	, $\phi 20 \times L63\text{mm}$
SL.Sti4068	, $\phi 25 \times L70\text{mm}$

<L1>

Ea. Price, ₩

(2) Cowie® Premium-type, <UK-Made>

* -200℃+280℃

CW.001.610	Oval, $\phi 5 \times L10\text{mm}$
CW.001.615	, $\phi 6 \times L15\text{mm}$
CW.001.620	, $\phi 10 \times L20\text{mm}$
CW.001.625.10	, $\phi 10 \times L25\text{mm}$
CW.001.625	, $\phi 12 \times L25\text{mm}$
CW.001.630	, $\phi 16 \times L30\text{mm}$
CW.001.635	, $\phi 16 \times L35\text{mm}$
CW.001.650.17	, $\phi 17 \times L50\text{mm}$ New
CW.001.640	, $\phi 20 \times L40\text{mm}$
CW.001.650	, $\phi 20 \times L50\text{mm}$
CW.001.664	for Production, $\phi 20 \times L64\text{mm}$
CW.001.670	, $\phi 20 \times L70\text{mm}$
CW.001.670.27	, $\phi 27 \times L70\text{mm}$
CW.001.6100	, $\phi 30 \times L100\text{mm}$
CW.001.6150	, $\phi 40 \times L150\text{mm}$

Ea. Price, ₩

PTFE Oval- & Octagonal-type Mixed Stirrer Bar-Set, for Lab & Industry, L10~75mm, 18pcs/set

Excellent for Chemical and Corrosion Resistance, -200℃+260℃, PTFE 타원/달걀형+육타고날형 마그네틱바 세트



■ Supplied by SciLab Korea

■ Useful Assortment of the Most Popular Size of PTFE Magnetic Stirrer Bars

■ Component ; Oval - ($\phi 5 \times L10\text{mm}$) · ($\phi 7 \times L15\text{mm}$) · ($\phi 15 \times L35\text{mm}$) · ($\phi 16 \times L48\text{mm}$) × 2pcs,
($\phi 25 \times L70\text{mm}$) × 1pcs

Octagonal - ($\phi 8 \times L15\text{mm}$) · ($\phi 8 \times L20\text{mm}$) · ($\phi 8 \times L30\text{mm}$) · ($\phi 8 \times L50\text{mm}$) × 2pcs,
($\phi 8 \times L75\text{mm}$) × 1pcs

* For Universal Application * Ideal for Storage of Small Parts and Accessory.

* Sectional Utility Box, with Assembly Partition

* Including Hanging Loop & Hinged Locking Lid



Spare Box

<L1>

SL.Sti4107 Set of 18 PTFE Oval- & Octagonal-type Mixed Stirrer Bars

B3.802 Spare Box, Divider & Hinged Locking Lid, 10-divid. type ; 200×130×h35mm

Set. Price, ₩

18pcs

1ea



PTFE Cross-type Magnetic Stirrer Bar, for Lab & Industry, L5~50mm, PTFE Cross-type 마그네틱바

* Excellent for Chemical and Corrosion Resistance

* Very Stable General Purpose Stirring

(1) Standard-type

* -200℃+260℃

SL.Sti4079	Cross, $\phi 8 \times L20\text{mm}$
SL.Sti4080	, $\phi 9 \times L25\text{mm}$
SL.Sti4081	, $\phi 10 \times L30\text{mm}$
SL.Sti4082	, $\phi 10 \times L40\text{mm}$
SL.Sti4083	, $\phi 15 \times L50\text{mm}$

<L1>

Ea. Price, ₩

(2) Cowie® Premium-type, <UK-Made>

* -200℃+280℃

CW.001.2401	Cross, $\phi 10\text{mm} \times h5\text{mm}$
CW.001.2402	, $\phi 20\text{mm} \times h8\text{mm}$
CW.001.2403	, $\phi 25\text{mm} \times h9\text{mm}$
CW.001.2404	, $\phi 30\text{mm} \times h10\text{mm}$
CW.001.2405	, $\phi 38\text{mm} \times h11\text{mm}$
CW.001.2406	, $\phi 50\text{mm} \times h15\text{mm}$
CW.001.2407	, $\phi 60\text{mm} \times h20\text{mm}$

Ea. Price, ₩



PTFE Elliptical-type Magnetic Stirrer Bar, for Round Bottom Flasks, -200℃+260℃, L20~41mm

Excellent for Chemical and Corrosion Resistance, for Lab & Industry, PTFE Elliptical-type 마그네틱바

■ Supplied by SciLab Korea

* Ideal for Round Bottom Flasks

* Rounded pivot point to reduce friction

<L1>

SL.Sti4071	Stirrer Bar, Elliptical, $\phi 10 \times L20\text{mm}$
SL.Sti4072	, $\phi 12 \times L26\text{mm}$
SL.Sti4073	, $\phi 14 \times L34\text{mm}$
SL.Sti4074	, $\phi 14 \times L41\text{mm}$
SL.Sti4075	, $\phi 16 \times L40\text{mm}$

Ea. Price, ₩

Articles		Cat. No	Description
 			Cowie® PTFE Wing-type Stirrer Bars, for Test-Tubes & V-Vials Excellent for Chemical and Corrosion resistance, $-200^{\circ}\text{C} + 280^{\circ}\text{C}$, <UK-Made> PTFE 날개형 마그네틱바
(1)	(2)		
 			For V-Vials & Test-Tubes, for Tube $\phi 6\sim 12\text{mm}$
(3)	(4)		
		CW.001.2201.1	(1) Stirrer bar, $\phi 5.5\times h9\text{mm}$, Tube id $\phi 6\text{-}7\text{mm}$
		CW.001.2201	(2) , $\phi 8\times h11\text{mm}$, Tube id $\phi 9\text{-}10\text{mm}$
		CW.001.2202	(3) , $\phi 10\times h16\text{mm}$, Tube id $\phi 11\text{-}12\text{mm}$
			For Test-Tubes & Centrifuge Tubes, for Tube $\phi 14\sim 15\text{mm}$
(5)			
		CW.001.2203	(4) Stirrer bar, Rounded, Tube id $\phi 14\text{-}15\text{mm}$
		CW.001.2204	(5) , Tapered , Tube id $\phi 14\text{-}15\text{mm}$




SciLab™ PTFE Magnetic Retriever, $\phi 9$ & 10 , L250~450mm $-200^{\circ}\text{C} + 260^{\circ}\text{C}$ Stable, PTFE 마그네틱 리트리버/홀더

- Supplied by SciLab Korea
- 실험 용기로부터 마그네틱 바를 꺼내는 집게 역할의 홀더임.
- Can be used to Remove the Stirrer Bar from the Vessel

SL.Sti4087	Retriever, PTFE, $\phi 9 \times L 250\text{mm}$
SL.Sti4088	, $\phi 9 \times L 300\text{mm}$
SL.Sti4089	, $\phi 9 \times L 340\text{mm}$
SL.Sti4090	, $\phi 10 \times L 450\text{mm}$



<L1>	
Ea.	Price, ₩



(1), (2) PTFE & PP Retriever

Cowie® Premium Magnetic Retriever, PP & PTFE, L150~450mm, <UK-Made> 마그네틱 리트리버

(1) PP Magnetic Retrievers with Hanging Ring, $-10^{\circ}\text{C} + 120/140^{\circ}\text{C}$

- * Made of polypropylene(PP)
- * Hanging Ring

CW.004.1300	Retriever, PP, $\phi 10 \times L 300\text{mm}$
CW.004.1350	, $\phi 10 \times L 350\text{mm}$
CW.004.1450	, $\phi 10 \times L 450\text{mm}$



<L1>	
Ea.	Price, ₩

(2) PTFE Magnetic Retrievers, $-200^{\circ}\text{C} + 280^{\circ}\text{C}$

- * Straight shaft

CW.004.150	Retriever, PTFE, $\phi 10 \times L 150\text{mm}$
CW.004.250	, $\phi 10 \times L 250\text{mm}$
CW.004.350	, $\phi 10 \times L 350\text{mm}$
CW.004.450	, $\phi 10 \times L 450\text{mm}$

Ea.	Price, ₩



Timer

Multi-Timer/Alarm Clock

DAIHAN™ Pocket size Digital Timer, LCD Display, Alarm, Count-Up/Down

Compact, 5 × 4 × h2cm, 포켓형 타이머, 알람, 카운트 업/다운

- * Maximum Time Count ; 99hr : 59min. : 59sec.
- * Count Up and Count Down Function.
- * CE Approval.
- * Compact Design
- * Strong Magnet at the Back for Easy Attachment.
- * Dimension ; 53 × 42 × 21 mm (LCD Size ; 33 × 12 mm)
- * Battery ; 1 × 1.5Volts "A76" Size (LR44, G13) or Equivalent

<L1>

A1.TM190 Timer, Count-Down/-Up, Timer up to 99:59:59

Ea. Price, ₩



Timer

SciLab™ Digital Timer, Count-Up/Down, LCD Display, 8 × 7 × h2cm

With 10 Setting Buttons, Alert Function, 디지털 타이머

- * Maximum Time Count ; 99min. : 59sec.
- * When the Count-down Time Finished, there will be Alert Sound for 60 Seconds and LED Indicator Light Flicking
- * Large LCD Display, 53 × 24 mm
- * 10 Setting Buttons for Setting the Desire Time
- * Flip Open Easel Stand, Spring Fastener Slip
- * Easy to Attach with Strong Magnet at the Back

<L1>

SL.Tim3001 Digital Timer, Count-up & Count-down, w/10 Setting Buttons

Ea. Price, ₩



Timer

size ;
69 × 74 × 23 mm

3-row Display Timer/Alarm

DAIHAN™ Handy 3-row Display Multi Timer, wide LCD, Alarm/Snooze set

Count-Up/Down, Compact, 3 Timers, 12/24hr Selectable, 1sec Divi, 3열 표시 멀티 타이머

- * 3 (T1/T2/T3) Timers in One Unit.
- * Maximum Count ; 23:59:59 (Hr. : Min. : Sec.)
- * Count Up and Count Down Function.
- * Pocket Size with Large LCD Display.
- * Flip Out Stand and Strong Magnet at the Back for Easy Attachment.
- * Triple Row Display for Instant Reading.
- * CE Approval.
- * Compact Design.
- * Different Alarm Sounds for Timers, Easy to Identify which Timer is Beeping.
- * Dimension ; 69 × 74 × 23 mm (LCD Size ; 37 × 34 mm)
- * Ambient Temperature ; 0°C~50°C or 32°F~122°F
- * Weight (g) / Color / Resolution ; 65 / Black / 1 Second
- * Battery / Package ; 1 × 1.5Volts "A76" Size (LR44, G13) or Equivalent / Sliding Blister with Print Card

<L1>

A1.TM150 3-row Display Digital Multi-Timer, Handy/Pocket, 3(T1/T2/T3) Timers in One Unit

Ea. Price, ₩



Timer

4-channel Multi-Timer/Alarm Clock

DAIHAN™ 4-Channel Multi-Timer, with Alarm/Snooze set, Count-Up/Down, Large LCD

With LCD, 4 Timers in One Unit, 99:59:59 hr, 12/24hr Selectable, 1sec Division, 디지털 멀티 타이머

- * 4 (T1/T2/T3/T4) Timers in One Unit.
- * Maximum Count ; 99:59:59 (Hr. : Min. : Sec.)
- * Alarm Clock Function
- * 4 Identified Alarm Sounds
- * Strong Magnet at the Back for Easy Attachment.
- * Dimension / Resolution ; 73 × 66 × 25 mm (LCD Size ; 50 × 23 mm) / 1 Second
- * CE Approval.
- * Battery ; 1 × 1.5Volts "A76" Size (LR44, G13) or Equivalent
- * Ambient Temperature ; 0°C~50°C or 32°F~122°F

<L1>

A1.TM044 4-channel Multi-Timer, Count Up & Down, 4 Timers in One Unit

Ea. Price, ₩



Stopwatch
-Timer

ETL® Digital Timer/Stopwatch, 100min 39sec, 1sec~99min 99sec

With Count-Up/Down, 타이머/스탑워치

- * 99min 99sec
- * Alarm Clock
- * 1sec Stopwatch (Timer and Stopwatch)
- * 1.5V Battery : about 2 years life-time
- * Size : w62 × d58 × h10 mm
- * Mounting way : Attach to the Table, Cloth, wall Using a Magnet Clip-Stand

<B3>

SL.Tim3003 Timer/Stopwatch, Count-Down/-Up, 99min. 99sec.

Ea. Price, ₩

DAIHAN™ Analog Ultrasonic Cleaner “WUC-A”, Timer/Temp Output Controller, with Certi. & Traceability, 1.2~22Lit
With Stainless-steel Flat Lid, Highly Effective Cleaning, Up to 85℃, 0~30min, 40kHz Frequency, without Basket
초음파 세척기, 온도 및 시간 설정, 고효율, 다용도, 리드 포함




3-Year Warranty

10Lit Analog Ultrasonic Cleaner with Flat Lid
“WUC-A10”



3-Year Warranty

1.2 & 6Lit Analog Ultrasonic Cleaner with Flat Lid
“WUC-A01 & WUC-A06”



Stainless-steel Wire Basket (Optional)



(6) Fine Mesh Strainer

- Supplied by DAIHAN Korea
- CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Made of Stainless-steel, Finished Surface
- Analog Controller, Timer, Temperature and Output Control
- PL(Product Liability) Insurance

- * Useful for Cleaning to Remove Dirt, Grease, Waxes and Oil from all Lab Instruments by Ultrasonic Frequency
- * Made of Stainless-steel, Finished Surface
- * Analog Controller, Timer, Temperature and Output Control
- * Temperature Range from Ambient to 85℃. (1.2 & 1.8 Lit ; 65℃ Fixed)
- * Convenient Preset Function ; Heating and Timer
- * HF-Frequency ; 40kHz
- * Bottom-Mounted Drain (WUC-A10 and WUC-A22)
- * Maximum 30 minutes of Cleaning Time (Except WUC-A10)

Specifications & Ordering information

<L1>

Spec. Model	“WUC-A01” Max. Capa. 1.2Lit	“WUC-A02” Max. Capa. 1.8Lit	“WUC-A03” Max. Capa. 3.3Lit	“WUC-A06” Max. Capa. 6Lit	“WUC-A10” Max. Capa. 10Lit	“WUC-A22” Max. Capa. 22Lit
Capacity	1.2Lit	1.8Lit	3.3Lit	6Lit	10Lit	22Lit
Frequency	40KHz					
HF-Peak Out	72W	75W	124W	168W	265W	379W
Heating Temp	65℃ Fixed			ambient to 85℃		
Heating Power	55W		110W		240W	320W
Timer	None	0~30min.	0~40min.			
Drain-Valve	None				1/2" NPT	
Tank Size inner (w×d×h) (mm)	150×135×60	150×135×100	240×140×100	290×150×150	290×240×150	500×300×150
Weight	1.7kg	2.7kg	4.4kg	5.7kg	7.8kg	12kg
Packing Dimension (w×d×h) (mm)	200×200×190	200×200×240	320×200×290	370×220×350	430×320×360	640×380×370
Gross Weight	2.1kg	3.1kg	5kg	6.7kg	9.4kg	14.3kg
Power Consumption	117W	117W	230W	280W	528W	739W
Electrical Reqts.	1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz					
Cat. No.	230V	DH.WUC.A01H	DH.WUC.A02H	DH.WUC.A03H	DH.WUC.A06H	DH.WUC.A10H
	120V	DH.WUC.A101H	DH.WUC.A102H	DH.WUC.A103H	DH.WUC.A106H	DH.WUC.A110H
Price, Each. ₩						

* Other Specifications are available upon Customer's Request.

Optional Stainless-steel Wire Basket

- DH.WUC.B02 (1) Basket, Wire, “B02”, for 1.8 Lit, 144×131×h92mm
- DH.WUC.B03 (2) , “B03”, for 3.3 Lit, 224×128×h83mm
- DH.WUC.B06 (3) , “B06”, for 6 Lit, 267×136×h126mm
- DH.WUC.B10 (4) , “B10”, for 10 Lit, 276×218×h132mm
- DH.WUC.B22 (5) , “B22”, for 22 Lit, 467×275×h131mm
- DH.WUC.B100 (6) Fine Mesh Strainer, “B100”, φ 100mm

Ea. Price, ₩

DAIHAN™ Digital Ultrasonic Cleaner Set "WUC-D", Microprocessor Control, 3.3~22Lit
With Stainless-steel Wire Basket & Flat Lid, Highly Effective Cleaning, Up to 80°C, 0~60min, 40kHz Frequency
디지털 초음파 세척기 세트, 와이어 바스켓 및 리드 포함, 마이크로프로세서 제어




10Lit Digital Ultrasonic Cleaner Set, "WUC-D10H"
with Flat Lid and Wire Basket

3-Year Warranty



3.3Lit Digital Ultrasonic Cleaner Set, "WUC-D03H"
with Flat Lid and Wire Basket

3-Year Warranty



Stainless-steel Wire Basket
(included)



(5) Fine Mesh Strainer

- Supplied by DAIHAN Korea
- CE Certified
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Useful for Cleaning to Remove Dirt, Grease, Waxes and Oil from all Lab Instruments by Ultrasonic Frequency
- Microprocessor Control ; Timer, Temperature and Output Control
- PL (Product Liability) Insurance
- * Useful for Cleaning to Remove Dirt, Grease, Waxes and Oil from all Lab instruments by Ultrasonic Frequency
- * Microprocessor Control ; Timer, Temperature and Output Control
- * Multi-Frequency Sonication Body
- * HF-Frequency ; 40kHz
- * Temperature Range from Ambient to 80°C
- * Convenient Preset Function ; Heating and Timer
- * Drain is mounted in the bottom (WUC-D10H, WUC-D22H)
- * Maximum 60 minutes of Cleaning Time
- * Very Strong Cleansing Degassing & Sweep Function

Specifications & Ordering information

<L1>

Spec.	Model	“WUC-D03H” Max. Capa. 3.3Lit	“WUC-D06H” Max. Capa. 6Lit	“WUC-D10H” Max. Capa. 10Lit	“WUC-D22H” Max. Capa. 22Lit
Capacity		3.3Lit	6Lit	10Lit	22Lit
Frequency		40kHz			
HF-peak out		100W	253W	407W	528W
Heat Temp		ambient to 80℃			
Heating Power		172W	172W	258W	344W
Timer		0~30min.	0~60min.		
Display		FND			
Drain-valve		None		1/2" NPT	
Tank Size inner (w×d×h) (mm)		240×140×100	290×150×150	290×240×150	500×300×150
Weight net		3.4kg	5.8kg	8kg	12.9kg
Packing Dimension (w×d×h) (mm)		320×200×290	360×220×350	440×320×360	640×380×370
Gross Weight		3.7kg	6.3kg	8.5kg	13.4kg
Power Consumption		290W	425W	665W	872W
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz			
Cat. No.	230V	DH.WUC.D03H	DH.WUC.D06H	DH.WUC.D10H	DH.WUC.D22H
	120V	DH.WUC.D103H	DH.WUC.D106H	DH.WUC.D110H	DH.WUC.D122H
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

Spare Wire Basket, Stainless-steel

- DH.WUC.B03 (1) Basket, Wire, "B03", for 3.3 Lit, 224×128×h83mm
 DH.WUC.B06 (2) , "B06", for 6 Lit, 267×136×h126mm
 DH.WUC.B10 (3) , "B10", for 10 Lit, 276×218×h132mm
 DH.WUC.B22 (4) , "B22", for 22 Lit, 467×275×h131mm
 DH.WUC.B100 (5) Fine Mesh Strainer, "B100", ∅ 100mm

Ea. Price, ₩

DAIHAN™ Large Capacity Ultrasonic Cleaner Set “WUC-N”, Remote Control System, 30~74Lit

With Digital Remote Generator · Stainless-steel Flat Lid · Drain Valve, Highly Effective Cleaning, without Basket, Up to 105°C, 40kHz
대용량 초음파 세척기 세트, 리모트 제너레이터 · 리드 · 드레인 밸브 포함



3-Year Warranty



47Lit Large Capacity Ultrasonic Cleaner Set
Without Basket “WUC-N47H”



Stainless-steel Wire Basket
(Optional)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL(Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Useful for Cleaning to Remove Dirt, Grease, Waxes and Oil from all Lab Instruments by Ultrasonic Frequency
- Remote Generator ; Power Control and Timer Function with LCD
- * BTL Sonicator Installed ; Robust and Stable Running
- * Sweep Function can Realize Complete Cleaning by the Dispersion of Ultrasound
- * Auto Optimization of Frequency as various Water Level and Temperature.
- * Supply Full Set with Lid, Foot Switch and Generator (Basket excluded)
- * Drain is mounted at the bottom

Specifications & Ordering information

<L1>

Spec.	Model	“WUC-N30H” Max. Capa. 30Lit	“WUC-N47H” Max. Capa. 47Lit	“WUC-N60H” Max. Capa. 60Lit	“WUC-N74H” Max. Capa. 74Lit
Capacity		30Lit	47Lit	60Lit	74Lit
HF-frequency		40kHz			
Heat Temp		ambient to 105°C			
Heating Power		990W			
Timer (min)		0~99min.			
Tank Size w×d×h(mm)	Inner	400×300×250	450×350×300	500×400×300	550×450×300
	External	480×380×400	530×430×450	530×430×450	580×480×450
Net Weight		26kg	31kg	36kg	41kg
Packing Size w×d×h(mm)		620×520×640	670×570×710	720×620×710	770×670×710
Gross Weight		32kg	37kg	42kg	47kg
Power Consumption		1,694W	1,980W	2,200W	2,420W
Remote Generator	Display	LCD			
	Power	600W	900W	1,200W	1,500W
	Weight	9kg	10kg		
Cat. No.	230V	DH.WUC.N30H	DH.WUC.N47H	DH.WUC.N60H	DH.WUC.N74H
	120V	DH.WUC.N130H	DH.WUC.N147H	DH.WUC.N160H	DH.WUC.N174H
Price, Each. ₩					

* Other Specifications are available upon Customer's Request.

Optional Wire Basket, Stainless-steel

Ea. Price, ₩

SL.Bas7021	Basket, Wire, “Bas7021”, for 30Lit, 340×240×h150mm
SL.Bas7022	“Bas7022”, for 47Lit, 390×290×h200mm
SL.Bas7023	“Bas7023”, for 60Lit, 440×340×h200mm
SL.Bas7024	“Bas7024”, for 74Lit, 490×390×h200mm



1Lit

5Lit

deconex® Ultrasonic Cleaning Special Detergent, 1 & 5Lit

Ideal for Max. Effect of the Ultrasonic Cleaning, Optimum Cleaning Results, 초음파 세척기 전용 세척제

■ Made in Switzerland

* Specially Developed for Ultrasonic Applications, Running 5 to 15 minutes for Ultrasonic Cleaning

<L1>

Model	(1) “11 Universal” Cleaner	(2) “15 PF-x” Cleaner	(3) “16 PLUS” Cleaner
Spec.			
pH	11.8 at 1%, Mildly Alkaline	11.9 at 1%, Alkaline	Neutral
Ingredients	<With Activated Chlorine> Alkalis, Surfactants, Inorganic Chlorine Carrier	<Phosphate-free> Alkalis, Complexing Agents, Surfactants, Solubilizer	<Not Contain Alkalis> Surfactants, Complexing Agents, Sequestering Agents
Used for	General Contamination, Oily and Greasy Residues, Soot, Ink, and Drilling · Grinding · Polishing · Lapping Residues	General Contamination, Oily and Greasy Residues, Dried Blood, Proteins	Corrosion-Sensitive Laboratory Material and Instruments
Not Suitable for	Aluminium, zinc, tin	Aluminium, zinc, tin	-
Dosage	1~5%	1~5%	3~5%
Cat. No.	1Lit	5Lit	5Lit
	DH.Det001	DH.Det002	DH.Det005
Pk. Price, ₩			

DAIHAN™ Ultrasonic Homogenizer Full Sets, "UH-300" • "UH-500" • "UH-1000", Up to 1,000ml, 19.5-25kHz
With Standard Probe, Automatic Data Storage, Timer(0.1sec~999sec), 디지털 초음파 호모지나이저, 사운드 프루프 박스 & 프로브 포함, 저소음




Fig. Example of Use
※ Beaker is not Included



Spare & Optional Titanium Alloy Probe

- Supplied by DAIHAN Korea
- 본체, 사운드 프루프 박스, 프로브 세트가 구성되어 경제적이고 편리성이 높음
- 사운드 프루프 박스 포함 - 샘플 처리 시 발생하는 소음을 최소화하여 작업 효율성을 높임
- 4.3인치 TFT 정전식 터치스크린 제어를 사용하여 쉽게 조작 가능
- Component of Set ;
 - (1) Ultrasonic Generator,
 - (2) Sonic System with Standard Horn ($\phi 6\text{mm}$ for UH-300 & UH-500, $\phi 20\text{mm}$ for UH-1000),
 - (3) Sound Proof Box, (4) Lab Support Jack, (5) Temp. Probe, (6) Power Cable, (7) Fuse, (8) User Manual
- Multi-Functional & Multi-Purpose Instrument that Uses Strong Ultrasound to Produce Cavitation Effect in Liquid and Ultrasonically Treat Substances
- Applications ; Cell Pulverization, Emulsification, Dispersion, Homogenization, Extraction, Defoaming, Degassing, Acceleration of Chemical Reactions
- * Large LCD Screen & Touch Operation Panel
- * Power Adjustable
- * Realize the Data Display of Man-Machine Friendly Operation Interface at a Glance
- * Intelligent Control ; Central Computer is used to Control the Total Working Time Centrally
- * Strong Sound Insulation - Made of Special Materials, which Greatly Reduces the Noise and Effectively Insulates the Noise
- * Durable Probe made by Titanium Alloy Material
- * Streamlined Appearance Design
- * Multiple Models Optional Matching Different Horn with Different Hosts
- * Data Storage Function - Can store 20sets of Data, Experimental Data and Convenient Recording
- * Liftable Operating Table Set in the Sound Proof Box - Easy to Adjust the Height of the Samples and Take Out
- * Can Realize the Voice Alarm Function caused by Malfunction, Over-Temperature or Operation Error

Specifications & Ordering information

<L1>

Spec.	Model	"UH-300"	"UH-500"	"UH-1000"
Mixing Capacity		0.1-300ml	0.1-500ml	1-1000ml
Power (Continuously Adjustable)		25-250W	10-650W	25-1200W
Frequency Range		20-25kHz		19.5-20.5kHz
Probe Included		$\phi 6\text{mm}$		$\phi 20\text{mm}$
Optional Probe		$\phi 2\text{mm}$	$\phi 2$ & $\phi 12\text{mm}$	$\phi 6$ & $\phi 15$ & $\phi 22\text{mm}$
Single Ultrasound Time		0.1-999s		
Single Gap Time		0.1-999s		
Total Time (Ultrasound + Gap)		1-99h 59m 59s		
Temperature Control Range		0-100 °C		
Alarm Function		Temp. Time, Overload, No Load, Over Temperature		
Voice Alarm & Prompt Function		Yes		
Display Mode		Touch Screen		
Data Storage		20 Groups		
Power Supply		220V / 60Hz		
Dimension & Weight		410 × 225 × h290mm, 12kg		410 × 225 × h290mm, 15kg
Packing Size & Gross Weight		470 × 305 × h430mm, 16kg		470 × 305 × h430mm, 19kg
Sound Box Packing Size & Gross Weight		350 × 350 × h550mm, 10kg		
Cat. No.		DH.Hom1001	DH.Hom1002	DH.Hom1003
Price, Set. ₩				

Spare & Optional Titanium Alloy Probe

Cat. No.	Probe Size (mm)	Power Range Suitable	Capacity Range	Power Rate	Ea. Price, ₩
DH.Pro7001	$\phi 2$	Power Less than 1000W	0.2-5ml	1-35%	
DH.Pro7002	$\phi 6$		10-100ml	1-60%	
DH.Pro7003	$\phi 12$		50-500ml	20-95%	
DH.Pro7004	$\phi 15$	Power Higher than 1200W	100-600ml	1-80%	
DH.Pro7005	$\phi 20$		100-1000ml	10-90%	
DH.Pro7006	$\phi 22$		200-1000ml	20-95%	



(1) FA0217, Single-Branch, 300ml STS Funnel



(2) FA0211, 3-Branch, 300ml STS Funnel



(3) FA0212, 6-Branch, 300ml STS Funnel



(4) FA0214, 3-Branch, 300ml Glass Funnel



(5) FA0215, 6-Branch, 300ml Glass Funnel

TJV[®] Vacuum Manifolds & Filter Holder Funnel, Ideal for Strong Solvents & HPLC Media With All Stainless-steel & Glass Funnel, Single · 3 · 6-Branch, 진공 매니 폴드 & 필터 홀더

- Ideal for Chemistry Analyze, Biochemical, Pharmaceutical, Sanitation Test, Environment Test, Water Quality Analyze, Beverage and Science Research
- Particular Design and Produce for Several Samples Filtrate at the Same Time
- Single or Manifolds Vacuum Holder Operate Together
- Stainless-steel Nature of the Manifolds Vacuum Holder can Filtrate at 180°C
- High-quality Extra Hard Glass is Elegant in Workmanship, and be to Bear up 200°C
- Type 316L Stainless-steel Resistant to Acids Alkalies, and be Able to bear Corrosion, Easy to High-Temperature Sterilization
- Complete Set with ; (1) Glass / Stainless-steel(316L) Funnel of 300ml.
- (2) PTFE Support Screen (Pore Size ; 1.5mm). (3) PTFE Valve.
- (4) Glass or Stainless-steel Holder(316L). (5) Clamp of Aluminum Alloy.

* Each Holder have Individual Control Valve, Convenience in Antisepsis, High Efficiency
 * Test Three or Six Samples at the Same Time, High Sensitivity, Low False-positive, and Easy to Operation
 * Suitable Membrane; The Choosing of the Membrane Filter is the Key to the Best Filtration. The Nature of Membrane Filter must Suit for Different Chemical Performance of the Filter Medium. At Present Mostly are PVDF(organic), Synthetic Fabric Colophony Membrane(Aqueous), and Nylon Membrane(both Aqueous and Organic).

Specification

Model	(1) FA0217	(2) FA0211	(3) FA0212	(4) FA0214	(5) FA0215
Specification	1-Branch	3-Branch	6-Branch	3-Branch	6-Branch
Funnel	Stainless-steel 316L of 300ml			Glass of 300ml	
Valve	PTFE			PTFE	
Support Screen	Stainless-steel 316L (Pore Size 100 μ m)			PTFE (Pore Size 1.5mm)	
Connection	Silicone Stopper			Glass Joint	
Holder	Stainless-steel 316L			Stainless-steel 316L	

<L1>

(1)~(3) Stainless-steel Manifold, Stainless-steel Base

- TJ.FA0217 (1) STS Manifolds, 1-Branch, 300ml STS Funnel, PTFE Valve, Screen 100 μ m, STS Connection with Silicone Stopper
 TJ.FA0211 (2) STS Manifolds, 3-Branch, 300ml STS Funnel, PTFE Valve, Screen 100 μ m, STS Connection with Silicone Stopper
 TJ.FA0212 (3) STS Manifolds, 6-Branch, 300ml STS Funnel, PTFE Valve, Screen 100 μ m, STS Connection with Silicone Stopper

Set.	Price, ₩

(4) & (5) Glass Manifold, Stainless-steel Base

- TJ.FA0214 (4) Glass Manifolds, 3-Branch, 300ml Glass Funnel, PTFE Support Screen 1.5mm, PTFE Valve, Glass 50/40 Joint Connection
 TJ.FA0215 (5) Glass Manifolds, 6-Branch, 300ml Glass Funnel, PTFE Support Screen 1.5mm, PTFE Valve, Glass 50/40 Joint Connection

Set.	Price, ₩

Components / Spare Parts

- TJ.FA0207 Stainless-steel 316L Funnel, 300ml, without Lid
 TJ.FA0208 Glass Funnel 300ml
 TJ.FA0221 Anodized Aluminium Clamp
 TJ.FA0544 Lid, for 300ml Stainless-steel 316L Funnel

Ea.	Price, ₩

SciLab[™] Premium Vacuum Manifold, DURAN Glass Stopcock, 3~5 Places Single- or Double-Line, 유리콕형 진공매니폴드

- Supplied by SciLab Korea
- Made of DURAN[®] Glass Stopcock & Borosilicate Glass 3.3
- * Stopcock Pass Bore ; up to 2.5 & 4 mm
- * Hose Connection Tube ; od ϕ 10mm Olive-Type

<L1>

(1) Single-Line, with Glass Stopcock (SM.28640480), 3~5 Places

- SL.Man2020 Vacu. Manifold, Single, GI-Stopcock, 1-way, 3 $\times$ Places
 SL.Man2021 , 4 $\times$ Places
 SL.Man2022 , 5 $\times$ Places

Ea.	Price, ₩


(2) Double-Line, with 2 Bore Stopcock (SM.28642500), 3~5 Places

- SL.Man2024 Vacu. Manifold, Double, GI-Stopcock, 2-way, 3 $\times$ Places
 SL.Man2025 , 4 $\times$ Places
 SL.Man2026 , 5 $\times$ Places

Ea.	Price, ₩



(1) Single-Line ; 3~5 Places, PTFE Valved



(2) Double-Line ; 3~5 Places, PTFE Valved

SciLab™ Premium 10⁻⁹ Torr High Vacuum Manifold, PTFE Valve, 3~5 Places
Single-or Double-Lines, with PTFE Rotaflo™ Valve, PTFE 밸브형 고진공 매니폴드

- Supplied by SciLab Korea ■ Made of Pyrex® Rotaflo® PTFE Valve & Boro-glass 3.3
- Custom-Made Available / 규격외 Items 자유롭게 주문하세요!
- "Pyg.HP3ST.7" PTFE Valve Bore 3mm * Hose Connection Tube ; od ϕ10mm Olive-Type

<L1>

(1) Single-Line, With PTFE Needle Valve, 3~5 Places

SL.Man2011	Vacu. Manifold. Single, 1-Way PTFE Valve, 3×Places
SL.Man2012	, 4×Places
SL.Man2013	, 5×Places

Ea. Price, ₩

(2) Double-Line, With PTFE Needle Valve, 3~5 Places

SL.Man2015	Vacu. Manifold. Double, 2-Way PTFE Valve, 3×Places
SL.Man2016	, 4×Places
SL.Man2017	, 5×Places

Ea. Price, ₩

Solid-phase Extraction Vacuum Manifold, for Chromatography, Used with SPE Cartridges, 12~24 port, 270×90×h160mm
Ideal for Sample Preparation in Laboratories, High Recovery Rate, Effective Removal, Concentrate Analytes

고체상추출(SPE) 매니폴드, SPE컬럼 & 진공펌프 별도



Solid-phase Extraction Vacuum Manifold, 19 port

- Supplied by SciLab Korea
- 액체 시료에서 SPE 카트리지의 흡착층을 통과시켜 표적 화합물 추출 하는데 사용되는 실험실 장비
- HPLC, GC의 의한 미량 유기물 분석 전 분석시료의 전처리(정제, 농축, 분획/분리, 용매변경)목적
- 액-액 추출법(Liquid-liquid Extraction, LLE) 단점 보완, 높은 회수율, 실험시간 단축
- 일반적으로 환경분석, 제약 및 식품 안전과 같은 다양한 응용 분야에서 사용
- Solid Phase Extraction is a Chromatographic Technique Used to Prepare Samples for Subsequent Analysis
- The General Procedure includes Sorbent Activation, Loading, Washing and Eluting
- * Increase Rear Valve Height by 10mm to Minimize Interference with Front Valve
- * Allows to Observe the Flow of all Possible Ports * Individual Flow Control Stopcocks
- * Sufficient Clearance between Ports
- * Includes Vacuum Gauge and Vacuum Control Valve * Made of Special Solid Glass
- * Can Withstand more than -0.085Mpa of the High Negative Pressure
- * 12 ~ 24 Columns Maximum * Processing ; Vacuum Manifold

※ SPE Cartridges ; Visit DAIHAN-Sci.com Shopping Mall

<L1>

SL.SpeM102	SPE Manifold, 12 port
SL.SpeM103	, 19 port
SL.SpeM104	, 24 port
DH.Pu0300	Vacuum Pump

Set. Price, ₩

DAIHAN™ 105mbar Oil Free Vacuum Pump "PU-300", with Fine Vacuum Regulator and Moisture Trap, 20Lit/min
With Prote™ Overflow Protection, Built-in Thermal Protection, Protection Class IP30, 진공 펌프


105mbar Oil Free Vacuum Pump
"PU-300"

3-Year Warranty


- 진공 처리가 필요한 다양한 실험에서 사용이 가능 ■ CE, CSA, and RoHS Certified
- Used for General Filtration, Microbiological Test, Suspended Solid Test, Solid-Phase Extraction, Microplate Vacuum Manifold, Air Sampling, etc.

- * Driven by Piston without Lubricant, Regular Oil Change and Maintenance ; No Oil Pollution
- * Prote™ Overflow Protection is Installed to Prevent a Large Amount of Liquid from being Sucked into the Pump.
- * Equipped with Filter Cartridge in Air Inlet to Filter Particle and Moisture to Prolong the Life of Pump.


<L1>

Spec.	Model	"PU-300"	
Speed of Evacuation		20 Lit/min	23 Lit/min
Max. Vacuum		105 mbar, 78.8 torr	
Motor Power & Max. Current		60 W, 0.3A	65 W, 0.7A
Motor Rotation & Noise Level		1450 rpm, 50 dB	1750 rpm, 50 dB
Hose Barb		8mm	
Size & Weight		268×135×h204mm, 4.1kg	
Electric Rqts.		220V, 50Hz	110V, 60Hz
Cat. No.		DH.Pu0300	DH.Pu0301
Price, Each. ₩			

“witeg” Cannon-Fenske Routine Flow Capillary Viscometer, for Transparent, Forming & Liquids Mixture, ASTM/ISO
With Individual Certificate of Calibration, Constant K-value, with Ring Mark, Constant, 기본형 캐논-펜스케 점도계, 투명 액상용



With Certificate,
Cannon-Fenske
Viscometer,
Routine-type



■ Supplied by “witeg” in Germany ■ With Certificated, Constant K-value.

* ASTM D445 and 446, ISO 3104 and 3105. Calibrated, Adjusted at 40 °C

* Transparent Liquids, Forming Liquids, Liquids Mixture with Highly Volatile Components

* Filling Quantity ; Approx. 7~10mℓ, Overall Approx. 245mm

<L1>


Cat. No.	Capillary No.	Range (mm ² / S)	Constant K (Approx.)	Ea. Price ₩
WI.8589101C	25	0.4~2	0.002	
WI.8589102C	50	0.8~4	0.004	
WI.8589103C	75	1.6~8	0.008	
WI.8589104C	100	3~15	0.015	
WI.8589105C	150	7~35	0.035	
WI.8589106C	200	20~100	0.1	
WI.8589107C	300	50~200	0.25	
WI.8589108C	350	100~500	0.5	
WI.8589109C	400	240~1200	1.2	
WI.8589110C	450	500~2500	2.5	
WI.8589111C	500	1600~8000	8	
WI.8589112C	600	4000~20000	20	

“witeg” Cannon-Fenske Reverse Flow Capillary Viscometer, for Opaque Liquids, ASTM/ISO

With Individual Certificate of Calibration, Constant K-value, with 3 Ring Marks, 2 Constants, 역류형 캐논-펜스케 점도계, 불투명 액상용



With Certificate,
Cannon-Fenske
Viscometer,
Reverse-type



■ Supplied by “witeg” in Germany ■ With Certificated, Constant K-value.

* ASTM D445 and 446, ISO 3104 and 3105. Calibrated, Adjusted at 40 °C

* for Opaque Liquids

* Filling Quantity ; Approx. 12mℓ, Overall Approx. 295mm

<L1>


Cat. No.	Capillary No.	Range (mm ² / S)	Constant K (Approx.)	Ea. Price ₩
WI.8587001C	25	0.4~2	0.002	
WI.8587002C	50	0.8~4	0.004	
WI.8587003C	75	1.6~8	0.008	
WI.8587004C	100	3~15	0.015	
WI.8587005C	150	7~35	0.035	
WI.8587006C	200	20~100	0.1	
WI.8587007C	300	50~200	0.25	
WI.8587008C	350	100~500	0.5	
WI.8587009C	400	240~1200	1.2	
WI.8587010C	450	500~2500	2.5	
WI.8587011C	500	1600~8000	8	
WI.8587012C	600	4000~20000	20	

“witeg” Ubbelohde Capillary Viscometer, Routine-Type, for Transparent Liquids & High-Temp or Low-Temp Measurements

With Individual Certificate of Calibration, Constant K-value, ASTM/ISO, 우베로데 점도계, 투명 액상용에 적합



With Certificate,
Ubbelohde
Viscometer,
Routine-type



■ Supplied by “witeg” in Germany ■ With Certificated, Constant K-value.

* ASTM D445 and 446, ISO 3104 and 3105. Calibrated, Adjusted at 40 °C

* Transparent Liquids, High-Temperature or Low-Temperature Measurements


* Filling Quantity ; Approx. 15~20mℓ, Overall Approx. 283mm

<L1>

Cat. No.	Capillary No.	Range (mm ² / S)	Constant K (Approx.)	Ea. Price ₩
WI.8589101	0	0.3-1	0.001	
WI.8589102	0C	0.6-3	0.003	
WI.8589103	0B	1-5	0.005	
WI.8589104	1	2-10	0.01	
WI.8589105	1C	6-30	0.03	
WI.8589106	1B	10-50	0.05	
WI.8589107	2	20-100	0.1	
WI.8589108	2C	60-300	0.3	
WI.8589109	2B	100-500	0.5	
WI.8589110	3	200-1000	1	
WI.8589111	3C	600-3000	3	
WI.8589112	3B	1000-5000	5	
WI.8589113	4	2000-10000	10	
WI.8589114	4C	6000-30000	30	
WI.8589115	4B	10000-50000	50	
WI.8589116	5	20000-100000	100	

Cannon-Fenske Routine Flow Capillary Viscometer, for Transparent, Forming & Liquids Mixture With Calibration Certificate, Constant K-value, with Ring Mark, ASTM/ISO, 기본형 캐논-펜스케 점도계, 투명 액상용

NEW



With Certificate,
Cannon-Fenske
Viscometer,
Routine-type



- Supplied by SciLab Korea ■ With Certificated, Constant K-value.
- 중력 흐름만을 이용하여 동점도 측정에 적합한 직관형 모세관 점도계
- 투명한 액체, 발포성 샘플, 휘발성 혼합물 등 일반 점도 측정에 사용
- * ASTM D445 and 446, ISO 3105. Calibrated, Adjusted at 40 °C
- * Transparent Liquids, Forming Liquids, Liquids Mixture with Highly Volatile Components
- * Filling Quantity ; Approx. 7~10ml, Overall Approx. 250mm

<L1>


Cat. No.	Capillary No.	Capillary Inner Diameter (mm)	Range (mm ² / S)	Constant K (Approx.)	Ea. Price ₩
SL.Met1001	25	0.3	0.4~2	0.002	
SL.Met1002	50	0.44	0.8~4	0.004	
SL.Met1003	75	0.54	1.6~8	0.008	
SL.Met1004	100	0.63	3~15	0.015	
SL.Met1005	150	0.78	7~35	0.035	
SL.Met1006	200	1.01	20~100	0.1	
SL.Met1007	300	1.27	50~200	0.25	
SL.Met1008	350	1.52	100~500	0.5	
SL.Met1009	400	1.92	240~1200	1.2	
SL.Met1010	450	2.35	500~2500	2.5	
SL.Met1011	500	3.2	1600~8000	8	
SL.Met1012	600	4.2	4000~20000	20	

Ubbelohde Capillary Viscometer, Routine-Type, for Transparent Liquids & High-Temp or Low-Temp Measurements With Calibration Certificate, Constant K-value, ASTM/ISO, 우베로데 점도계, 투명 액상용에 적합

NEW



With Certificate,
Ubbelohde
Viscometer,
Routine-type



- Supplied by SciLab Korea ■ With Certificated, Constant K-value.
- 중력 보정 가능한 구조로 외부 압력의 영향을 최소화한 고정밀 모세관 점도계
- 반복 측정에 우수하며, 고온 또는 저온 환경에서 정밀 측정에 적합
- * ASTM D445 and 446, ISO 3105. Calibrated, Adjusted at 40 °C
- * Transparent Liquids, High-Temperature or Low-Temperature Measurements
- * Filling Quantity ; Approx. 15~20ml, Overall Approx. 283mm

<L1>

Cat. No.	Capillary No.	Capillary Inner Diameter (mm)	Range (mm ² / S)	Constant K (Approx.)	Ea. Price ₩
SL.Met1031	0	0.24	0.3-1	0.001	
SL.Met1032	0C	0.36	0.6-3	0.003	
SL.Met1033	0B	0.46	1-5	0.005	
SL.Met1034	1	0.58	2-10	0.01	
SL.Met1035	1C	0.78	6-30	0.03	
SL.Met1036	1B	0.88	10-50	0.05	
SL.Met1037	2	1.03	20-100	0.1	
SL.Met1038	2C	1.36	60-300	0.3	
SL.Met1039	2B	1.55	100-500	0.5	
SL.Met1040	3	1.83	200-1000	1	
SL.Met1041	3C	2.43	600-3000	3	
SL.Met1042	3B	2.75	1000-5000	5	
SL.Met1043	4	3.27	2000-10000	10	
SL.Met1044	4C	4.32	6000-30000	30	
SL.Met1045	4B	5.2	10000-50000	50	
SL.Met1046	5	6.25	20000-100000	100	

DAIHAN™ Digital Precise Viscosity Bath “WVB”, for 5 Viscometers, Max. 30Lit/min, Up to 100℃, ±0.1℃
Stainless-steel Lid with 5 Holes for Viscometer Holders, Available Reverse & Routine-type Viscometer, Transparent Window
투시형 정밀 점도 항온수조, 5×점도계 사용가능, 디지털 퍼지 제어 시스템, 5-구 스텐리스 스틸 리드포함



Max. 30Lit/min Viscosity Bath “WVB-30”
Includes Flat Lid for Viscometer Holder.
(Viscometer Holder ; Option)



Fig1. Flat Lid with 5 Viscometer Holders



(1) Viscometer Holder “WVB02010” (Optional)

- Supplied by DAIHAN Korea ■ CE Certified ■ PL (Product Liability) Insurance
- with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Used for Precise Measurement with Capillary Viscometers
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Powerful Circulation Pump ensures the Highest Uniformity of Temperature
- Viscometer Holder ; Available Reverse & Routine-type Viscometer
- Patented Jog-Shuttle Control System and Digital LCD with Back-Light Function.
- * Stainless-steel Lid Contains Five Holes for Insertion of Viscometer Holders
- * User's Self-compensation Function to Control Difference between Real Temperature and Indicated Temperature
- * Storage Function for the Set Values of Temperature and Time
- * Locking Mode Supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Includes Flat Lid for Viscometer(without Viscometer Holder)
- * Over Temp., Over Current Protection and Sensor Error Detection
- * Digital Timer Function ; 'Delay Time' and 'Operation Time' ; Can be Controlled

Specifications & Ordering information

<L1>

Spec.	Model	“WVB-30” Max. Flow 30Lit/min
Range & Accuracy		Room Temp +5℃ to 100℃, ±0.1℃
Temp. Resolution		0.1℃-Display, 0.1℃-Control
Temp. Uniformity & Probe		±0.2℃, PT100
Heating Power		2 kW
Pump		Magnetic Drive, Max. Flow 30 Lit/min
Timer & Alarm		99hr 59min (Delay & Continuous run), Error status & Timer-end
Display		Digital LCD with Back-Light Function
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn+Push)
Safety Circuit		Over-Temp./-Current protection, Sensor Error Detection
Lid		Included 5 Holes Stainless-steel Lid for Viscometer Holder
Drain Valve		Mounted
Viewing Window Size(mm)		230×180
Material		Internal ; Stainless-steel (#304), External ; Powder Coated steel Window ; Vacuum Sealed and Tempered Safety Pair Glasses
Dimension (w×d×h)	Effective Space (mm)	330×300×h260mm
	Internal (mm)	350×300×h300mm
	External (mm)	585×375×h420mm
Packing Size (mm) & Gross Weight (kg)		640×450×h500mm, 24kg
Electric Reqs.		1 Phase, AC 120V, 60Hz or 230V, 50 / 60Hz
Cat. No.	230V	DH.WVB01030
	120V	DH.WVB02030
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

※ 운임 별도

Spare Stainless-steel Flat Lid for Viscometer Holder

DH.WVB11030 Flat Lid, for 30 Lit, w380×d330, 5×φ50mm Holes, “VBF030”

Ea. Price, ₩

Optional Viscometer Holder

DH.WVB02010 (1) Viscometer Holder, “WVB02010”
for Ostwald, Ubbelohde
DH.WVB02011 Viscometer Holder, “WVB02011”
for “witeg” Ubbelohde, Cannon-Fenske
DH.WVB02012 Viscometer Holder, “WVB02012”
for Thread-type Ubbelohde

Ea. Price, ₩

Viscosity Bath Validation (IQ, OQ) Service

DH.WVB31030 Validation Service(IQ, OQ), “VBF030”, with “WVB-30”

Ea. Price, ₩

DAIHAN™ Standard Rotary Viscometer-full Set "WVS-0.1M" & "WVS-2M", with Calibration Certificate, 1~2,000,000 cP
With Standard Spindle-kit(LV1~4), Lifting Stand, Hand Handling Case, 0.3~60rpm, 표준형 디지털 회전 점도계-풀세트

~ 2,000,000cP

3-Year Warranty




with Calibration Certificate



Hand Handling Case

Standard Rotary Viscometer Full Set "WVS"



Spindle-Rack / Stand

(1) Standard LV Spindle kit, LV1 ~ 4
Included for "WVS-0.1M" & "WVS-2M" Meter



(2) Ultra-Low Viscosity Adaptor



(3) RTD Temp. Probe



(4) Micro Printer, with Cable

- Supplied by DAIHAN Korea
- Ideal for measuring the Viscosity of Various Liquid such as Grease, Painting, Foodstuff, Dope, Paper making, Cosmetics, Chemical industry, Capsule Stickiness Agent and Medicines.
- Component of Set ; Viscometer Unit "WVS", Spindle Kit(LV 1~4), Lifting Stand, Hand Handling Case, Calibration Certificate
- Included ASTM Standard Spindle-kit of 4 Ea. ; LV1 ; 15~20K, LV2 ; 50~100K, LV3 ; 200~400K, LV4 ; 1K~2M
- Direct Readout of all Measurement Parameters on an illuminated LCD.
- The Viscometer is equipped with many Pre-sets speeds ; ("WVS-0.1M" ; 6- / 12- / 30- / 60- RPM, "WVS-2M" ; 0.3- / 0.6- / 1.5- / 3- / 6- / 12- / 30- / 60-RPM)
- Auto Range Function to determine the Maximum Viscosity with each Spindle / Speed Combination.
- Safety Protection Against Sudden Rises in the Power Supply.

- * Data on screen ; • Speed Selected ; RPM • Spindle Selected ; S.P.
- Viscosity Reading ; cP(mPa.s) • Percentage of Full Scale ; %
- Sample Temperature ; °C
- * Viscosity Reading ; Dynamic Viscosity (cP or mPa.s)
- * The Viscometer and the PT100 can be calibrated by the User.
- * Ultra-sensitive Front Keyboard that allows Easy Input of all Parameters.
- * Error Message in the Screen with Sound Alarm when the Viscosity of the Test Material exceeds the Full Scale of the selected Combination(SP / rpm.)
- * Sound Alarm if the Equipment is Working at Under 20% of Selected full Scale
- * Direct reading in cP(mPa.s) * without Limits
- * Can be Use with Various Optional Accessories ; Ultra-Low Viscosity Adaptor RTD Temp. Probe, Micro Printer

Specifications & Ordering information

<L1>

Spec.	Model	A. "WVS-0.1M" Max. Viscosity. 100,000 cP	B. "WVS-2M" Max. Viscosity. 2,000,000 cP
Viscosity Range, cP(mPa.s)		1~100,000 cP	1~2,000,000 cP
Speed Range (RPM)		6, 12, 30, 60	0.3, 0.6, 1.5, 3, 6, 12, 30, 60
Certificate		with Individual Calibration Certificate	
ASTM Standard Spindle-kit & Measure Viscosity Range ; cP		Included Standard Spindle-kit(LV1~LV4), Stainless-steel ① LV1 ; 15~20,000 cP ② LV2 ; 50~100,000 cP ③ LV3 ; 200~400,000 cP ④ LV4 ; 1,000~2,000,000 cP ※ LV0 Spindle for 1~20 cP(Optional)	
Measurement Error		± 5% (Newton Liquid)	
Accuracy		± 1.0% of range	
Repeatability		± 0.5%	
Ambient Temperature		5 °C ~ 35 °C Relative Humidity ; ≤80%	
Weight		4 kg	
Electrical Reqts.		AC 230V, 50/60Hz	
Cat. No.	230V	DH.WV1011	DH.WV1012
Price, Each. ₩			

* K=1,000 cP(mPa.s), M =1,000,000 cP(mPa.s)

Spare Spindle-kit

DH.WV2011 (1) Standard LV Spindle kit, LV1 ~ 4,
LV1 ; 15~20K, LV2 ; 50~100K, LV3 ; 200~400K, LV4 ; 1K~2M

Ea. Price, ₩

Optional Accessories

DH.WV2011.2 (2) Ultra-Low Viscosity Adaptor for 1~20cP (LV0 Spindle)
* Please order with the product!, Calibration is required because there is no Software and Data of LV0 in the Viscometer
DH.WV2041 (3) RTD Temp. Probe, 0.0~100.0 °C, Accu. ± 0.1, Resolution 0.1 °C
DH.WV2043 (4) Micro Printer, with RS-232, 1.5m Cable
(Print Information ; Temp., RPM, Spindle, Torque, Viscosity)

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ Multi-speed Rotary Viscometer-full Set "WVM-0.6M" & "WVM-6M", 1~6,000,000 cP


With Cal. Certi. & Standard Spindle kit(LV1~4), Lifting Stand, Hand Handling Case, 0.1~60/100-rpm, 멀티 스피드 디지털 회전 점도계

~ 6,000,000cP

3-Year Warranty



Multi-speed Rotary Viscometer Full Set "WVM"

(1) Standard LV Spindle kit, LV1~4
Included for "WVM-0.6M" & "WVM-6M" Meter

(2) Ultra-Low Viscosity Adaptor



(3) RTD Temp. Probe



(4) Small Sample Adapter



(5) Micro Printer, with Cable

- Supplied by DAIHAN Korea
- Ideal for measuring the Viscosity of Various Liquid such as Grease, Painting, Foodstuff, Dope, Paper making, Cosmetics, Chemical industry, Capsule Stickiness Agent and Medicines.
- Component of Set ; Viscometer Unit "WVM", Spindle Kit(LV 1~4), Lifting Stand, Hand Handling Case, Calibration Certificate
- Included ASTM Standard Spindle-kit of 4 Ea. ; LV1 ; 15~20K, LV2 ; 50~100K, LV3 ; 200~400K, LV4 ; 1K~2M
- Direct Readout of all Measurement Parameters on an illuminated LCD.
- Available Pre-sets speeds, 0.1~100 rpm
- Auto Range Function to determine the Maximum Viscosity with each Spindle / Speed Combination.
- Safety Protection Against Sudden Rises in the Power Supply.

- * Data on screen ; • Speed Selected ; RPM • Spindle Selected ; S.P.
- Viscosity Reading ; cP(mPa.s) • Percentage of Full Scale ; %
- Sample Temperature ; °C
- * Viscosity Reading ; Dynamic Viscosity (cP or mPa.s)
- * The Viscometer and the PT100 can be calibrated by the User.
- * Ultra-sensitive Front Keyboard that allows Easy Input of all Parameters.
- * Error Message in the Screen with Sound Alarm when the Viscosity of the Test Material exceeds the Full Scale of the selected Combination(SP / rpm.)
- * Sound Alarm if the Equipment is Working at Under 20% of Selected full Scale
- * Direct reading in cP(mPa.s) * without Limits
- * Can be Use with Various Optional Accessories ; Ultra-Low Viscosity Adaptor RTD Temp. Probe, Micro Printer, Small Sample Adapter

Specifications & Ordering information

<L1>

Spec.	Model	A. "WVM-0.6M" Max. Viscosity. 600,000cP	B. "WVM-6M" Max. Viscosity. 6,000,000cP
Viscosity Range cP(mPa.s)		1~600,000cP	1~6,000,000cP
Speed Range (RPM)		1~60	0.1~100
Certificate		with Individual Calibration Certificate	
Certificate		Included Standard Spindle-kit(LV1~LV4), Stainless-steel ① LV1 ; 15~20,000cP ② LV2 ; 50~100,000cP ③ LV3 ; 200~400,000 cP ④ LV4 ; 1,000~2,000,000cP ※ LV0 Spindle for 1~20cP(Optional)	
Measurement Error		± 5% (Newton Liquid)	
Accuracy		± 1.0% of range	
Repeatability		± 0.5%	
Ambient Temperature		5°C~35°C Relative Humidity ; ≤ 80%	
Weight		4 kg	
Electrical Reqts.		AC 230V, 50/60Hz	
Cat. No.		DH.WV1021	DH.WV1022
Price, Each. ₩			

* K=1,000 cP(mPa.s), M=1,000,000 cP(mPa.s)

Spare Spindle-kit

DH.WV2011 (1) Standard LV Spindle kit, LV1~4,
LV1 ; 15~20K, LV2 ; 50~100K, LV3 ; 200~400K, LV4 ; 1K~2M

Ea. Price, ₩

Optional Accessories

DH.WV2011.2 (2) Ultra-Low Viscosity Adaptor for 1~20cP (LV0 Spindle)
* Please order with the product!, Calibration is required
because there is no Software and Data of LV0 in the Viscometer

DH.WV2041 (3) RTD Temp. Probe, 0.0~100.0°C, Accu. ±0.1, Resolution 0.1°C

DH.WV2042 (4) Small Sample Adapter, Stainless-steel, 2~16ml
* Please order with the product!, Calibration is required
because there is no Software and Data of SSA in the Viscometer

DH.WV2043 (5) Micro Printer, with RS-232, 1.5m Cable
(Print Information ; Temp., RPM, Spindle, Torque, Viscosity)

Ea. Price, ₩

Ea. Price, ₩

Ea. Price, ₩

DAIHAN™ Glass/Labware Washer "WSH400.L"

Ideal for Laboratory, etc.

Glass/Labware Washer

HIGH
PERFORMANCE
CLEANING

USER
SAFETY

PROGRAM
CYCLE

High-Performance Cleaning

- Outstanding Fast Cycle and High-Performance Cleaning results from Efficient Direct Injection Washing.
- High performance Pump ensure High Flow Rate of Water.
- Efficient Water Spray Pressure bring about Better and Faster Cleaning.
- Water, Liquid Detergent and Distilled Water Automatically Injected by Program Washing Cycle.

Safety Door Lock

- Door is Automatically Locked for User Safety when Washing Cycle Starts.

Compact, Robust, and Efficient Construction

- Compact design for Space Saving and Maximize Working Capacity in Laboratory.
- Double Layer Structure ensures Soundproof and Heatproof during Operation.
- Built-in Steam Condenser Prevents Steam to Enter into the Laboratory.
- Telescopic Sliding Rails ensures Safe and Easy to Load and Unload.

Efficient Washing Cycle Program

- Provide 7 Program Washing Cycles ; Pre-set 6 Programs & Customized 1 Program.

Pre-set 6 Program Cycles Setting Value

Pre-set 6 Washing Cycle	PreWash	Light Wash	Standard Wash	Light Demi	Standard Demi	Intense Wash
Tap Water Wash	3min	-	-	-	-	-
Detergent Wash	-	60℃ / 8min	75℃ / 8min	60℃ / 5min	75℃ / 5min	80℃ / 10min
Rinse Wash	-	75℃ / 5min	80℃ / 5min	75℃ / 3min	80℃ / 3min	80℃ / 5min
Distilled Water Wash	-	-	-	2 cycle 3min / cycle	2 cycle 3min / cycle	2 cycle 2min / cycle



Irrigation Nozzle



Spray Arm

Various Interchangeable Racks & Baskets

- Includes 1 Stainless-steel Wire Basket.
- Various Custom Racks ; Basket, Pipette Rack, Volumetric Flask Rack, Direct Injection Kit, etc.



DAIHAN™ Glass/Labware Washer “WSH400.L”, Max 85 °C, Washing Surface 0.52m², 400Lit/min, Max 5bar
With 1 Stainless-steel Wire Basket, Water · Detergent · Distilled Water-Auto Injection, Program Washing Cycle, Double Wall, Telescopic Rail, 초자 세척기



3-Year Warranty

Glass/Labware Washer “WSH400.L”



3-Year Warranty

Glass/Labware Washer “WSH400.L”
with Customizing Rack (Option)

(1) STS Wire Basket



(2) STS Injection Rack



(3) STS Rack, for Vessel



(4) STS Rack, for Test Tube



(5) STS Rack, Mesh-type



(6) Cleaner

(7) Rinse

Optional Liquid Detergent

- Supplied by DAIHAN Korea ■ PL(Product Liability) Insurance.
 - with Certificate & Traceability ; Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System.
 - Ideal for Clinical, Pharmaceutical, Bioresearch, Chemical Research, etc.
 - Outstanding Fast Washing Cycle and High-performance Cleaning Results from Efficient Direct Injection Washing.
 - Water-/Detergent-/Distilled Water- Automatically Injected into Wash Chamber via Pumps.
 - Double Layer Wall ensures Soundproof and Heatproof.
 - Automatic Door Lock system for User Safety.
 - Special Steam Condensing Mechanism for User Safety.
- * Built-in Telescopic Sliding Rail for Safety and Easy to Use.
* Provides Efficient Program Washing Cycle ; Pre-set 6Programs and Customizable 1Program.
* Built-in Steam Condenser Prevents Steam to enter into the Laboratory.
* Customizable Various Rack, Basket, Pipette Rack, Volumetric Flask Rack, Direct Injection Kit, etc.

Specifications & Ordering information

<L1>

Spec.	Model	“WSH400.L”
Total Washing Surface		0.52m ²
Temp Range		Room Temp +5 °C ~85 °C
Display		Digital LCD with Back-Light Function
Controller		Digital Fuzzy-Control System (Fully Automatic), Touch Button UI
Programmable Function		Pre-set 6 Programs & Customizable 1 Programs
Wash Pump		Recirculation ; 400Lit/min, Water Pressure ; 3~5bar, Water Consumption ; 8~10Lit
Rack		Included 1ea Stainless-steel Standard Rack
Safety Circuit		Over-Temp/-Pressure/-Current protection, Sensor Error Detection
Material		Stainless-steel (#304)-Internal, Stainless-steel (#304)-External
Dimension (w × d × h)	Interior	508 × 515 × h407 mm
	Exterior	600 × 625 × h870 mm
Power Consumption		5 kW
Net Weight		95 kg
Electric Reqts.		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 or 60 Hz
Cat. No.	60Hz	230V DH.WGW0400
		120V DH.WGW0401
	50Hz	230V DH.WGW5400
Price, Each. ₩		

* Other Specifications are available upon Customer's Request.

※ Electrical Work (Earth Leakage Circuit Breaker 30A) and Water Drainage System should be Prepared in Advance for Installation.

※ Note! Freight and Installation Charges are Extra.

설치 운임 - 서울 · 경기 ; ₩ 300,000 충청 · 강원 영서 ; ₩ 350,000

대구 · 경북 · 전북 · 강원 영동 ; ₩ 400,000 포항 · 경남 · 전남 ; ₩ 500,000

(※ 공사비가 아닌 설치 · 운임비용입니다.)

Spare Stainless-steel Standard Rack

DH.WGW01001 (1) STS Wire Basket “WB495”, for Rail, 495 × 510mm

Ea. Price, ₩

Optional Stainless-steel Injection Rack for Various Vessel

DH.WGW01002 (2) STS Injection Rack “IR40” for Various Vessel, 495 × 510mm

Ea. Price, ₩

Optional Stainless-steel Wire Rack, Use for Included Standard Rack

DH.WGW01003 (3) Wire Rack “WRV12”, for Various Vessel, 248 × 460 × h180mm

DH.WGW01004 (4) “WRT56”, for Test Tube, 248 × 460 × h160mm

DH.WGW01005 (5) “WRM248”, Mesh-type, 248 × 460 × h65mm

Ea. Price, ₩

Optional Liquid Cleaning Detergent & Rinse

DH.WGW02010 (6) Potassium Hydroxide Liquid Cleaning Detergent, 18Lit

DH.WGW02011 (7) Nonionic Surfactant Liquid Rinse, 18Lit

Ea. Price, ₩

Optional Tube & Bottle for Distilling Water Supply

SL.Bot3105 Bottle, Aspiration, White, 50L, Neck id-φ95, φ380 × h700mm

TUB.Si020 Tubing, 10-size ; id/od-φ10/13 × thick-1.5mm

0338-403 Clamp, Width 8, Grip-φ6 ~ 16mm

Ea. Price, ₩

DAIHAN™ Electric Classic Water Still “WWS”, All Stainless-steel, Built-in Prefilter & Auto On/Off system

For Laboratory Water, Pyrogen-Free, $0.3\text{M}\Omega\cdot\text{cm}$ (resistivity), pH5.4 to 7.3, Capa. 4 · 8 · 12 Lit/hr

(1) with Stainless-steel 100 Mesh Prefilter, (2) with PP $1\mu\text{m}$ Pore Cartridge Prefilter, 필터부착형 증류수 제조기

3-Year Warranty



(1) with Stainless-steel 100 Mesh Prefilter, “WWS-S”; 4 · 8 · 12 Lit/hr

- Supplied by DAIHAN Korea
- 3-Year Warranty
- High-performance / All Stainless-steel (STS-430/304) models ; Designed for the Durability & Consistent Quality
- “WWS-S” models, Built-in Stainless-steel(STS-304) Prefilter Kit, interchangeable 100 mesh ($150\mu\text{m}$) ; Remove Particles & Scale-forming impurity.
- “WWS-P” models, Built-in PP(Polypropylene) $1\mu\text{m}$ Prefilter / Cartridge Kit, interchangeable ; Remove Rust, Particles & Scale-forming impurity.(Filter Replace Period ; 20~40 Day)
- All Models, Complete Set with ;
 - ① Stainless-steel Unit Kit,
 - ② External Consistent-Level Device,
 - ③ Drain Valve, ④ Prefilter Kit,
 - ⑤ 2m Feeding Water Tubing (Silicone) Assembly for Connection of Water Tap $\phi 15$ or $\phi 19$ mm,
 - ⑥ PP Un-equal Tapered Connector $\phi 8, 10, 12 / 14, 16\text{mm}$ -for Different Diameter Tubings
- Auto On/Off with Water Level Sensor/Micro Level Sensor.

- * This Space-Saving Horizontal Condenser was Designed for Durability and Consistent Quality
- * Energy Saving Features include a Condenser that Preheats the Feedwater and a Double-Walled Evaporator.
- * A Low Water Cutoff protects Heating Elements from Burn-Out when Water Supply is Interrupted
- * Stainless-steel Heater Built-in 0.4mm Ni/Cr. Element, Boiler, Efficient Energy and Water Consumption.
- * Sealing ; Silicone O-ring packing.

3-Year Warranty



(2) with PP $1\mu\text{m}$ Pore Cartridge Prefilter, “WWS-P”; 4 · 8 · 12 Lit/hr



Specifications & Ordering information

<L1>

Spec. / Model		(1) With Stainless steel 100 Mesh Prefilter			(2) With PP 1 μ m Pore Cartridge Prefilter		
		“WWS-S4”	“WWS-S8”	“WWS-S12”	“WWS-P4”	“WWS-P8”	“WWS-P12”
Capacity(Water)		4 Lit/hr	8 Lit/hr	12 Lit/hr	4 Lit/hr	8 Lit/hr	12 Lit/hr
Filter	Function	Remove Particles & Scale-forming impurity			Remove Rust, Particles & Scale-forming impurity		
	Material	Stainless-steel, STS304			PP(Polypropylene)		
	Housing	POM(Polyoxymethylene)			PE(Polyethylene)		
	Pore size	150 μ m(100 mesh)			1 μ m		
Heat-resisting		-29℃ + 85℃			85℃		
Water Purity		0.3M Ω ·cm (resistivity), pH5.4 to 7.3					
Level Sensor		Auto On/Off with Water Level Sensor / Micro Level Sensor					
Material	Boiler	STS-430	STS-430		STS-430	STS-430	
	Condenser	STS-304	STS-304		STS-304	STS-304	
	Base	STS-430	Aluminum		STS-304	Aluminum	
Size	Condenser(only)	$\phi 85\times\text{L}375$	$\phi 120\times\text{L}430$	$\phi 120\times\text{L}490$	$\phi 85\times\text{L}375$	$\phi 120\times\text{L}430$	$\phi 120\times\text{L}490$
	Body	$\phi 350\times\text{h}670$	$\phi 350\times\text{h}760$	$\phi 350\times\text{h}760$	$\phi 400\times\text{h}670$	$\phi 400\times\text{h}760$	$\phi 400\times\text{h}760$
Heater		3kW	6kW (3kW $\times$ 2Ea)	9kW (3kW $\times$ 3Ea)	3kW	6kW (3kW $\times$ 2Ea)	9kW (3kW $\times$ 3Ea)
Plug Style		Plug Code “C”	Plug Code “F”		Plug Code “C”	Plug Code “F”	
Electrical Reqts.		1 Phase, AC 220V, 60Hz					
Cat. No.	220V	DH.WWS01004	DH.WWS01008	DH.WWS01012	DH.WWS02004	DH.WWS02008	DH.WWS02012
Price, Each. ₩							

Spare Prefilter & Optional Storage Bottle

- DH.WWS10101 (a) Spare, Stainless-steel 100 mesh Prefilter Kit, “SF101”, POM Housing $\phi 36 \times \text{h}61\text{mm}$, STS-304 Filter only $\phi 22 \times \text{h}46\text{mm}$
- DH.WWS10102 (b) , PP $1\mu\text{m}$ pore Cartridge Prefilter & Housing Kit, “PF102”, PE Housing $\phi 110 \times \text{h}300\text{mm}$, PP Filter only $\phi 60 \times \text{h}252\text{mm}$
- DH.WWP.10PF Spare, PP $1\mu\text{m}$ pore Cartridge Prefilter, “PF103”, ϕ (od.) PP Filter $60 \times \phi$ (id.) $30 \times \text{h}252\text{mm}$
- DH.WWS10104 , PE Housing only, “PF104”, $\phi 110 \times \text{h}300\text{mm}$
- SL.Bot3102 Optional, Bottle, HDPE, Aspiration, White, 5L, Neck id- $\phi 70$, $\phi 200 \times \text{h}360\text{mm}$
- SL.Bot3103 , 10L, - $\phi 70$, $\phi 250 \times \text{h}460\text{mm}$
- SL.Bot3104 , 20L, - $\phi 95$, $\phi 300 \times \text{h}555\text{mm}$
- SL.Bot3105 , 50L, - $\phi 95$, $\phi 380 \times \text{h}700\text{mm}$

Ea. Price, ₩

DAIHAN™ 3.5 & 7.5 Lit/hr Electric Chamber Water Still "CWS", All Stainless-steel, Auto On/Off System

Good for Laboratory Water, Double Wall Still Chamber, Built-in PP 1 μ m pore Cartridge Prefilter, 스텐레스 사각 증류수 제조기, 필터부착형

3.5 & 7.5 Lit/hr Cabinet-type Water Still "CWS"



Plug Code "C"



Plug Code "F"



(1) Cartridge Prefilter & PE Housing



(2) HDPE Storage Bottle 5~50L

- Supplied by DAIHAN Korea
- All Stainless-steel for Corrosion Resistance
- Double-walled Chamber Structure for Protects Burn Injuries
- High-performance & Compact Design

- * Designed for the Durability & Consistent Quality
- * Built-in Interchangeable PP(Polypropylene) 1 μ m Prefilter Cartridge Kit
- * Condenser and Boiler are Separable for Easy Cleaning and Maintenance
- * Auto Heater On/Off with Water Level Sensor & OPT Sensor
- * Pre-filter Removes Rust, Particles & Scale-forming Impurity
- * Economical Energy and Water Consumption
- * Built in Low water Sensor for Protects Heating Elements
- * Sealing ; Silicone Gasket

Specifications & Ordering information

<L2>

Spec.	Model	"CWS-4"	"CWS-8"
Capacity		Max. 3.5 Lit / hr	Max. 7.5 Lit / hr
Cooling Water		30~50 Lit	60~80 Lit
Filter	Function	Remove Rust, Particles & Scale-forming Impurity	
	Material	PP (Polypropylene)	
	Housing	PC (Polycarbonate) Bottle + PP (Polypropylene) Lid	
	Pore size	1 μ m	
	Heat-resisting	85 $^{\circ}$ C	
Water Purity		0.3M Ω -cm (resistivity), pH5.4 to 7.3	
Level Sensor		Auto On / Off with Water Level Sensor / Micro Level Sensor	
Material	Interior	Stainless-steel	
	Exterior	Stainless-steel Hair Lined with Powder Coating	
	Heater	Stainless-steel Pipe	
Safety		Low Water & Over Temperature Protection, Electric Leakage Circuit Breaker	
Dimension(w×d×h)		430×315×h570 mm	445×325×h630 mm
Net Weight		12 kg	18 kg
Packing Size(w×d×h)		440×370×h610 mm	450×380×h670 mm
Gross Weight		13 kg	19 kg
Power Consumption		3 kW	6 kW
Plug Style		Plug Code "C"	Plug Code "F"
Electric Reqts.		220 V, 50 / 60 Hz, 1 Phase	
Cat. No.		DH.WatS8001	DH.WatS8002
Price, Each. ₩			

* Other Specifications are available upon Customer's Request.

※ 운임 별도

(1) Spare Prefilter

- DH.WatA1001 PP 1 μ m pore Cartridge Prefilter Housing Kit, "CWS-A1001"
 DH.WatA1002 PP 1 μ m pore Cartridge Prefilter Only, "CWS-A1002"
 DH.WatA1003 Housing Only, PC Bottle & PP Lid "CWS-A1003"

Ea. Price, ₩

(2) Optional HDPE Storage Bottle


- SL.Bot3102 Bottle, Aspiration, White, 5L, Neck id- ϕ 70, ϕ 200×h360 mm
 SL.Bot3103, 10L, - ϕ 70, ϕ 250×h460 mm
 SL.Bot3104, 20L, - ϕ 95, ϕ 300×h555 mm
 SL.Bot3105, 50L, - ϕ 95, ϕ 380×h700 mm

Ea. Price, ₩

(3) Optional Inlet/Outlet Tube & Clamp

- TUB.PV034 (A) Tubing, PVC, id/od- ϕ 10/14×thick-2mm, for Feeding Water
 TUB.Si020 (B) Tubing, 10-size; id.od- ϕ 10/13 thick-1.5mm, for Distilled & Cooling Water
 0338-405 (C) Clamp, width 8. Grip- ϕ 6~22mm

	m/dh.pk. Price, ₩	Roll. Price, ₩
(A) Tubing, PVC, id/od- ϕ 10/14×thick-2mm, for Feeding Water	1m	90m
(B) Tubing, 10-size; id.od- ϕ 10/13 thick-1.5mm, for Distilled & Cooling Water	1m	10m
(C) Clamp, width 8. Grip- ϕ 6~22mm		Pk. Price, ₩
		10ea



DAIHAN™ Water Purification System “New-P.NIX RO(Reverse Osmosis)”, 15 · 25 · 35L/hr
With Pretreatment System, 2-Steps of Filter Exchange Indicator, (RO) Product 0.2~30 μS/cm, 순수 제조기



3-Year Warranty



Included
Pretreatment
System
'HMC-PT1'



(1) White, HDPE
Aspiration Bottle

(2) Amber, HDPE
Aspiration Bottle



(3) 14" Super A/C Filter



(4) 10" Prefilter



(5) P.NIX RO Pack I Filter
for "New-P.NIX RO 180"



(6), (7) P.NIX RO Pack II, III Filter
for "New-P.NIX RO 280", "New-P.NIX RO 380"

- Supplied by DAIHAN Korea ■ Certified CE, ISO 9001, ISO 14001
- Built-In RS 232C Port and Volumetric Dispensing Function by Timer, Included Level sensor
- Auto Regular Recycling to Maintain the Optimum water Quality Level
- Programming by 2 Steps for the Filter Exchange
- Filter Replace Message Indicating by LCD Display & Blinking LED Light

- * Microprocess One Touch Technology Controller
- * 16 × 4 Characters Backlight LCD Digital Display
- * Awarded the Korean Excellent Products
- * Included Pretreatment System HMC-PT1, 14" Super A/C Filter, RO Pack

Specifications & Ordering information

<L2>

Spec.	Model	"New-P.NIX RO 180" Max. Daily Volume. (Max. 60L/day)	"New-P.NIX RO 280" Max. Daily Volume. (Max. 100L/day)	"New-P.NIX RO 380" Max. Daily Volume. (Max. 120L/day)
Certified Mark	CE, ISO 9001, ISO 14001			
Product Flow Rate(RO)		Max. 15L/hr	Max. 25L/hr	Max. 35L/hr
Daily Volume		60L/day	100L/day	120L/day
Water Quality Display		• RO Product ; 0.2 to 250.0 μS/cm • Feed water ; 10 to 1000 μS/cm • Temperature ; 0°C to 99.9°C		
Display		16 × 4 Characters Backlight LCD Digital Display		
Controller		Microprocess One Touch Technology		
Dimensions (w × d × h)		300 × 460 × h470mm		
Power		210~230V, 50/60 Hz, 60W		
Cat. No.		DH.WWP.RO180	DH.WWP.RO280	DH.WWP.RO380
Price, Each. ₩				

Optional Accessory

SL.Bot3102	(1) Bottle, HDPE, Aspiration, White,	5L, Neck id-φ70, φ200 × h360mm
SL.Bot3103		, 10L, -φ70, φ250 × h460mm
SL.Bot3104		, 20L, -φ95, φ300 × h555mm
SL.Bot3105		, 50L, -φ95, φ380 × h700mm
SL.Bot3106	(2) Bottle, HDPE, Aspiration, Amber,	5L, Neck id-φ70, φ200 × h360mm
SL.Bot3107		, 10L, -φ70, φ250 × h460mm
SL.Bot3108		, 20L, -φ95, φ300 × h555mm
SL.Bot3109		, 50L, -φ95, φ380 × h700mm

Ea.	Price, ₩

Spare Filter

Cat. No.	Description	Filter Replace Period	Ea. Price, ₩
DH.WWP.PR14AC	(3) 14" Super A/C Filter	Approx. 3~4 Month	
DH.WWP.10PF	(4) 10" Prefilter for HMC-PT1	20 ~ 40 Day	
DH.WWP.ROM1	(5) P.NIX RO Pack I for DH.WWP.RO180	Approx. 1 Year	
DH.WWP.ROM2	(6) P.NIX RO Pack II for DH.WWP.RO280		
DH.WWP.ROM3	(7) P.NIX RO Pack III for DH.WWP.RO380		

- The above mentioned Filter replacement period can be changed according to various factors such as feed water condition, producing volume and so on.

DAIHAN™ Ultra Pure(UP) Water Purification System “New-PNIX UP 900”, Max 1.5L/min

With Total Organic Carbon(TOC), 2-Steps of Filter Exchange Indicator, (UP) Up to 18.3MΩ·cm, 초순수 제조기



3-Year Warranty



(1) White, HDPE Aspiration Bottle

(2) Amber, HDPE Aspiration Bottle



(3) 0.2μm Final Capsule Filter



(4), (5) P.NIX UP Pack Filter

- Supplied by DAIHAN Korea
- Certified CE, ISO 9001, ISO 14001
- RS 232C Port and Volumetric Dispensing Function by Timer
- Auto Regular Recycling to Maintain the Optimum water Quality Level
- Programming by 2 Steps for the Filter Exchange
- Filter Replace Message Indicating by LCD Display & Blinking LED Light

- * Microprocess One Touch Technology Controller
 - * 16×4 Characters Backlight LCD Digital Display
 - * Obtain GD, QD Mark, Awarded the Korean Excellent Products
 - * Included 0.2μm Final Capsule Filter, UP Membrane Pack
 - * This Device needs Pure Water entering.
- (Please see the Preceding Pages for the Water Purification System.)

Specifications & Ordering information

<L2>

Spec.	Model	“New-P.NIX UP 900”			
		Scholar-type	UV-type	Bio-type	Integrate-type
Certified Mark		CE, ISO 9001, ISO 14001			
Product Flow Rate(UP)		Max. 1.5L/min			
Water Quality Display		• UP Product ; Up to 18.3MΩ·cm • Feed water ; 0 to 250μS/cm • Temperature ; 0℃ to 99.9℃			
Display		16×4 Characters Backlight LCD Digital Display			
Included Filter Composition		• 0.2μm Final Capsule Filter • UP Pack (DH.WWP.UPS)	• 0.2μm Final Capsule Filter • UP Pack (DH.WWP.UPU) • 254/185UV Lamp	• 0.2μm Final Capsule Filter • UP Pack (DH.WWP.UPS) • MW 5000 UF Filter(L)	• 0.2μm Final Capsule Filter • UP Pack(DH.WWP.UPU) • 254/185UV Lamp • MW 5000 UF Filter(L)
Controller		Microprocess One Touch Technology			
Dimensions (w×d×h)		300×500×h470mm			
Power		210~230V, 50/60 Hz			
Cat. No.		DH.WWP.UP.S	DH.WWP.UP.U	DH.WWP.UP.B	DH.WWP.UP.I
Price, Each. ₩					

- TOC monitoring system is also available for the Integrate-type on request. (TOC monitoring display range ; 0~200ppb)

Optional Accessory

SL.Bot3102	(1) Bottle, HDPE, Aspiration, White, 5L, Neck id-φ70, φ200×h360mm
SL.Bot3103	, 10L, -φ70, φ250×h460mm
SL.Bot3104	, 20L, -φ95, φ300×h555mm
SL.Bot3105	, 50L, -φ95, φ380×h700mm
SL.Bot3106	(2) Bottle, HDPE, Aspiration, Amber, 5L, Neck id-φ70, φ200×h360mm
SL.Bot3107	, 10L, -φ70, φ250×h460mm
SL.Bot3108	, 20L, -φ95, φ300×h555mm
SL.Bot3109	, 50L, -φ95, φ380×h700mm

Ea.	Price, ₩

Spare Filter & UV Lamp

Cat. No.	Description	Filter Replace Period	Ea. Price, ₩
DH.WWP.0.2CF	(3) 0.2μm Final Capsule Filter	Approx. 6~8 Month	
DH.WWP.UPS	(4) P.NIX UP Pack, Scholar / Bio	Approx. 1 Year	
DH.WWP.UPU	(5) P.NIX UP Pack, Scholar-UV / Integrate		
DH.WWP.MWL	MW 5000 UF Filter(L)	Approx. 6~8 Month	
DH.WWP.UV	254/185 UV Lamp	Approx. 3 Year	

- The above mentioned Filter replacement period can be changed according to various factors such as feed water condition, producing volume and so on.

DAIHAN™ Ultra Pure(UP) & Pure(RO) Water Purification All-in-One System “New-PNIX Power I”, Max 15L/hr
With Pretreatment System, 2-Steps of Filter Exchange Indicator, (RO) 0.2~250 μ S/cm, (UP) Up to 18.3M Ω ·cm, 초순수/순수 제조기



3-Year Warranty



Included Pretreatment System 'HMC-PT-3'



(3) 10" A/C Filter (4) 10" Prefilter



(5) P.NIX RO Pack I Filter



(6) 0.2 μ m Final Capsule Filter



(7), (8) P.NIX UP Pack Filter

- Supplied by DAIHAN Korea ■ Certified CE, ISO 9001, ISO 14001
- Built-In RS 232C Port and Volumetric Dispensing Function by Timer, Included Level sensor
- Auto Regular Recycling to Maintain the Optimum water Quality Level
- Programming by 2 Steps for the Filter Exchange
- Filter Replace Message Indicating by LCD Display & Blinking LED Light
- Self-Diagnosis for Ext- & Int-System, Auto Cleaning, Auto Sensor Calibration
- * Microprocess One Touch Technology Controller * 16×4 Characters Backlight LCD Digital Display
- * Awarded the Korean Excellent Products * Daily Volume ; 60L/day
- * Included Pretreatment System HMC-PT-3, RO Pack, UP Pack, 0.2 μ m Final Filter

Specifications & Ordering information

<L2>

Spec.	Model	“New-P.NIX Power I”			
		Scholar-type	UV-type	Bio-type	Integrate-type
Certified Mark		CE, ISO 9001, ISO 14001			
Product Flow Rate (RO & UP)		Max. 15L/hr			
Water Quality Display		• Pure water(RO) ; 0.2 ~ 250 μ S/cm • Ultra Pure water(UP) ; Up to $\leq$ 18.3M Ω ·cm • Feed water ; 0 to 990 μ S/cm • Temperature ; 0℃ to 99.9℃			
Included Filter Composition		• Pretreatment System HMC-PT-3 • RO Pack • 0.2 μ m Final Filter • UP Pack(Scholar)	• Pretreatment System HMC-PT-3 • RO Pack • 0.2 μ m Final Filter • UP Pack(S-UV) • 254/185 UV Lamp	• Pretreatment System HMC-PT-3 • RO Pack • 0.2 μ m Final Filter • UP Pack(BIO) • MW5000UF filter(M)	• Pretreatment System HMC-PT-3 • RO Pack • 0.2 μ m Final Filter • UP Pack(INT) • 254/185 UV Lamp • MW5000UF filter(M)
Display		16×4 Characters Backlight LCD Digital Display			
Controller		Microprocess One Touch Technology			
Dimensions (w×d×h)		395×500×h470mm			
Power		210~230V, 50/60 Hz			
Cat. No.		DH.WWP.POWER.S	DH.WWP.POWER.U	DH.WWP.POWER.B	DH.WWP.POWER.I
Price, Each. ₩					

• Product Flow Rate 25L/hr and 35L/hr are also available on Request.

Optional Accessory

SL.Bot3102	(1) Bottle, HDPE, Aspiration, White, 5L, Neck id- ϕ 70, ϕ 200×h360mm
SL.Bot3103	, 10L, - ϕ 70, ϕ 250×h460mm
SL.Bot3104	, 20L, - ϕ 95, ϕ 300×h555mm
SL.Bot3105	, 50L, - ϕ 95, ϕ 380×h700mm
SL.Bot3106	(2) Bottle, HDPE, Aspiration, Amber, 5L, Neck id- ϕ 70, ϕ 200×h360mm
SL.Bot3107	, 10L, - ϕ 70, ϕ 250×h460mm
SL.Bot3108	, 20L, - ϕ 95, ϕ 300×h555mm
SL.Bot3109	, 50L, - ϕ 95, ϕ 380×h700mm

Ea.	Price, ₩

Spare Filter & UV Lamp

Cat. No.	Description	Filter Replace Period	Ea. Price, ₩
DH.WWP.PR10AC	(3) 10" A/C Filter	Approx. 3~4 Month	
DH.WWP.10PF	(4) 10" Prefilter for HMC-PT-3	20 ~ 40 Day	
DH.WWP.ROM1	(5) P.NIX RO Pack I	Approx. 1 Year	
DH.WWP.0.2CF	(6) 0.2 μ m Final Capsule Filter	Approx. 6~8 Month	
DH.WWP.UPS	(7) P.NIX UP Pack, Scholar/Bio	Approx. 1 Year	
DH.WWP.UPU	(8) P.NIX UP Pack, Scholar-UV/Integrate	Approx. 1 Year	
DH.WWP.MWM	MW5000UF Filter (M)	Approx. 6~8 Month	
DH.WWP.UV	UV Lamp	Using time ; 2000 Hrs	

• The above mentioned Filter replacement period can be changed according to various factors such as feed water condition, producing volume and so on.

Handwriting practice area with horizontal dashed lines.



Handwriting practice lines consisting of 20 sets of three horizontal dashed lines.

CERTIFICATES

- ✓ Quality
- ✓ Safety
- ✓ Environment Management

- CE** EMC & LVD Compliance
The Certified DAIHAN Laboratory Instruments are Specified in this Catalog
- UL** UL Certification
Laboratory-use Electrical Equipment
for General-type Incubators & Oven, B.O.D Incubators, Vacuum Ovens, Hotplate Stirrers, ULT Freezers
- GD (Good Design) Mark**
- GMP** GMP Certification
Good Manufacturing Practices of Medical Device
- ISO 9001** System Certification
Design, Development, Manufacture & Servicing
for All DAIHAN Laboratory Instruments
- ISO 13485** System Certification
Design, Development, Manufacture & Servicing
for All DAIHAN Laboratory Instruments
- SQ (Service Qualification) Certification**



DAIHAN-Sci.com Shopping-mall
ALLforLAB.com Open-market

DAIHANmedical.com
DAIHANused.com
LabRental.co.kr
Lab-Service.co.kr
Pipettermall.com

Calibration.kr (ml & g)

검 색 대 한 과 학



Since 1980



Your Neighborhood Friend
SOUTH KOREA

 **DAIHAN Scientific Co., Ltd.**
Manufacturer and General Supplier

Shopping mall : DAIHAN-Sci.com and ALLforLAB.com

Head Office, Factory, LAB, and Logistics Center :
326, Sinpyeongseokhwa-ro, Jijeong-myeon, Wonju,
Gangwon State, South Korea

Tel : 080-008-3000 or 070-5020-1516
E-mail : dh@daihansci.co.kr
Fax : 033-737-7575

International Tel : +82-70-5020-9000
E-mail : info@daihan-sci.com
Fax : +82-2-963-5231



Cat 424 (23 Mar. 2026)

Distributor